**MACRODUCT<sup>®</sup>**  
Sweat Collection System


U.S. Patent Number 4,542,751  
U.K. Patent Number 2,116,850  
Other Patents Pending

**ABBREVIATED INSTRUCTIONS FOR USE**

**DO NOT TOUCH OR OTHERWISE CONTAMINATE THE CONICAL SURFACE OF MACRODUCT**

- 1. Attach the appropriately-sized strap to the MACRODUCT collector with proper orientation for subsequent attachment to limb.
- 2. Thoroughly clean and dry the stimulated skin area.
- 3. Apply the conical surface of MACRODUCT precisely over the stimulated area.
- 4. Fix MACRODUCT very tightly in position using the strap.  
**MACRODUCT REMAINS ATTACHED DURING STEPS 5 THROUGH 8.**
- 5. The protective transparent cover over the spiral microbore tube may be removed by inserting a pointed tool under one of the cut-out sections and prying upward. (The nippers may be conveniently employed for this operation.) It can be removed immediately after attachment of the MACRODUCT unless an elastic bandage or other overwrap is to be used.
- 6. When sufficient sweat has accumulated,\* connect the free end of the microbore tube to the needle of the sweat dispenser or blunt end needle on a syringe. Grasp the end of the microbore tube (not the dispenser or syringe) and carefully lift the tube away from the adhesive surface until its entire length is free.
- 7. Use the nippers to sever the attached end of the microbore tube as closely as possible to the upper surface of the MACRODUCT collector body.
- 8. Immediately transfer the sweat to the sealable microsample cup by gently and continuously squeezing the dispenser to expel the sweat from the microbore tube, or connect the tubing to the nipple of the conductivity cell.
- 9. Remove the MACRODUCT collector body from the patient's limb; retain the attachment strap and discard remainder.

\*The bore of the collecting tube is controlled to contain at least 2.7 µL/cm, corresponding to the volumes indicated in the diagram below.

**INFORMATION FOR PARENTS:**

**SWEAT TESTING POSES A REMOTE RISK OF MINOR SKIN BURNS**

There is an element of risk inherent in all medical procedures, no matter how simple. The sweat test has been an important laboratory tool since the 1950s. It provides a quantitative test result to confirm or exclude a clinical diagnosis of cystic fibrosis. Unfortunately, the test has been accompanied by occasional minor burns.

The sweat test consists of three sequential procedures: (1) sweat stimulation, (2) sweat collection, and (3) sweat analysis. The first procedure is known as pilocarpine iontophoresis. It is universally accepted by medical authorities as a safe and effective method of stimulating sweat glands. A sweat-inducing drug, pilocarpine, is delivered from the surface of the skin through the watery pathways of the sweat ducts into the sweat glands by a small electric current that is made to flow through the dermal layers. The electric current is supplied by a battery-powered device through a pair of electrodes fitted to the limb of the patient.

Minor skin burns have been an unwelcome, adverse side-effect of pilocarpine iontophoresis from the beginning. Some types of iontophoretic apparatus are prone to cause burns, particularly if there is procedural error. Fortunately, such burns are extremely rare with the Wescor iontophoretic system. It uses a sophisticated microprocessor current controller and a very low delivery current of only 1.5 milliamperes. Pilocarpine is contained in unique Pilogel gel drug reservoirs that are 96% water. These features substantially reduce, but do not totally eliminate, the possibility of skin burns.

Burn descriptions vary from "tiny black pinholes in the skin" to "crater-like, third-degree burns two to three millimeters in diameter." In most of the incidents reported, the children have exhibited no sign of pain or discomfort during iontophoresis, and the burn was not discovered until the electrodes were removed.

Most individuals exhibit a sensitivity to pilocarpine that is typically manifested as mild erythema (redness) of the skin at the electrode locations. In some cases, one or more blister-like welts may also form. These are often mistaken as burns, but they are simply the reaction of the skin to pilocarpine. Such "blisters" invariably disappear within 2 to 3 hours, leaving no after-effects.

Based on current data and reported events, the apparent burn rate is less than 1 in 50,000. Wescor prescribes proper test procedures which minimize the risk of burns from its equipment. It is highly unlikely that your child will suffer a burn during the sweat stimulation phase of the sweat test.

We realize these statistics will be of scant comfort to the parents of a child who has the misfortune of suffering the "one burn in 50,000." However, experience has shown that when burns do occur, the injuries are minor and there are no lasting effects. The burns usually heal completely within one to two weeks with little or no scarring.

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**Manufacturer's Name:** Wescor, Inc  
**Manufacturer's Address:** 459 South Main Street  
Logan, Utah 84321-5294  
USA

**DECLARATION  
of  
CONFORMITY**

**declares, the product**

**Product Name:** Macroduct® Sweat Collection System  
**Model Number:** 3700-SYS  
**Type of Equipment:** Battery-powered device used to induce sweat by iontophoresis.

**conforms to the applicable sections of the European Union Council Directive concerning medical devices — 93/42/EEC and as a result complies with:**

EN55011/CISPR11(1991) (Emissions) RF Emissions Group 1, Class A  
EN50082-1(1992) (Immunity)  
EN61000-4-2(1991) ESD — 8 kV Air Discharge, 4 kV Contact Discharge  
EN61000-4-3(1984) RF Immunity — 3 V/m  
and  
EN60601-1 (1990) (Safety)

**therefore bears the CE Marking as a Class IIA medical device. Certificate issued by Notified Body BSI, 0086.**

This declaration applies to the instrument itself and to the accessories and supplies listed hereafter.

Burt  
Name  
Director, QA/RA  
Title  
6/10/04  
Date  
Logan, Utah USA  
Place

## Macroduct Sweat Collection System Model 3700-SYS

### Accessories and Supplies

SS-032	Macroduct Supply Kit (for 6 sweat tests; contains Pilogel® Disc (12 each), Macroduct Sweat Collector (6 each), Small Sealable Container (6 each)).
SS-032-ND	Macroduct Supply Kit, no dye (for 6 sweat tests; contains Pilogel® Disc (12 each), Macroduct Sweat Collector with no dye (6 each), Small Sealable Container (6 each)).
SS-107	Small Sealable Containers (6 each)
SS-128	Macroduct Strap, Small, 14 cm (replaces RP-083)
SS-129	Macroduct Strap, Medium, 18 cm (replaces RP-084)
SS-130	Macroduct Strap, Large, 25 cm (replaces RP-085)
SS-131	Macroduct Strap, Extra Large, 39 cm (replaces RP-086)
SS-132	Macroduct Strap Set (4 straps, 1 of each size; replaces RP-087)
SS-142	Macroduct (collector only, 1 each; replaces SS-042)
SS-142-ND	Macroduct, no dye (collector only, 1 each; replaces SS-042)

### SECTION 1: Identification of the substance/mixture and of the company/undertaking

#### 1.1. Product identifier

Product form	: Mixture
Product name	: MACRODUCT® Supply Kit or Pilogel® Discs
Product code	: SS-032, SS-032G, SS-032IT, SS-032-ND, SS-032G-ND, SS-023, SS-023G, or SS-023IT
Product group	: Trade product

#### 1.2. Relevant identified uses of the substance or mixture and uses advised against

##### 1.2.1. Relevant identified uses

Industrial/Professional use spec	: For professional use only
Use of the substance/mixture	: Sweat Induction for Diagnosis of Cystic Fibrosis

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#### 1.3. Details of the supplier of the safety data sheet

ELITechGroup Inc.  
370 West 1700 South  
84321 Logan, UT - USA  
T +1 (435) 752-6011 - F +1 (435) 752-4127  
[qara\\_ebs@elitechgroup.com](mailto:qara_ebs@elitechgroup.com) - [www.elitechgroup.com](http://www.elitechgroup.com)

#### 1.4. Emergency telephone number

Emergency number	: Contact your distributor or poison control center in your country. InfoTrac Emergency Response: Calls within the USA, phone: 1-800-535-5053. Calls outside the USA, phone: +1 352-323-3500 (call collect) Customer ID: #90104 (NOTE: this number is required when a customer calls into either phone number above).
------------------	---

### SECTION 2: Hazards identification

#### 2.1. Classification of the substance or mixture

##### Classification according to Regulation (EC) No. 1272/2008 [CLP]

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#### 2.2. Label elements

##### Labeling according to Regulation (EC) No. 1272/2008 [CLP]

Hazard pictograms (CLP)

GHS07

Signal word (CLP)	: Warning
Hazardous ingredients	: Pilocarpine nitrate
Hazard statements (CLP)	: H302 - Harmful if swallowed
Precautionary statements (CLP)	: P264 - Wash hands, forearms and face thoroughly after handling P102 - Keep out of reach of children P270 - Do not eat, drink or smoke when using this product P301+P312 - IF SWALLOWED: Call a POISON CENTER, a doctor if you feel unwell P330 - Rinse mouth P501 - Dispose of contents/container to an authorized waste collection point
Security closing plug for children	: No

# MACRODUCT® Supply Kit or Pilogel® Discs

## Safety Data Sheet

according to Regulation (EU) 2015/830

Tactile warning : No

### 2.3. Other hazards

Adverse physicochemical, human health and environmental effects : Harmful if swallowed.

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Not applicable

### 3.2. Mixture

Name	Product identifier	%	Classification according to Directive 67/548/EEC	Classification according to Regulation (EC) No. 1272/2008 [CLP]
Pilocarpine nitrate	(CAS No) 148-72-1	< 1	Not classified	Acute Tox. 2 (Oral), H300 Acute Tox. 2 (Inhalation), H330

Full text of R- and H- phrases: see section 16

## SECTION 4: First aid measures

### 4.1. Description of first aid measures

First-aid measures general : Call a poison center/doctor/physician if you feel unwell.  
First-aid measures after inhalation : Remove person to fresh air and keep comfortable for breathing.  
First-aid measures after skin contact : Wash skin with plenty of water.  
First-aid measures after eye contact : Rinse eyes with water as a precaution.  
First-aid measures after ingestion : Rinse mouth. Call a poison center/doctor/physician if you feel unwell.

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## SECTION 5: Firefighting measures

### 5.1. Extinguishing media

Suitable extinguishing media : Water spray. Dry powder. Foam.

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Methods for cleaning up : Recover mechanically the product.

# MACRODUCT® Supply Kit or Pilogel® Discs

## Safety Data Sheet

according to Regulation (EU) 2015/830

Other information : Dispose of materials or solid residues at an authorized site.

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### SECTION 11: Toxicological information

#### 11.1. Information on toxicological effects

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## SECTION 1: Identification of the substance/mixture and of the company/undertaking

### 1.1. Product identifier

Product form : Mixture  
Trade name : CALIBRATOR for Sweat-Chek™  
Product code : SS-140  
Product group : Trade product

#44.2

### 1.2. Relevant identified uses of the substance or mixture and uses advised against

#### 1.2.1. Relevant identified uses

Industrial/Professional use spec : For professional use only  
Use of the substance/mixture : Sodium Chloride Solution, Saline

#### 1.2.2. Uses advised against

No additional information available

### 1.3. Details of the supplier of the safety data sheet

ELITechGroup Inc.  
370 West 1700 South  
84321 Logan, UT - USA  
T +1 (435) 752-6011 - F +1 (435) 752-4127  
[qara\\_ebs@elitechgroup.com](mailto:qara_ebs@elitechgroup.com) - [www.elitechgroup.com](http://www.elitechgroup.com)

### 1.4. Emergency telephone number

Emergency number : Contact your distributor or poison control center in your country.  
InfoTrac Emergency Response: Calls within the USA, phone: 1-800-535-5053. Calls outside the USA, phone: +1 352-323-3500 (call collect)  
Customer ID: #90104 (NOTE: this number is required when a customer calls into either phone number above).

## SECTION 2: Hazards identification

### 2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]

Not classified

### 2.2. Label elements

Labeling according to Regulation (EC) No. 1272/2008 [CLP]

Security closing plug for children : No  
Tactile warning : No

### 2.3. Other hazards

No additional information available

## SECTION 3: Composition/Information on ingredients

### 3.1. Substance

Not applicable

### 3.2. Mixture

This mixture does not contain any substances to be mentioned according to the criteria of section 3.2 of REACH annex II

# CALIBRATOR for Sweat-Chek™

## Safety Data Sheet

according to Regulation (EU) 2015/830

### SECTION 4: First aid measures

#### 4.1. Description of first aid measures

- First-aid measures general : Never give anything by mouth to an unconscious person. If you feel unwell, seek medical advice (show the label where possible).
- First-aid measures after inhalation : Allow victim to breathe fresh air. Allow the victim to rest.
- First-aid measures after skin contact : Remove affected clothing and wash all exposed skin area with mild soap and water, followed by warm water rinse.
- First-aid measures after eye contact : Rinse immediately with plenty of water. Obtain medical attention if pain, blinking or redness persist.
- First-aid measures after ingestion : Rinse mouth. Do NOT induce vomiting. Obtain emergency medical attention.

#### 4.2. Most important symptoms and effects, both acute and delayed

- Symptoms/injuries : Not expected to present a significant hazard under anticipated conditions of normal use.

#### 4.3. Indication of any immediate medical attention and special treatment needed

No additional information available

### SECTION 5: Firefighting measures

#### 5.1. Extinguishing media

- Suitable extinguishing media : Foam. Dry powder. Carbon dioxide. Water spray. Sand.
- Unsuitable extinguishing media : Do not use a heavy water stream.

#### 5.2. Special hazards arising from the substance or mixture

No additional information available

#### 5.3. Advice for firefighters

- Firefighting instructions : Use water spray or fog for cooling exposed containers. Exercise caution when fighting any chemical fire. Prevent fire-fighting water from entering environment.
- Protection during firefighting : Do not enter fire area without proper protective equipment, including respiratory protection.

### SECTION 6: Accidental release measures

#### 6.1. Personal precautions, protective equipment and emergency procedures

##### 6.1.1. For non-emergency personnel

- Emergency procedures : Evacuate unnecessary personnel.

##### 6.1.2. For emergency responders

- Protective equipment : Equip cleanup crew with proper protection.
- Emergency procedures : Ventilate area.

#### 6.2. Environmental precautions

No additional information available

#### 6.3. Methods and material for containment and cleaning up

No additional information available

#### 6.4. Reference to other sections

See Heading 8. Exposure controls and personal protection.

### SECTION 7: Handling and storage

#### 7.1. Precautions for safe handling

- Precautions for safe handling : Wash hands and other exposed areas with mild soap and water before eating, drinking or smoking and when leaving work.

#### 7.2. Conditions for safe storage, including any incompatibilities

- Incompatible products : Strong bases. Strong acids.

#### 7.3. Specific end use(s)

No additional information available

# CALIBRATOR for Sweat-Chek™

## Safety Data Sheet

according to Regulation (EU) 2015/830

### SECTION 8: Exposure controls/personal protection

#### 8.1. Control parameters

No additional information available

#### 8.2. Exposure controls

Personal protective equipment : Avoid all unnecessary exposure.

Other information : Do not eat, drink or smoke during use.

### SECTION 9: Physical and chemical properties

#### 9.1. Information on basic physical and chemical properties

Physical state	: Liquid
Color	: Colorless.
Odor	: characteristic.
Odor threshold	: No data available
pH	: No data available
Relative evaporation rate (butyl acetate=1)	: No data available
Melting point	: No data available
Freezing point	: No data available
Boiling point	: 100 °C
Flash point	: No data available
Auto-ignition temperature	: No data available
Decomposition temperature	: No data available
Flammability (solid, gas)	: Non flammable
Vapor pressure	: No data available
Relative vapor density at 20 °C	: No data available
Relative density	: No data available
Solubility	: Water: No available data
Log Pow	: No data available
Viscosity, kinematic	: No data available
Viscosity, dynamic	: No data available
Explosive properties	: No data available
Oxidizing properties	: No data available
Explosion limits	: No data available

#### 9.2. Other information

No additional information available

### SECTION 10: Stability and reactivity

#### 10.1. Reactivity

No additional information available

#### 10.2. Chemical stability

Not established.

#### 10.3. Possibility of hazardous reactions

Not established.

#### 10.4. Conditions to avoid

No additional information available

#### 10.5. Incompatible materials

Strong acids. Strong bases.

# CALIBRATOR for Sweat-Chek™

## Safety Data Sheet

according to Regulation (EU) 2015/830

### 10.6. Hazardous decomposition products

No additional information available

## SECTION 11: Toxicological information

### 11.1. Information on toxicological effects

Acute toxicity	: Not classified
Skin corrosion/irritation	: Not classified Based on available data, the classification criteria are not met
Serious eye damage/irritation	: Not classified Based on available data, the classification criteria are not met
Respiratory or skin sensitization	: Not classified Based on available data, the classification criteria are not met
Germ cell mutagenicity	: Not classified Based on available data, the classification criteria are not met
Carcinogenicity	: Not classified Based on available data, the classification criteria are not met
Reproductive toxicity	: Not classified Based on available data, the classification criteria are not met
Specific target organ toxicity (single exposure)	: Not classified Based on available data, the classification criteria are not met
Specific target organ toxicity (repeated exposure)	: Not classified Based on available data, the classification criteria are not met
Aspiration hazard	: Not classified Based on available data, the classification criteria are not met
Potential Adverse human health effects and symptoms	: Based on available data, the classification criteria are not met.

## SECTION 12: Ecological information

### 12.1. Toxicity

No additional information available

### 12.2. Persistence and degradability

CALIBRATOR for Sweat-Chek	
Persistence and degradability	Not established.

### 12.3. Bioaccumulative potential

CALIBRATOR for Sweat-Chek	
Bioaccumulative potential	Not established.

### 12.4. Mobility in soil

No additional information available

### 12.5. Results of PBT and vPvB assessment

No additional information available

### 12.6. Other adverse effects

Additional information : Avoid release to the environment

## SECTION 13: Disposal considerations

### 13.1. Waste treatment methods

Waste disposal recommendations : Dispose in a safe manner in accordance with local/national regulations.

# CALIBRATOR for Sweat-Chek™

## Safety Data Sheet

according to Regulation (EU) 2015/830

### SECTION 14: Transport information

In accordance with ADR / RID / IMDG / IATA / ADN

#### 14.1. UN number

Not regulated for transport

#### 14.2. UN proper shipping name

Proper Shipping Name (ADR) : Not applicable

Proper Shipping Name (IMDG) : Not applicable

Proper Shipping Name (IATA) : Not applicable

Proper Shipping Name (ADN) : Not applicable

Proper Shipping Name (RID) : Not applicable

#### 14.3. Transport hazard class(es)

##### ADR

Transport hazard class(es) (ADR) : Not applicable

##### IMDG

Transport hazard class(es) (IMDG) : Not applicable

##### IATA

Transport hazard class(es) (IATA) : Not applicable

##### ADN

Transport hazard class(es) (ADN) : Not applicable

##### RID

Transport hazard class(es) (RID) : Not applicable

#### 14.4. Packing group

Packing group (ADR) : Not applicable

Packing group (IMDG) : Not applicable

Packing group (IATA) : Not applicable

Packing group (ADN) : Not applicable

Packing group (RID) : Not applicable

#### 14.5. Environmental hazards

Dangerous for the environment : No

Marine pollutant : No

Other information : No supplementary information available

#### 14.6. Special precautions for user

##### - Overland transport

No data available

##### - Transport by sea

No data available

##### - Air transport

No data available

##### - Inland waterway transport

No data available

# CALIBRATOR for Sweat-Chek™

## Safety Data Sheet

according to Regulation (EU) 2015/830

### - Rail transport

No data available

### 14.7. Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not applicable

## SECTION 15: Regulatory information

### 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

#### 15.1.1. EU-Regulations

Contains no REACH substances with Annex XVII restrictions

Contains no REACH candidate substance

Contains no REACH Annex XIV substances.

### 15.2. Chemical safety assessment

No chemical safety assessment has been carried out

## SECTION 16: Other information

Other information : None.

SDS EU Custom - EBS

*This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product.*

Reason for change: Format change to meet GHS compliance.

## SECTION 1: Identification of the substance/mixture and of the company/undertaking

### 1.1. Product identifier

Product form : Mixture #44.3  
 Product name : Sweat Controls for Cystic Fibrosis Testing  
 Product code : SS-150  
 Product group : Trade product

### 1.2. Relevant identified uses of the substance or mixture and uses advised against

#### 1.2.1. Relevant identified uses

Industrial/Professional use spec : For professional use only  
 Use of the substance/mixture : Sodium and Chloride in Simulated Human Sweat Matrix. Sweat Control Levels: 1, 2, and 3.

#### 1.2.2. Uses advised against

No additional information available

### 1.3. Details of the supplier of the safety data sheet

ELITechGroup Inc.  
 370 West 1700 South  
 84321 Logan, UT - USA  
 T +1 (435) 752-6011 - F +1 (435) 752-4127  
[qara\\_ebs@elitechgroup.com](mailto:qara_ebs@elitechgroup.com) - [www.elitechgroup.com](http://www.elitechgroup.com)

### 1.4. Emergency telephone number

Emergency number : Contact your distributor or poison control center in your country.  
 InfoTrac Emergency Response: Calls within the USA, phone: 1-800-535-5053. Calls outside the USA, phone: +1 352-323-3500 (call collect)  
 Customer ID: #90104 (NOTE: this number is required when a customer calls into either phone number above).

## SECTION 2: Hazards identification

### 2.1. Classification of the substance or mixture

#### Classification according to Regulation (EC) No. 1272/2008 [CLP]

Hazardous to the aquatic environment - Chronic Hazard Category 3 H412

Full text of classification categories and H statements : see section 16

### 2.2. Label elements

#### Labeling according to Regulation (EC) No. 1272/2008 [CLP]

Hazard pictograms (CLP) :

GHS07

Signal word (CLP) : Warning  
 Hazard statements (CLP) : H317 - May cause an allergic skin reaction  
 H412 - Harmful to aquatic life with long lasting effects  
 Precautionary statements (CLP) : P273 - Avoid release to the environment  
 P302+P352 - IF ON SKIN: Wash with plenty of water  
 EUH phrases : EUH208 - Contains METHYLCHLOROISOTHIAZOLINONE(26172-55-4), METHYLISOTHIAZOLINONE(2682-20-4). May produce an allergic reaction  
 Security closing plug for children : No  
 Tactile warning : No

# Sweat Controls for Cystic Fibrosis Testing

## Safety Data Sheet

according to Regulation (EU) 2015/830

### 2.3. Other hazards

Adverse physicochemical, human health and environmental effects : To our knowledge, this product does not present any particular risk, provided it is handled in accordance with good occupational hygiene and safety practice.

## SECTION 3: Composition/Information on ingredients

### 3.1. Substance

Not applicable

### 3.2. Mixture

Name	Product identifier	%	GHS-US classification
5-chloro-2-methyl-2H-isothiazol-3-one	(CAS No) 26172-55-4 (EC Number) 247-500-7	< 0.5	Acute Tox. 3 (Oral), H301 Acute Tox. 3 (Dermal), H311 Skin Corr. 1B, H314 Eye Dam. 1, H318 Skin Sens. 1, H317 STOT SE 3, H335 Aquatic Acute 1, H400
2-methyl-2H-isothiazol-3-one	(CAS No) 2682-20-4 (EC Number) 220-239-6	< 0.5	Acute Tox. 3 (Oral), H301 Acute Tox. 3 (Dermal), H311 Skin Corr. 1B, H314 Eye Dam. 1, H318 Skin Sens. 1, H317 STOT SE 3, H335 Aquatic Acute 1, H400

## SECTION 4: First aid measures

### 4.1. Description of first aid measures

First-aid measures after inhalation : Remove person to fresh air and keep comfortable for breathing.  
First-aid measures after skin contact : Wash skin with plenty of water.  
First-aid measures after eye contact : Rinse eyes with water as a precaution.  
First-aid measures after ingestion : Call a poison center/doctor/physician if you feel unwell.

### 4.2. Most important symptoms and effects, both acute and delayed

No additional information available

### 4.3. Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

## SECTION 5: Firefighting measures

### 5.1. Extinguishing media

Suitable extinguishing media : Water spray. Dry powder. Foam. Carbon dioxide.

### 5.2. Special hazards arising from the substance or mixture

Hazardous decomposition products in case of fire : Toxic fumes may be released.

### 5.3. Advice for firefighters

Protection during firefighting : Do not attempt to take action without suitable protective equipment. Self-contained breathing apparatus. Complete protective clothing.

## SECTION 6: Accidental release measures

### 6.1. Personal precautions, protective equipment and emergency procedures

#### 6.1.1. For non-emergency personnel

Emergency procedures : Ventilate spillage area.

#### 6.1.2. For emergency responders

Protective equipment : Do not attempt to take action without suitable protective equipment. For further information refer to section 8: "Exposure controls/personal protection".

### 6.2. Environmental precautions

Avoid release to the environment.

# Sweat Controls for Cystic Fibrosis Testing

## Safety Data Sheet

according to Regulation (EU) 2015/830

### 6.3. Methods and material for containment and cleaning up

Methods for cleaning up : Take up liquid spill into absorbent material.  
Other information : Dispose of materials or solid residues at an authorized site.

### 6.4. Reference to other sections

For further information refer to section 13.

## SECTION 7: Handling and storage

### 7.1. Precautions for safe handling

Precautions for safe handling : Ensure good ventilation of the work station. Wear personal protective equipment.  
Hygiene measures : Do not eat, drink or smoke when using this product. Always wash hands after handling the product.

### 7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Store in a well-ventilated place. Keep cool.

### 7.3. Specific end use(s)

No additional information available

## SECTION 8: Exposure controls/personal protection

### 8.1. Control parameters

No additional information available

### 8.2. Exposure controls

Appropriate engineering controls : Ensure good ventilation of the work station.  
Hand protection : Protective gloves  
Environmental exposure controls : Avoid release to the environment.

## SECTION 9: Physical and chemical properties

### 9.1. Information on basic physical and chemical properties

Physical state : Liquid  
Color : No data available  
Odor : No data available  
Odor threshold : No data available  
pH : No data available  
Relative evaporation rate (butyl acetate=1) : No data available  
Melting point : Not applicable  
Freezing point : No data available  
Boiling point : No data available  
Flash point : No data available  
Auto-ignition temperature : No data available  
Decomposition temperature : No data available  
Flammability (solid, gas) : Not applicable  
Vapor pressure : No data available  
Relative vapor density at 20 °C : No data available  
Relative density : No data available  
Solubility : Water: No data available  
Log Pow : No data available  
Viscosity, kinematic : No data available  
Viscosity, dynamic : No data available  
Explosive properties : No data available  
Oxidizing properties : No data available

# Sweat Controls for Cystic Fibrosis Testing

## Safety Data Sheet

according to Regulation (EU) 2015/830

Explosion limits : No data available

### 9.2. Other information

No additional information available

## SECTION 10: Stability and reactivity

### 10.1. Reactivity

The product is non-reactive under normal conditions of use, storage and transport.

### 10.2. Chemical stability

Stable under normal conditions.

### 10.3. Possibility of hazardous reactions

No dangerous reactions known under normal conditions of use.

### 10.4. Conditions to avoid

None under recommended storage and handling conditions (see section 7).

### 10.5. Incompatible materials

No additional information available

### 10.6. Hazardous decomposition products

Under normal conditions of storage and use, hazardous decomposition products should not be produced.

## SECTION 11: Toxicological information

### 11.1. Information on toxicological effects

Acute toxicity : Not classified

5-chloro-2-methyl-2H-isothiazol-3-one (2682-20-4)	
ATE US (oral)	100.000 mg/kg body weight
ATE US (dermal)	300.000 mg/kg body weight
2-methyl-2H-isothiazol-3-one (26172-55-4)	
ATE US (oral)	100.000 mg/kg body weight
ATE US (dermal)	300.000 mg/kg body weight

Skin corrosion/irritation : Not classified

Serious eye damage/irritation : Not classified

Respiratory or skin sensitization : Not classified

Germ cell mutagenicity : Not classified

Carcinogenicity : Not classified

Reproductive toxicity : Not classified

Specific target organ toxicity (single exposure) : Not classified

Specific target organ toxicity (repeated exposure) : Not classified

Aspiration hazard : Not classified

## SECTION 12: Ecological information

### 12.1. Toxicity

Ecology - general : The product is not considered harmful to aquatic organisms or to cause long-term adverse effects in the environment.

### 12.2. Persistence and degradability

No additional information available

### 12.3. Bioaccumulative potential

No additional information available

# Sweat Controls for Cystic Fibrosis Testing

## Safety Data Sheet

according to Regulation (EU) 2015/830

### 12.4. Mobility in soil

No additional information available

### 12.5. Results of PBT and vPvB assessment

No additional information available

### 12.6. Other adverse effects

No additional information available

## SECTION 13: Disposal considerations

### 13.1. Waste treatment methods

Waste treatment methods : Dispose of contents/container in accordance with licensed collector's sorting instructions.

## SECTION 14: Transport information

In accordance with ADR / RID / IMDG / IATA / ADN

### 14.1. UN number

Not regulated for transport

### 14.2. UN proper shipping name

Proper Shipping Name (ADR) : Not applicable

Proper Shipping Name (IMDG) : Not applicable

Proper Shipping Name (IATA) : Not applicable

Proper Shipping Name (ADN) : Not applicable

Proper Shipping Name (RID) : Not applicable

### 14.3. Transport hazard class(es)

#### ADR

Transport hazard class(es) (ADR) : Not applicable

#### IMDG

Transport hazard class(es) (IMDG) : Not applicable

#### IATA

Transport hazard class(es) (IATA) : Not applicable

#### ADN

Transport hazard class(es) (ADN) : Not applicable

#### RID

Transport hazard class(es) (RID) : Not applicable

### 14.4. Packing group

Packing group (ADR) : Not applicable

Packing group (IMDG) : Not applicable

Packing group (IATA) : Not applicable

Packing group (ADN) : Not applicable

Packing group (RID) : Not applicable

### 14.5. Environmental hazards

Dangerous for the environment : No

Marine pollutant : No

Other information : No supplementary information available

# Sweat Controls for Cystic Fibrosis Testing

## Safety Data Sheet

according to Regulation (EU) 2015/830

### 14.6. Special precautions for user

#### - Overland transport

No data available

#### - Transport by sea

No data available

#### - Air transport

No data available

#### - Inland waterway transport

No data available

#### - Rail transport

No data available

### 14.7. Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not applicable

## SECTION 15: Regulatory information

### 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

#### 15.1.1. EU-Regulations

Contains no REACH substances with Annex XVII restrictions

Contains no REACH candidate substance

Contains no REACH Annex XIV substances.

### 15.2. Chemical safety assessment

No chemical safety assessment has been carried out

## SECTION 16: Other information

Full text of H-phrases:

Acute Tox. 3 (Dermal)	Acute toxicity (dermal) Category 3
Acute Tox. 3 (Oral)	Acute toxicity (oral) Category 3
Aquatic Acute 1	Hazardous to the aquatic environment - Acute Hazard Category 1
Eye Dam. 1	Serious eye damage/eye irritation Category 1
Skin Corr. 1B	Skin corrosion/irritation Category 1B
Skin Sens. 1	Skin sensitization Category 1
STOT SE 3	Specific target organ toxicity (single exposure) Category 3
H301	Toxic if swallowed
H311	Toxic in contact with skin
H314	Causes severe skin burns and eye damage
H317	May cause an allergic skin reaction
H318	Causes serious eye damage
H335	May cause respiratory irritation
H400	Very toxic to aquatic life

SDS EU Custom – EBS

*This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product.*

Reason for change: Format change to meet GHS compliance.

#42.2

**GLASS ENCAPSULATED CALIBRATOR SOLUTION  
FOR SWEAT-CHEK™ CONDUCTIVITY ANALYZER  
For In Vitro Diagnostic Use**

**Intended Use**

Wescor's CALIBRATOR solution is for calibration of the Sweat-Chek Sweat Conductivity Analyzer.

**Product Description**

Wescor's CALIBRATOR solution is liquid, ready to use, and requires no reconstitution or dilution. Each ampule contains 0.75 mL of NaCl solution, at 90 mmol/L molarity.

**Warnings and Precautions**

No special precautions are required in handling the product other than those routinely used in the laboratory.

**Storage and Stability**

1. CALIBRATOR should be stored at room temperature.
2. Each ampule contains 0.75 mL of solution. This volume adequately mitigates evaporative concentration for a few hours after the ampule is broken. Ampules are intended for one-time use only.

**Instructions for Use**

1. Invert gently to ensure homogeneity of the contents.
2. Before opening, flip the neck of the ampule or tap the ampule bottom against a hard surface to clear solution from the neck.
3. Place the ampule in the ampule breaking position of the organizer.
4. Place the safety sleeve over the neck of the ampule.
5. Holding ampule organizer down against the counter, grasp the neck of the ampule firmly between thumb and forefinger and snap the neck off.
6. Sample directly from the ampule, using a fresh Sweat-Chek take up tube to avoid contamination of the solution.
7. Follow procedures in User's Manual to assay CALIBRATOR solution.

**Expected Values**

The Sweat-Chek conductivity cell has an extremely stable response characteristic that is virtually linear through the critical range from 75 to 110 mmol/L. The main purpose of calibration is to ascertain that the conductivity cell temperature is within the correct range and that the electronics are otherwise functioning normally. If the conductivity cell is free of residual salt or water, a reading within  $\pm 2$  of the labeled value of 90 mmol/L indicates overall performance to a very high level of confidence. If the reading falls outside of this range, follow the instructions listed in the User's Manual for calibrating the instrument.

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IS20A46

Glass encapsulated calibrator solution for sweat-check™ conductivity analyzer

## SECTION 1: Identification of the substance/mixture and of the company/undertaking

### 1.1. Product identifier

Product form : Mixture  
Trade name : CALIBRATOR for Sweat-Chek™  
Product code : SS-140  
Product group : Trade product

### 1.2. Relevant identified uses of the substance or mixture and uses advised against

#### 1.2.1. Relevant identified uses

Industrial/Professional use spec : For professional use only  
Use of the substance/mixture : Sodium Chloride Solution, Saline

#### 1.2.2. Uses advised against

No additional information available

### 1.3. Details of the supplier of the safety data sheet

ELITechGroup Inc.  
370 West 1700 South  
84321 Logan, UT - USA  
T +1 (435) 752-6011 - F +1 (435) 752-4127  
[qara\\_ebs@elitechgroup.com](mailto:qara_ebs@elitechgroup.com) - [www.elitechgroup.com](http://www.elitechgroup.com)

### 1.4. Emergency telephone number

Emergency number : Contact your distributor or poison control center in your country.  
InfoTrac Emergency Response: Calls within the USA, phone: 1-800-535-5053. Calls outside the USA, phone: +1 352-323-3500 (call collect)  
Customer ID: #90104 (NOTE: this number is required when a customer calls into either phone number above).

## SECTION 2: Hazards identification

### 2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]  
Not classified

### 2.2. Label elements

#### Labeling according to Regulation (EC) No. 1272/2008 [CLP]

Security closing plug for children : No  
Tactile warning : No

### 2.3. Other hazards

No additional information available

## SECTION 3: Composition/Information on ingredients

### 3.1. Substance

Not applicable

### 3.2. Mixture

This mixture does not contain any substances to be mentioned according to the criteria of section 3.2 of REACH annex II

# CALIBRATOR for Sweat-Chek™

## Safety Data Sheet

according to Regulation (EU) 2015/830

### SECTION 4: First aid measures

#### 4.1. Description of first aid measures

- First-aid measures general : Never give anything by mouth to an unconscious person. If you feel unwell, seek medical advice (show the label where possible).
- First-aid measures after inhalation : Allow victim to breathe fresh air. Allow the victim to rest.
- First-aid measures after skin contact : Remove affected clothing and wash all exposed skin area with mild soap and water, followed by warm water rinse.
- First-aid measures after eye contact : Rinse immediately with plenty of water. Obtain medical attention if pain, blinking or redness persist.
- First-aid measures after ingestion : Rinse mouth. Do NOT induce vomiting. Obtain emergency medical attention.

#### 4.2. Most important symptoms and effects, both acute and delayed

- Symptoms/injuries : Not expected to present a significant hazard under anticipated conditions of normal use.

#### 4.3. Indication of any immediate medical attention and special treatment needed

No additional information available

### SECTION 5: Firefighting measures

#### 5.1. Extinguishing media

- Suitable extinguishing media : Foam. Dry powder. Carbon dioxide. Water spray. Sand.
- Unsuitable extinguishing media : Do not use a heavy water stream.

#### 5.2. Special hazards arising from the substance or mixture

No additional information available

#### 5.3. Advice for firefighters

- Firefighting instructions : Use water spray or fog for cooling exposed containers. Exercise caution when fighting any chemical fire. Prevent fire-fighting water from entering environment.
- Protection during firefighting : Do not enter fire area without proper protective equipment, including respiratory protection.

### SECTION 6: Accidental release measures

#### 6.1. Personal precautions, protective equipment and emergency procedures

##### 6.1.1. For non-emergency personnel

- Emergency procedures : Evacuate unnecessary personnel.

##### 6.1.2. For emergency responders

- Protective equipment : Equip cleanup crew with proper protection.
- Emergency procedures : Ventilate area.

#### 6.2. Environmental precautions

No additional information available

#### 6.3. Methods and material for containment and cleaning up

No additional information available

#### 6.4. Reference to other sections

See Heading 8. Exposure controls and personal protection.

### SECTION 7: Handling and storage

#### 7.1. Precautions for safe handling

- Precautions for safe handling : Wash hands and other exposed areas with mild soap and water before eating, drinking or smoking and when leaving work.

#### 7.2. Conditions for safe storage, including any incompatibilities

- Incompatible products : Strong bases. Strong acids.

#### 7.3. Specific end use(s)

No additional information available

# CALIBRATOR for Sweat-Chek™

## Safety Data Sheet

according to Regulation (EU) 2015/830

### SECTION 8: Exposure controls/personal protection

#### 8.1. Control parameters

No additional information available

#### 8.2. Exposure controls

Personal protective equipment : Avoid all unnecessary exposure.

Other information : Do not eat, drink or smoke during use.

### SECTION 9: Physical and chemical properties

#### 9.1. Information on basic physical and chemical properties

Physical state	: Liquid
Color	: Colorless.
Odor	: characteristic.
Odor threshold	: No data available
pH	: No data available
Relative evaporation rate (butyl acetate=1)	: No data available
Melting point	: No data available
Freezing point	: No data available
Boiling point	: 100 °C
Flash point	: No data available
Auto-ignition temperature	: No data available
Decomposition temperature	: No data available
Flammability (solid, gas)	: Non flammable
Vapor pressure	: No data available
Relative vapor density at 20 °C	: No data available
Relative density	: No data available
Solubility	: Water: No available data
Log Pow	: No data available
Viscosity, kinematic	: No data available
Viscosity, dynamic	: No data available
Explosive properties	: No data available
Oxidizing properties	: No data available
Explosion limits	: No data available

#### 9.2. Other information

No additional information available

### SECTION 10: Stability and reactivity

#### 10.1. Reactivity

No additional information available

#### 10.2. Chemical stability

Not established.

#### 10.3. Possibility of hazardous reactions

Not established.

#### 10.4. Conditions to avoid

No additional information available

#### 10.5. Incompatible materials

Strong acids. Strong bases.

# CALIBRATOR for Sweat-Chek™

## Safety Data Sheet

according to Regulation (EU) 2015/830

### 10.6. Hazardous decomposition products

No additional information available

## SECTION 11: Toxicological information

### 11.1. Information on toxicological effects

Acute toxicity	: Not classified
Skin corrosion/irritation	: Not classified Based on available data, the classification criteria are not met
Serious eye damage/irritation	: Not classified Based on available data, the classification criteria are not met
Respiratory or skin sensitization	: Not classified Based on available data, the classification criteria are not met
Germ cell mutagenicity	: Not classified Based on available data, the classification criteria are not met
Carcinogenicity	: Not classified Based on available data, the classification criteria are not met
Reproductive toxicity	: Not classified Based on available data, the classification criteria are not met
Specific target organ toxicity (single exposure)	: Not classified Based on available data, the classification criteria are not met
Specific target organ toxicity (repeated exposure)	: Not classified Based on available data, the classification criteria are not met
Aspiration hazard	: Not classified Based on available data, the classification criteria are not met
Potential Adverse human health effects and symptoms	: Based on available data, the classification criteria are not met.

## SECTION 12: Ecological information

### 12.1. Toxicity

No additional information available

### 12.2. Persistence and degradability

CALIBRATOR for Sweat-Chek	
Persistence and degradability	Not established.

### 12.3. Bioaccumulative potential

CALIBRATOR for Sweat-Chek	
Bioaccumulative potential	Not established.

### 12.4. Mobility in soil

No additional information available

### 12.5. Results of PBT and vPvB assessment

No additional information available

### 12.6. Other adverse effects

Additional information : Avoid release to the environment

## SECTION 13: Disposal considerations

### 13.1. Waste treatment methods

Waste disposal recommendations : Dispose in a safe manner in accordance with local/national regulations.

# CALIBRATOR for Sweat-Chek™

## Safety Data Sheet

according to Regulation (EU) 2015/830

### SECTION 14: Transport information

In accordance with ADR / RID / IMDG / IATA / ADN

#### 14.1. UN number

Not regulated for transport

#### 14.2. UN proper shipping name

Proper Shipping Name (ADR) : Not applicable

Proper Shipping Name (IMDG) : Not applicable

Proper Shipping Name (IATA) : Not applicable

Proper Shipping Name (ADN) : Not applicable

Proper Shipping Name (RID) : Not applicable

#### 14.3. Transport hazard class(es)

##### ADR

Transport hazard class(es) (ADR) : Not applicable

##### IMDG

Transport hazard class(es) (IMDG) : Not applicable

##### IATA

Transport hazard class(es) (IATA) : Not applicable

##### ADN

Transport hazard class(es) (ADN) : Not applicable

##### RID

Transport hazard class(es) (RID) : Not applicable

#### 14.4. Packing group

Packing group (ADR) : Not applicable

Packing group (IMDG) : Not applicable

Packing group (IATA) : Not applicable

Packing group (ADN) : Not applicable

Packing group (RID) : Not applicable

#### 14.5. Environmental hazards

Dangerous for the environment : No

Marine pollutant : No

Other information : No supplementary information available

#### 14.6. Special precautions for user

##### - Overland transport

No data available

##### - Transport by sea

No data available

##### - Air transport

No data available

##### - Inland waterway transport

No data available

# CALIBRATOR for Sweat-Chek™

## Safety Data Sheet

according to Regulation (EU) 2015/830

### - Rail transport

No data available

### 14.7. Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not applicable

## SECTION 15: Regulatory information

### 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

#### 15.1.1. EU-Regulations

Contains no REACH substances with Annex XVII restrictions

Contains no REACH candidate substance

Contains no REACH Annex XIV substances.

### 15.2. Chemical safety assessment

No chemical safety assessment has been carried out

## SECTION 16: Other information

Other information : None.

SDS EU Custom - EBS

*This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product.*

Reason for change: Format change to meet GHS compliance.

## SECTION 1: Identification of the substance/mixture and of the company/undertaking

### 1.1. Product identifier

Product form : Mixture  
Product name : Sweat Controls for Cystic Fibrosis Testing  
Product code : SS-150  
Product group : Trade product

### 1.2. Relevant identified uses of the substance or mixture and uses advised against

#### 1.2.1. Relevant identified uses

Industrial/Professional use spec : For professional use only  
Use of the substance/mixture : Sodium and Chloride in Simulated Human Sweat Matrix. Sweat Control Levels: 1, 2, and 3.

#### 1.2.2. Uses advised against

No additional information available

### 1.3. Details of the supplier of the safety data sheet

ELITechGroup Inc.  
370 West 1700 South  
84321 Logan, UT - USA  
T +1 (435) 752-6011 - F +1 (435) 752-4127  
[qara\\_ebs@elitechgroup.com](mailto:qara_ebs@elitechgroup.com) - [www.elitechgroup.com](http://www.elitechgroup.com)

### 1.4. Emergency telephone number

Emergency number : Contact your distributor or poison control center in your country.  
InfoTrac Emergency Response: Calls within the USA, phone: 1-800-535-5053. Calls outside the USA, phone: +1 352-323-3500 (call collect)  
Customer ID: #90104 (NOTE: this number is required when a customer calls into either phone number above).

## SECTION 2: Hazards identification

### 2.1. Classification of the substance or mixture

#### Classification according to Regulation (EC) No. 1272/2008 [CLP]

Hazardous to the aquatic environment - Chronic Hazard Category 3 H412

Full text of classification categories and H statements : see section 16

### 2.2. Label elements

#### Labeling according to Regulation (EC) No. 1272/2008 [CLP]

Hazard pictograms (CLP) :

GHS07

Signal word (CLP) : Warning  
Hazard statements (CLP) : H317 - May cause an allergic skin reaction  
H412 - Harmful to aquatic life with long lasting effects  
Precautionary statements (CLP) : P273 - Avoid release to the environment  
P302+P352 - IF ON SKIN: Wash with plenty of water  
EUH phrases : EUH208 - Contains METHYLCHLORISOTHIAZOLINONE(26172-55-4), METHYLISOTHIAZOLINONE(2682-20-4). May produce an allergic reaction  
Security closing plug for children : No  
Tactile warning : No

# Sweat Controls for Cystic Fibrosis Testing

## Safety Data Sheet

according to Regulation (EU) 2015/830

### 2.3. Other hazards

Adverse physicochemical, human health and environmental effects : To our knowledge, this product does not present any particular risk, provided it is handled in accordance with good occupational hygiene and safety practice.

## SECTION 3: Composition/Information on ingredients

### 3.1. Substance

Not applicable

### 3.2. Mixture

Name	Product identifier	%	GHS-US classification
5-chloro-2-methyl-2H-isothiazol-3-one	(CAS No) 26172-55-4 (EC Number) 247-500-7	< 0.5	Acute Tox. 3 (Oral), H301 Acute Tox. 3 (Dermal), H311 Skin Corr. 1B, H314 Eye Dam. 1, H318 Skin Sens. 1, H317 STOT SE 3, H335 Aquatic Acute 1, H400
2-methyl-2H-isothiazol-3-one	(CAS No) 2682-20-4 (EC Number) 220-239-6	< 0.5	Acute Tox. 3 (Oral), H301 Acute Tox. 3 (Dermal), H311 Skin Corr. 1B, H314 Eye Dam. 1, H318 Skin Sens. 1, H317 STOT SE 3, H335 Aquatic Acute 1, H400

## SECTION 4: First aid measures

### 4.1. Description of first aid measures

First-aid measures after inhalation : Remove person to fresh air and keep comfortable for breathing.  
First-aid measures after skin contact : Wash skin with plenty of water.  
First-aid measures after eye contact : Rinse eyes with water as a precaution.  
First-aid measures after ingestion : Call a poison center/doctor/physician if you feel unwell.

### 4.2. Most important symptoms and effects, both acute and delayed

No additional information available

### 4.3. Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

## SECTION 5: Firefighting measures

### 5.1. Extinguishing media

Suitable extinguishing media : Water spray. Dry powder. Foam. Carbon dioxide.

### 5.2. Special hazards arising from the substance or mixture

Hazardous decomposition products in case of fire : Toxic fumes may be released.

### 5.3. Advice for firefighters

Protection during firefighting : Do not attempt to take action without suitable protective equipment. Self-contained breathing apparatus. Complete protective clothing.

## SECTION 6: Accidental release measures

### 6.1. Personal precautions, protective equipment and emergency procedures

#### 6.1.1. For non-emergency personnel

Emergency procedures : Ventilate spillage area.

#### 6.1.2. For emergency responders

Protective equipment : Do not attempt to take action without suitable protective equipment. For further information refer to section 8: "Exposure controls/personal protection".

### 6.2. Environmental precautions

Avoid release to the environment.

# Sweat Controls for Cystic Fibrosis Testing

## Safety Data Sheet

according to Regulation (EU) 2015/830

### 6.3. Methods and material for containment and cleaning up

Methods for cleaning up : Take up liquid spill into absorbent material.  
Other information : Dispose of materials or solid residues at an authorized site.

### 6.4. Reference to other sections

For further information refer to section 13.

## SECTION 7: Handling and storage

### 7.1. Precautions for safe handling

Precautions for safe handling : Ensure good ventilation of the work station. Wear personal protective equipment.  
Hygiene measures : Do not eat, drink or smoke when using this product. Always wash hands after handling the product.

### 7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Store in a well-ventilated place. Keep cool.

### 7.3. Specific end use(s)

No additional information available

## SECTION 8: Exposure controls/personal protection

### 8.1. Control parameters

No additional information available

### 8.2. Exposure controls

Appropriate engineering controls : Ensure good ventilation of the work station.  
Hand protection : Protective gloves  
Environmental exposure controls : Avoid release to the environment.

## SECTION 9: Physical and chemical properties

### 9.1. Information on basic physical and chemical properties

Physical state : Liquid  
Color : No data available  
Odor : No data available  
Odor threshold : No data available  
pH : No data available  
Relative evaporation rate (butyl acetate=1) : No data available  
Melting point : Not applicable  
Freezing point : No data available  
Boiling point : No data available  
Flash point : No data available  
Auto-ignition temperature : No data available  
Decomposition temperature : No data available  
Flammability (solid, gas) : Not applicable  
Vapor pressure : No data available  
Relative vapor density at 20 °C : No data available  
Relative density : No data available  
Solubility : Water: No data available  
Log Pow : No data available  
Viscosity, kinematic : No data available  
Viscosity, dynamic : No data available  
Explosive properties : No data available  
Oxidizing properties : No data available

# Sweat Controls for Cystic Fibrosis Testing

## Safety Data Sheet

according to Regulation (EU) 2015/830

Explosion limits : No data available

### 9.2. Other information

No additional information available

## SECTION 10: Stability and reactivity

### 10.1. Reactivity

The product is non-reactive under normal conditions of use, storage and transport.

### 10.2. Chemical stability

Stable under normal conditions.

### 10.3. Possibility of hazardous reactions

No dangerous reactions known under normal conditions of use.

### 10.4. Conditions to avoid

None under recommended storage and handling conditions (see section 7).

### 10.5. Incompatible materials

No additional information available

### 10.6. Hazardous decomposition products

Under normal conditions of storage and use, hazardous decomposition products should not be produced.

## SECTION 11: Toxicological information

### 11.1. Information on toxicological effects

Acute toxicity : Not classified

5-chloro-2-methyl-2H-isothiazol-3-one (2682-20-4)	
ATE US (oral)	100.000 mg/kg body weight
ATE US (dermal)	300.000 mg/kg body weight
2-methyl-2H-isothiazol-3-one (26172-55-4)	
ATE US (oral)	100.000 mg/kg body weight
ATE US (dermal)	300.000 mg/kg body weight

Skin corrosion/irritation : Not classified

Serious eye damage/irritation : Not classified

Respiratory or skin sensitization : Not classified

Germ cell mutagenicity : Not classified

Carcinogenicity : Not classified

Reproductive toxicity : Not classified

Specific target organ toxicity (single exposure) : Not classified

Specific target organ toxicity (repeated exposure) : Not classified

Aspiration hazard : Not classified

## SECTION 12: Ecological information

### 12.1. Toxicity

Ecology - general : The product is not considered harmful to aquatic organisms or to cause long-term adverse effects in the environment.

### 12.2. Persistence and degradability

No additional information available

### 12.3. Bioaccumulative potential

No additional information available

# Sweat Controls for Cystic Fibrosis Testing

## Safety Data Sheet

according to Regulation (EU) 2015/830

### 12.4. Mobility in soil

No additional information available

### 12.5. Results of PBT and vPvB assessment

No additional information available

### 12.6. Other adverse effects

No additional information available

## SECTION 13: Disposal considerations

### 13.1. Waste treatment methods

Waste treatment methods : Dispose of contents/container in accordance with licensed collector's sorting instructions.

## SECTION 14: Transport information

In accordance with ADR / RID / IMDG / IATA / ADN

### 14.1. UN number

Not regulated for transport

### 14.2. UN proper shipping name

Proper Shipping Name (ADR) : Not applicable

Proper Shipping Name (IMDG) : Not applicable

Proper Shipping Name (IATA) : Not applicable

Proper Shipping Name (ADN) : Not applicable

Proper Shipping Name (RID) : Not applicable

### 14.3. Transport hazard class(es)

#### ADR

Transport hazard class(es) (ADR) : Not applicable

#### IMDG

Transport hazard class(es) (IMDG) : Not applicable

#### IATA

Transport hazard class(es) (IATA) : Not applicable

#### ADN

Transport hazard class(es) (ADN) : Not applicable

#### RID

Transport hazard class(es) (RID) : Not applicable

### 14.4. Packing group

Packing group (ADR) : Not applicable

Packing group (IMDG) : Not applicable

Packing group (IATA) : Not applicable

Packing group (ADN) : Not applicable

Packing group (RID) : Not applicable

### 14.5. Environmental hazards

Dangerous for the environment : No

Marine pollutant : No

Other information : No supplementary information available

# Sweat Controls for Cystic Fibrosis Testing

## Safety Data Sheet

according to Regulation (EU) 2015/830

### 14.6. Special precautions for user

#### - Overland transport

No data available

#### - Transport by sea

No data available

#### - Air transport

No data available

#### - Inland waterway transport

No data available

#### - Rail transport

No data available

### 14.7. Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not applicable

## SECTION 15: Regulatory information

### 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

#### 15.1.1. EU-Regulations

Contains no REACH substances with Annex XVII restrictions

Contains no REACH candidate substance

Contains no REACH Annex XIV substances.

### 15.2. Chemical safety assessment

No chemical safety assessment has been carried out

## SECTION 16: Other information

Full text of H-phrases:

Acute Tox. 3 (Dermal)	Acute toxicity (dermal) Category 3
Acute Tox. 3 (Oral)	Acute toxicity (oral) Category 3
Aquatic Acute 1	Hazardous to the aquatic environment - Acute Hazard Category 1
Eye Dam. 1	Serious eye damage/eye irritation Category 1
Skin Corr. 1B	Skin corrosion/irritation Category 1B
Skin Sens. 1	Skin sensitization Category 1
STOT SE 3	Specific target organ toxicity (single exposure) Category 3
H301	Toxic if swallowed
H311	Toxic in contact with skin
H314	Causes severe skin burns and eye damage
H317	May cause an allergic skin reaction
H318	Causes serious eye damage
H335	May cause respiratory irritation
H400	Very toxic to aquatic life

SDS EU Custom – EBS

*This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product.*

Reason for change: Format change to meet GHS compliance.

#42.3

**Catalog Number SS-150**  
GLASS ENCAPSULATED SWEAT CONTROLS FOR CYSTIC FIBROSIS TESTING  
For In Vitro Diagnostic Use

Lot No:  
Expiration Date: 28 February 2009

<b>Analyte / Instrument</b>	<b>LEVEL 1</b> (Source lot 19141) Normal	<b>LEVEL 2</b> (Source lot 19142) High Normal / Equivocal	<b>LEVEL 3</b> (Source lot 19143) Abnormal
<b>Chloride</b> Chloridometer <sup>1</sup>	34 ± 6 mmol/L	68 ± 11 mmol/L	137 ± 21 mmol/L
<b>Conductivity</b> Wescor <sup>®</sup> Sweat•Chek™ 3100 and 3120	42 ± 5 mmol/L	77 ± 6 mmol/L	142 ± 12 mmol/L
<b>Osmolality</b> Wescor <sup>®</sup> Model 5520 Vapor Pressure Osmometer	74 ± 12 mmol/kg	139 ± 12 mmol/kg	281 ± 22 mmol/kg

<sup>1</sup> Mean calculated using data from Labconco Digital Chloridometer (Labconco Corp., Kansas City MO 64132).

### **Intended Use**

Wescor Sweat Controls for Cystic Fibrosis Testing are assayed controls intended as a means of validating the measurement of conductivity<sup>A</sup>, osmolality, or chloride in patient samples.

### **Product Description**

Wescor Sweat Controls for Cystic Fibrosis Testing are liquid, ready to use, and require no reconstitution or dilution. They are supplied in three levels, representing normal, high normal/equivocal, and abnormal electrolyte concentrations. There are 36 ampules per box, 12 x 0.75 mL of each level. They are prepared in an aqueous simulated human sweat matrix<sup>B</sup> to which preservatives, including 0.05% sodium azide have been added.

### **Warnings and Precautions**

The Wescor Sweat Controls for Cystic Fibrosis Testing do not contain any human source material. No special precautions are required in handling the product other than those routinely used in the laboratory. Dispose of properly. Sodium azide may form metal azides in plumbing and pose a threat of explosion.

### **Storage and Stability**

1. The Sweat Controls should be stored at room temperature.
- Each ampule contains 0.75 mL of solution. This volume adequately mitigates evaporative concentration for a few hours after the ampule is broken. Ampules are intended for one-time use only.

### **Instructions for Use**

1. Invert gently to ensure homogeneity of the contents.
2. Before opening, flip the neck of the ampule or tap the bottom of the ampule against a hard surface to clear solution from the neck.
3. Place the ampule in the ampule organizer in the position marked for breaking ampules.
4. Place the safety sleeve over the neck of the ampule.
5. Holding the ampule organizer down against the counter, grasp the neck of the ampule firmly between thumb and forefinger and snap the neck off.
6. Sample directly from the ampule, use a fresh Sweat Chek™ take up tube, micropipettor tip, or other sample loading device each time to avoid contamination of the solution.
7. Treat the controls as you would a patient sample in accordance with the manufacturer's requirements of the test method for each type of analysis: Chloride, Sodium, Conductivity, Osmolality.

### **Expected Values**

The expected values have been established in the laboratory and from interlaboratory data using the listed manufacturer's instruments. Use of control materials having known component concentrations are an integral part of diagnostic procedures. Monitoring of control values establishes intralaboratory parameters for accuracy and precision of the test method. Each laboratory should establish its own quality control standards for the test method used.

### **Limitations**

Any future changes made by the manufacturer of a test method may give a different value from the one indicated. Limitations of the test method are included in the information provided by the instrument manufacturer.

### **References**

<sup>A</sup> Hammond KB, Turcios NL, Gibson LE. Clinical Evaluation of the Macroduct Sweat Collection System and Conductivity Analyzer in the Diagnosis of Cystic Fibrosis. *Pediatrics* 1994; 124: 255-260.

<sup>B</sup> Emrich HM, Stoll E, Colombo JP, Richterrich R, Rossi E: Sweat Composition in Relation to Rate of Sweating in Patients with Cystic Fibrosis of the Pancreas, *Pediatric Research*, 2: 464-478 (1968).

800 453 2725  
435 752 6011  
Fax: 435 752 4127  
www.wescor.com

Wescor, Inc  
459 South Main Street  
Logan, Utah 84321 USA  
DOC4-00014G.doc

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**Priedai ir tiekiamos medžiagos**

AC-071

Ampule Organizer for Sweat Controls

SS-006

Deionized Water, 60 mL bottle, with dropper

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Vaidas Vilmantas (MB „Beikeris“, jm .k. 304539005)

EC DECLARATION OF CONFORMITY

Manufacturer's Name: ELITechGroup Inc. (d.b.a. Wescor, Inc.)
Manufacturer's Address: 370 West 1700 South Logan, Utah 84321 USA
European Authorized Representative: MT Promedt Consulting, GmbH Altenhofstr. 80 D-66386 St. Ingbert / Germany Tel.: +49 6894 - 58 10 20 Fax: +49 6894 - 58 10 21 www.mt-procons.com

Declares, under sole responsibility, that the follow product,

Product Name: Sweat-Chek™
Model Number: 3120
Type of Equipment: AC line-powered, Conductivity-based Sweat Analyzer
GMDN Code: 60268 (Device Only)
IVDD Classification: Other/General...not under List A or B or Self-Testing

Conforms to the requirements of RoHS 2 Directive (2011/65/EU), the applicable section of the European Union Council Directive for InVitro Medical Devices, IVDD 98/79/EC, and the associated Essential Requirements under Annex I. The following standards have been applied:

- EN 1040: 2008
EN 55011:2006 +A1: 2007
EN 61326-1 2006
IEC 61000-4-2:2001, 2008
IEC 61000-4-3:2008
IEC 61000-4-4:2004
IEC 60601-1-2: 2007
IEC 60601-1: 1988, 1991, 1995
EN 55022: 2006 + A1: 2007
IEC 61000-4-6:2006
IEC 61000-4-8:2001
IEC 61000-4-11:2004
IEC 61000-3-2:2005
IEC 61000-3-3:2005
EN 61010-1: 1993, 2001

Additionally, the Accessories and Supplies listed hereafter conform to the requirements in the applicable section of the European Union Council Directive for InVitro Medical Devices, IVDD 98/79/EC, and the associated Essential Requirements under Annex I.

This declaration applies to the instrument itself and to the accessories and supplies listed hereafter.
Accessories and Supplies

- AC-071 Ampule Organizer for Sweat Controls
SS-006 Deionized Water, 60 mL bottle, with dropper
SS-044 Take-Up Tubes, 14 inch, package of 100
SS-045 Syringe/Needle Set, 3 ea: 1 cc & #22 blunt needle
SS-107 Small sealable containers, 6 each
SS-140 Sweat-Chek Calibrator, 0.75 mL, 90 mmol/L, pk/60
SS-150 Sweat Controls; L1, L2, L3 12ea; 36 0.75 mL ampules

And as such bears the CE and IVD markings.

A technical file is maintained at ELITechGroup Inc. Biomedical System Division located at 370 West 1700 South, Logan, Utah, 84321, USA and a copy is provided to our European Authorized Representative.

Handwritten signature of Dawn T. Perdue, Date: 22 Dec 2016, Location: Logan, UT U.S.A.
Dawn T. Perdue
Director of Quality Assurance

EC DECLARATION OF CONFORMITY

BIOMEDICAL SYSTEMS
370 West 1700 South
Logan, Utah 84321 USA
435 752 6011
800 453 2725
fax 435 752 4127
web www.elitechgroup.com

Manufacturer's Name: ELITechGroup Inc. (d.b.a. Wescor, Inc.)
Manufacturer's Address: 370 West 1700 South Logan, Utah 84321 USA
European Authorized Representative: MT Promedt Consulting, GmbH
Altenhofstr. 80
D-66386 St. Ingbert / Germany
Tel.: +49 6894 - 58 10 20
Fax: +49 6894 - 58 10 21
www.mt-procons.com

Declares, under sole responsibility, that the following Class I items that are part of the Macroduct Supply Kit (SS-032),

- Item Number: SS-107 Small Sealable Containers
GMDN code: 16291
Item Number: SS-142 Macroduct Collector
GMDN Code: 60298

Type of Equipment: Supplies used to collect and transport sweat samples for the Cystic Fibrosis (CF) analysis with Model 3700-SYS Macroduct®.

Conforms to the applicable sections of the European Union Council Directive concerning medical devices, MDD 93/42/EEC, and is classified as a Class I medical device and therefore bears the CE Marking as a Class I medical device.

As such, this item bears the CE marking.

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