A, B, DVI-, cti/A₁, B**Determination of the ABO/Rh blood groups combined with reverse grouping****Product-Identification: 50092****INTRODUCTION**

ABO blood group typing, using anti-A and anti-B test sera, is known as direct or forward grouping test.

Reverse grouping uses red cell reagents of known ABO antigen specificity to indicate the presence or absence of anti-A and anti-B isoagglutinins, the results of which determine the reverse group. Discrepancies between forward and reverse grouping require further investigation [1].

Classification of blood groups must be based on both forward and reverse grouping.

The ID-Card "DiaClon ABO/D + Reverse Grouping" allows combined testing of forward and reverse grouping as well as RhD determination.

REAGENTS

ID-Card "DiaClon ABO/D + Reverse Grouping" contains monoclonal anti-A [cell line A5], anti-B [cell line G½] and anti-D [cell lines LHM 59 / 20 (LDM3) + 175-2] within the gel matrix. The microtube cti is the negative control. Two microtubes with "neutral" gel serve for reverse grouping with A₁ and B cells. Preservative: < 0.1% NaN₃.

Caution: All reagents should be treated as potentially infectious.

Do not store near any heat, air conditioning sources or ventilation outlets.

Stability: see expiry date on label.

ADDITIONAL REAGENTS REQUIRED

- ID-Diluent 2: modified LISS for red cell suspensions.
- Test cell Reagents: ID-DiaCell A₁ and B in a 0.8% suspension, in 10 ml vials, ready-to-use.

(see related package insert)

FURTHER MATERIALS REQUIRED

- ID-Dispenser
- ID-Pipetor
- ID-Tips (pipetor tips)
- Tubes for suspensions
- ID-Working table
- ID-Centrifuge 6, 12 or 24

SAMPLE MATERIAL

For optimal results, the determination should be performed using a freshly drawn sample, or in accordance with local laboratory procedures for sample acceptance criteria. Preferably, blood samples should be drawn into citrate, EDTA or CPD-A anticoagulant. Samples drawn into plain tubes (no anticoagulant) may also be used.

When the use of serum instead of plasma is required, the serum must be well cleared, by centrifugation at 1,500 g for 10 minutes, before use avoid fibrin residues, which may interfere with the reaction pattern.

PREPARATION OF BLOOD SAMPLE**a) Red cell suspension (for ABO/D determination)**

Prepare a 5% red cell suspension in ID-Diluent 2 as follows:
Allow the diluent to reach room temperature before use.

1. Dispense 0.5 ml of ID-Diluent 2 into a clean tube.
2. Add 50 µl of whole blood or 25 µl of packed cells, mix gently.
The cell suspension may be used immediately.

b) Plasma or Serum for reverse grouping

Where samples are not for immediate testing they should be stored at 2–8 °C after separation for a maximum of 48 hours, thereafter at -20 °C, or in accordance with local / national policies / guidelines.

CONTROLS

Known positive and negative samples should be included in accordance with the relevant guidelines of quality assurance.

TEST PROCEDURE

Do not use ID-Cards which show signs of drying, have bubbles, damaged seals, drops of gel or supernatant in the upper part of the microtubes or on the underside of the aluminium foil.

Allow the test cell reagent to reach room temperature before use.

1. Identify the ID-Card with the unique patient or donor number / details as appropriate.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Pipette 50 µl of ID-DiaCell A₁ to microtube 5 (A₁).
4. Pipette 50 µl of ID-DiaCell B to microtube 6 (B).
5. Pipette 50 µl of the patient serum or plasma to both microtubes 5 and 6. An incubation for 10 minutes at room temperature is recommended (see remarks, point 4).
6. Pipette 10 or 12.5 µl of the patient's red cell suspension to the microtubes 1-4 (A, B, D, ctl).
7. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
8. Read and record the reactions.

INTERPRETATION OF THE RESULTS

A) Principle [2]

Positive: Agglutinated cells forming a red line on the surface of the gel or agglutinates dispersed in the gel.
Negative: Compact button of cells on the bottom of the microtube.

B) Reactions for blood groups ABO

Anti-A	Anti-B	Blood group
+++ to +++++	negative	A
negative	+++ to +++++	B
+++ to +++++	+++ to +++++	AB
negative	negative	O

Weaker reactions than +++ may indicate A or B subgroups. For correct interpretation, a complete grouping test should be performed (anti-A, anti-B, anti-AB). In the presence of weak or very weakly expressed antigens the reaction can be negative. The anti-B of monoclonal origin does not react with the acquired B antigen.

Important: The microtube ctl must show a negative reaction. If the ctl is positive, the ABO determination is not valid. Repeat the test as described under "Remarks 1."

C) Reactions for reverse grouping

A ₁	B	Blood group
+ to +++++	negative	B
negative	+ to +++++	A
+ to +++++	+ to +++++	O
negative	negative	AB

If questionable reactions are obtained, repeat reverse grouping with 4 red cell reagents (A₁, A₂, B and O).

D) Reactions for RhD

+++ to +++++	± to ++*	negative
RhD positive	RhD weak positive	RhD negative

* ±, trace or weak reactions should be subject to further investigations to distinguish between weak and partial D types as appropriate for the category of sample being tested.

The anti-D used has been selected so as not to react with DVI variants.

Weak D may give a negative reaction. If all weak / partial D's are required to be detected, all D negative results must be retested. Note that most guidelines do not recommend further testing for weak or partial-D in patients.

Important: The microtube ctl must show a negative reaction. If the ctl is positive, the RhD determination is not valid. Repeat the test as described under "Remarks 1."

REMARKS

1. The negative control must always show a negative reaction.
 - If the negative control is positive, wash the red cells first in warm isotonic saline solution or ID-Diluent 2, before preparing the red cell suspension.
 - Then proceed as under "Preparation of blood sample" and "Test procedure".
 - If the negative control subsequently shows a negative result, the reactions can be interpreted as described in sections B, C and D.
 - If the negative control remains positive, the results of the ABO/Rh determination should be considered invalid and further investigations following recommended techniques should be undertaken to ascertain the reason, before valid antigen typing can be assured.
2. Where a discrepancy occurs between the results of ABO typing and reverse grouping, consult the "DiaMed ABO discrepancy chart" for appropriate information.
3. Full forward and reverse grouping requires the use of anti-A, -B, -AB, and A₁, A₂, B and O cells. The ID-Card "DiaClon ABO/D + Reverse Grouping" does not contain anti-AB and allows the use of A₁ and B cells only. It should be used only for confirmation of previously fully tested samples with established blood group status or in accordance with local or national guidelines / recommendations for blood grouping.
4. For reverse grouping, an incubation of at least 10 minutes at 18–25 °C prior to centrifugation will enhance the reactions and minimise repeat testing due to weak isoagglutinins.

LIMITATIONS

- a) ID-Cards which show air bubbles in the gel or drops in the upper part of the microtubes and / or the seal, must be centrifuged before use.
- b) Bacterial or other contamination of materials used can cause false positive or false negative results.
- c) Fibrin residues in the serum or the red cell suspension may trap non-agglutinated cells presenting a fine pink line on top of the gel while most of the cells are on the bottom of the microtube after centrifugation.
- d) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.
- e) Use of suspension solutions other than ID-Diluent 2 may modify the reactions.
- f) Too heavy or too weak red cell suspensions can cause aberrant results.

BIBLIOGRAPHY

1. Mollison, P.L., Engelfriet, C.P. and Contreras, M.: Blood Transfusion in Clinical Medicine, 10th ed. 1997; Blackwell Scientific Publications, Oxford.
2. Lapiere, Y., Rigal, D., Adam, J. et al.: The gel test; A new way to detect red cell antigen-antibody reactions. Transfusion 1990; 30: 109–113.

PRODUCTS

ID-Card "DiaClon ABO/D + Reverse Grouping"	4 x 12	REF 001234
	24 x 12	REF 001237
	60 x 12	REF 001236
	112 x 12	REF 001235

These products are guaranteed to perform as described on the label and in the instruction sheet. The manufacturer declines all responsibility arising out of the use or sale of these products in any way or for any purpose other than those described therein.

Changes to the version 07.11 are shaded grey.

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Declaration of conformity

MANUFACTURER: **DiaMed GmbH**
ADDRESS: **Pra Rond 23**
1785 Cressier FR
Switzerland
Phone: **+41(0)26 674 51 11**

PRODUCT NAME: **DiaClon ABO/D + Reverse Grouping**
Id-n°: **50092**
REF: **001234 / 001237 / 001236 / 001235**

We hereby declare that the above mentioned product meets the provisions of the following Directives:

APPLICABLE DIRECTIVE: Directive 98/79/EC of the European Parliament and of the Council of 27 October 1998 on in vitro Diagnostic medical devices

CLASSIFICATION: Annex II List A

CONFORMITY ROUTE: Annex IV

GMDN Code: 45308

Generic Device Group Term: ABO/Rh(D) multiple blood grouping IVD, kit, agglutination

NOTIFIED BODY: TÜV Product Service GmbH
Ridlerstrasse 65-80339 München (Germany)
CE-N° 0123

Name:	Function:	Issued in:	Date:	Signature
Galéa Diane	Site Quality Management Representative	Cressier FR	11.12.2017	

DiaMed GmbH
Pra Rond 23
1785 Cressier FR / Šveicarija
Telefonas: +41 (0)26 674 51 11
Faksas: +41 (0)26 674 54 45

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Version number 3

Revision: 21.08.2018

SECTION 1: Identification of the substance/mixture and of the company/ undertaking

- **1.1 Product identifier**
 - Trade name: **DiaMed GmbH IVD MD**
 - IHD Id-N°
10270 / 10291 / 10301 / 10600 / 10610 / 10620 / 10630 / 10640 / 10650 / 11270 /
11280 / 12300 / 12310 / 12320 / 12330 / 12340 / 12350 / 14060 / 14070 / 17030 /
17141 / 17230 / 17323 / 17420 / 17620 / 17720 / 45470 / 45480 / 46130 / 50012 /
50053 / 50057 / 50071 / 50092 / 50093 / 50110 / 50115 / 50171 / 50200 / 50212 /
50221 / 50242 / 50250 / 50411 / 50481 / 50492 / 50601 / 50682 / 50710 / 50741 /
50742 / 50850 / 50870 / 50961 / 50981 / 51011 / 51051 / 51090 / 51160 / 51190 /
51210 / 51241 / 51250 / 60300 / 60400 / 61470 / 61480 / 62380 / 63390 / 64370 /
65570 / 08610 / 08620 / 09410 / 52000 / 52010 / 52020 / 52030 / 52040 / 52050
 - Bio-Rad MSDS Number: 2704M
- **1.2 Relevant identified uses of the substance or mixture and uses advised against**
No further relevant information available.
 - **Application of the substance / the mixture**
In vitro diagnostic medical device or component.
- **1.3 Details of the supplier of the safety data sheet**
 - **Manufacturer/Supplier:**
DiaMed GmbH
Pra Rond 23
CH-1785 Cressier FR
(Switzerland/Schweiz/Suisse/Svizzera)
Tel: +41 (0)26 674 51 11
Fax: +41 (0)26 674 51 45
 - **Further information obtainable from:** fds-msds.ch@bio-rad.com
- **1.4 Emergency telephone number:** CHEMTREC UK (London) +(44)-870-8200418

SECTION 2: Hazards identification

- **2.1 Classification of the substance or mixture**
 - **Classification according to Regulation (EC) No 1272/2008**
The product is not classified according to the CLP regulation.
- **2.2 Label elements**
 - **Labelling according to Regulation (EC) No 1272/2008** Void
 - **Hazard pictograms** Void
 - **Signal word** Void
 - **Hazard statements** Void
- **2.3 Other hazards**
 - **Results of PBT and vPvB assessment**
 - **PBT:** Not applicable.
 - **vPvB:** Not applicable.

SECTION 3: Composition/information on ingredients

- **3.2 Chemical characterisation: Mixtures**
 - **Description:** Mixture of substances listed below with nonhazardous additions.
 - **Dangerous components:** Void
 - **Additional information:**
Caution: All reagents should be treated as potentially infectious.
For the wording of the listed hazard phrases refer to section 16.

GB

(Contd. on page 2)

Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 21.08.2018

Version number 3

Revision: 21.08.2018

Trade name: DiaMed GmbH IVD MD

(Contd. of page 1)

SECTION 4: First aid measures

- **4.1 Description of first aid measures**
 - **General information:** No special measures required.
 - **After inhalation:** Supply fresh air; consult doctor in case of complaints.
 - **After skin contact:** Generally the product does not irritate the skin.
 - **After eye contact:** Rinse opened eye for several minutes under running water.
 - **After swallowing:** If symptoms persist consult doctor.
- **4.2 Most important symptoms and effects, both acute and delayed**
No further relevant information available.
- **4.3 Indication of any immediate medical attention and special treatment needed**
No further relevant information available.

SECTION 5: Firefighting measures

- **5.1 Extinguishing media**
 - **Suitable extinguishing agents:**
CO₂, powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
Use fire extinguishing methods suitable to surrounding conditions.
- **5.2 Special hazards arising from the substance or mixture**
No further relevant information available.
- **5.3 Advice for firefighters**
 - **Protective equipment:** No special measures required.

SECTION 6: Accidental release measures

- **6.1 Personal precautions, protective equipment and emergency procedures** Not required.
- **6.2 Environmental precautions:**
Dilute with plenty of water.
Do not allow product to reach sewage system or any water course.
- **6.3 Methods and material for containment and cleaning up:**
Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
- **6.4 Reference to other sections**
See Section 7 for information on safe handling.
See Section 8 for information on personal protection equipment.
See Section 13 for disposal information.

SECTION 7: Handling and storage

- **7.1 Precautions for safe handling** No special measures required.
 - **Information about fire - and explosion protection:** No special measures required.
- **7.2 Conditions for safe storage, including any incompatibilities**
 - **Storage:**
 - **Requirements to be met by storerooms and receptacles:** No special requirements.
 - **Information about storage in one common storage facility:** Not required.
 - **Further information about storage conditions:** see related package insert.
- **7.3 Specific end use(s)** No further relevant information available.

SECTION 8: Exposure controls/personal protection

- **Additional information about design of technical facilities:** No further data; see item 7.

(Contd. on page 3)

GB

Safety data sheet
according to 1907/2006/EC, Article 31

Printing date 21.08.2018

Version number 3

Revision: 21.08.2018

Trade name: DiaMed GmbH IVD MD

(Contd. of page 2)

- **8.1 Control parameters**
 - **Ingredients with limit values that require monitoring at the workplace:**
 The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
 - **Additional information:** The lists valid during the making were used as basis.
- **8.2 Exposure controls**
 - **Personal protective equipment:**
 - **General protective and hygienic measures:**
 The usual precautionary measures are to be adhered to when handling chemicals.
 - **Respiratory protection:** Not required.
 - **Protection of hands:**

Protective gloves

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.

- **Material of gloves**
 The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.
- **Penetration time of glove material**
 The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

· **Eye protection:**

Safety glasses

· **Body protection:** Protective work clothing

SECTION 9: Physical and chemical properties

- **9.1 Information on basic physical and chemical properties**
 - **General Information**
 - **Appearance:**
 - **Form:** Fluid
 - **Colour:** Colourless
 - **Odour:** Odourless
 - **Odour threshold:** Not determined.
 - **pH-value:** Not determined.
 - **Change in condition**
 - **Melting point/freezing point:** Undetermined.
 - **Initial boiling point and boiling range:** 100 °C
 - **Flash point:** Not applicable.
 - **Flammability (solid, gas):** Not applicable.

(Contd. on page 4)

GB

Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 21.08.2018

Version number 3

Revision: 21.08.2018

Trade name: DiaMed GmbH IVD MD

(Contd. of page 3)

· Ignition temperature:	
· Decomposition temperature:	Not determined.
· Auto-ignition temperature:	Product is not selfigniting.
· Explosive properties:	Product does not present an explosion hazard.
· Explosion limits:	
· Lower:	Not determined.
· Upper:	Not determined.
· Vapour pressure at 20 °C:	23hPa
· Density:	Not determined.
· Relative density	Not determined.
· Vapour density	Not determined.
· Evaporation rate	Not determined.
· Solubility in / Miscibility with	
· water:	Fully miscible.
· Partition coefficient: n-octanol/water:	Not determined.
· Viscosity:	
· Dynamic:	Not determined.
· Kinematic:	Not determined.
· Solvent content:	
· Water:	50.0%
· Solids content:	0.0%
· 9.2 Other information	No further relevant information available.

SECTION 10: Stability and reactivity

- **10.1 Reactivity** No further relevant information available.
- **10.2 Chemical stability**
 - **Thermal decomposition / conditions to be avoided:**
No decomposition if used according to specifications.
- **10.3 Possibility of hazardous reactions** No dangerous reactions known.
- **10.4 Conditions to avoid** No further relevant information available.
- **10.5 Incompatible materials:** No further relevant information available.
- **10.6 Hazardous decomposition products:** No dangerous decomposition products known.

SECTION 11: Toxicological information

- **11.1 Information on toxicological effects**
 - **Acute toxicity** Based on available data, the classification criteria are not met.
 - **LD/LC50 values relevant for classification:**
Sodium Azid is used as a preservative. (Concentration <0.1%)

26628-22-8 Sodium Azide

Oral	LD50	27 mg/kg (rat)
Dermal	LD50	20 mg/kg (rabbit)

- **Primary irritant effect:**
 - **Skin corrosion/irritation**
Based on available data, the classification criteria are not met.

(Contd. on page 5)

GB

Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 21.08.2018

Version number 3

Revision: 21.08.2018

Trade name: **DiaMed GmbH IVD MD**

(Contd. of page 4)

- **Serious eye damage/irritation**
Based on available data, the classification criteria are not met.
- **Respiratory or skin sensitisation**
Based on available data, the classification criteria are not met.
- **Additional toxicological information:** No further relevant information available.
- **CMR effects (carcinogenity, mutagenicity and toxicity for reproduction)**
 - **Germ cell mutagenicity**
Based on available data, the classification criteria are not met.
 - **Carcinogenicity** Based on available data, the classification criteria are not met.
 - **Reproductive toxicity**
Based on available data, the classification criteria are not met.
- **STOT-single exposure** Based on available data, the classification criteria are not met.
- **STOT-repeated exposure**
Based on available data, the classification criteria are not met.
- **Aspiration hazard** Based on available data, the classification criteria are not met.

SECTION 12: Ecological information

- **12.1 Toxicity**
 - **Aquatic toxicity:** No further relevant information available.
- **12.2 Persistence and degradability** No further relevant information available.
- **12.3 Bioaccumulative potential** No further relevant information available.
- **12.4 Mobility in soil** No further relevant information available.
 - **Additional ecological information:**
 - **General notes:** Generally not hazardous for water
- **12.5 Results of PBT and vPvB assessment**
 - **PBT:** Not applicable.
 - **vPvB:** Not applicable.
- **12.6 Other adverse effects** No further relevant information available.

SECTION 13: Disposal considerations

- **13.1 Waste treatment methods**
 - **Recommendation** Smaller quantities can be disposed of with household waste.
- **Uncleaned packaging:**
 - **Recommendation:** Disposal must be made according to official regulations.
 - **Recommended cleansing agents:**
Water, if necessary together with cleansing agents.

SECTION 14: Transport information

- | | | |
|--|---|------|
| · 14.1 UN-Number | · ADR, ADN, IMDG, IATA | Void |
| · 14.2 UN proper shipping name | · ADR, ADN, IMDG, IATA | Void |
| · 14.3 Transport hazard class(es) | · ADR, ADN, IMDG, IATA
· Class | Void |

(Contd. on page 6)

GB

Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 21.08.2018

Version number 3

Revision: 21.08.2018

Trade name: DiaMed GmbH IVD MD

(Contd. of page 5)

· 14.4 Packing group · ADR, IMDG, IATA	Void
· 14.5 Environmental hazards: · Marine pollutant:	No
· 14.6 Special precautions for user	Not applicable.
· 14.7 Transport in bulk according to Annex II of Marpol and the IBC Code	Not applicable.
· Transport/Additional information:	Not dangerous according to the above specifications.
· UN "Model Regulation":	Void

SECTION 15: Regulatory information

- **15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture**
 - **Labelling according to Regulation (EC) No 1272/2008** Void
 - **Hazard pictograms** Void
 - **Signal word** Void
 - **Hazard statements** Void
 - **Directive 2012/18/EU**
 - **Named dangerous substances - ANNEX I** None of the ingredients is listed.
- **15.2 Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

SECTION 16: Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

- **Department issuing SDS:**
DiaMed GmbH
Pra Rond 23
CH-1785 Cressier FR
(Switzerland/Schweiz/Suisse/Svizzera)
 - **Abbreviations and acronyms:**
ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)
IMDG: International Maritime Code for Dangerous Goods
IATA: International Air Transport Association
GHS: Globally Harmonised System of Classification and Labelling of Chemicals
EINECS: European Inventory of Existing Commercial Chemical Substances
ELINCS: European List of Notified Chemical Substances
CAS: Chemical Abstracts Service (division of the American Chemical Society)
LC50: Lethal concentration, 50 percent
LD50: Lethal dose, 50 percent
PBT: Persistent, Bioaccumulative and Toxic
vPvB: very Persistent and very Bioaccumulative
 - *** Data compared to the previous version altered.**

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#46.1.2

B009290 03.13

Modified LISS for red cell suspensions

Product-Identification: 05761

INTRODUCTION

Low ionic strength solution (LISS) increases the rate of antibody association and thus enhances antigen/antibody reactions.

ID-Diluent 2[®] is a modified low ionic strength solution made for the ID-System, for preparing 5% red cell suspensions for blood grouping as well as 0.8% red cell suspensions for crossmatching, autocontrol, direct anti-human globulin test, blood grouping of newborns and test cells prepared in the laboratory.


REAGENTS

ID-Diluent 2[®]: modified LISS for Red Cell Suspension, in 100 and 500 ml bottles or racks of 60 x 700 µl for the IH-Analyzers.

Preservatives: the antibiotics trimethoprim and sulfamethoxazole.

Once opened and if handled in accordance with good laboratory practice (GLP) principles and stored correctly, the 100 and 500 ml bottles can be used until the date of expiry or for up to 6 months afterwards, whichever is the sooner.

Stability: see expiry date on label.

FURTHER MATERIALS REQUIRED (FOR 100 AND 500 ML BOTTLES)

- ID-Dispenser
- Pipette 25 and 50 µl
- ID-Tips (pipetor tips)
- Suspension Tubes
- ID-Working Table

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PREPARATION OF BLOOD SAMPLE (FOR 100 AND 500 ML BOTTLES)

A) For blood group determinations

Prepare a 5% red cell suspension of the patient's red cells in ID-Diluent 2[®] as follows:

Allow the diluent to reach room temperature before use.

1. Dispense 0.5 ml of ID-Diluent 2[®] into a clean tube
2. Add 50 µl of whole blood or 25 µl of packed cells, mix gently.

The cell suspension may be used immediately.

B) For crossmatching, direct anti-human globulin test, autocontrol in antibody screening and other tests as specified in the relevant box inserts

Prepare a 0.8% red cell suspension in ID-Diluent 2[®] as follows:

Allow the diluent to reach room temperature before use.

1. Dispense 1.0 ml of ID-Diluent 2[®] into a clean tube.
2. Add 10 µl of packed red cells, mix gently.

The cell suspension may be used immediately.

C) For preparation of test cell reagents in the laboratory

1. Wash the red cells with isotonic saline solution or with ID-Diluent 2[®] three times or more, until the supernatant is clear.
2. Decant the supernatant and resuspend the packed cells to a 0.8% suspension with ID-Diluent 2[®].
3. Under aseptic conditions place the suspension into sterile glass vials or tubes, close with sterile caps.

Red cell reagents thus prepared and stored at 2–8 °C are stable for one day.

Important: consult the related instruction sheet of the ID-Cards for precise working procedures.

PREPARATION OF BLOOD SAMPLE (RACKS WITH 60 X 700 µL)

Please refer to the user manual of the IH-Analyzers.

LIMITATIONS

- a) Bacterial or other contamination of materials used can cause false positive or false negative results.
- b) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.
- c) "ID-Diluent 2" has been formulated in order to reduce the incidence of detection of non-significant red cell bound complement components.
- d) Some antibodies in the Kell system may react more weakly in LISS techniques.
- e) LISS solutions may affect autoantibody activity and therefore may cause problems with certain samples.

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Changes to the version 08.12 are shaded grey.

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#46.1.2

Declaration of conformity

MANUFACTURER: **DiaMed GmbH**
ADDRESS: **Pra Rond 23**
1785 Cressier FR
Switzerland
Phone: +41(0)26 674 51 11

PRODUCT NAME: **ID-Diluent 2**
Id-n°: **05761**
REF: **009260 / 009280 / 009290**

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NOTIFIED BODY: /

Name:	Function:	Issued in:	Date:	Signature
Galéa Diane	Site Quality Management Representative	Cressier FR	01.02.2017	

DiaMed GmbH
Pra Rond 23
1785 Cressier FR / Šveicarija
Telefonas: +41 (0)26 674 51 11
Faksas: +41 (0)26 674 54 45

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- **1.1 Product identifier**
 - Trade name: **DiaMed GmbH IVD MD**
 - IHD Id-N°
12060 / 12140 / 20100 / 40230 / 40240 / 40560 / 45270 / 46140 / 50560 / 50620 /
50830 / 50900 / 05710 / 05730 / 05740 / 05751 / 05761 / 05780 / 06311 / 06760 /
06940 / 09550
 - Bio-Rad MSDS Number: 2702M
- **1.2 Relevant identified uses of the substance or mixture and uses advised against**
No further relevant information available.
 - **Application of the substance / the mixture**
In vitro diagnostic medical device or component.
- **1.3 Details of the supplier of the safety data sheet**
 - **Manufacturer/Supplier:**
DiaMed GmbH
Pra Rond 23
CH-1785 Cressier FR
(Switzerland/Schweiz/Suisse/Svizzera)
Tel: +41 (0)26 674 51 11
Fax: +41 (0)26 674 51 45
 - **Further information obtainable from:** fds-msds.ch@bio-rad.com
- **1.4 Emergency telephone number:** CHEMTREC UK (London) +(44)-870-8200418

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- **3.2 Chemical characterisation: Mixtures**
 - **Description:** Mixture of substances listed below with nonhazardous additions.
 - **Dangerous components:** Void
 - **Additional information:** For the wording of the listed hazard phrases refer to section 16.

SECTION 4: First aid measures

- **4.1 Description of first aid measures**
 - **General information:** No special measures required.
 - **After inhalation:** Supply fresh air; consult doctor in case of complaints.
 - **After skin contact:** Generally the product does not irritate the skin.
 - **After eye contact:** Rinse opened eye for several minutes under running water.

(Contd. on page 2)

GB

Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 14.02.2018

Version number 2

Revision: 14.02.2018

Trade name: DiaMed GmbH IVD MD

(Contd. of page 1)

- **After swallowing:** If symptoms persist consult doctor.
- **4.2 Most important symptoms and effects, both acute and delayed**
No further relevant information available.
- **4.3 Indication of any immediate medical attention and special treatment needed**
No further relevant information available.

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- **Additional information about design of technical facilities:** No further data; see item 7.
- **8.1 Control parameters**
 - **Ingredients with limit values that require monitoring at the workplace:**
The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
 - **Additional information:** The lists valid during the making were used as basis.

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· **8.2 Exposure controls**

· **Personal protective equipment:**

· **General protective and hygienic measures:**

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· **Penetration time of glove material**

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· **9.1 Information on basic physical and chemical properties**

· **General Information**

· **Appearance:**

· **Form:**

Not determined.

· **Colour:**

Not determined.

· **Odour:**

Odourless

· **Odour threshold:**

Not determined.

· **pH-value:**

Not determined.

· **Change in condition**

· **Melting point/freezing point:**

Undetermined.

· **Initial boiling point and boiling range:**

Undetermined.

· **Flash point:**

Not applicable.

· **Flammability (solid, gas):**

Not applicable.

· **Ignition temperature:**

· **Decomposition temperature:**

Not determined.

· **Auto-ignition temperature:**

Product is not selfigniting.

· **Explosive properties:**

Product does not present an explosion hazard.

(Contd. on page 4)

GB

Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 14.02.2018

Version number 2

Revision: 14.02.2018

Trade name: DiaMed GmbH IVD MD

(Contd. of page 3)

<ul style="list-style-type: none"> · Explosion limits: · Lower: Not determined. · Upper: Not determined.
<ul style="list-style-type: none"> · Vapour pressure: Not determined.
<ul style="list-style-type: none"> · Density: Not determined. · Relative density: Not determined. · Vapour density: Not determined. · Evaporation rate: Not determined.
<ul style="list-style-type: none"> · Solubility in / Miscibility with · water: Fully miscible.
<ul style="list-style-type: none"> · Partition coefficient: n-octanol/water: Not determined.
<ul style="list-style-type: none"> · Viscosity: · Dynamic: Not determined. · Kinematic: Not determined.
<ul style="list-style-type: none"> · Solvent content: · Water: 1.0%
<ul style="list-style-type: none"> · Solids content: 0.0%
<ul style="list-style-type: none"> · 9.2 Other information No further relevant information available.

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SECTION 11: Toxicological information

- **11.1 Information on toxicological effects**
 - **Acute toxicity** Based on available data, the classification criteria are not met.
 - **Primary irritant effect:**
 - **Skin corrosion/irritation**
Based on available data, the classification criteria are not met.
 - **Serious eye damage/irritation**
Based on available data, the classification criteria are not met.
 - **Respiratory or skin sensitisation**
Based on available data, the classification criteria are not met.
- **Additional toxicological information:** No further relevant information available.
 - **CMR effects (carcinogenicity, mutagenicity and toxicity for reproduction)**
 - **Germ cell mutagenicity**
Based on available data, the classification criteria are not met.
 - **Carcinogenicity** Based on available data, the classification criteria are not met.
 - **Reproductive toxicity**
Based on available data, the classification criteria are not met.
 - **STOT-single exposure** Based on available data, the classification criteria are not met.

(Contd. on page 5)

GB

Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 14.02.2018

Version number 2

Revision: 14.02.2018

Trade name: DiaMed GmbH IVD MD

(Contd. of page 4)

- **STOT-repeated exposure**
Based on available data, the classification criteria are not met.
- **Aspiration hazard** Based on available data, the classification criteria are not met.

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- **Department issuing SDS:**
 DiaMed GmbH
 Pra Rond 23
 CH-1785 Cressier FR
 (Switzerland/Schweiz/Suisse/Svizzera)
- **Abbreviations and acronyms:**
 ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)
 IMDG: International Maritime Code for Dangerous Goods
 IATA: International Air Transport Association
 GHS: Globally Harmonised System of Classification and Labelling of Chemicals
 EINECS: European Inventory of Existing Commercial Chemical Substances
 ELINCS: European List of Notified Chemical Substances
 CAS: Chemical Abstracts Service (division of the American Chemical Society)
 PBT: Persistent, Bioaccumulative and Toxic
 vPvB: very Persistent and very Bioaccumulative
- *** Data compared to the previous version altered.**

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Test cells for reverse grouping

INTRODUCTION

Test cell reagents are routinely used in blood group serology to detect the presence or absence of anti-A and anti-B isoagglutinins (in reverse grouping).

For reverse grouping, test cell reagents of group A₁ and B / A₁, A₂ and B / A₁, B and O or A₁, A₂, B and O are used, following various requirements and guidelines.

The test cell reagents are specially designed for the ID-System.

REAGENTS


IVD

All test cell reagents are of human origin, in a buffered suspension medium at 0.8% (\pm 0.1%).
Preservatives: the antibiotics trimethoprim and sulfamethoxazole.

ID-DiaCell ABO: A₁, A₂, B, O / A₁, A₂, B / A₁, B, O / A₁, B / A₁ / A₂ / B / O

Shipment on standing order, every 4 weeks.

Caution: The source materials from which these products were manufactured, were found non-reactive for HBsAg, HCV and HIV (1+2) when tested with licensed reagents. However, no known test method can assure that infectious agents are absent. Products from human blood should be considered potentially infectious.

Stability: see expiry date on label.

ADDITIONAL REAGENTS REQUIRED

- ID-Card "NaCl, Enzyme Test and Cold Agglutinins" with 6 microtubes containing neutral gel (Id-n°: 50520).
- ID-Card "Reverse Grouping with Antibody Screening" with 3 microtubes containing neutral gel and 3 microtubes containing polyspecific anti-human globulin (AHG) serum (Id-n°: 50510).
- ID-Card "ABO/D + Reverse Grouping", A, B, D, ctl / A₁, B (Id-n°: 50081).
- ID-Card "DiaClon ABO/D + Reverse Grouping", A, B, DVI-, ctl / A₁, B (Id-n°: 50092).

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CONTROLS

Controls should be included in accordance with the relevant guidelines of quality assurance.

USE OF THE ID-TEST CELL REAGENTS

- All Test cell reagents are for use with the ID-Cards of the **ID-System** only.
- **Strictly follow the test procedures as described in the specific package inserts of the ID-Cards to be used.**
- Always gently resuspend the red cells, by inverting the vial several times before use and also before placing the vials into a pipetting automat.
- Make sure that the test cells are at room temperature (18–25 °C) when in use.
- During the working procedures, check that the test cell reagents remain in suspension. If there is settling of the cells, resuspend again.
- **For the ID-System, precise pipetting is of importance.** Use the ID-Pipetors for serial pipetting.
- Avoid contamination of the test cell reagents.
- After use, close the vials and replace them in the refrigerator.

INTERPRETATION OF THE RESULTS

A) Principle

Agglutination indicate the presence of the agglutinins.

No agglutination indicate the absence of agglutinins.

B) Reaction with ID-DiaCell ABO

A ₁	A ₂	B	O	Isoagglutinins	Blood groups
-	-	+	-	anti-B	A
+	+	-	-	anti-A	B
-	-	-	-	none	AB
+	+	+	-	anti-A and anti-B	O

Note: The results of the reverse grouping must agree with the blood group observed from forward (antigen) typing. Discrepant results should be subjected to further investigation, preferably with a new sample.

REMARKS

1. The optimal reaction temperature for isoagglutinins is 4 °C. If weak or doubtful reactions are observed, repeat the test with an incubation at 2–8 °C for 15 minutes.
2. Atypical reactions require further studies.
3. In general, anti-A and anti-B iso-antibodies are not detectable in newborns and infants. Older people and patients with agammaglobulinaemia may also lack these iso-antibodies.

LIMITATIONS

- a) Bacterial or other contamination of materials used can cause false positive or false negative results.
- b) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.

BIBLIOGRAPHY

1. Technical Manual of the American Association of Blood Banks, 13th edition, 1999.
2. Lapiere, Y., Rigal, D., Adam, J. et al.: The gel test; A new way to detect red cell antigen-antibody reaction. Transfusion 1990; 30:109-113.

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LITERATŪRA

- Technical Manual of the American Association of Blood Banks, 13th edition, 1999.
- Lapierre, Y., Rigal, D., Adam, J. et al.: The gel test; A new way to detect red cell antigen-antibody reaction. Transfusion 1990; 30:109-113.

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#46.1.3

Declaration of conformity

MANUFACTURER: DiaMed GmbH
ADDRESS: Pra Rond 23
1785 Cressier FR
Switzerland
Phone: +41(0)26 674 51 11

PRODUCT NAME: ID-DiaCell ABO A₁, B
Id-n°: 45092
REF: 003624

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CONFORMITY ROUTE: Annex IV

GMDN Code: 52691

Generic Device Group Term: Reverse ABO grouping reagent red blood cell (A1/A2/B/AB/O) IVD, kit, agglutination

NOTIFIED BODY: TÜV Product Service GmbH
Ridlerstrasse 65-80339 München (Germany)
CE-N° 0123

<i>Name:</i>	<i>Function:</i>	<i>Issued in:</i>	<i>Date:</i>	<i>Signature</i>
Galéa Diane	Site Quality Management Representative	Cressier FR	31.01.2017	

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Bio-Rad MSDS Number: 2703M

1.2 Relevant identified uses of the substance or mixture and uses advised against

No further relevant information available.

Application of the substance / the mixture

In vitro diagnostic medical device or component.

1.3 Details of the supplier of the safety data sheet

Manufacturer/Supplier:

DiaMed GmbH
Pra Rond 23
CH-1785 Cressier FR
(Switzerland/Schweiz/Suisse/Svizzera)
Tel: +41 (0)26 674 51 11
Fax: +41 (0)26 674 51 45

Further information obtainable from: fds-msds.ch@bio-rad.com

1.4 Emergency telephone number: CHEMTREC Netherlands: +(31)-858880596

SECTION 2: Hazards identification

2.1 Classification of the substance or mixture

Classification according to Regulation (EC) No 1272/2008

The product is not classified according to the CLP regulation.

2.2 Label elements

Labelling according to Regulation (EC) No 1272/2008 Void

- Hazard pictograms Void
- Signal word Void
- Hazard statements Void

2.3 Other hazards

Results of PBT and vPvB assessment

- PBT: Not applicable.
- vPvB: Not applicable.

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3.2 Chemical characterisation: Mixtures

· **Description:** Mixture of substances listed below with nonhazardous additions.

· **Dangerous components:** Void

· **Additional information:** For the wording of the listed hazard phrases refer to section 16.

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Trade name: DiaMed GmbH IVD MD

(Contd. of page 2)

· 8.1 Control parameters
· Ingredients with limit values that require monitoring at the workplace:

The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.

· **Additional information:** The lists valid during the making were used as basis.

· 8.2 Exposure controls
· Personal protective equipment:
· General protective and hygienic measures:

The usual precautionary measures are to be adhered to when handling chemicals.

· **Respiratory protection:** Not required.

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· Material of gloves

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer.

As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

· Penetration time of glove material

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

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· 9.1 Information on basic physical and chemical properties
· General Information
· Appearance:

· **Form:**

Not determined.

· **Colour:**

Not determined.

· **Odour:**

Odourless

· **Odour threshold:**

Not determined.

· **pH-value:**

Not determined.

· Change in condition

· **Melting point/freezing point:**

Undetermined.

· **Initial boiling point and boiling range:** Undetermined.

· **Flash point:**

Not applicable.

· **Flammability (solid, gas):**

Not applicable.

(Contd. on page 4)

NL

Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 14.02.2018

Version number 2

Revision: 14.02.2018

Trade name: DiaMed GmbH IVD MD

(Contd. of page 3)

· Ignition temperature:	
· Decomposition temperature:	Not determined.
· Auto-ignition temperature:	Product is not selfigniting.
· Explosive properties:	Product does not present an explosion hazard.
· Explosion limits:	
· Lower:	Not determined.
· Upper:	Not determined.
· Vapour pressure:	Not determined.
· Density:	Not determined.
· Relative density	Not determined.
· Vapour density	Not determined.
· Evaporation rate	Not determined.
· Solubility in / Miscibility with	
· water:	Fully miscible.
· Partition coefficient: n-octanol/water:	Not determined.
· Viscosity:	
· Dynamic:	Not determined.
· Kinematic:	Not determined.
· Solvent content:	
· Water:	1.0%
· Solids content:	0.0%
· 9.2 Other information	No further relevant information available.

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SECTION 11: Toxicological information

- **11.1 Information on toxicological effects**
 - **Acute toxicity** Based on available data, the classification criteria are not met.
 - **Primary irritant effect:**
 - **Skin corrosion/irritation**
Based on available data, the classification criteria are not met.
 - **Serious eye damage/irritation**
Based on available data, the classification criteria are not met.
 - **Respiratory or skin sensitisation**
Based on available data, the classification criteria are not met.
 - **Additional toxicological information:** No further relevant information available.

(Contd. on page 5)

NL

Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 14.02.2018

Version number 2

Revision: 14.02.2018

Trade name: **DiaMed GmbH IVD MD**

(Contd. of page 4)

- **CMR effects (carcinogenity, mutagenicity and toxicity for reproduction)**
 - **Germ cell mutagenicity**
Based on available data, the classification criteria are not met.
 - **Carcinogenicity** Based on available data, the classification criteria are not met.
 - **Reproductive toxicity**
Based on available data, the classification criteria are not met.
- **STOT-single exposure** Based on available data, the classification criteria are not met.
- **STOT-repeated exposure**
Based on available data, the classification criteria are not met.
- **Aspiration hazard** Based on available data, the classification criteria are not met.

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SECTION 13: Disposal considerations

- **13.1 Waste treatment methods**
 - **Recommendation** Smaller quantities can be disposed of with household waste.
 - **Uncleaned packaging:**
 - **Recommendation:** Disposal must be made according to official regulations.
 - **Recommended cleansing agents:**
Water, if necessary together with cleansing agents.

SECTION 14: Transport information

- | | |
|---|------|
| <ul style="list-style-type: none"> · 14.1 UN-Number <ul style="list-style-type: none"> · ADR, ADN, IMDG, IATA | Void |
| <ul style="list-style-type: none"> · 14.2 UN proper shipping name <ul style="list-style-type: none"> · ADR, ADN, IMDG, IATA | Void |
| <ul style="list-style-type: none"> · 14.3 Transport hazard class(es) <ul style="list-style-type: none"> · ADR, ADN, IMDG, IATA · Class | Void |
| <ul style="list-style-type: none"> · 14.4 Packing group <ul style="list-style-type: none"> · ADR, IMDG, IATA | Void |
| <ul style="list-style-type: none"> · 14.5 Environmental hazards: <ul style="list-style-type: none"> · Marine pollutant: | No |

(Contd. on page 6)

NL

Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 14.02.2018

Version number 2

Revision: 14.02.2018

Trade name: DiaMed GmbH IVD MD

(Contd. of page 5)

· 14.6 Special precautions for user	Not applicable.
· 14.7 Transport in bulk according to Annex II of Marpol and the IBC Code	Not applicable.
· Transport/Additional information:	Not dangerous according to the above specifications.
· UN "Model Regulation":	Void

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- **Abbreviations and acronyms:**

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)
 IMDG: International Maritime Code for Dangerous Goods
 IATA: International Air Transport Association
 GHS: Globally Harmonised System of Classification and Labelling of Chemicals
 EINECS: European Inventory of Existing Commercial Chemical Substances
 ELINCS: European List of Notified Chemical Substances
 CAS: Chemical Abstracts Service (division of the American Chemical Society)
 PBT: Persistent, Bioaccumulative and Toxic
 vPvB: very Persistent and very Bioaccumulative

- *** Data compared to the previous version altered.**

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Anti-D monoclonal, Cell-Line MS 201

Anti-D monoclonal, Cell-Line RUM-1

Reagent for specific detection of D-antigen

Bloodgroup-Testreagent for tube-, slide/plate- and microplate techniques

For in vitro diagnostic use only, store at + 2 - 8 °C when not in use.

Product description:

Anti-D (cell line MS 201 and cell line Rum-1) are monoclonal human IgM blood grouping reagents which detect the D antigen when tested according to the recommended techniques. It has been obtained from an IgM antibody producing clone.

The Rh Blood Group system:

The observations of Levine and Statson in 1939 and of Landsteiner and Weiner in 1940 provided the basis for current understanding of the clinical significance and laboratory detection of Anti-D.

Approximately 15% of Caucasians lack the RhD antigen and are easily stimulated by a RhD positive pregnancy of blood transfusion to produce Anti-D. This may cause haemolytic disease of the newborn or severe haemolytic transfusion reaction.

Weakened Expression of the RhD antigen / Limitations:

The collective term D_w is widely used to describe red cells which have a weaker expression of the D antigen than normal. The term D weak denotes individuals with a reduced number of entire D antigen sites per red cell. The term D partial denotes individuals with missing D epitopes. D category VI is the D partial category which lacks most D epitopes. Both cell lines detect most examples of D weak and partial D red cells by direct agglutination, but it will not detect D category VI. The reagent are recommended as particular suitable for grouping patients. It is recommended that the tube technique is used for the detection of D weak or variant cells. The slide and microplate techniques are not recommended for the detection of weak D or variant cells.

Principle of the reagent:

When used by the recommended techniques this reagent will cause agglutination of red cells carrying the specific antigen. (positive test) Lack of agglutination of the red cells demonstrates the absence of the specific antigen. (negative test)

These reagents have been optimised for use by the recommended techniques without further dilution or additions.

The product is supplied sterile filtered to 0,2 µm.

Materials:

46.2 Talpa 10 mL

Anti-D cell line MS 201 as well as Anti-D cell line RUM-1 are composed of monoclonal IgM antibodies in a buffer solution containing macromolecular chemical potentiators. Each vial (10ml) contains sufficient material for approximately 250 tests. These reagents both contain 0,1% (w/v) sodium azide and Bovine material obtained from a USDA approved source free of TSES.

Precautions:

- All blood products should be treated as potential infectious. The human donor of the cell line used to produce this reagent has been tested and found to be negative for Anti-HIV, Anti-HCV, Anti-HbsAg, EBV and Mouse Antibody Production (MAP) viruses. No known regime of testing can completely guarantee that any product derived from human blood is incapable of transmitting infectious agents. Care should be exercised in the use and disposal of each container and its contents.
- Sodium azide (0,1% w/w final concentration) is added as a preservative. Sodium azide can cause high explosive metal azide combinations with lead and cooper. When pouring rinse with a lot of water.
- This product should be clear. Cloudiness may indicate bacterial contamination. It should not be used if a precipitate, fibrin-gel or particles are present.
- This reagent is for professional in vitro diagnostic use only.
- False positive or negative results may occur through contamination of test materials or any deviation from the recommended techniques.

Advice to users:

It is recommended that a positive control and a negative control should be tested in parallel with each batch of tests. Tests must be considered invalid if controls do not show the expected reactions.

It is not required to use a reagent control in parallel with all tests using this reagent. Only in typing the red cells of patients known to have auto antibodies or protein abnormalities is the use of a reagent control recommended. This should be tested in parallel with the reagent.

The reagent has been characterised by the procedure recommended in this package insert, its suitability or use in other techniques must be determined by the user.

In the event of changes in the analytical performance of the device or damage to the packaging please contact the Quality Assurance Department of the producer.

Storage:

Store the opened/unopened product at 2-8°C until the expiry date detailed on the product label.

Failure to store the product at the correct temperature, for example, storage at higher temperatures or repeated freezing and thawing may result in accelerated loss of reagent activity.

CE-IMMUNDIAGNOSTIKA GmbH

Am Seerain 13, D-74927 Eschelbronn Telefon: +49 6226 42100 Fax: +49 6226 42012

Specimen Collection:

Blood should be collected by an approved technique. The specimen should be tested as soon as possible following the collection. If a delay in testing should occur, store at 2-8°C. Specimen displaying gross haemolysis or microbial contamination should not be tested with this reagent. Failure to store the specimen at the correct temperature, for example, storage at higher temperatures or repeated freezing and thawing may result in false positive or false negative results.

Samples may be drawn into EDTA, ACD, others or may be drawn without anticoagulant.

Material required but not provided:

Tube technique: Test tube, Centrifuge, Isotonic saline, Timer, 37°C Incubator

Plate technique: Plates, Timer, Isotonic saline

Slide technique: Slides, Timer, Isotonic saline

Microplate technique: Microplates, Microplates shaker, Centrifuge, Isotonic saline, Timer

Recommended techniques:**Tube test:**

1. Prepare a 2-3 suspension of red cells in isotonic saline.
2. Add 1 drop of anti-D and one drop of red cells to the appropriately labelled tube and mix.
3. Centrifuge for 1 minute at 1500 UpM, or about 20 seconds at 900 – 1000 xg resp. a time ,appropriate to produce the strongest reaction.
4. Gently agitate each tube to resuspend the cell buttons and examine macroscopically for agglutination.
5. Record results and reactivity strengths.
Incubate all negative or weakly positive tests at 37°C for 5 minutes and repeat step 3 and 4. This may enhance the reaction strength in typing red cells of rare phenotypes.

Slide-test / plate test

1. Slide-tests are performed with whole blood, plate-tests with washed Erythrocytes or whole blood.
2. Place one drop of the reagent (appr. 50 µl) on a clean glass-, plastic slide or plate.
3. Add one drop of whole blood (resp. 35-45 suspension of red cells) to the slides or 10 red-cell suspension in 0,9 saline solution resp. whole blood to the plates using a transfer pipette or applicator stick.
4. Mix blood and reagent. On glass slides, use a separate clean applicator stick to mix each reagent/cell mixture over an area approximately 20 mm diameter. On plastic slides follow the manufacturer's insert.
5. Read and record results. This is achieved on glass slides by slow rotation over a period up to 2 minutes and on plates after an incubation time of 5 – 10 minutes. Incubation time for whole blood testings is limited to 5 min. maximum.
Do not place the slides or plates on a heated illuminated surface.
6. Observe for macroscopic agglutination and record results. Care should be taken not to mistake peripheral drying or fibrin strands as agglutination.

Microplate technique:

1. Prepare a 3-5 suspension of test red cells in Isotonic saline.
2. Add one drop of anti-D reagent to one appropriate test wells of a microplate.
3. Add an equal volume of the test cell suspension to the appropriate test wells.
4. Mix the contents of each well using manual means or a microplate shaker.
5. Incubate the microplate at room-temperatures for 15-20 minutes.
6. Centrifuge the microplate at 100 rd for 40 seconds.
7. Resuspend the red cells using microplate shaker (as no.4)
8. Read tests macroscopically or with an automated reader. The use of an automated plate reader must be validated by the customer.

Performance characteristic:

Anti-D reagents both have been tested by each of the recommended techniques with donor, clinical and neonatal specimen. These were collected in either ESTA, CPDA, Citrate or as clotted sample. The sample population represented all major Rh D phenotypes. Based on the number of tests the SENSITIVIT and the SPECIFICIT for each techniques was calculated with >99,5 resp. 100 . Definition from the Common technical Specifications (CTS):

Diagnostic Sensitivity: The probability that the device gives positive results in the presence of the target marker.

Detection of weak D-variants (D weak):

1. Add and mix 1 drop of Anti-D incomplete or Anti-D monoclonal with 1 drop of 3 red cell suspension.
2. Incubate this mixture for 5 minutes at 37 C.
3. Washing 3 times with isotonic saline.
4. Add 2 drops of Anti-Humanglobuline
5. Centrifuge for 1 minute at 160 rcf or 20 seconds at 1.000rcf.
6. Gently agitate each tube to resuspend the cell buttons and examine for agglutination.
7. Record results and reactivity strength. In case of agglutination a weak D-antigen is present. All negative reactions must be checked with IgG sensitibilized red cells for testing the reactivity of Anti-Humanglobulin.

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#46.2

EC Certificate

mdc medical device certification GmbH

Notified Body 0483
herewith certifies that

**CE-Immundiagnostika GmbH
Karl-Landsteiner-Str. 6
69151 Neckargemünd
Germany**

for the scope

**Reagents for blood typing:
ABO system, Rhesus (C, c, D, E, e), Kell system, Duffy system,
Kidd system and determination of anti-irregular erythrocyte antibodies**

has introduced and applies a

Quality System

for the design, manufacture and final inspection.

The mdc audit has proven that this quality system
meets all requirements according to

**Annex IV – excluding Section 4 and 6
of the Council Directive 98/79/EC**

of the European Parliament and of the Council of
27 October 1998 on in vitro diagnostic medical devices.

The surveillance will be held as specified in Annex IV, Section 5.

Valid from	2021-05-10
Valid until	2024-05-26
Registration no.	D1415300016
Report no.	P20-01441-184463
Stuttgart	2021-02-10

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CE-IMMUNDIAGNOSTIKA GmbH Karl-Landsteiner-Straße 6 D-69151 Neckargemünd Phone: +49 6223 80094 00/ Fax: +49 6223 80094 99 www.ce-immundiagnostika.com			
Instruction for use Use by professionals only Tests/ml: 25: with drop size 40µl when using separate volumetric pipettes			
Revision:		29/07-2019	
Product-Name:	Product-Code:	Product-Name:	Product-Code:
Anti-A A-11H5	A-mono-11H5	Anti-B B-6F9	B-mono-6F9
Anti-A BIRMA-1	A-mono-BIRMA	Anti-B LB-2	B-mono-LB2
Anti-AB A-5E10-B-2D7	AB-mono-5E10	(Mouse IgM)	
Reagent for specific detection of the corresponding antigen. Bloodgroup testreagent for microplate-, tube-, slide- and plate tests. All described test methods are only valid for manual applications as recommended in this instruction. The user must determine their suitability for use in other techniques (automates, semi-automates, gel-cards, others) according to recognized techniques and hints of the machine- or cards manufacturer in individual responsibility. Only for in-vitro diagnostic laboratory use. Store at + 2 - 8 °C when not in use.			

Product description:	Anti-A, Anti-B and Anti-AB are monoclonal <i>Mouse IgM</i> reagents which detect the corresponding red cell antigens in a direct agglutination reaction. Lack of agglutination indicates the absence of the corresponding antigen. The respective monoclonal antibodies derived from cultured hybridoma cell lines (<i>murine</i>). The antibodies (IgM) are diluted in a buffered isotonic saline solution containing bovine albumin (without stabilizers), EDTA and ingredients to facilitate the resuspension of red cell buttons following centrifugation. Sodium azide at a final concentration of < 0,1% w/w is added as a preservative.										
Clone:	Anti-A: A-11H5 and BIRMA-1, Anti-B: B-6F9 and LB-2, Anti-AB: A-5E10-B-2D7 Reagents Anti-A and Anti-B are coloured (blue respectively yellow) in order to avoid confusion of the reagents and to allow better control of the preparations.										
Note/Caution:	Sodium azide can cause high explosive metal azide combinations with lead and cooper. When pouring rinse with a lot of water. All blood products should be treated as potentially infectious. No known regime of testing can completely guarantee that any product derived from human blood is incapable of transmitting infectious agents. Care should be exercised in the use and the disposal of the container and its contents.										
Shelf life	Test-reagents can be used until the end of the product shelf life indicated on the label if stored correctly between 2 – 8° C when not in use. After the first opening the reagents must be stored properly closed at 2 – 8° C. Bacterial contamination are to be avoided. Clearly clouded reagents should not be used any longer.										
Test methods:	Samples may be drawn aseptic into the common anticoagulants (EDTA, ACD,). Testing should be performed as soon as possible after blood withdrawal to minimise the chance that falsely positive or falsely negative reactions will be encountered due to contamination or improper storage of a specimen. Samples that cannot be tested immediately should be stored at +2 – 8 °C..										
Additional materials required:	Isotonic saline, transfer pipettes, glass slides, applicator slicks, slides or plates, test tubes and test tube racks, validated serological centrifuge, cell panel, timer. Microplate tests: microplates, microplate shaker (optional), validated serological centrifuge, when used with reading machine or on automates it is the users responsibility to validate any accessory device for the intended use, NaCl solution, timer, transfer pipette, if necessary Bovine Albumin.										
Microplate test:	<p>MTP from different suppliers show different characteristics which might have non specific reaction of the red blood cells as a consequence. It is recommended to pre-treat new MTP before its first use in order to minimize the fastening of the red cells. Recommended are MTP with U-profile out of plastic.</p> <ol style="list-style-type: none"> Add 1 drop (30-50µl) of Bovine Albumin 22% to each well. Through careful movements or on a shaker mix well so that all wells uniformly are coated. Incubate 10-15 min. at RT (18-25°C). Pour off Bovine Albumin and give the contents of the wells in suitable waste containers Rinse MTP at least 10 x with tap water. Rinse MTP subsequently 2 x with distilled water. MTP tip over and dap away in order to remove surplus water. Dry MTP before use at the air. <p>Alternative methods possible as far as validated by the users.</p> <ol style="list-style-type: none"> Prepare a 2-4 % suspension of test red cells in isotonic solution. (Recommendation 2% suspension) Add one drop of the respective reagent (30-50µl) to the appropriate test wells of a U well microplate. Add an equal volume of the cell suspension to the appropriate test wells. Mix the contents of each well using manual means or a microplate shaker. (30 sec.) No incubation time necessary except for titrations or to the strengthening of weak Phenotypes. Centrifuge the microplates at 1.500 UpM until 60 sec. or other appropriate time and UpM. Re-suspend the red cells using the microplate shaker. (as in no. 4.) Read tests macroscopically or with an automated plate reader. The use of an automated plate reader must be validated by the customer. The use of additional visual remedies as mirror or magnifier can ease the reading. 										
Tube test:	<p>In order to improve the test-results it is recommended to wash the cells before re-suspension at least one time in 0,9% isotonic saline.</p> <table border="1" style="width: 100%;"> <tr> <td style="width: 5%;">1.</td> <td>Prepare a 2 – 3 % suspension of the washed red cells in 0,9% isotonic saline, plasma or serum.</td> </tr> <tr> <td>2.</td> <td>Add 1 drop of the respective test serum Anti-A or Anti-B or Anti-AB and one drop of red cells to the appropriately labelled tubes and mix.</td> </tr> <tr> <td>3.</td> <td>Centrifuge for 1 Min. at 400g (1.500 UpM) or with alternative UpM at appropriate time.</td> </tr> <tr> <td>4.</td> <td>Control supernatant on absence of Hemolysis which can have its origin on bacteriological contamination. Gently agitate each tube to re-suspend the cell buttons and examine for agglutination just after centrifugation.</td> </tr> <tr> <td>5.</td> <td>Record results and reactivity strengths.</td> </tr> </table>	1.	Prepare a 2 – 3 % suspension of the washed red cells in 0,9% isotonic saline, plasma or serum.	2.	Add 1 drop of the respective test serum Anti-A or Anti-B or Anti-AB and one drop of red cells to the appropriately labelled tubes and mix.	3.	Centrifuge for 1 Min. at 400g (1.500 UpM) or with alternative UpM at appropriate time.	4.	Control supernatant on absence of Hemolysis which can have its origin on bacteriological contamination. Gently agitate each tube to re-suspend the cell buttons and examine for agglutination just after centrifugation.	5.	Record results and reactivity strengths.
1.	Prepare a 2 – 3 % suspension of the washed red cells in 0,9% isotonic saline, plasma or serum.										
2.	Add 1 drop of the respective test serum Anti-A or Anti-B or Anti-AB and one drop of red cells to the appropriately labelled tubes and mix.										
3.	Centrifuge for 1 Min. at 400g (1.500 UpM) or with alternative UpM at appropriate time.										
4.	Control supernatant on absence of Hemolysis which can have its origin on bacteriological contamination. Gently agitate each tube to re-suspend the cell buttons and examine for agglutination just after centrifugation.										
5.	Record results and reactivity strengths.										

Instruction for use

Use by professionals only
 Tests/ml: 25: with drop size 40µl when using separate volumetric pipettes

Revision:

29/07-2019

Product-Name:	Product-Code:	Product-Name:	Product-Code:
Anti-A A-11H5	A-mono-11H5	Anti-B B-6F9	B-mono-6F9
Anti-A BIRMA-1	A-mono-BIRMA	Anti-B LB-2	B-mono-LB2
Anti-AB A-5E10-B-2D7	AB-mono-5E10	(Mouse IgM)	

Reagent for specific detection of the corresponding antigen. Bloodgroup testreagent for microplate-, tube-, slide- and plate tests.
 All described test methods are only valid for manual applications as recommended in this instruction. The user must determine their suitability for use in other techniques (automates, semi-automates, gel-cards, others) according to recognized techniques and hints of the machine- or cards manufacturer in individual responsibility.
 Only for in-vitro diagnostic laboratory use. Store at + 2 - 8 °C when not in use.

Slide-test/plate-test:	<ol style="list-style-type: none"> Slide-tests are performed with whole blood, plate-tests with washed Erythrocytes or whole blood. Place one drop of the reagent on a clean glass-, plastic slide or plate. Add one drop of whole blood (respectively 35-45% suspension of red cells) to the slides or one drop of whole blood in 0,9% saline solution (respectively 10% red-cell suspension) to the plate using a transfer pipette or applicator stick. Mix blood and reagent. This is achieved by slow rotation over a period up to 2 minutes (slides) and on plates after an incubation time of 5 – 10 minutes. Incubation time for whole blood testing is limited to 5 min. maximum. Observe macroscopically for agglutination and record results. Care should be taken not to mistake peripheral drying as agglutination. Do not place slides on or before a heated illuminated surface.
Interpretation of test results:	There must be agreement between the results of the antigen determination (cell grouping) and the determination of the alloantibodies (serum grouping) of a blood specimen. The interpretation of reactions obtained when testing infant bloods may be complicated by the fact that the infant's serum does not necessarily contain antibody for any antigen absent from the cells, and passive anti-a and/or anti-B antibodies from the mother's circulation may yield conflicting results when tests are performed on cord blood specimens. Cord blood specimens may also give weaker-than-normal reactions in the cell grouping, as the ABH antigens are imperfectly developed at birth.
Advice to users:	<ul style="list-style-type: none"> - Control: On each test, positive and negative control red cells have to be tested in parallel, - On plates, blood samples may occasionally react by rouleaux formation, which can be mistaken for a weak agglutination and may incorrectly be read as a positive result. This phenomenon has non-immunological causes. Rouleaux formation also occurs in heparin blood, blood from patients treated with plasma expanders (e.g. dextran) as well as in blood from patients with plasmacytoma (high protein content, changes in protein composition), oncological disease (abnormal blood count) or coagulation dysfunctions. For testing these patients, use the tube test, which usually avoids this phenomenon. - The reagent agglutinates weakly expressed antigens with normal or weaker agglutination strength (A3, Bweak) or weaker to negative reaction (Ax). To determine weakly expressed antigens the tube test should be used because of its higher sensitivity, if necessary with incubation for 30 minutes. - Unusually weak reactions until non reactivity may occur from A- and B subgroups. The red cells of people with some disease states may give falsely positive or falsely negative reactions with Anti-A or Anti-B. - Cord cells contaminated with Wharton's jelly may give falsely positive reactions. - Do not use monoclonal reagents with mouse antibodies in indirect Anti-Human-Globulin tests with AHG-reagent. - Certain subgroups of A and B may produce reactions that are weaker or even no reactions than those obtained with A or B cells of most random donors. The presence of such antibodies cannot be predicted. When sufficiently strong they can cause the non-specific agglutination of reagent A1 and B cells in serum (reverse) grouping tests. They can also produce non-specific agglutination in cell (forward) tests with anti-A, and -B and -AB when unwashed, plasma- or cells respectively in serum suspended cells or plasma samples are determined. Consequently for the determination of blood groups always ABO and Isoagglutinin determinations are mandatory. - Discrepancies between forward and reverse results should be investigated thoroughly before an ABO group is assigned, regardless of the strength of the reactions obtained in any cell or serum test.
Performance data:	The reagent fulfils the common technical specifications' requirements according to Annex II, List A der Directive 98/79/EC for in vitro diagnostics. It has the same or a better performance characteristics as comparable reagents in use. It was tested on more than 1000 samples with sensitivity and specificity of 100%.
Limitations:	Strengths of the test results are depending from the age of the blood. Falsely positive or falsely negative test results can occur from insufficient cell-concentration, inadequate incubation time or temperature, improper centrifugation, improper storage of materials or non- consideration of the instructions for the different test methods. They can occur as well from bacterial or chemical contamination of the anti-serum, cells or the saline solution. When less sensitive test methods than the tube test are used negative results can occur with Anti-AB and weak Ax cells. Than it is recommended to repeat the test with the tube method. <i>Clone BIRMA-1 has not been tested by the company with regard to recognizing Ax-cells.</i> The use of the anti-sera in machines or on gel cards may require dilutions. The use of such manipulated sera asks for re-validation under the responsibility of the user. This is valid for all manipulations as for example the cold freezing of the sera on microplates.
References:	<ol style="list-style-type: none"> Kohler C. & Milstein C. (1975), Continuous cultures of fused cells secreting antibody of predefined specificity. Nature, 256, 495-497. Lee H.H., Rouger P., Germain C., Muller A. & Salmon C. (1983). The production and standardisation of monoclonal antibodies as AB blood group typing reagents. Symposium of International Association of Biological Standardisation on monoclonal antibodies. Human Blood Groups, by Geoff Daniels, 1st Ed., Blackwell Science, Oxford 1995. HMSO, Guidelines for Blood Transfusion Services., 2nd Ed., 1994.

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EC Certificate

#46.3.1; 46.4; 46.5

mdc medical device certification GmbH

Notified Body 0483
herewith grants

**CE-Immundiagnostika GmbH
Karl-Landsteiner-Str. 6
69151 Neckargemünd
Germany**

for the scope

**Reagents for blood typing: AB0 system
(see attachment)**

the

EC Design Examination Certificate

The examination of the design of the product by mdc has proven
that the design meets the requirements according to

Annex IV – Section 4 of the Council Directive 98/79/EC

of the European Parliament and of the Council of
27 October 1998 on in vitro diagnostic medical devices.

This certificate is only valid in connection with a valid mdc certificate
according to Annex IV – excluding section 4 and 6 for the above mentioned products.

Valid from	2021-07-01
Valid until	2024-05-26
Registration no.	D1415300017
Report no.	P20-01451-184566
Stuttgart	2021-06-18

Head of Certification Body

Attachment of the certificate

No. D1415300017

Date 2021-06-18

Page 1 of 1

Product category	Product	Class
Monoclonal antisera: AB0-System	Anti-A, Klon A-11H5: A-11H5-0010-01 / -05 / -10 / -20 / -50 Anti-A, Klon BIRMA-1: A-BIRMA-0010-01 / -05 / -10 / -20 / -50 Anti-B, Klon B-6F9: B-6F9-0010-01 / -05 / -10 / -20 / -50 Anti-B, Klon LB-2: B-LB2-0010-01 / -05 / -10 / -20 / -50 Anti-AB, Klon A-5E10-B-2D7: AB-5E10-0010-01 / -05 / -10 / -20 / -50	List A, Annex II

Head of Certification Body

#46.3.1; 46.3.2; 46.4; 46.5 buteliukų talpa 10 ml

P.O. Confirmation No. 9152

Thank you for your order received, please find below our order a business accessible at www.ce-immundiagnostika.com

Item	Description
	Unfortunately can we only provide one LOT for each of the ordered reagents.
1	Anti - A monoclonal / 10 ml REF:A-mono-11H5-10 LOT: OAM155-1 Titre: A1: 1:512, A2B: 1:128 Clone: A-11 H 5 Exp. date: 2023-04 IFU: Rev. 29/07-2019
2	Anti - B monoclonal / 10ml REF:B-mono-6F9-10 LOT: OBM161-3 Titre : 1:512 Clone: B-6F9 Exp. date: 2023-02 IFU: Rev. 29/07-2019
3	Anti - D / 10 ml REF: D-mono-RUM-10 LOT: ODM320-1 Titer : 1:128 Clone: RUM1 Exp. date: 2023-09 IFU: Rev. 29/07-2019
4	Anti - Kell / 5 ml REF: Kell-mono-MS56-05 LOT: OKM206-2 Titre : 1:32 Clone: MS-56 Exp. date: 2023-06 IFU: Rev. 29/07-2019
5	Control reagent monoclonal/10 ml REF: KontRea-mono-10 LOT: MK-033-1

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1) SUBSTANCE/PREPARATION AND COMPANY IDENTIFICATION

Test reagents for use in the identification of the mentioned blood grouping antigens by the recommended techniques in the product inserts.

Product name:	Anti-A A-11H5	Anti-A BIRMA-1	Anti-B B-6F9	Anti-B LB-2	Anti-AB A-5E10-B-2D7
Product code:	A-mono-11H5	A-mono-BIRMA	B-mono-6F9	B-mono-LB2	AB-mono-5E10

Product name:	Anti-D MS-201	Anti-D RUM-1	Anti-D blend TH-28/MS-26	Anti-D blend 175 2-415 1E4
Product code:	D-mono-MS201	D-mono-RUM1	D-mono-blend-TH	D-mono-blend-175

2) COMPOSITION, INFORMATION ON INGREDIENTS

Chemical characterisation of the product:

Test reagents based on monoclonal antibodies.

Dangerous components:	Ingredients	CAS Number	Concentration
	Sodium Azide	26628-22-8	< 0,1 % (w/w)

Sodium azide at a product concentration below 0,1% by weight is not subject to labeling requirements.

3) HAZARDS IDENTIFICATION

Not a hazardous product according to Regulation (EC) No. 1272 / 2008.

4) FIRST AID MEASURES

Skin contact:	Wipe off affected area and flush with plenty of soap and water.
Eye contact:	Immediately flush eyes with plenty of water the eyelids open.
Inhalation:	Not applicable
Ingestion:	Rinse mouth and drink plenty of water. In case of indisposition seek medical advice.

5) FIRE FIGHTING MEASURES

Suitable extinguishing media:	Select according to surrounding materials.
Special hazards:	None

6) ACCIDENTAL RELEASE MEASURES

Personal precaution:	Avoid skin and eye contact. For personal protective equipment see heading 8.
Cleaning methods:	Collect spilled product with absorbing material (paper towels, cellulose), remove for disposal. Wipe the area until clean and dry.

7) HANDLING AND STORAGE

Handling:	For In Vitro Diagnostic Use according to the instructions for use, Treatment as potentially infectious.
Storage:	Store in tightly closed containers at 2...8°C.

8) EXPOSURE CONTROLS, PERSONAL PROTECTION**General protective and hygiene measures**

The usual safety rules and precautionary measures for handling chemical and potentially infectious material must be adhered to. Protective work clothing should be worn. Do not eat or drink while working with the product. After working with the product, disinfect hands, wash hands carefully and remove protective work clothing.

Personal protective equipment:

Respiratory protection:	Not necessary.
Hand protection:	Disposable laboratory gloves.
Eye protection:	Not necessary.

Safety Data Sheet	
According to Regulation (EC) No. 1907 / 2006	
Manufacturer / Supplier:	
CE-IMMUNDIAGNOSTIKA GmbH	
Germany, 74927 Eschelbronn, Am Seerain 13	
Phone +49 6226-42100 Fax +49 6226-42012	
Revision Date:	27/06-2016

Product name:	Anti-A A-11H5	Anti-A BIRMA-1	Anti-B B-6F9	Anti-B LB-2	Anti-AB A-5E10-B-2D7
Product code:	A-mono-11H5	A-mono-BIRMA	B-mono-6F9	B-mono-LB2	AB-mono-5E10

Product name:	Anti-D MS-201	Anti-D RUM-1	Anti-D blend TH-28/MS-26	Anti-D blend 175 2-415 1E4
Product code:	D-mono-MS201	D-mono-RUM1	D-mono-blend-TH	D-mono-blend-175

9) PHYSICAL AND CHEMICAL PROPERTIES	
Physical state:	Liquid
Color:	Anti-A: blue, Anti-B: yellow; Anti-AB: yellowish, Anti-D and Anti-D blend: clear
Odeur:	Odeurless
pH-value:	6,6 – 7,7 (25°C)
Boiling point:	100°C (212° F) aqueous solution)
Flash point:	Not applicable
Flammable limits:	Not applicable
Vapour density:	Not available
Specific gravity:	Not available

10) STABILITY AND REACTIVITY	
Conditions to avoid:	Store below 2°C and/or above 8°C, fire and excessive heat, freezing and thawing may decompose the product.
Substances to avoid:	Lead and copper salt.
Hazardous decomposition products:	Sodium azide when heated to decomposition liberates nitrogen gas and sodium, which is explosive.

11) TOXICOLOGICAL INFORMATION	
ACUTE TOXICITY:	Toxicological dates for the products are not available. The product is classified according to the Regulation (EC) No. 1272 / 2008 as not hazardous product .
CHRONIC EFFECTS.	No information found.
LISTED CARCINOGENS.	None
Nevertheless the product should be handled with the usual caution whilst handling biological products.	

12) ECOLOGICAL INFORMATION	
Data regarding ecological effects of this product are not available. If handled and used appropriately, ecological problems are not expected to occur.	

13) DISPOSAL CONSIDERATIONS	
Product:	In case there are residues to be disposed of (the quantities involved may only be a few Millilitres), these should be placed in disinfectant solution before disposal.
Packaging:	Disposal is to be made in accordance with official local regulations. Uncontaminated and completely emptied packaging may be treated as domestic waste or recycled

14) TRANSPORTATION INFORMATION	
Cooled at +2.. 8°C. No other special measures need be taken for transporting the products.	

15) REGULATORY INFORMATION	
This MSDS meets the requirements of Regulation (EC) Nr. 1907/2006.	

16) OTHER INFORMATION	
The information supplied here is based on data considered accurate and on our current state of knowledge. No warranty is expressed or implied regarding the accuracy of this data. Liability is expressly disclaimed for loss or injury arising out of use of this information or the use of any materials designated. It does not establish any contractual relationship.	
H300	Fatal if swallowed.
H400	Very toxic to aquatic life.
H410	Very toxic to aquatic life with long lasting effects.
EUH032	Contact with acids liberates very toxic gas.
EUH210	Safety data sheet available on request.

CE-IMMUNDIAGNOSTIKA GmbH Am Seerain 13 74927 Eschelbronn, Germany Tel: 06226-42100 / Fax: 06226-42012 www.ce-immundiagnostika.com	
Product information	CE
<i>Control-serum for monoclonal sera</i>	
Negative control for blood typing tests with monoclonal sera. For in vitro diagnostic use. Store at + 2 – 8 °C when not in use.	
Product description:	Parallel to blood typing tests it is mandatory to perform control tests to detect unspecific agglutination. Control-serum with monoclonal testsera is prepared from human serum without blood group specific antibodies and a supplement that is used for the preparation of monoclonal testsera.
Note/Caution:	Sodium azide as conservative (0,1) can cause high explosive metal azide combinations with lead and cooper. When pouring rinse with a lot of water. All sera derived from human blood should be treated as potential infectious. The source of material used to produce this reagent has been tested and found to be negative. No known regime of testing can completely guarantee that any testsera is incapable of transmitting infectious agents. Care should be exercised in the use and the disposal of the container and its contents.
Direction for use:	Control-serum for monoclonal sera is applied as the corresponding tests serum.
Test methods:	Samples may be drawn into EDTA, ACD, others or may be drawn without anticoagulants. Testing should be performed as soon as possible to minimise the chance that falsely positive or falsely negative reactions will be encountered due to contamination or improper storage of a specimen. Samples that cannot be tested immediately should be stored at +2 – 8 °C..
Tube test:	<p>In order to improve the test results it is recommended to wash the blood cells at least one time in 0,9 saline solution.</p> <ol style="list-style-type: none"> 1. Prepare a 2 - 3 suspension of the washed red cells in isotonic saline. 2. Add 1 drop of the respective test serum and one drop of red cells to the appropriately labelled tube and mix. 3. Incubate 10 – 15 minutes at room temperature. Centrifuge for 1 minute at 400g (1.500 UpM) or alternative UpM with appropriate time. 4. Gently agitate each tube to re-suspend the cell buttons and examine for agglutination. 5. Record results and reactivity strengths. Do not forget a negative and positive control test.
Slide-test/plate-test:	<ol style="list-style-type: none"> 1. It is recommended to wash patients or donors cells. 2. Place one drop of the reagent (appr. 50µl) on a clean glass-, plastic slide or plate. 3. Add one drop of whole blood (resp. 35-45 suspension of red cells) or 10 red cell suspension in 0,9 saline solution from the sample using a transfer pipette or applicator stick. 4. Do not place the slides or plates on a heated illuminated surface. 5. Mix the blood and reagent. On glass slides, use a separate clean applicator stick to mix each reagent/cell mixture over an area approximately 20mm diameter. On plastic slides follow the manufacturer's insert. Read and record results. This is achieved on glass slides by slow rotation over a period up to 2 minutes and on plates after an incubation time of 5-10 minutes. 6. Observe für macroscopic agglutination and record results. Care should be taken not to mistake peripheral drying or fibrin strands as agglutination.
Results and interpretation:	In case of a positive result (the recommended method has to be followed carefully) the test with the specific antiserum cannot be read. Further investigations should clarify whether the positive reaction is due to rouleaux formation, cold agglutination or autoantibodies. Then antisera for the saline test should be applied.
Limitations:	<p>Tube tests should be examined directly after the centrifugation, slide tests within 2 minutes and plate tests after 5 - 15 minutes in order to avoid falsely positive reaction due to dry up occurrences.</p> <p>The use of the anti-sera in machines may require dilutions. The use of such manipulated sera asks for re-validation under the responsibility of the user. This is valid for all manipulations as for example the cold freezing of the sera for mikroplates.</p> <p>Falsely positive or falsely negative test results can occur from bacterial or chemical contamination of test materials, inadequate incubation time or temperature, improper centrifugation, improper storage of materials or non consideration of the instructions of the different test methods. Strengths of the results are also depending from the age of the blood-cells.</p> <p>Do not freeze the testsera and use sera only until the expiry date indicated on the label / package.</p> <p>Care should be exercised in the use and the disposal of the container and is contents. Bovine serum albumin and the corresponding raw material were derived from bovine resources exclusively free of BSE.</p>

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
Revision: 01/07-2021

Valid until: 09/05-2024

Product Name:	Control reagent
Product-Code:	33110
Clones:	monoclonal
EU-Classification: List A or List B (Annex II Directive 98/79/EC)	IVD other

As exclusively responsible manufacturer we herewith declare that the above-mentioned products meet the provisions of the Council Directive Annex I 98/79/EC on in vitro diagnostic Medical Devices. All supporting documentation is retained under the Premises of the Manufacturer.

Supervising authority :	Regierungspräsidium Karlsruhe
Notified Body: Identification Number: 0483	mdc medical device certification GmbH Kriegerstraße 6 70191 Stuttgart Germany
Conformity Assessment According to:	Annex III.4 Directive 98/79/EC
Certificate No.:	D1415300015 valid till 2024-05-09
Place, Date of Issue:	Neckargemuend, 10/05.2021
Signature of QMB: (Angela Grajek)	 <p>CE-IMMUNDIAGNOSTIKA GmbH Karl-Landsteiner-Str. 6 69151 Neckargemünd Tel. +49 (0)6223 - 80094 00 Fax +49 (0)6223 - 80094 99 www.ce-immundiagnostika.com</p>

DNr.:252EN 07.2021

#46.3.2

Safety Data Sheet According to Regulation (EC) No. 1907 / 2006 Manufacturer / Supplier: CE-IMMUNDIAGNOSTIKA GmbH Germany, 74927 Eschelbronn, Am Seerain 13 Phone +49 6226-42100 Fax +49 6226-42012	
Revision Date:	27/06-2016

1) SUBSTANCE/PREPARATION AND COMPANY IDENTIFICATION

Test reagents for use in the identification of the mentioned blood grouping antigens by the recommended techniques in the product inserts.

Product name:	Control reagent monoclonal
Product code:	KontRea-mono

2) COMPOSITION, INFORMATION ON INGREDIENTS

Chemical characterisation of the product:

Test reagents based on human sera containing antibodies of the relevant specificity.

Dangerous components:	Ingredients	CAS Number	Concentration
	Sodium Azide	26628-22-8	< 0,1 % (w/w)

Sodium azide at a product concentration below 0,1% by weight is not subject to labeling requirements.

3) HAZARDS IDENTIFICATION

Not a hazardous product according to Regulation (EC) No. 1272 / 2008.

4) FIRST AID MEASURES

Skin contact:	Wipe off affected area and flush with plenty of soap and water.
Eye contact:	Immediately flush eyes with plenty of water the eyelids open.
Inhalation:	Not applicable
Ingestion:	Rinse mouth and drink plenty of water. In case of indisposition seek medical advice.

5) FIRE FIGHTING MEASURES

Suitable extinguishing media:	Select according to surrounding materials.
Special hazards:	None

6) ACCIDENTAL RELEASE MEASURES

Personal precaution:	Avoid skin and eye contact. For personal protective equipment see heading 8.
Cleaning methods:	Collect spilled product with absorbing material (paper towels, cellulose), remove for disposal. Wipe the area until clean and dry.

7) HANDLING AND STORAGE

Handling:	For In Vitro Diagnostic Use according to the instructions for use, Treatment as potentially infectious.
Storage:	Store in tightly closed containers at 2...8°C.

8) EXPOSURE CONTROLS, PERSONAL PROTECTION

General protective and hygiene measures

The usual safety rules and precautionary measures for handling chemical and potentially infectious material must be adhered to. Protective work clothing should be worn. Do not eat or drink while working with the product. After working with the product, disinfect hands, wash hands carefully and remove protective work clothing.

Personal protective equipment:

Respiratory protection:	Not necessary.
Hand protection:	Disposable laboratory gloves.
Eye protection:	Not necessary.

Safety Data Sheet	
According to Regulation (EC) No. 1907 / 2006	
Manufacturer / Supplier:	
CE-IMMUNDIAGNOSTIKA GmbH	
Germany, 74927 Eschelbronn, Am Seerain 13	
Phone +49 6226-42100 Fax +49 6226-42012	
Revision Date:	27/06-2016

Product name:	Control reagent monoclonal
Product code:	KontRea-mono

9) PHYSICAL AND CHEMICAL PROPERTIES	
Physical state:	Liquid
Color:	Yellowish
Odeur:	Odeurless
pH-value:	7,1 +/- 0,3 (25°C)
Boiling point:	100°C (212° F) aqueous solution)
Flash point:	Not applicable
Flammable limits:	Not applicable
Vapour density:	Not available
Specific gravity:	Not available

10) STABILITY AND REACTIVITY	
Conditions to avoid:	Store below 2°C and/or above 8°C, fire and excessive heat, freezing and thawing may decompose the product.
Substances to avoid:	Lead and copper salt.
Hazardous decomposition products:	Sodium azide when heated to decomposition liberates nitrogen gas and sodium, which is explosive.

11) TOXICOLOGICAL INFORMATION	
ACUTE TOXICITY:	Toxicological dates for the products are not available. The product is classified according to Regulation (EC) No. 1272 / 2008 as not hazardous product .
CHRONIC EFFECTS.	No information found.
LISTED CARCINOGENS.	None
Nevertheless the product should be handled with the usual caution whilst handling biological products.	

12) ECOLOGICAL INFORMATION	
Data regarding ecological effects of this product are not available. If handled and used appropriately, ecological problems are not expected to occur.	

13) DISPOSAL CONSIDERATIONS	
Product:	In case there are residues to be disposed of (the quantities involved may only be a few Millilitres), these should be placed in disinfectant solution before disposal.
Packaging:	Disposal is to be made in accordance with official local regulations. Uncontaminated and completely emptied packaging may be treated as domestic waste or recycled

14) TRANSPORTATION INFORMATION	
Cooled at +2 .. 8°C. No other special measures need be taken for transporting the products.	

15) REGULATORY INFORMATION	
This MSDS meets the requirements of Regulation (EC) Nr. 1907/2006.	

16) OTHER INFORMATION	
The information supplied here is based on data considered accurate and on our current state of knowledge. No warranty is expressed or implied regarding the accuracy of this data. Liability is expressly disclaimed for loss or injury arising out of use of this information or the use of any materials designated. It does not establish any contractual relationship.	
H300	Fatal if swallowed.
H400	Very toxic to aquatic life.
H410	Very toxic to aquatic life with long lasting effects.
EUH032	Contact with acids liberates very toxic gas.
EUH210	Safety data sheet available on request.

Biotestcell-A₁, -A₂, -B, -0

Reagent red blood cells for reverse typing and control of AB0- and Rh-blood grouping reagents

REF	816 009	46.6 pirkimo pozicijai siūlomi reagentai yra pakuotėje 4x10 ml (A₁, A₂, B, O).
VOL	4 x 10 ml	
PRES	Preservative: 0.01 % Neomycin sulfate, 0.033 % Chloramphenicol, 5 ppm Amphotericin B	

IVD Diagnostic reagent for in vitro use only
To be used by trained laboratory personnel only. ◀

Test purpose

- For the determination of serum characteristics in blood typing on slide, tile or tube test. According to the national guidelines the determination of AB0 blood types also requires the determination of serum markers (Anti-A and/or Anti-B). A blood type determination is regarded safe if the serum characteristics correspond to the red cell characteristics.
- As control red blood cells for AB0 blood typing reagents
- As control red blood cells for anti-A₁- and anti-H-reagents
- As control red blood cells for Rhesus typing reagents
- For the detection of irregular anti-A₁ and cold antibodies

Test principle

The test principle is a hemagglutination test. The antigens of the reagent red cells react with the respective antibodies in the serum or plasma to be tested. The existence or lack of Anti-A and/or Anti-B antibodies must correspond with the lack or existence of A and/or B antigens on the reagent red cells.

Reagent

The reagent red cells are available ready for use suspended approx. 4 % in modified Alsevers solution and can be used immediately following careful resuspension.

Biotestcell-A₁, -A₂, -B, -0 have the following antigen combinations:

Biotestcell-A ₁	A ₁ Rh positive (D positive)	(CcD.Ee)
Biotestcell-A ₂	A ₂ Rh negative (D negative)	(ccddee)
Biotestcell-B	B Rh negative (D negative)	(ccddee)
Biotestcell-0	0 Rh positive (D positive)	(CcD.Ee)

Materials required but not supplied

- slides, tiles, Bioplate
- glass tubes
- pipettes (drop volume 40-50µl)
- laboratory centrifuge
- isotonic saline
- Antihumaglobulin (e.g. Anti-Human-Globulin **REF** 804 020 or Anti-Human-Globulin Color **REF** 804 120 / 804 115 / 804130)
- Control Set QC (**REF** 816034)

Sample material

The specimen should be tested as soon as possible after collection. If testing is delayed, EDTA or citrate anticoagulated and clotted whole blood specimens should be stored at 2...8°C. Use of samples older than 14 days should be avoided unless there is no other alternative since antibody reactivity has been shown to decrease in older samples. ◀

For the determination of serum characteristics the sera to be tested do not need to be inactivated, because the stabilizer inhibits hemolysis. For reverse typing and as control red blood cells Biotestcell-A₁, -A₂, -B, -0 is used without pretreatment. For demonstration of hemolysis Biotestcell-A₁, -A₂, -B, -0 must be washed once with isotonic saline prior to use. Blood samples should be taken following general blood sampling guidelines.

Test procedure

Resuspend reagent red cells prior to use and bring up to room temperature.

Tile / slide Test

- Place 2 drops of serum to be tested on tile or Bioplate in the respective position.
- Add 1 drop of Biotestcell-A₁, -A₂, -B, -0 to the corresponding positions and mix well.
- Incubate for 15-30 minutes at room temperature.
- Slightly rock the tile back and forth and observe for agglutination.

Tube Test

- Place 2 drops of sample material to be tested into properly marked tubes.
- Add to each tube 1 drop of Biotestcell-A₁, -A₂, -B, -0 and mix. ◀
- Centrifuge for 2 minutes at 150-200 x g or 20 seconds at 800-1000 x g.
- Gently dislodge the cell button and observe for agglutination.

Detection of immune anti-A and anti-B

Anti-A and anti-B antibodies of the IgG class may cause immune hemolytic anemia, in newborns mhn. For the detection of such antibodies, Biotestcell-A₁, -A₂, -B, -0 must be washed once with isotonic saline solution prior to use. It is then used in an indirect antiglobulin test. Hemolysis of A or B cells indicates the presence of immune anti-A or anti-B antibodies.

Interpretation of results**Reaction**

test sera + patient red cells			Reagent red cells + serum/plasma				blood-group
Anti-A	Anti-B	Anti-AB	A ₁	A ₂	B	O	
+	-	+	-	-	+	-	A
-	+	+	+	+	-	-	B
-	-	-	+	+	+	-	0
+	+	+	-	-	-	-	AB

+ = agglutination

- = no agglutination

Evaluation of the reaction strength is carried out according to the Technical Manual (1):

Reaction strength	Agglutination
4+	One single agglutinate
3+	Several large agglutinates
2+	Medium-size agglutinates, clear background
1+	Small agglutinates, turbid background
+/-	Barely visible agglutination, turbid background
-	No agglutination

Performance characteristics and limitations of the method

- In very rare cases weak reactions (reaction strength under 3+) or hemolysis may occur.
- Isoantibodies may be absent in newborns, infants, elderly patients and patients with agammaglobulinemia.
- Since serum characteristics may react at difference strengths, incubation for 15 to 30 minutes at room temperature may be performed.
- The reactivity of the product may decrease during the dating period.
- The rate of decrease in reactivity is partially dependent on individual donor characteristics that are neither controlled nor predicted by the manufacturer.

In case of questionable results with unknown causes, our Bio-Rad Service (Tel. +49-6103-3130-611) will assist you.

Shelf life

After opening the vial the product can be stored under proper storage conditions (2...8°C) until the expiry date. Since antigenicity may decrease, the reagent red cells should not be used after the expiry date. Do not use damaged vials.

Warning and precautions

- The reagent cells should not be used if the cells are darkly discolored, spontaneously agglutinated or show considerable hemolysis or turbidity.
- Manual techniques are to be performed according to the manufacturer's instructions. Each deviation from these instructions is the sole responsibility of the user. Biotestcell-A₁, -A₂, -B, -0 is suited for use in automated systems. The use in automated systems is to be validated by the user.
- Used test material must be discarded as hazardous material. Waste management information can be found in the safety data sheet.
- All materials of human origin used for this product were tested for HbsAg, Anti-HCV, and Anti-HIV-1/2 and proved to be non-reactive. Nevertheless, all products of human origin must be regarded as potential transmitters of hepatitis, HIV, or other infectious diseases. Appropriate safety precautions are recommended.
- This product contains natural rubber latex which may cause allergic reactions.
- Internal quality controls according to national guidelines are recommended at regular intervals.

Biotestcell-A₁, -A₂, -B, -0 is produced every 4 weeks.

References

(1) Technical Manual, 17th edition, Section 1, American Association of Blood Banks

Key: Underline = Addition or significant change ◀ = Deletion of text

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#46.6

EC Design-Examination Certificate
Directive 98/79/EC Annex IV, Section 4
In Vitro Diagnostic Medical Devices

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Product

Identification: Reagent red blood cells for determination of serum characteristics of the AB0 system; control reagents for AB0, Rh and Kell reagents

(see attachment for products included)
Replaces Certificate, Registration No.: IL 60091479 0001

The Notified Body hereby declares that an examination of the design dossier relating to the listed products has been performed according to Annex IV, section 4 of the directive 98/79/EC and that the design of the devices conforms to the requirements of the abovementioned directive.

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TÜV Rheinland LGA Products GmbH is a Notified Body according to Directive 98/79/EC concerning in vitro diagnostic medical devices with the identification number 0197.

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Dr. H. Lüdemann

**SAFETY DATA SHEET**

according to Regulation (EC) No. 1907/2006 (REACH) and Regulation (EU) No. 2015/830

Revision date: 31/5/2017

Version: 3

Language: en-GB,IE

Date of print: 13/12/2017

Erytypecell

Material number 186018

Page: 1 of 7

SECTION 1: Identification of the substance/mixture and of the company/undertaking**1.1 Product identifier**

Trade name: Erytypecell

This safety data sheet pertains to the following products:

816006 Erytypecell-A1, -A2, -B, -O, 4x10 mL

816009 Biotestcell-A1, -A2, -B, -O, 4x10 mL

816041 Erytypecell-A1, 10 mL

816042 Erytypecell-A2, 10 mL

816043 Erytypecell-B, 10 mL

816044 Erytypecell-O, 10 mL

1.2 Relevant identified uses of the substance or mixture and uses advised againstGeneral use: Use as laboratory reagent.
Only for industrial users.**1.3 Details of the supplier of the safety data sheet**

Company name: Bio-Rad Medical Diagnostics GmbH

Street/POB-No.: Industriestr. 1

Postal Code, city: 63303 Dreieich
GermanyWWW: www.medizinische-diagnostik-dreieich.deE-mail: contact.bmd@bio-rad.com

Telephone: +49 (0)6103-3130-0

Telefax: +49 (0)6103-3130-646

Dept. responsible for information:

Produktmanagement Transfusion

Telephone: 06103 3130-611

Telefax: 06103 3130-724

1.4 Emergency telephone number**GIZ-Nord, Göttingen, Germany,
Telephone: +49 551-19240****SECTION 2: Hazards identification****2.1 Classification of the substance or mixture****Classification according to EC regulation 1272/2008 (CLP)**

This mixture is classified as not hazardous.

2.2 Label elements**Labelling (CLP)**

Hazard statements: not applicable

Precautionary statements: not applicable

SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006 (REACH) and Regulation (EU) No. 2015/830

Revision date: 31/5/2017
Version: 3
Language: en-GB,IE
Date of print: 13/12/2017

Erytypecell

Material number 186018

Page: 2 of 7

2.3 Other hazards

Results of PBT and vPvB assessment:
No risks worthy of mention.
No data available

SECTION 3: Composition / information on ingredients

3.1 Substances: not applicable

3.2 Mixtures

Chemical characterisation: Aqueous solution of anorganic salts and organic compounds.
Contains human source proteins.

SECTION 4: First aid measures

4.1 Description of first aid measures

In case of inhalation: Provide fresh air. Seek medical treatment in case of troubles.
Following skin contact: Remove residues with water. Change contaminated clothing. In case of skin reactions, consult a physician.
After eye contact: Immediately flush eyes with plenty of flowing water for 10 to 15 minutes holding eyelids apart. Remove contact lenses, if present and easy to do. Continue rinsing. In case of eye irritation consult an ophthalmologist.
After swallowing: Rinse mouth and drink large quantities of water.
Induce vomiting when the affected person is not unconscious. Observe risk of aspiration if vomiting occurs.
Seek medical attention. Never give anything by mouth to an unconscious person.

4.2 Most important symptoms and effects, both acute and delayed

No data available

4.3 Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

SECTION 5: Firefighting measures

5.1 Extinguishing media

Suitable extinguishing media:
Product is non-combustible. Extinguishing materials should therefore be selected according to surroundings.

5.2 Special hazards arising from the substance or mixture

Fires in the immediate vicinity may cause the development of dangerous vapours. In the event of a fire, the following may be produced when the water evaporates: Carbon monoxide and carbon dioxide.

5.3 Advice for firefighters

Special protective equipment for firefighters:
Wear self-contained breathing apparatus. Wear appropriate protective equipment.

SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006 (REACH) and Regulation (EU) No. 2015/830

Revision date: 31/5/2017
Version: 3
Language: en-GB,IE
Date of print: 13/12/2017

Erytypecell

Material number 186018

Page: 3 of 7

Additional information: Hazchem-Code: -
Do not allow fire water to penetrate into surface or ground water.

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6.1 Personal precautions, protective equipment and emergency procedures

Avoid contact with the substance. Provide adequate ventilation. Wear personal protection equipment. Do not breathe vapour/aerosol.

6.2 Environmental precautions

Do not allow to penetrate into soil, waterbodies or drains.

6.3 Methods and material for containment and cleaning up

Soak up with absorbent materials such as sand, siliceous earth, acid- or universal binder. Store in special closed containers and dispose of according to ordinance. Final cleaning.

6.4 Reference to other sections

Refer additionally to section 8 and 13.

SECTION 7: Handling and storage

7.1 Precautions for safe handling

Advices on safe handling: Avoid contact with skin and eyes. Do not breathe vapours. Wear appropriate protective equipment.
Keep all containers, equipment and working place clean.

7.2 Conditions for safe storage, including any incompatibilities

Requirements for storerooms and containers:

Keep containers tightly closed and at a temperature between 2 °C and 8 °C.

Hints on joint storage: Do not store together with acids/alkalies and oxidation agents.

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8.1 Control parameters

Additional information: Contains no substances with occupational exposure limit values.

8.2 Exposure controls

Provide good ventilation and/or an exhaust system in the work area.

Personal protection equipment

Occupational exposure controls

Respiratory protection: If vapours form, use respiratory protection. The filter class must be suitable for the maximum contaminant concentration (gas/vapour/aerosol/particulates) that may arise when handling the product.

SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006 (REACH) and Regulation (EU) No. 2015/830

Revision date: 31/5/2017
Version: 3
Language: en-GB,IE
Date of print: 13/12/2017

Erytypecell

Material number 186018

Page: 4 of 7

Hand protection: Protective gloves according to EN 374.
Glove material: Nitrile rubber - Breakthrough time: > 480 min.
Observe glove manufacturer's instructions concerning penetrability and breakthrough time.

Eye protection: Tightly sealed goggles according to EN 166.

Body protection: Wear suitable protective clothing.

General protection and hygiene measures:
Avoid contact with skin and eyes. Change contaminated clothing. After work, wash hands and face.

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9.1 Information on basic physical and chemical properties

Appearance: Form: liquid
Colour: reddish

Odour: odourless

Odour threshold: No data available

pH value: No data available

Melting point/freezing point: No data available

Initial boiling point and boiling range: approx. 100 °C

Flash point/flash point range: No data available

Evaporation rate: No data available

Flammability: not applicable

Explosion limits: LEL (Lower Explosion Limit): not applicable
UEL (Upper Explosive Limit): not applicable

Vapour pressure: No data available

Vapour density: No data available

Density: No data available

Solubility: No data available

Partition coefficient: n-octanol/water: No data available

Auto-ignition temperature: not applicable

Decomposition temperature: No data available

Viscosity, kinematic: No data available

Explosive properties: not applicable

Oxidizing characteristics: No data available

9.2 Other information

Additional information: No data available

SECTION 10: Stability and reactivity

10.1 Reactivity

refer to 10.3

10.2 Chemical stability

Stable under recommended storage conditions.

SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006 (REACH) and Regulation (EU) No. 2015/830

Revision date: 31/5/2017
Version: 3
Language: en-GB,IE
Date of print: 13/12/2017

Erytypecell

Material number 186018

Page: 5 of 7

10.3 Possibility of hazardous reactions

No hazardous reactions known.

10.4 Conditions to avoid

Protect against heat /sun rays.

10.5 Incompatible materials

Strong oxidizing agents, strong acids and alkalis.

10.6 Hazardous decomposition products

Thermal decomposition: No data available

SECTION 11: Toxicological information

11.1 Information on toxicological effects

Toxicological effects: Acute toxicity (oral): Lack of data.
Acute toxicity (dermal): Lack of data.
Acute toxicity (inhalative): Lack of data.
Skin corrosion/irritation: Lack of data.
Serious eye damage/irritation: Lack of data.
Sensitisation to the respiratory tract: Lack of data.
Skin sensitisation: Lack of data.
Germ cell mutagenicity/Genotoxicity: Lack of data.
Carcinogenicity: Lack of data.
Reproductive toxicity: Lack of data.
Effects on or via lactation: Lack of data.
Specific target organ toxicity (single exposure): Lack of data.
Specific target organ toxicity (repeated exposure): Lack of data.
Aspiration hazard: Lack of data.

SECTION 12: Ecological information

12.1 Toxicity

Further details: No data available

12.2 Persistence and degradability

Further details: No data available

12.3 Bioaccumulative potential

Partition coefficient: n-octanol/water:
No data available

12.4 Mobility in soil

No data available

SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006 (REACH) and Regulation (EU) No. 2015/830

Revision date: 31/5/2017
Version: 3
Language: en-GB,IE
Date of print: 13/12/2017

Erytypecell

Material number 186018

Page: 6 of 7

12.5 Results of PBT and vPvB assessment

No data available

12.6 Other adverse effects

General information: Do not allow to penetrate into soil, waterbodies or drains.

SECTION 13: Disposal considerations

13.1 Waste treatment methods

Product

Waste key number: 16 10 02 = aqueous liquid wastes other than those mentioned in 16 10 01
Recommendation: Small quantities can be disposed of with the domestic waste or burnt, with due observance of regulations of local authorities. Large quantities are hazardous waste and must be disposed of accordingly.

Contaminated packaging

Recommendation: Dispose of waste according to applicable legislation.
Handle contaminated packages in the same way as the substance itself.
Non-contaminated packages may be recycled.

SECTION 14: Transport information

14.1 UN number

ADR/RID, IMDG, IATA-DGR:
not applicable

14.2 UN proper shipping name

ADR/RID, IMDG, IATA-DGR:
Not restricted

14.3 Transport hazard class(es)

ADR/RID, IMDG, IATA-DGR:
not applicable

14.4 Packing group

ADR/RID, IMDG, IATA-DGR:
not applicable

14.5 Environmental hazards

Marine pollutant: no

14.6 Special precautions for user

No dangerous good in sense of these transport regulations.

14.7 Transport in bulk according to Annex II of Marpol and the IBC Code

No data available

SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006 (REACH) and Regulation (EU) No. 2015/830

Revision date: 31/5/2017
Version: 3
Language: en-GB,IE
Date of print: 13/12/2017

Erytypecell

Material number 186018

Page: 7 of 7

SECTION 15: Regulatory information

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

National regulations - Great Britain

Hazchem-Code: -
No data available

15.2 Chemical Safety Assessment

For this mixture a chemical safety assessment is not required.

SECTION 16: Other information

Further information

Reason of change: General revision
Date of first version: 8/9/2014

Department issuing data sheet

Contact person: see section 1: Dept. responsible for information

For abbreviations and acronyms, see: ECHA Guidance on information requirements and chemical safety assessment, chapter R.20 (Table of terms and abbreviations).

The information in this data sheet has been established to our best knowledge and was up-to-date at time of revision. It does not represent a guarantee for the properties of the product described in terms of the legal warranty regulations.