

Integruoti
moduliai

1.



B40 Patient Monitor

B40 Paciento monitorius

Affordable clinical excellence

Patients with acute, life-threatening conditions need the best possible care. The B40 Monitor from GE Healthcare provides a continuous flow of quality information to enhance clinical decision-making for adult, pediatric and neonatal¹ patients in various care areas. 1.

Advanced clinical parameters

The B40 Monitor is designed with advanced measurement technologies for accurate and reliable patient monitoring

- GE EK-Pro arrhythmia analysis
- GE DINAMAP* SuperSTAT non-invasive blood pressure
- GE TruSignal enhanced SpO₂ saturation monitoring options available: Nellcor** OxiMax** SpO₂ and Masimo** SET** SpO₂ algorithms
- GE Datex-Ohmeda* gas technology to support non-invasive monitoring in anesthesia and critical care areas
- GE Entropy* monitoring that provides information on the patient's central nervous system during general anesthesia
- Comprehensive package of neonatal¹ measurements

Ease of use for fast decision-making

The B40 Monitor makes it easy to acquire accurate patient data to support timely decision-making.

- 12.1" crystal-clear monitor displays up to six waveforms simultaneously
- Intuitive menus and one-button access to commonly-used functions
- 72-hour trend display with graphical and numerical data to review patient progress
- HL7** direct output and connectivity with the CARESCAPE* Gateway enables communications to EMR systems
- Capability to work in CARESCAPE Network and S/5 Network environments

Performance and reliability

With its streamlined design, the portable B40 Monitor fits into crowded spaces and is easily moved to different care areas as needed. The system's rugged design stands up to harsh environments and the everyday wear-and-tear of busy care areas. It will provide the performance and accuracy that you expect of GE equipment-so you can provide the care that your patients expect.

5. Galimyb jungti monitori Carescape arba S/5 informacin tinkl

Programin ranga HL7 skirta monitoriaus

8. duomen perdavimui ligonin s informacin sistem



¹ CO₂ measurement through E-miniC is intended for use with patients weighing over 5kg (11 lb) only. Entropy is intended for use with adult and pediatric patients older than 2 years.

Parameters and modules 9.1–9.8

Parameters	E-Modules*	Impedance Respiration
ECG		Range
Resp Kv pavimas		Adult/pediatric: 4 to 120 resp/min
SpO ₂		Neonate ¹ : 4 to 180 resp/min
NIBP Neinvazinis kraujosp dis	Configured Hemo module	Accuracy
Temp Temperat ra		±5% or ±5 resp/min, whichever is greater
2 channel InvBP Tiesioginis arterinis kraujospaudimas 2 kanalai		Gain range
Entropy 9.8	E-ENTROPY S mon s b kl s monitoravimas	0.1 to 5 cm/Ohm
Sidestream CO ₂	E-miniC	GE TruSignal SpO ₂
Sidestream CO ₂ , O ₂ , and N ₂ O	E-sCO	Pulse oximetry
Sidestream CO ₂ , O ₂ Agents and N ₂ O Kv puojam j ir anestetini duj monitoravimas	E-sCAiO	1 to 100% 1–100%
		Pulse rate
		30 to 250 bpm
		Measurement accuracy
		Saturation
		Matavimo tikslumas
		Without motion-adult/pediatric
		Finger Sensor: 70 to 100% ±2%
		Ear sensor: 70 to 100% ±3%
		Without motion-neonate ¹
		70 to 100% ±3%
		With motion-adult/pediatric/
		neonate ¹ 70 to 100% ±3%
		Low perfusion adult/pediatric 70 to
		100% ±3% (1~69% unspecified)
		Without motion: ±2 bpm
		(Adult/Pediatric/Neonate ¹)
		With motion: ±3 bpm
		(Adult/Pediatric/Neonate ¹)
		Low Perfusion: ±5 bpm
		(Adult/Pediatric)

NOTE: The monitor also is compatible with the E-sCOV and E-sCAiOV modules without Spirometry function.

NOTE: When monitoring neonatal¹ or other patients that have high respiration rate or low tidal volume, the CARESCAPE Respiratory Module E-sCO or E-sCAiO shall be used within the limits of respiration rates and tidal volumes to ensure specified measurement accuracy.

ECG EKG multi-derivacinis kanalas 9.1

Leads available	3-lead configuration: I, II, III	Pulse Rate
	5-lead configuration: I, II, III, aVR, aVL, aVF and V	Without motion: ±2 bpm
Sweep speed	12.5, 25 or 50mm/s	(Adult/Pediatric/Neonate ¹)
Gain range	0.2 to 5.0 cm/mV	With motion: ±3 bpm
Heart rate accuracy	30 to 300 bpm, ±5 % or ±5 bpm, whichever is greater	(Adult/Pediatric/Neonate ¹)
Pulsas		Low Perfusion: ±5 bpm
Bandwidth		(Adult/Pediatric)
50/60 Hz power supply	Monitor: 0.5 to 40 Hz	Nellcor OxiMax SpO ₂
	ST: 0.05 to 40 Hz	Measurement range
	Diagnostic: 0.05 to 150 Hz	Pulse oximetry
Pacemaker detection	Range 2 to 700 mV	1 to 100%
	Pulse width 0.5 to 2 ms	Pulse rate
Arrhythmia analysis	Asystole, Bradycardia, Tachycardia, Ventricular fibrillation, Ventricular tachycardia	20 to 250 bpm
ST segment analysis	ST numeric range -9 to +9 mm (-0.9 to +0.9 mV)	Measurement accuracy
	ST accuracy -8 mm to +8 mm ±0.2 mm or ±10%, whichever is greater	Saturation
	ST numeric resolution 0.1 mm (0.01 mV)	Adult 70 to 100% ±2%
	ST Trends Up to 72 hours	Neo 70 to 100% ±3%
		Low perfusion 70 to 100% ±2%
		Pulse Rate
		±3 bpm

¹ CO₂ measurement through E-miniC is intended for use with patients weighing over 5kg (11 lb) only.

Entropy is intended for use with adult and pediatric patients older than 2 years.

Masimo SET SpO₂

Measurement range

Pulse oximetry	1 to 100%
Pulse rate	25 to 240 bpm

Measurement accuracy

Saturation	Without motion-adult/pediatric 70 to 100% ±2% Without motion-neonate ¹ 70 to 100% ±3% With motion-adult/pediatric/ neonate ¹ 70 to 100% ±3% Low perfusion 70 to 100% ±2% (0~69% unspecified)
Pulse Rate	Without motion ±3 bpm With motion ±5 bpm

NIBP

Measurement technique Oscillometric with step deflation

Modes Manual, automatic and stat

NIBP Measurement ranges

Systolic Adult/Pediatric: 30 to 290 mmHg

**Matavimo ribos-
MAP diapazonas** Neonate¹: 30 to 140 mmHg
Adult/Pediatric: 20 to 260 mmHg

Neonate¹: 20 to 125 mmHg

Diastolic Adult/Pediatric: 10 to 220 mmHg

Neonate¹: 10 to 110 mmHg

Accuracy Meets AAMI SP10

**Default initial
inflation pressure** Adult/Pediatric: 135 ±15 mmHg
Neonate¹: 100 ±15 mmHg

**Maximum
determination time** Adult/Pediatric 2 minutes
Neonate¹: 85 s

Over pressure monitor Adult/Pediatric: 300 ±6 to
330 mmHg
Neonate¹: 150 ±3 to 165 mmHg

Invasive blood pressure 13.1–13.2

Measurement range -40 to 320 mmHg (-5.3 to 42.7 kPa)

Measurement accuracy ±5% or ±2 mmHg,
whichever is greater

Frequency response 4 to 22 Hz

Transducer sensitivity 5 μV/V/mmHg

13.1 invazinio arterinio kraujo spaudimo matavimo diapazonas

13.2 invazinio arterinio kraujo spaudimo matavimo paklaida

¹ CO₂ measurement through E-miniC is intended for use with patients weighing over 5kg (11 lb) only.

Entropy is intended for use with adult and pediatric patients older than 2 years.

Temperature

Numerical display	T1, T2, T2-T1
Scale	° Fahrenheit (F) ° Celsius (C)

Measurement range 10 to 45°C (50 to 113°F)

Measurement accuracy ±0.1°C without probe

Display Resolution ± 0.1°C at 25 to 45°C with
reusable probes

Probe YSI probes recommended by
GE Healthcare

Carbon dioxide (CO₂) E-sCO or E-sCAiO

GE infrared absorption sensor technology

CO₂ waveform

EtCO₂ End-tidal CO₂ concentration

FiCO₂ Inspired CO₂ concentration

Measurement range 0 to 15 vol%
(0 to 15 kPa, 0 to 113 mmHg)

Accuracy ±(0.2 vol% +2% of reading)

Rise time <260 ms

Adjustable low and high alarm limits for EtCO₂ or FiCO₂

Carbon dioxide (CO₂) – E-miniC

Measurement range 0 to 20 vol%

Accuracy 0 to 15 vol%
±(0.2 vol% +2% of reading)

15 to 20 vol%
±(0.7 vol% +2% of reading)

Resolution 0.1%

Rise time <300 ms with nominal flow

Warm-up time 1 minute for operation
30 minutes for full specification

NOTE: CO₂ measurement is intended for use with patients weighing over 5kg (11 lb) only.

Respiration rate (RR) – E-sCO or E-sCAiO

Measurement range 4 to 100 breaths/min

Detection criteria 1 vol% change in CO₂ level

Alarm note sent to host device if no breath detected in
20 seconds

Respiration rate (RR) – E-miniC

Breath detection	1% variation in CO ₂ level
Measurement range	4 to 80 bpm
Accuracy	±1 bpm in the range 4 to 20 bpm
Resolution	1 bpm

NOTE: CO₂ measurement is intended for use with patients weighing over 5kg (11 lb) only.

Patient oxygen (O₂)

GE differential paramagnetic sensor

O₂ waveform

FiO ₂	Inspired O ₂ concentration
EtO ₂	End-tidal O ₂ concentration
FiO ₂ -EtO ₂	Inspired-expired difference
Measurement range	0 to 100 vol%
Accuracy	±(1 vol% +2% of reading)
Rise time	<260 ms

Nitrous oxide (N₂O)

GE infrared absorption sensor

FiN ₂ O	Inspired N ₂ O concentration
EtN ₂ O	End-tidal N ₂ O concentration
Measurement range	0 to 100 vol%
Accuracy	±(2 vol% +2% of reading) N ₂ O < 85%

Anesthetic agent (AA)

GE infrared absorption sensor

Anesthetic agent waveform

FiAA	Inspired anesthetic agent concentration
EtAA	End-tidal anesthetic agent concentration

MAC or MACage value options

Agent mixture detection

Measurement range

Sevoflurane	0 to 8 vol%
Desflurane	0 to 20 vol%
Isoflurane, enflurane, halothane	0 to 6 vol%
Accuracy	±(0.15 vol% +5% of reading)

Agent identification

Identification threshold	0.20 (+0.15/-0.05) vol%
Detection time	<20 sec

Entropy

Display range	RE 0 to 100 SE 0 to 91 BSR 0 to 100%
Display resolution	RE/SE 1 BSR 1%

Thermal Recorder

Operating Principle	Thermal array
Waveforms	Selectable 1, 2, or 3 waveforms
Numerics	HR, SpO ₂ , NIBP, IBP1, IBP2, ETCO ₂ , T1, T2, Resp, O ₂ , AA
Tabular trend printout	HR, NIBP, IBP1, IBP2, T1, T2, Et/FiCO ₂ , RR, Et/Fi O ₂ , Et/Fi AA
Graphical trend printout	HR, ST, IBP1, IBP2, NIBP, SpO ₂ , Pleth, CO ₂ , N ₂ O, O ₂ , AA, Resp, T1+T2, Entropy

Resolution

Vertical	8 dots/mm (200 dots/inch)
Horizontal	24 dots/mm (600 dots/inch)
Paper width	50 mm, printing width 48 mm
Paper speed	1, 6.25, 12.5, 25 mm/s

I/O connectors

RS-232 computer serial output, Defibrillation synch, Nurse call

Networking

Work in CARESCAPE Network and S/5 Network environments

Printing

Network laser printer supported in S/5 network

Mounting

GCX compatible

Integrated transport handle

Electrical specifications

AC input	100 to 240V ±10%, 50/60 Hz, 150VA
Protection	Class I
Battery	Exchangeable lithium-ion, 2 pcs max
Charging time	2 hours per battery pack
Run time	Up to 4.5 hours

Physical specifications

Dimensions (HxWxD)

Without extension rack 312 mm x 312 mm x 158 mm
(±5 mm)

With extension rack 312 mm x 352 mm x 178 mm
(±5 mm)

Weight ≤7Kg

Ingress protection IP21

Warranty

One year

Certifications

IEC 60601-1 passed

CE marking according to Directive 93/42/ EEC

Technical specifications

Display Ekranas 2

Display size 2 **12.1 inch** Ekranas dydis

Resolution 2 **SVGA resolution, 800 x 600** Ekranas rezoliucija
Snapshot

Number of traces 3 **Up to 6** 6 kreivės

Display layout and colors User-configurable

Controls Monitoriaus Valdymas 4

TrimKnob* control and hard keys

Standard

Alarms

Priority High, Medium, Low and Message

Notification Audible and visual

Setting Default and individual

Visual alarm notification Red, yellow, cyan

Audio silence message

General alarm message

Audio pause 2 minutes

Adjustment Central alarm display and adjustment page

Trending 10 minutes graphical trends referenced to set alarm limits

Trends

Graphical **All parameters**, selectable time scales from 20 minutes to 72 hours

Numerical **All parameter**s, every 5 minutes sampling or after NIBP determination

Up to 10 snapshots

Manual or alarm triggered

OCRG trend Real time or snapshot

Neonate mode only

Trend cursor In both graphical and numerical trends

Minitrends 5 or 30 minute minitrends can be displayed for a continuous historical view

About GE Healthcare

GE Healthcare provides transformational medical technologies and services that are shaping a new age of patient care. Our broad expertise in medical imaging and information technologies, medical diagnostics, patient monitoring systems, drug discovery, biopharmaceutical manufacturing technologies, performance improvement and performance solutions services help our customers to deliver better care to more people around the world at a lower cost. In addition, we partner with healthcare leaders, striving to leverage the global policy change necessary to implement a successful shift to sustainable healthcare systems.

Our “healthymagination” vision for the future invites the world to join us on our journey as we continuously develop innovations focused on reducing costs, increasing access and improving quality around the world. Headquartered in the United Kingdom, GE Healthcare is a unit of General Electric Company (NYSE: GE). Worldwide, GE Healthcare employees are committed to serving healthcare professionals and their patients in more than 100 countries. For more information about GE Healthcare, visit our website at www.gehealthcare.com.

GE Healthcare
P.O. Box 900, FIN-00031 GE, Finland
GE Direct United Kingdom: +44 (0)800 0329201

www.gehealthcare.com



GE imagination at work

© 2013 General Electric Company – All rights reserved.

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE representative for the most current information.

Trademark of General Electric Company is the preferred form of trademark notice. This form of notice can be used with registered and unregistered marks, and marks that are the subject of a pending application.

GE and the GE Monogram are trademarks of the General Electric Company.

* CARESCAPE, Entropy, DINAMAP, Ohmeda and TrimKnob are trademarks of GE Medical Systems Information Technologies, Inc.

** Masimo and SET are trademarks of Masimo Corporation.

** Nellcor and OxiMax are trademarks of Nellcor Puritan Bennett, Inc.

** HL7 is a registered trademark of Health Level Seven, Inc

General Electric Company, by and through its GE Healthcare division.

EMEA DOC1290796 3/13
(Global version DOC1266491)

Neuromuscular Transmission Module, E-NMT

14.



For integrated measurement of the level of neuromuscular block

The Neuromuscular transmission module provides a unique, integrated measurement of the level of neuromuscular block.

Features

- Provides quantitative, automatic measurement of muscle response to an electrical stimulus
- All common stimulation modes: train of four (TOF), single twitch (ST), double burst stimulation (DBS) and tetanic stimulation
- Covers the whole range of neuromuscular blockade
- Two measurement sensor options: the easy-to-apply MechanoSensor (KMG) for routine clinical use, and ElectroSensor (EMG) for research use
- Two MechanoSensor options: for adults and for pediatric patients
- Automatic setting of supramaximal current
- Automatic measurement at user-selectable interval
- Recovery Note can be set to notify when the block is wearing off
- Recall function enables moving the module with the patient without losing the reference values and the measured supramaximal current
- In addition to the integrated neuromuscular block measurement, the module enables nerve location for regional block by utilizing a comfortable 40 μ s pulse



Technical Data

Direct function keys

Start-up	Automatically initiates the measurement by setting supramaximal current, reference value and starting cycle
Stop/Continue	Stops / continues measurement with same patient

Dimension and weight

Module size (WxDxH)	37 x 186 x 112 mm/1.5 x 7.3 x 4.4 in
Module weight	0.35 kg/0.8 lb

NMT

Stimulation modes	Train of four, TOF Double burst, DBS (3,3) Single twitch, ST 50 Hz tetanic & post tetanic count, PTC
14.1	
Numeric display	TOF%/DBS%, Count, T1%, PTC
Measurement intervals for TOF/DBS	Manual, 10 sec, 12 sec, 15 sec, 20 sec, 1 min, 5 min, 15 min in Anesthesia Monitor and Compact Anesthesia Monitor Manual, 20 sec, 1 min, 5 min, 15 min, 30 min, 60 min, 120 min in Critical Care Monitor and Compact Critical Care Monitor
Measurement intervals for ST	Manual, 1 sec, 10 sec, 20 sec
Stimulus pulse	Square wave, constant current
Pulse width 14.2	100, 200 or 300 μs
Stimulus current range	supramax 10 to 70 mA manual 10 to 70 mA with 5 mA steps

Stimulus current accuracy	10% or ± 3 mA, whichever is greater
Max load for 70mA	3 k Ω
Max voltage	300 V

Regional Block mode

Stimulation mode	Single twitch
Stimulation intervals	1, 2, 3 sec
Stimulus pulse	Square wave, constant current
Pulse width	40 μ s
Stimulus current range	0 to 5.0 mA with 0.1 mA steps
Stimulus current accuracy	20% or ± 0.3 mA, whichever is greater

Ordering information

Description	Order code
Neuromuscular Transmission Module	E-NMT
NMT Sensor Cable, 3.3 m/11 ft	888414
1.5 m/5 ft	888415
NMT MechanoSensor	888418
NMT MechanoSensor, pediatric	897439
NMT ElectroSensor	888416
NMT Regional Block Adapter	888417
NMT Electrodes 30 pcs	57268

GE Healthcare
P.O. Box 900, FIN-00031 GE, Finland
Tel. +358 10 394 11 • Fax +358 9 146 3310

www.gehealthcare.com

© 2005 General Electric Company – All rights reserved.
GE and GE Monogram are trademarks of
General Electric Company.

GE Healthcare Finland Oy, a General Electric company,
going to market as GE Healthcare.

www.datex-ohmeda.com



Neuroraumeninio laidumo modulis, E-NMT

Integruotam neuroraumeninės blokados lygio matavimui

Neuroraumeninio laidumo modulis suteikia išskirtinę integruoto neuroraumeninės blokados lygio matavimo galimybę.

Savybės

- Atlieka kiekybinį, automatinį raumens atsako į elektrinį dirgiklį matavimą
- Visi įprastiniai stimuliacijos būdai: 4 impulsų sekos (angl. *train of four*, TOF), pavienio impulso (angl. *single twitch*, ST), dvigubo pliūpsnio (angl. *double burst stimulation*, DBS) ir tetaninė stimuliacija
- Apima visą neuroraumeninės blokados gylio skalę
- Dviejų rūšių matavimo davikliai: patogus naudoti „MechanoSensor“ (KMG) kasdieniniam klinikiniam vartojimui ir „ElektroSensor“ (EMG) moksliniam naudojimui
- Dvi „MechanoSensor“ pasirinktys: suaugusiems ir vaikams
- Automatinis supramaksimalios srovės nustatymas
- Automatinis matavimas vartotojo pasirinktame intervale
- Gali būti nustatytas „atsigavimo“ signalas, kad praneštų apie besibaigiančią blokadą
- Atminties funkcija leidžia transportuoti modulį kartu su pacientu neprarandant pradinių reikšmių ir nustatytos supramaksimalios srovės
- Be integruoto neuroraumeninės blokados matavimo, modulis leidžia nustatyti nervo lokalizaciją regioninei blokadai, panaudojant 40 μ s dažnio impulsus

Techniniai duomenys

Tiesioginių funkcijų mygtukai

„Start-up“	Automatiškai pradeda matavimą nustatydamas supramaksimalią srovę, pradinę reikšmę ir paleidimo ciklą
„Stop/Continue“	Sustabdo / tęsia matavimą tam pačiam pacientui

Matmenys ir svoris

Modulio dydis (plotis, gylis, aukštis)	37 x 186 x 112 mm/1.5 x 7.3 x 4.4 colio
Modulio svoris	0.35 kg/0.8 svaro

NMT

Stimuliacijos būdai	4 impulsų seka (TOF) Dvigubas pliūpsnis (DBS) (3,3) Pavienis impulsas (ST) 50 Hz tetaninis ir potetaninis stimulus (PTC)
Skaitmeninė išraiška	TOF%/DBS%, Count, T1%, PTC
TOF/DBS matavimo intervalai	Rankinis, kas 10 s, 12 s, 15 s, 20 s, 1 min, 5 min, 15 min anestezijos aparato monitoriuje ir nešiojamame anestezijos monitoriuje Rankinis, kas 20 s, 1 min, 5 min, 15 min, 30 min, 60 min, 120 min intensyvios terapijos skyriaus monitoriuje ir nešiojamame intensyvios terapijos skyriaus monitoriuje
ST matavimo intervalai	Rankinis, 1 s, 10 s, 20 s
Impulsų ritmas	„Laiptuota“ banga, pastovi srovė
Impulsų pasikartojimo dažnis	100, 200 ar 300 μ s
Impulsų srovės diapazonas	supramaksimali, nuo 10 iki 70 mA rankinis: nuo 10 iki 70 mA, kas 5 mA
Impulsų srovės tikslumas	10% arba ± 3 mA, kuris didesnis
70mA maksimali varža	3 k Ω
Maksimali įtampa	300 V

Regioninės blokados režimas

Stimuliacijos būdas	Pavienis impulsas
Stimuliacijos intervalai	1, 2, 3 s
Impulsų ritmas	“Laiptuota” banga, pastovi srovė
Impulsų pasikartojimo dažnis	40 μ s
Impulsų srovės diapazonas	nuo 0 iki 5.0 mA, kas 0.1 mA
Impulsų srovės tikslumas	20% arba ± 0.3 mA, kuris didesnis

Informacija užsakymui

Aprašas	Užsakymo kodas
Neuroraumeninio laidumo modulis	E-NMT
NMT daviklio kabelis,	
3.3 m/11 pėdų	888414
1.5 m/5 pėdų	888415
NMT „MechanoSensor“	888418
NMT „MechanoSensor“, vaikiškas	897439
NMT „ElectroSensor“	888416
NMT regioninės blokados adapteris	888417
NMT elektrodai, 30 vnt.	57268