



CHROMABOND® PS-RP / PS-OH⁻ / PS-H⁺ / PS-Mix / PS-Ag⁺ / PS-Ba²⁺

phases for RP and ion chromatography

★ Key features

- Very low degree of swelling, thus very well suited for chromatography, reliable function over the whole pH range from 0–14

🔧 Technical characteristics

- Base material high purity polystyrene-divinylbenzene copolymers (PS/DVB), pore size 100 Å, particle size 100 µm
- Different modifications for different applications from the elimination of nonpolar compounds up to the removal of specific polar components

✓ Recommended application

- Removal of interfering compounds
- Improves chromatographic separation, if the interfering components overlap with the analyte in the chromatogram
- Improves lifetime of the chromatographic column, since interfering components can irreversibly block the column packing
- Enrichment of the analytes

Properties of the individual modifications

| | | |
|---------------------|--|---|
| PS-RP | hydrophobic PS/DVB copolymer | removal of organic interfering components from water |
| PS-OH ⁻ | strong PS/DVB anion exchanger, OH ⁻ form capacity 0.6 meq/g | removal or concentration of anions from water increasing the pH value in acidic samples |
| PS-H ⁺ | strong PS/DVB cation exchanger, H ⁺ form capacity 2.9 meq/g | removal or concentration of cations from water decreasing the pH value of basic samples |
| PS-Mix | mixture of PS-OH ⁻ and PS-H ⁺ | desalting of water |
| PS-Ag ⁺ | strong PS/DVB cation exchanger, Ag ⁺ form | removal of halide ions from water |
| PS-Ba ²⁺ | strong PS/DVB cation exchanger, Ba ²⁺ form | removal of sulfate ions from water |

Removal of halides from aqueous samples shown for the trace analysis of nitrate besides an excess of chloride or bromide

MN Appl. No. 301930 / 302750

Compounds investigated:

20 ppm nitrate besides 2500 ppm chloride or 500 ppm bromide

Column type:

CHROMAFIX® PS-Ag⁺ (M) 0.8 mL, min. 250 mg
REF 731865

Column conditioning: 1 mL dist. water

Sample application and Elution:

apply 4 x 1 mL sample fractions to the cartridge, discard 1st mL, collect 2nd, 3rd and 4th mL separately

Further analysis: HPLC with column 250 x 4 mm NUCLEOSIL® Anion II; eluent 2 mmol/L potassium hydrogen phthalate pH 6, 2 mL/min; detection: indirect UV, 280 nm (see applications 110440 and 110450 at www.mn-net.com/apps)

| Phases | Adsorbent weight → 3 mL / 200 mg | 3 mL / 500 mg | 6 mL / 500 mg | 6 mL / 900 mg | Pack of | | |
|--------------------------------------|-------------------------------------|----------------------------|------------------|----------------------------|---------|----------------------------|---------|
| CHROMABOND® PS polypropylene columns | | | | | | | |
| PS-OH ⁻ | | | 730378 | | 30 | | |
| PS-H ⁺ | 730690 | 730376 | 730377 | | 30 | | |
| PS-Mix | | 730394 | | 730310 | 30 | | |
| Phases | Size S | Minimum adsorbent weight → | Size M | Minimum adsorbent weight → | Size L | Minimum adsorbent weight → | Pack of |
| CHROMAFIX® PS cartridges | | | | | | | |
| PS-RP | 731877 | 60 mg | 731875 | 160 mg | | | 50 |
| PS-OH ⁻ | 731868 | 70 mg | 731860 | 180 mg | 731862 | 510 mg | 50 |
| PS-H ⁺ | 731867 | 90 mg | 731861 | 220 mg | 731863 | 620 mg | 50 |
| PS-Mix | 731909 | 70 mg | | | | | 50 |
| PS-Ag ⁺ | 731866 | 100 mg | 731865 | 250 mg | | | 50 |
| PS-Ba ²⁺ | 731871 | 100 mg | 731870 | 250 mg | | | 50 |