

| REF          |     | SYSTEM   |
|--------------|-----|--|
| 05950929 190 | 100 | Elecsys 2010<br>MODULAR ANALYTICS E170<br><b>cobas e 411</b><br><b>cobas e 601</b><br><b>cobas e 602</b> |

**English****Please note**

The measured HE4 value of a patient's sample can vary depending on the testing procedure used. The laboratory finding must therefore always contain a statement on the HE4 assay method used. HE4 values determined on patient samples by different testing procedures cannot be directly compared with one another and could be the cause of erroneous medical interpretations. If there is a change in the HE4 assay procedure used while monitoring therapy, then the HE4 values obtained upon changing over to the new procedure must be confirmed by parallel measurements with both methods.

**Intended use**

Immunoassay for the quantitative determination of HE4 in human serum and plasma. The assay is used as an aid in monitoring recurrence or progressive disease in patients with epithelial ovarian cancer. Serial testing for patient HE4 values should be used in conjunction with other clinical findings used for monitoring ovarian cancer.

It is further intended to be used in conjunction with the Elecsys CA 125 II assay as an aid in estimating the risk of epithelial ovarian cancer in premenopausal and postmenopausal women presenting with pelvic mass. The results must be interpreted in conjunction with other methods in accordance with standard clinical management guidelines.

The electrochemiluminescence immunoassay "ECLIA" is intended for use on Elecsys and **cobas e** immunoassay analyzers.

**Summary**

The human epididymal protein 4 (HE4, also known as WFDC2) belongs to the family of whey acidic four-disulfide core (WFDC) proteins with suspected trypsin inhibitor properties.<sup>1,2</sup> The corresponding gene codes for a 13 kD protein. In its mature glycosylated form the protein has a molecular weight of approximately 20-25 kD and consists of a single peptide chain containing two WFDC domains.<sup>3</sup>

HE4 was first determined in the epithelium of the distal epididymis.<sup>4</sup> It shows low expression in epithelia of respiratory and reproductive tissues including ovary, but high expression in ovarian cancer tissue.<sup>5</sup> High secreted levels can also be found in the serum of ovarian cancer patients.<sup>6</sup> HE4 is also expected to help in the risk assessment of epithelial ovarian cancer.

Ovarian cancer is the fourth most common cause of cancer-related death in women worldwide. It is the most lethal form of gynecological cancer, and potentially curable if diagnosed early<sup>7,8</sup> and treated by surgeons familiar with the management of ovarian cancer.<sup>9,10</sup> However, the symptoms of ovarian cancer are related to the presence of adnexal masses and are often vague and unspecific. Thus, 70-75 % of ovarian cancers are detected at a late stage. According to the International Agency for Research on Cancer, the 5 year survival rate of ovarian cancer patients is 46 %. However, if the disease is diagnosed early, the survival rate increases up to 94 %.

As a single tumor marker, HE4 had the highest sensitivity for detecting ovarian cancer, especially in stage I disease, the early non-symptomatic stage. Combined, CA 125 and HE4 yielded the highest sensitivity with 76.4 % at 95 % specificity. Additionally, HE4 is more sensitive in early-stage endometrial cancer compared to CA 125. Elevated serum HE4 with normal CA 125 would suggest the presence of either ovarian or other type of cancer, for example endometrial cancer.<sup>11,12</sup>

Combined with other markers such as CA 125, HE4 can help determining whether a pelvic mass is benign or malignant in pre- and post-menopausal women. The dual marker combination CA 125 and HE4 is a more accurate predictor of malignancy than either alone.<sup>11</sup> Huhtinen et al. reported a 78.6 % sensitivity at 95 % specificity in ovarian carcinoma vs. endometriotic cysts.<sup>13</sup> HE4 levels correlate with clinical response to therapy or recurrence status in women with diagnosis of ovarian carcinoma as determined by CT

imaging. HE4 could thus be an important early indicator for disease recurrence.<sup>14</sup>

**Test principle**

Sandwich principle. Total duration of assay: 18 minutes.

- 1st incubation: 10 µL of sample, a biotinylated monoclonal HE4-specific antibody, and a monoclonal HE4-specific antibody labeled with a ruthenium complex<sup>a)</sup> form a sandwich complex.
- 2nd incubation: After addition of streptavidin-coated microparticles, the complex becomes bound to the solid phase via interaction of biotin and streptavidin.
- The reaction mixture is aspirated into the measuring cell where the microparticles are magnetically captured onto the surface of the electrode. Unbound substances are then removed with ProCell/ProCell M. Application of a voltage to the electrode then induces chemiluminescent emission which is measured by a photomultiplier.
- Results are determined via a calibration curve which is instrument-specifically generated by 2-point calibration and a master curve provided via the reagent barcode.

a) Tris(2,2'-bipyridyl)ruthenium(II)-complex (Ru(bpy)<sub>3</sub><sup>2+</sup>)

**Reagents - working solutions**

The reagent rackpack is labeled as HE4.

- M Streptavidin-coated microparticles (transparent cap), 1 bottle, 6.5 mL:  
Streptavidin-coated microparticles 0.72 mg/mL; preservative.
- R1 Anti-HE4-Ab~biotin (gray cap), 1 bottle, 10 mL:  
Biotinylated monoclonal anti-HE4 antibody (mouse) 0.75 mg/L;  
phosphate buffer 100 mmol/L, pH 6.5; preservative.
- R2 Anti-HE4-Ab~Ru(bpy)<sub>3</sub><sup>2+</sup> (black cap), 1 bottle, 10 mL:  
Monoclonal anti-HE4 antibody (mouse) labeled with ruthenium complex 1.5 mg/L; phosphate buffer 100 mmol/L, pH 7.4;  
preservative.

**Precautions and warnings**

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents.

Disposal of all waste material should be in accordance with local guidelines. Safety data sheet available for professional user on request.

Avoid foam formation in all reagents and sample types (specimens, calibrators and controls).

**Reagent handling**

The reagents in the kit have been assembled into a ready-for-use unit that cannot be separated.

All information required for correct operation is read in from the respective reagent barcodes.

**Storage and stability**

Store at 2-8 °C.

Do not freeze.

Store the Elecsys reagent kit **upright** in order to ensure complete availability of the microparticles during automatic mixing prior to use.

| Stability:              |                                  |
|-------------------------|----------------------------------|
| unopened at 2-8 °C      | up to the stated expiration date |
| after opening at 2-8 °C | 12 weeks                         |
| on the analyzers        | 28 days                          |

**Specimen collection and preparation**

Only the specimens listed below were tested and found acceptable.

Serum collected using standard sampling tubes or tubes containing separating gel.

Li-heparin, K<sub>2</sub>- and K<sub>3</sub>-EDTA plasma as well as Li-heparin plasma tubes containing separating gel.

Criterion: Slope 0.9-1.1 + intercept  $\leq \pm 10$  pmol/L, coefficient of correlation  $\geq 0.95$ .

Stable for 48 hours at 2-8 °C, 5 hours at 15-25 °C, 12 weeks at -20 °C. The samples may be frozen twice.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

Do not use heat-inactivated samples.

Do not use samples and controls stabilized with azide.

Ensure the samples, calibrators and controls are at 20-25 °C prior to measurement.

Due to possible evaporation effects, samples, calibrators and controls on the analyzers should be analyzed/measured within 2 hours.

**Materials provided**

See "Reagents – working solutions" section for reagents.

**Materials required (but not provided)**

- [REF] 05950945190, HE4 CalSet, for 4 x 1 mL
- [REF] 05950953190, PreciControl HE4, for 2 x 1 mL each of PreciControl HE4 1 and 2
- [REF] 03609987190, Diluent MultiAssay, 2 x 16 mL sample diluent
- General laboratory equipment
- Elecsys 2010, MODULAR ANALYTICS E170 or **cobas e** analyzer

For epithelial ovarian cancer risk assessment with ROMA (Risk of Ovarian Malignancy Algorithm):

- [REF] 11776223322, CA 125 II, 100 tests
- [REF] 11776240322, CA 125 II CalSet, 4 x 1 mL
- [REF] 11776452122, PreciControl Tumor Marker, for 2 x 3 mL each of PreciControl Tumor Marker 1 and 2
- [REF] 11732277122, Diluent Universal, 2 x 16 mL sample diluent or [REF] 03183971122, Diluent Universal, 2 x 36 mL sample diluent

Accessories for Elecsys 2010 and **cobas e** 411 analyzers:

- [REF] 11662988122, ProCell, 6 x 380 mL system buffer
- [REF] 11662970122, CleanCell, 6 x 380 mL measuring cell cleaning solution
- [REF] 11930346122, Elecsys SysWash, 1 x 500 mL washwater additive
- [REF] 11933159001, Adapter for SysClean
- [REF] 11706802001, Elecsys 2010 AssayCup, 60 x 60 reaction vessels
- [REF] 11706799001, Elecsys 2010 AssayTip, 30 x 120 pipette tips

Accessories for MODULAR ANALYTICS E170, **cobas e** 601 and **cobas e** 602 analyzers:

- [REF] 04880340190, ProCell M, 2 x 2 L system buffer
- [REF] 04880293190, CleanCell M, 2 x 2 L measuring cell cleaning solution
- [REF] 03023141001, PC/CC-Cups, 12 cups to prewarm ProCell M and CleanCell M before use
- [REF] 03005712190, ProbeWash M, 12 x 70 mL cleaning solution for run finalization and rinsing during reagent change
- [REF] 03004899190, PreClean M, 5 x 600 mL detection cleaning solution
- [REF] 12102137001, AssayTip/AssayCup Combimagazine M, 48 magazines x 84 reaction vessels or pipette tips, waste bags

- [REF] 03023150001, WasteLiner, waste bags

- [REF] 03027651001, SysClean Adapter M

Accessories for all analyzers:

- [REF] 11298500316, Elecsys SysClean, 5 x 100 mL system cleaning solution

**Assay**

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

Resuspension of the microparticles takes place automatically prior to use. Read in the test-specific parameters via the reagent barcode. If in exceptional cases the barcode cannot be read, enter the 15-digit sequence of numbers.

MODULAR ANALYTICS E170, **cobas e** 601 and **cobas e** 602 analyzers: PreClean M solution is necessary.

Bring the cooled reagents to approximately 20 °C and place on the reagent disk (20 °C) of the analyzer. Avoid foam formation. The system automatically regulates the temperature of the reagents and the opening/closing of the bottles.

**Calibration**

Traceability: This method has been standardized against the HE4 EIA method from Fujirebio Diagnostics, Inc.

Every Elecsys reagent set has a barcoded label containing specific information for calibration of the particular reagent lot. The predefined master curve is adapted to the analyzer using the relevant CalSet.

*Calibration frequency:* Calibration must be performed once per reagent lot using fresh reagent (i.e. not more than 24 hours since the reagent kit was registered on the analyzer). Renewed calibration is recommended as follows:

- after 12 weeks when using the same reagent lot
- after 7 days (when using the same reagent kit on the analyzer)
- as required: e.g. quality control findings outside the defined limits

**Quality control**

For quality control, use PreciControl HE4.

In addition, other suitable control material can be used.

Controls for the various concentration ranges should be run individually at least once every 24 hours when the test is in use, once per reagent kit, and following each calibration.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

**Calculation**

The analyzer automatically calculates the analyte concentration of each sample in pmol/L.

**Limitations - interference**

The assay is unaffected by icterus (bilirubin < 1130 µmol/L or < 66 mg/dL), hemolysis (Hb < 0.621 mmol/L or < 1.0 g/dL), lipemia (Intralipid < 2000 mg/dL) and biotin (< 205 nmol/L or < 50 ng/mL).

Criterion: Recovery within  $\pm 10$  % of initial value.

Samples should not be taken from patients receiving therapy with high biotin doses (i.e. > 5 mg/day) until at least 8 hours following the last biotin administration.

No interference was observed from rheumatoid factors up to a concentration of 1500 IU/mL.

There is no high-dose hook effect at HE4 concentrations up to 40000 pmol/L.

In vitro tests were performed on 18 commonly used pharmaceuticals. No interference with the assay was found.

Special cancer drugs were tested with concentrations shown in the table below. No interference with the assay was found.

Criterion: Recovery within  $\pm 10$  % of initial value.

| Drug                    | Concentration (µg/mL) |
|-------------------------|-----------------------|
| Carboplatin             | 600                   |
| Cisplatin               | 180                   |
| Cyclophosphamide        | 500                   |
| Dexamethasone           | 20                    |
| Doxorubicin             | 120                   |
| Leucovorin              | 750                   |
| Melphalan               | 15                    |
| Methotrexat-disodium    | 150                   |
| Paclitaxel              | 265                   |
| Fluorouracil            | 900                   |
| Bevacizumab (Avastin)   | 750                   |
| Erlotinib (Tarceva)     | 150                   |
| Rituximab (MabThera)    | 750                   |
| Trastuzumab (Herceptin) | 600                   |

In rare cases, interference due to extremely high titers of antibodies to analyte-specific antibodies, streptavidin or ruthenium can occur. These effects are minimized by suitable test design.

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

### Limits and ranges

#### Measuring range

15.0-1500 pmol/L (defined by the Limit of Detection and the maximum of the master curve). Values below the Limit of Detection are reported as < 15.0 pmol/L. Values above the measuring range are reported as > 1500 pmol/L (or up to 30000 pmol/L for 20-fold diluted samples).

#### Lower limits of measurement

*Limit of Blank (LoB), Limit of Detection (LoD) and Limit of Quantitation (LoQ)*

Limit of Blank = 5.00 pmol/L

Limit of Detection = 15.0 pmol/L

Limit of Quantitation = 20.0 pmol/L with a total allowable error of 30 %

The Limit of Blank, Limit of Detection and Limit of Quantitation were determined in accordance with the CLSI (Clinical and Laboratory Standards Institute) EP17-A requirements.

The Limit of Blank is the 95<sup>th</sup> percentile value from  $n \geq 60$  measurements of analyte-free samples over several independent series. The Limit of Blank corresponds to the concentration below which analyte-free samples are found with a probability of 95 %.

The Limit of Detection is determined based on the Limit of Blank and the standard deviation of low concentration samples. The Limit of Detection corresponds to the lowest analyte concentration which can be detected (value above the Limit of Blank with a probability of 95 %).

The Limit of Quantitation is defined as the lowest amount of analyte in a sample that can be accurately quantitated with a total allowable error of  $\leq 30$  %.

A study was performed using 5 equine serum samples and diluted human serum samples each for LoB and LoD respectively. The samples were tested in 6 runs over 3 days on 2 analyzers resulting in  $n = 60$  values. LoB and LoD were calculated to be 0.358 pmol/L and 0.661 pmol/L respectively. For LoQ 3 human serum samples were diluted and measured in 6 runs over 3 days on 2 analyzers. At a total allowable error of  $\leq 30$  % the LoQ was 4.42 pmol/L.

#### Linearity

The Elecsys HE4 assay is linear across the measuring range from 15.0-1500 pmol/L. Samples were prepared according to CLSI EP6-A, by diluting 3 serum and 3 plasma samples each with Diluent MultiAssay in multiple steps ranging from > 1500 pmol/L down to LoB.

#### Dilution

Samples with HE4 concentrations above the measuring range can be diluted with Diluent MultiAssay. The recommended dilution is 1:20 (either

automatically by the MODULAR ANALYTICS E170, Elecsys 2010 or **cobas e** analyzers or manually). The concentration of the diluted sample must be > 75 pmol/L.

After manual dilution, multiply the result by the dilution factor.

After dilution by the analyzers, the MODULAR ANALYTICS E170, Elecsys 2010 and **cobas e** software automatically takes the dilution into account when calculating the sample concentration.

#### Expected values

A study in one clinical center in Germany with the Elecsys HE4 assay on sera from 358 apparently healthy women yielded the following results:

|             |     | HE4 (pmol/L) |                             |
|-------------|-----|--------------|-----------------------------|
| Age (years) | N   | Median       | 95 <sup>th</sup> percentile |
| < 40        | 127 | 42.0         | 60.5                        |
| 40-49       | 65  | 44.3         | 76.2                        |
| 50-59       | 60  | 47.9         | 74.3                        |
| 60-69       | 60  | 55.0         | 82.9                        |
| $\geq 70$   | 46  | 62.1         | 104                         |

The distribution in percentage (%) of HE4 assay values determined in two clinical centers in Spain and Germany with the Elecsys HE4 assay in 896 female specimens is summarized in the table below:

|                                  |     | Elecsys HE4 values (pmol/L) |                |                |                |              |
|----------------------------------|-----|-----------------------------|----------------|----------------|----------------|--------------|
|                                  |     | 0.0-70.0                    | 70.1-140       | 140.1-500      | 500.1-1500     | > 1500       |
| N (percentage distribution)      |     |                             |                |                |                |              |
| Apparently healthy               |     |                             |                |                |                |              |
| Premenopausal                    | 90  | 76<br>(84.4 %)              | 13<br>(14.4 %) | 1<br>(1.1 %)   | 0<br>(0.0 %)   | 0<br>(0.0 %) |
| Postmenopausal                   | 106 | 63<br>(59.4 %)              | 40<br>(37.7 %) | 3<br>(2.8 %)   | 0<br>(0.0 %)   | 0<br>(0.0 %) |
| Benign conditions                |     |                             |                |                |                |              |
| Premenopausal                    | 177 | 160<br>(90.4 %)             | 16<br>(9.0 %)  | 1<br>(0.6 %)   | 0<br>(0.0 %)   | 0<br>(0.0 %) |
| Postmenopausal                   | 102 | 62<br>(60.8 %)              | 31<br>(30.4 %) | 9<br>(8.8 %)   | 0<br>(0.0 %)   | 0<br>(0.0 %) |
| Pregnancy                        | 50  | 50<br>(100 %)               | 0<br>(0.0 %)   | 0<br>(0.0 %)   | 0<br>(0.0 %)   | 0<br>(0.0 %) |
| Non-gynecological disease        | 35  | 16<br>(45.7 %)              | 6<br>(17.1 %)  | 6<br>(17.1 %)  | 7<br>(20.0 %)  | 0<br>(0.0 %) |
| CHF <sup>b)</sup>                | 23  | 9<br>(39.1 %)               | 11<br>(47.8 %) | 3<br>(13.0 %)  | 0<br>(0.0 %)   | 0<br>(0.0 %) |
| Cancer                           |     |                             |                |                |                |              |
| Ovarian cancer<br>Premenopausal  | 39  | 12<br>(30.8 %)              | 7<br>(17.9 %)  | 13<br>(33.3 %) | 5<br>(12.8 %)  | 2<br>(5.1 %) |
| Ovarian cancer<br>Postmenopausal | 97  | 10<br>(10.3 %)              | 19<br>(19.6 %) | 34<br>(35.1 %) | 28<br>(28.9 %) | 6<br>(6.2 %) |
| Endometrial cancer               | 49  | 18<br>(36.7 %)              | 20<br>(40.8 %) | 9<br>(18.4 %)  | 1<br>(2.0 %)   | 1<br>(2.0 %) |
| Breast cancer                    | 47  | 22<br>(46.8 %)              | 19<br>(40.4 %) | 5<br>(10.6 %)  | 1<br>(2.1 %)   | 0<br>(0.0 %) |
| Gastrointestinal cancer          | 46  | 19<br>(41.3 %)              | 20<br>(43.5 %) | 6<br>(13.0 %)  | 1<br>(2.2 %)   | 0<br>(0.0 %) |
| Lung cancer                      | 23  | 5<br>(21.7 %)               | 7<br>(30.4 %)  | 10<br>(43.5 %) | 1<br>(4.3 %)   | 0<br>(0.0 %) |
| Bladder cancer                   | 12  | 3<br>(25.0 %)               | 4<br>(33.3 %)  | 4<br>(33.3 %)  | 1<br>(8.3 %)   | 0<br>(0.0 %) |

b) CHF = congestive heart failure

In this study 84 % of the apparently healthy premenopausal women had an Elecsys HE4 assay value at or below 70 pmol/L and 97 % of the apparently healthy postmenopausal women had an Elecsys HE4 assay value at or below 140 pmol/L.

In this study the 95<sup>th</sup> percentiles for the apparently healthy pre- and postmenopausal women (all ages) were 92.1 pmol/L and 121 pmol/L, respectively.

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

#### Monitoring of disease status in patients diagnosed with ovarian cancer

The effectiveness of the Elecsys HE4 assay as an aid in monitoring of disease status in ovarian cancer patients was determined by assessing changes in HE4 levels in serial serum samples from 100 patients compared to changes in disease status. A study involving a total of 375 pairs of observations was undertaken with  $\geq 3$  blood withdrawals per patient. A positive change HE4 was defined as an increase in the value that was at least 20 % greater than the previous value of the test. 58.0 % (29/50) of the patient samples with a positive change correlated with the disease progression while 84.0 % (273/325) of the patient serial samples with no significant change in HE4 value correlated with no progression. The total concordance was 80.5 % (302/375). The following table presents the data in a 2 x 2 format.

| Change in disease state per sequential pair |             |                |       |
|---|-------------|----------------|-------|
| Increase in HE4 concentration               | Progression | No progression | Total |
| > 20 %                                      | 29          | 52             | 81    |
| $\leq 20$ %                                 | 21          | 273            | 294   |
| Total                                       | 50          | 325            | 375   |

#### Risk estimation in patients with pelvic mass

The effectiveness of the Elecsys HE4 assay in combination with the Elecsys CA 125 II assay for risk estimation of epithelial ovarian cancer of patients presenting with pelvic mass was determined in an international multi-center clinical trial using repository samples. An algorithm (ROMA = Risk of Ovarian Malignancy Algorithm) was developed for estimation of the risk of epithelial ovarian cancer. The algorithm takes into account the HE4 and CA 125 values as well as the menopausal status of the patient. The algorithm calculates a predictive probability of finding epithelial ovarian cancer on surgery.

#### Calculation of Predictive Index (PI)<sup>15</sup>

A Predictive Index is calculated for premenopausal and postmenopausal patients separately using equations (1) and (2) below. To calculate the PI, the assay values obtained from the Elecsys HE4 assay and the Elecsys CA 125 II assay are inserted into the equations below, depending on the menopausal status of the women.

(1) Premenopausal:

$$PI = -12.0 + 2.38 \cdot \ln[HE4] + 0.0626 \cdot \ln[CA125]$$

(2) Postmenopausal:

$$PI = -8.09 + 1.04 \cdot \ln[HE4] + 0.732 \cdot \ln[CA125]$$

where, LN = Natural Logarithm. Do not use LOG = Log<sub>10</sub>.

#### Calculation of ROMA value<sup>15</sup>

To calculate the ROMA value (i.e. predictive probability), insert the calculated value for PI into equation (3):

$$(3) \text{ ROMA value (\%)} = \exp(PI) / [1 + \exp(PI)] \cdot 100 \text{ where, } \exp(PI) = e^{PI}$$

NOTE: These equations were used for the calculation of ROMA values with the Elecsys HE4 assay from 28.8-3847 pmol/L and with the Elecsys CA 125 II assay from 6.42-5000 U/mL.

The examples below should be used in order to validate calculations of PI and ROMA before reporting patient results:

| Menopausal status | Elecsys values |                  | PI calculation                        | PI       | ROMA % |
|-------------------|----------------|------------------|---------------------------------------|----------|--------|
|                   | HE4 (pmol/L)   | CA 125 II (U/mL) |                                       |          |        |
| Premenopausal     | 37.5           | 74.9             | -12.0 + (2.38*3.624) + (0.0626*4.316) | -3.10388 | 4.29   |
|                   | 387            | 21.8             | -12.0 + (2.38*5.957) + (0.0626*3.082) | 2.371517 | 91.5   |
| Postmenopausal    | 66.7           | 11.3             | -8.09 + (1.04*4.200) + (0.732*2.425)  | -1.94683 | 12.5   |
|                   | 383            | 22.7             | -8.09 + (1.04*5.948) + (0.732*3.122)  | 0.381799 | 59.4   |

#### Stratification into low risk and high risk groups

In a study a total of 384 repository patient samples were included and the predictive probability for ovarian cancer as well as the ability for separation into a low and a high risk group based on ROMA values were determined.

The risk of ovarian malignancy algorithm was used to stratify women into risk groups for finding epithelial ovarian cancer. The following cut-points were used in order to provide a specificity level of 75 % for the Elecsys HE4 and Elecsys CA 125 II assay combination:

##### Premenopausal women

ROMA value  $\geq 11.4$  % = high risk of finding epithelial ovarian cancer

ROMA value < 11.4 % = low risk of finding epithelial ovarian cancer

##### Postmenopausal women

ROMA value  $\geq 29.9$  % = high risk of finding epithelial ovarian cancer

ROMA value < 29.9 % = low risk of finding epithelial ovarian cancer

The risk stratification of all 384 patients (194 pre- and 190 postmenopausal) presenting with pelvic mass using the ROMA values for the Elecsys HE4 and Elecsys CA 125 II assay combination is shown in the following table:

| Patient groups presenting with pelvic mass | Premenopausal patients |               |                    | Postmenopausal patients |               |                    |
|--|------------------------|---------------|--------------------|-------------------------|---------------|--------------------|
|  | N                      | ROMA < 11.4 % | ROMA $\geq 11.4$ % | N                       | ROMA < 29.9 % | ROMA $\geq 29.9$ % |
| Stage I-II EOC <sup>c)</sup>               | 16                     | 6 (37.5 %)    | 10 (62.5 %)        | 16                      | 6 (37.5 %)    | 10 (62.5 %)        |
| Stage I-IIIC EOC <sup>d)</sup>             | 21                     | 7 (33.3 %)    | 14 (66.7 %)        | 34                      | 9 (26.5 %)    | 25 (73.5 %)        |
| Stage I-IV EOC                             | 25                     | 7 (28.0 %)    | 18 (72.0 %)        | 53                      | 10 (18.9 %)   | 43 (81.1 %)        |
| Stage III-IV EOC                           | 9                      | 1 (11.1 %)    | 8 (88.9 %)         | 37                      | 4 (10.8 %)    | 33 (89.2 %)        |
| Unstaged EOC                               | 12                     | 2 (16.7 %)    | 10 (83.3 %)        | 44                      | 2 (4.5 %)     | 42 (95.5 %)        |
| Benign                                     | 157                    | 118 (75.2 %)  | 39 (24.8 %)        | 93                      | 71 (76.3 %)   | 22 (23.7 %)        |

c) EOC = epithelial ovarian cancer

d) Stage I-IIIB and Stage I-IIIC (omentum negative, lymph node positive) EOC

The sensitivity for stratifying patients with stage I-IV epithelial ovarian cancer into the high risk group was 84.3 % at the set specificity of 75 %, such that 75.6 % of women with benign pelvic mass were classified into the low risk group. The positive and negative predictive values were 64.9 % and 90 % respectively.

AUC (95 % CI):

Premenopausal women = 0.858 (0.779-0.937)

Postmenopausal women = 0.923 (0.885-0.962)

**The following has to be taken into consideration**

- The level of HE4 cannot be used as absolute evidence for the presence or absence of malignant disease and the Elecsys HE4 assay should not be used as a cancer screening test.
- Elecsys HE4 results should be used in conjunction with other clinical data; e.g. symptoms, medical history, etc.
- If the Elecsys HE4 results are inconsistent with clinical evidence, additional testing is suggested to confirm the result.
- Elecsys HE4 results should not be used interchangeably with other manufacturers' methods for HE4 determinations.
- Elecsys CA 125 II results should not be used interchangeably with other manufacturers' methods for CA 125 determinations in the ROMA calculation.
- Patients with confirmed ovarian cancer may have Elecsys HE4 assay values in the same range as healthy women. Certain histological types of ovarian cancer (e.g. mucinous or germ cell tumors) rarely express HE4, therefore the use of the Elecsys HE4 assay is not recommended for monitoring of patients with known mucinous or germ cell ovarian cancer.<sup>5</sup> Conversely, elevated levels of HE4 antigen may be present in individuals with renal, liver and non-malignant diseases.
- The ROMA has not been validated for the following patient groups: patients previously treated for malignancy, patients currently being treated with chemotherapy, and patients less than 18 years of age. Failure of the Elecsys HE4 assay and/or the Elecsys CA 125 II assay to perform as indicated, or error in the calculation of results, could lead to inaccurate risk assessment and improper management of the patient. Specifically, a falsely low result of the assay(s) could result in a determination that the patient is at lower risk of having epithelial ovarian cancer, which could triage the patient to a less specialized level of care.

**Specific performance data**

Representative performance data on the analyzers are given below. Results obtained in individual laboratories may differ.

**Precision**

Precision was determined using Elecsys reagents, samples and controls in a protocol (EP5-A2) of the CLSI (Clinical and Laboratory Standards Institute); 2 runs per day in duplication each for 21 days (n = 84). The following results were obtained:

| Elecsys 2010 and <b>cobas e 411</b> analyzers |             |               |      |                        |      |
|---|-------------|---------------|------|------------------------|------|
|   |             | Repeatability |      | Intermediate precision |      |
| Sample  | Mean pmol/L | SD pmol/L     | CV % | SD pmol/L              | CV % |
| Human serum 1                                 | 25.3        | 0.450         | 1.8  | 0.945                  | 3.7  |
| Human serum 2                                 | 53.7        | 0.988         | 1.8  | 2.28                   | 4.2  |
| Human serum 3                                 | 142         | 2.33          | 1.6  | 6.11                   | 4.3  |
| Human serum 4                                 | 779         | 11.3          | 1.5  | 32.6                   | 4.2  |
| Human serum 5                                 | 1437        | 18.9          | 1.3  | 39.4                   | 2.7  |
| PreciControl HE4 1                            | 45.7        | 0.661         | 1.4  | 1.92                   | 4.2  |
| PreciControl HE4 2                            | 345         | 5.67          | 1.6  | 11.8                   | 3.4  |

| MODULAR ANALYTICS E170, <b>cobas e 601</b> and <b>cobas e 602</b> analyzers |             |               |      |                        |      |
|---|-------------|---------------|------|------------------------|------|
|   |             | Repeatability |      | Intermediate precision |      |
| Sample  | Mean pmol/L | SD pmol/L     | CV % | SD pmol/L              | CV % |
| Human serum 1   | 27.4        | 0.481         | 1.8  | 0.798                  | 2.9  |
| Human serum 2   | 57.7        | 1.06          | 1.8  | 1.59                   | 2.8  |
| Human serum 3   | 155         | 2.28          | 1.5  | 3.95                   | 2.6  |
| Human serum 4   | 852         | 13.5          | 1.6  | 25.2                   | 3.0  |
| Human serum 5   | 1390        | 26.9          | 1.9  | 44.9                   | 3.2  |

| MODULAR ANALYTICS E170, <b>cobas e 601</b> and <b>cobas e 602</b> analyzers |             |               |      |                        |      |
|---|-------------|---------------|------|------------------------|------|
|   |             | Repeatability |      | Intermediate precision |      |
| Sample  | Mean pmol/L | SD pmol/L     | CV % | SD pmol/L              | CV % |
| PreciControl HE4 1  | 45.2        | 0.670         | 1.5  | 1.52                   | 3.4  |
| PreciControl HE4 2  | 345         | 6.06          | 1.8  | 11.2                   | 3.2  |

**Method comparison**

A comparison of the Elecsys HE4 assay (y) with a manual HE4 method (x) using clinical samples gave the following correlations:

Number of samples measured: 1502

|                              |                    |
|------------------------------|--------------------|
| Passing/Bablok <sup>16</sup> | Linear regression  |
| $y = 1.018x + 5.52$          | $y = 1.10x + 1.64$ |
| $r = 0.808$                  | $r = 0.978$        |

The sample concentrations were between approximately 22 and approximately 1150 pmol/L.

A comparison of the Elecsys HE4 assay (y) with an automated HE4 method (x) using clinical samples gave the following correlations:

Number of samples measured: 703

|                              |                     |
|------------------------------|---------------------|
| Passing/Bablok <sup>16</sup> | Linear regression   |
| $y = 0.900x + 6.41$          | $y = 0.875x + 11.0$ |
| $r = 0.836$                  | $r = 0.985$         |

The sample concentrations were between approximately 20 and approximately 1400 pmol/L.

**Analytical specificity**

The following cross-reactivities were found, tested with HE4 concentrations of 65 and 155 pmol/L:

| Proteins (WFDC family)                    | Concentration tested pmol/L | Cross-reactivity % |
|---|-----------------------------|--------------------|
| Elafin <sup>e)</sup> /SKALP <sup>f)</sup> | 54500                       | 0.025              |
| SLPI <sup>g)</sup>                        | 20833                       | 0.088              |

e) Elafin = elastase-specific inhibitor

f) SKALP = skin-derived antileukoproteinase

g) SLPI = secretory leucocyte protease inhibitor

**References**

- 1 Israeli O, Goldring-Avram A, Rienstein S, et al. In silico chromosomal clustering of genes displaying altered expression patterns in ovarian cancer. *Cancer Genet Cytogenet* 2005;160:35-42.
- 2 Bouchard D, Morisset D, Bourbonnais Y, et al. Proteins with whey-acidic-protein motifs and cancer. *Lancet Oncol* 2006;7:167-174.
- 3 Bingle L, Singleton V, Bingle CD. The putative ovarian tumour marker gene HE4 (WFDC2), is expressed in normal tissues and undergoes complex alternative splicing to yield multiple protein isoforms. *Oncogene* 2002;21:2768-2773.
- 4 Kirchhoff C. Molecular characterization of epididymal proteins. *Rev Reprod* 1998;3:86-95.
- 5 Drapkin R, von Horsten HH, Lin Y, et al. Human epididymis protein 4 (HE4) is a secreted glycoprotein that is overexpressed by serous and endometrioid ovarian carcinomas. *Cancer Res* 2005;65:2162-2169.
- 6 Hellström I, Raycraft J, Hayden-Ledbetter M, et al. The HE4 (WFDC2) protein is a biomarker for ovarian carcinoma. *Cancer Res* 2003;63:3695-3700.
- 7 Pickle LW, Hao Y, Jemal A, et al. A new method of estimating United States and state-level cancer incidence counts for the current calendar year. *CA Cancer J Clin* 2007;57:30-42.
- 8 Cancer Mondial Statistical Information System: <http://www-dep.iarc.fr/>

- 9 Earle CC, Schrag D, Neville BA, et al. Effect of surgeon specialty on processes of care and outcomes for ovarian cancer patients. J Natl Cancer Inst 2006;98:172-180.
- 10 Engelen MJ, Kos HE, Willemse PH, et al. Surgery by consultant gynecologic oncologists improves survival in patients with ovarian carcinoma. Cancer 2006;106:589-598.
- 11 Moore RG, Brown AK, Miller MC, et al. The use of multiple novel tumor biomarkers for the detection of ovarian carcinoma in patients with a pelvic mass. Gynecol Oncol 2008;108:402-408.
- 12 Moore RG, Brown AK, Miller MC, et al. Utility of a novel serum tumor biomarker HE4 in patients with endometrioid adenocarcinoma of the uterus. Gynecol Oncol 2008;110:196-201.
- 13 Huhtinen K, Suvitie P, Hiissa J, et al. Serum HE4 concentration differentiates malignant ovarian tumours from ovarian endometriotic cysts. Br J Cancer 2009;100:1315-1319.
- 14 Anastasi E, Marchei GG, Viggiani V, et al. HE4: a new potential early biomarker for the recurrence of ovarian cancer. Tumor Biol 2010;31:113-119.
- 15 Moore RG, McMeekin DS, Brown AK, et al. A novel multiple marker bioassay utilizing HE4 and CA125 for the prediction of ovarian cancer in patients with a pelvic mass. Gynecol Oncol 2009;111:40-46.
- 16 Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.

For further information, please refer to the appropriate operator's manual for the analyzer concerned, the respective application sheets, the product information and the Method Sheets of all necessary components (if available in your country).

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.



**FUJIREBIO®**  
Diagnostics, Inc.

ROMA is a registered trademark of Fujirebio  
Diagnostics, Inc.

HE4 antigens and antibodies used in the  
Roche HE4 products are licenced by  
Fujirebio Diagnostics, Inc.

## Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard.

|  |   |
|--|---|
|  | Contents of kit                                     |
|  | Analyzers/Instruments on which reagents can be used |
|  | Reagent   |
|  | Calibrator  |
|  | Volume after reconstitution or mixing               |

COBAS, COBAS E, ELECSYS, MODULAR and PRECICONTROL are trademarks of Roche. INTRALIPID is a trademark of Fresenius Kabi AB.

All other product names and trademarks are the property of their respective owners.

Significant additions or changes are indicated by a change bar in the margin.

© 2013, Roche Diagnostics



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim  
www.roche.com

