

N-MID Osteocalcin

Osteocalcin (OCN)

cobas®

REF		SYSTEM
12149133 122	100	Elecsys 2010 MODULAR ANALYTICS E170 cobas e 411 cobas e 601 cobas e 602

English

Intended use

Immunoassay for the in vitro quantitative determination of N-MID osteocalcin in human serum and plasma. The determination is used for the control of antiresorptives therapeutic efficiency, e.g. for patients with osteoporosis or hypercalcemia.

The electrochemiluminescence immunoassay "ECLIA" is intended for use on Elecsys and **cobas e** immunoassay analyzers.

Summary

References^{1,2,3,4,5,6,7,8,9,10}

Osteocalcin, the most important non-collagen protein in bone matrix, is a bone-specific, calcium-binding protein which is dependent on vitamin K. It contains 49 amino acids and has a molecular weight of approx. 5800 daltons. It contains up to three γ-carboxyglutamic acid residues (bone-GLA-protein, BGP).

During bone synthesis osteocalcin is produced by the osteoblasts. Its production is dependent upon vitamin K (formation of γ-carboxyglutamic acid residues) and is stimulated by vitamin D3. After release from the osteoblasts, osteocalcin is not only assimilated into the bone matrix but also secreted into the blood stream. Accordingly, the serum (plasma) osteocalcin level is related to the rate of bone turnover in various disorders of bone metabolism, e.g. osteoporosis in particular, but also in primary and secondary hyperparathyroidism or Paget's disease.

Osteocalcin is therefore termed a bone turnover marker and is used for this purpose. By means of osteocalcin measurements it is possible to monitor therapy with antiresorptive agents (bisphosphonates or hormone replacement therapy, HRT) in, for example, patients with osteoporosis or hypercalcemia.

Both intact osteocalcin (amino acids 1-49) and the large N-MID fragment (amino acids 1-43) occur in blood. Intact osteocalcin is unstable due to protease cleavage between amino acids 43 and 44. The N-MID-fragment resulting from cleavage is considerably more stable.

The Elecsys N-MID Osteocalcin assay uses two monoclonal antibodies specifically directed against epitopes on the N-MID-fragment and the N-terminal-fragment. The assay hence detects the stable N-MID-fragment as well as the (still) intact osteocalcin. The test is non-dependent on the unstable C-terminal-fragment (amino acids 43-49) of the osteocalcin molecule and thus ensures constant measurement results under routine conditions in the laboratory.

Test principle

Sandwich principle. Total duration of assay: 18 minutes.

- 1st incubation: 20 µL of sample, a biotinylated monoclonal N-MID osteocalcin-specific antibody, and a monoclonal N-MID osteocalcin-specific antibody labeled with a ruthenium complex^{a)} react to form a sandwich complex.
- 2nd incubation: After addition of streptavidin-coated microparticles, the complex becomes bound to the solid phase via interaction of biotin and streptavidin.
- The reaction mixture is aspirated into the measuring cell where the microparticles are magnetically captured onto the surface of the electrode. Unbound substances are then removed with ProCell/ProCell M. Application of a voltage to the electrode then induces chemiluminescent emission which is measured by a photomultiplier.
- Results are determined via a calibration curve which is instrument-specifically generated by 2-point calibration and a master curve provided via the reagent barcode.

a) Tris(2,2'-bipyridyl)ruthenium(II)-complex (Ru(bpy)₃)²⁺

Reagents - working solutions

The reagent rackpack is labeled as OSTEOC.

- M Streptavidin-coated microparticles (transparent cap), 1 bottle, 6.5 mL: Streptavidin-coated microparticles 0.72 mg/mL; preservative.
- R1 Anti-N-MID Osteocalcin-Ab-biotin (gray cap), 1 bottle, 10 mL: Biotinylated monoclonal anti-N-MID Osteocalcin antibody (mouse) 1.5 mg/L; phosphate buffer 100 mmol/L, pH 6.0; preservative.
- R2 Anti-N-MID Osteocalcin-Ab-Ru(bpy)₃²⁺ (black cap), 1 bottle, 10 mL: Monoclonal anti-N-MID Osteocalcin antibody (mouse) labeled with ruthenium complex 1.3 mg/L; phosphate buffer 100 mmol/L, pH 6.0; preservative.

Precautions and warnings

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents.

Disposal of all waste material should be in accordance with local guidelines. Safety data sheet available for professional user on request.

Avoid foam formation in all reagents and sample types (specimens, calibrators and controls).

Reagent handling

The reagents in the kit have been assembled into a ready-for-use unit that cannot be separated.

All information required for correct operation is read in from the respective reagent barcodes.

Storage and stability

Store at 2-8 °C.

Do not freeze.

Store the Elecsys reagent kit **upright** in order to ensure complete availability of the microparticles during automatic mixing prior to use.

Stability:	
unopened at 2-8 °C	up to the stated expiration date
after opening at 2-8 °C	12 weeks
on the analyzers	8 weeks

Specimen collection and preparation

Only the specimens listed below were tested and found acceptable.

Serum collected using standard sampling tubes.

Li-heparin and K₃-EDTA plasma.

Criterion: Recovery within 90-110 % of serum value or slope 0.9-1.1 + intercept within < ± 2x analytical sensitivity (LDL) + coefficient of correlation > 0.95.

Note: Avoid hemolysis! Erythrocytes contain proteases which degrade osteocalcin. It is recommended that blood be centrifuged immediately.

Stability of serum and heparinized plasma: 8 hours at 15-25 °C, 3 days at 2-8 °C, 3 months at -20 °C. Freeze once only.

Stability of EDTA-plasma: 2 days at 15-25 °C, 3 days at 2-8 °C, 3 months at -20 °C. Freeze only once.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

N-MID Osteocalcin

Osteocalcin (OCN)



Do not use heat-inactivated samples.

Do not use samples and controls stabilized with azide.

Ensure the samples, calibrators and controls are at 20-25 °C prior to measurement.

Due to possible evaporation effects, samples, calibrators and controls on the analyzers should be analyzed/measured within 2 hours.

Materials provided

See "Reagents – working solutions" section for reagents.

Materials required (but not provided)

- [REF] 11972111122, N-MID Osteocalcin CalSet, for 4 x 1 mL
- [REF] 05618860190, PreciControl Varia, for 2 x 3 mL each of PreciControl Varia 1 and 2
- [REF] 11732277122, Diluent Universal, 2 x 16 mL sample diluent or [REF] 03183971122, Diluent Universal, 2 x 36 mL sample diluent
- General laboratory equipment
- Elecsys 2010, MODULAR ANALYTICS E170 or **cobas e** analyzer

Accessories for Elecsys 2010 and **cobas e** 411 analyzers:

- [REF] 11662988122, ProCell, 6 x 380 mL system buffer
- [REF] 11662970122, CleanCell, 6 x 380 mL measuring cell cleaning solution
- [REF] 11930346122, Elecsys SysWash, 1 x 500 mL washwater additive
- [REF] 11933159001, Adapter for SysClean
- [REF] 11706802001, Elecsys 2010 AssayCup, 60 x 60 reaction vessels
- [REF] 11706799001, Elecsys 2010 AssayTip, 30 x 120 pipette tips

Accessories for MODULAR ANALYTICS E170, **cobas e** 601 and **cobas e** 602 analyzers:

- [REF] 04880340190, ProCell M, 2 x 2 L system buffer
- [REF] 04880293190, CleanCell M, 2 x 2 L measuring cell cleaning solution
- [REF] 03023141001, PC/CC-Cups, 12 cups to prewarm ProCell M and CleanCell M before use
- [REF] 03005712190, ProbeWash M, 12 x 70 mL cleaning solution for run finalization and rinsing during reagent change
- [REF] 03004899190, PreClean M, 5 x 600 mL detection cleaning solution
- [REF] 12102137001, AssayTip/AssayCup Combimagazine M, 48 magazines x 84 reaction vessels or pipette tips, waste bags
- [REF] 03023150001, WasteLiner, waste bags
- [REF] 03027651001, SysClean Adapter M

Accessories for all analyzers:

- [REF] 11298500316, Elecsys SysClean, 5 x 100 mL system cleaning solution

Assay

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

Resuspension of the microparticles takes place automatically prior to use. Read in the test-specific parameters via the reagent barcode. If in exceptional cases the barcode cannot be read, enter the 15-digit sequence of numbers.

MODULAR ANALYTICS E170, **cobas e** 601 and **cobas e** 602 analyzers: PreClean M solution is necessary.

Bring the cooled reagents to approx. 20 °C and place on the reagent disk (20 °C) of the analyzer. Avoid foam formation. The system automatically regulates the temperature of the reagents and the opening/closing of the bottles.

Calibration

Traceability: This method has been standardized against in-house reference standards: osteocalcin in analyte-free human serum matrix.

Every Elecsys reagent set has a barcoded label containing specific information for calibration of the particular reagent lot. The predefined master curve is adapted to the analyzer using the relevant CalSet.

Calibration frequency: Calibration must be performed once per reagent lot using fresh reagent (i.e. not more than 24 hours since the reagent kit was registered on the analyzer). Renewed calibration is recommended as follows:

- after 12 weeks when using the same reagent lot
- after 7 days (when using the same reagent kit on the analyzer)
- as required: e.g. quality control findings outside the defined limits

Quality control

For quality control, use PreciControl Varia.

In addition, other suitable control material can be used.

Controls for the various concentration ranges should be run individually at least once every 24 hours when the test is in use, once per reagent kit, and following each calibration.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

Calculation

The analyzer automatically calculates the analyte concentration of each sample (either in ng/mL or µg/L).

Limitations - interference

The assay is unaffected by icterus (bilirubin < 1112 µmol/L or < 65 mg/dL), lipemia (Intralipid < 1500 mg/dL) and biotin (< 205 nmol/L or < 50 ng/mL).

Criterion: Recovery within ± 10 % of initial value.

Hemolysis interferes. Erythrocytes contain proteases which degrade osteocalcin.

Samples should not be taken from patients receiving therapy with high biotin doses (i.e. > 5 mg/day) until at least 8 hours following the last biotin administration.

No interference was observed from rheumatoid factors up to a concentration of 2200 U/mL.

There is no high-dose hook effect at N-MID osteocalcin concentrations up to 4200 ng/mL.

In vitro tests were performed on 16 commonly used pharmaceuticals. No interference with the assay was found.

In rare cases, interference due to extremely high titers of antibodies to analyte-specific antibodies, streptavidin or ruthenium can occur. These effects are minimized by suitable test design.

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

Limits and ranges

Measuring range

0.500-300 ng/mL (defined by the lower detection limit and the maximum of the master curve). Values below the lower detection limit are reported as < 0.500 ng/mL. Values above the measuring range are reported as > 300 ng/mL (or up to 1500 ng/mL for 5-fold diluted samples).

Lower limits of measurement

Lower detection limit of the test

Lower detection limit: < 0.500 ng/mL

The lower detection limit represents the lowest measurable analyte level that can be distinguished from zero. It is calculated as the value lying two standard deviations above that of the lowest standard (master calibrator, standard 1 + 2 SD, repeatability study, n = 21).

Dilution

Samples with N-MID osteocalcin concentrations above the measuring range can be diluted with Diluent Universal. The recommended dilution is 1:5 (either automatically by the MODULAR ANALYTICS E170, Elecsys 2010 or **cobas e** analyzers or manually). The concentration of the diluted sample must be > 60 ng/mL).

After manual dilution, multiply the result by the dilution factor.

N-MID Osteocalcin

Osteocalcin (OCN)



After dilution by the analyzers, the MODULAR ANALYTICS E170, Elecsys 2010 and **cobas e** software automatically takes the dilution into account when calculating the sample concentration.

Expected values

The reference ranges are test-dependent. Completed studies (study protocol No. 9905 - 8/2000) with the Elecsys N-MID Osteocalcin assay have revealed the following ranges in ng/mL:

	Number	N-MID osteocalcin	
		50 th Perc.	5-95 th Perc.
Healthy women			
• Premenopausal, > 20 yrs.	108	23	11-43
• Postmenopausal (no HRT)	102	27	15-46
Osteoporosis patients	120	27	13-48
Healthy men			
• 18- < 30 yrs.	183	40	24-70
• 30-50 yrs.	179	25	14-42
• > 50-70 yrs.	125	24	14-46

In patients with renal failure the osteocalcin values can be elevated, both directly, due to impaired clearance and indirectly, due to renal osteodystrophy.¹¹

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

Specific performance data

Representative performance data on the analyzers are given below. Results obtained in individual laboratories may differ.

Precision

Precision was determined using Elecsys reagents, pooled human sera and controls in a protocol (EP5-A2) of the CLSI (Clinical and Laboratory Standards Institute): 2 runs per day in duplication each for 21 days (n = 84). The following results were obtained:

Elecsys 2010 and cobas e 411 analyzers					
		Repeatability		Intermediate precision	
Sample	Mean ng/mL	SD ng/mL	CV %	SD ng/mL	CV %
Human serum 1	6.01	0.085	1.4	0.186	3.1
Human serum 2	12.2	0.135	1.1	0.373	3.1
Human serum 3	35.6	0.601	1.7	1.06	3.0
Human serum 4	169	3.12	1.8	5.56	3.3
Human serum 5	8.11	0.091	1.1	0.159	2.0
PreciControl Varia 1	19.3	0.164	0.8	0.267	1.4
PreciControl Varia 2	93.2	1.01	1.1	1.65	1.8

MODULAR ANALYTICS E170, cobas e 601 and cobas e 602 analyzers					
		Repeatability		Intermediate precision	
Sample	Mean ng/mL	SD ng/mL	CV %	SD ng/mL	CV %
Human serum 1	6.11	0.056	0.9	0.120	2.0
Human serum 2	12.0	0.126	1.1	0.240	2.0
Human serum 3	34.5	0.361	1.0	0.677	2.0
Human serum 4	160	2.03	1.3	3.65	2.3
Human serum 5	7.49	0.066	0.9	0.107	1.4
PreciControl Varia 1	17.9	0.166	0.9	0.207	1.2

MODULAR ANALYTICS E170, **cobas e** 601 and **cobas e** 602 analyzers

		Repeatability		Intermediate precision	
Sample	Mean ng/mL	SD ng/mL	CV %	SD ng/mL	CV %
PreciControl Varia 2	85.9	0.755	0.9	1.12	1.3

Method comparison

A comparison of the Elecsys N-MID Osteocalcin assay (y) with a commercially available N-MID osteocalcin test (x) using clinical samples gave the following correlations (ng/mL):

Number of samples measured: 185

Passing/Bablok¹² Linear regression

$y = 1.29x - 2.79$ $y = 1.43x - 6.24$

$r = 0.866$ $r = 0.987$

The sample concentrations were between approx. 10 and 210 ng/mL.

Analytical specificity

For the monoclonal antibodies used, the following cross-reactivities were found:

No cross-reactivity detectable for β -CrossLaps, parathyroid hormone, and bone-specific alkaline phosphatase.

References

- Rosenquist C, Qvist P, Bjarnason N, et al. Measurement of a More Stable Region of Osteocalcin in Serum by ELISA with Two Monoclonal Antibodies. Clin Chem 1995;41(10):1439-1445.
- Delmas PD, Wahner HW, Mann KG, et al. Assessment of bone turnover in postmenopausal osteoporosis by measurement of serum bone Gla-protein. J Lab Clin Med 1983;102(4):470-476.
- Delmas PD, Wilson DM, Mann KG, et al. Effects of Renal Function on Plasma Levels of Bone Gla Protein. J Clin Endocrinol Metab 1983;57(5):1028-1030.
- Garnero P, Grimaux M, Seguin P, et al. Characterisation of Immunoreactive Forms of Human Osteocalcin Generated In Vivo and In Vitro. J Bone Miner Res 1994;9(2):255-264.
- Chen J-T, Hosoda K, Hasumi K, et al. Serum N-Terminal Osteocalcin is a Good Indicator for Estimating Responders to Hormone Replacement Therapy in Postmenopausal Women. J Bone Miner Res 1996;11(11):1784-1792.
- Garnero P, Delmas PD. New Developments in Biological Markers for Osteoporosis. Calcif Tissue Int 1996;59(1):2-9.
- Epstein S. Bone-Derived Proteins. Trends Endocrinol Metab 1989;1:9-14.
- Brown JP, Malaval L, Chapuy MC, et al. Serum bone GLA protein: A specific marker for bone formation in postmenopausal osteoporosis. Lancet 1984;1091-1093.
- Ravn P, Christensen JO, Baumann M, et al. Changes in Biochemical Markers and Bone Mass After Withdrawal of Ibendronate Treatment: Prediction of Bone Mass Changes During Treatment. Bone 1998;22(5):559-564.
- Gundberg CM. Biology, Physiology, and Clinical Chemistry of Osteocalcin. J Clin Ligand Assay 1998;21(2):128-138.
- Masters PW, Jones RG, Purves DA, et al. Commercial assays for serum osteocalcin give clinically discordant results. Clin Chem 1994;40(3):358-363.
- Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.

For further information, please refer to the appropriate operator's manual for the analyzer concerned, the respective application sheets, the product information and the Method Sheets of all necessary components (if available in your country).

N-MID Osteocalcin






Osteocalcin (OCN)



A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard.

 CONTENT	Contents of kit
 SYSTEM	Analyzers/Instruments on which reagents can be used
 REAGENT	Reagent
 CALIBRATOR	Calibrator
	Volume after reconstitution or mixing

COBAS, COBAS E, ELECSYS, MODULAR and PRECICONTROL are trademarks of Roche. INTRALIPID is a trademark of Fresenius Kabi AB.

All other product names and trademarks are the property of their respective owners.

Significant additions or changes are indicated by a change bar in the margin.

© 2013, Roche Diagnostics



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim
www.roche.com

