
Baxter

ST60/ST100/ST150 SETS WITH AN 69

Membrane for CRRT Powered By
PrisMax and **Prismaflex**

ST60 SET/ST100 SET/ST150 SET

The **Prismaflex** set is indicated for use only with the **Prismaflex** control unit or with the **PrisMax** control unit (in countries where **PrisMax** is cleared or registered) in providing continuous fluid management and renal replacement therapies. The system is intended for patients who have acute renal failure, fluid overload, or both. This set is intended for use in the following veno-venous therapies: SCUF, CVVH, CVVHD, CVVHDF.

PHYSICAL CHARACTERISTICS⁽¹⁾

	ST60 set	ST100 set	ST150 set
Membrane effective surface area	0.6 m ²	1 m ²	1.5 m ²
Fiber internal diameter (wet)	240 µm		
Fiber wall thickness	50 µm		
Blood volume in set	97 mL	155 mL	193 mL
Overall dimensions	27 x 22 x 9 cm	27 x 22 x 9 cm	27 x 22 x 9 cm
Weight	772 g	828 g	894 g
Minimal patient weight	11 kg	30 kg	30 kg

MATERIALS

AN 69 ST hollow fiber:
 - Acrylonitrile and sodium methallyl sulfonate copolymer
 - Surface treatment agent: Polyethylene imine

Filter housing and headers: Polycarbonate
 Filter potting compound: Polyurethane
 Tubing material: Plasticized polyvinyl chloride (PVC)
 Cartridge: PETG
 Sterilization mode: EtO (ethylene oxide)

OPERATING PARAMETERS

Maximum TMP (mmHg/kPa)	450/60		
Maximum blood pressure (mmHg/kPa)	500/66.6		
	ST60 set	ST100 set	ST150 set
Minimum blood flow rate	50 mL/min	75 mL/min	100 mL/min
Maximum blood flow rate	180 mL/min	400 mL/min	450 mL/min

PERFORMANCE SPECIFICATIONS⁽²⁾

Maximum ultrafiltration rate (mL/min)⁽³⁾
 (bovine blood, Hct 32%, Cp 60 g/L, 37°C)

	ST60 set			ST100 set			ST150 set			
QB (mL/min)	100	180	100	200	300	400	100	200	300	450
Max.QUF (± 20%)	39	56	45	70	91	109	52	82	106	136

Sieving coefficient
 (bovine plasma, Cp 60 g/L, 37°C)
 QB = 100 mL/min, QUF = 20 mL/min

• Urea	1									
• Vitamin B ₁₂	1									
• Inulin	0.96									

Sieving coefficient
 (human plasma, Cp 60 g/L, 37°C)
 QB = 100 mL/min, QUF = 20 mL/min

• Myoglobin	0.70									
• Albumin	<0.0045									

CLEARANCE

Clearance (mL/min) (saline solution ; 37°C)

Parameters:	ST60 set			ST100 set				ST150 set			
QB/QS	100 mL/min			150 mL/min				200 mL/min			
QUF	0 mL/min			0 mL/min				0 mL/min			
QD L/h	1	2.5	4	1	2.5	4	8	1	2.5	4	8
QD mL/min	17	42	67	17	42	67	133	17	42	67	133
Urea (±10%)	17	40	56	17	41	63	97	17	42	66	117
Vitamin B ₁₂ (±20%)	15	26	30	16	32	41	50	17	38	51	68
Inulin(±20%)	13	19	22	15	26	30	35	16	33	40	49

(1) Nominal values – given for indication

(2) Typical mean values obtained from laboratory testing of post-sterilization sample lots. Results may vary depending on patient and clinical conditions.

(3) Ultrafiltration is controlled by the control unit and is independent of the ultrafiltration coefficient (KUF).

ACRONYMS

TMP: Transmembrane pressure
 QB/QS: Arterial blood flow rate
 QUF: Ultrafiltration flow rate (fluid removal + replacement flow rate + pre blood pump flow rate)
 QD: Dialysate flow rate
 Hct: Hematocrit
 Cp: Protein concentration

ORDERING INFORMATION

	Code N°	N° units/box
ST60 set	107643	4
ST60 set CKT	115308	
ST100 set	107636	4
ST100 set CKT	115309	
ST150 set	107640	4
ST150 set CKT	115310	
5-liter effluent bag	115959 (SP414) or 114423 (A6001)	40 (SP414) 50 (A6001)
9-liter effluent bag	107650 (SP418)	30

For safe and proper use of the devices mentioned herein, please refer to the Instructions for Use or Operator's Manual.

www.baxter.com

Baxter Healthcare/Gambro Lundia AB, Magistratsvägen 16, SE-220 10 Lund, SWEDEN
 t: 46 46 16 90 00

Baxter, AN 69, Prismaflex and Prismax are trademarks of Baxter International Inc. or its subsidiaries.
 GBU-AT1-240001 (V1.0) 02/2024

Baxter

M60 / M100 / M150 AND ST60 / ST100 / ST150 SETS WITH AN 69

Membrane for CRRT Powered By
PrisMax and Prismaflex

The Prismaflex ST Set is Authorized only for the duration of the declaration that circumstances exist justifying the authorization of the emergency use of the Prismaflex ST Set under section 564(b)(1) of the Act, 21 U.S.C. § 360bbb-3(b)(1).

M60/M100/M150 AND ST60/ST100/ST150 SETS

The Prismaflex set is indicated for use only with the Prismaflex control unit or with the PrismaMax control unit (in countries where PrismaMax is cleared or registered) in providing continuous fluid management and renal replacement therapies. The system is intended for patients who have acute renal failure, fluid overload, or both. These sets are intended for use in the following veno-venous therapies: SCUF, CWH, CWHHD, CWHDF.

PHYSICAL CHARACTERISTICS ⁽¹⁾						
	M60 set	M100 set	M150 set	ST60 set	ST100 set	ST150 set
Membrane effective surface area	0.6 m ²	0.9 m ²	1.5 m ²	0.6 m ²	1 m ²	1.5 m ²
Fiber internal diameter (wet)	240 µm			240 µm		
Fiber wall thickness	50 µm			50 µm		
Blood volume in set	97 mL	155 mL	193 mL	97 mL	155 mL	193 mL
Overall dimensions	27 x 22 x 9 cm			27 x 22 x 9 cm		
Weight	758 g	814 g	874 g	772 g	828 g	894 g
Minimal patient weight	11 kg	30 kg	30 kg	11 kg	30 kg	30 kg

MATERIALS						
	M60 set	M100 set	M150 set	ST60 set	ST100 set	ST150 set
Hollow fiber	AN 69 HF hollow fiber: Acrylonitrile and sodium methallyl sulfonate copolymer			AN 69 ST hollow fiber: Acrylonitrile and sodium methallyl sulfonate copolymer Surface treatment agent: Polyethylene imine		
Filter housing and headers	Polycarbonate			Polycarbonate		
Filter potting compound	Polyurethane			Polyurethane		
Tubing material	Plasticized polyvinyl chloride (PVC)			Plasticized polyvinyl chloride (PVC)		
Cartridge	PETG			PETG		
Sterilization mode	EtO (ethylene oxide)			EtO (ethylene oxide)		

OPERATING PARAMETERS						
	M60 set	M100 set	M150 set	ST60 set	ST100 set	ST150 set
Maximum TMP (mmHg/kPa)	450/60			450/60		
Maximum blood pressure (mmHg/kPa)	500/66.6			500/66.6		
Minimum blood flow rate	50 mL/min	75 mL/min	100 mL/min	50 mL/min	75 mL/min	100 mL/min
Maximum blood flow rate	180 mL/min	400 mL/min	450 mL/min	180 mL/min	400 mL/min	450 mL/min

ORDERING INFORMATION		
	Code N°	N° units/box
M60 set	106696	4
M100 set	106697	4
M150 set	109990	4
ST60 set	107643	4
ST100 set	107636	4
ST150 set	107640	4
5-liter effluent bag	114423 (A6001)	50 (A6001)
9-liter effluent bag	107650 (SP418)	30

PERFORMANCE SPECIFICATIONS ⁽²⁾																				
Maximum ultrafiltration rate (mL/min) ⁽³⁾ (bovine blood; Hct 32%; Cp 60 g/L, 37°C)																				
	M60 set			M100 set			M150 set			ST60 set			ST100 set			ST150 set				
QB (mL/min)	100	180	100	200	300	400	100	200	300	450	100	180	100	200	300	400	100	200	300	450
Max.QUF (± 15%)	38	55	44	68	89	107	52	82	106	138	39	56	45	70	91	109	52	82	106	136
Sieving coefficient (bovine plasma, Cp 60 g/L, 37°C); QB = 100 mL/min; QUF = 20 mL/min																				
Urea	1						1						1							
Creatinine	1						None						None							
Vitamin B ₁₂	1						1						1							
Inulin	0.95						0.95						0.96							
Sieving coefficient (human plasma, Cp 60 g/L, 37°C); QB = 100 mL/min; QUF = 20 mL/min																				
Myoglobin	0.70						0.70						0.70							
Albumin	<0.0045						<0.0045						<0.0045							

CLEARANCE (mL/min) (saline solution; 37°C)																					
Parameters:	M60 set			M100 set			M150 set			ST60 set			ST100 set			ST150 set					
QB/QS	100 mL/min			150 mL/min			200 mL/min			100 mL/min			150 mL/min			200 mL/min					
QUF	0 mL/min			0 mL/min			0 mL/min			0 mL/min			0 mL/min			0 mL/min					
QD L/h	1	2.5	4	1	2.5	4	8	1	2.5	4	8	1	2.5	4	1	2.5	4	8	1	2.5	4
QD mL/min	17	42	67	17	42	67	133	17	42	67	133	17	42	67	17	42	67	133	17	42	67
Urea (±10%)	17	39	54	17	41	63	95	17	42	66	117	17	40	56	17	41	63	97	17	42	66
Vitamin B ₁₂ (±20%)	14	23	28	16	30	37	45	17	37	49	64	15	26	30	16	32	41	50	17	38	51
Inulin (±20%)	12	17	19	14	23	26	30	16	31	37	45	13	19	22	15	26	30	35	16	33	40

ACRONYMS	
TMP:	Transmembrane pressure
QB/QS:	Arterial blood flow rate
QUF:	Ultrafiltration flow rate (fluid removal + replacement flow rate + pre blood pump flow rate)
QD:	Dialysate flow rate
Hct:	Hematocrit
Cp:	Protein concentration

Rx Only. For the safe and proper use of the devices mentioned herein, please refer to the Instructions for Use, or Operator's Manual.

(1) Nominal values – given for indication (2) Typical mean values obtained from laboratory testing of post-sterilization sample lots. Results may vary depending on patient and clinical conditions. (3) Ultrafiltration is controlled by the control unit and is independent of the ultrafiltration coefficient (KUF).

PRISMAFLEX ST SET

Emergency Use Authorization for the United States

The **Prismaflex** ST Set has been Authorized by the FDA to provide continuous renal replacement therapy (CRRT) to treat patients in an acute care environment during the Coronavirus Disease 2019 (COVID-19) pandemic.

The **Prismaflex** ST Set has neither been cleared or approved to provide CRRT in an acute care environment.

The **Prismaflex** ST Set has been authorized by FDA under EUA200704.

The **Prismaflex** ST Set is Authorized only for the duration of the declaration that circumstances exist justifying the authorization of the emergency use of the **Prismaflex** ST Set under section 564(b) (1) of the Act, 21 U.S.C. § 360bbb-3(b)(1), unless the authorization is terminated or revoked sooner.

Intended Use for Patients with COVID-19

The **Prismaflex** ST Set is indicated for use only with the **Prismaflex** control unit or with the **PrisMax** control unit in providing continuous fluid management and renal replacement therapies in an acute care environment during the Coronavirus Disease 2019 (COVID-19) pandemic. The system is intended for patients who have acute renal failure, fluid overload, or both.

Relative contraindications (individual risk/benefit to be determined by treating physician) for the use of **Prismaflex** ST Sets include:

- The inability to establish vascular access
- Severe hemodynamic instability
- Known hypersensitivity to any component of the **Prismaflex** ST

This set is intended for use in the following veno-venous therapies: SCUF; CVVH; CVVHD; CVVHDF.

All treatments administered with the **Prismaflex** ST Set must be prescribed by a physician. The size, weight, metabolic and fluid balance, cardiac status, and general clinical condition of the patient must be carefully evaluated by the prescribing physician before each treatment.

Prismaflex ST Set

Emergency Use Authorization for the United States

The Prismaflex ST Set has been Authorized by the FDA to provide continuous renal replacement therapy (CRRT) to treat patients in an acute care environment during the Coronavirus Disease 2019 (COVID-19) pandemic.

The Prismaflex ST Set has neither been cleared or approved to provide CRRT in an acute care environment.

The Prismaflex ST Set has been authorized by FDA under EUA200704.

The Prismaflex ST Set is Authorized only for the duration of the declaration that circumstances exist justifying the authorization of the emergency use of the Prismaflex ST Set under section 564(b)(1) of the Act, 21 U.S.C. § 360bbb-3(b)(1), unless the authorization is terminated or revoked sooner.

Intended Use for Patients with COVID-19

The Prismaflex ST Set is indicated for use only with the Prismaflex control unit or with the PrisMax control unit in providing continuous fluid management and renal replacement therapies in an acute care environment during the Coronavirus Disease 2019 (COVID-19) pandemic. The system is intended for patients who have acute renal failure, fluid overload, or both.

Relative contraindications (individual risk/benefit to be determined by treating physician) for the use of Prismaflex ST Sets include:

- The inability to establish vascular access
- Severe hemodynamic instability
- Known hypersensitivity to any component of the Prismaflex ST Set

This set is intended for use in the following veno-venous therapies: SCUF; CVVH; CVVHD; CVVHDF.

All treatments administered with the Prismaflex ST Set must be prescribed by a physician. The size, weight, metabolic and fluid balance, cardiac status, and general clinical condition of the patient must be carefully evaluated by the prescribing physician before each treatment.

Additional Product Information for the United States

The tubing contained in the following product codes are labeled as DEHP-free in the IFU, this means that the product was not made with DEHP plasticizer: 107643, 115308, 107636, 115309, 107640 and 115310.

The tubing contained in the following product codes were made using DEHP plasticizer: 955468 and 955596.

The following codes do not include reference to the PrisMax control unit; however, all Prismaflex ST Sets may be used with the Prismaflex control unit or with the PrisMax control unit: 115308, 115309, 115310, 955468 and 955596.

To access COVID-19 Resources, product details, product use information, and the comprehensive Prismaflex Control Unit Operator's Manual and PrisMax Control Unit Operator's Manual please visit the Baxter Healthcare Acute Therapies website at <http://www.renalacute.com>

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
The **Prismaflex ST60/ST100/ST150** set is manufactured by **GAMBRO Industries, 7 avenue Lionel Terray, BP 126, 69883 MEYZIEU CEDEX, FRANCE.**

⚠ **Caution:** Federal law (USA) restricts this device to sale by or on the order of a physician.

DEFINITION OF EXPRESSIONS USED IN THIS MANUAL

In this document :

⚠ **“Warning”** indicates a hazardous situation which, if not avoided, could result in death or serious injury.

⚠ **“Caution”** indicates a hazardous situation which, if not avoided, could result in minor or moderate injury;

“Note” to give additional information.

SCUF: Slow Continuous UltraFiltration.

CVVH: Continuous Venous-Hemofiltration.

CVVHD: Continuous Venous-Hemodialysis.

CVVHDF: Continuous Venous-Hemodiafiltration.

Prdilution: addition of replacement fluid to the blood stream upstream to the filter.

Postdilution: addition of replacement fluid to the blood stream downstream to the filter.

"Control unit" refers to the **PrismaFlex** control unit, or to the **PrisMax** control unit (in countries where **PrisMax** is cleared or registered).

PRODUCT DESCRIPTION

- The **Prismaflex ST60 / ST100 / ST150** set is a disposable, extracorporeal circuit for use with the **PrismaFlex** control unit or with the **PrisMax** control unit (in countries where **PrisMax** is cleared or registered).
- The **Prismaflex ST60 / ST100 / ST150** set consists of a AN69 ST hollow fiber hemofilter/dialyzer* and tubing lines ; refer to the control unit operator manual drawing for details.
- The blood return line (blue-striped) is equipped with a Luer-lock connection** near the deaeration chamber, dedicated to the connection of authorized devices and accessories described in the control unit operator's manual.
- All line connectors are compatible with the ISO 594-1 and ISO 594-2 international standards concerning conical fittings.
- The fluid pathways of the **Prismaflex** set are guaranteed sterile and non pyrogenic.
- The **Prismaflex ST60 / ST100 / ST150** set is sterilized by ethylene oxide (EtO). Deaeration is such that EtO residuals comply with the standards in ISO 10993.
- Expiration date: please refer to product label.

* In this document the hemofilter/dialyzer will be referred to as "filter".

INTENDED USE / INDICATIONS

The **Prismaflex** set is indicated for use only with the **PrismaFlex** control unit or with the **PrisMax** control unit (in countries where **PrisMax** is cleared or registered) in providing continuous fluid management and renal replacement therapies. The system is intended for patients who have acute renal failure, fluid overload, or both.

This set is intended for use in the following veno-venous therapies : SCUF ; CVVH ; CVVHD ; CVVHDF.

All treatments administered via the **Prismaflex** set must be prescribed by a physician. The size, weight, state of uremia, cardiac status, and general physical condition of the patient must be carefully evaluated by the prescribing physician before each treatment.

CONTRAINDICATIONS

There are no known absolute contraindications to continuous renal replacement therapies.

For the following conditions a careful assessment of the individual risk/benefit ratio has to be made by the treating physician (relative contraindications):

- inability to establish vascular access,
- severe hemodynamic instability,
- known hypersensitivity to any component of the **Prismaflex** set.

CAUTIONS AND WARNINGS

Note: refer to the control unit user interface and operator's manual for additional cautions and warnings.

⚠ Cautions

- Particular attention must be paid to extra corporeal blood volume with respect to patient size. Consider the sum of the **Prismaflex** set blood volume (refer to "Specifications") plus the blood volume of any accessory or device if used. The **Prismaflex ST60** set should be restricted to patients with a body weight greater than 11kg (24lb). The **Prismaflex ST100** set and **ST150** set should be restricted to patients with a body weight greater than 30kg (66lb).
- If the patient is not immediately connected to the **Prismaflex** set after priming is complete, flush the set with at least 1 000 mL priming solution [saline or alkaline solution (pH ≥ 7.3), with or without heparin added according to usual institutional practice] prior to connecting the patient. This requires use of a new bag of priming solution.
- When not using the pre blood pump infusion line, it is recommended to clamp this line close to its connection to the access blood tubing line; this will prevent the sedimentation of blood into the pre-blood pump infusion line.
- Using the **Prismaflex** set with blood flow rates lower than the recommended minimum values (see "Operating Parameters" section) may impair filter performance due to hemoconcentration, or to increased risk of coagulation.
- Since drugs can be removed by the membrane of the filter, the dosage of associated drug treatments may need to be adjusted for patients on continuous renal replacement therapy. Monitoring of blood drug levels of relevant compounds should be performed. The removal of other water-soluble compounds (e.g. vitamins, trace elements) during therapy also requires clinical consideration.

⚠ Warnings

- Carefully read these instructions for use and the control unit operator's manual before using this product.
- The use of operating procedures other than those published by the manufacturer or the use of accessory devices not recommended by the manufacturer can result in patient injury or death.
- Store the **Prismaflex** set in a dry place, between 0° C (32° F) and 30° C (86° F).
- Do not use this set if the packaging is damaged, if the sterilization caps are missing or loose, or if any of the lines in the set are kinked.
- To prevent contamination, this **Prismaflex** set must be used as soon as its packaging and sterilization caps are removed.
- Do not try to remove the filter from the cartridge plate.
- Use aseptic techniques when handling all blood and fluid lines in the set.
- Prismaflex** sets are compatible with the usual disinfection agents used for aseptic setup; however solvents and other chemicals, if used in contact with the product, could damage the set.
- During priming and operation, observe closely for leakage at joints within the set, and connections to other approved accessories and bags. Leakage can cause blood loss, fluid imbalance or air embolism. If a leakage is detected at a Luer connection and cannot be stopped by tightening the connections, or if leakage occurs at any other location, replace the set.
- Tightening Luer connections with an excessive force can damage the connectors.
- Once priming is complete the set's blood circuit will still contain heparinised saline solution. Depending on the level of the patient's bleeding risk the physician must decide if an additional priming using 500 mL **non heparinised** saline solution is necessary.
- Do not allow air to enter the blood compartment of the filter after priming is started. If a large amount of air enters, the set must be replaced.
- Should acute allergic reactions (first-use syndrome) occur in patients receiving treatment, **immediately stop the treatment** and administer appropriate intervention. Pay special attention to patients receiving ACE inhibitors and/or having already shown similar allergic reactions (see "Hypersensitivity Reactions" section).
- Use a 21-gauge or smaller needle to obtain blood/fluid samples or remove trapped air from the **Prismaflex** set. Use of larger needles can cause holes in the sample sites, resulting in external leak or air intake.
- External blood leakage may not be immediately identified by monitoring equipment and could result in significant blood loss. Check the filter and all connections of the disposable tubings during treatment to minimize the risk of leakage.

16. To assure adequate filter performance, it is recommended that the set be changed every 24 hours of use. However, the set must be changed after 3 days (72 hours). Continued use beyond this limit could result in rupture of the pump segments, with risk of patient injury or death.
17. Destroy this set after single use, using aseptic technique for potentially contaminated equipment and following local regulation for disposal. Do not re-sterilize. The Prismaflex set is intended for single use only. Re-using the Prismaflex set may cause serious damage to the product resulting in patient injury or death.
18. Use only drugs compatible with plastics listed in the specifications section. Some plastics can be incompatible with drugs when in contact with solutions with pH > 10.

SPECIFICATIONS

See Tables at end of document.

SET MATERIALS

AN69 ST hollow fiber	: Acrylonitrile and sodium methallyl sulfonate copolymer + Polyethylene Imine (surface treatment agent)
Housing and headers	: Polycarbonate
Potting compound	: Polyurethane
Tubing material	: Plasticized polyvinyl chloride (PVC)
Cartridge	: Polyethyleneterephthalate Glycol

Note : the following information is available from the manufacturer upon request :

- information about test methods used to obtain performance characteristics,
- the number and range of particles in the effluent from the dialyzer prepared as recommended for clinical use,
- the types and amounts of residue from the sterilization process.

Note: the **Prismaflex** set is not made with rubber natural latex.

Note: all fluid pathways in direct or indirect blood contact are DEHP-free.

INSTRUCTIONS FOR USE

Note: use the set by following the detailed on-line instructions provided by the control unit. Additional information is available in the control unit operator's manual.

Load Set

Install the set onto the control unit using the photographs on the inside cover as a guide - the same procedure applies for both **Prismaflex** and **PrisMax** control units (in countries where **PrisMax** is cleared or registered).

Prepare and Connect Solutions

In order to gain full benefit from the AN69ST in terms of improvement of hemocompatibility, it is recommended to add 5000IU of unfractionated heparin per liter of priming/rinsing solution. This procedure allows the adsorption of active heparin onto the AN69ST before the start of extracorporeal circulation.

Consequently, the systemic anticoagulation strategy during treatment will be adapted with respect to patient specificity. In the cases where priming/rinsing without addition of unfractionated heparin, we recommend infusing the loading dose of heparin to the patient 2 to 5 minutes before connection to the filter.

Hang bag of priming solution [saline or alkaline solution (pH ≥ 7.3) with added 5000IU unfractionated heparin/liter according to usual institutional practice, correctly homogenised] on priming hook. Connect access (red)/effluent (yellow) Y-line to priming solution bag.

SPECIAL PROCEDURES IN CASE OF COMPLICATION

External Blood Leaks

Note: see Warning no. 15.

If an external blood leakage is observed, immediately stop the blood pump. Initiate corrective action by securing connections or replacing the Prismaflex set.

If necessary, administer adequate replacement solution to the patient to compensate for blood loss.

Hypersensitivity Reactions

Note: see Warning no. 13.

Should acute allergic reactions (first use syndrome) occur within the first few minutes of the treatment, it is important to react immediately by discontinuing the session and administering appropriate treatment.

Adverse reactions may occur due to the complex interaction between blood and the artificial surfaces of the entire extracorporeal circuit. These reactions may also be precipitated and/or exacerbated by other external factors involved with the individual patient's specific disease process and the treatment of renal insufficiency. Certain types of adverse reaction may occur due to operational factors associated with the treatment. Therefore, proper management of the fluid removal, electrolyte balance, anticoagulation and blood flow rate as well as monitoring of the overall treatment parameters are essential to avoid side-effects which may be associated with hemodialysis/hemofiltration therapies.

Hypersensitivity reactions have been observed during dialysis. Symptoms of a hypersensitivity reaction may be gastrointestinal, mucocutaneous, respiratory, cardiovascular or systemic in nature and range from very mild to severe. Such symptoms have been described as anaphylactic-like reactions within the first few minutes. Manifestations include nausea, malaise, weakness, a sensation of burning or heat throughout the body, profuse perspiration, respiratory distress and in some instances hypotension and cardiopulmonary arrest. Should a combination of such symptoms appear, particularly at the start of the treatment session, it is important to react immediately by discontinuing the session and administering appropriate treatment. Blood in the extracorporeal circuit must not be returned to the patient.

Extra care must be taken when treating patients who have exhibited possible hypersensitivity symptoms during previous treatments, or patients who have a history of being highly sensitive and allergic to a variety of substances. A physician must be consulted to evaluate the risk and prescribe the appropriate precautions if a possible sensitivity is suspected.

The following factors are considered essential to minimize the risk of hypersensitivity reaction and other side effects:

- Strict adherence to the set-up, priming and rinsing procedures detailed in the manufacturer's instructions for use.
- Setting up and monitoring the treatment operating parameters according to the manufacturer's recommendations specified for each type of **Prismaflex** set and to the patient's needs and tolerance.
- Strict adherence to all WARNINGS and CAUTIONS given by the manufacturer in the instructions for use.

Patients receiving angiotensin converting enzyme (ACE) inhibitors as medication can develop, within the first few minutes of a treatment, symptoms similar to acute allergic reactions i.e bronchospasm, edema of airways or larynx, dyspnea, angioedema, urticaria, nausea, vomiting, diarrhea, respiratory arrest, abdominal cramping, hypotension, hypovolemic shock and death.

However, for these patients, administration of antihistamines often does not alleviate the symptoms. In this case, treatment must be stopped and a more aggressive first-line therapy for an anaphylactoid reaction should be initiated immediately after the onset of symptoms.

Therefore, the prescribing physician must pay special attention to patients receiving ACE inhibitors and/or having already shown similar reactions.

WARRANTY AND LIMITATION OF LIABILITY

- a) The manufacturer warrants that the **Prismaflex** set has been manufactured in accordance with its specifications and in compliance with good manufacturing practices, other applicable industry standards and regulatory requirements.
If provided with the lot/serial number of the defective product, the manufacturer will, by replacement or credit, remedy manufacturing defects in the **Prismaflex** set becoming apparent before the expiration date.
- b) The warranty under paragraph a) above is in lieu of, and to the exclusion of, any other warranty, whether written or oral, express or implied, statutory or otherwise, and there are no warranties of merchantability or other warranties, which extend beyond those described in paragraph a) above. The remedy set out above for manufacturing defects is the sole remedy available to any person due to defects in the **Prismaflex** set and the manufacturer shall not be liable for any consequential or incidental loss, damage, injury or expense arising directly or indirectly from the use of the **Prismaflex** set, whether as a result of any defect therein or otherwise.
- c) The manufacturer shall not be liable for any misuse, improper handling, non-compliance with warnings and instructions, damage arising from events after the manufacturer's release of the **Prismaflex** set, failure or omission to inspect the **Prismaflex** set before use in order to ensure that the **Prismaflex** set is in proper condition, or any warranty given by independent distributors or dealers.
- d) The manufacturer is GAMBRO Industries, 7 avenue Lionel Terray, BP 126, 69883 MEYZIEU CEDEX, FRANCE.

Indonesia - Distributor : PT. Tawada Healthcare, Rukan Permata Senayan Blok A No. 18-19, Jl. Tentara Pelajar No. 5, Grogol Utara, Jakarta Selatan - Indonesia

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⁽¹⁾ Nominal values - given for indication

⁽²⁾ Typical mean values obtained from laboratory testing of post-sterilization sample lots.
Results may vary depending on patient and clinical conditions.


⁽³⁾ Ultrafiltration is controlled by the control unit and is independent of the ultrafiltration coefficient (KUF).
Note : a TMP > 40 kPa (300 mmHg) does not allow a higher ultrafiltration.




Prismatherm

SP420 line connexion

INSTALL SET DIRECTLY ONTO THE CONTROL UNIT

Patented ready to use PrismaFlex set packaging.