Baxter

PrisMax

ACUTE CARE SYSTEM

MAXIMISE NURSING **EFFICIENCY** WITH THE AUTO-EFFLUENT DRAIN ACCESSORY


powered by the **PrisMax** system

The **PrisMax** System with the optional Auto-Effluent Drain Accessory eliminates the need to change effluent bags during CRRT, freeing up nursing time and enabling more time for patient care.¹

Save Valuable Time

Without the optional Auto-Effluent Drain Accessory, as often as every 1-2 hours, an ICU nurse will have to change effluent bags for each patient on CRRT:

Making Treatment Delivery More Efficient

The unique Auto-Effluent Drain Accessory is designed to reduce treatment interventions, freeing up valuable nursing time.¹


By eliminating the manual task of changing effluent bags as often as every 1-2 hours, users may **save an average of more than 30 minutes** per patient per day in nursing time.^{2,*}

Gain up to 30 minutes additional therapy time for achieving prescription targets by decreasing therapy downtime due to effluent bag changes.^{2,*}

Bag changes may be reduced by 60% compared to similar CRRT systems.²

*Based on data from a pilot evaluation of 305 treatments with the **PrisMax** system at eight ICUs in six countries

"The benefit goes beyond the actual number of minutes saved because the nurse can focus on more important things in a complex patient scenario."³


The Auto-Effluent Drain Accessory

The Auto-Effluent Drain Accessory continuously pumps effluent collected during therapy into a drain, while providing the accuracy of a scale-based measurement system. Eliminating effluent bag changes minimises user exposure to effluent.¹

Using a hook and line system having up to an 18 m (60 ft) range,^{*} the Auto-Effluent Drain Accessory works with a range of room layouts and drains.^{1,†}

The Auto-Effluent Drain Accessory can remain in place for up to six days, even when changing CRRT therapy modes or blood filter sets.¹

8. **An extension line may be needed depending on the distance between the device and drain**
 †Disposal should be managed in accordance with your hospital's infection control policy

- 7.
- | | |
|---------------------|--------------------------|
| 1. Cartridge | 5. Effluent tube |
| 2. Drain hook | 6. Pinch clamps |
| 3. Effluent bag, 1L | 7. Female luer connector |
| 4. Effluent bag, 5L | 8. Check valve |

- The Auto-Effluent Drain Accessory eliminates the need to change effluent bags during treatment by alternating effluent fill and drain between the 1 litre and 5 litre bags.
- Effluent flows into one bag while the other bag drains.
- When a bag fills to 1 litre, the pinch valves change position automatically to begin draining that bag and filling the other.
- Effluent draining can be paused as needed to allow for changes in patient location or drain point.
- The one-way check valve at the proximal end of the drain line prevents effluent backflow.

For more information, please speak with a Baxter representative or visit www.baxter.com.

For safe and proper use of products mentioned herein, please refer to the appropriate Operators Manual or Instructions for Use.

References:

1. Baxter – PrisMax Operators Manual. AW8005 Rev A. 2017; Sep. 2. Baxter – PrisMax Limited Controlled Distribution Report #ER6877 Rev A. Data on file. 2018. 3. Broman M, Bell M, Joannes-Boyau O, Ronco C, The Novel PrisMax Continuous Renal Replacement Therapy System in a Multinational, Multicentre Pilot Setting. Blood Purif 2018;;220-227..

Baxter.dk, Baxter.fi, Baxter.no, Baxter.se

Baxter, Prismaflex, Prismax and Thermax are trademarks of Baxter International Inc. or its subsidiaries.


Baxter Confidential - INTERNAL USE ONLY		
Part Number: AW6254	Date: 19OCT23	Proofread No.: P03
Designer: R.B	Page: 1 of 28	
Colour Reference:	Black	

prisMax


REF 115370  A6003 - Auto Effluent Accessory

REF 115485 A6005 - Auto Effluent Extension

eng	Instructions For Use
fre	Mode d'emploi
ger	Gebrauchsanweisung
ita	Istruzioni per l'uso
swe	Bruksanvisning
dut	Gebruiksaanwijzing
dan	Brugervejledning
fin	Käyttöohjeet
nor	Bruksanvisning
spa	Instrucciones de uso
br-pt	Instruções de uso
gre	Οδηγίες χρήσης
tur	Kullanma Talimatları
rus	Инструкции по эксплуатации
pol	Instrukcja użytkowania
por	Instruções de Utilização
cze	Návod k použití
rom	Instrucțiuni de utilizare
vie	Hướng dẫn sử dụng
hrv	Upute za uporabu
hun	Használati útmutató
ind	Petunjuk Penggunaan
lit	Naudojimo instrukcijos

A6003

A6005

Baxter Confidential - INTERNAL USE ONLY		
Part Number: AW6254	Date: 190CT23	Proofread No.: P03
Designer: R.B	Page: 3 of 28	
Colour Reference:	Black	

ENGLISH 

For more information please contact:

 Gambro UF Solutions, Inc.
7601 Northland Drive, Suite 170
Brooklyn Park MN 55428
USA
800-525-2623

www.baxter.com

See product label for country of origin

 Gambro Industries
7, Avenue Lionel Terray
69883 Meyzieu Cedex
FRANCE

Definitions of Symbols Used on Labeling of Product

Symbol	Title of Symbol	Description of Symbol
	Do not re-use	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
	Use-by date	Indicates the date after which the medical device is not to be used.
	Date of manufacture	Indicates the date when the medical device was manufactured.
	Batch Code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified.
	Temperature limit	Indicates the temperature limits to which the medical device can be safely exposed.
	Do not use if package is damaged	Indicates a medical device that should not be used if the package has been damaged or opened.
	Sterile fluid path. Sterilized using ethylene oxide	Indicates the presence of a sterile fluid path within the medical device. The method of sterilization is indicated within the box. (Auto Effluent Accessory)
	Non-Sterile	Indicates a medical device that has not been subjected to a sterilization process. (Auto Effluent Extension)
	Green dot certification	European trademark symbol that indicates that the manufacturer contributes to the cost of recovery and recycling.
	Manufacturer	Indicates the medical device manufacturer.
	Authorized representative in the European Community	Indicates the authorized representative in the European community.
	Prescription only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Caution	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
	Fragile, handle with care	Indicates a medical device that can be broken or damaged if not handled carefully.
	Keep dry	Indicates a medical device that needs to be protected from moisture.
	Humidity limitation	Indicates the range of humidity to which the medical device can be safely exposed.
	Biological risks	Indicates that there are potential biological risks associated with the medical device.
	Consult instructions for use	Indicates the need for the user to consult the instructions for use.
	Medical device	Indicates the item is a medical device.
	Unique Device Identifier	Indicates the Unique Device Identifier information.
	This way up	To indicate correct upright position of the transport package.

Definitions of Expressions Used in this Manual

 **“Warning”** indicates a hazardous situation which, if not avoided, could result in death or serious injury.

 **“Caution”** indicates a hazardous situation which, if not avoided, could result in minor or moderate injury.

“Note” is used as a reminder to the user/operator on normal treatment activity and on what is suitable action in a particular situation.

Product Description

- The Auto Effluent Accessory is a disposable circuit for use with the PrisMax monitor in combination with other Prismaflex or PrisMax disposable sets.
- The Auto Effluent Accessory consists of a cartridge, drain attachment, 2 bags, and several tubing lines and connectors.
- The Auto Effluent Extension consists of a tubing line with connectors.
- The Auto Effluent Accessory is a sterile disposable product free of toxics.
- The Auto Effluent Extension is non-sterile.
- All line Connectors are compatible with the international standards concerning conical fittings.
- The fluid pathways of the Auto Effluent Accessory are guaranteed sterile and non-pyrogenic.
- The Auto Effluent Accessory is sterilized by ethylene oxide (EO).
- Expiration date: Please refer to product label.
- Not made with natural rubber latex.

Use of the Auto Effluent Accessory in the PrisMax system requires a specific panel on the lower left side of the control unit; contact your local Baxter representative.

Intended Use and Indications for Use

The Auto Effluent Accessory is indicated for use only with the PrisMax control unit in providing automatic effluent draining during Continuous Renal Replacement Therapy (CRRT).

The Auto Effluent Accessory is intended for application up to 560,000 pump revolutions (approx. 6 days at 3,500mL/hr)

Auto Effluent Extension is a non-sterile accessory to the Auto Effluent Accessory. It is intended to extend the length of the drain line of the Auto Effluent Accessory in cases where a drain is longer than twenty feet (6 meters) away.

Contraindications

There are no known contraindications when used as intended and within the specified operating conditions.

Cautions and Warnings

Note: Additional warnings and cautions pertaining to the PrisMax system are included in the PrisMax Monitor or Baxter Control Unit Operator's Manual.

Cautions

- Read the PrisMax Operator's Manual carefully prior to use. 
- Read the IFU of the CRRT Therapy Set carefully prior to use. 
- Read the IFU of the Auto Effluent Accessory and Auto Effluent Extension carefully prior to use.
- Use aseptic technique according to hospital policy.
- Dispose of the product after use according to clinical and legal regulations.
- Do not expose to heat or sun.
- The product should be stored in a cool, dry and dark place.
- Do not use solvents on the product such as alcohol, ether, acetone, liquid inhalation anaesthesia (e.g. halothan, enfuran, etc.). They can damage the products from the inside and outside.
- In order to avoid product damage, do not use lubricants, disinfectants, etc. for tube connections.
- Tubes and Luer Lock connections have to be fixed and fastened by hand only. Other aids can damage the connectors.
- Once the product has been primed it should be used within hours to minimize the danger of microbiological contamination and/or precipitation in priming fluid.
- Rx Only.
- Auto Effluent Accessory and Auto Effluent Extension is intended for single use only. Reprocessing this extension tubing may result in patient or user hazard due to contamination. 
- Ensure there is an air gap between the Auto Effluent outlet and the fluid of the sink/drain.

Warnings

- The use of operating procedures, other than those published by the manufacturer, or the use of other devices not recommended by the manufacturer can result in patient injury or death.
- Sterility is only guaranteed if the packaging is undamaged and unopened.
- Pay attention to the expiration date.
- The product must be carefully checked by the user before application for damages and intactness. Do not use if package is damaged, sterilization caps are missing, or lines are kinked.
- For single use only. Do not resterilize.
- All products have to be used under aseptic conditions.
- This product should be used only for medical applications and should only be used under medical supervision and by suitably trained operators.
- This product must be used as soon as the packaging is removed.

Specifications

Accessory Materials

Component	Material
1-Cartridge	Polyethylene Terephthalate Glycol (PETG)
2-Drain Hook	High-Density Polyethylene (HDPE)
3-Bags, 1L(a) & 5L(b)	Polyvinyl Chloride (PVC)
4-Tubing	Polyvinyl Chloride (PVC)
5-Pinch Clamps	Polypropylene (PP)
6-Check Valve	Acrylonitrile Butadiene Styrene (ABS)
7-Fittings	PVC, PETG, HDPE, PP

Instructions For Use

Note: Use the Auto Effluent Accessory and the Auto Effluent Extension by following the detailed on-screen instructions provided by the PrisMax monitor. Additional information is available in the PrisMax monitor Operator's Manual.

Warranty and Limitation of Liability

- A) The manufacturer warrants that the PrisMax accessory has been manufactured in accordance with its specifications and in compliance with good manufacturing practices, other applicable industry standards and regulatory requirements. If provided with the lot/serial number of the defective product, the manufacturer will, by replacement or credit, remedy manufacturing defects in the PrisMax accessory becoming apparent before the expiration date.
- B) The warranty under paragraph A) above is in lieu of, and to the exclusion of, any other warranty whether written or oral, express or implied, statutory or otherwise, and there are no warranties of merchantability or other warranties, which extend beyond those described in paragraph A) above. The remedy set out above for manufacturing defects is the sole remedy available to any person due to defects in the PrisMax accessory, and the manufacturer shall not be liable for any consequential or incidental loss, damage, injury, or expense arising directly, or indirectly, from the use of the PrisMax accessory, whether as a result of any defect therein or otherwise.
- C) The manufacturer shall not be liable for any misuse, improper handling, non-compliance with warnings and instructions, damage arising from events after the manufacturer's release of the PrisMax accessory, failure or omission to inspect the PrisMax accessory before use in order to ensure that the PrisMax accessory is in proper condition, or any warranty given by independent distributors or dealers.

Serious Adverse Event Reporting:

Per Regulation (EU) 2017/745 on medical devices: For a user and/or patient if, during the use of this device, or as a result of its use, a serious incident has occurred, please report this incident to the manufacturer, and/or its authorized representative, and to the competent authority of the Member State in which the user and/or patient is established.

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Baxter Confidential - INTERNAL USE ONLY	
Part Number: AW6254	Date: 19OCT23 Proofread No.: P03
Designer: R.B	Page: 26 of 28
Colour Reference:	Black

This page intentionally left blank.

Baxter Confidential - INTERNAL USE ONLY	
Part Number: AW6254	Date: 19OCT23 Proofread No.: P03
Designer: R.B	Page: 27 of 28
Colour Reference:	Black

This page intentionally left blank.

Baxter Confidential - INTERNAL USE ONLY	
Part Number: AW6254	Date: 19OCT23 Proofread No.: P03
Designer: R.B	Page: 28 of 28
Colour Reference:	Black