
mont blanc[®]

mont blanc surgical technique

YOU LEAD
WE FIND **the way**

The mont blanc[®] spinal system is a top-loading polyaxial screw system with a secured locking mechanism, a lower functional profile and a smart, intuitive and customisable instrumentation.


mont blanc[®] is a comprehensive system for thoracic, lumbar and sacral spine that has been designed for treatment of degenerative disease, deformity, trauma and tumor indications.

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
Pedicle screw

Buttress sriegis ir torx formos atsuktuvus

Buttress thread and torx footprint

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
Top loading design
Vis "Top loading"
Top loading design

Color coding to identify
screw diameter
Code couleur pour
l'identification du
diamètre de la vis
Código de color para
identificar el
diámetro del tornillo

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
Pre-operatively, mont blanc[®] pre-bent rods template can be a support to confirm the lordosis of the instrumented segment(s).

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


Instrument overview

Screwdriver

4 screwdriver shafts for 4 screws design.

Screwdriver Assembly/Disassembly can be obtained by pushing the knob button.

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

Patient positioning

Position of the patient

The patient is generally positioned in ventral decubitus avoiding compression points.

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
Screw entry points

Landmarks

Determine entry point and position of pedicle screws

• Thoracic spine

Entry points are just below the rim of the upper facet joint. The screws should be inserted at an angle of 7–10° towards the midline and 10–20° caudally.

• Lumbar spine

Entry points are at the intersection of a vertical line tangential to the lateral border of the superior articular process and the horizontal line bisecting the transverse processes.

The screws should be inserted at an angle of 5–10° towards the midline at the thoracolumbar junction. They should converge by 10° at L2 and increase to 15° at L5.

• Sacrum

At S1 level, the entry point is at the infero-lateral margin of the basis of the superior articular process of the sacrum. The direction of the screws follow the axis of the pedicle. It is slightly oblique towards the midline (on a average 10–15°).

• Iliac

Entry point for iliac are located 1cm Cephalad to Distal Posterior Superior Iliac Spine (PSIS) and 1cm inferior to the superior PSIS on the medial slope.

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
After having determined the entry point of the pedicle, create a cortical breach into the base of the pedicle by using the sharp tip and its depth adjustable ring.

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Note 1: Sharp tip with tap (optional) can be used to create a cortical breach into the base of the pedicle and prepare the pedicle thread.

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Sharp tip with tap is not available for Ø4 and Ø5 pedicle screws

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Note 2: Sharp tip with stop (optional) can be used to create a cortical breach into the base of the pedicle with a fix depth: 12mm.

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
Introduce suitable probe with care to avoid perforation of the medial cortex, to gently create a path in the pedicle. Graduation shaft allowing for more accurate depth control.

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
Pedicle finder
Palpeur pédiculaire
Palpador pedicular
003PF3

Use the pedicle finder to confirm five distinct bony borders: floor and four walls (medial, lateral, superior and inferior).

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
Connect the selected tap to the removable handle and tap the pedicle and the vertebral body.

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Screw placement

Screwdriver connection with screw

Engage the screwdriver into the saddle of the screw head.

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

Fully seat the inner shaft into the screw head.

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Screw insertion

Insert the screws into the vertebral bodies. Screws should be checked radiographically at this time to ensure intraosseous screw placement.

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
Disengage the screwdriver from the screw by counterclockwise rotation of the handle.

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Screw placement

Changing screwdriver shafts

In standard kit, screwdrivers are delivered assembled with a shaft for polyaxial screw.

Shafts specific for different types of screws may be delivered separately.

To assemble a new screwdriver:

- 1 Push the "Push only for cleaning" button off the knob.

- 2 Slide it through the shaft, choose the appropriate inner shaft and assemble it within outer tube and black sleeve. Then assemble the knob within screwdriver.

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

Screw adjustment

Specific screwdrivers for polyaxial screws or for monobloc screws may be used for correction, in the case where the screw position is not satisfactory and requires **unscrewing**.

Assemble the screwdriver using a quick connect handle. Insert the tip in the screw and unscrew it.

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
Rod holder
Porte tige
Porta-barra
018HR2

Screw head alignment

In order to facilitate the rod placement and locking screws insertion, rod holder can be used to align polyaxial screw heads and screwdriver monobloc for monobloc screws.

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Trial rod
Tige malléable
Barra maleable
250 mm 339D250TR1
500 mm 351D500TR1

Rod selection

Trial rods can be used to estimate the desired rod length and curvature.

Trial rods are malleable material (titanium shape memory). When heated (cleaning and sterilization process), the rod will recover its original straight shape.

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Use rod cutter with definitive rod only.

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Rod bender
Cintreuse de tige
Doblador de barra
035BR1

Rod contouring

To fit the desired spinal contours, rod bending is performed with the rod bender.

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Rod placement and locking screw insertion

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Rod placement and locking screw positioning


After having chosen the appropriate length with the desired lordosis (cf. implants references p.42), rod is placed into the screws with the rod holder.

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After picking-up locking screw with nut holder, engage it. **Do not tighten completely** to keep the rod free if further correction is necessary.

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
Open rod pusher
Pousse-tige ouvert
Empujador de barra abierto
028IOP1

Option 1 : Open rod pusher

When a minimal amount of reduction is required, the open rod pusher can be used to seat the rod into the head of the pedicle screw.

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
Closed rod pusher
Pousse-tige fermé
Empujador serrado de barra
029ICP1

Option 2 : Closed rod pusher

When a more important reduction is required, the closed rod pusher can be used.

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
Option 3 : Vertical rod pusher

Use the vertical rod pusher when additional force is needed to bring the rod to the implant. (Refer to next page for more details).

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Reduction procedures

Position of the laser marking indicates the position of the rod compared to the pedicle screw head:

1 Vertical rod pusher can be connected to the pedicle screw. To achieve the position, rotate the black handle counterclockwise.

2 At this position, the vertical rod pusher is locked to the pedicle screw and the rod can be pushed into the screw.

3 Then rotate the black handle clockwise until the indication line is in the closer handle position. At this position, the rod is fully seated and the locking nut can be inserted using the screwdriver shaft T30 retaining.

Do not tighten completely the locking screw to keep the rod free if further correction is necessary.

Screwdriver shaft, T30 retaining is mandatory to engage the locking screw with vertical rod pusher.

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


Vertical rod pusher removal

Dismantling

To remove the vertical rod pusher, rotate the black handle counterclockwise completely. Then pull up the grey part and twist the instrument to disengage from the screw.

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
Pedicle traction screws

mont blanc[®] pedicle traction screws can be used during a reduction procedure. The dedicated screwdriver shaft is used to insert traction screws into the pedicles.

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
Rod holder
Porte tige
Porta-barra
018HR2


Wing breaker

Remove the wings from the pedicle screw heads using the wing breaker (optional) or the rod holder.

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Deformity correction

The mont blanc[®] system was designed to offer solutions that accommodate various surgical philosophies.

Deformity correction may be obtained using one of three different reduction procedures :

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
Correction procedures

Traditionnal rod derotation

To start the rotational correction maneuver, rod need to be inserted into all of screws and locking screws engaged but not tightened.

The rod may be rotated using a rod clamp and the derotation wrench.

Once the rod has been fully rotated, locking screws are provisionally tightened.

Additional deformity correction may be obtained by further distraction/compression maneuvers.

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

Correction procedures

Distraction/Compression

Spinal deformities can be further effected by creating a distraction on the concavity of the deformity with the distraction clamp and compression on the convexity of the deformity with the compression clamp.

The compression or distraction maneuvers should be performed with all locking screws inserted but not finally tightened.

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
Use bender in-situ to bend rod in the sagittal plan and coronal bender to bend rod in the frontal plan.

Care should be taken to not wake extreme bends, so as to avoid stress concentration and notching of the rod.

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Iliac fixation


mont blanc[®] Iliac Fixation solution creates a strong anchorage point and stable rigid fixation. After having determined entry point of the iliac wing, create a cortical breach by using the sharp tip and its depth adjustable ring. As with pedicle screw procedure, follow same steps to insert the iliac screw. Screw should be checked at this time to ensure intraosseous screw placement.

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Iliac lateral connector and rod placement

First, pass the lateral connector through the rod, then seat the connector onto the head of the iliac screw. When all implants are securely in place, final tightening can be performed.

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
The torque limiter allows to transmit the optimal torque force for final tightening (10 Nm). Three "click" validate the optimal tightening. Use the counter torque lever, in way to avoid a torque effort on the vertebrae.

In case of using the counter torque lever XL, use the screwdriver shaft T30 retaining quick connect.

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
Fixation

Position the transverse connector onto the rods:

- 1- Unscrew the nut by a quarter of a turn
- 2- Apply pressure with the screwdriver until you feel a "CLICK".
- 3- Proceed in the same way on the opposite side.

At this point the connector remains adjustable.

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
Transverse connector setting

Positioning

The modularity of the transverse connector allows a positioning adapted to the patient's anatomy without screw positioning compromise.

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Final locking

Proceed with definitive tightening of nuts with the screwdriver for transverse connector.

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

Choice of transverse connector dimension

Select the size of the connector with the measuring caliper (optional).

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Hook Implantation site preparation

Pedicle elevator
Rugine pédiculaire
Legra pedicular
040EP2


Pedicle

Use the pedicle elevator to prepare the insertion of a pedicle hook.

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Transverse elevator
Rugine transverse
Legra transversal
038ET1

Transverse

Use the transverse elevator to prepare the insertion of a transverse hook.

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Laminar elevator
Rugine lamaire
Legra laminar
039EL1


Laminar

Use the laminar elevator to prepare the insertion of a laminar hook. This elevator is only used for lumbar vertebrae.

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
Hook holder
Pince porte-crochet
Pinza portaganchos
021HH1


Option 1: Hook holders

According to accessibility, use the hook holder or the offset hook holder to set hooks

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
Hook pusher
Pousse-crochet
Empujador de ganchos
030IHP2

Option 2: Hook holder + hook pusher

Use the hook pusher when you need more strength to set the hook. Especially to set pedicle hooks.

Hook pusher cannot be used with the offset hook holder.

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

Closed rod pusher
Pousse-tige fermé
Empujador serrado de barra
029ICP1

Rod pusher, vertical
Pousse-tige vertical
Empujador de barra, vertical
034IRP1

Rod reduction with hooks follows strictly the same process as rod reduction with screws.

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Using C-rings


After setting up mont blanc hooks, the installation of C-Ring may be necessary when a movement of distraction or derotation is needed.

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

C ring
Anneau en C
Anillo en C
080HCR1

Tighten the C-ring on the rod with the screwdriver.

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The C-ring offers a support zone for the instrument and allows to keep the hook in place. Distraction with hooks follows strictly the same process as distraction with screws.

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Remove the C-ring after using. (Although made out of implantable stainless steel 316L).

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
Hooks final tightening

«CLICS»

Final tightening with hooks follows strictly the same steps as final tightening with screws.

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Removal or revision

To remove or revise a mont blanc[®] previous surgery, a mont blanc[®] ablation set (C1STD02) is required.

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Transverse connector removal

If transverse connector removal is required, use the screwdriver for transverse connector (398SNS1)

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Locking screw removal

Slip the T-screwdriver (425ST1 or 009ST1 if Hex. 4.5mm) through the counter torque (074HAR2), then untighten the locking screws.

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Rod removal

If needed, loosen all the locking screws, holding the rod with the rod holder (018HR2).

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Pedicle screw removal

In the case where revising or extracting the material is necessary, use the type of screwdriver that corresponds to the screw (259SP3D1 for polyaxial screws or 258SM3D1 for monobloc screws)

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Hook removal

In the case where hook removal is required, remove the locking screws (see Locking screw removal), then use the hook holder (021HH1) to remove the hook.

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Implants references

All the implants are in titanium alloy TA6V, except cobalt-chromium rods.

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All measurements are in millimeter

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Pedicle polyaxial screw

Vis pédiculaire polyaxiale
Tornillo pedicular poliaxial

L	Ø4	Ø5	Ø6	Ø7	Ø8
25	MB1PS425	-	-	-	-
30	MB1PS430	MB1PS530	MB1PS630	MB1PS730	MB1PS830
35	MB1PS435	MB1PS535	MB1PS635	MB1PS735	MB1PS835
40	MB1PS440	MB1PS540	MB1PS640	MB1PS740	MB1PS840
45	-	MB1PS545	MB1PS645	MB1PS745	MB1PS845
50	-	MB1PS550	MB1PS650	MB1PS750	MB1PS850
55	-	-	MB1PS655	-	-

Pedicle polyaxial traction screw

Vis pédiculaire de traction polyaxiale
Tornillo pedicular poliaxial de traccion

L	Ø5	Ø6	Ø7
30	-	MB1PSS630	-
35	-	MB1PSS635	MB1PSS735
40	MB1PSS540	MB1PSS640	MB1PSS740
45	MB1PSS545	MB1PSS645	MB1PSS745
50	MB1PSS550	MB1PSS650	MB1PSS750
55	-	MB1PSS655	-

Iliac polyaxial screw

Vis polyaxiale iliaque
Tornillo poliaxial iliaco

L	Ø7	Ø8
55	MB1PS755	MB1PS855
60	MB1PS760	MB1PS860
65	MB1PS765	MB1PS865
70	MB1PS770	MB1PS870
75	MB1PS775	MB1PS875
80	MB1PS780	MB1PS880
85	MB1PS785	MB1PS885
90	MB1PS790	MB1PS890
95	MB1PS795	MB1PS895
100	MB1PS7100	MB1PS8100

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Removable handle,
quick connect, Ratchet
Manche amovible
encliquetable, Ratchet
Mango amovible,
conectable, Ratchet
313HRMR1

Counter torque lever, XL
Lever anti-rotation, XL
Palanca anti rotación, XL
424HARL1

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FDA cleared

Spinal implantable medical devices class IIb / class I Instruments, manufacturer Spineway. For a good use of these medical devices, please read the IFU and the surgical technique. DM registered on the LPPR list referred to Article L.165-1 of the French Code of Social Security : See terms on ameli.fr website. First EC marking: 2006.

SPINEWAY, mont blanc are registered trademarks of Spineway company.

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Instructions for use

INDICATIONS FOR USE:

The MONT BLANC and MONT BLANC MIS Spinal Systems are intended to provide immobilization and stabilization of spinal segments in skeletally mature patients as an adjunct to fusion in the treatment of acute and chronic instabilities or deformities of the thoracic, lumbar and sacral spine and sacral/iliac screw fixation.

The MONT BLANC and MONT BLANC MIS Spinal Systems are intended for noncervical pedicle fixation and nonpedicle fixation for the following indications: degenerative disc disease (defined as back pain of discogenic origin with degeneration of the disc confirmed by history and radiographic studies); spondylolisthesis; trauma (i.e., fracture or dislocation); spinal stenosis; curvatures (i.e., scoliosis, kyphosis, and/or lordosis); tumor, pseudoarthrosis; and failed previous fusion in skeletally mature patients.

When used in a posterior percutaneous approach with MIS instrumentation, the MONT BLANC MIS System is intended for noncervical pedicle fixation and nonpedicle fixation for the following indications: degenerative disc disease (defined as back pain of discogenic origin with degeneration of the disc confirmed by history and radiographic studies); spondylolisthesis; trauma (i.e., fracture or dislocation); spinal stenosis; curvatures (i.e., scoliosis, kyphosis, and/or lordosis); tumor, pseudoarthrosis; and failed previous fusion in skeletally mature patients.

When used for posterior non-cervical pedicle screw fixation in pediatric patients, the MONT BLANC and MONT BLANC MIS Spinal Systems metallic implants are indicated as an adjunct to fusion to treat adolescent idiopathic scoliosis. The MONT BLANC and MONT BLANC MIS Spinal Systems are intended to be used with autograft and/or allograft. Pediatric pedicle screw fixation is limited to a posterior approach.

CONTRAINDICATIONS:

The Mont Blanc Spinal Systems are contraindicated for:

- Any medical or surgical condition precluding the potential benefit of spinal surgery
- Any medical or mental condition which would exclude the patient at high risk from surgery of this severity;
- Acute or chronic systemic, spinal or localized infections
- Systemic and metabolic diseases
- Obesity
- Pregnancy
- Dependency on pharmaceutical drugs, drug abuse, or alcoholism
- Lack of patient cooperation
- Foreign body sensitivity to the implant material
- Significant osteopenia
- Soft tissue deficit not allowing wound closure;
- Congenital abnormalities, tumors or other conditions that would prevent secure component fixation that has the potential to decrease the useful life of the device;
- Inadequate pedicles of the thoracic, lumbar, and sacral vertebrae.
- Pars defect
- Clinically compromised vertebral bodies at affected level due to current or past trauma

Relative contraindications include disease conditions that have been shown to be safely and predictably managed without the use of internal fixation devices. Severe osteoporosis is a relative contraindication because it may prevent adequate fixation of spinal anchors and thus preclude the use of this or any other spinal instrumentation system. Any entity or condition that totally precludes the possibility of fusion, i.e., cancer, kidney dialysis, or osteopenia is a relative contraindication. Other relative contraindications include obesity, certain degenerative diseases, and foreign body sensitivity. In addition, the patient's occupation, activity level, or mental capacity may be relative contraindications to this surgery. Specifically, patients who may place undue stresses on the implant during bony healing and may be at higher risk for implant failure because of their occupation or lifestyle, or because of conditions such as mental illness, alcoholism, or drug abuse.

WARNINGS:

Warnings, Precautions and Possible Adverse Effects concerning Temporary Metallic Internal Fixation Devices
These warnings, precautions and possible adverse effects

do not include all possible adverse surgical effects, but are particular to metallic internal fixation devices. Explain general surgical risks to the patient before surgery.

Potential risks identified with the use of this device system, which may require additional surgery, include: device component fracture, loss of fixation, nonunion, vertebral fracture, neurological injury, and vascular or visceral injury.

1. **CORRECT IMPLANT SELECTION IS VITAL.** Selecting the proper implant size, shape, and design increases the potential for satisfactory fixation. While proper selection can help minimize risks, the size and shape of human bones present implant size, shape, and strength limitations. Metallic internal fixation devices cannot withstand activity levels equal to those placed on normal healthy bone. No implant can be expected to withstand indefinitely the unsupported stress of full weight bearing.
2. **IMPLANTS CAN BREAK WHEN SUBJECTED TO THE INCREASED LOADING ASSOCIATED WITH DELAYED UNION OR NONUNION.** Internal fixation appliances are load-sharing devices that are used to obtain alignment until normal healing occurs. If healing is delayed or does not occur, the implant may eventually break due to metal fatigue. The degree or success of union, loads produced by weight bearing, and activity levels among other conditions will dictate implant longevity. Notches, scratches or implant bending during the surgery may also contribute to early failure. Fully inform patients of the implant failure risks.
3. **MIXING METALS CAN CAUSE CORROSION.** There are many forms of corrosion damage, and several of these occur on metals surgically implanted in humans. General or uniform corrosion is present on all implanted metals and alloys. The rate of corrosive attack on metal implant devices is usually very low due to the presence of passive surface films. Dissimilar metals in contact, such as titanium and stainless steel, accelerate the corrosion process of stainless steel, and more rapid attack occurs. The presence of corrosion often accelerates fatigue fracture of implants. The amount of metal compounds released into the body system will also increase. Internal fixation devices, such as rods, hooks, wires, etc., that come into contact with other metal objects, must be made from like or compatible materials.
4. **PATIENT SELECTION:** The following factors can be extremely important to the eventual success of the procedure:
 - A. The patient's weight. An overweight or obese patient can produce loads on the device that can lead to failure of the appliance and the operation.
 - B. The patient's occupation or activity. If the patient is involved in an occupation or activity that includes heavy lifting, muscle strain, twisting, repetitive bending, stooping, running, substantial walking, or manual labor, he/she should not return to these activities until the bone is fully healed. Even with full healing, the patient may not be able to return to these activities successfully.
 - C. Senility, mental illness, alcoholism, or drug abuse. These conditions, among others, may cause the patient to ignore certain necessary limitations and precautions in the device use, leading to implant failure or other complications.
 - D. Certain degenerative diseases. In some cases, degenerative disease progression may be so advanced at implantation that it may substantially decrease the device's expected useful life. For such cases, orthopedic devices can only be considered a delaying technique or temporary remedy.
 - E. Foreign body sensitivity. No pre-operative test can completely exclude the possibility of sensitivity or allergic reaction. Patients can develop sensitivity or allergy after implants have been in the body for a period of time.
 - F. Smoking. Patients who smoke have been observed to experience higher rates of pseudoarthrosis following surgical procedures where bone graft is used. Additionally, smoking has been shown to cause diffuse degeneration of intervertebral discs. Progressive degeneration of adjacent segments caused by smoking can lead to late clinical failure (recurring pain) even after successful fusion and initial clinical improvement.

The safety and effectiveness of these devices has not been established for use as part of a pediatric growing rod construct. These devices are only intended to be used when definitive

fusion is being performed at all instrumented levels.

The use of pedicle screw fixation in the pediatric population may present additional risks when patients are of smaller stature and skeletally immature. Pediatric patients may have smaller spinal structures (pedicle diameter or length) that may preclude the use of pedicle screws or increase the risk of pedicle screw malpositioning and neurological or vascular injury. Patients not skeletally mature that undergo spinal fusion procedures may have longitudinal spinal growth, or may be at risk for rotational spinal deformities (the "crankshaft phenomenon") due to continued differential growth of the anterior spine. Other adverse events related to pedicle screw fixation, such as screw or rod bending, breakage, or loosening, may also occur in pediatric patients. Pediatric patients may be at increased risk for device related injury because of their smaller stature.

PRECAUTIONS:

1. **SURGICAL IMPLANTS MUST NEVER BE REUSED.** Even though the device appears undamaged, it may have small defects and internal stress patterns which may lead to early breakage.
2. **CORRECT IMPLANT HANDLING IS VITAL.** Only contour metal implants with proper equipment. Avoid any notching. Scratching or reverse bending of the devices when contouring. Alterations will produce defects in surface finish and internal stresses that may become the focal point for eventual breakage. Do not use the implant if damage is suspected. Bending of screws will significantly decrease the fatigue life and may cause failure. Excessive torque applied to the screws when seating the rod may cause bone failure resulting in stripped threads and/or compromised screw purchase.
3. **BENDING THE CONSTRUCT.** Titanium alloy components should never be bent sharply or reverse bent. If a construct is over-contoured contour a new construct correctly rather than reverse bending the over-contoured construct.
4. **IMPLANT REMOVAL AFTER HEALING.** If the device is not removed after the completion of its intended use, any of the following complications may occur:
 - (1) Corrosion, with localized tissue reaction or pain;
 - (2) Implant migration resulting in injury;
 - (3) Risk of additional injury from postoperative trauma;
 - (4) Bending, loosening, and/or breakage, which could make removal impractical or difficult;
 - (5) Pain, discomfort, or abnormal sensations due to device presence;
 - (6) Possible increased risk of infection;
 - (7) Bone loss due to stress shielding. Carefully weigh the risks versus benefits when deciding whether to remove the implant. Implant removal should be followed by adequate postoperative management to avoid refracture or deformity. If the patient is older and has a low activity level, the surgeon may choose not to remove implant thus eliminating the risks involved in a second surgery.
5. **ADEQUATELY INSTRUCT THE PATIENT.** Postoperative care and the patient's ability and willingness to follow instructions are among the most important aspects of successful bone healing. Inform the patient about the implant limitations, and to limit physical activities, especially lifting and twisting motions and participating in any type of sports. Tell the patient that a metallic implant is not as strong as normal healthy bone and could loosen, bend and/or break if excessive demands are placed on it, especially in the absence of complete bone healing. Implants displaced or damaged by improper activities may migrate and damage the nerves or blood vessels. An active, debilitated, or demented patient who cannot properly use weight-supporting devices may be particularly at risk during postoperative rehabilitation.

The implantation of pedicle screw spinal systems in pediatric patients should be performed only by experienced spinal surgeons with specific training in the use of this pedicle screw spinal system in pediatric patients because this is a technically demanding procedure presenting a risk of serious injury to the patient.

Preoperative and operating procedures, including knowledge of surgical techniques, good reduction, and proper selection and placement of the implants are important considerations in the successful utilization of the system in pediatric patients. **PHYSICIAN NOTE:** Although the physician is the learned intermediary between the company and the patient, the important medical information given in this document should be conveyed to the patient.

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