

OPERATOR'S MANUAL

Defibrillator/Monitor D100

EU representative

EC REP OBELIS S.A

Bd. Général Wahis, 53, 1030 Brussels, Belgium

Manufacturer

 **Mediana Co., Ltd.**

132, Donghwagongdan-ro, Munmak-eup,
Wonju-si, Gangwon-do, Korea

Tel: (82) 2 542 3375 (82) 33 742 5400

Fax: (82) 2 542 7447 (82) 33 742 5483

D100 Operator's Manual

Part Number: A7741-2

Revised Date: 2022-10



Copyright © 2022 All rights reserved.

Directive

- Copyright law allows no part of this instruction manual to be reproduced without permission.
- The content of this manual are subject to change without notice.
- The contents of this manual should be correct. If for some reason, there are any questionable points, please do not hesitate to contact our service center.
- The manual will be replaced if any pages are missing or collation is incorrect.

Warranty

- Please contact your local distributor about the warranty period.
- Device failure or damage related to the following situations during the guarantee period is not covered by this warranty:
 - Installation, transfer installation, maintenance and repairs by any person other than an authorized employee or technician by Mediana.
 - Damage sustained to the Mediana product(s) caused by product(s) from another company excluding products delivered by Mediana.
 - Damage – caused by mishandling and/or misuse – is the responsibility of the user.
 - Maintenance and repairs utilizing maintenance components that are not specified by Mediana.
 - Device modifications or use of accessories not recommended by Mediana.
 - Damage caused by accidents or natural disasters (earthquakes, flooding, etc.).
 - Damage resulting from usage where caution statements and operating instructions shown in this manual have not been followed.
 - Damage due to neglect of specified maintenance checks.
- This warranty only covers the hardware of the D100. The warranty does not cover the following selections:
 - Whatever damage or loss results from the attachment of accessories or their operation.
 - In the event of a defect in the product, contact our sales outlet or EU representative as noted on the back cover.
- The D100 conforms to the EMC standard IEC60601-1-2.

Note that mobile phones should not be used in the vicinity of the D100.

Note, however, any device not complying to the EMC standard that is used with the D100 renders the D100 as non-compliant to the EMC standard.

Trademark

Product brand names shown in this manual are likely to be the trademark or registered trademark of the company concerned.

CONTENTS

CONTENTS.....	iii
SAFETY INFORMATION	9
General Safety Information	9
Warning	9
Cautions	14
INTRODUCTION.....	15
Intended Use for the D100	15
Indications for Use.....	15
About This Manual.....	16
Identifying the D100 Configurations	17
DESCRIPTION OF THE DEFIBRILLATOR/MONITOR.....	19
Front Panel Components.....	19
Top Panel Components	21
Rear Panel Components	22
Left Panel Components	23
Right Panel Components.....	24
SETTING UP THE DEFIBRILLATOR/MONITOR.....	27
Unpacking and Inspection	27
List of Components.....	28
Power Cable Connections	29
Battery Operation	32
USING THE DEFIBRILLATOR/MONITOR	37
Turning On and Off the Defibrillator/Monitor	38
Displays	40
Operating the Defibrillator/Monitor	42
Software Menu Area	42
Configure	44
ALARMS AND LIMITS.....	51
General	51
Alarm Priority and Messages	51
Visual Alarm Indication	53
Audible Alarm Indication	54
Changing Alarm Limits.....	55
Alarm Audio Acknowledged (Silence).....	56
Audio Paused and OFF	57
AED(Automated External Defibrillator) MODE	59
AHA/ERC Guidelines (Rescue protocol)	61
Preparing for Defibrillation	61
Performing CPR	62
Description of AED Mode Menu Functions	64
MANUAL MODE	65
General	65
Preparing for Defibrillation	66
Operating the Manual Mode of Defibrillator/Monitor	67
Defibrillating (Async mode)	68
Synchronized Cardioversion (Sync mode).....	69
Description of the Manual Mode Menu Functions.....	70
MONITOR MODE.....	71
General	71
Description of Monitor Mode Menu Functions.....	71
Freeze	72
PACING MODE	73
General	74

Description of Pacing Mode Menu Functions	74
Demand Mode and Fixed Mode	75
Preparing for Pacing	75
Operating the Pacing Mode of Defibrillator/Monitor	75
ECG MONITORING	77
General	78
Setup Connections	78
Description of ECG Waveform Menu Functions	81
SpO₂ MONITORING	83
General	84
Setup Connections	85
Description of SpO ₂ Menu Functions	86
Description of SpO ₂ Waveform Menu Functions	87
Theory of Operation	88
Clinical Studies	90
RESPIRATION MONITORING	93
General	93
Setup Connections	93
Description of Respiration Menu Functions	94
Description of Respiration Waveform Menu Functions	95
Theory of Operation	96
SELF-TEST FUNCTION	97
General	97
External Shock Test	98
Trouble Shooting	98
EVENT REPORT	99
General	99
TRENDS	101
General	101
DATA MANAGEMENT	103
General	103
Data Management Mode	104
PRINTING	107
General	107
Type of Print-Out	109
Replace printer paper	112
EXTERNAL COMMUNICATIONS	113
General	113
SD Memory Card	113
Central System Communication	113
MENU STRUCTURE	115
MAINTENANCE	121
Recycling and Disposal	121
Returning the Defibrillator/Monitor and System Components	122
Service	122
Periodic Safety Checks	122
Cleaning	122
Battery Maintenance	123
TROUBLESHOOTING	125
General	125
Obtaining Technical Assistance	125
EMI (Electromagnetic Interference)	126
FACTORY DEFAULTS	129
General	129
Parameter Ranges and Default Settings	129
SPECIFICATION	133

Display	133
Controls	133
Alarms	133
Physical Characteristics and Printer	133
Electrical	134
Environmental Conditions	135
Tone Definition	135
Measurement Parameters	136
Event/Trend	140
Defibrillator (Technical Specification)	141
Compliance	155
Manufacturer's EMC Declaration	157

Figures

Figure 1. Front panel components	19
Figure 2. Top panel components	21
Figure 3. Defibrillator/Patient monitor paddle controller	21
Figure 4. Rear panel components	22
Figure 5. Left panel components	23
Figure 6. Right panel components	24
Figure 7. Right panel components (Option)	25
Figure 8. AC power connection	29
Figure 9. Cradle installation	30
Figure 10. Battery placement	33
Figure 11. Mode select knob	39
Figure 12. Displays (Manual mode)	40
Figure 13. Software menu area (Manual mode)	42
Figure 14. Configure window	44
Figure 15. 4mV waveform	47
Figure 16. AED mode menu	64
Figure 17. Manual mode screen	70
Figure 18. Monitor mode menu – Full color mode	71
Figure 19. Monitor mode menu – Grayscale mode	72
Figure 20. Monitor mode menu – Freeze	72
Figure 21. Pacing mode menu	74
Figure 22. Pacing display (Fixed mode)	75
Figure 23. Standard 3 electrode placement	78
Figure 24. 5 electrode placement	79
Figure 25. HR/PR display	80
Figure 26. ECG waveform display	81
Figure 27. SpO ₂ display	86
Figure 28. SpO ₂ waveform display	87
Figure 29. Oxyhemoglobin dissociation curve	89
Figure 30. Modified bland-altman plot	92
Figure 31. Respiration display	94
Figure 32. Respiration waveform display	95
Figure 33. Event report	99
Figure 34. Event review	100
Figure 35. Trend – Tabular trend	101
Figure 36. Data management	104
Figure 37. Data management – Export	105
Figure 38. Data management – Review	106
Figure 39. Measurement data printing	109
Figure 40. Trend data printing	109
Figure 41. Event data printing	109
Figure 42. Event summary review data printing	110
Figure 43. Setting value printing	111
Figure 44. Self-test result printing	111
Figure 45. Replace print paper	112
Figure 46. Biphasic waveforms at 200 joules	145
Figure 47. Biphasic waveforms at 175 joules	145
Figure 48. Biphasic waveforms at 150 joules	146
Figure 49. Biphasic waveforms at 125 joules	146
Figure 50. Biphasic waveforms at 100 joules	147
Figure 51. Biphasic waveforms at 75 joules	147
Figure 52. Biphasic waveforms at 50 joules	148
Figure 53. Biphasic waveforms at 40 joules	148
Figure 54. Biphasic waveforms at 30 joules	149
Figure 55. Biphasic waveforms at 20 joules	149

Figure 56. Biphasic waveforms at 10 joules.....	150
Figure 57. Biphasic waveforms at 9 joules.....	150
Figure 58. Biphasic waveforms at 8 joules.....	151
Figure 59. Biphasic waveforms at 7 joules.....	151
Figure 60. Biphasic waveforms at 6 joules.....	152
Figure 61. Biphasic waveforms at 5 joules.....	152
Figure 62. Biphasic waveforms at 4 joules.....	153
Figure 63. Biphasic waveforms at 3 joules.....	153
Figure 64. Biphasic waveforms at 2 joules.....	154
Figure 65. Biphasic waveforms at 1 joule	154

Tables

Table 1. D100 controls.....	20
Table 2. Panel and label symbols.....	26
Table 3. Standard accessories	28
Table 4. Optional accessories	28
Table 5. Front panel indications for power source.....	34
Table 6. The defibrillator/monitor battery status icon.....	35
Table 7. Display Symbols	41
Table 8. Display menu	45
Table 9. Sound menu	45
Table 10. Waveform menu	46
Table 11. The number of ECG 4mV waveform (Full option).....	46
Table 12. Parameter menu	47
Table 13. Alarm setup menu.....	48
Table 14. Patient menu.....	48
Table 15. Printer menu	48
Table 16. Review menu	49
Table 17. Service setting menu	49
Table 18. Others menu	49
Table 19. Alarm priority condition	51
Table 20. Visual alarm characteristics.....	53
Table 21. Audible alarm characteristics.....	54
Table 22. Alarm setup menu.....	55
Table 23. Alarm limits ranges	56
Table 24. ECG lead colors.....	79
Table 25. ECG lead pairs	79
Table 26. HR/PR menu.....	80
Table 27. ECG waveform menu	81
Table 28. SpO ₂ sensors.....	85
Table 29. SpO ₂ menu	86
Table 30. SpO ₂ waveform menu.....	87
Table 31. Demographic data	90
Table 32. Demographic data	91
Table 33. SpO ₂ accuracy for Medtronic™ sensors	92
Table 34. Respiration menu.....	94
Table 35. Respiration waveform menu	95
Table 36. Printer Menu	107
Table 37. Printing waveform menu.....	108
Table 38. Parameter ranges and factory defaults	129
Table 39. ECG analysis algorithm accuracy: Adult	143
Table 40. ECG analysis algorithm accuracy: Pediatric	143
Table 41. Delivered energy at every defibrillator settings into a range of loads.....	144
Table 42. Electromagnetic emissions (IEC60601-1-2).....	157
Table 43. Electromagnetic immunity (IEC60601-1-2)	157
Table 44. Electromagnetic immunity (IEC60601-2-4)	159
Table 45. Electro-surgical unit interference (IEC60601-2-27, IEC-60601-2-30)	160
Table 46. Recommended separation distances	160
Table 47. Immunity to proximity fields from RF wireless communications equipment (IEC60601-1-2)	161
Table 48. Cables (IEC60601-1-2).....	161

SAFETY INFORMATION

General Safety Information

This section contains important safety information related to general use of the D100 defibrillator/monitor. Other important safety information appears throughout the manual. The D100 defibrillator/monitor will be referred to as the defibrillator/monitor throughout this manual.

Important! Before use, read carefully the manual, accessory directions for use, all precautionary information and specifications.

Warning

Warnings alert you to potential serious outcomes (death, injury, or adverse events) to the patient or user.

 WARNING	When using the defibrillator/monitor with anesthetics, nitrous oxide or high concentrations of oxygen, connect the gas outlet to a scavenger system.
 WARNING	When using the defibrillator/monitor with a commercial electric power source, use the defibrillator/monitor with an electric power wall socket with a grounding wire for medical use. Not doing so could cause electric shock.
 WARNING	Do not connect grounding wire to gas pipes. This could cause fire.
 WARNING	Only doctors and officially certified personnel should use this defibrillator/monitor. Do not allow patients to touch the defibrillator/monitor. Allowing patients to touch the defibrillator/monitor could cause accidents.
 WARNING	The defibrillator/monitor cannot be used when MRI is in progress. If MRI is in use, keep patient attachments away from patients to prevent accidents.
 WARNING	The defibrillator/monitor conforms to the requirements of the EMC standard (IEC60601-1-2), and may therefore be used simultaneously with pacemakers and other electrical simulators. It should, however, be noted that the defibrillator/monitor may be affected by electrical scalpels and microwave therapeutic apparatus. Please check operation of the defibrillator/monitor during and after use of such equipment.
 WARNING	Do not take mobile phones or transceivers into a room where this defibrillator/monitor is installed, as such devices may cause accidents.
 WARNING	In order to avoid accidents, do not use any unauthorized accessories or options.
 WARNING	Thoroughly read the instruction manuals supplied with accessories and options to ensure correct use. This instruction manual does not carry the caution selections for such equipment.
 WARNING	Do not open cover or disassemble this defibrillator/monitor. Doing so could cause electric shock or fire. It is prohibited by law to modify the defibrillator/monitor without authorization.
 WARNING	Do not use power source other than the specified voltage, (100-240V~50/60Hz) as this may cause fire or electric shock.

 WARNING	Pre-use inspection and preventive maintenance must be performed for safe use.
 WARNING	The defibrillator/monitor may be used with electrical surgical equipment.
 WARNING	Follow the instruction manuals for medical instruments – notably electrosurgical and diathermy instruments – when used, as their high-frequency energy units may cause burns to patients via attachments.
 WARNING	This defibrillator/monitor is protected against the discharge of a defibrillator. However, do not touch the defibrillator/monitor when a defibrillator is being discharged (electrified), as doing so may cause electric shock.
 WARNING	The following cautions apply when connecting the defibrillator/monitor with other equipment. 1. Ensure that the connected equipment is in accordance with the IEC60601-1 or IEC safety standards, so that the system complies with IEC60601-1. 2. Employ additional protective measures (e.g. additional protective earthing) as necessary.
 WARNING	Do not connect devices that do not meet medical safety standards (such as commercial PCs), as they may cause electric shock. This defibrillator/monitor meets the restricted level of leakage current required for medical devices. Therefore, this defibrillator/monitor cannot be connected to a device that would give a combined total of leakage current beyond the restricted level.
 WARNING	Avoid connecting the patient to several devices at once. Leakage current limits may be exceeded. Do not use a second defibrillator on the patient while pacing with the defibrillator/monitor.
 WARNING	Do not place anything on top of this defibrillator/monitor. If something is spilled over the defibrillator/monitor or gets into it, such spillage may cause fire or electric shock. If fluid spills on the defibrillator/monitor accidentally, disconnect power cord, wipe dry immediately, and have the defibrillator/monitor serviced to make sure that no hazard exists.
 WARNING	Do not place heavy objects on the power cord, as doing so may cause fire or electric shock.
 WARNING	Before conducting maintenance work, turn the power off and unplug the power cord from the wall socket to prevent electric shock.
 WARNING	When the following occur, turn the power off immediately and unplug the power cord from the wall socket. Continued use in such situations may cause fire or electric shock. • There is smoke or a strange odor leaking out of the defibrillator/monitor. • The defibrillator/monitor has been dropped or impacted by an object. • Liquid or foreign matter gets inside the defibrillator/monitor. • Defibrillator/monitor failure has occurred. Also, when any of the above occurs, promptly do the following: 1. Check to see that the power cord has been unplugged from the wall socket. 2. Place an 'Out of Order' sign on the defibrillator/monitor and do not use it. 3. Have the defibrillator/monitor serviced to make sure that no hazard exists.

 WARNING	Do not connect more than one patient to the defibrillator/monitor. Do not connect more than one defibrillator/monitor to a patient.
 WARNING	The defibrillator/monitor is a prescription device and is to be operated by qualified personnel only.
 WARNING	As with any medical equipment, carefully route patient cabling to reduce the possibility of patient entanglement or strangulation.
 WARNING	Never lift the defibrillator/monitor by the sensor cable, power cord, or any other accessory. Such accessories could detach, causing the defibrillator/monitor to fall on the patient.
 WARNING	Do not make any clinical judgments based on this defibrillator/monitor's measurement only.
 WARNING	Emergency defibrillation should be performed only by appropriately trained, skilled and qualified personnel who are aware of the protocol for handling a patient in medical emergency such as cardiac arrest and have been certified in Advanced Cardiac Life Support (ACLS) or Basic Life Support (BLS).
 WARNING	Synchronized electrical cardioversion should be performed only by skilled personnel trained in Advanced Cardiac Life Support (ACLS) and practiced in equipment operation. The precise cardiac arrhythmia must be determined prior to performing defibrillation.
 WARNING	The defibrillator delivers up to 200 joules of electrical energy. Unless properly used as described in this manual, this electrical energy may cause serious injury or death. Do not attempt to operate this defibrillator/monitor unless thoroughly familiar with this manual and the function of all controls, indicators, connectors and accessories.
 WARNING	Do not discharge standard paddles on top of pads or ECG electrodes. Do not allow standard paddles (or pads) to touch each other, ECG electrodes, lead wires, dressings, transdermal patches, etc. Such contact can cause electrical arcing and patient skin burns during defibrillation and may divert defibrillating energy away from the heart muscle.
 WARNING	Discharging the defibrillator/monitor with the standard paddle surfaces shorted together can pit or damage the paddle electrode surface. Pitted or damaged paddle surfaces may cause patient skin burns during defibrillation. Discharge the defibrillator only as described in these operating instructions.
 WARNING	If a person is touching the patient, bed, or any conductive material in contact with the patient during defibrillation, the delivered energy may be partially discharged through that person. Clear everything away from contact with the patient, bed, and other conductive material before discharging the defibrillator/monitor.
 WARNING	Do not discharge the defibrillator into the open air. To remove an unwanted charge, change the energy selection, or change the mode, or turn off the defibrillator/monitor.
 WARNING	Conductive gel on the paddle handles can allow the electrical energy to be discharged through the operator during defibrillation. Completely clean the paddle electrode surfaces, handles, and storage area after defibrillation.
 WARNING	A gel pathway on the skin between the standard paddles will cause defibrillating energy to arc between paddles and delivery energy away from the heart muscle. Do not allow conductive gel (wet or dry) to

	become continuous between paddle sites.
 WARNING	During defibrillation checks, the discharged energy passes through the cable connectors. Securely attach cable connectors to the simulator.
 WARNING	Do not touch the patient or any equipment connected to the patient during defibrillation. Warn all persons around patient to DO NOT TOUCH THE PATIENT prior to defibrillation.
 WARNING	The defibrillator/monitor should be out of contact with water (puddles or water spray). It may cause electrical shocks and defibrillator/monitor failure. Electrical safety of the defibrillator/monitor may not work properly when wet.
 WARNING	Defibrillation may cause implanted devices to malfunction. Place standard paddles or pads away from implanted devices if possible. Check implanted device function after defibrillation.
 WARNING	Implanted pacemakers may cause the heart rate meter to count the pacemaker rate during incidents of cardiac arrest or other arrhythmias. Pacemaker patients should be carefully observed. Check the patient's pulse; do not rely solely on heart rate meters. Dedicated pacemaker detection circuitry may not detect all implanted pacemaker spikes. Patient history and physical exam are important in determining the presence of an implanted pacemaker.
 WARNING	To ensure patient electrical isolation, connect only to other equipment with circuits that are electrically isolated.
 WARNING	To avoid risk of electric shock, the defibrillator/monitor must only be connected to a supply mains with protective earth.
 WARNING	The use of ACCESSORY equipment not complying with the equivalent safety requirements of this equipment may lead to a reduced level of safety of the resulting system. Consideration relating to the choice shall include: Use of the accessory in the PATIENT VICINITY . Evidence that the safety certification of the ACCESSORY has been performed in accordance with the appropriate IEC (EN) 60601-1 and/or IEC (EN) 60601-2-XX particular standards.
 WARNING	Check leakage levels prior to use. Leakage current may be excessive if more than one defibrillator/monitor or other piece of equipment is connected to the patient.
 WARNING	Do not place or use the defibrillator/monitor when flammable gas presents in atmospheres or other flammable material exist near/around. Eg. Oxygenrich area, flammable anesthetic gases, gasoline and any flammable material or gas. Do not use the defibrillator/monitor near the place of a gasoline or other volatile substances spill may cause an explosion.
 WARNING	Do not disassemble the defibrillator. It contains no operator serviceable components and dangerous high voltages may be present. Contact authorized service personnel for repair.
 WARNING	Do not immerse any portion of this defibrillator in water or other fluids. Avoid spilling any fluids on defibrillator or accessories. Spilled liquids may cause the defibrillator and accessories to perform inaccurately or fail. Do not clean with ketones or other flammable agents. Do not autoclave or sterilize this defibrillator or accessories unless otherwise specified.
 WARNING	Use care when operating this device close to oxygen sources (such as bag-valve-mask devices or ventilator tubing). Turn off gas source or

	move source away from patient during defibrillation.
 WARNING	If you are monitoring a patient and using the system connector, all equipment connected to the system connector must be battery powered or electrically isolated from AC power according to EN 60601-1. If in doubt, disconnect the patient from the defibrillator before using the system connector. Only use Mediana recommended data transmission cables. For more information, contact Mediana Technical Support.
 WARNING	The defibrillator/monitor delivers up to 200 joules of electrical energy. When discharging the defibrillator, do not touch the paddle electrode surfaces or disposable paddle/pads electrodes. Do not attempt to perform this test unless you are qualified by training and experience and are thoroughly familiar with these operating instructions.
 WARNING	Electric shock hazards exist internally. Do not remove assembly screws. Refer servicing to qualified personnel.
 WARNING	Monitors, defibrillators, and their accessories (including electrodes and cables) contain ferromagnetic materials. As with all ferromagnetic equipment, these products must not be used in the presence of the high magnetic field created by a Magnetic Resonance Imaging (MRI) device. The high magnetic field created by an MRI device will attract the equipment with a force sufficient to cause death or serious personal injury to persons between the equipment and the MRI device. This magnetic attraction may also damage the equipment. Skin burns will also occur due to heating of electrically conductive materials, such as patient leads and pulse oximeter sensors. Consult the MRI manufacturer for more information.
 WARNING	Medical electrical equipment which does not incorporate defibrillator protection should be disconnected during defibrillation.
 WARNING	Operating the defibrillator/monitor or its accessories in conditions outside the environmental specifications can result in device or accessory malfunction. The defibrillator/monitor should be allowed to stabilize within the operating temperature range prior to operation.
 WARNING	Do not position the defibrillator/monitor so it is difficult to operate the disconnection device when a separable plug is used as isolation.

Cautions

Caution statements identify conditions or practices that could result in damage to the equipment or other property.

 CAUTION	Always check that the defibrillator/monitor functions properly and is in proper condition before use.
 CAUTION	Federal law restricts this device to sale by or on the order of a physician.
 CAUTION	The defibrillator/monitor may not operate properly if it is operated or stored at conditions outside the ranges stated in this manual, or subjected to excessive shock or dropping.
 CAUTION	When connecting the defibrillator/monitor to any instrument, verify proper operation before clinical use. Both the defibrillator/monitor and the instrument connected to it must be connected to a grounded outlet.
 CAUTION	Accessory equipment connected to the defibrillator/monitor's data interface must be certified according to IEC60950 for data-processing equipment or IEC60601-1 for electromedical equipment. All combinations of equipment must be in compliance with IEC60601-1-1 system requirements. Anyone who connects additional equipment to the signal input or signal output port configures a medical system and is therefore responsible that the system complies with the requirements of IEC 60601-1-1 and the electromagnetic compatibility system standard IEC60601-1-2. If in doubt, consult Mediana Technical Support Representative.
 CAUTION	Where the integrity of the external protective conductor in the installation or its arrangement is in doubt, equipment shall be operated from its internal electrical power source.
 CAUTION	Risk of explosion if battery is replaced by an incorrect type.
 CAUTION	Software alone does not encrypt data and should be used only in closed hospital networks, and use a PC equipped with a security system such as a firewall or vaccine that conforms to the hospital's policy.
 CAUTION	In the event of an accident involving cybersecurity, contact qualified service personnel or your local supplier for assistance and wait until the administrator takes action.

INTRODUCTION



Patient conditions may result in erroneous readings. If the measurements are suspect, verify the reading using another clinically accepted measurement method.

Intended Use for the D100

1 The defibrillator/monitor is intended for use by trained medical technician, doctor, nurse or medical specialist in outdoor and indoor emergency care settings including ground ambulances within the environmental conditions specified. Manual and Automated external defibrillation, External pacing, Diagnostic electrocardiography are intended for use on adult and pediatric patients.

The other monitoring functions (Electrocardiography (ECG, Heart Rate), Functional arterial oxygen saturation (SpO₂), Respiration (RESP)) are intended for use on adult, pediatric and neonatal patients.

Indications for Use

	Indications	Contraindications
Manual Defibrillation 2.3.1 2.3.2	Manual defibrillation is indicated for the termination of certain potentially fatal arrhythmias, such as ventricular fibrillation and ventricular tachycardia. Delivery of this energy in the synchronized mode is a method for treating atrial fibrillation, atrial flutter.	Defibrillation is contraindicated in the treatment of Pulseless Electrical Activity (PEA), such as idioventricular or ventricular escape rhythms, and in the treatment if asystole. Defibrillation is not intended for patient who is conscious and has a palpable pulse and is contraindicated for use on neonates.
Automated External Defibrillation 2.3.3	AED mode is to be used only on patients in cardiopulmonary arrest. The patient who must be unconscious, pulseless and not breathing normally before using the defibrillator to analyze the patient's ECG rhythm.	AED is not intended for patient who is conscious and has a palpable pulse and is contraindicated for use on neonates.
External Pacing	External pacing is indicated for symptomatic bradycardia in patients with a pulse. Increased heart rates in response to external pacing often suppress ventricular ectopic activity and may prevent tachycardia.	External pacing is contraindicated for the treatment of ventricular fibrillation and asystole. External pacing in the presence of severe hypothermia may be contraindicated and is contraindicated for use on neonates.

	Indications	Contraindications
Electrocardiography	The electrocardiogram is used to identify, diagnose and treat patients with cardiac disorders and is useful in the early detection.	No known contraindications
Pulse Oximetry	Pulse oximetry monitoring is intended to be used to monitor functional arterial oxygen saturation and pulse rate.	Pulse Oximetry is not intended for use with severe peripheral vascular disease and severe anemia (decreased Hemoglobin).

About This Manual

This manual explains how to set up and use the defibrillator/monitor.

Read the entire manual including the ***Safety Information*** section, before you operate the defibrillator/monitor.

Identifying the D100 Configurations

The following table identifies D100 defibrillator/monitor configurations and how they are indicated. The reference number and serial number are located on the back of the product.

All information in this manual, including the illustrations, is based on a device configured with the Battery, Pacing module, Mediana or Medtronic SpO₂ module, Infoware X module and Cradle/External adaptor.

Product Code	Description
Basic	0: Manual Defibrillation, Automated External Defibrillation, ECG (3/5 lead), Respiration, Printer
Top case	H: with paddles
	X: without paddles
Pacing module	P: Mediana Pacing module
	X: Not installed
SpO ₂ module	S: Medtronic SpO ₂ module
	D: Mediana SpO ₂ module
	X: Not installed
Infoware X module	W: Mediana InfoX-W(Wi-Fi) module
	E: Mediana InfoX-E(mobile) module
	X: Not installed
Cradle / External adaptor	A: Cradle / External AC Adaptor
	C: Cradle / External DC Adaptor
	X: Not installed

Note: 'X' in the product code is not indicated from the real product code.

This page is intentionally left blank.

DESCRIPTION OF THE DEFIBRILLATOR/MONITOR

Front Panel Components






Figure 1. Front panel components

1	Paddle	13	Service LED
2	Paddle energy level button (+,-)	14	Battery charging status LED
3	Patient type button	15	Mode selection knob
4	ECG Lead button	16	Energy level button(+,-)
5	ECG Size button	17	Charge button
6	LCD	18	Analyze button
7	Print button	19	Shock button
8	Soft key (1 st from left)	20	Sync button
9	Speaker	21	Event button
10	Paddle charge button	22	Alarm button
11	Handle	23	Multi function knob
12	AC/DC in LED		

Table 1. D100 controls

Symbol	Function
	Mode selection knob selects five modes of operation. (AED, Manual, Monitor, OFF, Pacing)
	Energy level button(+,-) selects the defibrillation energy level.
	Analyze button operates the advisory mode from the manual mode. analyzes the patient's ECG in the AED mode.
	Sync button selects sync or async mode.
	Charge button charges defibrillation energy.
	Shock button delivers a shock.
	Patient type button selects patient type. (Adult, Pediatric, Neonatal)
	ECG Lead button selects lead to display the first ECG waveform on the screen.
	ECG Size button selects to display the amplitude of an ECG waveform.
	Print button starts and stops the printing on operating.
	Alarm button pauses the audible alarm temporarily. turns off the audible alarm.
	Event button selects and records event data.
	Soft key x 4EA operates the assigned function according to the operation mode.
	Multi function knob selects the function for change, search, select, enter, add, cancel, etc..

Top Panel Components




Figure 2. Top panel components

- 1 Paddle
- 2 Shock Button
- 3 Paddle Cradle




Figure 3. Defibrillator/Patient monitor paddle controller

- 1 Shock Button
- 2 Energy Level Button(+)
- 3 Energy Level Button(-)
- 4 Shock Button
- 5 Charge Button

Rear Panel Components




Figure 4. Rear panel components

- | | | | |
|---|--------------------|---|----------------|
| 1 | Cradle Hook | 4 | AC inlet cover |
| 2 | DC Power Connector | 5 | Battery Hook |
| 3 | Battery | 6 | Rubber Foot |

Left Panel Components




Figure 5. Left panel components

- 1 Printer
- 2 Option Case
- 3 ECG Connector
- 4 SpO₂ Connector

Right Panel Components




Figure 6. Right panel components

- 1 Side Cover
- 8. 2 USB & SD Card Cover
- 3 Paddle/Pads Connector





Figure 7. Right panel components (Option)

- 1 Wireless Module
- 2 SD Card Insertion Hole
- 3 USB port
- 4 Paddle/Pads Connector

Table 2. Panel and label symbols

Symbols	Description	Symbols	Description
	AC indicator		Follow instructions for use
	Direct current		Separate Collection
	Battery charging status indicator		EU representative
	Service indicator		Manufacturer
	Type BF- Defibrillator Proof		Date of manufacture
	Type CF- Defibrillator Proof		Reference number
ECG	ECG connector		Serial number
SpO₂	SpO ₂ connector		Environmental shipping/storage altitude limitations
	AC power input rating		Environmental shipping/storage humidity limitations
	Equipotential terminal		Environmental shipping/storage temperature limitations
Rx ONLY	Prescription only device		Fragile
IP44	Dust and water resistance (D100)		Keep dry
IP22	Dust and water resistance (Cradle/External adapter)		This way up
	CE mark		Single patient use only

SETTING UP THE DEFIBRILLATOR/MONITOR

 WARNING	To ensure accurate performance and prevent defibrillator/monitor failure, do not expose the defibrillator/monitor to extreme moisture, including direct exposure to rain. Such exposure may cause inaccurate performance or device failure. Refer to Specification section.
 WARNING	The defibrillator/monitor should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the defibrillator/monitor should be observed to verify normal operation in the configuration it is to be used.
 WARNING	Make sure that the defibrillator/monitor speaker is not obstructed. Failure to do so could result in an inaudible alarm tone.
 CAUTION	Recharging the battery is strongly recommended when the battery has not been recharged for 6 or more months.
 CAUTION	Follow local government ordinances and recycling instructions regarding disposal or recycling of defibrillator/monitor components, including batteries.
 CAUTION	The use of accessories, cables, transducers and sensors sourced from manufacturers, which Mediana has not recommended may cause incorrect analysis.
 CAUTION	Electrical installation of the room or the building in which the defibrillator/monitor is to be used must comply with regulations specified by the country in which the equipment is to be used.

Unpacking and Inspection

The defibrillator/monitor is shipped in one carton. Examine the carton carefully for evidence of damage. Contact Mediana Technical Support Representative immediately if any damage is discovered. Refer to the **Maintenance** section for instructions on returning damaged items.

Note: Refer to the Performance Verification section in the service manual for detailed information.

Set the defibrillator/monitor to the user's intended position where the user can easily recognize the visual and audible monitoring conditions. Normally it is recommended to set at a distance of 1m from the user. Also the viewpoint is at any point within the base of a cone by an angle of 30° to the center of the monitoring display.

List of Components

The following items are standard in the package.

Table 3. Standard accessories

Items	Qty
D100 defibrillator/monitor	1
Operator's manual	1
AC Power cord	1
Print paper	2
Shock/Pacing Pads	1
Pads extension cable	1
ECG 3 leads Cable (SNAP)	1
For Mediana SpO ₂ Module SpO ₂ reusable sensor YM-1 SpO ₂ extension cable MEX03 * Only when Mediana SpO ₂ option is installed	1
For Medtronic SpO ₂ Module SpO ₂ reusable sensor DS-100A SpO ₂ extension cable DOC-10 * Only when Medtronic SpO ₂ option is installed	1
Rechargeable Li-ion Battery Pack	1

Table 4. Optional accessories

Items	Qty
External Paddle	-
Shock only pads	-
ECG 3 leads cable (SNAP)	-
ECG 3 leads cable (GRAB)	-
ECG 5 leads cable (SNAP)	-
ECG 5 leads cable (GRAB)	-
ECG 3 leads trunk cable	-
ECG 5 leads trunk cable	-
ECG 3 lead set (SNAP)	-
ECG 3 lead set (GRAB)	-
ECG 5 lead set (SNAP)	-
ECG 5 lead set (GRAB)	-
Service manual (English)	-
Cradle / External AC Adaptor	-
Cradle / External DC Adaptor	-
Mediana InfoX-W(Wi-Fi) module	-
Mediana InfoX-E(mobile) module	-
Carry bag	-

Optional items may be ordered if needed. Contact your local supplier for pricing and ordering information.

Power Cable Connections

⚠ WARNING	Do not connect to an electrical outlet controlled by a wall switch or dimmer because the defibrillator/monitor may be accidentally turned off.
⚠ WARNING	Do not turn off the Mediana defibrillator/monitor by using Mode selection knob or disconnecting the power cord or battery. The data stored will be damaged.
⚠ CAUTION	If the integrity of the AC power source is in doubt, the defibrillator/monitor must be operated from its internal battery.

AC Power

Make sure that the AC outlet is properly grounded and supplies the specified voltage and frequency (100-240V~ 50/60 Hz).




Figure 8. AC power connection

1. Connect the AC power cord to AC power connector on the rear of the Defibrillator/Monitor.
2. Connector the power cord to a grounded AC socket.
3. Connect a grounding wire if needed. After connection a ground wire socket to the equipotential grounding terminal on the rear, and then attach the clip end of the grounding wire to the medical device grounding terminal on the outer wall.
4. After connection, check if the **Battery charging status LED** and **AC/DC in LED** on the front of the Defibrillator/Monitor.

*Note: **Battery charging status LED** indicates battery is installed and being charged by the AC power.*

*Note: Even if the defibrillator/monitor is not turned on, the **Battery charging status LED** is lit when the AC power cord is connected into a mains outlet.*

Note: Do not place the defibrillator/monitor so that it is difficult to disconnect the AC power cord.

Cradle / External adapter

The DC power source can be used, when the defibrillator/monitor is used for the emergency condition in the moving car. Cradle uses AC or DC power.




Figure 9. Cradle installation

-
1. Connect an AC or DC power to cradle and verify that the **AC/DC Power Input LED** is lit.
 2. Set the defibrillator/monitor as Figure 9.
 3. Verify that the DC power input icon appears on the screen, **Battery charging status LED** and **AC/DC in LED** on the defibrillator/monitor's front panel is lit.

*Note: **Battery charging status LED** turns on in orange and **AC/DC in LED** turns on in green when battery is installed and being charged by the AC or DC power.*

*Note: Even if the defibrillator/monitor is not turned on, the **Battery charging status LED** is lit when a Cradle/External DC adaptor is connected into AC or DC input connector.*

Note: If the Battery charging status LED is not lit, check:

- the power cord
- the AC or DC input connector
- the Cradle/External DC adaptor
- the power/ mains outlet
- No Battery

Note: If the Battery charging status LED still is not lit although no problem is found, contact qualified service personnel or your local supplier for assistance.

Battery Operation

 WARNING	Do not disassemble, puncture, crush, heat above 100°C (212°F), or incinerate the battery. Be careful not to short the battery terminals because this could result in a fire hazard.
 WARNING	Mediana has no information regarding the performance or effectiveness of its Mediana defibrillator/monitors if other manufacturers' batteries or battery chargers are used. Using other manufacturers' batteries or battery chargers may cause the defibrillator/monitor to perform improperly and invalidate the safety agency certifications. Use only Mediana's battery or battery charger.
 WARNING	Grounding reliability can only be achieved when the equipment is connected to an equivalent receptacle marked 'HOSPITAL ONLY' or 'HOSPITAL GRADE'. If the grounding integrity of the line cord or AC receptacle is in doubt, operate on battery only.
 WARNING	Do not turn off the Mediana defibrillator/monitor by using Mode selection knob or disconnecting the power cord or battery. The data stored will be damaged.
 CAUTION	Recharging the battery is strongly recommended when it has not been fully recharged for 6 or more months.
 CAUTION	Do not install the battery into the defibrillator/monitor when storage may exceed 90 days. Battery damage may occur.
 CAUTION	When the voltage of the battery is very low, it is a possibility of not operating.
 CAUTION	Check for the battery-in-use indication when the defibrillator/monitor is operating on mains and apply corrective action.
 CAUTION	Do not operate the defibrillator/monitor without a battery. Keep a fully charged spare battery pack with the defibrillator/monitor at all times.
 CAUTION	Partial charge of battery results in a shortened service life.
 CAUTION	Storing at temperatures above 40°C (104°F) for extended periods of time will significantly reduce a battery's life-expectancy.
 CAUTION	Using an improperly maintained battery to power the defibrillator/monitor may cause power failure without warning. Use the appropriate battery charger to charge batteries.
 CAUTION	Battery pins in the defibrillator/monitor may be damaged if batteries are dropped or forced into battery wells. Inspect pins routinely for signs of damage. Keep batteries installed at all times except when defibrillator/monitor is removed from service for storage.
 CAUTION	When storing the defibrillator/monitor for an extended period of time, the battery should be removed from the defibrillator/monitor.
 CAUTION	Stored batteries lose charge. Failure to charge a stored battery before use may cause device power failure without warning. Always charge a stored battery before placing it in active use.

Note: It is recommended that keep AC power source connected to the defibrillator/monitor when not in use. This will ensure a fully charged battery whenever it is needed.

Note: As the battery is used and recharged over a period of time, the amount of time between the onset of the low battery alarm and the defibrillator/monitor shut-off may become shorter. It is recommended for service personnel to check periodically or replace the battery if necessary.

Operating the defibrillator/monitor on Battery Power

The defibrillator/monitor has the rechargeable battery that can be used to power the defibrillator/monitor when AC power source is not available. The battery status icon appears on the screen when the defibrillator/monitor is on battery power.




Figure 10. Battery placement

1. Turn off the defibrillator/monitor.
2. Push the **Battery hook**.
3. Insert the battery into the defibrillator/monitor carefully.

Table 5. Front panel indications for power source

Power Connections	Front Panel Indications
AC source	AC power input icon appears on the screen. AC/DC in LED is lit.
DC source	DC power input icon appears on the screen. AC/DC in LED is lit.
Battery	Battery status icon appears on the screen.

The defibrillator/monitor cannot operate with a completely discharged battery. Before turning on the defibrillator/monitor with a battery that has been completely discharged, plug the defibrillator/monitor into an AC outlet to charge the battery for a minimum of 3 minutes. The defibrillator/monitor may then be powered on.

A new, fully charged battery will provide 5 hours monitoring operation in *Monitor mode* or at least 200 times of shock delivery in 200J selected energy in *Manual mode or AED mode* or 2 hours monitoring operation in *Pacing mode* (pacing rate: 80 ppm, pacing output: 60mA, load resistance: 50Ω) under the following conditions:

- No audible alarm sound
- No data Output, No communication (Communication system is installed)
- No printing out (Printer module is installed)
- All monitoring parameters are active
- LCD Backlight: default
- Ambient temperature at 25°C

Battery Status Indication

When operating on batteries, the battery status icon in the lower part of the display indicates the battery charge condition. See Table 6.

Table 6. The defibrillator/monitor battery status icon

Battery Status Icons	Battery Status
	Normal, Fully charged
	More than 4/5 from fully charged
	More than 3/5 from fully charged
7. 	More than 2/5 from fully charged
	More than 1/5 from fully charged
	Battery charging
	Battery charging error
	Low
	Critical low

Note: The battery status icon is displayed in 5 levels depending on the amount of battery. User can check the approximate amount of battery by the number of bars inside the battery status icon.

A low priority alarm occurs when the remaining battery power is only enough for 15 minutes of operation. The **'System: Low battery condition'** message is displayed on the screen and visual alarm window is lit with yellow. Connect the AC power or replace the battery with a fully charged battery immediately when alarm occurs.

This alarm audio cannot be paused while running on battery power. Connecting the defibrillator/monitor to AC power will stop the alarm.

A high priority alarm occurs for about 5 minutes before the defibrillator/monitor shuts off. The **'System: Critically low battery condition'** message is displayed on the screen and the visual alarm window is lit with red. After that, the monitor will automatically shut down. Connect the defibrillator/monitor to an AC power source immediately to avoid any loss of trend data or settings.

User cannot use printer while the device is in low battery and critically low battery condition.

Note: The battery will not be charged for safety if the operating temperature exceeds 40°C.

Note: In critically low battery condition, the monitor is not turned on once the monitor is turned off.

Measurement Cable Connections

 WARNING	For best product performance and measurement accuracy, use only accessories supplied or recommended by Mediana. Use accessories according to the manufacturer's directions for use and your facility's standards. Use only accessories that have passed the recommended biocompatibility testing in compliance with ISO10993-1.
 WARNING	Connect the sensor, probe, cuff/hose or transducer firmly into socket and do not use a damaged sensor, probe or cuff/hose.
 WARNING	To avoid damage to the cable, always hold by the connector rather than the cable, when connecting or disconnecting either end.
 WARNING	The sensor connector should not be connected to anything other than a sensor.

Note: Both frequent checks by the operator on a daily basis and more comprehensive technical checks less frequently are covered by this requirement in order to detect mechanical damage and damage to cables, etc.

ECG Cables and Leads

1. Connect an ECG cable to the ECG connector making sure that the connector arrow is pointing panel (see Figure 5)
2. Connect an ECG leadwire cable to the ECG cable.

Defibrillator Paddle and Pads

1. Connect a paddle or pads to paddle / pads connector on the defibrillator/monitor's right panel. (see Figure 6)
2. Use **ECG lead button** to select pads. (If paddle is connected, the mode is automatically changed to paddle.)

SpO₂ Cables and Sensors (if configured with SpO₂ option)

1. Select an appropriate sensor for the patient and desired application. (Refer to the **SpO₂ Monitoring** section.)
2. Connect the extension cable to the SpO₂ connector on the defibrillator/monitor's front panel. (see Figure 5)
3. Attach the sensor to the end of the cable.

Note: If lead wire cable, sensor and transducer are not connected firmly, the defibrillator/monitor could lose signal from patient.

USING THE DEFIBRILLATOR/MONITOR

 WARNING	Each time the defibrillator/monitor is used, check alarm limits to make sure that they are appropriate for the patient being monitored.
 WARNING	If different alarm presets are used for the same or similar equipment in any single area, e.g. an intensive care unit or cardiac operating room, a potential hazard can exist.
 WARNING	Keep patients under close surveillance when monitoring. It is possible, although unlikely, that radiated electromagnetic signals from sources external to the patient and the defibrillator/monitor can cause inaccurate measurement readings. Do not rely entirely on the defibrillator/monitor readings for patient assessment.
 WARNING	Be aware of patient cables, including ECG monitoring equipment when used with high frequency surgical equipment.
 WARNING	When inaccurate analysis is suspected because of motion artifact do not analyze in a moving vehicle. Motion artifact may affect the ECG signal resulting in an inappropriate SHOCK ADVISED or NO SHOCK ADVISED message. Motion detection may delay analysis. Stop vehicle and stand clear of patient during analysis.
 WARNING	To avoid risk of electrical shock, do not touch gelled area of the pads while pacing or shock. When defibrillating with paddles use your thumbs to operate the Shock button in order to avoid inadvertent operator shock. No portion of the hand should be near the paddle plates.
 WARNING	Pediatric defibrillation energy levels should be set based on site-specific clinical protocols.
 WARNING	The metronome sounds do not indicate information regarding the patient's condition. Because patient status can change in a short time, the patient should be assessed at all times. Do not perform CPR on a patient who is responsive or is breathing normally.
 WARNING	Place the patient on a firm surface before performing CPR.
 WARNING	Ventricular fibrillation may be induced with improper synchronization. Do not use the ECG from another monitor (slaving) to synchronize the defibrillator/monitor's discharge. Always monitor the patient's ECG directly through the defibrillator/monitor's ECG cable or pad cable.
 WARNING	Pitted or damaged paddles may cause patient skin burns during defibrillation.
 WARNING	The ECG rhythm analysis function does not warn the operator of patient asystole, as it is not a shockable rhythm.
 WARNING	Determination of electrical capture should only be performed by viewing the ECG on the screen with its ECG cable directly attached to the patient.
 WARNING	Connect the device only to a three-wire, grounded, hospital grade receptacle. The three-conductor plug must be inserted into a properly wired three-wire receptacle; if a three-wire receptacle is not available, a qualified electrician must install one in accordance with the governing electrical code. Do not under any circumstances remove the grounding connector from the power plug. Do not use extension cords or adapters of any type. The power cord and plug must be intact and undamaged.
 CAUTION	Only use thumbs to depress the paddle Shock button. Failure to do so could result in the inadvertent depression of the energy select buttons, causing the defibrillator/monitor to disarm itself.

⚠ CAUTION	Do not discharge the defibrillator/monitor except as indicated in the instructions.
⚠ CAUTION	Do not permit gel to accumulate between the paddle electrodes on the chest wall (gel bridge). This could cause burns and reduce the amount of energy delivered to the heart.
⚠ CAUTION	Changing the selected energy while the defibrillator/monitor is charging or charged will cause the defibrillator/monitor to disarm itself. Press the Charge button again to charge the defibrillator/monitor.
⚠ CAUTION	If using defibrillation pads, make sure that the size of the pad is large enough to cover the entire paddle electrode area.
⚠ CAUTION	Pads should be replaced after 8 hours of continuous pacing to ensure maximum patient benefit.
⚠ CAUTION	Be sure to safely discharge external paddles.
⚠ CAUTION	Do not use high conductivity gel manufactured by other than manufacturer recommended for defibrillation using paddle. It may cause burns to the patients or device failure.

Note: Check if the defibrillator/monitor correctly operates according to the displays or signal sounds during the measurement.

Note: If the power has some changes, the buzzer sound of the monitor will occur.

Note: To reduce the risks due to the software error, use the new software.

Note: If there is a software error, please contact qualified service personnel or your local supplier.

Note: If the software or hardware is not restored after rebooting the device, the defibrillator / monitor will beep. Forcibly shut down the equipment by disconnecting the power and contact your local dealer or headquarter service.

Note: After the power turns on, check if the battery status indicator, display and other indicators operate normally.

Turning On and Off the Defibrillator/Monitor

Before using the defibrillator/monitor, confirm that the defibrillator/monitor is working properly and is safe to use as described below.

⚠ CAUTION	Look for display motion before accepting any displayed data as a current measurement.
------------------	--

Note: If unusable sound like buzzer can be heard, do not use the defibrillator/monitor. Instead, please contact qualified service personnel or your local supplier.

Turning On

The defibrillator/monitor can start to be used by rotating the Mode select knob of the device that turned off.

1. Rotate **Mode select knob** to select one mode from AED, manual, monitor or pacing mode.
2. While the defibrillator/monitor turns on, the LEDs of Alarm button and Shock button are turned on.
3. The defibrillator/monitor automatically starts the Power On Self Test (POST) that does the test for the function of the device.
4. When the Power On Self Test (POST) is completed successfully, the defibrillator/monitor will display the selected mode screen.

Note: If internal problem is detected or the final self test result is failed during the Power

On Self Test (POST), the alarm of the defibrillator/monitor will occur and the error code is displayed. When the error code is displayed, please contact qualified service personnel or your local supplier.




Figure 11. Mode select knob

- Pacing mode - to enable pacing mode for Demand or Fixed mode pacing.
- Manual mode - to enable manual mode for asynchronous or synchronous defibrillation.
- Monitor mode - to enable monitor mode for 3- or 5-lead ECG monitoring and vital signs trending, or monitoring of optional parameters.
- AED mode - to enable AED mode for semi-automated external defibrillation.

Turning Off

The defibrillator/monitor can start to be used by rotating the Mode select knob of the device that turned off.

1. To turn off the defibrillator/monitor, rotate **Mode select knob** to select Off mode. After selecting Off mode, the settings will be stored and the defibrillator/monitor will be turned off.

Displays

The screen below is a default screen in manual mode. The screen may be different according to other modes. For screens in other modes, refer to the **AED mode**, **manual mode**, **monitor mode** and **pacing mode** section.

A typical screen in Manual Mode is shown below.




Figure 12. Displays (Manual mode)

1	Patient type	8	Alarm volume
2	Waveform area	9	Connection type
3	Operating information	10	CMS status
4	Soft key	11	Patient information
5	Runtime	12	Physiological/Technical alarm message
6	Numeric area	13	Date/Time display
7	Power status		

Note: CMS (Central Monitoring System) bed number is not displayed on the screen if there is no set bed number.