
en	Instructions for use/Technical description AESCULAP® Reusable trocar system 5,5 mm – 13 mm
USA	Note for U.S. users This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusaifus.com . If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.
de	Gebrauchsanweisung/Technische Beschreibung AESCULAP® Wiederverwendbares Trokarsystem 5,5 mm, 13 mm
fr	Mode d'emploi/Description technique AESCULAP® Système de trocart réutilisable 5,5-13 mm
es	Instrucciones de manejo/Descripción técnica AESCULAP® Sistema de trocar reutilizable 5,5 mm – 13 mm
it	Istruzioni per l'uso/Descrizione tecnica AESCULAP® Sistema trocar riutilizzabile da 5,5 mm-13 mm
pt	Instruções de utilização/Descrição técnica AESCULAP® Sistema de trocarte reutilizável de 5,5 mm, 13 mm
nl	Gebruiksaanwijzing/Technische beschrijving AESCULAP® Herbruikbaar trocarsysteem 5,5 mm – 13 mm
da	Brugsanvisning/Teknisk beskrivelse AESCULAP® Trokarsystem til flergangsbrug 5,5 mm – 13 mm
no	Bruksanvisning/ beskrivelse AESCULAP® Gjenbrukbart trokarsystem 5,5 mm-13 mm
sv	Bruksanvisning/Teknisk beskrivning AESCULAP® Återanvändbart trocarsystem 5,5 mm – 13 mm
fi	Käyttöohjeet/Tekninen kuvaus AESCULAP® Uudelleenkäytettävä trokaarijärjestelmä 5,5 mm – 13 mm
et	Kasutusjuhend/Tehniline kirjeldus AESCULAP® Korduvkasutatav trokaarisüsteem 5,5-13 mm
lv	Lietošanas instrukcija/Tehniskais apraksts AESCULAP® Atkārtoti lietojama trokāra sistēma 5,5-13 mm
lt	Naudojimo instrukcija/Tekninis aprašymas AESCULAP® Daugkartinė trokaro sistema, 5,5-13 mm
ru	Инструкция по применению/Техническое описание AESCULAP® Система многоразовых трокаров 5,5 мм – 13 мм
cs	Návod k použití/Technický popis AESCULAP® Systém trokaru 5,5 mm – 13 mm k opakovanému použití
pl	Instrukcja użytkowania/Opis techniczny AESCULAP® System trokarów wielokrotnego użytku 5,5-13 mm
sk	Návod na používanie/Technický opis AESCULAP® Trokárový systém na opakované použitie 5,5 mm – 13 mm
hu	Használati utasítás/Műszaki leírás AESCULAP® Újrahasznátható trokárrendszer, 5,5-13 mm
sl	Navodila za uporabo/Tehnični opis AESCULAP® Trokarski sistem za večkratno uporabo, 5,5-13 mm
hr	Upute za uotrebu/Tehnički opis AESCULAP® Sustav trokara za višekratnu uporabu 5,5 mm – 13 mm
ro	Manual de utilizare/Descriere tehnică AESCULAP® Sistem de trocar reutilizabil 5,5 mm – 13 mm
bg	Инструкции за употреба/техническо описание AESCULAP® Троакарна система за многократна употреба 5,5 – 13 mm
tr	Kullanım Kılavuzu/Teknik Açıklama AESCULAP® Tekrar kullanılabilir trokar sistemi 5,5 mm – 13 mm
el	Οδηγίες χρήσης/Τεχνική περιγραφή AESCULAP® Επαναχρησιμοποιήσιμο σύστημα τροκάρ 5,5 mm – 13 mm

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 (0) 7461 95-0 | Fax +49 (0) 7461 95-26 00 | www.bbraun.com

AESCULAP® – a B. Braun brand

TA011883 2023-08 Change No. AE0063131

Legend

- 1 Sealing cap
- 2 Trocar body
- 3 Silicone lip valve
- 4 Ring (with or without Luer Lock insufflation connector)
- 5 Sealing cap for Luer Lock connector
- 6 Spacer with Inflation Stopcock
- 7 Spring cap
- 8 Stopcock
- 9 Coding ring
- 10 Trocar sleeve
- 11 Fixation screw
- 12 Obturator
- 13 Trocar
- 14 Suture material fixation
- 15 Cone attachment

1 About this document

Note

General risk factors associated with surgical procedures are not described in these instructions for use.

1.1 Scope

These instructions for use apply for trocars.

Note

The applicable CE mark for the product can be seen on the label or packaging of the product.

Note

Instructions for use and further information about B. Braun/AESCULAP products can be found on the B. Braun eFU website at eifu.bbraun.com

1.2 Safety messages

Safety messages make clear the dangers to patient, user and/or product that could arise during the use of the product. Safety messages are labeled as follows:

⚠ WARNING

Indicates a possible threat of danger. If not avoided, minor or moderate injury may result.

⚠ CAUTION

Indicates a possible threat of material damage. If not avoided, the product may be damaged.

2 Clinical use

2.1 Available sizes

The reusable trocar system of the EJ series is available in the following variants:

- Working lengths:
 - 60 mm
 - 110 mm
 - 150 mm
- Diameter:
 - 5.5 mm (red)
 - 10 mm (green)
 - 13 mm (yellow)

For a complete overview of our trocar systems, please see the Endoscopy Online Catalog under www.endoscopy-catalog.com

2.2 Areas of use and limitations of use

2.2.1 Intended use

The reusable trocar system is used in operations in laparoscopic general surgery, gynecology and urology. It serves to create and maintain an approach to the operating field for instruments and endoscopes in laparoscopic operations. Endoscopic instruments of 5 mm (5 mm trocar), 5–10 mm (10 mm trocar) and of 5–12 mm (13 mm trocar) can be inserted.

It serves to create and maintain an approach to the operating field for instruments and endoscopes in laparoscopic operations using the HASSON technique. Endoscopic instruments of up to 5–10 mm (10 mm trocar) diameter can be inserted.

They may only be used in conjunction with the corresponding silicone flap valve.

2.2.2 Indications

Note

The manufacturer is not responsible for any use of the product against the specified indications and/or the described applications.

For indications, see Intended use.

2.2.3 Contraindications

Do not use if laparoscopic operating techniques are contraindicated.

2.3 Safety information

2.3.1 Clinical user

General safety information

To prevent damage caused by improper setup or operation, and to not compromise the manufacturer warranty and liability:

- ▶ Use the product only according to these instructions for use.
- ▶ Follow the safety and maintenance instructions.
- ▶ Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge and experience.
- ▶ Store any new or unused products in a dry, clean, and safe place.
- ▶ Prior to use, check that the product is in good working order.
- ▶ Keep the instructions for use accessible for the user.

Note

The user is obligated to report all severe events in connection with the product to the manufacturer and the responsible authorities of the state in which the user is located.

Notes on surgical procedures

It is the user's responsibility to ensure that the surgical procedure is performed correctly.

Appropriate clinical training as well as a theoretical and practical proficiency of all the required operating techniques, including the use of this product, are prerequisites for the successful use of this product.

The user is required to obtain information from the manufacturer if there is an unclear preoperative situation regarding the use of the product.

2.3.2 Product-specific safety information

- ▶ Only use the product with visual control.
- ▶ Only combine Aesculap trocar components of the EJ series.
- ▶ Only combine trocar sleeves and trocar obturators with the same color code (diameter)/label and the same working length.
- ▶ Use cone attachment 15 only with smooth trocar shaft EJ432R.
- ▶ Use trocars with silicone flap valve.
- ▶ Use suture material fixation 14 only with smooth trocar shaft 10.
- ▶ Observe handling instructions in the instrument-specific instructions for use.

2.3.3 Sterility

Non-sterile packaged products

The product is supplied non-sterile and intended to be used in sterile condition.

- ▶ Clean the new product after removing its transport packaging and prior to its initial sterilization.

2.4 Application

⚠ WARNING

Risk of injury and/or malfunction!

- ▶ Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.
- ▶ Always carry out a function test prior to each use of the product.
- ▶ To prevent damage to the sealing element, apply appropriate care when inserting any instruments.
- ▶ If possible, insert instruments in their closed position, straight and central through the sealing element.
- ▶ Prior to inserting the trocar into the patient, prepare an abdominal pneumoperitoneum e.g. with a Verress cannula.
- ▶ Apply skin incisions in order to avoid excessive application of force.
- ▶ Position trocars under intra-abdominal visual control.

⚠ CAUTION

Malfunction due to incompatible instruments!

▶ Check for mutual compatibility of the trocar system and instruments. To do this, carefully insert the instrument into the trocar and check for patency.

- ▶ Mounting the trocar, see Assembly.
- ▶ When using the trocar sleeve and the spacer with inflation stopcock 6: close stopcock 8. For this, turn the valve 90° to the right.
- ▶ Insert the trocar mandrel 12 into the trocar 13.
- ▶ Insert the trocar into the patient by alternating left/right rotating movements, applying even and controlled pressure.
- ▶ For trocar sleeves 10 with thread: turn in a clockwise direction until the trocar is in the desired position.
- ▶ Pull out the obturator 12. Trocar sleeve 10 remains positioned in the patient.
- ▶ When using the suture material fixation 14: fixate trocar shaft 10 in the incision with previously placed fixation suture. Wind the threads of the suture around suture material fixation 14 at least twice, into the slots.
- ▶ When using the cone attachment 15:
 - Slide cone attachment onto smooth trocar sleeve and tighten locking screw.
 - Fixate cone attachment with prepared fixation sutures in the incision.
 - Wind the threads of the fixation sutures around the thread fixation at least twice, into the slots.
- ▶ When connecting the inflation to the inflation stopcock 8: connect, start inflation and open stopcock 8.
- ▶ Carry out intra-abdominal visual control.
- ▶ Post application, unwind the suture thread from the suture material fixation 14.
- ▶ Post application, unwind the suture thread from the cone attachment 15.
- ▶ After use, with trocar sleeve 10 with thread: turn trocar sleeve counterclockwise until the thread is fully released.

Note

Apply rotary movement to remove hook-shaped instruments.

3 Validated processing procedure

3.1 General safety instructions

Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the processing of products.

Note

Mechanical processing should be favored over manual cleaning as it gives better and more reliable results.

Note

Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.

Note

Sterilization should always be favored over disinfection under all circumstances. Only if final sterilization is not possible should disinfection with a virucidal agent be preferred.

Note

Up-to-date information about processing and material compatibility can be found on the B. Braun eFU site at eifu.bbraun.com

The validated steam sterilization procedure was carried out in the AESCULAP sterile container system.

3.2 General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.

Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel.

Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying.

Additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.
- ▶ Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.
- ▶ Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at www.a-k-i.org, link to "AKI-Brochures", "Red brochure".

3.3 Service life

Influences of the processing using the validated procedure which lead to damage to the product are not known.
Careful visual and functional inspection before each use is the best opportunity to recognize a product that is no longer functional, see Inspection.

3.4 Preparations at the place of use

- If applicable, rinse non-visible surfaces preferably with deionized water, with a disposable syringe for example.
- Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

3.5 Preparing for cleaning

- Disassemble the product prior to cleaning, see Disassembly.
- Remove the sealing cap from the Luer lock connector.
- If there is a spacer with inflation stopcock 6 open stopcock 8.

3.6 Disassembly

- Remove sealing cap 1 from trocar body 2.
- Unscrew the body of the trocar 2 from the trocar shaft 10.
- Remove silicone flap valve 3 from trocar body 2.
- If necessary, remove the spacer ring 4/6 from the trocar sleeve 10.
- Remove the sealing cap for Luer Lock connector 5 if applicable.
- If necessary, remove the code ring 9 from the trocar sleeve 10.
- If necessary, slide suture material fixation 14 from trocar sleeve 10 and completely unscrew fixation screw 11.
- If necessary, loosen fixation screw of cone attachment 15 and remove cone attachment from trocar sleeve 10.

Spacer with Inflation Stopcock

- Unscrew the spring cap 7 from the spacer with the inflation stopcock 6.
- Remove stopcock 8.

3.7 Cleaning/Disinfection

3.7.1 Product-specific safety information on the processing method

Damage to or destruction of the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Following the manufacturer's instructions, use cleaning and disinfecting agents,
 - be approved for plastic material (thermoplastics, silicon) and high-grade steel,
 - that do not attack softeners (e.g., in silicone).
- Observe specifications regarding concentration, temperature and exposure time.
- Do not exceed the maximum allowable disinfection temperature of 95 °C.

3.7.2 Validated cleaning and disinfection procedure

Validated procedure	Specific requirements	Reference
Manual cleaning with immersion disinfection	<ul style="list-style-type: none">► Use a suitable cleaning brush, e.g. TA007747► Use a disposable syringe 20 ml.► When cleaning instruments with movable hinges, ensure that these are in an open position and, if applicable, move the hinge while cleaning.► Drying phase: Use a lint-free cloth or medical compressed air.	see Manual cleaning/disinfecting and subsection: ■ see Manual cleaning with immersion disinfection
Manual cleaning with ultrasound and immersion disinfection ■ EJ430R	<ul style="list-style-type: none">► Use a suitable cleaning brush, e.g. TA007747► Use a disposable syringe 20 ml.► When cleaning instruments with movable hinges, ensure that these are in an open position and, if applicable, move the hinge while cleaning.► Drying phase: Use a lint-free cloth or medical compressed air.	see Manual cleaning/disinfecting and subsection: ■ see Manual cleaning with ultrasound and immersion disinfection
Manual pre-cleaning with brush and subsequent mechanical alkaline cleaning and thermal disinfection	<ul style="list-style-type: none">► Use a suitable cleaning brush, e.g. TA007747► Use a disposable syringe 20 ml.► Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots).► Connect components with lumens and channels directly to the rinsing port of the injector carriage.► Place the product on the tray with all product links and joints open.► Use jetting lances or irrigation sleeves in the injector vehicle to rinse out the lumens/ducts.	see Mechanical cleaning/disinfection with manual pre-cleaning and subsection: ■ see Manual pre-cleaning with a brush ■ see Mechanical alkaline cleaning and thermal disinfecting
Manual pre-cleaning with ultrasound and brush, and subsequent mechanical alkaline cleaning and thermal disinfection ■ EJ430R	<ul style="list-style-type: none">► Use a suitable cleaning brush, e.g. TA007747► Use a disposable syringe 20 ml.► Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots).► Connect components with lumens and channels directly to the rinsing port of the injector carriage.► Place the product on the tray with all product links and joints open.► Use jetting lances or irrigation sleeves in the injector vehicle to rinse out the lumens/ducts.	see Mechanical cleaning/disinfection with manual pre-cleaning and subsection: ■ see Manual pre-cleaning with ultrasound and brush ■ see Mechanical alkaline cleaning and thermal disinfecting

3.8 Manual cleaning/disinfecting

- Prior to manual disinfecting, allow water to drip off for a sufficient length of time to prevent dilution of the disinfectant.
- After manual cleaning/disinfection, check visible surfaces visually for residues.
- Repeat the cleaning/disinfecting process if necessary.

3.8.1 Manual cleaning with immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ≈ 9*
II	Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfection	RT (cold)	5	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ≈ 9*
IV	Final rinse	RT (cold)	1	-	FD-W	-
V	Drying	RT	-	-	-	-

D-W: Drinking water
FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)
RT: Room temperature
*Recommended: B. Braun Stabimed fresh
► Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- Drain any remaining water fully.

Phase III

- Fully immerse the product in the disinfectant solution.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- Rinse lumens at least 5 times at the beginning of the exposure time using an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

Phase IV

- Rinse/flush the product thoroughly (all accessible surfaces).
- Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- Rinse lumens with an appropriate disposable syringe at least five times.
- Drain any remaining water fully.

Phase V

- Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

3.8.2 Manual cleaning with ultrasound and immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ≈ 9*
II	Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfection	RT (cold)	5	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ≈ 9*
IV	Final rinse	RT (cold)	1	-	FD-W	-
V	Drying	RT	-	-	-	-

D-W: Drinking water
FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)
RT: Room temperature
*Recommended: B. Braun Stabimed fresh
► Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- Drain any remaining water fully.

Phase III

- Fully immerse the product in the disinfectant solution.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- Rinse lumens at least five times at the beginning of the exposure time with an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

Phase IV

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- Rinse lumens with an appropriate disposable syringe at least five times.
- Drain any remaining water fully.

Phase V

- Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

3.9 Mechanical cleaning/disinfection with manual pre-cleaning

Note

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. compliance with EN ISO 15883).

Note

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

3.9.1 Manual pre-cleaning with a brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ≈ 9*
II	Rinsing	RT (cold)	1	-	D-W	-

D-W: Drinking water

RT: Room temperature

*Recommended: B. Braun Stabimed fresh

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

3.9.2 Manual pre-cleaning with ultrasound and brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ≈ 9*
II	Rinsing	RT (cold)	1	-	D-W	-

D-W: Drinking water

RT: Room temperature

*Recommended: B. Braun Stabimed fresh

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

3.9.3 Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	<div><div></div>Concentrate, alkaline:<ul style="list-style-type: none">- pH = 13- < 5 % anionic surfactant<div><div></div>0.5% working solution<ul style="list-style-type: none">- pH = 11*</div></div>
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfecting	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Recommended: B. Braun Helimatic Cleaner alkaline

- Check visible surfaces for residues after mechanical cleaning/disinfecting.

3.10 Inspection

- Allow the product to cool down to room temperature.
- Dry the product if it is wet or damp.

3.10.1 Visual inspection

- Ensure that all soiling has been removed. In particular, pay attention to mating surfaces, hinges, shafts, recessed areas, drill grooves and the sides of the teeth on rasps.
- If the product is dirty: repeat the cleaning and disinfection process.
- Check the product for damage, e.g. insulation or corroded, loose, bent, broken, cracked, worn or severely scratched and fractured components.
- Check the product for missing or faded labels.
- Check the surfaces for rough spots.
- Check the product for burrs that could damage tissue or surgical gloves.
- Check the product for loose or missing parts.
- Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical service.

3.10.2 Assembly

⚠ CAUTION

Leaks caused by damaged or defective sealing elements (sealing cap, silicone flap valve, sealing cap for Luer Lock connector, code ring!)

- **Replace damaged/defective sealing elements.**
- Put code ring 9 on trocar shaft if applicable 10.
- Put spacer ring 4/6 on trocar shaft if applicable 10.
- Attach the seal cap for Luer Lock connection 5 if applicable.
- Insert silicone flap valve 3 in trocar body 2. Make certain that silicone flap valve 3 is pushed down to its position below the thread.
- Screw trocar body 2 onto trocar shaft 10.
- Place sealing cap 1 on the trocar body 2.
- Slide suture material fixation 14 onto trocar shaft 10 and screw in fixation screw 11.

Spacer with Inflation Stopcock

- Put stopcock 8 into the spacer with inflation stopcock 6.
- Screw on spring cap 7.

3.10.3 Functional test

- Assemble disassembled products, see Assembly.
- Check that the product functions correctly.
- Check that all moving parts are working properly (e.g. hinges, locks/latches, sliding parts etc.).
- Check for compatibility with associated products.
- Immediately put aside inoperative products and send them to Aesculap Technical Service, see Technical service.

3.11 Packaging

- Place the product in its holder or on a suitable tray. Ensure that sharp edges are covered.
- Package trays appropriately for the sterilization process (e.g. in AESCULAP sterile containers).
- Ensure that the packaging provides sufficient protection against contamination of the product during storage.

3.12 Steam sterilization

⚠ CAUTION

Damage or leaking of the silicon flap valve caused by incorrect handling during sterilization!

- **Sterilize trocar mandrels separately (not inserted in trocar sleeve or valve unit).**
- Remove sealing cap 1 before sterilization, see Disassembly
- Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g., by opening any valves and faucets).
- Validated sterilization process
 - Steam sterilization using fractional vacuum process
 - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
 - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min
- If several devices are sterilized at the same time in the same steam sterilizer: Ensure that the maximum permitted load according to the manufacturers' specifications is not exceeded.

3.13 Storage

- Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

4 Maintenance and service

4.1 Technical service

⚠ CAUTION

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

- **Do not modify the product.**
- **For service and repairs, please contact your national B. Braun/AESCULAP agency.**

Service addresses

Aesculap Technischer Service

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78532 Tuttlingen / Germany

Phone: +49 7461 95-1601

Fax: +49 7461 16-2887

E-Mail: ats@aesculap.de

Other service addresses can be obtained from the address indicated above.

4.2 Accessories / spare parts

Arl. no.	Designation
EJ570P	Silicone flap seal, 5,5/7,0 mm (pack of 20)
EJ571P	Silicone flap seal, 10/12,5 mm (pack of 20)
EJ572P	Silicone flap seal, 13 mm (pack of 20)
EJ700251	Spacer ring for Luer Lock inflation connector, 5,5/7,0 mm
EJ750251	Spacer ring for Luer Lock inflation connector, 10/12,5 mm
EJ651R	Spacer ring with inflation stopcock, 5,5 mm
EJ650R	Spacer ring with inflation stopcock, 10/12,5 mm
EJ700250	Spacer ring without inflation connector 5,5/7,0 mm
EJ750250	Spacer ring without inflation connector 10/12,5 mm
SK271C	Metal stopcock with Luer Lock attachment
EJ751251	Sealing cap for Luer Lock connector (pack of 20)
EJ702250	Coding ring, red, 5,5 mm (pack of 20)
EJ752250	Coding ring, green, 10 mm (pack of 20)
EJ772250	Coding ring, black, 12,5 mm (pack of 20)
EJ644P	Coding ring, yellow, 13 mm (pack of 20)
EJ446255	Sealing cap, red, 5,5 mm (pack of 20)
EJ449255	Sealing cap, green, 10 mm (pack of 20)
EJ450255	Sealing cap, black, 12,5 mm (pack of 20)
EJ643P	Sealing cap, yellow, 13 mm (pack of 20)
EJ640P	Sealing cap set, green, 5 mm/10 mm (pack of 5)
EJ641P	Sealing cap set, black, 5 mm/10 mm/12,5 mm (pack of 5)
EJ642P	Sealing cap set, yellow, 5 mm/10 mm/13 mm (pack of 5)

See Endoscopy online catalogue www.endoscopy-catalog.com

5 Disposal

WARNING

Risk of infection due to contaminated products!

► Adhere to national regulations when disposing of or recycling the product, its components and its packaging.

WARNING

Risk of injury due to sharp-edged and/or pointed products!

► When disposing of or recycling the product, ensure that the packaging prevents injury by the product.

Note

The user institution is obliged to process the product before its disposal, see Validated processing procedure.

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