EVIS X1 Video System Center

CV-1500

A Unified Platform with 5 LED Spectrum Technology

Eil. 2.16

A Unified Platform with 5 LED Spectrum Technology

By integrating the LED light source with the video processor, Olympus has developed a powerful system that is much more compact and lightweight than the predecessors*1.

Broad Compatibility

The CV-1500 can be connected to many different types of endoscopes, providing access to a wide variety of endoscopy-supporting functions.

Enhanced Observations Eil. Nr. 2.2.1 Eil. Nr. 2.2.2

In addition to conventional white light and NBI (Narrow Band Imaging) and AFI (Auto Fluorescence Imaging) observation, the CV-1500 offers three other powerful enhanced observations to improve diagnostic and therapeutic capability:

Eil. Nr. 2.3.1 TXI (Texture & Color Enhancement Imaging) optimizes the structure, color tone and brightness of the mucosal surface.

Eil. Nr. 2.2.3 RDI (Red Dichromatic Imaging) improves visibility of deep blood vessels and bleeding points.

Eil. Nr. 2.3.2 BAI-MAC (Brightness Adjustment Imaging with MAintenance of Contrast) improves brightness in darker portions.

Intuitive, User-friendly Functions

With One-Touch Connector for quick, easy connection and no need for white balance adjustment*2, setup is simplified, with the aim of streamlining workflow and accelerating procedure time. Touch-sensitive panel facilitates intuitive operation, while convenient functions like Pre-freeze and MyCV mode ensure user-friendly working environment. Downtime is reduced thanks to the use of LED bulbs that last years without needing replacement.

*1 Combination of EVIS EXERA III/EVIS LUCERA ELITE series light source and processor *2 Olympus 1100/1200/1500 series endoscopes only

Specifications		
Power Supply	Rated voltage	100-240 V AC; Within $\pm 10\%$
	Frequency	50/60 Hz; within ± 3 Hz
	Rated input	600 VA
Size	Dimensions (W x H x D)	370 x 198 x 488 mm; 398 x 218 x 580 mm (maximum)
	Weight	19.4 kg
Classification (Medical Electrical Equipment)	Type of protection against electric shock	Class I
	Degree of protection against electric shock of applied part	Depend on applied part. (The degree of protection against electric shock of this product is BF type if the mounting part to be connected to this product is BF type. However CF type is not subject to combination in this product.)
	Degree of protection against explosion	The video system center should be kept away from flammable gases.
	Analog signal output	VBS composite and Y/C; simultaneous outputs possible.
Eil. Nr. 2.9	Digital signal output	12G-SDI (SMPTE ST 2082), 3G-SDI (SMPTE424M), HD-SDI (SMPTE292M), SD-SDI (SMPTE259M)
	User settings	The function settings for up to 20 users can be stored.
Eil. Nr. 2.7	Color tone adjustment	Adjust the color tone of each endoscopic image for White light observation mode, NBI observation mode, and RDI observation mode. · Red adjustment : ± 8 steps · Blue adjustment : ± 8 steps · Chroma adjustment : ± 8 steps
	Automatic gain control (AGC)	The image can be electronically amplified when the light is inadequate due to the distal end of the endoscope being too far from the object.
	Contrast	· H (High) : Darkens the dark part and brightens the bright part. · L (Low) : Brightens the dark part and darkens the bright part.
Eil. Nr. 2.8	BAI-MAC	Brightness adjustment with maintenance of contrast
	Iris	The iris modes can be switched. · Auto : The brightness is adjusted based on the brightest part of the central part and the average brightness of the periphery part. · Peak : The brightness is adjusted based on the brightest part of the endoscopic image. · Average : The brightness is adjusted based on the average brightness of the endoscopic image.
	Image enhancement settings	Fine patterns or edges in the endoscopic images can be enhanced electrically to increase the image sharpness. · Enhancement type A : Emphasizes the pattern and contour of the endoscopic image. · Enhancement type B : Emphasizes the finer parts than structure emphasis type A.
Observation Eil. Nr. 2.13 Eil. Nr. 2.5 Eil. Nr. 2.4 Eil. Nr. 2.10	Switching the enhancement modes	The enhancement level can be selected from 3 levels (OFF, 1, 2, and 3)
	Image size selection	The size of the endoscopic image can be selected from 2 modes. (Except SDTV)
	Electric zoom	Switch between mode 1, mode 2, and mode 3.
	PIP/POP	Switch between PIP and POP.
	Aspect ratio	Switch between 16:9 and 4:3. (Except SDTV)
Eil. Nr. 2.6	Freeze	Freeze the endoscopic image.
	Pre-freeze	The image with the least blur is selected from the images captured in the set time period before freeze operation and displayed.
	Optical-digital observation	The optical-digital observation can be performed. The endoscope compatible with the optical-digital observation is required. · NBI observation : This observation mode uses the narrow band light. · RDI observation : This observation mode uses the red dichromatic lights. · AFI observation : This observation mode uses the blue light. · TXI observation : This observation mode enhances color, texture and brightness.
	Beginning and ending examination	Beginning and ending examination timing can be set interlock with the particular operation.
	Custom switch	Assign specific functions to the following buttons. · Remote switches (Up to 5) · Foot switches (Up to 2) · Keyboard custom key (Up to 4) · Touch panel custom button of basic functions screen (Up to 3) · Touch panel custom button of custom functions screen (Up to 10)
Documentation Eil. Nr. 2.11 Eil. Nr. 2.14.2	MyCV mode	Switch setting values of multiple functions at once.
	Remote control	The following peripheral device can be controlled (specified models only). · Portable memory · Video recorder · Color video printer · Image filing system · Server
	Patient information	The following data can be displayed on the monitor. · Patient ID · Patient name · Gender · Age · Date of birth · Comment
	Displaying the record state	The recording state of the following peripheral device can be displayed on the monitor. · Portable memory : Remaining capacity · Video recorder : Number of shots / Recording status · Color video printer : Number of shots · Image filing system : Number of shots
	Displaying the image information	The following data can be displayed on the monitor. · Image enhancement · Electric zoom ratio · Color mode · Focus · Observation mode
Memory Backup	Advanced registration of patient information	Up to 50 patient information can be registered. · Patient ID · Patient name · Gender · Age · Date of birth
	Recording format	Standard image quality: TIFF; Low image quality: JPEG
Memory Backup	Memorization of user settings	The settings are held in memory even after the video system center is turned OFF.
	White balance	The white balance that is once set is held in memory (only when using the compatible endoscope).

EVIS X1 VIDEO SYSTEM CENTER OLYMPUS CV-1500

Specifications, design, and accessories are subject to change without any notice or obligation on the part of the manufacturer.

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4K UHD LCD Monitor

OEV321UH

32-Inch LCD Monitor Optimized for Olympus Endoscopy Systems

OEV321UH

Eil.Nr. 2.4

High-Performance Display

This 32-inch LCD monitor with Ultra High Definition picture quality brings out the full potential of Olympus endoscopy systems. The 4K upscaling function accommodates existing HD imaging systems.

A.I.M.E.™ Technology

The Advanced Image Multiple Enhancer (A.I.M.E.™) produces sharp, vivid images by enhancing both structure and color to support a more detailed observation.

Versatile Signal Routing

The OEV321UH allows 4K/HD video signals to be routed via a single 12G-SDI output. Various input and output terminals offer impressive connectivity suited to user requirements.

Multiple Display Modes

Picture-in-Picture (PIP) and Picture-out-Picture (POP) display modes offer optimal support during procedures.


Convenient CLONE OUT Functionality

Users are able to duplicate the 4K/HD video signals as displayed on the screen, including PIP/POP to a second monitor or recording device.

User-Friendly Design


The front panel is easy to clean thanks to its flat surface including the control buttons. Due to downward-facing connectors, a detachable cable cover and a quick access window, cabling of this monitor is simplified.

Input/Output Port

OEV321UH Specifications			
Display	Display Size	31.5 in (800.8 mm)	Eil.Nr. 2.1
	Panel Technology	TFT Active Matrix LCD	
	Resolution / Aspect Ratio	3840 × 2160 / 16:9	Eil.Nr. 2.5
	Luminance	450 cd/m²	
	Contrast	1000:1	
	Color Space	BT.2020 / BT.709	
	Viewing Angle	178°/178°	
	Colors	1.07 billion	Eil.Nr.2.6
	Backlight	LED	
	Image Enhancement	A.I.M.E.™	
Functions	Multiple Image Display	PIP/POP	
	Flip Pattern	Rotation	
	4K	in 12G-SDI ×2*, Display Port ×1, HDMI ×1	
	out	12G-SDI ×2*	
Input Output	2K	in 3G-SDI ×1, DVI-D ×1	
	out	3G-SDI ×1	
	CLONE OUT / AUX IN	12G-SDI ×1* / Any port	
	Power Supply	AC 100V-240V, 50/60Hz, 1.7A-0.8A DC supplied by the optional AC adaptor AC-300MD	
General	Size	753.9 × 476.3 × 79.2 mm	
	Weight	11.8 kg	
	Remote	RS-232C connector D-sub 9 pin	

* When a 12G-SDI cable longer than MAJ-2429 (8.5 m) is used, please contact Olympus.

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MEDICAL DEVICE Medicinos priemonė

WM-NP3

Mobile Workstation

WM-NP3 Workstation

Features

The WM-NP3 has been designed and manufactured to enhance the user experience and add value to the Olympus Imaging Platform systems. Using a slim line profile and offering increased loading capacity while taking up less floor space and lowering running costs, the WM-NP3 will support Olympus systems without dominating the clinical environment, making it a valuable successor to WM-NP2.

Customization

Dedicated accessories are available to optimize system and procedural efficiency.

Electrical Safety Eil. Nr.3.3

The WM-NP3 is supplied with a separation transformer as standard. Output from the transformer is controlled using the central on/off switch, which allows all equipment to be powered up simultaneously. Eil. Nr.3.4

Inrush Current Control and Low-Power Standby

The WM-NP3 uses a level detector to detect mains distortion and protect the transformer from inrush current. The inrush current protection system is powered from a medically approved AC-DC power supply. To reduce the environmental impact of the WM-NP3, the low-power standby mode reduces standby power consumption giving greater efficiency and energy savings.

Easy to Move

Ergonomically designed handles and twin-wheeled castors allow the workstation to be moved into the ideal position in the clinical environment.

Cable Management

Optimized cable management has been achieved by locating the cable management modules within the workstation hoop, increasing the cable capacity and providing access to route and remove cables easily.

Imaging

The WM-NP3 workstation is offered with a choice of monitor arms that optimize the range of compatible monitors and enhance the procedural view with both swivel and tilt functionality. The MAJ-2216 has a weight range between 6.5 and 12 kg, and the MAJ-2217 has a weight range between 12 and 14 kg. Using a VESA fitting, both arms will accommodate a monitor up to 32".

WM-NP3 Workstation

Accessories

MAJ-2154
Consumable Storage Holder

Article number — K10030423

MAJ-2158
Suction Jar Holder

Article number — K10028139

MAJ-2159
Scope Pole Kit

Article number — K10028140

MAJ-2160
Irrigation Tube Holder

Article number — K10028141

MAJ-2165
Drawer Unit

Article number — K10030178

MAJ-2166
Sterile Water Holder

Article number — K10030036

MAJ-2167
Nurses Control Arm

Article number — K10030210

MAJ-2173
Side Shelf


Article number — K10030389

MAJ-2211
Side Handles

Article number — K10035108

MAJ-2216
LCD Monitor Arm 6.5 kg - 12 kg

Article number — K10035789

Standard Sets Include:

WM-NP3 Workstation plus items marked (●)

Article Number	Part Description	Region	Basic Set	GI Standard	URO Set
		EUROPE ROW	K10035360 K10035363	K10035365 K10035368	K10035370 K10035373
	Workstation Fitting Kit		●	●	●
	Transformer (market-specific)		●	●	●
K10028141	MAJ-2160 Irrigation Tube Holder		●	●	●
K10027573	MAJ-2146 Sliding Keyboard Tray			●	●
K10028140	MAJ-2159 Scope Pole Kit			●	
K10035789	MAJ-2216 LCD Monitor Arm 6.5-12 kg			●	
K10035108	MAJ-2211 Side Handles			●	●
K10009210	IEC Lead Set		●	●	●

Optional Compatible Accessories

Article Number	Part Description
K10021041	MAJ-1639 CO ₂ Cylinder Holder-140
K10016952	MAJ-1642 IV Pole
K10021042	MAJ-1650 CO ₂ Cylinder Holder-205
K10021043	MAJ-1653 CO ₂ Cylinder Holder — Double
K10021352	MAJ-1654 Equipotential Terminal Strip
K10021791	MAJ-1657 Keyboard Arm Side-Mounted
K10021795	MAJ-1661 Side-Mounted LCD Monitor Arm
K10021797	MAJ-1663 EUS Arm-Mount Kit
K10022056	MAJ-1665 Camera Head Holder
K10027573	MAJ-2146 Sliding Keyboard Tray
K10027575	MAJ-2149 Dual Monitor Arm

Specifications**WM-NP3**

Dimensions	1400 mm (H) × 675 mm (D) × 665 mm (W)
Weight	84 kg (unladen) including fitted transformer
Maximum Load Capacity	Maximum 20 kg on top tray
	Maximum 31 kg per shelf
	Maximum 35 kg on base panel
Castors	125 mm diameter
	2× braked conductive twin-wheeled castors
	2× braked twin-wheeled castors

Voltage**The WM-NP3 is supplied with a separation transformer based on one of the following voltage options:**

MAJ-2155	WM-T3 Transformer 100 V
MAJ-2156	WM-T3 Transformer 110-120 V
MAJ-2157	WM-T3 Transformer 220-240 V

Construction**Workstation**

Hoop	Extruded profiled aluminum
Shelves	Mild steel
Top Tray	Injection-molded ABS

Olympus reserves the right of errors, modification and changes of the service and/or product offerings.

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
SUCTION PUMP

KV-6

The KV-6 suction pump is the successor to the KV-5 pump. It is designed for endoscopic aspiration and for non-endoscopy applications. The unique flow change buttons allow a fast vacuum build-up with a single touch. It has a “whisper mode” for quiet operation and “turbo mode” for faster suction. The pump connects to an extensive fluid collection system.

Key Benefits

- **High Suction Performance**

The KV-6 suction pump has vacuum capabilities of up to **-95 kPa**. There are **three pre-set speed** settings (**40 l/min**, **50 l/min**, and **60 l/min**), which enables the operator the flexibility to control the rate of vacuum build-up at the touch of a button.

Eil. Nr. 4.2

Eil. Nr. 4.3

Eil. Nr. 4.3.1

Eil. Nr. 4.3.2

Eil. Nr. 4.3.3

- **Cleaning and Hygiene**

- Single-use collection system allows easy and safe disposal of clinical waste.
- Touch-sensitive buttons are easy to clean.

- **Compatibility**

Compatible with existing KV-5 accessories and Olympus workstations.

The KV-6 suction pump is ideal for use either mounted on a procedure trolley or as a stand-alone unit. The variable suction control allows adjustment to meet individual requirements. A high-performance pump and motor combination ensure that the unit is “whisper quiet” for any working environment.

The KV-6 may be adapted to suit individual requirements by the use of optional components, to allow a single-use or reusable suction system. The single-use liners are conveniently supplied with solidifier for safe and simple use and disposal.

Product Details

Compact Design for Convenient Use

The KV-6's compact and easy-to-use design means that it is suitable for use in office space or the operating theater.

Use with 1.5- or 2.5-Liter Jars and Liners

The KV-6 can be used with either 1.5- or 2.5-liter jars, enabling its use in a variety of procedures.

Simple and Easy to Operate with Clearly Labelled Controls

A clear control panel ensures that the KV-6 suction pump can be used safely and efficiently by all staff.


Suction Jar Mountable on Olympus Workstations

As the suction jar can be mounted on Olympus workstations, it is compatible with the Olympus range of devices and equipment for the office or operating theater.

Trolley-Mountable for Convenience

The KV-6 is ideal for use either mounted on a procedure trolley or as a stand-alone unit for maximum flexibility of use.

Overview of Accessories

Suction Tubing  K10026130 with 2 coupling pieces (1 x 90°), 7 x 12 mm, 60 cm 14921  K10026804 with 2 coupling pieces, 7 x 12 mm, 60 cm 14903	Patient Connecting Tubing  K7503432 50 pack 14919
Disposable Liners  K10021894 K10021895 1.5 liter 14908  K10021897 K10021898 2.5 liter 14909	Reusable Lids  K10026127 0.25 liter 14912  K7503760 1 & 2 liter 14911
Jars for Liners  K10021893 1.5 liter 14906  K10021896 2.5 liter 14907	Reusable Jars  K10026128 0.25 liter 14943  K7503762 1 liter  K7503763 2 liter Eil. Nr. 4.1
Filters  K10027321 Disposable bacteria/ odor filter 10 pack 14945  K10027322 Disposable bacteria filter 10 pack 14946	 14913  14921  14915 K10026130 K10026127 K10026128
Clamp Holder  K10002819 MAJ-191 14901  K10021793 MAJ-1659 K10001115 MAJ-190 14922  K7503556 WM-60 / D-60 14944	Two or More Jars Disposable collection system  K10021885 MAJ-1673 14900  K10021188 10 Pack 14899  K10021878 MAJ-1672 14902 Disposable/reusable collection system  K10021882 MAJ-1674 14920

Specifications

Classification (Electromedical Equipment)	Standards Compliance	Complies with EN ISO 10079-1 and EN IEC 60601-1 and UL 60601-1.		
	Electromagnetic Compatibility	This product complies with the requirements of EN IEC 60601-1-2 for emissions and immunity, and as such its operation is unlikely to be affected by, or cause interference with, equipment meeting appropriate EMC standards. As a precaution, equipment which may be sensitive to interference outside the limits specified by EN IEC 60601-1-2 should not be placed in close proximity to the KV-6.		
	Type of Protection against Electrical Shock	Classification according to EN IEC 60601-1 and UL 60601-1: Class 1 with Type CF applied part.		
	Degree of Protection against Electrical Shock	In accordance with EN IEC 60601-1 and UL 60601-1, the KV-6 is marked with the symbol to indicate the provision of an adequate degree of protection against electrical shock and that it has an applied part isolated from all other parts of the equipment.		
	Degree of Protection against Explosion	None: the Olympus KV-6 must NOT be used in the zone of risk of flammable anesthetic gases.		
	Mode of Operation	Continuous.		
Regulatory Status	European Economic Area (EEA)	CE: This mark on the product indicates compliance with Directive 93/42/EEC relating to medical devices, Class IIa.		
Dimensions and Weight	Dimensions	Height: 210 mm	Width: 305 mm	Depth: 375 mm
	Weight	9.3 kg		
Power Switch	Marking	The power switch is marked. A white LED in the center signalizes it is turned on.		
Power Requirements	Power Supply	100–240 V~	Fusing	2 x T1.6 AH
	Frequency	50/60 Hz	Power Rating	120 W
Environmental Conditions	Ambient Temperature	Operational:	+5°C to +40°C (+41°F to +104°F)	
		Storage:	–20°C to +50°C (–4°F to +122°F)	
	Relative Humidity	Operational:	30–75%	
		Storage:	20–95%	
	Atmospheric Pressure	Operational/ Storage:	70.0 to 106 kPa	
Fluid Ingress	In accordance with EN IEC 60601-1 and UL 60601-1, the product is marked IP21 to indicate that it is provided with an enclosure that prevents entry of such an amount of falling liquid as might interfere with its safe and satisfactory operation when correctly positioned.			
Resistance to Chemicals	The external surfaces of the Olympus KV-6 are resistant to: 2% aqueous neutral detergent, 70% ethyl alcohol, isopropyl alcohol, water.			
Pump Specification	Nominal Vacuum: 95 kPa		Nominal Free Air Flow Rate: 50 l/min, adjustable to 40 or 60 l/min.	
	In accordance with ISO 10079-1, the unit is marked with the words “HIGH VACUUM, HIGH FLOW” to indicate it has a vacuum of at least 60 kPa in 10 seconds and a free air flow rate of 20 l/min or greater.			
Thermal Cut-Out	The temperature is observed by a sensor in the housing. If the temperature in the housing exceeds 80°C (176°F), the software stops the motor and turns on the yellow LED. It is possible to restart the motor as soon as the temperature in the housing drops below 70°C (158°F).			
Single-Use Suction Liner	To prevent fluid from being drawn up to the pump the single-use suction liner contains a mechanism to shut off the flow when the liner is overfilled. This mechanism is activated when contact is made with aspirated fluid.			
Suction Jar	The suction jars can be autoclaved up to 135°C or disinfected. The jar is impact-resistant.			
Suction Jar Lid	The suction jar lid clips and where applicable the float mechanism are autoclavable up to 135°C. The inlet connection (from the endoscope) is marked with the word PATIENT, and the outlet connection (to the pump) is marked with the word VACUUM.			
Tubing	Use only tubing specified in the KV-6 Instructions for Use, otherwise equipment malfunction may result.			
Vacuum Gauge	Range: 0 to 100 kPa	Accuracy: ±2.5%		

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E0428423 · 2.000 · 11/15 · PR · ED

OLYMPUS **FLUSHING** **PUMP** OFP-2

With its powerful but controllable peristaltic pump the OFP-2 will irrigate fluid down either the instrumentation or auxiliary water channel of Olympus Gastrointestinal EVIS ACTERA, EXERA, LUCERA, SMARTAGE and EUS endoscopes.

The OFP-2 facilitates two main functions:

- The washing of gastric and colonic mucosa resulting in improved visualisation, diagnosis and therapy during endoscopic procedures
- To aid the use of trans-endoscopic ultrasound probes by providing rapid filling of the organ to be examined

Main Unit:

- Designed to prevent over pressurisation of Olympus endoscopes
- Microprocessor control ensures accurate and precise, flow rate delivery
- Ergonomically designed touch buttons for simple use
- Bright LED display for easy verification of selected settings

Practicality:

- Use of the MAJ-1606 irrigation adaptor allows either irrigation or a combination of irrigation and the passage of an EndoTherapy device through the endoscopes instrumentation channel

Eil. Nr. 5.2 Automatinio išjungimo funkcija

Safety Features:

- Pump Head automatically cuts out if opened during operation preventing accidental user injury
- Standby mode allows quick and safe exchange of water channel tubes without the need to operate the mains switch
- A 20 second automatic cut-off timer allows controlled irrigation

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

Water Container:

- Holds two litres with markings to indicate fluid levels
- A specially designed lid guides the water tubing to optimise fluid usage
- The water container and lid are autoclavable for effective sterilisation
- Washer disinfector compatible

Size and functionality:

- Improved functionality allows use with either the instrumentation or auxiliary water channel, including EUS ultrasound endoscopes for the first time
- Compact design for essential space saving on workstations
- Operates via either remote controlled* scope operation or foot switch

Eil. Nr. 5.3 Nuotolinio valdymo funkcija

* requires the use of MAJ-920 remote control cable

Consumables:

- The MAJ-1607 and MAJ-1608 channel water tubes and the MAJ-1606 adaptor are supplied sterile for immediate use
- The MAJ-1608 auxiliary channel water tube is designed for one single days use
- Instrument channel adaptor and MAJ-1607 water tube are single patient use items

MAJ-1606

Single patient use


Supplied sterile

Instrument channel adaptor

Pack 100: Part number - K10007072

Pack 10: Part number - K10016091

MAJ-1607

Single patient use

Supplied sterile

Instrument channel water tube

Pack 50: Part number - K10001146

Pack 10: Part number - K10016136

MAJ-1608

Single days use

Supplied sterile

Auxiliary channel water tube

Pack 50: Part number - K10001147

Pack 10: Part number - K10016135

Accessories:

MAJ-1603

Pack of three

Water container two litre

Part number - K10007071

MAJ-855

Reusable

Supplied non sterile

Auxiliary water tube

Part number - MAJ-855

MAJ-920

Remote control cable

Part number - MAJ-920

Specification:

Dimensions: 200(W) x 173(H) x 385(D) mm

Weight: 4 kg (with empty water container)

Maximum Flow Rate: Instrument channel: 750 ml/min

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Safety: Complies with EN/UL/IEC60601-1

Regulatory: CE marked and ETL Listed

Power

Requirements: Power Supply: 100-240 V ~

Frequency: 50/60 Hz,

Fusing: 2 x T2.5 AH 250 V

Power rating: 100 VA

Olympus reserves the right to alter the specification without notice.

OLYMPUS

OLYMPUS KEYMED GROUP
KeyMed House, Stock Road, Southend-on-Sea,
Essex SS2 5QH, United Kingdom
Telephone: +44 (0)1702 616333 Fax: +44 (0)1702 465677
e-mail: info@olympus.co.uk www.olympus.co.uk

OLYMPUS MEDICAL SYSTEMS CORP.
Shinjuku Monolith, 3-1, Nishi Shinjuku 2-chome, Shinjuku-ku, Tokyo 163-0914,
OLYMPUS EUROPA HOLDING GMBH
(Premises/Goods delivery) Wendenstraße 14-18, 20097 Hamburg, Germany
(Letters) Postfach 10 49 08, D-20034 Hamburg, Germany

OLYMPUS AMERICA INC.
3500 Corporate Parkway, PO Box 610, Center Valley, PA 18034-0610, USA

OLYMPUS SINGAPORE PTE LTD.
491B River Valley Road #12-01/04, Valley Point Office Tower, Singapore 248373, Singapore

UCR CO₂ Insufflator

Exploiting CO₂'s Ability to Minimize Luminal Distension, the UCR Features Easy Operation, a Selectable Flow Rate, and a Compact Design


Main Features

- About 150 times faster than conventional air, carbon dioxide's rapid absorption properties keep abdominal distension and related pain to a minimum during the procedure and speed up recovery afterwards.
- Easy, reliable one-button start/stop operation, pressure display, and timer function to automate CO₂ insufflation shutoff.
- Compact design allows the UCR to fit easily on an endoscopic workstation.
- Three controlled flow rate settings are available, by using the appropriate corresponding tubing. Reusable and disposable tubes are optional.

Eil. Nr. 6.4

Rapid absorption of CO₂ gas

Carbon dioxide is absorbed by human tissue about 150 times faster than conventional air.


Specifications		
Power Requirements	Voltage (AC)	100-120V (USA), 100-240V (Other areas)
	Frequency	50/60Hz
	Input	40VA
	Voltage fluctuation	Within ±10%
Dimensions	125 (W) x 300 (D) x 150 (H) mm (housing dimensions)	
Weight	4.9 kg	
Applicable Gas	CO ₂ gas for medical use (Connectable to CO ₂ cylinder or medical gas pipeline)	
Air Feeding Pressure	Maximum pressure feed	45 kPa

Eil. Nr. 6.3

System Chart

Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.

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HF Generator ESG-150

SIMPLY SMART ...

ESG-150 Is User-Friendly and Offers Outstanding Effectiveness

Olympus understands the specific demands of modern electrosurgery, so it knows how important it is to have a safe, dedicated energy source to maintain the precision physicians require in their work.

... At the Simple Touch of a Button

A simple, intuitive touch-screen interface means users can easily set up the system to meet their personal needs in just a few steps.

Eil. Nr. 7.3

... Supports User Level

The user can choose between two user levels: **advanced and basic**. The advanced level is best for users that want to precisely adjust their HF settings and have maximum control. The basic level is ideal for users that prefer ease of use and the ability to quickly adapt the HF settings required for the procedure.

... Offers You Guidance


The smart assistance messaging system provides simple and easily understandable feedback that empowers users to confidently operate and troubleshoot the ESG-150.

... Talks Your Language

A choice of 29 display languages means users can understand and operate the system with ease.

... Provides Advanced Safety

The power range within the basic user level is limited to a reasonable amount to avoid unintended excessive power output. This results in additional safety for the user and the patient.

... Built On Experience

Fast Spark Monitor – measures spark intensity during cutting procedures, automatically adjusting power output to achieve reproducible tissue effects and enhanced protection of EndoTherapy instruments.

High-Power Cut Support – immediate spark ignition enables cutting without initial delays, which means superior cutting performance.

Contact Quality Monitor – analyses contact resistance of the neutral electrode, reducing the risk of thermal injuries to the patient.

Leakage-Protection Sensor – measures and controls leakage current, reducing the risk of unintended harm to the patient or user.

Monopolar HF Modes

Cut Modes

PureCut

Continuous cutting mode with low coagulation effect

PulseCut Slow/Fast

Intermittent cutting mode which is commonly used for various endoscopic procedures

Coagulation Modes

SoftCoag

Slow and deep coagulation – optimal for coagulation of thick blood vessels

ForcedCoag

Fast and effective coagulation

High Power Cut Support (HPCS) – Immediate Cutting Fast Spark Monitor (FSM) – Optimized Cutting Power

Versatile Coagulation Technology

Technical data	
Olympus ESG-150	Specifications
High frequency	356 kHz \pm 20%
Type, protection class	CF, Class I
Power supply	AC 100 V to 120 V/AC – 220 V to 240 V, 50/60 Hz, 650 VA
Size (W x D x H)	295 x 400 x 115 mm
Weight, without foot switch	6.6 kg

Monopolar output		
Sockets	3-pin plugs: Pin Ø 4 mm (international/Valleylab standard)	
	1-pin plugs: Pin Ø 8 mm (Bovie standard)	
	1-pin plugs: Pin Ø 4 mm	
Split neutral electrode	10 mm 2-pins	
PureCut (Effect 1–3)	120 W @ 500 Ohms	Eil. Nr. 7.4.1
PulseCut slow/fast (Effect 1–5)	120 W @ 500 Ohms	Eil. Nr. 7.4.2
SoftCoag (Effect 1–5)	120 W @ 500 Ohms	Eil. Nr. 7.4.3
ForcedCoag (Effect 1–5)	120 W @ 500 Ohms	Eil. Nr. 7.4.4
Bipolar output		
Sockets	28.8 mm 2-pins and Ø 4/8 mm coaxial	
BipolarCut (Effect 1–3)	100 W @ 500 Ohms	Eil. Nr. 7.4.5
BiSoftCoag (Effect 1–3)	120 W @ 100 Ohms	Eil. Nr. 7.4.6

Electrosurgical generator
(ESG-150)

Ordering Information				
Order No.	Description	Patient Weight	Contact Area	Unit
WA90150W	Electrosurgical generator "ESG-150"			
E0420539S	ESG-150 set incl. WA99516C			
WA99512C	Neutral Electrode, for newborns, with cable	< 5 kg	40 cm ²	15 pcs
WA99513A	Neutral Electrode, for children	5 – 15 kg	70 cm ²	100 pcs
WA99513C	Neutral Electrode, for children, with cable	5 – 15 kg	71 cm ²	15 pcs
WA99514A	Eil. Nr. 7.6.1 Neutral Electrode, for adults	> 15 kg	140 cm ²	100 pcs
WA99514C	Neutral Electrode, for adults, with cable	> 15 kg	140 cm ²	15 pcs
WA99516A	Neutral Electrode, universal size	> 5 kg	110 cm ²	100 pcs
WA99516C	Neutral Electrode, universal size, with cable	> 5 kg	110 cm ²	15 pcs
WA99510A	Eil. Nr. 7.6.2 Cable, for neutral electrode			
A02915A	Eil. Nr. 7.6.3 HF-cable "MAJ-860", monopolar, 2 m, UES-30, Erbe Intl. and Valleylab (new) HF-unit to SD handle			

As medical knowledge is constantly growing, technical modifications or changes of the product design, product specifications and accessories may be required.

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