HDLC3

DTL cholesterolis, 3-ia karta

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cobas®

Substratai

REF	CONTENT	Analizatorius (-iai), su kuriuo (-iais) galima naudoti cobas c pakuotę (-es)
04399803 190	HDL-Cholesterol plus 3rd generation (200 tyrimų)	Sistemos-ID 07 6833 2 COBAS INTEGRA 400 plus COBAS INTEGRA 800
12172623 122	C.f.a.s. Lipids (3 × 1 mL)	Sistemos-ID 07 6570 8
12172623 160	C.f.a.s. Lipids (3 × 1 mL, skirta JAV)	Sistemos-ID 07 6570 8
10781827 122	Precinorm L (4 × 3 mL)	Sistemos-ID 07 9026 5
11778552 122	Precipath HDL/LDL-C (4 × 3 mL)	Sistemos-ID 07 9028 1
05117003 190	PreciControl ClinChem Multi 1 (20 × 5 mL)	Sistemos-ID 07 7469 3
05947626 190	PreciControl ClinChem Multi 1 (4 × 5 mL)	Sistemos-ID 07 7469 3
05947626 160	PreciControl ClinChem Multi 1 (4 × 5 mL, skirta JAV)	Sistemos-ID 07 7469 3
05117216 190	PreciControl ClinChem Multi 2 (20 × 5 mL)	Sistemos-ID 07 7470 7
05947774 190	PreciControl ClinChem Multi 2 (4 × 5 mL)	Sistemos-ID 07 7470 7
05947774 160	PreciControl ClinChem Multi 2 (4 × 5 mL, skirta JAV)	Sistemos-ID 07 7470 7
20756350 322	NaCl Diluent 9 % (6 × 22 mL)	Sistemos-ID 07 5635 0

Lietuvių

Sistemos informacija

Tyrimas HDLC3, tyrimo ID 0-331 COBAS INTEGRA 400 plus sistemose; tyrimo ID 0-333 COBAS INTEGRA 800 sistemose

Paskirtis

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Nuorodos

- Dominiczak M, McNamara J. The system of Cardiovascular prevention. 103-125; Nauk M, Wiebe D, Warnick G. Measurement of High-Density-Lipoprotein Cholesterol. 221-244. In: Handbook of Lipoprotein Testing (eds. Rifai, Warnick and Dominiczak), 2nd edition.
- Linsel-Nitschke P, Tall AR. HDL as a target in the treatment of atherosclerotic cardiovascular disease. Nature Reviews 2005;4:193-205.
- Ng DS. Treating low HDL - From bench to bedside. Clinical Biochemistry 2004;37:649-659.
- Sugiuchi H, Uji Y, Okabe H, et al. Direct Measurement of High-Density Lipoprotein Cholesterol in Serum with Polyethylene Glycol-Modified Enzymes and Sulfated α -Cyclodextrin. Clin Chem 1995;41:717-723.
- Matsuzaki Y, Kawaguchi E, Morita Y, et al. Evaluation of Two Kinds of Reagents for Direct Determination of HDL-Cholesterol. J Anal Bio-Sc 1996;19:419-427.
- Nauck M, März W, Jarausch J, et al. Multicenter evaluation of a homogeneous assay for HDL-cholesterol without sample pretreatment. Clin Chem 1997;43:1622-1629.
- Zawta B, Klüber J. Brochure "Wissenswertes zu Apolipoproteinen". Fragen/Antworten (Boehringer Mannheim 1991). In: Henry JB, ed. Clinical Diagnosis and Management by Laboratory Methods. 17th ed. Philadelphia: WB Saunders 1984;251-282.
- AVP Fettstoffwechselstörungen. Therapieempfehlungen 1, 1st ed. 1996;2-16.
- Hatch FT, Lees RS. Practical methods for plasma lipoprotein analysis. Adv Lipid Res 1968;6:1-68.
- Narayan KA, Kummerow FA. Disk electrophoresis of human serum lipoprotein. Nature 1965;205:246-248.
- Okazaki M, Shiraishi K, Ohno Y, et al. Heterogeneity of human high density lipoproteins on high performance liquid chromatography. J Biochem 1982;92:517-524.
- Burstein M, Scholnick HR, Morfix R. Rapid method for the isolation of lipoproteins from human serum by precipitation with polyanions. J Lipid Res 1970;11:583-595.
- Musto J, Lawlor JF. HDL-cholesterol: online separation and analysis utilizing an automated chemistry analyzer [Abstract]. Clin Chem 1993;39:1125.
- Kakuyama T, Kimura S, Hashiguchi Y. Fully automated determination of HDL-cholesterol from human serum with Hitachi 911 [Abstract]. Clin Chem 1994;40:1104.
- Harris N, Galpichian V, Rifai N. Three routine methods for measuring high-density lipoprotein cholesterol compared with the Reference method. Clin Chem 1996;42:738-743.

- Sugiuchi H. History of development and technical details of the homogeneous assay for HDL and LDL cholesterol. The Fats of Life 2005;IX No. 1:4-11.
- Cohn JS, McNamara JR, Schaefer EJ. Lipoprotein Cholesterol Concentrations in the Plasma of Human Subjects as Measured in the Fed and Fasted States. Clin Chem 1988;34:2456-2459.
- Pisani T, Gebbski CP, Leary ET, et al. Accurate Direct Determination of Low-density Lipoprotein Cholesterol Using an Immunoseparation Reagent and Enzymatic Cholesterol Assay. Arch Pathol Lab Med 1995 Dec;119(12):1127-1135.
- Data on file at Roche Diagnostics.
- Kimberly M, Leary E, Cole T, et al. Selection, Validation, Standardization and Performance of a Designated Comparison Method for HDL-Cholesterol for Use in the Cholesterol Reference Method Laboratory Network. Clin Chem 1999;45:1803-1812.
- Tietz NW. Textbook of Clinical Chemistry, 3rd Edition 1999;842-843.
- Cooper GR, Myers GL, Smith SJ, et al. Standardization of Lipid, Lipoprotein, and Apolipoprotein Measurements. Clin Chem 1988;34/8B:B99.
- Kadri N, Douville P, Lachance P. Letter to editor. Clin Chem 2002;48:964.
- Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-475.
- Breuer J. Report on the Symposium "Drug effects in Clinical Chemistry Methods". Eur J Clin Chem Clin Biochem 1996;34:385-386.
- Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385.
- Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med 2007;45(9):1240-1243.
- Thomas L, ed. Labor und Diagnose, 4th ed. Marburg: Die Medizinische Verlagsgesellschaft 1992;208.
- Assmann G. At what levels of total low- or high-density lipoprotein cholesterol should diet/drug therapy be initiated? European guidelines. Amer J Cardiol 1990;65:11F.
- Assmann G, Schriewer H, Schmitz G, et al. Quantification of high-density-lipoprotein cholesterol by precipitation with phosphotungstic acid/MgCl₂. Clin Chem 1983;29(12):2026-2030.
- Third Report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). NIH Publication No 01-3670; May 2001.
- Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.

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Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim
www.roche.com

JAV platina:

Roche Diagnostics, Indianapolis, IN

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cobas[®]
Substratai

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© 2014, Roche Diagnostics