
BiCart® 1150g cartridge

Essential component of a system designed for ultrapure dialysis fluid to avoid treatment-related inflammation

The BiCart dry cartridge is central to the dialysis fluid delivery system:

- A hygienic method of bicarbonate supply from a dry-concentrate cartridge
- Minimizing the risk of contamination and bacterial growth
- Lighter and smaller container, removing strain associated with lifting heavy canisters
- Cleaner, safer disposal after treatment with cap closure
- Less storage space required

Tougher hygienic and ergonomic demands required an upgraded bicarbonate supply. With its introduction, the BiCart cartridge set a new global standard in bicarbonate dialysis.

BiCart® 1150g—specifications

General description

The BiCart cartridge is a polypropylene container sodium bicarbonate which enables on-line production of liquid bicarbonate concentrate.

When attached to a special holder, water passes through the BiCart cartridge, thus producing a saturated bicarbonate solution, ready for use.

This concentrated solution is then proportioned in the dialysis machine.

Product code: 113 779 (EU-Pallet)
113 882 (container-pallet)

Units per box: 6

Specifications

Contents: Each BiCart cartridge contains min. 1150 g of sodium bicarbonate powder.

Capacity: Together with a suitable acid concentrate each BiCart cartridge will yield enough sodium bicarbonate solution to produce 300 l of dialysis fluid with a bicarbonate concentration of 34 mmol/l.

Duration: Dialysis fluid flow rate
500 ml/min 10 h
600 ml/min 8 h 20 min
700 ml/min 7 h 10 min

Storing: Store below +40°C.

Shelf life: 24 months from date of manufacture

Quality control

Quality control as well as manufacturing procedures for the BiCart cartridge follow the current Good Manufacturing Practice (GMP) for pharmaceuticals. As well as European and US legislation for medical devices.

The information may be subject to change without further notice. For further information and operating instructions, please refer to the operator's manual.

CE 0086 The products meet the applicable provisions of Annex I (Essential Requirements) and Annex II (Full quality assurance system) of the Council Directive 93/42/EEC of 14 June 1993 concerning medical devices.

Gambro and BiCart are trademarks belonging to the Gambro Group.

The inspection of components, the process control and finished product testing are rigorous. All documentation and results from tests and inspections are checked and reviewed.

Requirements on acid concentrate

The acid concentrate must be designed for use with pure sodium bicarbonate concentrate and have a dilution ratio of 1:35 (1+34) or 1:45 (1+44). As an example, Gambro SoftPac™ can be used. When mixed 1:45 (1+44) this will yield (mmol/l):

Sodium	103.00
Potassium	2.00
Calcium	1.75
Magnesium	0.50
Chloride	109.50
Acetic acid	3.00

Together with the sodium bicarbonate concentrate from the BiCart cartridge the final dialysis fluid will have a sodium concentration of 140 mmol/l and a bicarbonate concentration of 34 mmol/l, at a normal machine setting.

Water purity

Water for diluting concentrated hemodialysis solutions should comply with local regulations and if no such regulations are available follow ISO 13959.

Top filter (A)

A filter, placed in the top of the BiCart 1150g cartridge, prevents powder from going backwards up into the dialysis machine.

Parking lot for cap (B)

Two caps come with the BiCart 1150g cartridge. The caps protect the connectors during transportation and by putting them back on the connectors after a treatment session, spillage and dripping can be minimized. During the treatment session the caps can be placed in the two parking lots.

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„Gambro Lundia AB“