


Vascular Intervention // **Coronary**  
Drug-Eluting Stent System

# Orsiro



Clinically proven

---



Highly deliverable



---



Ultrathin 60  $\mu$ m struts



**BIOTRONIK**  
excellence for life




**3.8%**  
restenosis rate  
at five years

## The new benchmark for DES

### BIOFLOW-V 12-month clinical outcomes compared to Xience

In a post-hoc analysis of pooled patient-level data from three RCTs, Orsiro achieved a 96.9% probability of superiority\* on TLF rate versus Xience.<sup>10</sup>


### BIOFLOW-V / -IV / -II Bayesian Population (n=2,208)



## Proven long term clinical outcomes

### All stents implanted from 2007 until January 11, 2017 unadjusted (SCAAR)<sup>11,12</sup>

Orsiro showed a lower restenosis rate than all DES out to five years.



# Orsiro

Clinically proven, highly deliverable with ultrathin 60 µm\* struts.

## Clinically proven

### Extensive clinical program\*\*

- >32,500 patients enrolled
- >50,500 patients planned in total
- >44 studies ongoing
- >55 studies planned in total

\*\*status as of Feb 2017

### Outstanding clinical results even in challenging subgroups

Orsiro has demonstrated consistently low target lesion failure (TLF) in all-comers trials compared to major modern drug-eluting stents (DES).


#### BIO-RESORT<sup>1,2</sup> (n=3,514 patients)



#### BIO-SCIENCE<sup>3,4</sup> (n=2,121 patients)



#### SORT-OUT VII<sup>5,6</sup> (n=2,314 patients)



STEMI  
5.4%  
10.8% Xience  
TLF<sup>4</sup> 24 months  
BIO-SCIENCE<sup>7</sup>

Diabetics  
0.0%  
ST<sup>8</sup> 60 months  
BIOFLOW-II<sup>9</sup>

Small vessels  
0.0%  
ST<sup>8</sup> 60 months  
BIOFLOW-II<sup>9</sup>



ST - Stent Thrombosis

\*ø 2.25 - 3.0 mm

>50,500  
patients planned  
in total

0.0%  
ST<sup>8</sup> at 5 years  
BIOFLOW-II<sup>9</sup>





**57%**  
better push


**13%**  
lower crossing  
profile

## Highly deliverable

Designed for challenging cases, the Orsiro stent system provides better push and easier cross with a lower crossing profile.


### Better push

Transmitting up to 57%<sup>13</sup> more force from hub to tip.<sup>14</sup>




### Easier cross

Up to 68%<sup>15,16</sup> less force<sup>15,16</sup> needed to successfully cross demanding anatomies.



### Lower crossing profile

Improved acute performance - up to 13% lower crossing profile.<sup>15</sup>



Strut thickness  
in perspective<sup>17</sup>

Orsiro  
BIOTRONIK  
CoCr-SES  
60 µm\*

Synergy  
Boston Scientific  
PtCr-EES  
74 µm

Ultimaster  
Terumo  
CoCr-SES  
80 µm

Resolute Onyx  
Medtronic  
CoNi-ZES  
81 µm

Xience Family  
Abbott  
CoCr-EES  
81 µm

Promus  
Boston Scientific  
PtCr-EES  
81 µm

BioMatrix  
Biosensors  
316L-BES  
120 µm

Ultrathin 60 µm struts

Thinner struts make  
the difference

Thinner struts create:

- Less disrupted flow<sup>18</sup>
- Less arterial injury<sup>18</sup>


Which leads to:

- Improved re-endothelialization<sup>18</sup>
- Reduced risk of restenosis  
and thrombosis<sup>18</sup>



The thinner the better, as long as the radial  
force can be maintained<sup>18</sup>

Up to 15% more radial strength<sup>19,20</sup> for stronger  
scaffolding once implanted.



15%  
more radial  
strength

\* ø 2.25 – 3.0 mm

# Orsiro

Vascular Intervention  
Coronary



Indicated for discrete de novo stenotic lesions and in-stent restenotic lesions.\*

| Technical Data             |  | Stent  |
|----------------------------|--|--|
| Stent material             |  | Cobalt chromium, L-605   |
| Passive coating            |  | proBIO (Amorphous Silicon Carbide)                                 |
| Active coating             |  | BIOLute bioabsorbable Poly-L-Lactide (PLLA) eluting a limus drug   |
| Drug dose                  |  | 1.4 µg/mm <sup>2</sup>   |
| Strut thickness            |  | ø 2.25 - 3.0 mm: 60 µm (0.0024"); ø 3.50 - 4.0 mm: 80 µm (0.0031") |
| Delivery system            |  |  |
| Catheter type              |  | Rapid exchange   |
| Recommended guide catheter |  | 5F (min. I.D. 0.056")  |
| Lesion entry profile       |  | 0.017"   |
| Guide wire diameter        |  | 0.014"   |
| Usable catheter length     |  | 140 cm   |
| Balloon material           |  | Semi crystalline polymer material                                  |
| Coating (distal shaft)     |  | Hydrophilic coating  |
| Marker bands               |  | Two swaged platinum-iridium markers                                |
| Proximal shaft diameter    |  | 2.0F   |
| Distal shaft diameter      |  | 2.6F: ø 2.25 - 3.5 mm; 2.8F: ø 4.0 mm                              |
| Nominal pressure (NP)      |  | 8 atm  |
| Rated burst pressure (RBP) |  | 16 atm   |

| Compliance Chart           |        | Balloon diameter x length (mm) |               |               |               |               |               |
|----------------------------|--------|--------------------------------|---------------|---------------|---------------|---------------|---------------|
|                            |        | ø 2.25 x 9-40                  | ø 2.50 x 9-40 | ø 2.75 x 9-40 | ø 3.00 x 9-40 | ø 3.50 x 9-40 | ø 4.00 x 9-40 |
| Nominal Pressure (NP)      | atm**  | 8                              | 8             | 8             | 8             | 8             | 8             |
|                            | ø (mm) | 2.25                           | 2.50          | 2.75          | 3.00          | 3.50          | 4.00          |
| Rated Burst Pressure (RBP) | atm**  | 16                             | 16            | 16            | 16            | 16            | 16            |
|                            | ø (mm) | 2.50                           | 2.77          | 3.05          | 3.33          | 3.88          | 4.44          |

\*\*1 atm = 1.013 bar

| Ordering Information | Stent ø (mm) | Catheter length 140 cm<br>Stent length (mm) |        |        |        |        |        |        |        |        |
|----------------------|--------------|---|--------|--------|--------|--------|--------|--------|--------|--------|
|                      |              | 9   | 13     | 15     | 18     | 22     | 26     | 30     | 35     | 40     |
|                      | 2.25         | 364469                                      | 364475 | 364481 | 364487 | 364499 | 364505 | 364511 | 391234 | 391238 |
|                      | 2.50         | 364470                                      | 364476 | 364482 | 364488 | 364500 | 364506 | 364512 | 391235 | 391239 |
|                      | 2.75         | 364471                                      | 364477 | 364483 | 364489 | 364501 | 364507 | 364513 | 391236 | 391240 |
|                      | 3.00         | 364472                                      | 364478 | 364484 | 364490 | 364502 | 364508 | 364514 | 391237 | 391241 |
|                      | 3.50         | 364473                                      | 364479 | 364485 | 364491 | 364503 | 364509 | 364515 | 391018 | 391020 |
|                      | 4.00         | 364474                                      | 364480 | 364486 | 364492 | 364504 | 364510 | 364516 | 391019 | 391021 |

1. von Birgelen et al. Very thin strut biodegradable polymer everolimus-eluting stents versus durable polymer zotarolimus-eluting stents in all-comers with coronary artery disease (BIO-RESORT): a three-arm, randomised, non-inferiority trial. The Lancet 2016. 10.1016/S0140-6736(16)31920-1 and presentation at TCT 2016; 2. TLF as a composite of cardiac death, target vessel-related myocardial infarction, or clinically indicated target lesion revascularization; 3. Pilgrim et al. Ultrathin strut biodegradable polymer sirolimus-eluting stent versus durable polymer everolimus-eluting stent for percutaneous coronary revascularization (BIO-SCIENCE): a randomised, single-blind, non-inferiority trial. The Lancet 2014.10.1016/S0140-6736(14)61038-2; 4. TLF as a composite of cardiac death, target vessel myocardial infarction, and clinically indicated target lesion revascularization; 5. Jensen et al. Randomized comparison of a sirolimus-eluting Orsiro stent with a biolimus-eluting Nobori stent in patients treated with percutaneous coronary intervention: Rationale and study design of the Scandinavian Organization of Randomized Trials with Clinical Outcome VII trial. 10.1016/j.ahj.2015.05.009; 6. Target Lesion Failure as a composite of cardiac death, myocardial infarction (not related to other than index lesion), or target lesion revascularization; 7. Piccolo R. Biodegradable polymer sirolimus-eluting stents vs. Durable polymer Everolimus-eluting stents in patients with STEMI: Two-year follow-up of the BIOSCIENCE oral presentation, EuroPCR 2016; 8. Definite or probable stent thrombosis per ARC definition; 9. Preliminary analysis based on non locked data – Ton Slagboom, poster presentation, presented at TCT, November 2016; 10. Kandzari et al. Ultrathin Bioresorbable Polymer Sirolimus-Eluting Stents versus thin durable Polymer Everolimus-eluting stents in patients Undergoing Coronary Revascularization (BIOFLOW-V): a randomized trial, The Lancet 2017; 11. Adapted from SCAAR data (January 11, 2017) <http://www.ucr.uu.se/swedeheart/99-scaar/forskning-scaar>; 12. Compared to other DES included in SCAAR at five years; 13. Compared to Resolute Onyx; 14. The stent system is advanced through a model, to a point of blockage (simulating a total occlusion). The force at the proximal hub and the blockage is measured. Pushability is the force transmitted along the length of the catheter. IIB(P)31/2015 – IIB(P)85/2014-2; 15[16]. Compared to Xience Xpedition; 16. The stent system is advanced through a stenosis model. Crossability is the mean resistance (mean force) registered by the stenosis during the complete passage of the stent delivery system. IIB(P)31/2015 – IIB(P)85/2014-2; 17. Stefanini GG, Taniwaki M, Windecker S. Coronary stents: novel development, Heart doi:10.1136/heartjnl-2012-303522; 18. Foin et al. Impact of stent strut design in metallic stents and biodegradable scaffolds. Int J Cardiol.2014 Dec 20;177(3):800-8; 19. Compared to Xience Expedition; 20. Expanded 3.0 mm diameter stents are radially compressed (15% of ø) along full length. The force required to compress the stent is radial strength. BIOTRONIK data on file.

Synergy and Promus are registered trademarks of Boston Scientific/Resolute, Integrity, Resolute Integrity and Resolute Onyx are registered trademarks of Medtronic/Xience, Xience Prime and Xience Xpedition are registered trademarks of Abbott Cardiovascular Systems. Nobori and Ultimaster are registered trademarks of Terumo/BioMatrix is a registered trademark of Biosensors.

\*Indication as per IFU.

BIOTRONIK AG  
Ackerstrasse 6  
8180 Bülach, Switzerland  
Tel +41 (0) 44 8645111  
Fax +41 (0) 44 8645005

info.vi@biotronik.com  
www.biotronik.com

© 2018 BIOTRONIK AG – All rights reserved.  
Specifications are subject to modification, revision and improvement.

 **BIOTRONIK**  
excellence for life