MINI VIDAS[®]

Simply reliable

PIONEERING DIAGNOSTICS

MINI VIDAS®

Simply reliable

With **MINI VIDAS®** your lab can carry out tests simply and quickly, in a self-contained, space-saving unit. Its user-friendly Load & Go concept means your teams can process up to 36 tests an hour. It also offers excellent reliability, with a MTBF (1) of over 1,100 days.

An improved Design for Increased User Comfort

DESIGN


- Streamline your workflow by optimizing test loading
- Increase user-friendliness for data entry and reporting patient results
- Compact solution with integrated computer & printer

11 punktaz

RELIABLE

- System robustness (MTBF ⁽¹⁾ > 1,100 days)
- Single dose concept, ready-to-use reagents
- Automated barcode identification
- Just Load & Go

1 punktaz

* SPR: Solid Phase Receptacle
(1) MTBF = Mean Time Between Failure
(2) Pending list of assay registered in a specific country

FLEXIBLE

- Two independent 6-test sections
- Processing of different parameters simultaneously
- Packaging suite to different test volumes (30 or 60 tests)
- 100 parameters available ⁽²⁾

5 punktaz

7 punktaz


	MINI VIDAS®
Weight	40 kg (88 lbs)
Voltage	100-240 VAC
Electric Consumption	1.5-0.8 A
Frequency	50-60 Hz
Power	150 Watts

Discover the VIDAS community and become a member by registering on www.myvidas.com

Vidas

Assay Solutions

VIDAS® - **PROTOCOLS - UNITS**

VIDAS® - **PARAMETERS**

Compatible assays

Parameter Code	Sample volume	Parameter Code	Sample volume	Parameter Code	Sample volume	Parameter Code	Sample volume	Parameter Code	Sample volume
TSH	200 µL	VWF	100 µL	HIV5	200 µL	HBS	150 µL	LYM	100 µL
T4	200 µL	PC	100 µL	HIV6	200 µL	AHBS	200 µL	LYG	100 µL
TSH3	200 µL	HCG	100 µL	P24	200 µL	HBS ^(a) confirmation test	150 µL	LYGS	100 µL
FT3	100 µL	FER	100 µL	P24 ^(a) confirmation test	200 µL	HBST	150 µL	VCAG	100 µL
FT4N	100 µL	IgE	100 µL	MSG	100 µL	HCV	100 µL	VCAM	100 µL
CORS	100 µL	AFP	100 µL	MPG	100 µL	HBCT	150 µL	EBNA	100 µL
ATPO	100 µL	E2II	200 µL	VZG	100 µL	HAVT	150 µL	VITD	100 µL
ATG	100 µL	AMH ^(a)	200 µL	RBM	100 µL	HBE	150 µL	TES2	100 µL
T3	100 µL	PRL	200 µL	TXC	125 µL	HBL	150 µL		
DIG	100 µL	LH	200 µL	TXGA ^(b)	2 x 100 µL	HBET ^(a) confirmation test	150 µL		
CKMB	250 µL	FSH	200 µL	TXG	100 µL	HAVM	100 µL		
TNHS	200 µL	PRG	200 µL	TXM	100 µL				
TNIU	200 µL	CEAS	200 µL	RBG	100 µL				
PBN2	200 µL	199	200 µL	B2M	100 µL				
GAL3	200 µL	125	200 µL	CMVG	100 µL				
PCT	200 µL	TPSA	200 µL	CMVA	2 x 100 µL				
MYO	150 µL	FPSA	200 µL	HPY	100 µL				
DEX2	200 µL	153	100 µL	CMVM	100 µL				
		GDH ^(a)	300 µL	HBCM	100 µL				
		CDAB ^(a)	300 µL						

(a) = After sample processing
(b) = After dilution to obtain 15 IU/mL

Units VIDAS® / mini VIDAS® / VIDAS 3®

Only parameters for which results can be expressed in different units are given in this table.

Parameter Code	Proposed VIDAS Unit	Conversion Factors	Parameter Code	Proposed VIDAS Unit	Conversion Factors	Parameter Code	Proposed VIDAS Unit	Conversion Factors	Parameter Code	Proposed VIDAS Unit	Conversion Factors
E2II	pmol/L nmol/L pg/mL	pg/mL x 3.67 → pmol/L pmol/L x 0.272 → pg/mL	PRL	µIU/mL µIU/mL ng/mL 1 IRP ng/mL 3 IS	1 ng 1 st IRP = 32 µ IU 1 ng 3 rd IS = 21 µ IU 1 ng/mL 3 rd IS = 1.524 ng/mL 1 st IRP	CEAS	ng/mL	ng/mL x 15.43 → mIU/mL	ST	kIU/L	
→ FT4	pg/mL ng/dL ng/L ng/100 mL pmol/L	pmol/L x 0.777 → pg/mL pg/mL x 1.29 → pmol/L	→ B2M	µg/mL UI/mL UI/mL kUI/L kIU/L mg/L	mg/L x 14 → UI/mL UI/mL x 0.071 → mg/L	→ AMH**	ng/mL pmol/L ng/dL	pmol/L x 0.14 → ng/mL ng/mL x 7.14 → pmol/L	→ TSH	mIU/L µIU/mL mIU/L µIU/mL	
→ T4	µg/L µg/100 mL µg/dL nmol/L	nmol/L x 0.777 → µg/L µg/L x 1.29 → nmol/L	→ AFP	µg/L ng/mL µg/mL IU/mL UI/mL	ng/mL x 0.826 → UI/mL UI/mL x 1.21 → ng/mL	→ FSH	IU/L mIU/mL UI/L mUI/mL		→ TSH3	µIU/mL µIU/mL	
→ FT3	pg/mL ng/L pmol/L	pmol/L x 0.651 → pg/mL pg/mL x 1.54 → pmol/L	→ TNHS	ng/L pg/mL		→ HAVT	IU/L mIU/mL UI/L mUI/mL		→ HBST	mIU/mL	
→ T3	ng/mL µg/L µg/100 mL µg/dL nmol/L	nmol/L x 0.651 → µg/L µg/L x 1.54 → nmol/L	→ TNIU	ng/mL µg/L		→ LH	IU/L kUI/L kIU/L mUI/mL		→ AHBS	mIU/mL	
→ FER	µg/L ng/mL		→ PBN2	pg/mL		→ HCG	UI/L mIU/mL IU/L kUI/L kIU/L mUI/mL		→ P24	pg/mL Ag P24 pg/mL Ag VIH	1pg/mL of P24 = 3.65 pg/mL
→ CORS	nmol/L µg/dL µg/L ng/mL	nmol/L x 0.362 → ng/mL ng/mL x 2.76 → nmol/L	→ PCT	µg/L ng/mL		→ DIG	nmol/L µg/L ng/mL	nmol/L x 0.781 → ng/mL ng/mL x 1.28 → nmol/L	→ CMVG	AU/mL UA/mL	
			→ MYO	µg/L ng/mL		→ PRG	nmol/L ng/mL	nmol/L x 0.3145 → ng/mL ng/mL x 3.1796 → nmol/L	→ RBG	UI/mL IU/mL	
			→ VITD	ng/mL nmol/L		→ TES2	nmol/L ng/mL ng/dL	nmol/L x 0.288 → ng/mL ng/mL x 3.47 → nmol/L	→ TXG	IU/mL UI/mL	
			→ GAL3	ng/mL		→ IgE	UI/mL kIU/L IU/mL kUI/L		→ vWF	UI/mL	
									→ PC	IU/mL %	
									→ DEX2	µg/mL ng/mL	(FEU) (FEU)
									→ HBCM	PEIU/mL	

FEU: Fibrinogen Equivalent Unit
→ Default Unit

Range	Parameter	Reference	Code	Sample Volume	Calibration-Control		Sample Pre-Treatment	Test validated using	Test time (minutes)	Measurement Range or Interpretation	
					Test	Frequency					
EMERGENCY	BACTERIAL INFECTION	B.R.A.H.M.S PCT™	30450	PCT	200 µL	S1 S1 S2 S2 C1 C2	28 days	no	Serum, Plasma (Hep.) ⁽²⁾	20	0.05 - 200 ng/mL
	3 punktás	NT-proBNP2	30458	PBN2	200 µL	S1 S1 S2 S2 C1 C2	28 days	no	Serum, Plasma (Hep.)	20	15 - 25,000 pg/mL
		Troponin I Ultra	30448	TNIU	200 µL	S1 S1 S2 S2 C1 C2	28 days	no	Serum, Plasma (Hep.)	20	0.01 - 30 µg/L
		Hs Troponin I ^(C)	415386	TNHS	200 µl	S1 S1 S2 S2 C1 C2	28 days	no	Serum, Plasma (Hep.)	20	
	CARDIAC	Galectin-3	411191	GAL3	200 µL	S1 S1 C1	28 days	no	Serum, Plasma (Hep., EDTA)	20	3.3 - 100 ng/mL
		Myoglobin	30446	MYO	150 µL	S1 S1 C1	14 days	no	Serum, Plasma (Hep.)	17	5 - 1,000 µg/L
		CK-MB	30421	CKMB	250 µL	S1 S1 C1	14 days	no	Serum, Plasma (Hep., EDTA)	30	0.8 - 300 ng/mL
		Digoxin	30603	DIG	100 µL	S1 S1 S1 C1	14 days	no	Serum, Plasma (Hep., EDTA)	20	0.2 - 5 ng/mL
	VENOUS THROMBO-EMBOLISM / COAGULATION	D-Dimer Exclusion™ II	30455	DEX2	200 µL	S1 S1 C1 C2	28 days	no	Plasma (Cit.)	20	45 - 10,000 ng/mL (FEU)
		Protein C	30115	PC	100 µL	S1 S1 C1	14 days	no	Plasma (Cit.)	35	1 - 120 %
		vWF	30436	vWF	100 µL	S1 S1 C1	14 days	no	Plasma (Cit.)	35	1 - 120 %
INFECTIOUS DISEASES	HEPATITIS	HAV IgM	30307	HAVM	100 µL	S1 S1 C1 C2	14 days	no	Serum, Plasma (EDTA, Hep.)	60	Qualitative test: Negative, Positive or Equivocal
		Anti-HAV Total	30312	HAVT	150 µL	S1 S1 C1 C2	14 days	no	Serum, Plasma (EDTA, Hep., Cit.)	90	15 - 400 mIU/mL
		HBs Ag Ultra ⁽¹⁾	30315	HBS• HBL	150 µL	S1 S1 C1 C2	14 days	no	Serum, Plasma (Hep.)	60 / 90	Qualitative test: Negative or Positive
		Anti-HBsT II	30318	AHBS	200 µL	S1 S1 C1 C2	28 days	yes ⁽²⁾	Serum, Plasma (Hep.)	60	3 - 500 mIU/mL
		Anti-HBc Total II	30314	HBCT	150 µL	S1 S1 S1 C1 C2	14 days	no	Serum, Plasma (EDTA, Hep., Cit.)	90	Qualitative test: Negative, Positive or Equivocal
		HBc IgM II	30439	HBCM	100 µL	S1 S1 C1 C2	14 days	no	Serum, Plasma (EDTA, Hep., Cit.)	55	0 - 200 PEIU/mL
		HBe Ag	30305	HBE	150 µL	S1 S1 C1 C2	14 days	no	Serum, Plasma (EDTA, Hep., Cit.)	90	Qualitative test: Negative or Positive
		Anti-HCV	30308	HCV	100 µL	S1 S1 C1 C2	28 days	no	Serum, Plasma (Hep.)	40	Qualitative test: Negative or Positive
	HIV	HIV DUO Ultra	30443	HIV5	200 µL	S1 S1 S2 S2 C1 C2 C3	14 days	no	Serum, Plasma (EDTA, Hep.) ⁽²⁾	120	Ag and Ab qualitative test: Negative or Positive
		HIV DUO Quick	30447	HIV6	200 µL	S1 S1 C1 C2 C3	14 days	no	Serum, Plasma (EDTA, Hep.) ⁽²⁾	80	Qualitative test: Negative or Positive
		HIV P24 II ⁽¹⁾	30117	P24	200 µL	S1 S1 C1 C2	14 days	no	Serum, Plasma (EDTA, Hep., Ok., Thrombin)	90	0 - 400 pg/mL of Ag P24
	ToRC	Toxo IgM	30202	TXM	100 µL	S1 S1 C1 C2	14 days	no	Serum ⁽³⁾	40	Qualitative test: Negative, Positive or Equivocal
		Toxo IgG II	30210	TXG	100 µL	S1 S1 C1 C2	14 days	no	Serum ⁽³⁾ , Plasma (EDTA, Hep.)	40	0 - 300 IU/mL
		Toxo Competition	30211	TXC	125 µL	S1 S1 C1 C2	14 days	no	Serum, Plasma (EDTA, Hep.)	40	Qualitative test: Negative or Positive
		Toxo IgG Avidity ⁽⁴⁾	30222	TXGA	2 x 100 µL	C1 C2	14 days	no	Serum ⁽³⁾ , Plasma (EDTA, Hep., Cit)	40	Avidity index (AI) An AI ≥ 0.3 is a strong indication of a primary infection dating back more than 4 months. An AI < than 0.3 does not enable a recent infection to be differentiated from a former infection.
		Rub IgM	30214	RBM	100 µL	S1 S1 C1 C2	14 days	no	Serum	60	Qualitative test: Negative, Positive or Equivocal
		Rub IgG II	30221	RBG	100 µL	S1 S1 C1 C2	14 days	no	Serum, Plasma (EDTA, Hep.)	40	0 - 400 IU/mL
		CMV IgM	30205	CMVM	100 µL	S1 S1 C1 C2	14 days	no	Serum	60	Qualitative test: Negative, Positive or Equivocal
		CMV IgG	30204	CMVG	100 µL	S1 S1 C1 C2	14 days	no	Serum	40	0 - 400 AU/mL
		CMV IgG Avidity II	413557	CMVA	2 x 100 µL	CH CL	14 days	no	Serum	40	Avidity index (AI) An AI ≥ 0.65 is a strong indication of a primary infection dating back more than 3 months. An AI < than 0.40 is a strong indication of a primary infection dating back less than 3 months. An 0.40 < AI < 0.65 does not enable to distinguish a recent infection from a former infection.
	OTHER SEROLOGIES	EBV VCA IgM	30237	VCAM	100 µL	S1 S1 S1 C1 C2	28 days	no	Serum	40	Qualitative test: Negative, Positive or Equivocal
		EBV VCA/EA IgG	30236	VCAG	100 µL	S1 S1 S1 C1 C2	28 days	no	Serum	40	Qualitative test: Negative, Positive or Equivocal
		EBV EBNA IgG	30235	EBNA	100 µL	S1 S1 S1 C1 C2	28 days	no	Serum	40	Qualitative test: Negative, Positive or Equivocal
Lyme IgM		30319	LYM	100 µL	S1 S1 C1 C2	28 days	no	Serum, Plasma (Hep.)	27	Qualitative test: Negative, Positive or Equivocal	
Lyme IgG		30320	LYG	100 µL	S1 S1 C1 C2	28 days	no	Serum, Plasma (Hep.), /CSF	27	Qualitative test: Negative or Positive	
Measles IgG		30219	MSG	100 µL	S1 S1 C1 C2	14 days	no	Serum	40	Qualitative test: Negative, Positive or Equivocal	
Mumps IgG		30218	MPG	100 µL	S1 S1 C1 C2	14 days	no	Serum	40	Qualitative test: Negative, Positive or Equivocal	
Varicella-Zoster IgG	30217	VZG	100 µL	S1 S1 C1 C2	14 days	no	Serum	40	Qualitative test: Negative, Positive or Equivocal		
H. pylori IgG	30192	HPY	100 µL	S1 S1 C1 C2	14 days	no	Serum, Plasma (EDTA)	40	Qualitative test: Negative, Positive or Equivocal		
ANTIGEN DETECTION	C. difficile GDH	30125	GDH	300 µL	S1 S1 C1 C2	28 days	yes ⁽²⁾	Stool	50	Qualitative test: Negative or Positive	
	C. difficile Toxin A&B	30118	CDAB	300 µL	S1 S1 C1 C2 C3	14 days	yes ⁽²⁾	Stool	75	Qualitative test: Negative, Positive or Equivocal	
REPRODUCTION / FERTILITY	AMH ^(d)	417011	AMH	200 µL	S1 S1 C1	28 days	no	Serum, Plasma (Lith., Hep.)	40	0.02 – 9 ng/mL	
	Estradiol II	30431	E2II	200 µL	S1 S1 S1 C1	14 days	no	Serum, Plasma (Hep.)	60	9 - 3,000 pg/mL	
	FSH	30407	FSH	200 µL	S1 S1 C1	14 days	no	Serum, Plasma (Hep.)	40	0.1 - 110 mIU/mL	
	HCG	30405	HCG	100 µL	S1 S1 C1	14 days	no	Serum, Plasma (Hep., EDTA)	30	2 - 1,500 mIU/mL	
	LH	30406	LH	200 µL	S1 S1 C1	14 days	no	Serum, Plasma (Hep.)	40	0.1 - 100 mIU/mL	
	Prolactin	30410	PRL	200 µL	S1 S1 C1	14 days	no	Serum, Plasma (Hep.)	40	0.5 - 200 ng/mL	
	Progesterone	30409	PRG	200 µL	S1 S1 S1 C1	14 days	no	Serum, Plasma (Hep., EDTA, Silicone, Gel)	45	0.25 - 80 ng/mL	
	Testosterone II	414320	TES2	100 µl	S1 S1 C1	28 days	no	Serum, Plasma	40	0.05 to 13.50 ng/mL	
	TSH	30400	TSH	200 µL	S1 S1 C1	14 days	no	Serum, Plasma (Hep., Beads, Gel)	40	0.05 - 60 µIU/mL	
	TSH3	30441	TSH3	200 µL	S1 S1 S2 S2 C1 C2	14 days	no	Serum, Plasma (Hep., Silicone, Gel)	80	0.005 - 100 µIU/mL	
	FT3	30402	FT3	100 µL	S1 S1 S1 C1	14 days	no	Serum, Plasma (Hep.)	40	0.7 - 45 pmol/L	
	FT4	30459	FT4N	100 µL	S1 S1 C1	14 days	no	serum, plasma (Hep., Silicone, Gel)	40	1 - 100 pmol/mL	
	T3	30403	T3	100 µL	S1 S1 S1 C1	14 days	no	Serum, Plasma (Hep.)	40	0.4 - 9 nmol/L	
	IMMUNOCHEMISTRY	T4	30404	T4	200 µL	S1 S1 S1 C1	14 days	no	Serum, Plasma (Hep.)	40	6 - 320 nmol/L
		Anti-TPO	30461	ATPO	100 µL	S1 S1 C1	28 days	no	Serum, Plasma (Hep., EDTA, Silicone, Gel)	25	0.8 - 1,000 IU/mL
		Anti-Tg	30462	ATG	100 µL	S1 S1 C1	28 days	no	Serum, Plasma	25	6.4 - 800 IU/mL
		AFP	30413	AFP	100 µL	S1 S1 C1	14 days	no	Serum, Plasma (Hep., EDTA), Amniotic Fluid	30	0.5 - 400 IU/mL
		CA 15-3*	30429	153	100 µL	S1 S1 C1	14 days	no	Serum, Plasma (Hep, EDTA)	60	2 - 400 U/mL
CA 19-9™		30427	199	200 µL	S1 S1 C1	14 days	no	Serum, Plasma (Hep, EDTA)	60	3 - 500 U/mL	
CA 125 II™		30426	125	200 µL	S1 S1 C1	14 days	no	Serum, Plasma (Hep, EDTA)	60	4 - 600 U/mL	
CEA S		30453	CEAS	200 µL	S1 S1 C1	14 days	no	Serum, Plasma (Hep.)	60	0.5 - 200 ng/mL	
FPSA		30440	FPSA	200 µL	S1 S1 C1	14 days	no	Serum, Plasma (Hep., EDTA)	60	0.05 - 10 ng/mL	
TPSA		30428	TPSA	200 µL	S1 S1 C1	14 days	no	Serum, Plasma (Hep., EDTA)	60	0.07 - 100 ng/mL	
ALLERGY	Total IgE	30419	IgE	100 µL	S1 S1 C1	14 days	no	Serum, Plasma (Lith., Hep., EDTA)	30	0.5 - 1,000 kIU/L	
OTHERS	Cortisol S	30451	CORS	100 µL	S1 S1 S1 C1	14 days	yes (urine)	Serum, Plasma (Hep., EDTA), 24-hour Urine	40	2 - 650 ng/mL	
	Ferritin	30411	FER	100 µL	S1 S1 C1	14 days	no	Serum, Plasma (Hep., EDTA)	30	1.5 - 1,200 ng/mL	
	B2 Microglobulin	30420	B2M	100 µL	S1 S1 C1 C2	14 days	yes (urine)	Serum, Plasma (Hep., EDTA), Urine	40	0.004 - 4 mg/L	
	25 OH Vitamin D Total	30463	VITD	100 µL	S1 S1 C1	28 days	no	Serum, Plasma (Lith., Hep.)	40	8.1 - 126 ng/mL	

• HBL = long protocol

(c) = Some of these reagents have not yet obtained regulatory clearance in some countries. Please contact your local bioMérieux representative for further information and product availability.

Some references may vary according the country. Please contact your local bioMérieux representative for further information and product availability.

(d) = Not available in the following countries, regions and states: Armenia, Australia, Austria, Azerbaijan, Belarus, Belgium, Canada, Denmark, France, Germany, Hong-Kong, Ireland, Israel, Italy, Japan, Kazakhstan, Kyrgyzstan, Lichtenstein, Moldavia, New Zealand, Portugal, Russia, Spain, Switzerland, Tajikistan, The Netherlands, Turkey, Turkmenistan, United Kingdom.

(1) = Confirmation test available
(2) = See package insert

(3) = Inactivated or not 30 mins at 56°C
(4) = Double strips

Benefit from bioMérieux's experience with VIDAS® services, ensuring your lab's continued development

Instrument services

- Instant remote support service
- Extensive network of local subsidiaries and distributors:
skilled systems engineers and application specialists available wherever you are

Certification support and integrated quality control

- VIDAS validation package for your instrument and methodology
- Support in obtaining ISO 15189 certification

Online training

- Online training sessions complement VIDAS on-site service
- Continually develop your lab's skills and expertise

A worldwide VIDAS community

- The myvidas.com platform strengthens your expertise and makes
the most of your VIDAS solution
- Share information and experiences with your VIDAS peers

Vidas

Scan this code to get access to
myvidas.com and
the constant updated
information on
VIDAS assay solutions.

bioMérieux S.A.
69280 Marcy l'Etoile
France
Tél. : +33 (0)4 78 87 20 00
Fax : +33 (0)4 78 87 20 90

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Marcy l'Etoile, September 20th, 2013

Dear Customer,

You have chosen the VIDAS® system for your laboratory and we thank you for your confidence.

The aim of this letter is to give you the bioMérieux recommendations regarding calibration and controls.

The VIDAS® system processes samples in random access and batch mode using a single test format (strip and SPR) based on the Enzyme Linked Fluorescent Assay (ELFA) technology for all types of analysis: serology, immunochemistry, antigen detection and industrial microbiology.

For each parameter, a calibration should be done at the opening of the kit and a recalibration every 14 or 28 days (depending on the assay, please refer to package insert) to readjust the calibration curve.

The purpose of the controls included in the kit is to validate calibrations. They should also be run immediately after opening a new kit to ensure that reagent performances have not been altered.

Compliance with these recommendations allows your laboratory to be in line with bioMérieux guidelines.

Nevertheless, you may adapt your quality control frequency according to your laboratory organisation based on a risk analysis. bioMérieux recommends to use the Ishikawa diagram to identify causes likely to present a risk and find the appropriate answer (training, procedures, quality controls ...). We can help you in this approach.

Please distribute this information to all appropriate personnel in your laboratory, keep a copy in your files, and forward it to other parties to whom you may have transferred our product.

We thank you for applying these procedures in your laboratory. If you have any questions or need assistance, please contact your local bioMérieux representative.

We are committed to providing you with products, support and services of the highest quality. Thank you for your continued confidence in bioMérieux and our products.

Gabriel Pedone
Global Marketing Manager
ImmunoAssay Unit

Marie-Claude Bernard
Product Manager ImmunoAssays
Middle-East/ Africa/ Europe Distributors

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THE VIDAS CONCEPT

1 patient, 1 test, 1 result

can be used 24/7, with over 100 references in single test format

EMERGENCY

PBN2

DEX2

TNHS

VWF

CKMB

DIG

MYO

PC

GAL3

PCT

INFECTIOUS DISEASES

HBE

TXC

HAUT

CMVG

HIV6

CMVM

TXM

HBCM

HBET

AMH

TES2

ATG

125

IGE

PRG

IMMUNOCHEMISTRY

TSH

153

VITD

ATPO

COR5

T4

PRL

FT4N

LH

AFP

E2N

CEAS

FSH

B2M

TSH3

HCG

FT3

FER

YES

FPSA

199

T3

TPSA

TXGA

Hand icon

Reagent icons

Some of these reagents have not yet obtained regulatory clearance in some countries. Please contact your local bioMérieux representative for further information and product availability. Please refer to your local and/or international guidelines for the use of biomarkers.


Dear Customer,

The world of medicine is changing. Today, its focus has shifted from diagnosing a pathology to assessing what is unique about each patient.

The individual is always at the heart of the VIDAS® approach, yet its strength also lies in its exceptional network of over **25,000 customers**. Together, the VIDAS global teams and network make up an impressive force, ready to embrace the new medical, scientific and economic realities of healthcare today.

The VIDAS Solution will make an essential contribution to meeting the demands of a world where controlling health expenditure has become crucial and where, as a result, diagnostics have taken on critical importance.

Your bioMérieux team

* VIDAS® AMH is not available in the following countries, states and regions: Armenia, Australia, Austria, Azerbaijan, Belarus, Belgium, Canada, Denmark, France, Germany, Hong-Kong, Ireland, Israel, Italy, Japan, Kazakhstan, Kyrgyzstan, Liechtenstein, Moldavia, New Zealand, Portugal, Russia, Spain, Switzerland, Tajikistan, The Netherlands, Turkey, Turkmenistan, United Kingdom.

THE SINGLE TEST CONCEPT

The success of VIDAS is based on the robustness and reliability of its reagent design, combined with its ready-to-use single test concept with ELFA* technology.

HOW IT WORKS

The VIDAS principle is based on the interaction of two elements: the coated SPR® (Solid Phase Receptacle), containing antigens or antibodies, and the Strip, made up of a series of wells containing the correct amount of reagent necessary for the test.

CONFIDENCE IN RESULT REPORTING AND COST MANAGEMENT

High quality VIDAS reagents offer a dynamic measurement range **with proven accuracy and performance**, reducing unnecessary repeat tests.

All types of assays can be run simultaneously **without any risk of contamination** due to absence of external washing solution, drains and syringes for sample or reagent pipetting.

Each parameter has a single reference with reagent, standard, control and diluent included in the kit. Management is simple with minimum operator handling and intervention.

* Enzyme Linked Fluorescent Assay.

THE VIDAS® INSTRUMENTS RANGE

VIDAS instruments are bench-top automated immunoassay systems based on ELFA technology. They offer your laboratory high quality, on demand test results for small volume testing.

All VIDAS instruments are operational 24/7 and use the same reagent references. They allow you to process several different parameters simultaneously, in both single random access sample and batch test mode, for all types of analyses. Their unique design guarantees that there is no inter-reagent or inter-sample contamination.

mini VIDAS

SIMPLY RELIABLE

With mini VIDAS® your lab can carry out tests simply and quickly, in a self-contained, space-saving unit. Its user-friendly Load & Go concept means your teams can process up to 36 tests an hour. It also offers excellent reliability, with a Mean Time Between Failures (MTBF) of over 1,100 days.

VIDAS

THE PIONEER

VIDAS® offers your lab simplicity and productivity, with its loading capacity of 30 tests. It will enhance your workflow, allowing you to handle batch series while maintaining specific management of emergency tests, thanks to its five independent compartments.

VIDAS® also guarantees fast, reliable results, processing up to 80 tests per hour and offering a Mean Time Between Failures (MTBF) of over 700 days.

VIDAS 3

ENHANCED IN-LAB ORGANISATION

The VIDAS 3 system brings maximum peace of mind to your lab. Flexible and fully automated, it enables you to run up to 12 tests in 4 independent sections with complete security, guaranteed traceability and 24/7 access.

VIDAS 3 is designed with:

- Intuitive software,
- Fully automated pipettor action.

VIDAS 3 ensures test security thanks to:

- Complete traceability of all actions,
- Integrated quality control management,
- Bi-directional Laboratory Information System (LIS) connectivity.


VIDAS 3 is flexible with access at all times:

- Temperature-controlled loading bay, 27 samples, diluents and calibrators in three separate racks,
- 4 independent test compartments of 3 tests each - Up to 36 tests per hour.

ACUTE CARE

High sensitive Troponin I, NT-proBNP, D-Dimer, PCT, a unique test offer on a single device

Efficient and fast patient triage is essential in the Emergency Department, where common symptoms like chest pain, shortness of breath and fever require rapid decisions on hospitalization, investigations and treatment.¹

* Among the top-10 reasons for visiting an emergency department (National Hospital Ambulatory Medical Care Survey, 2011 Emergency Department Fact Sheet: www.cdc.gov/nchs/ahcd/factsheets.htm).
HF: Heart Failure - ACS: Acute Coronary Syndrome - PE: Pulmonary Embolism

VIDAS® and its unique, comprehensive and user-friendly panel allows you to manage urgent samples separately from routine activity and comply with the short turnaround time required by such complex clinical situations. That means improved patient management and less stress for doctors.

	VIDAS® B.R.A.H.M.S PCT™	VIDAS® D-Dimer Exclusion™ II	VIDAS® High sensitive Troponin I	VIDAS® NT-proBNP2	VIDAS® Galectin-3
Sample type	Plasma (heparin) Serum	Plasma (citrate)	Plasma (heparin) Serum	Plasma (heparin) Serum	Plasma (heparin, EDTA) Serum
Time to result	20 min (no warm-up period, load and start)	20 min (no warm-up period, load and start)	20 min (no warm-up period, load and start)	20 min (no warm-up period, load and start)	20 min (no warm-up period, load and start)

1. Wiler JL, Gentle C, Halfpenny JM, Heins A, Mehrotra A, Mikhail MG, Fite D. Optimizing emergency department front-end operations. *Ann Emerg Med*. 2010;55:142-160.e1.


SEPSIS

B.R.A.H.M.S PCT™ Guidance for Antibiotic therapy

When emergency teams suspect severe bacterial infection, sepsis or lower respiratory tract infection (LRTI), early diagnosis is essential.

Procalcitonin (PCT) blood levels rise rapidly during bacterial infections. It has been demonstrated that PCT helps to assess the severity of an infection and supports the early recognition of sepsis.¹ Procalcitonin has also been shown to be a valuable decision-support tool for the initiation, continuation or cessation of antibiotic treatment.

The use of PCT to guide antibiotic therapy has been described in various clinical settings. One example is the management of patients with community-acquired pneumonia (CAP) in emergency departments:

With PCT guidance, patients were treated for a mean of 7 days compared to 11.1 days in the control group, indicating a reduction in antibiotic exposure of around 40% (Schuetz *et al.*, 2012²).

1. Bouadma L, Luyt CE, Tubach F *et al.* Use of procalcitonin to reduce patients' exposure to antibiotics in intensive care units (PRORATA trial): a multicentre randomised controlled trial. *Lancet* 2010;375 (9713 463-74).
2. Schuetz P, Briel M, Christ-Crain M *et al.* Procalcitonin to guide initiation and duration of antibiotic treatment in acute respiratory infections: an individual patient data meta-analysis. *Clin Infect Dis*. 2012;55:651-662.

ACUTE CARE PANEL
**B.R.A.H.M.S PCT™, D-Dimer Exclusion™ II,
High sensitive Troponin I, NT-proBNP2,
Galectin-3, CK-MB, Myoglobin, Digoxin**

CARDIAC


A comprehensive approach to cardiac patient management

ACUTE CORONARY SYNDROME (ACS)

VIDAS® High sensitive Troponin I

This assay is an aid in the diagnosis of myocardial infarction and in the risk stratification of patients with symptoms suggestive of ACS with respect to relative risk of all-cause mortality and major adverse cardiac events (MACE) consisting of MI revascularization, at 30 days. VIDAS® High sensitive Troponin I offers analytical and clinical performances in accordance with guidelines¹:

- CV=7%
- accelerated and validated algorithm for delivering a diagnosis in only 2 hours - safe rule-out and accurate rule-in - of acute Myocardial Infarction (Non ST Elevation Myocardial Infarction) in the vast majority (70%) of suspected chest pain patients.

HEART FAILURE (HF)


VIDAS® NT-proBNP2

Early diagnosis of heart failure is key to improving patient outcome. Natriuretic peptides are recommended in international guidelines as a biomarker to rule out Heart Failure.¹

VIDAS® Galectin-3

Testing for elevated blood levels of galectin-3 allows an aid in diagnosis of chronic heart failure patients and identifies those with a 2 to 3 times higher risk of re-hospitalization or mortality.²⁻⁴

1. Eur Heart J. 2016; 37: 267-315.
2. Ponikowski P, et al. Eur Heart J. 2016 Jul 14;37(27):2129-200.
3. Lok DJ, et al. Clin Res Cardiol. 2010;99:323-8.
4. de Boer RA, et al. Ann Med. 2011;43:60-8.
5. Felker GM, et al. Circ Heart Fail. 2012;5:72-78.

CARDIAC PANEL

High sensitive Troponin I, NT-proBNP2, Galectin-3, CK-MB, Myoglobin, Digoxin

VENOUS THROMBOEMBOLISM (VTE)

VIDAS® D-Dimer EXCLUSION™ II The reference to exclude VTE

Two key steps are essential for cost-efficient and safe exclusion of VTE.^{1,2}

Clinical assessment using validated pre-test probability (PTP) classifies patient risk for VTE³⁻⁵ and should be followed by D-Dimer testing using a validated high-sensitivity assay (≥ 97%).⁶

The VIDAS® D-Dimer test safely excludes outpatients with suspected VTE.


VIDAS® D-Dimer Exclusion™ II with sensitivity of > 99% (NPV* > 99%) allows the safe exclusion of up to 50% of outpatients with suspected VTE and reduces missed diagnoses.

VIDAS® D-Dimer has the most published studies supporting it, including 7 prospective outcome studies (>8,000 patients enrolled)⁷. It is widely considered to be the Gold Standard and is FDA cleared for VTE exclusion.


* Negative Predictive Value.

1. Righini M, et al. J Thromb Haemost. 2007;5:1869-7.
2. Ten Cate-Hoek AJ, et al. J Thromb Haemost. 2005;3:2465-70.
3. Wells PS, et al. N Engl J Med. 2003;349:1227-35.
4. Wells PS, et al. Thromb Haemost. 2000;83:416-20.
5. Le Gal G, et al. Ann Intern Med. 2006;144:165-71.
6. Wayne PA. CLSI document H59-A. Clinical and Laboratory Standards Institute, PA, USA, 2011.
7. Carrier M, et al. Thromb Haemost. 2009;101:886-92.

Assay sensitivity and missed diagnoses per 100 D-dimer screens

PATIENTS: without VTE (grey), with VTE (red)
 MISSED DIAGNOSES: Assay A (Se 99%) (green), Assay B (Se 95%) (yellow), Assay C (Se 88%) (red)
D-Dimer is not recommended for exclusion in High PTP group

HEMOSTASIS PANEL


D-Dimer Exclusion™ II, Protein C, vWF

ToRC

Serenity throughout pregnancy


With the VIDAS® ToRC panel your lab has a complete solution for the screening and diagnosis of Toxoplasmosis, Rubella and Cytomegalovirus infections in pregnant women and newborns. This includes **high performance CMV and toxoplasmosis IgG Avidity assays** which are of utmost importance to help distinguish between recent and past infections if maternal IgM and IgG are both positive during pregnancy.

VIDAS® CMV IgG AVIDITY II

The VIDAS® CMV IgG Avidity assay correlates very well with the time from onset of infection.¹ A study conducted to assess VIDAS® CMV IgG Avidity II maturation kinetics on 135 sequential samples from 31 patients¹ confirmed a very good correlation to clinical diagnosis.

//The VIDAS assay was the most regular, showing the best correlation over time from onset of infection. No sample was wrongly misclassified.¹ //


1. C. Vauloup-Fellous *et al.* Re-evaluation of the VIDAS® cytomegalovirus (CMV) IgG avidity assay: determination of new cut-off values based on the study of kinetics of CMV-IgG maturation. J Clin Virol 2013; 56(2): 118-23.

ToRC PANEL


TOXO: **TOXO IgG II, TOXO IgM, TOXO Competition, TOXO IgG Avidity,**
 CMV: **CMV IgG, CMV IgM, CMV IgG Avidity II,**
 RUBELLA: **RUB IgG II, RUB IgM**

VIDAS® TOXO IgG AVIDITY

The most important parameter required for a Toxo IgG avidity assay is the ability to provide high avidity values for past toxoplasmosis infection.² VIDAS® TOXO IgG Avidity identifies more old infections than other commonly-used tests.

2. Murat J-B., *et al.* Comparison of the VIDAS System and Two Recent Fully Automated Assays for Diagnosis and Follow-Up of Toxoplasmosis in Pregnant Women and Newborns. Clin Vaccine Immunol 2013;20(8):1203-12.

//The VIDAS TOXO IgG Avidity had the best performance for the diagnosis of latent toxoplasmosis. //


Villard, O., *et al.* Comparison of Four Commercially Available Avidity Tests for Toxoplasma gondii-Specific IgG Antibodies. Clin Vaccine Immunol, 2013; 20(2): 197-204.

INFECTIOUS MONONUCLEOSIS EBV

Reassuring when symptoms are unclear

The VIDAS® EBV panel consists of 4 EBV-specific markers in 3 tests to establish the stage of infection:

- Absence of infection
- Primary infection
- Past infection

EBV PANEL

EBV VCA IgM, EBV VCA/EA IgG, EBV EBNA IgG


HIV

Ensuring early diagnosis

The VIDAS® HIV panel is a combination of advanced 4th generation tests based on a unique assay concept.

It offers one of the highest levels of sensitivity currently available, allowing earlier detection of HIV infection and the confirmation of the presence of p24 Ag - just two weeks after exposure to the virus. This rapid detection can make the difference in reducing transmission and allows clinicians to provide appropriate care as soon as possible.

HIV PANEL

HIV DUO Ultra, HIV DUO Quick, HIV P24 II, HIV P24 II Confirmation

HEPATITIS

Access to hepatitis profile made easy

Viral hepatitis infections are a major public health problem across the globe, affecting hundreds of millions worldwide and killing some 1.4 million every year.¹

The VIDAS® Hepatitis panel covers all the essential parameters for rapid and precise differential diagnosis of HAV, HBV and HCV infections, supporting appropriate treatment and enhanced patient care.

VIDAS® HBs Ag Ultra

Ultra sensitive test enabling to reduce the serological window and detect hidden infections.

High capacity for the detection of various mutants/variants in whole blood samples.²

VIDAS® Anti-HCV

Sensitive 3rd generation assay that detects all 6 HCV genotypes, ensuring clear and confident diagnosis of Hepatitis C, with results available in 40 minutes.

The solid performance of the VIDAS® anti-HCV assay in clinical trials confirms the high level of quality labs have come to expect from VIDAS, for complete peace of mind.

Diagnostic Sensitivity (95% Confidence Interval)	
HCV positive / HIV negative patients (n = 254)	100% [98.56% - 100%]
HCV positive / HIV positive patients (n = 61*)	98.36% [91.20% - 99.96%]
HCV positive / (HIV unknow) patients (n = 124)	100% [97.07% - 100%]
Total population (n = 439*)	99.77% [98.74% - 99.99%]

* 1 patient who was not detected using Anti-HCV either had a low antibody level or was not detected using equivalent method.

1. World Health Organisation, World Hepatitis Day Campaign 2013.

2. Servant-Delmas *et al.*, Variable capacity of 13 hepatitis B virus surface antigen assays for the detection of HBsAg mutants in blood samples, Journal of Clinical Virology, 53 (2012) 338-345.

HEPATITIS PANEL

HEPATITIS A: **HAV IgM, Anti-HAV Total**

HEPATITIS B: **HBs Ag Ultra, HBs Ag Ultra Confirmation, Anti-HBs Total II, Anti-HBc Total II, HBc IgM II, HBe/Anti-HBe**

HEPATITIS C: **Anti-HCV**

LYME DISEASE

Get the right diagnosis

Lyme Borreliosis is one of the world's most common and rapidly spreading tick-borne diseases.

The VIDAS® Lyme IgM and IgG panel allows your lab to establish a reliable serological profile in just 27 minutes. It is an aid for rapid differential diagnosis and easier classification of Lyme borreliosis. This ensures optimized patient treatment according to the stage of infection and rational use of antibiotics.

The good specificity and sensitivity of the VIDAS® Lyme IgM and IgG panel mean it detects all major pathogenic strains¹, reducing the need for additional testing. Studies have demonstrated the panel's good performance in areas with both high and low prevalence.

1. *Borrelia burgdorferi sensu stricto*, *Borrelia afzelii*, *Borrelia garinii*.

Excellent specificity observed with VIDAS Lyme IgG on sera from low-endemic area.

Dr. Olivier Péter, Valais Hospital – Central Institute Microbiology and Infectious Diseases, Sion Partner organization of the National Reference Center for Tick-Borne Diseases, Neuchâtel, Switzerland

H. PYLORI

Non-invasive detection of *H. pylori* infection

The *H. pylori* bacteria is one of the most common causes of stomach ulcers and may increase the risk of stomach cancer.

The VIDAS *H. pylori* IgG test is a qualitative, non-invasive test providing results in just 40 minutes to aid diagnosis of *H. pylori* infection in symptomatic adults. Easy to perform, with kits adapted to low volume testing, it detects anti-*Helicobacter pylori* IgG antibodies in human serum or plasma.


C. DIFFICILE

Simplify diagnosis and outbreak management through automation


Clostridium difficile accounts for up to 25% of cases of healthcare-associated diarrhea and is the primary cause of antibiotic-associated colitis.¹ Management of *C. difficile* infection requires patient isolation, extensive treatment and often leads to a longer hospital stay.

Part of the comprehensive bioMérieux *C. difficile* solution, the VIDAS® *C. difficile* automated panel rapidly gives you the information you need to help identify and contain the spread of the infection during sporadic or epidemic outbreaks.

The VIDAS® *C. difficile* panel proposes a two-step algorithm for the diagnosis of *C. difficile* infection: the VIDAS® *C. difficile* GDH assay provides highly sensitive screening of the *C. difficile* antigen, glutamate dehydrogenase (GDH) in stool specimens, followed by confirmation of positive GDH samples with VIDAS® *C. difficile* Toxin A & B.¹

1. Bartlett JG. Clinical Practice. Antibiotic-associated diarrhea. N Engl J Med. 2002;346:334-349.


OTHER INFECTIOUS DISEASES ASSAYS

***C. difficile* GDH, *C. difficile* Toxin A & B
Lyme IgM, Lyme IgG,
H. pylori IgG
Measles IgG, Mumps IgG, Varicella Zoster IgG**

INFERTILITY MANAGEMENT

When nature needs a little help

Hormone measurements are essential to help in the correct diagnosis of the cause of infertility and to monitor the progress of ovarian stimulation.

* *in vitro* fertilization
** Intracytoplasmic Sperm Injection

FERTILITY PANEL
AMH, LH, FSH
Estradiol II, Progesterone
Testosterone II, Prolactin, hCG

The VIDAS® Fertility solution makes testing simple for Fertility clinics

The VIDAS® Fertility panel is a convenient solution for on-site testing for the investigation and treatment of infertility. It includes the key parameters in this field - LH, FSH, Progesterone, Estradiol, Testosterone, Prolactin, hCG and now AMH*.

As an indicator of the ovarian follicle reserve, AMH is useful in different contexts all along the reproductive lifespan in women:

- Personalization of infertility management (IVF/ART)
- Detection of ovarian dysfunction (Polycystic Ovary Syndrome, Premature ovarian failure...)
- Monitoring of ovarian damage further to surgery or gonadotoxic therapies

Tests for infectious diseases can also be performed on the same instrument: HIV, Hepatitis, Toxoplasmosis, Rubella, Cytomegalovirus. This is important for prevention both before and during pregnancy and to allow gamete storage.

Clinical evaluation of VIDAS® Estradiol and VIDAS® Progesterone has confirmed the reliability of both assays for use in evaluating natural cycles and monitoring ovarian superovulation¹:

VIDAS E2 and PRG assays are precise and sensitive tools for evaluating natural cycles [...] and for monitoring ovarian superovulation.¹

1. E. Anckaert *et al.* Clinical Validation of a Fully Automated 17 β -Estradiol and Progesterone Assay (VIDAS®) for Use in Monitoring Assisted Reproduction Treatment. Clin Chem Lab Med 2002; 40(8):824-831.

* VIDAS® AMH is not available in the following countries, states and regions: Armenia, Australia, Austria, Azerbaijan, Belarus, Belgium, Canada, Denmark, France, Germany, Hong-Kong, Ireland, Israel, Italy, Japan, Kazakhstan, Kyrgyzstan, Liechtenstein, Moldavia, New Zealand, Portugal, Russia, Spain, Switzerland, Tajikistan, The Netherlands, Turkey, Turkmenistan, United Kingdom.


THYROID

Targeted treatment starts with accurate profiling

The VIDAS® Thyroid panel includes all the parameters your lab needs for more accurate, rapid diagnosis and follow-up of thyroid disorders. As well as the key TSH, FT4 and FT3 assays, the panel now also includes VIDAS® Anti-TPO and VIDAS® Anti-Tg tests, which are important when further investigation of thyroid dysfunction is required.


THYROID PANEL
**TSH, T3, T4, FT3, FT4,
Anti-TPO, Anti-Tg**

TUMOR MARKERS

Accurate monitoring for better patient management

The VIDAS® range of tumor markers offers your lab a rapid, automated, and reliable solution to meet your clinicians' main requirements.


TUMOR MARKER PANEL
**AFP, CEA (S), CA 15-3®, CA 19-9™,
CA 125 II™, FPSA, TPSA**

ALLERGY

Knowing the cause helps bring relief

The VIDAS® TOTAL IgE assay is the first global screening step in diagnosis of allergic diseases. A high concentration of IgE (adults > 150 kIU/L | children under 3 > 40 kIU/L, then increasing to reach adult values at about 8 years old) indicates a high probability of allergy.



ALLERGY PANEL
Total IgE

VITAMIN D TOTAL

Reach out for optimum wellness

VIDAS® 25 OH Vitamin D Total detects both vitamin D2 and D3 for reliable results and is easy to use, with no sample pretreatment. All reagents are included in a single strip and have very good correlation with LC-MS/MS.

OTHER ASSAYS
**Vitamin D, β2 Microglobulin,
Cortisol S, Ferritin**

SERVICES

As a VIDAS® user, you can also benefit from a range of services tailored to ensure your lab's continued development.

Skilled system engineers and application specialists available wherever your VIDAS® instrument is installed.

Our remote instrument support services provide you with immediate solutions for VIDAS® and VIDAS® 3 instruments.

myQC Quality Control program gives you a large statistical potential for QC plan assessment.

Our validation package includes all the documentation required to validate your instrument.

Support provided to accompany you on your **lab certification process**.

Online training sessions complement VIDAS® on-site service, and offer continuous improvement for your lab teams.

myQC SOLUTION!

Based on your lab risk analysis and quality control plan, internal Quality Control is essential to:

- sustain your quality management system,
- obtain and maintain ISO 15 189 certification

bioMérieux offers a complete solution: internal quality control and support with a single point of contact for greater efficiency.

myQC CONTROL TEST KITS

✓ myQC Routine + TM

✓ myQC Cardiac

✓ myQC Specialty

Name	myQC Routine + TM	myQC Cardiac	myQC Specialty
Common analytes	TSH, TSH3, T4, FT4, T3, FT3, Testosterone, Progesterone, HCG, Ferritin, 82 Microglobulin, AFP, LH, FSH, Estradiol, Cortisol, FPSA, TPSA, IgE, Prolactin, CEA, CA 125, CA 19-9, CA 15-3	Myoglobin, NT-proBNP, D-Dimer, High sensitive Troponin I	B.R.A.H.M.S PCT, Anti-Tg, Anti-TPO

myQC SERVICES

In addition to standard customer support, myQC services are part of a global offer designed by bioMérieux Performance Solutions™ with accreditation support through method verification, QC analysis and e-learning.

bioMérieux S.A.
69280 Marcy l'Etoile
France

Tel.: +33 (0)4 78 87 20 00
Fax: +33 (0)4 78 87 20 90

www.biomerieux.com
www.biomerieux-diagnostics.com

myvidas.com

EU DECLARATION OF CONFORMITY

1. Product model/product:
Ref. 410416 mini VIDAS®
Ref. 410417 VIDAS®
2. Name and address of the manufacturer or his authorised representative:
MANUFACTURER : bioMérieux S.A.
376 Chemin de l'Orme
69280 - MARCY L'ETOILE – FRANCE
3. This declaration of conformity is issued under the sole responsibility of the manufacturer.
4. Object of the declaration:
In Vitro Diagnostic Medical Device
5. The object of the declaration described above is in conformity with the relevant Union harmonisation legislation(s):
Directive 2011/65/EU – RoHS
6. References to the relevant harmonised standards used or references to the other technical specifications in relation to which conformity is declared:
EN 50581: 2012
7. Where applicable, the notified body:
-
8. Where applicable, description of accessories and components, including software, which allow the radio equipment to operate as intended and covered by the EU declaration of conformity:
-
9. Additional information:
VIDAS Family RoHS II Implementation Summary Report: SV-VD008

Signed for and on behalf of: **bioMérieux SA.**

Place : Craponne
Date : July 22th, 2016

Name & position :

Catherine Fritsch
Regulatory Affairs Director
Immunology Unit

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