# NADAL<sup>®</sup> Celiac Disease (tTG) (test cassette)

REF 790001

<b>DE</b>	Gebrauchsanweisung	p. 2	<b>PL</b>	Sposób użycia	p. 24
<b>EN</b>	Instruction for use	p. 7		Symbols	p. 32
<b>FR</b>	Instructions d'utilisation	p. 11		Our Teams	p. 32
<b>ES</b>	Instrucciones de uso	p. 15			
<b>IT</b>	Istruzioni per l'uso	p. 19			

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Carl-Zeiss-Strasse 12  
47445 Moers  
Germany

Tel: +49 (2841) 99820-0  
Fax: +49 (2841) 99820-1

Friedenstrasse 32  
93053 Regensburg  
Germany

Tel: +49 941 29010-0  
Fax: +49 941 29010-50

[www.nal-vonminden.com](http://www.nal-vonminden.com)  
[info@nal-vonminden.com](mailto:info@nal-vonminden.com)

Directors:  
Sandra von Minden  
Roland Meißner  
Thomas Zander

Commercial reg. Kleve  
HRB 5679  
Steuer-Nr. 244/133/00130  
UST-ID-Nr. DE 189 016 086

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#### 14. Referenzen

1. Zhang J. et al. "Modulation of the in situ activity of tissue transglutaminase by calcium and GTP". The Journal of Biological Chemistry (1998) 273 (4): 2288-2295.
2. Sárdy M. et al. "Recombinant human tissue transglutaminase ELISA for the diagnosis of gluten-sensitive enteropathy". Clinical Chemistry (1999) 45 (12): 2142-2149.
3. Sorell L. et al. "One-step immunochromatographic assay for screening of celiac disease". Lancet (2002) 359: 945-946.
4. Drago S. et al. "Recent developments in the pathogenesis, diagnosis and treatment of celiac disease". Expert Opinion in Therapy Patents (2002) 12 (1): 45-51.
5. Hansson T. et al. "Recombinant human tissue transglutaminase for diagnosis and follow-up of childhood coeliac disease". Pediatric Research (2002) 51 (6): 700-705.
6. Ferre-López, S. et al. "A reliable generation of immunochromatographic sticks for both t-TG and AGA screening in Coeliac Disease". Clinical Gastroenterology and Hepatology (2004) 2 (6): 480-484.
7. Jennings J. et al. "New developments in celiac disease". Current Opinion in Gastroenterology (2003) 19 (2): 118-129.

Rev. 1, 2012-01-23 SS

### 1. Intended Use

The **NADAL® Celiac Disease (tTG) Test** is an easy-to-use immunochromatographic test designed for the detection of antibodies (IgA/IgG/IgM) against human tissue transglutaminase (tTG), the main autoantigen recognized by anti-endomysial antibodies, in whole blood or serum.

### 2. Introduction and Clinical Significances

Celiac Disease (CD) is an intestinal disorder (enteropathy) caused by a permanent intolerance to gluten, a protein that is a component of wheat, barley and rye products. In genetically susceptible people the repeated ingestion of this protein may cause serious lesions of the intestinal mucosa. The diagnosis of CD is based on characteristic histological changes, such as crypt hyperplasia with complete or incomplete atrophy of intestinal villi, observed in intestinal biopsies. However, serological tests such as detection of circulating anti-gliadin, anti-tTG and anti-endomysial antibodies represent less invasive and cheaper analysis methods to detect this disease.

### 3. Principle of the Test

Anti-tTG antibodies present in the blood of patients having celiac disease react with latex particles conjugated with human recombinant tTG. These complexes formed of latex particles, tTG and anti-tTG antibodies, migrate through a chromatographic process towards the reaction area. This area contains human recombinant tTG, which reacts with the complex, generating a red line.

### 4. Reagents and Materials Supplied

- 10 NADAL® Celiac Disease (tTG) test cassettes
- 10 Sterile lancets
- 10 Collection pipettes (20 µl)
- 1 Dilution buffer (3 ml)

### 5. Additional Required Materials

- Centrifuge (for serum only)
- Timer

### 6. Storage & Stability

Store the components of the **NADAL® Celiac Disease (tTG) Test** kit cooled at 2°C to 8°C in the original packaging. The expiry date is indicated on the pouch. Do not use the components beyond the expiry date!

### 7. Warnings and Precautions

- For professional *in vitro* diagnostic only
- Patient samples may contain infectious agents and should be treated and discarded as potentially dangerous biological materials.
- Do not exchange components from kits with different lot numbers.
- Allow all kit components and samples to reach room temperature as cold reagents and samples can reduce the functionality of the test. About 20-30 minutes are usually sufficient for reaching room temperature.
- Do not use kit components beyond the expiration dates.
- If the package is broken, the product can still be used as long as none of the components have been damaged.
- It is important to add the correct amount of blood in each test. If less than indicated is used, the chromatography may

not work because the sample may not reach the reaction area; if there is too much, the reagent may be diluted and produce a weak line.


- The used product should be discarded as indicated by current legislation.
- Do not use the test if a coloured line appears in the result area before you start to use it.

### 8. Specimen Collection and Preparation

#### Sample Collection

- The test was designed for use with fresh, non-hemolyzed blood samples to be collected with the kit collection pipettes.
- The blood should be taken from the fingertip of the index finger, middle finger or ring finger.
- Clean the fingertip with alcohol and let it air dry.
- Push the yellow protector of the sterile lancet towards the lancet until you hear a "click" (1).
- Twist the protector to break it and take it off to free the lancet needle (2).
- Turn the hand palm up, place the open end of the lancet on the cleaned finger and press the activating pad (3).
- Discard the lancet in an appropriate container for sharp objects.

(1)

(2)

(3)

#### Specimen Storage

Serum samples can be stored in the fridge for one day.



### 9. Procedure of the Test


#### I. Preparation

Bring all components of the **NADAL® Celiac Disease (tTG) Test** kit to room temperature before testing. Take the test cassette out of the sealed pouch and place it on a clean and even surface. Take a blood sample as described in chapter 8 "sample collection".

## II. Test Procedure

### Whole blood with use of the collection pipette

1. Hold the finger lower than the elbow and rub the hand from wrist to palm to finger to form a round drop on the puncture site.
2. Wipe away the first drop of blood with a sterile cloth or cotton.
3. Collect the blood sample with a new sterile collection pipette (20 µl) by touching the blood drop with the tip of the collection pipette. Keep the pipette horizontally and don't close the ventilation hole. **Do not press the central body of the pipette to collect a blood sample!**

4. Fill the collection pipette to the line. Avoid air inclusion!
5. Immediately place the blood sample in the round sample well of the test cassette to prevent coagulation in the pipette. Position the collection pipette carefully in the center of the sample well (marked with an arrow), then cover the pipette ventilation hole with a finger and gently press the upper part of the pipette to release the blood sample.
 

6. Wait for **30 to 60 seconds** until the blood is absorbed, then add **two drops (100 µl) of dilution buffer** to the same window (hold the dropper bottle vertically)
 

7. Start the timer and read the result after **10 minutes**.

### Whole blood without use of the collection pipette

1. Put **20 µl** blood into the round sample well (marked with an arrow).
2. Wait for **30 to 60 seconds** until the blood is absorbed, then add **two drops (100 µl) of dilution buffer** to the same window (hold the dropper bottle vertically)
3. Start the timer and read the result after **10 minutes**.

### Serum sample

1. Prepare a **1:10 dilution of serum with the dilution buffer** (for example 15 µl serum + 135 µl dilution buffer) and mix thoroughly by repeated pipetting.
2. Add **125 µl of the serum dilution** drop by drop to the round sample well (marked with an arrow).
3. Start the timer and read the result after **10 minutes**.

## 10. Interpretation of the Results

### Positive result

#### Two distinct lines appear.

A blue line appears in the control region (C) and a red line in the test region (T). The colour intensity of the lines may vary due to the antibody concentration.

### Negative result

#### Only the blue line in the control region (C) appears.


There is no visible red line in the test region (T).

### Invalid result

#### No blue control line (C) appears.

In this case the result is invalid even if the red line appears in the test region (T). Deviations from the correct test procedure or an insufficient sample volume can be the cause. Please check the test procedure for possible errors and repeat the test with a new test cassette. If the problem continues, please contact us.

**Note:** Every line that appears after 10 minutes has no diagnostic value anymore.

**The final diagnosis should be based on more than this test result. You should be taking into account other test data and clinical symptoms of the patient.**

## 11. Quality Control

There is an internal control included in the **NADAL® Celiac Disease (tTG) Test**. The blue line which appears in the control line region (C) is considered as the internal procedural control. It confirms sufficient specimen volume and that the chromatography was completed. Control standards are not supplied with the **NADAL® Celiac Disease (tTG) Test kit**; yet, it is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

## 12. Limitations

- This test provides a qualitative result only. Don't do any quantitative interpretation due to the colour intensity of the test line.
- Test results are to be interpreted together with clinical symptoms of the patient and other test data available.
- To perform a correct testing it is important to use the correct amount of blood and to check the reaction time to don't receive any false results.

13. Performance Characteristics

Sensitivity and Specificity

**Evaluation 1:** A prospective study was carried out under the direction of Dr Carmen Ribes Koninckx, including all the patients on whom an intestinal biopsy was performed for diagnostic purposes in various units of Paediatric Gastroenterology. Anti-TG antibodies were determined simultaneously by ELISA/EIA and by **NADAL® Celiac Disease (tTG) Test**.

A total of 150 patients received analyses, 26 with normal histology (controls), plus 122 with Marsh 3 lesions and 2 with Marsh 2 lesions (patients having celiac disease). The blood was drawn by finger prick and used immediately.

The results are presented in the table below:

Patients with CD (n = 124)			
NADAL® Celiac Disease (tTG)/ELISA Concordance = 93.15%			
		NADAL® Celiac Disease (tTG)	
		Positive (110)	Negative (14)
ELISA	Positive (114)	108	6
	Negative (10)	2*	8
Control Patients (n = 26)			
NADAL® Celiac Disease (tTG)/ELISA Concordance = 80%			
		NADAL® Celiac Disease (tTG)	
		Positive (2)	Negative (24)
ELISA	Positive (5)	1	4
	Negative (21)	1	20

\*IgA Deficit

The efficacy of the different methods in diagnosing CD is presented in the table below:

	Sensitivity	Specificity	Efficacy
NADAL® Celiac Disease (tTG)	89%	92.5%	89.3%
IgA-tTG (EIA)	91.9%	80%	90%

**Evaluation 2:** A transversal population study was carried out on children from 5 to 9 years old, not known to have CD, who visited their pediatrician for any reasons. Under the direction of Dr Carmen Ribes Koninckx, the study was carried out by 34 pediatricians from 11 different centers in the Valencia metropolitan area.

After informed consent was received, some drops of blood were drawn by capillary puncture to determine anti-tTG

antibodies using the **NADAL® Celiac Disease (tTG) Test**. If the test was positive, the result was checked by ELIA-ImmunoCAP 100 (Phadia) after drawing blood by venipuncture.

A total of 850 tests using **NADAL® Celiac Disease (tTG) Test** were evaluated, obtaining 7 positive results:

- Four strong positives with values >100 U/ml in the Phadia test and with intestinal biopsies showing histological alterations compatible with CD.
- Three weak positives: Two results with Phadia test values of 38 U/ml and 36 U/ml. The first of these two cases presented very strong elevation (>100 U/ml) in a control 2 months later; an intestinal biopsy was therefore performed, the results also being compatible with CD. The second case has a first-degree relative diagnosed with CD and is pending developmental serologic control. In the third weak positive case, the venous blood determination was negative.

In the study sample (n=850), with 5 confirmed cases, a CD incidence of 1/170 was found; this could rise to 1/141 if the patient pending serologic control turns out to have the disease. In the population of 10,647 children of the same age range controlled by the study participants, there were 23 known cases of celiac disease; CD incidence was therefore 1/463.

For each case of previously-known celiac disease in this age range, 2.7 cases of silent/atypical CD have been detected. Extrapolating the results obtained (in the sample: 5/850 = 1/170; in the total non-CD population: 10624/170 = 62) and adding the known cases (23), we obtain a total prevalence of 85 cases in this population. This yields a CD incidence for this age of 1/125.

The results obtained in this study demonstrate that the **NADAL® Celiac Disease (tTG) Test** is highly useful for possible CD screening programs.

**Evaluation 3:** This evaluation covered 101 blood samples (analyzed with **NADAL® Celiac Disease (tTG) Test**) and the respective sera samples (Alphadia Immunodot Gliadin-tTG Dot Kit) from CD patients, relatives of patients with CD and patients with gastrointestinal symptoms (diarrhea, malnutrition, abdominal pain and so on).

The results are presented in the table below:

		Alphadia ImmunoDot	
		Positive (11)	Negative (90)
NADAL® Celiac Disease (tTG)	Positive (12)	11	1
	Negative (89)	0	89

Concordance percent: 99%

Sensitivity concordance: >99%

Specificity concordance: 98.9%

**Evaluation 4:** The **NADAL® Celiac Disease (tTG) Test** sensitivity and specificity were evaluated by analyzing 49 blood samples: 18 patients suspected of having CD, 31 negative patients (with CD and on gluten-free diet, or with other pathologies).

Organized by Eurospital, the study was carried out in a benchmark CD hospital in Italy.

The results are presented in the table below:

		Eurospital	
		Positive (18)	Negative (31)
<b>NADAL® Celiac Disease (tTG)</b>	Positive (17)	17	0
	Negative (32)	1	31
Sensitivity: 94.4%			
Specificity: 100%			
Positive Predictive Value: 100%			
Negative Predictive Value: 96.9%			

### Precision

#### Intra-Assay Precision

With one lot of **NADAL® Celiac Disease (tTG) Test**, 10 sensitivity curves (serial dilutions of a positive control serum) were performed. The analyses were performed by the same person on the same day. The results obtained in all the curves were similar, which demonstrates the high intra-assay precision of the **NADAL® Celiac Disease (tTG) Test**.

#### Inter-Day Precision

With one lot of **NADAL® Celiac Disease (tTG) Test**, a sensitivity curve (serial dilutions with a positive control serum) was performed for 10 consecutive days. The analyses were carried out by the same person, obtaining similar results each day, which demonstrates the high inter-day precision of the **NADAL® Celiac Disease (tTG) Test**.

#### Inter-Laboratory Precision

With one lot of **NADAL® Celiac Disease (tTG) Test**, 3 sensitivity curves (serial dilutions of a positive control serum) were performed. The analyses were performed by three technicians on the same day. The results obtained in the three cases were similar, which demonstrates the high inter-laboratory precision of the **NADAL® Celiac Disease (tTG) Test**.

#### Inter-Lot Precision

With three lots of **NADAL® Celiac Disease (tTG) Test**, a sensitivity curve (serial dilutions of a positive control serum) was performed. The analyses were carried out by the same technician, obtaining similar results with all three lots, which demonstrates the high inter-lot precision of the **NADAL® Celiac Disease (tTG) Test**.

### Hook Effect

The Hook effect is the effect produced when, in the presence of high concentrations of the analyte to detect, a decrease in the signal occurs, positive samples even turning negative.

Many positive samples (sera and blood) have been tested without observing a decrease in the intensity of the test line in case of higher concentrations of anti-tTG antibodies. Therefore, no Hook effect was observed.

### Cross reactions

With a lot of **NADAL® Celiac Disease (tTG) Tests**, several interference substances that are likely to be in blood were examined. The interference substances were added to the buffer in the final amount that is shown in the table and a sensitivity curve was performed. Two negative control sera were examined as well.

<b>Substance</b>	<b>Initial Concentration</b>	<b>Diluent</b>	<b>Concentration in buffer</b>
Aspirin	50 mg/ml	Ethanol	2,5 mg/ml
Caffeine	10 mg/ml	Water	0,5 mg/ml
Ibuprofen	8 mg/ml	Water	0,4 mg/ml
Ethanol	500 mg/ml	Water	50 mg/ml
Folic acid	0,4 mg/ml	Water	0,02 mg/ml
Nicotinic acid	20 mg/ml	Water	1 mg/ml
Haemoglobin	15 mg/ml	Water	0,75 mg/ml
hCG	250000 mUI/ml	Water	500 mUI/ml
EDTA	20 mg/ml	Water	1 mg/ml
Sodium citrate	0,2 M	Water	0,01M
Heparin	500 mg/ml	Water	50 mg/ml

None of the analyzed substances affected the proper performance of the **NADAL® Celiac Disease (tTG) Tests**.

### 14. References

- Zhang J. et al. "Modulation of the in situ activity of tissue transglutaminase by calcium and GTP". The Journal of Biological Chemistry (1998) 273 (4): 2288-2295.
- Sárdy M. et al. "Recombinant human tissue transglutaminase ELISA for the diagnosis of gluten-sensitive enteropathy". Clinical Chemistry (1999) 45 (12): 2142-2149.
- Sorell L. et al. "One-step immunochromatographic assay for screening of celiac disease". Lancet (2002) 359: 945-946.
- Drago S. et al. "Recent developments in the pathogenesis, diagnosis and treatment of celiac disease". Expert Opinion in Therapy Patents (2002) 12 (1): 45-51.
- Hansson T. et al. "Recombinant human tissue transglutaminase for diagnosis and follow-up of childhood coeliac disease". Pediatric Research (2002) 51 (6): 700-705.
- Ferre-López, S. et al. "A reliable generation of immunochromatographic sticks for both t-TG and AGA screening in Coeliac Disease". Clinical Gastroenterology and Hepatology (2004) 2 (6): 480-484.
- Jennings J. et al. "New developments in celiac disease". Current Opinion in Gastroenterology (2003) 19 (2): 118-129.

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## 14. Bibliographie

- 1) Zhang J. et al. "Modulation of the in situ activity of tissue transglutaminase by calcium and GTP". The Journal of Biological Chemistry (1998) 273 (4): 2288-2295.
- 2) Sárdy M. et al. "Recombinant human tissue transglutaminase ELISA for the diagnosis of gluten-sensitive enteropathy". Clinical Chemistry (1999) 45 (12): 2142-2149.
- 3) Sorell L. et al. "One-step immunochromatographic assay for screening of celiac disease". Lancet (2002) 359: 945-946.
- 4) Drago S. et al. "Recent developments in the pathogenesis, diagnosis and treatment of celiac disease". Expert Opinion in Therapy Patents (2002) 12 (1): 45-51.
- 5) Hansson T. et al. "Recombinant human tissue transglutaminase for diagnosis and follow-up of childhood coeliac disease". Pediatric Research (2002) 51 (6): 700-705.
- 6) Ferre-López, S. et al. "A reliable generation of immunochromatographic sticks for both t-TG and AGA screening in Coeliac Disease". Clinical Gastroenterology and Hepatology (in press)
- 7) Jennings J. et al. "New developments in celiac disease". Current Opinion in Gastroenterology (2003) 19 (2) : 117-129.

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14. Referencias

1. Zhang J. et al. "Modulation of the in situ activity of tissue transglutaminase by calcium and GTP". The Journal of Biological Chemistry (1998) 273 (4): 2288-2295.

2. Sárdy M. et al. "Recombinant human tissue transglutaminase ELISA for the diagnosis of gluten-sensitive enteropathy". Clinical Chemistry (1999) 45 (12): 2142-2149.

3. Sorrell L. et al. "One-step immunochromatographic assay for screening of celiac disease". Lancet (2002) 359: 945-946.

4. Drago S. et al. "Recent developments in the pathogenesis, diagnosis and treatment of celiac disease". Expert Opinion in Therapy Patents (2002) 12 (1): 45-51.

5. Hansson T. et al. "Recombinant human tissue transglutaminase for diagnosis and follow-up of childhood coeliac disease". Pediatric Research (2002) 51 (6): 700-705.

6. Ferre-López, S. et al. "A reliable generation of immunochromatographic sticks for both t-TG and AGA screening in Coeliac Disease". Clinical Gastroenterology and Hepatology (2004) 2 (6): 480-484.

7. Jennings J. et al. "New developments in celiac disease". Current Opinion in Gastroenterology (2003) 19 (2): 118-129.

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#### 14. Bibliografia

8. Zhang J. et al. "Modulation of the in situ activity of tissue transglutaminase by calcium and GTP". The Journal of Biological Chemistry (1998) 273 (4): 2288-2295.
9. Sárdy M. et al. "Recombinant human tissue transglutaminase ELISA for the diagnosis of gluten-sensitive enteropathy". Clinical Chemistry (1999) 45 (12): 2142-2149.
10. Sorell L. et al. "One-step immunochromatographic assay for screening of celiac disease". Lancet (2002) 359: 945-946.
11. Drago S. et al. "Recent developments in the pathogenesis, diagnosis and treatment of celiac disease". Expert Opinion in Therapy Patents (2002) 12 (1): 45-51.
12. Hansson T. et al. "Recombinant human tissue transglutaminase for diagnosis and follow-up of childhood coeliac disease". Pediatric Research (2002) 51 (6): 700-705.
13. Ferre-López, S. et al. "A reliable generation of immunochromatographic sticks for both t-TG and AGA screening in Coeliac Disease". Clinical Gastroenterology and Hepatology (2004) 2 (6): 480-484.
14. Jennings J. et al. "New developments in celiac disease". Current Opinion in Gastroenterology (2003) 19 (2): 118-129.

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nal von minden GmbH