FR

Nom du produit (et toxiques détectés)	Références
DRUGCHECK® 11 (AMP, BUP, BAR, BZO, COC, MTD, MET, MDMA, MOP, THC, TCA)	1140016
DRUGCHECK® 10 (AMP, BAR, BZO, COC, MTD, MET, MDMA, MOP, THC, TCA)	1140015
DRUGCHECK® 5 (AMP, COC, MDMA, MOP, THC)	1140014
DRUGCHECK® AMP	1140002
DRUGCHECK® BAR	1140003
DRUGCHECK® BUP	1140004
DRUGCHECK® BZO	1140005
DRUGCHECK® COC	1140006
DRUGCHECK® COT	1140001
DRUGCHECK® MTD	1140010
DRUGCHECK® MET	1140008
DRUGCHECK® MDMA	1140007
DRUGCHECK® MOP	1140009
DRUGCHECK® PCP	1140011
DRUGCHECK® THC	1140013
DRUGCHECK® TCA	1140012

1. OBJECTIF

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14. BIBLIOGRAPHIE

- 1) Tietz NW. Textbook of Clinical Chemistry. W.B. Saunders Company. 1986; 1735
- 2) Baselt RC. Disposition of Toxic Drugs and Chemicals in Man. 2nd Ed. Biomedical Publ., Davis, CA. 1982; 488
- 3) Hawks RL, CN Chiang. Urine Testing for Drugs of Abuse. National Institute for Drug Abuse (NIDA), Research Monograph 73, 1986

EN

Name of the product (and toxic detected)	Reference
DRUGCHECK® 11 (AMP, BUP, BAR, BZO, COC, MTD, MET, MDMA, MOP, THC, TCA)	1140016
DRUGCHECK® 10 (AMP, BAR, BZO, COC, MTD, MET, MDMA, MOP, THC, TCA)	1140015
DRUGCHECK® 5 (AMP, COC, MDMA, MOP, THC)	1140014
DRUGCHECK® AMP	1140002
DRUGCHECK® BAR	1140003
DRUGCHECK® BUP	1140004
DRUGCHECK® BZO	1140005
DRUGCHECK® COC	1140006
DRUGCHECK® COT	1140001
DRUGCHECK® MTD	1140010
DRUGCHECK® MET	1140008
DRUGCHECK® MDMA	1140007
DRUGCHECK® MOP	1140009
DRUGCHECK® PCP	1140011
DRUGCHECK® THC	1140013
DRUGCHECK® TCA	1140012

1. INTENDED USE

The DRUGCHECK® tests are lateral flow chromatographic immunoassays for the qualitative and simultaneous detection of multiple drugs and metabolites in urine at the following cut-off concentrations in urine:

Drug	Calibrator	Cut-off (ng/ml)
Amphetamine (AMP)	d-Amphétamine	1 000
Barbiturates (BAR)	Secobarbital	300
Benzodiazépines (BZO)	Oxazepam	300
Buprénorphine (BUP)	Buprénorphine	10
Cocaine (COC)	Benzoylcegonine	300
Cotinine (COT)	Cotinine	100
Marijuana (THC)	11-nor- Δ^9 -THC-9-COOH	50
Methadone (MTD)	Méthadone	300
Methamphetamine (MET)	d-Méthamphétamine	1 000
Methylenedioxymethamphetamine (MDMA)	d,l-Methylenedioxymethamphetamine	500
Morphine (MOP)	Morphine	300
Phencyclidine (PCP)	Phencyclidine	25
Tricyclic Antidepressants (TCA)	Nortriptyline	1 000

These tests detect also other related compounds (refer to chapter 13.3. Analytical specificity). The test is designed for professional in vitro diagnostic use only.

Urine based tests for multiple drugs of abuse range from simple immunoassay tests to complex analytical procedures. The speed and sensitivity of immunoassays have made them the most widely accepted method to screen urine for multiple drugs of abuse.

2. SUMMARY

The DRUGCHECK® test provide only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method.

Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used.

3. TEST PRINCIPLE

The DRUGCHECK® test are immunoassay based on the principle of competitive binding. It contains mouse monoclonal antibody-coupled particles and corresponding drug-protein conjugates. A goat antibody is employed in each control line. Drugs which may be present in the urine specimen compete against their respective drug conjugate for binding sites on their specific antibody.

During testing, a urine specimen migrates upward by capillary action. A drug, if present in the urine specimen below its cut-off concentration, will not saturate the binding sites of its specific antibody. The antibody will then react with the drug-protein conjugate and a visible colored line will show up in the test line region of the specific drug strip. The presence of drug above the cut-off concentration will saturate all the binding sites of the antibody. Therefore, the colored line will not form in the test line region.

A drug-positive urine specimen will not generate a colored line in the specific test line region of the strip because of drug competition, while a drug-negative urine specimen will generate a line in the test line region because of the absence of drug competition.

To serve as a procedural control, a colored line will always appear at the control line region, indicating that proper volume of specimen has been added and membrane wicking has occurred.

4. MATERIAL PROVIDED

DRUGCHECK® kit contains individually packed test devices (cassettes) each contains :

- 1 Device,
- 1 Pipette,
- 1 Desiccant,
- 1 bilingual instructions for use (French-English).

5. MATERIAL REQUIRED BUT NOT PROVIDED

- Timer with alarm.
- Specimen collection container.
- Positive and negative controls.

6. STORAGE AND STABILITY

- Store the kit in a dry area, at a temperature between 2°C and 30°C. The kit can be stored without the life and qualities of the product are altered, 12 days at 37°C, 8 days at 40°C, 5 days at 45°C. Do not freeze the kit.
- Use the kit before the expiration date.
- After opening the pouch, the device must be immediately used within one hour maximum.
- Take care to protect the components of this kit from contamination. Do not use if there is evidence of microbial contamination or precipitation. Biological contamination of dispensing equipment, containers or reagents can lead to false results.

7. PRECAUTIONS

- For *in vitro* diagnostic use by professionals only.
- Follow the product insert instruction carefully.
- Do not use after expiration date.
- For single use. Do not reuse tests.
- Do not use foreign reagents.
- Do not use test if its foil pouch has been damaged, due to humidity sensitive of the test.
- Handle all specimens as if they contain infectious agents. Do not eat, drink or smoke in the area where the specimens or kits are handled. Protective clothing such as laboratory coats, disposable gloves and eye protection are recommended.
- Observe established precautions against microbiological hazards throughout testing and follow standard procedures for proper disposal of specimens in accordance with local regulations.
- All materials used to perform the test should be treated in order to ensure the inactivation of pathogenic agents.
- Avoid splash and aerosol formation.
- Do not spill the specimens into the reaction zone.
- Bring all reagents to room temperature (15-30°C) before use.
- Humidity and temperature can adversely affect results.
- The components of the test (e.g. antibodies/chemicals) do not cause any danger if the test is used according to the instructions.

8. SPECIMEN COLLECTION AND STORAGE

- The urine specimen must be collected in a clean and dry container. Urine collected at any time of the day may be used. Urine specimens exhibiting visible precipitates should be centrifuged, filtered, or allowed to settle to obtain a clear supernatant for testing.
- Urine specimens may be stored for up to 24 hours at room temperature or for up to 48 hours at 2°C and 8°C. For prolonged storage, specimens may be stored below -20°C. Frozen specimens should be thawed and mixed well before testing.
- If specimens are to be shipped, pack them in compliance with all applicable regulations for transportation of etiologic agents.

9. PROCEDURE

1. Bring the kit and the samples to room temperature at least 30 minutes before testing.
2. Open the pouch of the device. Put the device on a flat surface. Once the pouch of the device is open, perform the test within one hour.

<p>PANEL TEST : Exp : DRUGCHECK 11</p> 	<p>SINGLE TEST DRUGCHECK</p> 	<p>3. Hold the pipette vertically and dispose 3 drops of urine (100 µl) in each specimen well of the device.</p> <p>Avoid trapping air bubbles in the specimen well.</p> <p>4. Read the result after 5 minutes. Do not interpret after 10 minutes.</p>
---	---	---

10. INTERPRETATION OF RESULTS

	<p>C control band and No T test band</p>	<p>The test is positive.</p> <p>Presence of the designated drug in sample at a concentration exceeding the cut-off level for that drug.</p>
	<p>C control band and T test band (even if weak intensity)</p>	<p>The test is negative.</p> <p>Absence of the designated drug in sample or presence at a concentration below the cut-off level for that drug.</p>
	<p>No C control band or No bands</p>	<p>The test is invalid: run the procedure again with another device.</p> <p>Insufficient specimen volume or incorrect procedure.</p>

11. QUALITY CONTROL

- An internal procedural control is included in the test. A coloured line appearing in the control line region (C) is an internal positive procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.
- External controls are not supplied with this kit. It is recommended to control on each new batch, or on each new delivery. Each laboratory has to set up its own planning of controls. (Reference controls available : Positive control (12 parameters) : 19001066 Negative control (13 parameters) : 19700000, Cotinine positive control : 19000102).

12. LIMITATIONS

- For the reliability of results, respect precautions strictly and follow the procedure. Refer to the results interpretation of instructions for use.
- As for each diagnostic test, confront results with clinical findings.
- Confirm the result of the test (qualitative, preliminary analytical result) with the results of a secondary analytical method. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method.
- Notice that there is a possibility that technical / procedural errors, as well as other interfering substances in the urine specimen may cause erroneous results.
- Consider that adulterants (bleach and/or alum) in urine specimens may produce erroneous results regardless of the analytical method used. If adulteration is suspected, the test should be repeated with another urine specimen.
- Keep in mind that a positive result does not indicate intoxication level, administration route or concentration in urine and that a negative result may not necessarily indicate drug-free urine. Indeed negative results can be obtained when drug is present but below the cut-off level of the test.
- Notice that the test does not distinguish between drugs of abuse and certain medications.
- Consider that a positive result might be obtained from certain foods or food supplements.
- The drug Subutex® (Buprenorphine) does not interfere with the MOP parameter.
- The drugs Subutex® (Buprenorphine), Prozac® (Fluoxetine) and Imovane® (Zopiclone) do not interfere neither with the BZO parameter nor with the TCA parameter. A false positive result can occur if the patient is treated simultaneously with the three of them.
- The drugs Zantac® and Zinetac® (Ranitidine) can interfere with the MET parameter and give false positive results.

13. PERFORMANCES

13.1. ACCURACY

A side-by-side comparison was conducted on 300 specimens (subjects presenting for Drug Screen Testing) using the DRUGCHECK® panel and single test and another available drug rapid test on urine. Presumptive positive results and 10% of negative specimens were confirmed by GC/MS.

• Percentage of agreement with the other test

Specimen panel / single	AMP	BAR	BZO	BUP*	COC	COT	THC
Positive	>99% / 96%	98% / >99%	98% / 90%	88%	>99%	>99%	>99%
Negative	>99%	>99% / 99%	>99% / 97%	>99%	99%	>99%	>99%
Total	>99% / 98%	99%	99% / 94%	97%	99%	>99%	>99%

* BUP test was compared to the self-reported use of Buprenorphine

Specimen panel / single	MTD	MET	MDMA	MOP	PCP	TCA
Positive	87% / >99%	>99% / 99%	98% / >99%	94% / >99%	>99% / 97%	92% / 95%
Negative	>99%	>99%	>99%	>99%	>99%	>99% / 96%
Total	94% / >99%	>99%	99%	97% / >99%	>99% / 99%	98% / 99%

• Percentage of agreement with the GC/MS method

Specimen panel / single	AMP	BAR	BZO	BUP*	COC	COT	THC
Positive	94% / 96%	92%	98% / 96%	98%	95% / 96%	>99%	95% / 97%
Negative	99% / 95%	99% / 98%	98% / 96%	99% / >99%	>99% / 90%	>99%	96%
Total	97% / 95%	96% / 95%	98% / 96%	99% / >99%	98% / 93%	>99%	95% / 97%

* BUP test was based on LC/MS data (liquid chromatography with mass spectrometry)

Specimen panel / single	MTD	MET	MDMA	MOP	PCP	TCA**
Positive	93% / 99%	90% / 99%	98% / >99%	98% / >99%	90% / >99%	>99%
Negative	>99% / 94%	>99% / 93%	98%	97% / 94%	99% / 97%	94% / 89%
Total	97% / 96%	95% / 96%	98% / 99%	97%	96% / 98%	95% / 91%

** TCA test was based on HPLC method (high performance liquid chromatography)

13.2. ANALYTICAL SENSITIVITY

A drug-free urine pool was spiked with drugs to the concentrations at $\pm 50\%$ cut-off and $\pm 25\%$ cut-off.

Concentration panel / single	AMP		BAR		BZO		BUP		COC		COT		THC	
	-	+	-	+	-	+	-	+	-	+	-	+	-	+
0% Cut-off	30	0	30	0	30	0	90	0	30	0	90	0	30	0
-50% Cut-off	30	0	30	0	30	0	90	0	30	0	90	0	30	0
-25% Cut-off	26 / 23	4 / 7	23 / 20	7 / 10	24 / 26	6 / 4	78 / 79	12 / 11	25 / 30	5 / 0	90	0	24 / 30	6 / 0
Cut-off	23 / 9	7 / 21	14 / 13	16 / 17	15 / 12	15 / 18	48 / 49	42 / 41	20 / 9	10 / 21	49	41	15 / 21	15 / 9
+25% Cut-off	7 / 1	23 / 1	7 / 8	23 / 22	6 / 3	24 / 27	24 / 10	66 / 80	5 / 7	25 / 23	4	86	6 / 17	24 / 13
+50% Cut-off	0	30	0	30	0	30	0	90	0	30	0	90	0	30

Concentration panel / single	MDMA		MOP		PCP		TCA		MTD		MET	
	-	+	-	+	-	+	-	+	-	+	-	+
0% Cut-off	30	0	30	0	30	0	30	0	30	0	30	0
-50% Cut-off	30	0	30	0	30	0	30	0	30	0	30	0
-25% Cut-off	27 / 24	3 / 7	20 / 28	10 / 2	26	4	25 / 26	5 / 4	26	4	25 / 24	5 / 6
Cut-off	17 / 18	13 / 15	18 / 20	12 / 10	19 / 11	11 / 19	13 / 14	17 / 16	18 / 16	17 / 14	23 / 18	7 / 12
+25% Cut-off	6 / 1	24	7 / 3	23 / 27	5 / 8	25 / 22	7 / 4	23 / 26	5 / 4	25 / 26	6 / 1	24 / 29
+50% Cut-off	0	30	0	30	0	30	0	30	0	30	0	30

13.3. ANALYTICAL SPECIFICITY

Here under the concentrations (ng/mL) detected positive in urine by the DRUGCHECK® test at 5 minutes.

AMPHETAMINE - AMP		MARIJUANA - THC	
d-Amphetamine	1.000	11-nor- Δ^9 -THC-9 COOH	50
d,l-Amphetamine	3.000	Cannabinol	20.000
l-Amphetamine	50.000	11-nor- Δ^8 -THC-9 COOH	30
Phentermine	3.000	Δ^8 -THC	15.000
3,4- Methylendioxyamphetamine	2.000	Δ^9 -THC	15.000
BARBITURATES - BAR		METHADONE - MTD	
Secobarbital	300	Methadone	300
Amobarbital	300	Doxylamine	50.000
Alphenol	150	METHAMPHETAMINE - MET	
Aprobarbital	200	d-Methamphetamine	1.000
Butabarbital	75	p-Hydroxymethamphetamine	30.000
Butethal	100	l-Methamphetamine	8.000
Butalbital	2.500	Mephentermine	50.000
Cyclopentobarbital	600	3,4-Methylenedioxy methamphetamine (MDMA)	2.000
Pentobarbital	300	METHYLENEDIOXYMETHAMPHETAMINE MDMA	
Phenobarbital	100	3,4-Methylenedioxy methamphetamine (MDMA)	500
BENZODIAZEPINES - BZO		3,4-Methylenedioxy amphetamine	3.000
Oxazepam	300	3,4-Methylenedioxy ethylamphetamine (MDEA)	300
Alprazolam	196	MORPHINE/OPIACE - MOP	
α -Hydroxylalprazolam	1.262	Morphine	300
Bromazepam	1.562	Codeine	300
Chlordiazepoxide	1.562	Ethylmorphine	6.250
Clonazepam	781	Hydrocodone	50.000
Clobazam	98	Hydromorphone	3.125
Clonazepam	781	Levorphanol	1.500
Clorazepate	195	6-Monoacetylmorphine	400
Delorazepam	1.562	Morphine 3- β -D-glucuronide	1.000
Desalkylflurazepam	390	Norcodeine	6.250
Diazepam	195	Normorphine	100.000
Estazolam	2.500	Oxycodone	30.000
Flunitrazepam	390	Oxymorphone	100.000
d,l Lorazepam	1.562	Procaine	15.000
RS-Lorazepam glucuronide	1.562	Thebaine	6.250
Midazolam	12.500	PHENCYCLIDINE - PCP	
Nitrazepam	98	Phencyclidine	25
Norchlordiazepoxide	195	4-Hydroxyphencyclidine	12.500
Nordiazepam	390	TRICYCLIC ANTIDEPRESSANTS - TCA	
Temazepam	98	Nortriptyline	1.000
Triazolam	2.500	Nordoxepin	1.000
BUPRENORPHINE - BUP		Trimipramine	3.000
Buprenorphine	10	Amitriptyline	1.500
Norbuprenorphine	20	Promazine	1.500
Buprenorphine 3-D-glucuronide	15	Desipramine	200
Norbuprenorphine 3-D-glucuronide	200	Imipramine	400
COCAINE - COC		Clomipramine	12.500
Benzoylcegonine	300	Doxepin	2.000
Cocaine	780	Maprotiline	2.000

Cocaethylene	12.500	Promethazine	25.000
Ecgonine	32.000	COTININE - COT	
		l-Cotinine	100
		S-l-Nicotine	12 500

13.4. CROSS REACTIVITY

A study was conducted to determine the cross-reactivity of the test with compounds in negative urine or containing: Amphetamine, Barbiturates, Benzodiazepines, Buprenorphine, Cocaine, Marijuana, Methadone, Methamphetamine, Methylenedioxyamphetamine, Morphine, Phencyclidine or Tricyclic Antidepressants.

Here under the different compounds showing no cross-reactivity when tested with the DRUGCHECK® test at a concentration of 100 μ g/ml.

Acetophenetidin	l-Cotinine	Ketamine	d-Pseudoephedrine
N-Acetylprocainamide	Creatinine	Ketoprofen	Quinidine
Acetylsalicylic acid	Deoxy corticosterone	Labetalol	Quinine
Aminopyrine	Dextromethorphan	Loperamide	Salicylic acid
Amoxicillin	Diclofenac	Meprobamate	Serotonin
Ampicillin	Diflunisal	Methoxyphenamine	Sulfamethazine
l-Ascorbic acid	Digoxin	Methylphenidate	Sulindac
Apomorphine	Diphenhydramine	Nalidixic acid	Tetracycline
Aspartame	Ethyl-p-aminobenzoate	Naproxen	Tetrahy drocortisone,
Atropine	β -Estradiol	Niacinamide	3-Acetate
Benzilic acid	Estrone-3-sulfate	Nifedipine	Cortisone
Benzoic acid	Erythromycin	Norethindrone	Tetrahydrozoline
Bilirubin	Fenoprofen	Noscapine	Thiamine
d,l-Brompheniramine	Furosemide	d,l-Octopamine	Thioridazine
Caffeine	Gentisic acid	Oxalic acid	d,l-Tyrosine
Cannabidiol	Hemoglobin	Oxolinic acid	Tolbutamide
Chloral hydrate	Hydralazine	Oxymetazoline	Triamterene
Chloramphenicol	Hydrochloro thiazide	Papaverine	Trifluoperazine
Chlorothiazide	Hydrocortisone	Penicillin-G	Trimethoprim
d,l-Chlorpheniramine	o-Hydroxyhippuric acid	Perphenazine	d,l-Tryptophan
Chlorpromazine	3-Hydroxytyramine	Phenelzine	Uric acid
Cholesterol	d,l-Isoproterenol	Prednisone	Verapamil
Clonidine	Isoxsuprine	d,l-Propanolol	Zomepirac

14. BIBLIOGRAPHY

- 1) Tietz NW. Textbook of Clinical Chemistry. W.B. Saunders Company. 1986; 1735
- 2) Baselt RC. Disposition of Toxic Drugs and Chemicals in Man. 2nd Ed. Biomedical Publ., Davis, CA. 1982; 488
- 3) Hawks RL, CN Chiang. Urine Testing for Drugs of Abuse. National Institute for Drug Abuse (NIDA), Research Monograph 73, 1986

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