PT (INR) test kit

for quantitative in vitro determination of prothrombin time on smart or CUBE laboratory photometer.

English

Order Information
Order number: ST0180
Order number: ST1800

Indication
PT (INR) test kit
PT (INR) control kit

Kit size
32 tests/kit
2 x 1 ml / N and AK

Test kit preparation: Allow single test min. 30 minutes to warm to room temperature by placing the test into the test kit rack. Put test kit package back into refrigerator.

Summary

PT (INR) test kit is used for the determination of prothrombin complex-activity in whole blood according to Owren. It is ideal for monitoring of oral anticoagulation with vitamin-K inhibitors like Warfarin or Marcumar. PT according to Owren measures the combined activity of the vitamin-K dependent coagulation factors II, VII and X. In analysis, the sample is mixed with thromboplastin, CaCl₂, fibrinogen and coagulation factor V. The clotting time is prolonged on the therapeutical level (<0.4 IU/mL) of factor II, V and X in the sample, but the assay is insensitive to heparin in therapeutical levels up to 1.0 IE/mL.

Assay principle

Photometrical clot detection.

Measurement Range: 1 - 8 INR
ISI: ~ 1,00 (Lot pending)

Testkit

ERS Cupette not pre-filled - see test preparation.
ERS Cap filled with buffer

Stability and Storage

Unopened test kit:
at 2 - 8°C: Until the expiration date stated on the label.
at 20 - 25°C: 14 days

R1 reconstituted (100µl prefilled in cuvette):
at 2 - 8°C: 45 days
at 20 - 25°C: 7 days

Warning

This testkit is for in vitro diagnostic use only. Do NOT INGEST. Avoid contact with skin and eyes. Contains sodium azide, which may react with lead or copper plumbin to form explosive compounds. Take the necessary precautions for the use of laboratory reagents.

Waste Management

Please refer to local legal requirements.

Sample Material

20µl fresh capillary finger blood OR ...
20µl well mixed citrated whole blood* (9+1 sodium citrate 3.2%) - not older than 4 hours! OR ...
20µl citrated plasma

* values may be slightly increased (~10%) due to the dilution factor of the sodium citrate added to the blood collection system.

Reference Ranges

Typical oral anticoagulant therapeutic range:
2,0 - 3,5 INR

Control Material

For internal quality assurance the PT (INR) control kit should be used; Order number: ST1800
Please make sure "Ctrl. Plasma" is chosen in the sample type menu when performing a measurement of a control.

EQUAS

For EQUAS procedures we solely can recommend Plasma Control from Technoclone®. Ask your EQUAS institute to source the right controls.

Precision for Plasma Controls

CV (INR 0,93) < 5,75 %
CV (INR 2,27) < 4,30 %

Precision for Citrated Whole Blood

CV (INR 1,31) < 7,05 %
CV (INR 2,41) < 1,28 %

Correlation for Citrated Plasma

N=39; R²=0,9337; y=1,0232x + 0,1333; x=ACL Top300 plasma INR; y=Eurolyser plasma INR

Correlation for Citrated Whole Blood

N=49; R²=0,9334; y=0,9273x + 0,174; x=ACL Top300 plasma INR; y=Eurolyser whole blood INR

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English

1. Reagent preparation

- 1.1 Open Reagent 1 (R1) and Diluent (D).
- 1.2 Pour Diluent into Reagent 1.
- 1.3 Close Reagent 1.
- 1.4 Gently rotate mixture and allow to warm to room temperature for 15 minutes. Meanwhile repeat mixing every 5min to achieve good mixture.

2. ERS Cuvette Preparation (Repeat steps 16x)

- 2.1 Take out cuvette from testkit box and remove cuvette cap.
- 2.2 Aspirate 100 µl from Reagent 1.
- 2.3 Dispense 100 µl into ERS cuvette. Close cuvette by use of cuvette cap.
- 2.4 Put back cuvette into testkit box.





3. Post-Processing

- 3.1 Insert date of mixing into the provided label on the test kit box.
- 3.2 Dispose Reagent 1 (remains) and Diluent vial.
- 3.3 Store filled cuvettes upright only (cap to top)

Note!

Testkit is stable for 45 days (stored at 2 - 8°C) or 7 days (stored at 20 - 25°C) after reconstitution.

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English

ATTENTION!

Allow single test min. 30 minutes to warm to room temperature!

Test can be processed on instruments with firmware version 1.33 or higher only.

Correct use of the capillary

- A Insert the plunger into the capillary opening at the wide, black marker. Then, hold capillary in a sloping position to blood drop until it is filled to the small, red marker completely.
- Aspirate**
- B WRONG: Capillary is not completely filled with blood to the small, red marker. Hold the capillary horizontally when aspirating to ensure it gets filled completely.
- C CORRECT: Capillary is filled with blood to the small, red marker with completely.
- Dispense**
- D WRONG: Plunger has not been pushed all the way into the capillary. Blood is left in capillary.
- E CORRECT: Plunger has been pushed completely into the capillary. White stopper is pressed to the end of the capillary. No blood is left.

1. Preparation of test system

- 1.1 Place RFID card
- 1.2 Place ERS cuvette in test kit rack
- 1.3 Place ERS cap in test kit rack
- 1.4 Press symbol on touch screen, required information is entered by use of the touch display of the analyser

2. Sample preparation for whole blood

- 2.1 Choose sampletype: Please make sure "Capil. Blood" is selected when performing a measurement from capillary finger blood or from citrated whole blood. When performing a measurement from citrated plasma or a control measurement please make sure "Ctrl. Plasma" is selected.
- 2.2 Aspirate 20 µl whole blood from fingertip (See "Correct use of capillary"!)
- 2.3 Attention: Entire reagent liquid has to be at the bottom of the cuvette. By pressing the plunger: Dispense 20 µl whole blood into ERS cuvette INTO LIQUID!

Carry out step 3 immediately!

3. Test preparation

- 3.1 Apply ERS Cap
- 3.2 Place ERS cartridge into laboratory photometer
- 3.3 Start automatic sample processing by pressing the [Start] button on the smart laboratory photometer or by closing the door of the CUBE laboratory photometer.

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