Troponinas T STAT

cobas®

Troponin T STAT (Short Turn Around Time), CARDIAC T (skubiai atliekamas)

04660307 190

100 tyrimų

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cobas®

Troponin T STAT (Short Turn Around Time), CARDIAC T (skubiai atliekamas)

4. Ohmann EM, et al. Risk stratification with admission cardiac troponin T levels in acute myocardial ischemia. *N Engl J Med* 1996;335:1333-1334.
5. Christenson RH, Duh SH, Newby LK, Ohman EM, Califf RM, Granger CB, et al. Cardiac troponin T and cardiac troponin I: relative values in short-term risk stratification of patients with acute coronary syndromes. *Clin Chem* 1998;44(3):494-501.
6. Lindahl B, Diderholm E, Lagerqvist B, Venge P, Wallentin L, and the FRISC II investigators. Mechanisms behind the prognostic value of Troponin T in unstable coronary artery disease: a FRISC II substudy. *J Am Coll Card* 2001;38:979-986.
7. Aviles RJ, Askari AT, Lindahl B, Wallentin L, Jia G, Ohman EM, Mahaffey KW, Newby LK, Califf RM, Simoons ML, Topol EJ, Lauer MS. Troponin T levels in patients with acute coronary syndromes, with or without renal dysfunction. *N Engl J Med* 2002;346:2047-52.
8. Lindahl B, Venge P, Wallentin L. Troponin T Identifies Patients With Unstable Coronary Artery Disease Who Benefit From Long-Term Antithrombotic Protection. *J Am Coll Cardiol* 1997;29(1):43-48.
9. Hamm CW, Heeschen C, Goldmann B, Vahanian A, Adgey J, Miguel MC, et al. Benefit of abciximab in patients with refractory unstable angina in relation to serum troponin T levels. *N Engl J Med* 1999;340(21):1623-1629.
10. Heeschen C, Hamm CW, Goldmann BU, et al. for PRISM study investigators. Troponin concentrations for stratification of patients with acute coronary syndromes in relation to therapeutic efficacy of tirofiban. *Lancet* 1999;354:1757-1762.
11. Lindahl B, Diderholm E, Lagerqvist B, Venge P, Wallentin L. Effects on mortality of long-term treatment with I.m.w. heparin in relation to troponin T level and ECG findings - a FRISC 2 substudy. *Eur Heart J* 2000;21(Suppl.):521.
12. Newby LK, Ohman EM, Christenson RH, Moliterno DJ, Harrington RA, White HD, et al. Benefit of Glycoprotein IIb/IIIa Inhibition in Patients With Acute Coronary Syndromes and Troponin T-Positive Status: The PARAGON-B Troponin T Substudy. *Circulation* 2001;103:2891-2896.
13. Alpert JS and Thygesen K. Myocardial Infarction redefined - A consensus document of the joint European Society of Cardiology/American College of Cardiology committee for the redefinition of myocardial infarction. *JACC* 2000;36:959-969.
14. Braunwald E, Antman EM, Beasley JW, Califf RM, Cheitlin MD, Hochman JS, et al. ACC/AHA guidelines for the management of patients with unstable angina and non-ST-elevation myocardial infarction: Executive summary and recommendations. *Circulation* 2000;102:1193-1209.
15. Jaffe AS, Ravkilde J, Roberts R, Naslund U, Apple FS, Galvani M, and Katus H. It's time for a change to a troponin standard. *Circulation* 2000;102:1216-1220.
16. Apple FS, Wu AHB. Myocardial infarction redefined: Role of cardiac troponin Testing. *Clin Chem* 2001;47:377-379.
17. Setsuta K, Seino Y, Takahashi N, Ogawa T, Sasaki K, Harada A, et al. Clinical significance of elevated levels of cardiac troponin T in patients with chronic heart failure. *Am J Cardiol* 1999;84:608-611.
18. Sato Y, Kataoka K, Matsumori A, Sasayama S, Yamada T, Ito H, Takatsu Y. Measuring serum aminoterminal type III procollagen peptide, 7S domain of type IV collagen, and cardiac troponin T in patients with idiopathic dilated cardiomyopathy and secondary cardiomyopathy. *Heart* 1997;78:505-508.
19. Lauer B, Niederau C, Kühl U, Schannwell M, Pauschinger M, Strauer BE, Schultheiss HP. Cardiac troponin T in patients with clinically suspected myocarditis. *JACC* 1997;30:1354-1359.
20. Swaanenburg JCJM, Klaase JM, DeJongste MJL, Zimmerman KW, ten Duis HJ. Troponin I, troponin T, CK-MB-activity and CK-MB mass as markers for the detection of myocardial contusion in patients who experienced blunt trauma. *Clin Chim Acta* 1998;272:171-181.
21. Dierkes J, Domrose U, Westphal S, Ambrosch A, Bosselmann HP, Neumann KH, Luley C. Cardiac troponin T predicts mortality in patients with renal disease. *Circulation* 2000;102:1964-1969.
22. Giannitsis E, Müller-Bardorff M, Kurowski V, Weidtmann B, Wiegand U, Kampmann M, Katus HA. Independent prognostic value of cardiac troponin T in patients with confirmed pulmonary embolism. *Circulation* 2000;102:211-217.
23. James P, Ellis CJ, Whitlock RML, McNeil AR, Henley J, Anderson NE. Relation between troponin T concentration and mortality in patients presenting with an acute stroke: observational study. *BMJ* 320:1502-1504.
24. ver Elst KM, Spapen HD, Nguyen DN, Garbar C, Huyghens LP, Gorus FK. Cardiac troponin I and T are biological markers of left ventricular dysfunction in septic shock. *Clin Chem* 2000;46:650-657.
25. Jenkins DP, Pugsley WB, Alkhulaifi AM, Kemp M, Hooper J, Yellon DM. Ischaemic preconditioning in patients undergoing coronary artery bypass surgery. *Heart* 1997;77:314-318.
26. Lopez-Jimenez F, Goldman L, Sacks DB, Thomas EJ, Johnson PA, Cook EF, Lee TH. Prognostic value of cardiac troponin T after non-cardiac surgery: 6 month follow up data. *JACC* 1997;29:1241-1245.
27. Johansen O, Brekke M, Stromme JH, Valen V, Seljeflot I, Skjaeggstad O, Arensen H. Myocardial damage during percutaneous transluminal coronary angioplasty as evidenced by troponin T measurements. *Eur Heart J* 1998;19:112-117.
28. Herman EH, Lipshultz SE, Rifai N, Zhang J, Papoian T, Yu ZX, et al. Use of cardiac troponin T levels as an indicator of doxorubicin-induced cardiotoxicity. *Cancer Res* 1998;58:195-197.
29. Ooi DS, Zimmerman D, Graham J, Wells GA. Cardiac troponin T predicts long-term outcomes in hemodialysis patients. *Clin Chem* 2001;47:412-417.
30. Apple FS, Murakami MM, Pearce LA, Herzog CA. Predictive value of cardiac troponin I and T for subsequent death in end-stage renal disease. *Circulation* 2002;106:2941-2945.

Troponinas T STAT

Troponin T STAT (Short Turn Around Time), CARDIAC T (skubiai atliekamas)

31. deFilippi C, Wasserman S, Rosanio S, Tiblier W, Sperger H, Tocchi M, Christenson R, Uretsky B, Smiley M, Gold J, Muniz H, Badalamenti J, Herzog C, Henrich W. Cardiac Troponin T and C-reactive protein for predicting prognosis, coronary atherosclerosis, and cardiomyopathy in patients undergoing long term hemodialysis. JAMA 2003;290:353-359.
32. Iliou MC, Fumeron C, Benoit MO, Tuppin P, Calonge VM, Moatti N, Buisson C, Jacquot C. Prognostic value of cardiac Markers in ESRD: Chronic Hemodialysis and New Cardiac markers Evaluation (CHANCE) study. Am J Kidney Dis 2003;42:513-523.
33. Choy JB, Armstrong PW, Ulan RA, Campbell PM, Gourishankar S, Prosser CI, Tymchak W. Do cardiac troponins provide prognostic insight in hemodialysis patients? Can J Cardiol 2003;19:907-911.
34. Scott B, Deman A, Peeters P, Van den Branden C, Stolear JC, Van Camp G, Verbeelen D. Cardiac Troponin T and malondialdehyde modified plasma lipids in haemodialysis patients. Nephrol Dial Transplant 2003;18:737-742.
35. Davis GK, Labugger R, Van Eyk JE, Apple FS. Cardiac troponin T is not detected in western blots of diseased renal tissue. Clin Chem 2001;47:782-783.
36. Ricchiuti V, Voss EM, Ney A, Odland M, Anderson PAW, Apple FS. Cardiac Troponin T isoforms expressed in renal diseased skeletal muscle will not cause false positive results by the second generation cTnT assay by Boehringer Mannheim; Clin Chem 1998;44(9):1919-1924.
37. Wu AHB, Feng Y-J, Moore R, et al. Characterisation of cardiac troponin subunit release into serum after acute myocardial infarction and comparison of assays for troponin T and I. Clin Chem 1998;44:1198-1208.
38. Hallermayer K, Klenner D and Vogel R. Use of recombinant human cardiac troponin T for standardization of third generation troponin T methods. Scand J Clin Invest 1999;59(Suppl 230):128-131.
39. Müller-Bardorff M, et al. Improved troponin T ELISA specific for cardiac troponin T isoform: assay development and analytical validation. Clin Chem 1997;43(3):458-466.
40. Kampmann M, Rauscher T, Müller-Bardorff M, Münch M, Klein G, Katus H. Clinical Evaluation of Troponin T and CK-MB Mass on the Elecsys 2010 analyzer. Poster presented at Medlab 97, 12th IFCC European Congress of Clinical Chemistry, August 17-22, 1997, Basel, Switzerland.
41. Kobayashii S, Tanaka M, Tamura N, Hashimoto H, Hirose S-I. Serum cardiac troponin T in polymyositis/dermatomyositis. Lancet 1992;340(8821):726.
42. Mair J, Dienstl F, Puschendorf B. Cardiac Troponin T in the Diagnosis of Myocardial Injury. Critical Reviews in Clin Lab Sci 1992;29(1):31-57.
43. Carrier M, Solymoss BC, Raymond C, Leclerc Y, Pelletier LC. Cardiac Troponin T and Creatine Kinase MB Isoenzyme as Biochemical Markers of Ischemia after Heart Preservation and Transplantation. J Heart Lung Transplant 1994;13(4):696-700.
44. Löfberg M, Tähtelä R, Härkönen M, Somer H. Myosin Heavy-Chain Fragments and Cardiac Troponins in the Serum in Rhabdomyolysis. Diagnostic Specificity of New Biochemical Markers. Arch Neurol 1995;52:1210-1214.
45. Anderson JR, Hossein-Nia M, Brown P, Holt DW, Murday A. Donor Cardiac Troponin T Predicts Subsequent Inotropic Requirements Following Cardiac Transplantation. Brief Communications in Transplant 1994;58(9):1056-1057.
46. Franz WM, Remppis A, Kandolf R, Kübler W, Katus HA. Serum Troponin T: Diagnostic Marker for Acute Myocarditis. Letter to the Editor in Clin Chem 1996;42(2):340-341.
47. Mair P, Mair J, Koller J, Wieser C, Artner-Dworzak E, Puschendorf B. Cardiac troponin T in the diagnosis of heart contusion. Lancet 1991;338:693.
48. Bablok W, et al. A General Regression Procedure for Method Transformation. J Clin Chem Clin Biochem 1988;26:783-790.

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©2007 Roche Diagnostics
Roche Diagnostics GmbH, D-68298 Mannheim