

Order information

REF	CONTENT	Analyzer(s) on which cobas c pack(s) can be used
0333795 190	Tina-quant α 1-Acid Glycoprotein Gen.2 (100 tests)	System-ID 07 6758 1 COBAS INTEGRA 400 plus COBAS INTEGRA 800
11355279 216	Calibrator f.a.s. Proteins (5 x 1 mL)	System-ID 07 6557 0
11355279 160	Calibrator f.a.s. Proteins (5 x 1 mL, for USA)	System-ID 07 6557 0
10557897 122	Precinorm Protein (3 x 1 mL)	System-ID 07 9105 9
10557897 160	Precinorm Protein (3 x 1 mL, for USA)	System-ID 07 9105 9
11333127 122	Precipath Protein (3 x 1 mL)	System-ID 07 9106 7
11333127 160	Precipath Protein (3 x 1 mL, for USA)	System-ID 07 9106 7
05117003 190	PreciControl ClinChem Multi 1 (20 x 5 mL)	System-ID 07 7469 3
05947626 190	PreciControl ClinChem Multi 1 (4 x 5 mL)	System-ID 07 7469 3
05947626 160	PreciControl ClinChem Multi 1 (4 x 5 mL, for USA)	System-ID 07 7469 3
05117216 190	PreciControl ClinChem Multi 2 (20 x 5 mL)	System-ID 07 7470 7
05947774 190	PreciControl ClinChem Multi 2 (4 x 5 mL)	System-ID 07 7470 7
05947774 160	PreciControl ClinChem Multi 2 (4 x 5 mL, for USA)	System-ID 07 7470 7
20756350 322	NaCl Diluent 9 % (6 x 22 mL)	System-ID 07 5635 0

English

System information

Test AAGP2, test ID 0-258

Intended use

In vitro test for the quantitative immunological determination of human α 1-acid glycoprotein in human serum and plasma on COBAS INTEGRA systems.

Summary^{1,2,3,4,5}

α 1-Acid glycoprotein is synthesized in hepatocytes and consists of a polypeptide chain having 5 carbohydrate chains N-glycosidically bonded to it (molar mass 41000 daltons). Structurally, it belongs to the lipocalin superfamily of secretory proteins (such as α 1-microglobulin and retinol-binding protein). α 1-Acid glycoprotein promotes fibroblast growth and interacts with collagen.

It is a sensitive acute phase reactant whose concentration can increase by a factor of 3 within 24-48 hours when inflammation occurs. α 1-Acid glycoprotein can also be used to differentiate between acute phase reactions (elevated serum level) and estrogen effects (normal or decreased serum level) whereas the serum level of other positive reactants such as ceruloplasmin and haptoglobin increases during such reactions. Along with haptoglobin it is perhaps the best protein for identifying slight in vivo hemolysis. An increased α 1-acid glycoprotein level and normal haptoglobin values indicate an acute phase reaction with concomitant slight in vivo hemolysis. Moderate and isolated increases occur when glomerular filtration is inhibited in the early stages of uremia. The determination is used in the assessment of the activity of acute and recurring inflammations as well as of tumors with cell necrosis.

Various assay methods for α 1-acid glycoprotein determination are available such as kinetic nephelometry, radial immunodiffusion (RID) and turbidimetry. The Roche α 1-acid glycoprotein assay is based on the principle of immunological agglutination.

Test principle⁵

Immunoturbidimetric assay

- Sample and addition of R1
- Addition of SR and start of reaction

Anti- α 1-acid glycoprotein antibodies react with antigen in the sample to form an antigen/antibody complex. Following agglutination, this is measured turbidimetrically.

Reagents - working solutions

R1 TRIS buffer: 50 mmol/L, pH 8.0; PEG: 7 %; NaCl: 300 mmol/L; preservative; stabilizer

SR Polyclonal anti-human α 1-acid glycoprotein antibody (goat): dependent on titer; TRIS buffer: 13 mmol/L, pH 7.5; NaCl: 198 mmol/L; preservative

R1 is in position B and SR is in position C.

Precautions and warnings

Pay attention to all precautions and warnings listed in Section 1 / Introduction of this Method Manual.

For USA: For prescription use only.

Reagent handling

Ready for use

Storage and stability

Shelf life at 2-8 °C See expiration date on **cobas c** pack label

COBAS INTEGRA 400 plus system

On-board in use at 10-15 °C 12 weeks

COBAS INTEGRA 800 system

On-board in use at 8 °C 12 weeks

Specimen collection and preparation

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable.

Serum

Plasma: Li-, Na-, NH₄⁺-heparin; K₂-, K₃-EDTA

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

Stability:⁶ < 72 hours at 4 °C
6 months at -20 °C

Materials provided

See "Reagents – working solutions" section for reagents.

Materials required (but not provided)

NaCl Diluent 9 %, Cat. No. 20756350 322, system-ID 07 5635 0 for automatic sample dilution and standard serial dilutions. NaCl Diluent 9 % is

AAGP2

Tina-quant α 1-Acid Glycoprotein Gen.2

placed in its predefined rack position and is stable for 4 weeks on-board COBAS INTEGRA 400 plus/800 analyzers.

Assay

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

Application for serum and plasma

COBAS INTEGRA 400 plus test definition

Measuring mode	Absorbance
Abs. calculation mode	Endpoint
Reaction mode	D-R1-S-SR
Reaction direction	Increase
Wavelength A/B	340/659 nm
Calc. first/last	33/69
Typical prozone effect	> 11.0 g/L (> 1100 mg/dL or > 275 μ mol/L)
Antigen excess check	No
Predilution factor	21
Unit	g/L

Pipetting parameters

		Diluent (H ₂ O)
R1	120 μ L	
Sample	12 μ L	8 μ L
SR	40 μ L	8 μ L
Total volume	188 μ L	

COBAS INTEGRA 800 test definition

Measuring mode	Absorbance
Abs. calculation mode	Endpoint
Reaction mode	D-R1-S-SR
Reaction direction	Increase
Wavelength A/B	340/659 nm
Calc. first/last	44/98
Typical prozone effect	> 11.0 g/L (> 1100 mg/dL or > 275 μ mol/L)
Antigen excess check	No
Predilution factor	21
Unit	g/L

Pipetting parameters

		Diluent (H ₂ O)
R1	120 μ L	
Sample	12 μ L	8 μ L
SR	40 μ L	8 μ L
Total volume	188 μ L	

Calibration

Calibrator	C.f.a.s. Proteins
Calibration dilution ratio	COBAS INTEGRA 400 plus analyzers: 1:7, 1:15, 1:30, 1:75, 1:150, and 0 g/L performed automatically by the instrument.

COBAS INTEGRA 800 analyzers:

1:7.5, 1:15, 1:30, 1:75, 1:150, and 0 g/L performed automatically by the instrument.

Calibration mode	Logit/log 4
Calibration replicate	Duplicate recommended
Calibration interval	Each lot and as required following quality control procedures

Enter the assigned lot-specific α 1-acid glycoprotein value of the undiluted calibrator indicated in the package insert for C.f.a.s. Proteins.

Traceability: This method is standardized against an internal method traceable to CRM 470.

Quality control

Reference range	Precinorm Protein or PreciControl ClinChem Multi 1
Pathological range	Precipath Protein or PreciControl ClinChem Multi 2
Control interval	24 hours recommended
Control sequence	User defined
Control after calibration	Recommended

For quality control, use control materials as listed in the "Order information" section. In addition, other suitable control material can be used.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

Calculation

COBAS INTEGRA analyzers automatically calculate the analyte concentration of each sample. For more details, please refer to Data Analysis in the Online Help (COBAS INTEGRA 400 plus/800 analyzers).

Conversion factors:	g/L \times 25 = μ mol/L
	g/L \times 100 = mg/dL
	mg/dL \times 0.25 = μ mol/L

Limitations - interference

Criterion: Recovery within \pm 10 % of initial value.

Icterus:⁷ No significant interference up to an I index of 60 for conjugated and unconjugated bilirubin (approximate conjugated and unconjugated bilirubin concentration: 60 mg/dL or 1026 μ mol/L).

Hemolysis:⁷ No significant interference up to an H index of 1000 (approximate hemoglobin concentration: 1000 mg/dL or 621 μ mol/L).

Lipemia:⁷ No significant interference up to an L index of 700. There is poor correlation between the L index (corresponds to turbidity) and triglycerides concentration.

Rheumatoid factors: Rheumatoid factors < 1200 IU/mL do not interfere.

Drugs: No interference was found at therapeutic concentrations using common drug panels.^{8,9}

In very rare cases, gammopathy, in particular type IgM (Waldenström's macroglobulinemia), may cause unreliable results.¹⁰

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

ACTION REQUIRED

Special Wash Programming: The use of special wash steps is mandatory when certain test combinations are run together on COBAS INTEGRA analyzers. Refer to the CLEAN Method Sheet for further instructions and for the latest version of the Extra wash cycle list.

Where required, special wash/carry-over evasion programming must be implemented prior to reporting results with this test.

Limits and ranges**Measuring range**0.25-4.0 g/L (6.25-100 μ mol/L or 25-400 mg/dL) (typical measuring range)

The upper limit of the measuring range depends on the actual calibrator value.

Determine samples having higher concentrations via the rerun function. Dilution of samples via the rerun function is a 1:2 dilution. Results from samples diluted using the rerun function are automatically multiplied by a factor of 2.

Lower limits of measurement

Lower detection limit of the test:

0.10 g/L (2.5 μ mol/L or 10 mg/dL)

The lower detection limit represents the lowest measurable analyte level that can be distinguished from zero. It is calculated as the value lying 3 standard deviations above that of a zero sample (zero sample + 3 SD, repeatability, n = 30).

Expected values*0.5-1.2 g/L (50-120 mg/dL or 12.5-30 μ mol/L)* Reference range according to CRM 470 protein standardization¹¹

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

Specific performance data

Representative performance data on the COBAS INTEGRA analyzers are given below. Results obtained in individual laboratories may differ.

Precision

Precision was determined using human samples and controls in an internal protocol with repeatability (n = 21) and intermediate precision (1 aliquot per run, 1 run per day, 21 days). The following results were obtained:

Repeatability	Mean g/L (μ mol/L)	SD g/L (μ mol/L)	CV %
Precinorm Protein	0.83 (20.8)	0.01 (0.3)	1.8
Precipath Protein	1.39 (34.8)	0.02 (0.5)	1.5
Human serum	0.87 (21.8)	0.02 (0.5)	2.0

Intermediate precision	Mean g/L (μ mol/L)	SD g/L (μ mol/L)	CV %
Precinorm Protein	0.80 (20.0)	0.02 (0.5)	2.3
Precipath Protein	1.34 (33.5)	0.02 (0.5)	1.7
Human serum	0.85 (21.3)	0.02 (0.5)	2.9

Method comparison α 1-Acid glycoprotein values for human serum and plasma samples obtained on a COBAS INTEGRA 700 analyzer with the COBAS INTEGRA AAGP Gen.2 reagent (y) were compared with those determined using the same reagent on a Roche/Hitachi 917 analyzer (x) and with the previous reagent (AAGP) on a COBAS INTEGRA 700 analyzer (x).**Roche/Hitachi 917 analyzer**

Sample size (n) = 55

Passing/Bablok¹²

Linear regression

y = 1.043x - 0.043 g/L

y = 1.027x - 0.025 g/L

r = 0.945

r = 0.996

SD (md 95) = 0.074

Sy.x = 0.041

The sample concentrations were between 0.4 and 2.85 g/L (40 and 285 mg/dL or 10 and 71.3 μ mol/L).**COBAS INTEGRA 700 analyzer**

Sample size (n) = 55

Passing/Bablok¹²

Linear regression

y = 1.005x - 0.005 g/L

y = 1.001x + 0.001 g/L

r = 0.964

r = 0.998

SD (md 95) = 0.056

Sy.x = 0.027

The sample concentrations were between 0.4 and 2.75 g/L (40 and 275 mg/dL or 10 and 68.8 μ mol/L).**References**

- Schmid K. α 1-Acid glycoprotein. In: The Plasma Proteins, 2nd ed. Putnam FW, ed. New York: Academic Press 1975;183-228.
- Greiling H, Gressner AM, eds. Lehrbuch der Klinischen Chemie und Pathobiochemie, 3rd ed. Stuttgart/New York: Schattauer Verlag 1995;234-236.
- Tietz NW, ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia, PA: WB Saunders Company 1995;66-67.
- Ganrot K. Plasma protein pattern in acute infectious disease. Scand J Clin Lab Invest 1974;34:75-81.
- Lievens M, Bienvenu J, Buitrago JMG, et al. Evaluation of four new Tina-quant assays for determination of α 1-acid glycoprotein, α 1-antitrypsin, haptoglobin and prealbumin. Clin Lab 1996;42:515-520.
- Wu AHB, ed. Tietz Clinical Guide to Laboratory Tests, 4th ed. St. Louis (MO): Saunders Elsevier 2006;42.
- Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-475.
- Breuer J. Report on the Symposium "Drug effects in Clinical Chemistry Methods". Eur J Clin Chem Clin Biochem 1996;34:385-386.
- Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385.
- Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med 2007;45(9):1240-1243.
- Consensus values of the Deutsche Gesellschaft für Laboratoriumsmedizin, the Deutsche Gesellschaft für Klinische Chemie und der Verband der Diagnostica-Industrie e.V. (VDGH). DG Klinische Chemie Mitteilungen 1995;26:119-122.
- Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard.

CONTENT

Contents of kit



Volume after reconstitution or mixing

GTIN

Global Trade Item Number

FOR US CUSTOMERS ONLY: LIMITED WARRANTY

Roche Diagnostics warrants that this product will meet the specifications stated in the labeling when used in accordance with such labeling and will be free from defects in material and workmanship until the expiration date printed on the label. THIS LIMITED WARRANTY IS IN LIEU OF ANY OTHER WARRANTY, EXPRESS OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. IN NO EVENT SHALL ROCHE DIAGNOSTICS BE LIABLE FOR INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES.

COBAS, COBAS INTEGRA, COBAS C, TINA-QUANT, PRECINORM, PRECIPATH and PRECICONTROL are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2015, Roche Diagnostics

0003333795190COINV7.0

AAGP2

Tina-quant α 1-Acid Glycoprotein Gen.2



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim
www.roche.com

Distribution in USA by:
Roche Diagnostics, Indianapolis, IN
US Customer Technical Support 1-800-428-2336



cobas[®]
Specific proteins