

REF	CONTENT	System-ID	Analyzers on which cobas c pack can be used
04512936 190	ONLINE DAT Amphetamines II (200 tests)	07 6835 9	COBAS INTEGRA 400 plus COBAS INTEGRA 800
03304671 190	Preciset DAT Plus I CAL 1-6 (6 × 5 mL)		
03304680 190	Preciset DAT Plus II CAL 1-6 (6 × 5 mL)		
03304698 190	C.f.a.s. DAT Qualitative Plus (6 × 5 mL)		
04590856 190	C.f.a.s. DAT Qualitative Plus Clinical (3 × 5 mL)		
03312968 190	Control Set DAT II (for 300 ng/mL assay) PreciPos DAT Set II (2 × 10 mL) PreciNeg DAT Set II (2 × 10 mL)		
03312950 190	Control Set DAT I (for 500 ng/mL assay) PreciPos DAT Set I (2 × 10 mL) PreciNeg DAT Set I (2 × 10 mL)		
04500873 190	Control Set DAT Clinical (for 500 ng/mL assay) PreciPos DAT Clinical (2 × 10 mL) PreciNeg DAT Clinical (2 × 10 mL)		
03312976 190	Control Set DAT III (for 1000 ng/mL assay) PreciPos DAT Set III (2 × 10 mL) PreciNeg DAT Set III (2 × 10 mL)		

English

System information

Test AM3S2, test-ID 0-358 for semiquantitative assay, 300 ng/mL

Test AM5S2, test-ID 0-359 for semiquantitative assay, 500 ng/mL

Test AM1S2, test-ID 0-360 for semiquantitative assay, 1000 ng/mL

Test AM3Q2, test-ID 0-330 for qualitative assay, 300 ng/mL

Test AM5Q2, test-ID 0-340 for qualitative assay, 500 ng/mL

Test AM1Q2, test-ID 0-350 for qualitative assay, 1000 ng/mL

Test AM5QC, test-ID 0-341 for qualitative assay, 500 ng/mL using

C.f.a.s. DAT Qualitative Plus Clinical

Intended use

Amphetamines II (AMPII) is an in vitro diagnostic test for the qualitative and semiquantitative detection of amphetamines and methamphetamines in human urine on COBAS INTEGRA systems at cutoff concentrations of 300 ng/mL, 500 ng/mL and 1000 ng/mL when calibrated with *d*-methamphetamine.

Semiquantitative test results may be obtained that permit laboratories to assess assay performance as part of a quality control program.

Semiquantitative assays are intended to determine an appropriate dilution of the specimen for confirmation by a confirmatory method such as gas chromatography/mass spectrometry (GC/MS).

Amphetamines II provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. GC/MS is the preferred confirmatory method.¹ Clinical consideration and professional judgement should be applied to any drug of abuse test result, particularly when preliminary positive results are used.

Summary

The amphetamines are known as the sympathomimetic amines as they mimic the effects of stimulation of the sympathetic nervous system. These small molecules, based on β-phenylethylamine, structurally resemble the bodies own catecholamines. A wide variety have been created via substitutions anywhere on the structure. The amphetamines are potent central nervous stimulants. As such they can increase wakefulness, physical activity, and decrease appetite. The amphetamines have some limited indications and approval for use in ADHD, narcolepsy, and obesity. However, because these CNS stimulants convey a sense of self-confidence, well being, and euphoria, they are highly addictive, widely abused, and consequently controlled substances.² Abuse can lead to medical, psychological, and social consequences. Adverse health effects include memory loss, aggression, psychotic behavior, heart damage, malnutrition, and severe dental problems.³ Amphetamine may be self-administered either orally or by intravenous injection in amounts of up to 2000 mg daily by tolerant addicts. It is a metabolite of a number of other

drugs including methamphetamine. Normally about 30 % is excreted unchanged in the 24 hour urine, but this may change to as much as 74 % in acid urine and may decrease to 1 % in alkaline urine.⁴

Amphetamines II is calibrated with *d*-methamphetamine and therefore the sensitivity towards amphetamines is different than *d*-methamphetamine, as indicated in the "Analytical specificity" section.

Test principle

The assay is based on the kinetic interaction of microparticles in a solution (KIMS)^{5,6} as measured by changes in light transmission. In the absence of amphetamines and/or methamphetamines, soluble drug conjugates bind to antibody-bound microparticles, causing the formation of particle aggregates. As the aggregation reaction proceeds in the absence of amphetamines and/or methamphetamines, the absorbance increases.

When a urine sample contains amphetamines and/or methamphetamines, this drug competes with the drug derivative conjugate for microparticle-bound antibody. Antibody bound to amphetamine and/or methamphetamine is no longer available to promote particle aggregation, and subsequent particle lattice formation is inhibited. The presence of amphetamines and/or methamphetamines diminishes the increasing absorbance in proportion to the concentration of drug in the sample. Amphetamines and/or methamphetamines content is determined relative to the value obtained for a known cutoff concentration of *d*-methamphetamine.

Reagents - working solutions

- R1** Conjugated amphetamine and methamphetamine derivatives; buffer; bovine serum albumin; 0.09 % sodium azide
- SR** Microparticles attached to amphetamine and methamphetamine antibodies (mouse monoclonal); buffer; bovine serum albumin; 0.09 % sodium azide

R1 is in position A and SR is in position C.

Precautions and warnings

Pay attention to all precautions and warnings listed in Section 1 / Introduction of this Method Manual.

For USA: For prescription use only.

Reagent handling

COBAS INTEGRA 400 plus analyzer

Mix all new (non-punctured) **cobas c** packs for 1 minute on a cassette mixer before loading on the analyzer. All in-use **cobas c** packs must also be mixed in the same manner at the beginning of each week (once a week).

COBAS INTEGRA 800 analyzer

Ready for use. After **cobas c** pack puncture, the analyzer automatically mixes the reagent for 1 minute and for half a minute during Begin of Day.

Storage and stability

Shelf life at 2-8 °C: See expiration date on
cobas c pack label

| **COBAS INTEGRA 400 plus analyzer**

On-board in use at 10-15 °C 12 weeks

COBAS INTEGRA 800 analyzer

On-board in use at 8 °C 12 weeks

- | Do not freeze reagents. Reagents that have been frozen should be discarded.

Specimen collection and preparation

Only the specimens listed below were tested and found acceptable.

Urine: Collect urine samples in clean glass or plastic containers. Fresh urine specimens do not require any special handling or pretreatment, but an effort should be made to keep pipetted samples free of gross debris. Samples should be within the normal physiological pH range of 5-8. No additives or preservatives are required. It is recommended that urine specimens be stored at 2-8 °C and tested within 5 days of collection.⁷

Centrifuge highly turbid specimens before testing.

Adulteration or dilution of the sample can cause erroneous results. If adulteration is suspected, another sample should be collected. Specimen validity testing is required for specimens collected under the *Mandatory Guidelines for Federal Workplace Drug Testing Programs*.⁸

Caution: Specimen dilutions should only be used as an estimation for GC/MS and are not intended for patient values. Dilution procedures, when used, should be validated.

Materials provided

See "Reagents – working solutions" section for reagents.

Assay

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

Application for urine| **COBAS INTEGRA 400 plus test definition**

	<i>Semiquantitative</i>	<i>Qualitative</i>
Measuring mode	Absorbance	Absorbance
Abs. calculation mode	Endpoint	Endpoint
Reaction mode	R1-S-SR	R1-S-SR
Reaction direction	Increase	Increase
Reaction start	SR	SR
Wavelength A	659 nm	659 nm
Test range		
	AM3S2 0-2000 ng/mL	0-2000 AM3Q2
	AM5S2, AM1S2 0-5000 ng/mL	0-5000 AM5Q2, AM5QC, AM1Q2
with postdilution		
	AM3S2 0-20000 ng/mL	
	AM5S2, AM1S2 0-50000 ng/mL	
Postdilution factor	10 recommended ^{a)}	No
Calc. first/last	51/70	51/70
Unit	ng/mL	

a) For use when estimating concentration in preparation for GC/MS analysis.

Pipetting parameters

AM3S2, AM3Q2		Diluent (H ₂ O)
R1	60 µL	20 µL
Sample	10 µL	5 µL

SR	47 µL	5 µL
Total volume	147 µL	
AM5S2, AM5Q2, AM5QC		Diluent (H ₂ O)
R1	60 µL	20 µL
Sample	7.5 µL	5 µL
SR	47 µL	5 µL
Total volume	144.5 µL	
AM1S2, AM1Q2		Diluent (H ₂ O)
R1	60 µL	20 µL
Sample	5 µL	5 µL
SR	47 µL	5 µL
Total volume	142 µL	

COBAS INTEGRA 800 test definition

	<i>Semiquantitative</i>	<i>Qualitative</i>
Measuring mode	Absorbance	Absorbance
Abs. calculation mode	Endpoint	Endpoint
Reaction mode	R1-S-SR	R1-S-SR
Reaction direction	Increase	Increase
Reaction start	SR	SR
Wavelength A	659 nm	659 nm
Test range		
	AM3S2 0-2000 ng/mL	0-2000
	AM5S2, AM1S2 0-5000 ng/mL	
with postdilution		
	AM3S2 0-20000 ng/mL	
	AM5S2, AM1S2 0-50000 ng/mL	
Postdilution factor	10 recommended ^{b)}	No
Calc. first/last	46/86	46/86
Unit	ng/mL	

b) For use when estimating concentration in preparation for GC/MS analysis.

Pipetting parameters

AM3S2, AM3Q2		Diluent (H ₂ O)
R1	60 µL	20 µL
Sample	10 µL	5 µL
SR	47 µL	3 µL
Total volume	145 µL	
AM5S2, AM5Q2, AM5QC		Diluent (H ₂ O)
R1	60 µL	20 µL
Sample	7.5 µL	5 µL
SR	47 µL	3 µL
Total volume	142.5 µL	
AM1S2, AM1Q2		Diluent (H ₂ O)
R1	60 µL	20 µL
Sample	5 µL	5 µL
SR	47 µL	3 µL
Total volume	140 µL	

Calibration

Calibrators

Semiquantitative applications

AM3S2, 0-358 Preciset DAT Plus II calibrators, CAL 1-6
0, 150, 300, 600, 1000, 2000 ng/mL
d-methamphetamine
(300 cutoff, DATS8, system-ID 07 6796 4)

AM5S2, 0-359;
AM1S2, 0-360 Preciset DAT Plus I calibrators, CAL 1-6
0, 250, 500, 1000, 3000, 5000 ng/mL
d-methamphetamine
(500 and 1000 cutoffs, DATS2,
system-ID 07 6764 6)

Qualitative applications

AM3Q2, 0-330 Preciset DAT Plus II calibrators, CAL 1
0 ng/mL or deionized water
and
Preciset DAT Plus II calibrators, CAL 3
300 ng/mL
(300 cutoff, DATQ3, system-ID 07 6770 0)
For qualitative applications, the cutoff of 300 ng/mL
is assigned a value of 1000.

AM5Q2, 0-340 (Preciset DAT Plus I calibrators, CAL 1)
0 ng/mL or deionized water
and
C.f.a.s. DAT Qualitative Plus
500 ng/mL
(500 cutoff, DATQ1, system-ID 07 6744 1)
For qualitative applications, the cutoff of 500 ng/mL
is assigned a value of 1000.

AM5QC, 0-341 Preciset DAT Plus I or II calibrators, CAL 1
0 ng/mL or deionized water
and
C.f.a.s. DAT Qualitative Plus Clinical
500 ng/mL
(500 cutoff, DATQ5, system-ID 07 6880 4)
For qualitative applications, the cutoff of 500 ng/mL
is assigned a value of 1000.

AM1Q2, 0-350 Preciset DAT Plus I calibrators, CAL 1
0 ng/mL or deionized water
and
Preciset DAT Plus I calibrators, CAL 4
1000 ng/mL
(1000 cutoff, DATQ2, system-ID 07 6768 9)
For qualitative applications, the cutoff of
1000 mg/dL is assigned a value of 1000.

Calibration mode

Semiquantitative applications

COBAS INTEGRA 400 plus analyzer:
Logit/log 5 (AM1S2, 0-360)
Logit/log 4 (AM3S2, 0-358; AM5S2, 0-359)
COBAS INTEGRA 800 analyzer:
Logit/log 4

Qualitative applications

Calibration replicate

Linear regression
Duplicate recommended

Calibration interval

COBAS INTEGRA 400 plus analyzer:
Each lot, every 4 weeks, and as required following
quality control procedures
COBAS INTEGRA 800 analyzer:
Each lot, every 4 weeks, and as required following
quality control procedures

A calibration curve is generated using the calibrators. Calibrators must be placed from the highest concentration first to the lowest last on the CAL/QC rack. This curve is retained in memory by the COBAS INTEGRA systems and recalled for later use.

Traceability: This method has been standardized against a primary reference method GC/MS.

Quality control

Quality control

300 ng/mL cutoff

Control Set DAT II

PreciPos DAT Set II
(DAT2P, system-ID 07 6771 9)

PreciNeg DAT Set II
(DAT2N, system-ID 07 6772 7)

500 ng/mL cutoff

Control Set DAT I

PreciPos DAT Set I
(DAT1P, system-ID 07 6753 0)

PreciNeg DAT Set I
(DAT1N, system-ID 07 6754 9)

or

Control Set DAT Clinical

PreciPos DAT Clinical
(DATCP, system-ID 07 6879 0)

PreciNeg DAT Clinical
(DATCN, system-ID 07 6878 2)

1000 ng/mL cutoff

Control Set DAT III

PreciPos DAT Set III
(DAT3P, system-ID 07 6773 5)

PreciNeg DAT Set III
(DAT3N, system-ID 07 6774 3)

Control sequence

User defined

Control after calibration

Recommended

For quality control, use control materials as listed in the "Order information" section. In addition, other suitable control material can be used.

Drug concentrations of Control Set DAT I, II, III, and Clinical have been verified by GC/MS.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

Results

COBAS INTEGRA systems report results with the following test flags:

*Semiquantitative result reporting***AM3S2 (300 ng/mL cutoff)**

Flag	COBAS INTEGRA	Value range
No flag	Negative	< 300 ng/mL

Flag	COBAS INTEGRA	Value range
<TEST RNG	Negative	< 0 ng/mL
>TEST RNG	Positive	> 2000 ng/mL
POS 300	Positive	≥ 300 ng/mL

Value ranges listed above are based on a cutoff value of 300 ng/mL.

Semiquantitative result reporting
AM5S2 (500 ng/mL cutoff)

Flag	COBAS INTEGRA	Value range
No flag	Negative	< 500 ng/mL
<TEST RNG	Negative	< 0 ng/mL
>TEST RNG	Positive	> 5000 ng/mL
POS 500	Positive	≥ 500 ng/mL

Value ranges listed above are based on a cutoff value of 500 ng/mL.

Semiquantitative result reporting
AM1S2 (1000 ng/mL cutoff)

Flag	COBAS INTEGRA	Value range
No flag	Negative	< 1000 ng/mL
<TEST RNG	Negative	< 0 ng/mL
>TEST RNG	Positive	> 5000 ng/mL
POS 1000	Positive	≥ 1000 ng/mL

Value ranges listed above are based on a cutoff value of 1000 ng/mL.

Qualitative result reporting
AM3Q2 (300 ng/mL cutoff)

Flag	COBAS INTEGRA	Value range
No flag	Negative	< 1000
<TEST RNG	Negative	< 0
>TEST RNG	Positive	> 2000
POS 1000	Positive	≥ 1000

Value ranges above are based on assigning the cutoff of 300 ng/mL a value of 1000.

Qualitative result reporting
AM5Q2, AM5QC (500 ng/mL cutoff)

Flag	COBAS INTEGRA	Value range
No flag	Negative	< 1000
<TEST RNG	Negative	< 0
>TEST RNG	Positive	> 2000
POS 1000	Positive	≥ 1000

Value ranges above are based on assigning the cutoff of 500 ng/mL a value of 1000.

Qualitative result reporting
AM1Q2 (1000 ng/mL cutoff)

Flag	COBAS INTEGRA	Value range
No flag	Negative	< 1000
<TEST RNG	Negative	< 0
>TEST RNG	Positive	> 2000
POS 1000	Positive	≥ 1000

Value ranges above are based on assigning the cutoff of 1000 ng/mL a value of 1000.

Semiquantitative result reporting

The semiquantitation of preliminary positive results should only be used by laboratories to determine an appropriate dilution of the specimen for confirmation by a confirmatory method such as GC/MS. It also permits the laboratory to establish quality control procedures and assess control performance.

Note: When using the post-dilution function (1:10 dilution), to ensure the sample was not over-diluted, the diluted result must be at least half the analyte cutoff value times 10. If the diluted result falls below half the analyte cutoff value times 10, repeat the sample with a smaller dilution. A dilution that produces a result closest to the analyte cutoff is the most accurate estimation.

To estimate the preliminary positive sample's concentration, multiply the result by the appropriate dilution factor. Dilutions should only be used as an estimation for GC/MS.

Limitations - interference

See the "Specific performance data" section of this document for information on substances tested with this assay. There is the possibility that other substances and/or factors may interfere with the test and cause erroneous results (e.g., technical or procedural errors).

A preliminary positive result with this assay indicates the presence of amphetamines and/or methamphetamines in urine. It does not measure the level of intoxication.

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

ACTION REQUIRED

Special Wash Programming: The use of special wash steps is mandatory when certain test combinations are run together on COBAS INTEGRA analyzers. Refer to the CLEAN Method Sheet for further instructions and for the latest version of the Extra wash cycle list.

Where required, special wash/carry-over evasion programming must be implemented prior to reporting results with this test.

Expected values

No drug should be present in individuals that have not ingested amphetamines and/or methamphetamines.

Specific performance data

Representative performance data on the COBAS INTEGRA analyzers are given below. Results obtained in individual laboratories may differ.

Precision

A *d*-methamphetamine (MAMP) solution (1 mg/mL) was added to 9 samples obtained from a human urine sample pool to achieve concentrations at approximately -100 %, -75 %, -50 %, -25 %, ±0 %, +25 %, +50 %, +75 %, and +100 % of the cutoff value. These samples were tested for precision in qualitative and semiquantitative modes. Following a CLSI (EP5-A2) precision protocol, samples were tested in 2 replicates per run, 2 runs per day for 21 days, total n = 84.

Qualitative - 300 ng/mL Cutoff

Drug	Concentration of Sample	Number of Determinations	Results # Neg / # Pos
MAMP	zero drug	84	84 Neg / 0 Pos
MAMP	-75 %	84	84 Neg / 0 Pos
MAMP	-50 %	84	84 Neg / 0 Pos
MAMP	-25 %	84	83 Neg / 1 Pos
MAMP	cutoff	84	27 Neg / 57 Pos
MAMP	+25 %	84	0 Neg / 84 Pos
MAMP	+50 %	84	0 Neg / 84 Pos
MAMP	+75 %	84	0 Neg / 84 Pos
MAMP	+100 %	84	0 Neg / 84 Pos

Qualitative - 500 ng/mL Cutoff

Drug	Concentration of Sample	Number of Determinations	Results # Neg / # Pos
MAMP	zero drug	84	84 Neg / 0 Pos
MAMP	-75 %	84	84 Neg / 0 Pos
MAMP	-50 %	84	84 Neg / 0 Pos
MAMP	-25 %	84	84 Neg / 0 Pos
MAMP	cutoff	84	28 Neg / 56 Pos
MAMP	+25 %	84	0 Neg / 84 Pos
MAMP	+50 %	84	0 Neg / 84 Pos
MAMP	+75 %	84	0 Neg / 84 Pos
MAMP	+100 %	84	0 Neg / 84 Pos

Qualitative - 1000 ng/mL Cutoff

Drug	Concentration of Sample	Number of Determinations	Results # Neg / # Pos
MAMP	zero drug	84	84 Neg / 0 Pos
MAMP	-75 %	84	84 Neg / 0 Pos
MAMP	-50 %	84	84 Neg / 0 Pos
MAMP	-25 %	84	84 Neg / 0 Pos
MAMP	cutoff	84	13 Neg / 71 Pos
MAMP	+25 %	84	0 Neg / 84 Pos
MAMP	+50 %	84	0 Neg / 84 Pos
MAMP	+75 %	84	0 Neg / 84 Pos
MAMP	+100 %	84	0 Neg / 84 Pos

Semiquantitative - 300 ng/mL Cutoff

Drug	Sample Conc.	Results # Neg / # Pos	Repeatability		Intermediate Precision	
			SD ng/mL	CV %	SD ng/mL	CV %
MAMP	zero drug	84 / 0	9.4	95.4	12.2	123.8
MAMP	-75 %	84 / 0	12.9	15.5	15.7	18.8
MAMP	-50 %	84 / 0	14.3	9.0	17.7	11.2
MAMP	-25 %	84 / 0	13.2	5.4	19.8	8.2
MAMP	cutoff	34 / 50	23.2	7.6	23.5	7.7
MAMP	+25 %	0 / 84	20.9	5.2	24.4	6.1
MAMP	+50 %	0 / 84	22.5	5.0	29.5	6.6
MAMP	+75 %	0 / 84	27.4	5.0	36.5	6.6
MAMP	+100 %	0 / 84	35.7	5.7	40.8	6.5

Semiquantitative - 500 ng/mL Cutoff

Drug	Sample Conc.	Results # Neg / # Pos	Repeatability		Intermediate Precision	
			SD ng/mL	CV %	SD ng/mL	CV %
MAMP	zero drug	84 / 0	9.6	117.3	11.0	134.1
MAMP	-75 %	84 / 0	14.6	10.5	20.6	14.9
MAMP	-50 %	84 / 0	20.0	7.9	21.9	8.6
MAMP	-25 %	84 / 0	24.1	6.2	27.2	7.0
MAMP	cutoff	39 / 45	28.2	5.6	39.4	7.8
MAMP	+25 %	0 / 84	31.5	4.9	48.5	7.5

Drug	Sample Conc.	Results # Neg / # Pos	Repeatability		Intermediate Precision	
			SD ng/mL	CV %	SD ng/mL	CV %
MAMP	+50 %	0 / 84	34.5	4.8	49.5	6.8
MAMP	+75 %	0 / 84	41.6	5.0	62.3	7.5
MAMP	+100 %	0 / 84	46.7	4.6	72.4	7.1

Semiquantitative - 1000 ng/mL Cutoff

Drug	Sample Conc.	Results # Neg / # Pos	Repeatability		Intermediate Precision	
			SD ng/mL	CV %	SD ng/mL	CV %
MAMP	zero drug	84 / 0	20.2	97.2	24.9	119.5
MAMP	-75 %	84 / 0	25.3	9.3	30.7	11.2
MAMP	-50 %	84 / 0	33.3	6.4	39.4	7.5
MAMP	-25 %	84 / 0	34.3	4.6	52.1	7.0
MAMP	cutoff	25 / 59	56.9	5.5	71.3	6.9
MAMP	+25 %	0 / 84	56.6	4.2	79.3	5.9
MAMP	+50 %	0 / 84	63.6	4.0	81.3	5.2
MAMP	+75 %	0 / 84	60.3	3.5	101.6	5.9
MAMP	+100 %	0 / 84	108.1	5.1	126.6	6.0

The same precision experiment was repeated utilizing *D*-amphetamine (AMP) as the target analyte instead of *D*-methamphetamine. The following tables show the results obtained on a COBAS INTEGRA 800 analyzer.

Qualitative - 300 ng/mL Cutoff

Drug	Concentration of Sample	Number of Determinations	Results # Neg / # Pos
AMP	zero drug	84	84 Neg / 0 Pos
AMP	-75 %	84	84 Neg / 0 Pos
AMP	-50 %	84	84 Neg / 0 Pos
AMP	-25 %	84	83 Neg / 1 Pos
AMP	cutoff	84	2 Neg / 82 Pos
AMP	+25 %	84	0 Neg / 84 Pos
AMP	+50 %	84	0 Neg / 84 Pos
AMP	+75 %	84	0 Neg / 84 Pos
AMP	+100 %	84	0 Neg / 84 Pos

Qualitative - 500 ng/mL Cutoff

Drug	Concentration of Sample	Number of Determinations	Results # Neg / # Pos
AMP	zero drug	84	84 Neg / 0 Pos
AMP	-75 %	84	84 Neg / 0 Pos
AMP	-50 %	84	84 Neg / 0 Pos
AMP	-25 %	84	82 Neg / 2 Pos
AMP	cutoff	84	6 Neg / 78 Pos
AMP	+25 %	84	0 Neg / 84 Pos
AMP	+50 %	84	0 Neg / 84 Pos
AMP	+75 %	84	0 Neg / 84 Pos
AMP	+100 %	84	0 Neg / 84 Pos

Qualitative - 1000 ng/mL Cutoff

Drug	Concentration of Sample	Number of Determinations	Results # Neg / # Pos
AMP	zero drug	84	84 Neg / 0 Pos
AMP	-75 %	84	84 Neg / 0 Pos
AMP	-50 %	84	84 Neg / 0 Pos
AMP	-25 %	84	82 Neg / 2 Pos
AMP	cutoff	84	7 Neg / 77 Pos
AMP	+25 %	84	1 Neg / 83 Pos
AMP	+50 %	84	0 Neg / 84 Pos
AMP	+75 %	84	0 Neg / 84 Pos
AMP	+100 %	84	0 Neg / 84 Pos

Semiquantitative - 300 ng/mL Cutoff

Drug	Sample Conc.	Results # Neg / # Pos	Repeatability		Intermediate Precision	
			SD ng/mL	CV %	SD ng/mL	CV %
AMP	zero drug	84 / 0	8.6	110.0	10.6	135.7
AMP	-75 %	84 / 0	10.4	16.8	13.6	21.9
AMP	-50 %	84 / 0	14.6	10.0	19.9	13.6
AMP	-25 %	84 / 0	15.7	6.6	20.6	8.7
AMP	cutoff	4 / 80	25.8	7.4	28.5	8.2
AMP	+25 %	0 / 84	30.7	7.1	32.2	7.4
AMP	+50 %	0 / 84	31.1	6.1	37.6	7.3
AMP	+75 %	0 / 84	29.1	4.9	41.2	7.0
AMP	+100 %	0 / 84	36.2	5.3	44.6	6.5

Semiquantitative - 500 ng/mL Cutoff

Drug	Sample Conc.	Results # Neg / # Pos	Repeatability		Intermediate Precision	
			SD ng/mL	CV %	SD ng/mL	CV %
AMP	zero drug	84 / 0	12.6	117.0	13.2	122.5
AMP	-75 %	84 / 0	21.4	16.9	25.2	20.0
AMP	-50 %	84 / 0	18.2	6.9	25.6	9.8
AMP	-25 %	83 / 1	25.6	5.7	29.7	6.6
AMP	cutoff	0 / 84	29.1	4.7	38.3	6.2
AMP	+25 %	0 / 84	42.0	5.5	51.7	6.8
AMP	+50 %	0 / 84	40.3	4.6	49.9	5.7
AMP	+75 %	0 / 84	58.7	5.9	59.7	6.0
AMP	+100 %	0 / 84	79.0	7.0	86.3	7.6

Semiquantitative - 1000 ng/mL Cutoff

Drug	Sample Conc.	Results # Neg / # Pos	Repeatability		Intermediate Precision	
			SD ng/mL	CV %	SD ng/mL	CV %
AMP	zero drug	84 / 0	15.9	111.1	19.0	132.6
AMP	-75 %	84 / 0	21.9	9.2	29.8	12.5
AMP	-50 %	84 / 0	30.3	5.4	40.9	7.3
AMP	-25 %	84 / 0	44.3	5.2	55.0	6.4

Drug	Sample Conc.	Results # Neg / # Pos	Repeatability		Intermediate Precision	
			SD ng/mL	CV %	SD ng/mL	CV %
AMP	cutoff	6 / 78	53.4	4.8	84.3	7.6
AMP	+25 %	0 / 84	72.5	5.6	87.5	6.7
AMP	+50 %	0 / 84	65.7	4.4	96.6	6.5
AMP	+75 %	0 / 84	74.3	4.5	105.2	6.4
AMP	+100 %	0 / 84	106.7	5.9	131.4	7.3

Accuracy

The accuracy of this assay was determined against *d*-methamphetamine or *d*-amphetamine GC/MS results. The evaluated cutoff concentrations for the GC/MS testing were the same as the assay screening cutoffs. For both *d*-methamphetamine and *d*-amphetamine, 36 urine samples, obtained from a clinical laboratory were confirmed to be negative by GC/MS and were evaluated with Amphetamines II. 100 % of these normal urines were negative with both the semiquantitative and qualitative assay relative to the 300 ng/mL, 500 ng/mL, and 1000 ng/mL cutoffs. 4 additional unaltered samples, containing *d*-methamphetamine between 85-204 ng/mL (GC/MS confirmed), were tested with the 500 ng/mL cutoff and gave negative results.

36 samples obtained from a clinical laboratory, where they screened preliminary positive with a commercially available enzyme immunoassay and were subsequently confirmed positive by GC/MS to contain *d*-methamphetamine or *d*-amphetamine, were evaluated with Amphetamines II. 100 % of these samples were positive relative to the 300 ng/mL, 500 ng/mL, and 1000 ng/mL cutoffs.

In addition, several unaltered near cutoff samples were run for each cutoff. These samples fell in the near cutoff negative range (between -50 % and cutoff) and the near cutoff positive range (between cutoff and +50 %). For each cutoff, 4 negative near cutoff samples and 4 positive near cutoff samples were assayed.

Data from the accuracy studies described above were combined with data generated from the unaltered near cutoff urine samples. The following results were obtained with the Amphetamines II assay on the COBAS INTEGRA 800 analyzer relative to the GC/MS values for both *d*-methamphetamine and *d*-amphetamine.

Amphetamines II Qualitative Assay Results (MAMP)

Roche ONLINE DAT AMPII assay	Low Neg	Near Cutoff Negative by GC/MS (between -50 % and cutoff)	Near Cutoff Positive by GC/MS (between cutoff and +50 %)	High Positive by GC/MS (greater than +50 %)	Percent Agreement with GC/MS (MAMP)
300 ng/mL Cutoff					
Positive	0	4	4	36	100 %
Negative	36	0	0	0	90 %
500 ng/mL Cutoff					
Positive	0	4	4	36	100 %
Negative	40	0	0	0	91 %
1000 ng/mL Cutoff					
Positive	0	4	4	36	100 %
Negative	36	0	0	0	90 %

Amphetamines II Semiquantitative Assay Results (MAMP)

Roche ONLINE DAT AMPII assay	Low Neg	Near Cutoff Negative by GC/MS (between -50 % and cutoff)	Near Cutoff Positive by GC/MS (between cutoff and +50 %)	High Positive by GC/MS (greater than +50 %)	Percent Agreement with GC/MS (MAMP)
300 ng/mL Cutoff					
Positive	0	4	4	36	100 %
Negative	36	0	0	0	90 %
500 ng/mL Cutoff					
Positive	0	4	4	36	100 %
Negative	40	0	0	0	91 %
1000 ng/mL Cutoff					
Positive	0	4	4	36	100 %
Negative	36	0	0	0	90 %

Accuracy samples were categorized based upon the *d*-methamphetamine GC/MS concentration only. The table below identifies those samples with a *d*-methamphetamine concentration below the cutoff, in which the observed result on a COBAS INTEGRA 800 analyzer was positive. The expected results column identifies the result expected with the Amphetamines II assay based upon the *d*-methamphetamine (MAMP) value relative to the cutoff.

GC/MS Summary of Discrepant Results (MAMP)

Cutoff Value (ng/mL)	Roche ONLINE DAT AMPII OBSERVED Result	Roche ONLINE DAT AMPII EXPECTED Result	GC/MS (ng/mL)	Drug / Metabolite
300 (SQ, Q)	Positive	Positive	173 181	MAMP AMP
300 (SQ, Q)	Positive	Positive	278 101	MAMP AMP
300 (SQ, Q)	Positive	Positive	171 220	MAMP AMP
300 (SQ, Q)	Positive	Positive	291 145	MAMP AMP
500 (SQ, Q)	Positive	Positive	488 466	MAMP AMP
500 (SQ, Q)	Positive	Negative	325 171	MAMP AMP
500 (SQ, Q)	Positive	Negative	291 145	MAMP AMP
500 (SQ, Q)	Positive	Positive	472 650	MAMP AMP
1000 (SQ, Q)	Positive	Positive	706 443	MAMP AMP
1000 (SQ, Q)	Positive	Positive	540 693	MAMP AMP
1000 (SQ, Q)	Positive	Positive	769 395	MAMP AMP
1000 (SQ, Q)	Positive	Positive	572 432	MAMP AMP

Amphetamines II Qualitative Assay Results (AMP)

Roche ONLINE DAT AMPII assay	Low Neg	Near Cutoff Negative by GC/MS (between -50 % and cutoff)	Near Cutoff Positive by GC/MS (between cutoff and +50 %)	High Positive by GC/MS (greater than +50 %)	Percent Agreement with GC/MS (AMP)
300 ng/mL Cutoff					
Positive	0	3	4	36	100 %
Negative	36	1	0	0	92.5 %
500 ng/mL Cutoff					
Positive	0	3	4	36	100 %
Negative	36	1	0	0	92.5 %
1000 ng/mL Cutoff					
Positive	0	3	4	36	100 %
Negative	36	1	0	0	97.5 %

Amphetamines II Semiquantitative Assay Results (AMP)

Roche ONLINE DAT AMPII assay	Low Neg	Near Cutoff Negative by GC/MS (between -50 % and cutoff)	Near Cutoff Positive by GC/MS (between cutoff and +50 %)	High Positive by GC/MS (greater than +50 %)	Percent Agreement with GC/MS (AMP)
300 ng/mL Cutoff					
Positive	0	3	4	36	100 %
Negative	36	1	0	0	92.5 %
500 ng/mL Cutoff					
Positive	0	3	4	36	100 %
Negative	36	1	0	0	92.5 %
1000 ng/mL Cutoff					
Positive	0	2	4	36	100 %
Negative	36	2	0	0	95 %

Accuracy samples were categorized based upon the *d*-amphetamine GC/MS concentration only. The table below identifies those samples with a *d*-amphetamine concentration below the cutoff, in which the observed result on a COBAS INTEGRA 800 analyzer was positive. The expected results column identifies the result expected with the Amphetamines II assay based upon the *d*-amphetamine (AMP) value relative to the cutoff.

GC/MS Summary of Discrepant Results (AMP)

Cutoff Value (ng/mL)	Roche ONLINE DAT AMPII OBSERVED Result	Roche ONLINE DAT AMPII EXPECTED Result	GC/MS (ng/mL)	Drug / Metabolite
300 (SQ, Q)	Positive	Positive	157 363	AMP MAMP
300 (SQ, Q)	Positive	Positive	181 173	AMP MAMP
300 (SQ, Q)	Positive	Positive	220 171	AMP MAMP
500 (SQ, Q)	Positive	Positive	438 121	AMP MAMP

Cutoff Value (ng/mL)	Roche ONLINE DAT AMPII OBSERVED Result	Roche ONLINE DAT AMPII EXPECTED Result	GC/MS (ng/mL)	Drug / Metabolite
500 (SQ, Q)	Positive	Positive	457 1152	AMP MAMP
500 (SQ, Q)	Positive	Positive	443 706	AMP MAMP
1000 (SQ, Q)	Positive	Negative	920	AMP
1000 (SQ, Q)	Positive	Positive	837 1163	AMP MAMP

Analytical specificity

The specificity of Amphetamines II, on the COBAS INTEGRA 800 analyzer, for various phenethylamines and structurally similar compounds was determined by generating inhibition curves for each of the compounds listed for both semiquantitative and qualitative modes and determining the approximate quantity of each compound that is equivalent in assay reactivity to the 300 ng/mL, 500 ng/mL, and 1000 ng/mL *d*-methamphetamine assay cutoff. The tables below show the semiquantitative results of the study for each assay cutoff. The same samples were run in the qualitative mode and all recovered appropriately negative or positive, based on the calculated cross-reactivity.

Compound	ng/mL Equivalent to 300 ng/mL MAMP	Approximate percent cross-reactivity
± MDMA ^{c)}	114	264
PMA ^{d)}	143	210
PMMA ^{e)}	181	166
± MDA ^{f)}	249	121
± MDEA ⁱ⁾	285	105
<i>d</i> -Amphetamine	311	96 ^{g)}
<i>d</i> -Methamphetamine	327	92
± MBDB HCl ^{h)}	339	88
± BDB HCl ^{j)}	648	46
meta-Chlorophenylpiperazine ^{k)}	1322	23
1-Methyl-3-phenylpropylamine ^{l)}	1792	17
<i>l</i> -Methamphetamine	2754	11
<i>l</i> -Amphetamine	6443	5
Dimethylamylamine ^{m)}	21111	1.4
Phendimetrazine	47473	0.63
Phentermine	66385	0.45
Tyramine	86271	0.35
<i>d</i> -Pseudoephedrine	87578	0.34
<i>l</i> -Ephedrine	94792	0.32
<i>d,l</i> -Phenylpropanolamine HCl	280374	0.11
<i>d</i> -Ephedrine	329670	0.09

Compound	ng/mL Equivalent to 500 ng/mL MAMP	Approximate percent cross-reactivity
± MDMA ^{c)}	173	289
PMA ^{d)}	231	216
PMMA ^{e)}	299	167
± MDA ^{f)}	433	115

Compound	ng/mL Equivalent to 500 ng/mL MAMP	Approximate percent cross-reactivity
<i>d</i> -Methamphetamine	444	113
<i>d</i> -Amphetamine	460	109 ^{g)}
± MDEA ⁱ⁾	494	101
± MBDB HCl ^{h)}	713	70
± BDB HCl ^{j)}	1209	41
meta-Chlorophenylpiperazine ^{k)}	2414	21
1-Methyl-3-phenylpropylamine ^{l)}	2925	17
<i>l</i> -Methamphetamine	4098	12
<i>l</i> -Amphetamine	11174	4
Phendimetrazine	72500	0.69
Phentermine	118483	0.42
Dimethylamylamine ^{m)}	37701	1.3
<i>d</i> -Pseudoephedrine	132275	0.38
Tyramine	141243	0.35
<i>l</i> -Ephedrine	154321	0.32
<i>d,l</i> -Phenylpropanolamine HCl	390625	0.13
<i>d</i> -Ephedrine	413223	0.12

Compound	ng/mL Equivalent to 1000 ng/mL MAMP	Approximate percent cross-reactivity
± MDMA ^{c)}	446	224
PMA ^{d)}	599	167
PMMA ^{e)}	695	144
± MDA ^{f)}	785	127
<i>d</i> -Methamphetamine	970	103
<i>d</i> -Amphetamine	1024	98 ^{g)}
± MBDB HCl ^{h)}	1194	84
± MDEA ⁱ⁾	1203	83
± BDB HCl ^{j)}	2262	44
meta-Chlorophenylpiperazine ^{k)}	4946	20
1-Methyl-3-phenylpropylamine ^{l)}	5479	18
<i>l</i> -Methamphetamine	9008	11
<i>l</i> -Amphetamine	23445	4
Dimethylamylamine ^{m)}	79780	1.3
Phendimetrazine	156740	0.64
<i>d</i> -Pseudoephedrine	269542	0.37
Phentermine	294118	0.34
<i>l</i> -Ephedrine	317460	0.32
Tyramine	323625	0.31
<i>d</i> -Ephedrine	793651	0.13
<i>d,l</i> -Phenylpropanolamine HCl	1111111	0.09

c) *d,l*-3,4-Methylenedioxymethamphetamine

d) para-Methoxyamphetamine

e) para-Methoxymethamphetamine

f) *d,l*-3,4-Methylenedioxymethamphetamine

g) Representative data from multiple lots demonstrate cross-reactivity in the range from approximately 75-125 %

h) *d,l*-N-Methyl-1-(3,4-methylenedioxyphenyl)-2-butanamine hydrochloridei) *d,l*-3,4-Methylenedioxyethylamphetamine

- j) *d,l*-3,4-Methylenedioxyphenyl-2-butanamine hydrochloride
 k) 1-(3-Chlorophenyl)piperazine (mCPP), major metabolite of Trazodone
 l) APB, metabolite of Labetalol
 m) 4-Methylhexan-2-amine, DMAA

Interference and cross-reactivity with unrelated drugs

The following compounds were added at the listed concentrations to a human urine pool spiked with either *d*-amphetamine or *d*-methamphetamine at approximately the negative and positive control concentrations for each cutoff ($\pm 25\%$ of assay cutoff). For each compound, the control level samples recovered properly for the 300 ng/mL, 500 ng/mL, and 1000 ng/mL cutoff in both semiquantitative and qualitative modes.

Compound	Conc. (ng/mL)	Semiquantitative All Cutoffs		Qualitative All Cutoffs	
		Low Ctrl	High Ctrl	Low Ctrl	High Ctrl
Acetaminophen	100000	NEG	POS	NEG	POS
Acetylsalicylic acid	100000	NEG	POS	NEG	POS
Amitriptyline	100000	NEG	POS	NEG	POS
Aspartame	40000	NEG	POS	NEG	POS
Benzocaine	100000	NEG	POS	NEG	POS
Benzoylcegonine	100000	NEG	POS	NEG	POS
Caffeine	100000	NEG	POS	NEG	POS
Cannabidiol	100000	NEG	POS	NEG	POS
Cocaine	100000	NEG	POS	NEG	POS
Codeine	100000	NEG	POS	NEG	POS
Desipramine HCl	100000	NEG	POS	NEG	POS
Dextromethorphan	100000	NEG	POS	NEG	POS
Dextropropoxyphene	100000	NEG	POS	NEG	POS
Diazepam	100000	NEG	POS	NEG	POS
Digoxin	100000	NEG	POS	NEG	POS
Diphenhydramine	100000	NEG	POS	NEG	POS
Diphenylhydantoin	100000	NEG	POS	NEG	POS
Doxepin	100000	NEG	POS	NEG	POS
Ecgonine	100000	NEG	POS	NEG	POS
Ecgonine methyl ester	100000	NEG	POS	NEG	POS
Erythromycin	100000	NEG	POS	NEG	POS
Furosemide	100000	NEG	POS	NEG	POS
Guaiacol glycerol ether	100000	NEG	POS	NEG	POS
Hydrochlorothiazide	100000	NEG	POS	NEG	POS
Ibuprofen	100000	NEG	POS	NEG	POS
Ketamine	100000	NEG	POS	NEG	POS
Levothyroxine	100000	NEG	POS	NEG	POS
LSD	2500	NEG	POS	NEG	POS
Meperidine	100000	NEG	POS	NEG	POS
Methadone	100000	NEG	POS	NEG	POS
Methaqualone	75000	NEG	POS	NEG	POS
Morphine	100000	NEG	POS	NEG	POS
Naloxone	100000	NEG	POS	NEG	POS
Naltrexone	100000	NEG	POS	NEG	POS
Naproxen	100000	NEG	POS	NEG	POS
Niacinamide	100000	NEG	POS	NEG	POS
Nicotine	100000	NEG	POS	NEG	POS

Nifedipine	100000	NEG	POS	NEG	POS
Nordiazepam	100000	NEG	POS	NEG	POS
Omeprazole	100000	NEG	POS	NEG	POS
Oxazepam	100000	NEG	POS	NEG	POS
Penicillin G	100000	NEG	POS	NEG	POS
Phencyclidine	40000	NEG	POS	NEG	POS
Phenobarbital	100000	NEG	POS	NEG	POS
Quinine	100000	NEG	POS	NEG	POS
Secobarbital	100000	NEG	POS	NEG	POS
Tetracycline	100000	NEG	POS	NEG	POS
Δ^9 -THC	10000	NEG	POS	NEG	POS

The compounds, including methylphenidate (Ritalin), were additionally added to aliquots of pooled drug-free human urine at a concentration of 100000 ng/mL. None of these compounds gave values in the assay that were equal to or greater than 0.19 % cross-reactivity and no results were greater than the assay cutoffs (300 ng/mL, 500 ng/mL, and 1000 ng/mL), with the following exception.

The compounds Labetalol HCl and Trazodone were additionally added to aliquots of pooled drug-free human urine at a concentration of 100000 ng/mL. Ranitidine was added to aliquots of pooled drug-free human urine at a concentration of 50000 ng/mL. The results obtained were between 0.21 % and 0.36 % for the 300 ng/mL, 500 ng/mL and the 1000 ng/mL assay cutoffs respectively.

The cross-reactivity for LSD was tested at a concentration of 2500 ng/mL. The results obtained were 0.32 % and 0.71 %, for the 300 ng/mL and 1000 ng/mL assay cutoffs respectively.

The cross-reactivity for procaine HCl was tested in drug-free urine at a concentration of 100000 ng/mL. The results obtained were 1.00 %, 0.99 %, and 0.90 %, for the 300 ng/mL, 500 ng/mL, and 1000 ng/mL assay cutoffs respectively.

Interference with other substances

Interfering substances were added to urine containing *d*-methamphetamine (MAMP) at -25% and $+25\%$ of the cutoff level at the concentration listed below. The same substances were additionally added to urine containing *d*-amphetamine (AMP) at -25% and $+25\%$ of the cutoff level at the concentration listed below. All samples were tested and the following results were obtained on a COBAS INTEGRA 800 analyzer. The value in the table indicates the level at which no interference was found for samples containing either *d*-methamphetamine or *d*-amphetamine.

Qualitative		300 ng/mL Cutoff		500 ng/mL Cutoff		1000 ng/mL Cutoff	
Compound	Cmpd. Conc.	Neg Level	Pos Level	Neg Level	Pos Level	Neg Level	Pos Level
Acetone	7.9 mg/mL	Neg	Pos	Neg	Pos	Neg	Pos
Ascorbic Acid	10 mg/mL	Neg	Pos	Neg	Pos	Neg	Pos
Conjugated Bilirubin	0.1 mg/mL	Neg	Pos	Neg	Pos	Neg	Pos
Creatinine	5 mg/mL	Neg	Pos	Neg	Pos	Neg	Pos
Ethanol	7.9 mg/mL	Neg	Pos	Neg	Pos	Neg	Pos
Glucose	12 mg/mL	Neg	Pos	Neg	Pos	Neg	Pos
Hemoglobin	1 mg/mL	Neg	Pos	Neg	Pos	Neg	Pos
Human serum albumin	3 mg/mL	Neg	Pos	Neg	Pos	Neg	Pos
Oxalic Acid	2 mg/mL	Neg	Pos	Neg	Pos	Neg	Pos
Sodium Chloride	23 mg/mL	Neg	Pos	Neg	Pos	Neg	Pos
Urea	60 mg/mL	Neg	Pos	Neg	Pos	Neg	Pos

The same experiment was performed in the semiquantitative mode for each cutoff. All negative and positive controls recovered properly in the presence of the interfering substance.

A protocol was executed in which samples containing either AMP or MAMP at control levels ($\pm 25\%$ of cutoff) with specific gravities ranging from 1.001 to 1.034 were tested. As with the other interferences, there were no control cross-overs on any of the 3 assay cutoffs at either extreme specific gravity level. Samples having a pH ranging from 4.5 to 8.0 and containing either AMP or MAMP at control levels ($\pm 25\%$ of cutoff) were also tested. There were $\leq 5\%$ cross-overs on any of the 3 assay cutoffs.




Any modification of the instrument as set forth in this labeling requires validation by the laboratory.

References

- 1 Karch SB, ed. Drug Abuse Handbook. Boca Raton, FL: CRC Press LLC 1998.
- 2 Borenstein M. Central Nervous System Stimulants. In: Troy D. ed. Remington: The Science and Practice of Pharmacy. 21st ed. Baltimore, MD: Lippincott Williams & Wilkins Co 2005;1551.
- 3 NIDA Research Report - Methamphetamine Abuse and Addiction: NIH Publication No. 06-4210. National Institute on Drug Abuse 6001 Executive Blvd. Room 5213, Bethesda, MD 20892-9561.
- 4 Baselt RC. Disposition of Toxic Drugs and Chemicals in Man. 7th ed. Foster City, CA: Biomedical Publications 2004;67.
- 5 Armbruster DA, Schwarzhoff RH, Pierce BL, et al. Method comparison of EMIT II and ONLINE with RIA for drug screening. J Forensic Sci 1993;38:1326-1341.
- 6 Armbruster DA, Schwarzhoff RH, Hubster EC, et al. Enzyme immunoassay, kinetic microparticle immunoassay, radioimmunoassay, and fluorescence polarization immunoassay compared for drugs-of-abuse screening. Clin Chem 1993;39:2137-2146.
- 7 Toxicology and Drug Testing in the Clinical Laboratory; Approved Guideline. 2nd ed. (C52-A2). Clinical and Laboratory Standards Institute 2007;27:33.
- 8 Mandatory Guidelines for Federal Workplace Drug Testing Programs. Fed Regist 2008 Nov 25;73:71858-71907.

Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard.

	Contents of kit
	Volume after reconstitution or mixing
	Global Trade Item Number

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

FOR US CUSTOMERS ONLY: LIMITED WARRANTY

Roche Diagnostics warrants that this product will meet the specifications stated in the labeling when used in accordance with such labeling and will be free from defects in material and workmanship until the expiration date printed on the label. THIS LIMITED WARRANTY IS IN LIEU OF ANY OTHER WARRANTY, EXPRESS OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. IN NO EVENT SHALL ROCHE DIAGNOSTICS BE LIABLE FOR INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES.

COBAS, COBAS C, COBAS INTEGRA, ONLINE DAT and PRECISET are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2015, Roche Diagnostics



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim
www.roche.com

Distribution in USA by:
Roche Diagnostics, Indianapolis, IN
US Customer Technical Support 1-800-428-2336

