

# Cleaner

## Order information

REF	CONTENT	Analyzer(s) on which reagent(s) can be used
20754765 322	Cleaner (1000 mL)	COBAS INTEGRA 400 plus COBAS INTEGRA 800 cobas c 111

## English

### Intended use

Cleaner is a wash solution for sample and reagent probes and tubing.

### Summary

For maintaining the integrity of sample and reagent probes and tubing, washing is required. Cleaner is used as a wash solution to prevent potential carry over from sample and reagent probes and tubing.

### Reagents - working solutions

HCl 0.3 mol/L

### Precautions and warnings

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents.

Disposal of all waste material should be in accordance with local guidelines. Safety data sheet available for professional user on request.

This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:

Hazardous components	Hydrogen chloride Isodecanol, ethoxylated
----------------------	--



Danger

H290	May be corrosive to metals.
H314	Causes severe skin burns and eye damage.

### Prevention

P280	Wear protective gloves/ protective clothing/ eye protection/ face protection.
------	---

### Response

P301 + P330 + P331	IF SWALLOWED: rinse mouth. Do NOT induce vomiting.
P303 + P361 + P353	IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower.
P305 + P351 + P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P310	Immediately call a POISON CENTER or doctor/physician.
P390	Absorb spillage to prevent material damage.

Contact phone: all countries: +49-621-7590, USA: +1-800-428-2336

### Reagent handling

Ready for use

### Storage and stability

Shelf life at 15-25 °C	See expiration date on bottle label.
------------------------	--------------------------------------

COBAS INTEGRA 400 plus/800 analyzers  
on-board in use 12 weeks

cobas c 111 analyzer  
in use 12 weeks

### Materials provided

See "Reagents – working solutions" section for reagents.

### Materials required (but not provided)

General laboratory equipment

Refer to the appropriate operator's manual for the additional required materials.

### Assay

Use Cleaner as specified in the operator's manuals of the respective analyzer.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

### Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard.

CONTENT	Contents of kit
REAGENT	Reagent
CALIBRATOR	Calibrator
→	Volume after reconstitution or mixing

### FOR US CUSTOMERS ONLY: LIMITED WARRANTY

Roche Diagnostics warrants that this product will meet the specifications stated in the labeling when used in accordance with such labeling and will be free from defects in material and workmanship until the expiration date printed on the label. THIS LIMITED WARRANTY IS IN LIEU OF ANY OTHER WARRANTY, EXPRESS OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. IN NO EVENT SHALL ROCHE DIAGNOSTICS BE LIABLE FOR INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES.

COBAS INTEGRA, COBAS and COBAS C are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Significant additions or changes are indicated by a change bar in the margin.

© 2013, Roche Diagnostics



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim  
www.roche.com

Distribution in USA by:  
Roche Diagnostics, Indianapolis, IN  
US Customer Technical Support 1-800-428-2336

