

ETOH2

Ethanol Gen.2

Order information

REF	CONTENT	Analyzer(s) on which cobas c pack(s) can be used
03183777 190	Ethanol Gen.2 (100 tests)	System-ID 07 6611 9 COBAS INTEGRA 400 plus COBAS INTEGRA 800
20751995 190	Ammonia/Ethanol/CO ₂ Calibrator (2 × 4 mL)	System-ID 07 5199 5
20752401 190	Ammonia/Ethanol/CO ₂ Control Normal (5 × 4 mL)	System-ID 07 5240 1
20753009 190	Ammonia/Ethanol/CO ₂ Control Abnormal (5 × 4 mL)	System-ID 07 5300 9

English

System information

Test ETOH2, test ID 0-611 (serum, plasma)

Test ETOU2, test ID 0-511 (urine)

Intended use

In vitro test for the quantitative determination of ethanol in human serum, plasma, and urine on COBAS INTEGRA systems.

Summary


Ethyl alcohol determinations are among the most frequent analyses required in the forensic and clinical toxicology laboratory. Ethyl alcohol measurements are used in the diagnosis and treatment of alcohol intoxication and poisoning.

Early techniques for blood alcohol determination used distillation, aeration, or diffusion to separate the alcohol from the plasma matrix. The distilled alcohol was then measured by oxidation of the alcohol by strong oxidizing agents. However, these methods lacked specificity, since other oxidizable compounds could also be distilled into and react in the reaction mixture.¹ While there are many acceptable published procedures, including gas chromatographic and osmometric methods, the enzymatic technique described below, based on the information given by Bucher and Redetzki², is specific and simple to perform.

Test principle

Enzymatic method with alcohol dehydrogenase

Ethyl alcohol and NAD are converted to acetaldehyde and NADH by ADH.



The NADH formed during the reaction, measured photometrically as a rate of change in absorbance, is directly proportional to the ethyl alcohol concentration. It is determined by measuring the increase in absorbance at 340 nm.

Reagents - working solutions

R1 Buffer; preservatives

SR NAD (yeast): ≥ 3 mmol/L; ADH (EC 1.1.1.1, yeast, 25 °C): ≥ 37 U/mL; stabilizers; preservatives

R1 is in position B and SR is in position C.

Precautions and warnings

Pay attention to all precautions and warnings listed in Section 1 / Introduction of this Method Manual.

For USA: Caution: Federal law restricts this device to sale by or on the order of a physician.

Reagent handling

Ready for use

Storage and stability

Shelf life at 2-8 °C	See expiration date on cobas c pack label
COBAS INTEGRA 400 plus system	
On-board in use at 10-15 °C	4 weeks
COBAS INTEGRA 800 system	
On-board in use at 8 °C	12 weeks

Specimen collection and preparation^{3,4}

Do not use alcohol or other volatile disinfectants at the site of venipuncture. Aqueous Zephiran (benzalkonium chloride), aqueous Merthiolate (thimerosal), or povidone-iodine may be used.

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable.

Serum, plasma: Li-, Na-, NH₄⁺-heparin and K₂-, K₃-EDTA

Stability: ⁵	2 days at 25 °C
	2 weeks at 5 °C
	4 weeks at -15 °C

Plasma: NaF/Na₂EDTA and NaF/K-Oxalate

Stability: ⁵	2 weeks at 25 °C
	3 months at 5 °C
	6 months at -15 °C

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Urine: Use random urine.

Stability: ⁶	30 days at 4 °C
-------------------------	-----------------

Storage: Samples must be tightly closed.

Centrifuge samples containing precipitates before performing the assay.

Materials provided

See "Reagents – working solutions" section for reagents.

Assay

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

Repeat assays must be performed on freshly poured cups, due to evaporation of alcohol.

When using Ammonia/Ethanol/CO₂ Calibrator: Do not leave calibrator cups open for longer than 30 minutes at 15-25 °C.

When using Ammonia/Ethanol/CO₂ Controls: Do not leave control cups open for longer than 1 hour at 15-25 °C.

Application for serum, plasma and urine

COBAS INTEGRA 400 plus test definition

Measuring mode	Absorbance
Abs. calculation mode	Kinetic
Reaction mode	R1-S-SR
Reaction direction	Increase
Wavelength A/B	340/659 nm
Calc. first/last	44/54
Unit	mmol/L

ETOH2

Ethanol Gen.2

cobas[®]
Drug abuse testing

Pipetting parameters

<i>Serum, plasma, urine</i>		Diluent (H ₂ O)
R1	50 µL	-
Sample	4 µL	16 µL
SR	50 µL	-
Total volume	120 µL	

COBAS INTEGRA 800 test definition

Measuring mode	Absorbance
Abs. calculation mode	Kinetic
Reaction mode	R1-S-SR
Reaction direction	Increase
Wavelength A/B	340/659 nm
Calc. first/last	62/79
Unit	mmol/L

Pipetting parameters

<i>Serum, plasma, urine</i>		Diluent (H ₂ O)
R1	50 µL	-
Sample	4 µL	16 µL
SR	50 µL	-
Total volume	120 µL	

Calibration

Calibrator	Roche Ammonia/Ethanol/CO ₂ Calibrator Use deionized water as zero calibrator.
Calibration mode	Linear regression
Calibration replicate	Duplicate recommended
Calibration interval	COBAS INTEGRA 400 plus system: Each cobas c pack and as required following quality control procedures COBAS INTEGRA 800 system: Each cobas c pack, every 6 weeks, and as required following quality control procedures

Calibration interval may be extended based on acceptable verification of calibration by the laboratory.

Traceability: This method has been standardized against NIST-traceable standard materials.

Quality control

Quality control	Ammonia/Ethanol/CO ₂ Control Normal and Abnormal
Control interval	8 hours recommended
Control sequence	User defined
Control after calibration	Recommended

For quality control, use control materials as listed in the "Order information" section. In addition, other suitable control material can be used.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

Calculation

COBAS INTEGRA analyzers automatically calculate the analyte concentration of each sample. For more details, please refer to Data Analysis in the Online Help (COBAS INTEGRA 400 plus/800 analyzers).

Conversion factor:⁷ mmol/L × 4.61 = mg/dL

Limitations - interference

Do not use volatile solvents in the work area when performing assays. Do not perform sample preparation (especially spiking of pools) in the immediate work area. Vapor contamination of reagents can impact calibration stability.

Criterion: Recovery within ± 10 % of initial value.

Serum/plasma

Icterus:⁸ No significant interference up to an I index of 60 for conjugated and unconjugated bilirubin (approximate conjugated and unconjugated bilirubin concentration: 1026 µmol/L or 60 mg/dL).

Hemolysis:⁸ No significant interference up to an H index of 1000 (approximate hemoglobin concentration: 621 µmol/L or 1000 mg/dL).

Lipemia (Intralipid):⁸ No significant interference up to an L index of 1200. There is poor correlation between the L index (corresponds to turbidity) and triglycerides concentration.

Drugs: No interference was found at therapeutic concentrations using common drug panels.^{9,10}

In very rare cases, gammopathy, in particular type IgM (Waldenström's macroglobulinemia), may cause unreliable results.¹¹

LDH/lactic acid (using a dose-response curve with purified LDH fractions added to 30 mmol/L lactic acid solution): No significant interference up to 2000 U/L LDH.

Urine

Drugs: No interference was found at therapeutic concentrations using common drug panels.¹⁰ Exception: Hypochlorite causes interference.

Urines containing sugars and contaminated with microorganisms may yield a false positive result due to fermentation of sugar to alcohol.

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

ACTION REQUIRED

Special Wash Programming: The use of special wash steps is mandatory when certain test combinations are run together on COBAS INTEGRA analyzers. Refer to the CLEAN Method Sheet for further instructions and for the latest version of the Extra wash cycle list.

Where required, special wash/carry-over evasion programming must be implemented prior to reporting results with this test.

Limits and ranges

Measuring range

2.17-108 mmol/L (10.0-498 mg/dL)

Determine samples having higher concentrations via the rerun function. Dilution of samples via the rerun function is a 1:3 dilution. Results from samples diluted using the rerun function are automatically multiplied by a factor of 3.

Lower limits of measurement

Lower detection limit of the test (serum, plasma, and urine): 2.17 mmol/L (10.0 mg/dL)

The lower detection limit represents the lowest measurable analyte level that can be distinguished from zero. It is calculated as the value lying 3 standard deviations above that of a zero sample (zero sample + 3 SD, repeatability, n = 21).

Expected values⁷

10.9-21.7 mmol/L (50.2-100 mg/dL)	Flushing, slowing of reflexes, impaired visual acuity
> 21.7 mmol/L (> 100 mg/dL)	Depression of CNS
> 86.8 mmol/L (> 400 mg/dL)	Fatalities reported

The legal definition of intoxication varies according to local law. Each laboratory should establish an acceptable reporting format and identify procedures for the reporting of abnormal results. Clinical consideration and

ETOH2

Ethanol Gen.2

professional judgment should be applied to the interpretation of any alcohol test results.

Specific performance data

Representative performance data on the COBAS INTEGRA analyzers are given below. Results obtained in individual laboratories may differ.

Precision

Precision was determined using human samples and controls in an internal protocol with repeatability (n = 21) and intermediate precision (1 aliquot per run, 1 run per day, 21 days). The following results were obtained:

Serum/plasma

Sample	Repeatability		Intermediate precision	
	Mean mmol/L (mg/dL)	CV %	Mean mmol/L (mg/dL)	CV %
Level 1	20.1 (93.0)	1.2	21.8 (100)	2.4
Level 2	42.0 (194)	1.1	42.8 (197)	3.9

Urine

Sample	Repeatability		Intermediate precision	
	Mean mmol/L (mg/dL)	CV %	Mean mmol/L (mg/dL)	CV %
Level 1	20.1 (93.0)	1.2	24.0 (111)	3.6
Level 2	31.9 (147)	1.7	30.7 (142)	3.3

Method comparison

Serum/plasma

Ethanol values for human serum and plasma samples obtained on a COBAS INTEGRA 700 analyzer using the COBAS INTEGRA Ethanol Gen.2 reagent (y) were compared with those determined using the corresponding reagent on a Roche/Hitachi 917 analyzer (x) and with those determined using the previous reagent (ETOH) on a COBAS INTEGRA 700 analyzer (x).

Roche/Hitachi 917 analyzer	Sample size (n) = 52
Passing/Bablok ¹²	Linear regression
$y = 0.958x + 0.242$ mmol/L	$y = 0.964x + 0.053$ mmol/L
$r = 0.970$	$r = 0.999$
SD (md 95) = 2.40	Sy.x = 1.06

The values were between 8.51 and 105 mmol/L (39.2 and 484 mg/dL).

COBAS INTEGRA 700 analyzer	Sample size (n) = 51
Passing/Bablok ¹²	Linear regression
$y = 0.957x - 0.474$ mmol/L	$y = 0.963x - 0.675$ mmol/L
$r = 0.969$	$r = 0.999$
SD (md 95) = 1.81	Sy.x = 0.818

The values were between 8.63 and 109 mmol/L (39.8 and 502 mg/dL).

Urine

Ethanol values for human urine samples obtained on a COBAS INTEGRA 700 analyzer using the COBAS INTEGRA Ethanol Gen.2 reagent (y) were compared with those determined using the corresponding reagent on a Roche/Hitachi 917 analyzer (x) and with those determined using the previous reagent (ETOH) on a COBAS INTEGRA 700 analyzer (x).

Roche/Hitachi 917 analyzer	Sample size (n) = 60
Passing/Bablok ¹²	Linear regression
$y = 0.964x - 0.217$ mmol/L	$y = 0.967x - 0.296$ mmol/L
$r = 0.978$	$r = 0.999$
SD (md 95) = 0.936	Sy.x = 0.779

The values were between 0.270 and 111 mmol/L (1.24 and 510 mg/dL).

COBAS INTEGRA 700 analyzer	Sample size (n) = 58
Passing/Bablok ¹²	Linear regression
$y = 0.997x - 0.235$ mmol/L	$y = 0.993x - 0.245$ mmol/L
$r = 0.979$	$r = 0.999$
SD (md 95) = 1.74	Sy.x = 0.699

The values were between 0.270 and 108 mmol/L (1.24 and 498 mg/dL).

Analytical specificity

COBAS INTEGRA Ethanol Gen.2 reagent is specific for ethyl alcohol. The following cross reactants were measured at 2000 mg/dL.

Compound	% Cross-reactivity (urine)
Acetaldehyde	-1.6
Acetone	0.0
N-butanol	0.1
Ethylene glycol	0.1
Isopropanol	0.3
Methanol	0.0
N-propanol	6.0

$$\text{Cross-reactivity (\%)} = \frac{100 \times (\text{analytical result} - \text{analyte concentration})}{\text{concentration of interferent}}$$

References

- Kaplan LA, Pesce AJ. Clinical Chemistry, Theory, Analysis and Correlation. Ladig D, Kasper R (ed), St Louis, CV Mosby Co 1984;1332-1334.
- Bucher T, Redetzki H. Specific photometric determination of ethyl alcohol based on an enzymatic reaction. Klin Wschr 1951;29:615-616.
- Proposed guidelines NCCLS: Blood Alcohol Testing in the Clinical Laboratory. NCCLS Vol. 8 No. 10. December 1988.
- Tietz NW. Fundamentals of Clinical Chemistry. 3rd ed. Philadelphia: WB Saunders 1987;890.
- Tietz NW, ed. Clinical Guide to Laboratory Tests, 4th ed. Philadelphia, PA: WB Saunders Company 2006;1344-1346.
- Levine B, Smith ML. Stability of drugs of abuse in biological specimens. Forensic Sci Rev 1990;2:148-156.
- Tietz NW. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia, PA: WB Saunders 1995;224-225.
- Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-475.
- Breuer J. Report on the Symposium "Drug effects in Clinical Chemistry Methods". Eur J Clin Chem Clin Biochem 1996;34:385-386.
- Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385.
- Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med 2007;45(9):1240-1243.
- Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard (for USA: see <https://usdiagnostics.roche.com> for definition of symbols used):

ETOH2

Ethanol Gen.2

cobas[®]
Drug abuse testing

CONTENT



GTIN

Contents of kit

Volume after reconstitution or mixing

Global Trade Item Number

FOR US CUSTOMERS ONLY: LIMITED WARRANTY

Roche Diagnostics warrants that this product will meet the specifications stated in the labeling when used in accordance with such labeling and will be free from defects in material and workmanship until the expiration date printed on the label. THIS LIMITED WARRANTY IS IN LIEU OF ANY OTHER WARRANTY, EXPRESS OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. IN NO EVENT SHALL ROCHE DIAGNOSTICS BE LIABLE FOR INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES.

COBAS, COBAS C and COBAS INTEGRA are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2016, Roche Diagnostics



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim
www.roche.com



Distribution in USA by:
Roche Diagnostics, Indianapolis, IN
US Customer Technical Support 1-800-428-2336