



08061505001V1.0

ISE Deproteinizer

cobas[®]

Order information

REF	CONTENT	Analyzer(s) on which kit(s) can be used
20763071 122	ISE Deproteinizer (6 × 21 mL)	COBAS INTEGRA 400 plus COBAS INTEGRA 800
04838181 001	ISE Deproteinizer (2 × 11 mL)	cobas c 111

English

Intended use

ISE Deproteinizer is a cleaning solution intended for use with the COBAS INTEGRA and **cobas c 111** ISE modules for cleaning the ion-selective electrodes.

Summary

ISE Deproteinizer is a cleaning solution used on the ISE module for the cleaning of the ion-selective electrodes, the mixing tower, and tubing during ISE maintenance. It is also used in the cleaning of the sample probe(s).

Reagents - working solutions

Sodium hypochlorite (NaOCl): approximately 1.2 %

Precautions and warnings

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents.

Disposal of all waste material should be in accordance with local guidelines. Safety data sheet available for professional user on request.

For USA: For prescription use only.

This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:



Warning

H315 Causes skin irritation.

H319 Causes serious eye irritation.

Prevention:

P264 Wash skin thoroughly after handling.

P280 Wear protective gloves/ eye protection/ face protection.

Response:

P302 + P352 IF ON SKIN: Wash with plenty of water.

P332 + P313 If skin irritation occurs: Get medical advice/attention.

P337 + P313 If eye irritation persists: Get medical advice/attention.

P362 + P364 Take off contaminated clothing and wash it before reuse.

Product safety labeling primarily follows EU GHS guidance.

Contact phone: all countries: +49-621-7590, USA: 1-800-428-2336

Reagent handling

Ready for use

COBAS INTEGRA 400 plus/800 analyzers:

Place the ISE Deproteinizer on the ISE rack.

cobas c 111 analyzer:

Place the aliquoted ISE Deproteinizer in the sample area as indicated on the touch screen of the analyzer.

Storage and stability

COBAS INTEGRA 400 plus/800 analyzer:

Shelf life at 2-8 °C: See expiration date on the label

On-board in use and refrigerated on the analyzer: 4 weeks

Replace the bottle when alerted to do so by the system. Discard any remaining solution.

cobas c 111 analyzer:

Shelf life at 2-8 °C: See expiration date on the label

after opening: 8 weeks at 2-8 °C or until the printed expiration date, whichever comes first

Materials provided

See "Reagents – working solutions" section for reagents.

Materials required (but not provided)

Refer to the appropriate instructions for use and to the operator's manual for additional required materials.

General laboratory equipment

Assay

For further information, please refer to the appropriate operator's manual for the analyzer concerned, the respective application sheets and method sheets of all necessary components.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard:

CONTENT	Contents of kit
REAGENT	Reagent
→	Volume after reconstitution or mixing
GTIN	Global Trade Item Number

FOR US CUSTOMERS ONLY: LIMITED WARRANTY

Roche Diagnostics warrants that this product will meet the specifications stated in the labeling when used in accordance with such labeling and will be free from defects in material and workmanship until the expiration date printed on the label. THIS LIMITED WARRANTY IS IN LIEU OF ANY OTHER WARRANTY, EXPRESS OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. IN NO EVENT SHALL ROCHE DIAGNOSTICS BE LIABLE FOR INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES.

COBAS, COBAS C and COBAS INTEGRA are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2016, Roche Diagnostics



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim
www.roche.com

Distribution in USA by:

Roche Diagnostics, Indianapolis, IN
US Customer Technical Support 1-800-428-2336

