

# ISE Etcher

## Order information

REF	CONTENT	Analyzer(s) on which kit(s) can be used
20763098 122	ISE Etcher (6 × 11 mL)	COBAS INTEGRA 400 plus COBAS INTEGRA 800
04838165 001	ISE Etcher (2 × 11 mL)	cobas c 111

## English

### Intended use

The ISE Etcher is a cleaning solution intended for use with the COBAS INTEGRA and **cobas c** 111 ISE modules for cleaning the ion-selective sodium electrode.

### Summary

The ISE Etcher is a cleaning solution used with the ISE module for cleaning the sodium electrode during ISE maintenance.

### Reagents - working solutions

Ammonium hydrogen difluoride  $[(\text{NH}_4)\text{HF}_2]$ : 100 mmol/L

### Precautions and warnings

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents.

Disposal of all waste material should be in accordance with local guidelines. Safety data sheet available for professional user on request.

For USA: For prescription use only.

This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:



### Warning

H315 Causes skin irritation.

H319 Causes serious eye irritation.

### Prevention:

P264 Wash skin thoroughly after handling.

P280 Wear protective gloves/ eye protection/ face protection.

### Response:

P302 + P352 IF ON SKIN: Wash with plenty of water.

P332 + P313 If skin irritation occurs: Get medical advice/attention.

P337 + P313 If eye irritation persists: Get medical advice/attention.

P362 + P364 Take off contaminated clothing and wash it before reuse.

Product safety labeling primarily follows EU GHS guidance.

Contact phone: all countries: +49-621-7590, USA: 1-800-428-2336

### Reagent handling

Ready for use

### COBAS INTEGRA 400 plus/800 analyzer:

Place the ISE Etcher on the ISE rack.

### cobas c 111 analyzer:

Place the ISE Etcher on the sample area as indicated on the touch screen of the analyzer.

### Storage and stability

#### COBAS INTEGRA 400 plus/800 analyzer:

Shelf life at 15-25 °C: See expiration date on the label

On-board in use and refrigerated on the analyzer: 8 weeks

Replace the bottle on alert from the system. Discard any remaining solution.

### cobas c 111 analyzer:

Shelf life at 2-8 °C: See expiration date on the label

On-board in use and refrigerated on the analyzer: 8 weeks

### Materials provided

See "Reagents – working solutions" section for reagents.

### Materials required (but not provided)

Refer to the appropriate instructions for use and to the operator's manual for additional required materials.

General laboratory equipment

### Assay

For further information, please refer to the appropriate operator's manual for the analyzer concerned, the respective application sheets and method sheets of all necessary components.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

### Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard:

CONTENT	Contents of kit
REAGENT	Reagent
→	Volume after reconstitution or mixing
GTIN	Global Trade Item Number

### FOR US CUSTOMERS ONLY: LIMITED WARRANTY

Roche Diagnostics warrants that this product will meet the specifications stated in the labeling when used in accordance with such labeling and will be free from defects in material and workmanship until the expiration date printed on the label. THIS LIMITED WARRANTY IS IN LIEU OF ANY OTHER WARRANTY, EXPRESS OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. IN NO EVENT SHALL ROCHE DIAGNOSTICS BE LIABLE FOR INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES.

COBAS, COBAS C and COBAS INTEGRA are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2016, Roche Diagnostics



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim  
www.roche.com

Distribution in USA by:  
Roche Diagnostics, Indianapolis, IN  
US Customer Technical Support 1-800-428-2336

