

# LDL\_C

LDL-Cholesterol plus 2nd generation

**cobas**<sup>®</sup>

Substrates

**Order information**

REF	CONTENT	Analyzer(s) on which <b>cobas c</b> pack(s) can be used
03038866 322	LDL-Cholesterol plus 2nd generation (175 tests)	System-ID 07 6627 5 COBAS INTEGRA 400 plus COBAS INTEGRA 800
12172623 122	Calibrator f.a.s. Lipids (3 x 1 mL)	System-ID 07 6570 8
12172623 160	Calibrator f.a.s. Lipids (3 x 1 mL, for USA)	System-ID 07 6570 8
10781827 122	Precinorm L (4 x 3 mL)	System-ID 07 9026 5
11778552 122	Precipath HDL/LDL-C (4 x 3 mL)	System-ID 07 9028 1
05117003 190	PreciControl ClinChem Multi 1 (20 x 5 mL)	System-ID 07 7469 3
05947626 190	PreciControl ClinChem Multi 1 (4 x 5 mL)	System-ID 07 7469 3
05947626 160	PreciControl ClinChem Multi 1 (4 x 5 mL, for USA)	System-ID 07 7469 3
05117216 190	PreciControl ClinChem Multi 2 (20 x 5 mL)	System-ID 07 7470 7
05947774 190	PreciControl ClinChem Multi 2 (4 x 5 mL)	System-ID 07 7470 7
05947774 160	PreciControl ClinChem Multi 2 (4 x 5 mL, for USA)	System-ID 07 7470 7
20756350 322	NaCl Diluent 9 % (6 x 22 mL)	System-ID 07 5635 0

**English****System information**

Test LDL\_C, test ID 0-301 on COBAS INTEGRA 400 plus analyzers;  
test ID 0-300 on COBAS INTEGRA 800 analyzers.

**Intended use**

In vitro test for the quantitative determination of LDL-cholesterol concentration in human serum and plasma on COBAS INTEGRA systems.

**Summary**

Low Density Lipoproteins (LDL) play a key role in causing and influencing the progression of atherosclerosis and, in particular, coronary sclerosis. The LDLs are derived from VLDLs (Very Low Density Lipoproteins) rich in triglycerides by the action of various lipolytic enzymes and are synthesized in the liver. The elimination of LDL from plasma takes place mainly by liver parenchymal cells via specific LDL receptors. Elevated LDL concentrations in blood and an increase in their residence time coupled with an increase in the biological modification rate results in the destruction of the endothelial function and a higher LDL-cholesterol uptake in the monocyte/macrophage system as well as by smooth muscle cells in vessel walls. The majority of cholesterol stored in atherosclerotic plaques originates from LDL.

The LDL-cholesterol value is the most powerful clinical predictor among all of the single parameters with respect to coronary atherosclerosis. Therefore, therapies focusing on lipid reduction primarily target the reduction of LDL-cholesterol which is then expressed in an improvement of the endothelial function, prevention of atherosclerosis and reducing its progression as well as preventing plaque rupture.

Various methods are available for the determination of LDL-cholesterol such as ultracentrifugation as the reference method, lipoprotein electrophoresis and precipitation methods. In the precipitation methods apolipoprotein B-containing LDL-cholesterol is, for example, precipitated using either polyvinyl sulfate, dextran sulfate or polycyclic anions. The LDL-cholesterol content is usually calculated from the difference between total cholesterol and cholesterol in the remainder (VLDL and HDL-cholesterol) in the supernatant after precipitation with polyvinyl sulfate and dextran sulfate. Lipid Research Clinics recommend a combination of ultracentrifugation and precipitation methods using polyanions in the presence of divalent cations. The precipitation methods are time-consuming, cannot be automated and are susceptible to interference by hyperlipidemic serum, particularly at high concentrations of free fatty acids. A more recent method is based on the determination of LDL-cholesterol after the sample is subjected to immunoadsorption and centrifugation.

The calculation of the LDL-cholesterol concentration according to Friedewald's formula is commonly practised. The formula is based on two cholesterol determinations, one triglyceride determination as well as precipitation of the HDL particles and presumes that a direct relationship exists between VLDL-cholesterol and triglycerides in fasting blood samples. Even in the presence of small amounts of chylomicrons or abnormal lipoproteins, the formula gives rise to falsely low LDL-cholesterol values. For this reason there is a great need for a simple and reliable method for

the determination of LDL-cholesterol without any preparatory steps or calculation.

This automated method for the direct determination of LDL-cholesterol takes advantage of the selective micellar solubilization of LDL-cholesterol by a nonionic detergent and the interaction of a sugar compound and lipoproteins (VLDL and chylomicrons). When a detergent is included in the enzymatic method for cholesterol determination (cholesterol esterase cholesterol oxidase coupling reaction), the relative reactivities of cholesterol in the lipoprotein fractions increase in this order:


HDL < chylomicrons < VLDL < LDL. In the presence of Mg<sup>++</sup>, a sugar compound markedly reduces the enzymatic reaction of the cholesterol measurement in VLDL and chylomicrons. The combination of a sugar compound with detergent enables the selective determination of LDL-cholesterol in serum.<sup>1,2,3,4,5,6,7,8</sup>

Non-fasting sample results are slightly lower than fasting results. Comparable non-fasting results were observed with the beta quantification method.<sup>9,10</sup>


This direct assay meets the 1995 NCEP goals of < 4 % total CV, bias ≤ 4 % versus reference method, and ≤ 12 % total analytical error.<sup>11</sup>

**Test principle**


Homogeneous enzymatic colorimetric assay.



Cholesterol esters are broken down quantitatively into free cholesterol and fatty acids by cholesterol esterase.



In the presence of oxygen, cholesterol is oxidized by cholesterol oxidase to Δ<sup>4</sup>-cholestenone and hydrogen peroxide.



a) Sodium N-(2-hydroxy-3-sulfoethyl)-3,5-dimethoxyaniline

In the presence of peroxidase, the hydrogen peroxide generated reacts with 4-aminoantipyrine and HSDA to form a purple-blue dye. The color intensity of this dye is directly proportional to the cholesterol concentration and is measured photometrically.

# LDL\_C

LDL-Cholesterol plus 2nd generation


  
Substrates

## Reagents - working solutions

- R1** MOPS<sup>b)</sup>: 20.1 mmol/L, pH 6.5; HSDA: 0.958 mmol/L; ascorbate oxidase (recombinant):  $\geq 50 \mu\text{kat/L}$ ; peroxidase (horseradish):  $\geq 167 \mu\text{kat/L}$ ; stabilizers, preservative
- SR** MOPS: 20.1 mmol/L, pH 6.8;  $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$ : 8.11 mmol/L; 4-aminoantipyrine: 2.46 mmol/L; cholesterol esterase (microbial):  $\geq 50 \mu\text{kat/L}$ ; cholesterol oxidase (microbial):  $\geq 33 \mu\text{kat/L}$ ; peroxidase (horseradish):  $\geq 334 \mu\text{kat/L}$ ; stabilizers; preservative

b) 3-morpholinopropane sulfonic acid buffer

R1 is in position B and SR is in position C.

## Precautions and warnings

Pay attention to all precautions and warnings listed in Section 1 / Introduction of this Method Manual.

For USA: For prescription use only.

## Reagent handling

Ready for use

## Storage and stability

Shelf life at 2-8 °C See expiration date on  
**cobas c** pack label

## COBAS INTEGRA 400 plus system

On-board in use at 10-15 °C 12 weeks

## COBAS INTEGRA 800 system

On-board in use at 8 °C 12 weeks

## Specimen collection and preparation

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable.

Serum

Plasma: Heparin (Li-, Na-,  $\text{NH}_4^+$ -) or  $\text{K}_3\text{-EDTA}$  plasma

EDTA plasma causes decreased values.

Fasting and non-fasting samples can be used.<sup>10</sup>

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

Stability:<sup>12</sup> 7 days at 2-8 °C  
30 days at -70 °C

It is reported that EDTA stabilizes lipoproteins.<sup>11</sup>

## Materials provided

See "Reagents – working solutions" section for reagents.

## Materials required (but not provided)

NaCl Diluent 9 %, Cat. No. 20756350322, system-ID 07 5635 0 for automatic sample postdilution. NaCl Diluent 9 % is placed in its predefined rack position and is stable for 4 weeks on-board

COBAS INTEGRA 400 plus/800 analyzers.

## Assay

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

## Application for serum and plasma

## COBAS INTEGRA 400 plus test definition

Measuring mode Absorbance  
Abs. calculation mode Endpoint

Reaction mode R1-S-SR  
Reaction direction Increase  
Wavelength A/B 583/659 nm  
Calc. first/last 33/69  
Unit mmol/L

## Pipetting parameters

R1 150  $\mu\text{L}$   
Sample 2  $\mu\text{L}$  7  $\mu\text{L}$   
SR 50  $\mu\text{L}$   
Total volume 209  $\mu\text{L}$

Diluent ( $\text{H}_2\text{O}$ )

## COBAS INTEGRA 800 test definition

Measuring mode Absorbance  
Abs. calculation mode Endpoint  
Reaction mode R1-S-SR  
Reaction direction Increase  
Wavelength A/B 583/659 nm  
Calc. first/last 44/98  
Unit mmol/L

## Pipetting parameters

R1 150  $\mu\text{L}$   
Sample 2  $\mu\text{L}$  7  $\mu\text{L}$   
SR 50  $\mu\text{L}$   
Total volume 209  $\mu\text{L}$

Diluent ( $\text{H}_2\text{O}$ )

## Calibration

Calibrator C.f.a.s. Lipids  
Use deionized water as zero calibrator.  
Calibration mode Linear regression  
Calibration replicate Duplicate recommended  
Calibration interval Each lot

Traceability: This method has been standardized against the beta quantification method as defined in the recommendations in the LDL Cholesterol Method Certification Protocol for Manufacturers.<sup>13</sup>

## Quality control

Quality control Precinorm L or PreciControl ClinChem Multi 1  
Precipath HDL/LDL-C or PreciControl ClinChem Multi 2  
Control interval 24 hours recommended  
Control sequence User defined  
Control after Recommended  
calibration

For quality control, use control materials as listed in the "Order information" section.

In addition, other suitable control material can be used.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

# LDL\_C

## LDL-Cholesterol plus 2nd generation

**cobas®**
**Substrates**

Follow the applicable government regulations and local guidelines for quality control.

### Calculation

COBAS INTEGRA analyzers automatically calculate the analyte concentration of each sample. For more details, please refer to Data Analysis in the Online Help (COBAS INTEGRA 400 plus/800 analyzers).

Conversion factor: mmol/L  $\times$  38.66 = mg/dL

### Limitations - interference

Criterion: Recovery within  $\pm$  10 % of initial value.

#### Serum, plasma

Icterus:<sup>14</sup> No significant interference up to an I index of 40 for conjugated and unconjugated bilirubin (approximate conjugated and unconjugated bilirubin concentration: 684  $\mu$ mol/L or 40 mg/dL).

Hemolysis:<sup>14</sup> No significant interference up to an H index of 1000 (approximate hemoglobin concentration: 0.62 mmol/L or 1000 mg/dL).

Lipemia (Intralipid):<sup>14</sup> No significant interference up to an L index of 400. No significant interference from native triglycerides up to 1200 mg/dL. There is poor correlation between the L index (corresponds to turbidity) and triglycerides concentration.

Drugs: No interference was found at therapeutic concentrations using common drug panels.<sup>15,16</sup>

Acetaminophen intoxications are frequently treated with N-Acetylcysteine. N-Acetylcysteine at the therapeutic concentration when used as an antidote and the Acetaminophen metabolite N-acetyl-p-benzoquinone imine (NAPQI) independently may cause falsely low results.

Venipuncture should be performed prior to the administration of Metamizole. Venipuncture immediately after or during the administration of Metamizole may lead to falsely low results.

Abnormal liver function affects lipid metabolism; consequently, HDL and LDL results are of limited diagnostic value. In some patients with abnormal liver function, the LDL\_C result is significantly negatively biased versus the reference method (beta quantification) result.

In very rare cases, gammopathy, in particular type IgM (Waldenström's macroglobulinemia), may cause unreliable results.<sup>17</sup>

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

### ACTION REQUIRED

**Special Wash Programming:** The use of special wash steps is mandatory when certain test combinations are run together on COBAS INTEGRA analyzers. Refer to the CLEAN Method Sheet for further instructions and for the latest version of the Extra wash cycle list.

**Where required, special wash/carry-over evasion programming must be implemented prior to reporting results with this test.**

### Limits and ranges

#### Measuring range

0.10-14.2 mmol/L (3.87-550 mg/dL)

Determine samples having higher concentrations via the rerun function. Dilution of samples via the rerun function is a 1:4 dilution. Results from samples diluted using the rerun function are automatically multiplied by a factor of 4.

The use of a dilution factor  $<$  4 is not allowed.

#### Lower limits of measurement

Lower detection limit of the test:

0.10 mmol/L (3.87 mg/dL)

The lower detection limit represents the lowest measurable analyte level that can be distinguished from zero. It is calculated as the value lying 3 standard deviations above that of a zero sample (zero sample + 3 SD, repeatability,  $n = 21$ ).

#### Expected values<sup>18</sup>

Levels in terms of risk for coronary heart disease.

Adult levels:

Optimal	$<$ 2.59 mmol/L	( $<$ 100 mg/dL)
Near optimal/above optimal	2.59-3.34 mmol/L	(100-129 mg/dL)
Borderline high	3.37-4.12 mmol/L	(130-159 mg/dL)

High 4.14-4.89 mmol/L (160-189 mg/dL)

Very high  $\geq$  4.92 mmol/L ( $\geq$  190 mg/dL)

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

### Specific performance data

Representative performance data on the analyzers are given below. Results obtained in individual laboratories may differ.

#### Precision

Precision was determined using human samples and controls in an internal protocol with repeatability ( $n = 21$ ) and intermediate precision (3 aliquots per run, 1 run per day, 21 days). The following results were obtained:

Repeatability	Level 1	Level 2
Mean	0.62 mmol/L (24.0 mg/dL)	4.29 mmol/L (166 mg/dL)
CV	1.5 %	1.1 %

Intermediate precision	Level 1	Level 2
Mean	2.03 mmol/L (78.5 mg/dL)	5.57 mmol/L (215 mg/dL)
CV	1.9 %	1.8 %

### Method comparison

LDL-cholesterol values for human serum samples obtained on a COBAS INTEGRA 700 analyzer using the COBAS INTEGRA LDL-Cholesterol plus 2nd generation reagent (y) were compared to those determined using the corresponding reagent on Roche/Hitachi 917 analyzer (x).

#### Roche/Hitachi 917 analyzer

Sample size (n)	55
Corr. coefficient (r)	0.996
Lin. regression	$y = 0.96x + 0.11$ mmol/L
Passing/Bablok <sup>19</sup>	$y = 1.01x + 0.01$ mmol/L

The sample concentrations were between 0.72 to 7.02 mmol/L (27.7 to 270 mg/dL).

### References

- Rifai N, Warnick GR, McNamara JR, et al. Measurement of Low-Density-Lipoprotein Cholesterol in Serum: a Status Report. Clin Chem 1992;38:150-160.
- Bachorik PS. Measurement of Low-Density Lipoprotein Cholesterol. In: Rifai N, Warnick GR, Dominiczak MH, eds. AACC Press 2000;12:245-263.
- National Cholesterol Education Program. Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel II). NIH Publication No. 93-3095 1995.
- Naito HK, Strong JP, Scott MG, et al. Atherogenesis: current topics on etiology and risk factors. Clin Chem 1995;41:132-133.
- Wieland H, Seidel D. Quantitative Lipoprotein Electrophoresis. In: Handbook of Electrophoresis, Vol III, ed. Lewis A, Boca Raton: CRC Press 1983;83-102.
- Armstrong V, Seidel D. Evaluation of a Commercial Kit for the Determination of LDL-Cholesterol in Serum Based on Precipitation of LDL with Dextran Sulfate. Arztl Lab 1985;31:325-330.
- Friedewald WF, Levy RI, Frederickson DS. Estimation of the Concentration of Low-Density Lipoprotein Cholesterol in Plasma, Without Use of the Preparative Ultracentrifuge. Clin Chem. 1972;18(6):499-502.
- Bachorik PS, Ross JW. National cholesterol education program recommendations for measurement of low-density lipoprotein cholesterol: executive summary. Clin Chem 1995;41:1414-1420.

# LDL\_C

## LDL-Cholesterol plus 2nd generation




**cobas®**  
Substrates

- 9 Cohn JS, McNamara JR, Schaefer EJ. Lipoprotein Cholesterol Concentrations in the Plasma of Human Subjects as Measured in the Fed and Fasted States. Clin Chem 1988;34:2456-2459.
- 10 Pisani T, Gebiski CP, Leary ET, et al. Accurate Direct Determination of Low-density Lipoprotein Cholesterol Using an Immunoseparation Reagent and Enzymatic Cholesterol Assay. Arch Pathol Lab Med 1995 Dec;119(12):1127-1135.
- 11 Cooper GR, Myers GL, Smith SJ, et al. Standardization of Lipid, Lipoprotein, and Apolipoprotein Measurements. Clin Chem 1988;34(8B):B95-B105.
- 12 Data on file at Roche Diagnostics.
- 13 LDL Cholesterol Method Certification Protocol for Manufacturers. National Reference System for Cholesterol. Cholesterol Reference Method Laboratory Network 1997, October.
- 14 Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-475.
- 15 Breuer J. Report on the Symposium "Drug Effects in Clinical Chemistry Methods". Eur J Clin Chem Clin Biochem 1996;34:385-386.
- 16 Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385.
- 17 Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med 2007;45(9):1240-1243.
- 18 Third Report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). NIH Publication No 01-3670; May 2001.
- 19 Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

### Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard.

	Contents of kit
	Volume after reconstitution or mixing
	Global Trade Item Number

### FOR US CUSTOMERS ONLY: LIMITED WARRANTY

Roche Diagnostics warrants that this product will meet the specifications stated in the labeling when used in accordance with such labeling and will be free from defects in material and workmanship until the expiration date printed on the label. THIS LIMITED WARRANTY IS IN LIEU OF ANY OTHER WARRANTY, EXPRESS OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. IN NO EVENT SHALL ROCHE DIAGNOSTICS BE LIABLE FOR INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES.

COBAS, COBAS INTEGRA, COBAS C, PRECINORM, PRECIPATH and PRECICONTROL are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2016, Roche Diagnostics



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim  
www.roche.com

Distribution in USA by:  
Roche Diagnostics, Indianapolis, IN  
US Customer Technical Support 1-800-428-2336

