

Order information

| REF | CONTENT | Analyzer(s) on which cobas c pack(s) can be used |
|--------------|---|--|
| 03333752 190 | Alkaline Phosphatase acc. to IFCC Gen.2 ALP2S 200 tests | System-ID 07 6761 1 Roche/Hitachi cobas c 311, cobas c 501/502 |
| 03333701 190 | Alkaline Phosphatase acc. to IFCC Gen.2 ALP2L 400 tests | System-ID 07 6760 3 Roche/Hitachi cobas c 311, cobas c 501/502 |
| 10759350 190 | Calibrator f.a.s. (12 x 3 mL) | Code 401 |
| 10759350 360 | Calibrator f.a.s. (12 x 3 mL, for USA) | Code 401 |
| 12149435 122 | Precinorm U plus (10 x 3 mL) | Code 300 |
| 12149435 160 | Precinorm U plus (10 x 3 mL, for USA) | Code 300 |
| 12149443 122 | Precipath U plus (10 x 3 mL) | Code 301 |
| 12149443 160 | Precipath U plus (10 x 3 mL, for USA) | Code 301 |
| 10171743 122 | Precinorm U (20 x 5 mL) | Code 300 |
| 10171735 122 | Precinorm U (4 x 5 mL) | Code 300 |
| 10171778 122 | Precipath U (20 x 5 mL) | Code 301 |
| 10171760 122 | Precipath U (4 x 5 mL) | Code 301 |
| 05117003 190 | PreciControl ClinChem Multi 1 (20 x 5 mL) | Code 391 |
| 05947626 190 | PreciControl ClinChem Multi 1 (4 x 5 mL) | Code 391 |
| 05947626 160 | PreciControl ClinChem Multi 1 (4 x 5 mL, for USA) | Code 391 |
| 05117216 190 | PreciControl ClinChem Multi 2 (20 x 5 mL) | Code 392 |
| 05947774 190 | PreciControl ClinChem Multi 2 (4 x 5 mL) | Code 392 |
| 05947774 160 | PreciControl ClinChem Multi 2 (4 x 5 mL, for USA) | Code 392 |
| 04489357 190 | Diluent NaCl 9 % (50 mL) | System-ID 07 6869 3 |

English

System information

For **cobas c** 311/501 analyzers:

ALP2S: ACN 158

ALP2L: ACN 683

For **cobas c** 502 analyzer:

ALP2S: ACN 8158

ALP2L: ACN 8683

Intended use

In vitro test for the quantitative determination of alkaline phosphatase in human serum and plasma on Roche/Hitachi **cobas c** systems.

Summary^{1,2,3,4,5,6}

Alkaline phosphatase in serum consists of four structural genotypes: the liver-bone-kidney type, the intestinal type, the placental type and the variant from the germ cells. It occurs in osteoblasts, hepatocytes, leukocytes, the kidneys, spleen, placenta, prostate and the small intestine. The liver-bone-kidney type is particularly important.

A rise in the alkaline phosphatase occurs with all forms of cholestasis, particularly with obstructive jaundice. It is also elevated in diseases of the skeletal system, such as Paget's disease, hyperparathyroidism, rickets and osteomalacia, as well as with fractures and malignant tumors. A considerable rise in the alkaline phosphatase activity is sometimes seen in children and juveniles. It is caused by increased osteoblast activity following accelerated bone growth.

The assay method was first described by King and Armstrong, modified by Ohmori, Bessey, Lowry and Brock and later improved by Hausamen et al. In 1983 the International Federation of Clinical Chemistry (IFCC) recommended a standardized method for the determination of alkaline phosphatase using an optimized substrate concentration and 2-amino-2-methyl-1-propanol as buffer plus the cations magnesium and zinc. The assay described here is based on this recommendation, but was optimized for performance and stability. The assay was standardized against the IFCC reference formulation proposed above.

Test principle⁶

Colorimetric assay in accordance with a standardized method. In the presence of magnesium and zinc ions, p-nitrophenyl phosphate is cleaved by phosphatases into phosphate and p-nitrophenol.



The p-nitrophenol released is directly proportional to the catalytic ALP activity. It is determined by measuring the increase in absorbance.

Reagents - working solutions

R1 2-amino-2-methyl-1-propanol: 1.724 mol/L, pH 10.44 (30 °C);
magnesium acetate: 3.83 mmol/L; zinc sulfate: 0.766 mmol/L;
N-(2-hydroxyethyl)-ethylenediamine triacetic acid: 3.83 mmol/L

R2 p-nitrophenyl phosphate: 132.8 mmol/L, pH 8.50 (25 °C);
preservatives

R1 is in position B and R2 is in position C.

Precautions and warnings

For in vitro diagnostic use.
Exercise the normal precautions required for handling all laboratory reagents.
Disposal of all waste material should be in accordance with local guidelines.
Safety data sheet available for professional user on request.

For USA: For prescription use only.

This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:



Warning

- | | |
|------|--|
| H315 | Causes skin irritation. |
| H319 | Causes serious eye irritation. |
| H412 | Harmful to aquatic life with long lasting effects. |

Prevention:

- | | |
|------|--------------------------------------|
| P264 | Wash skin thoroughly after handling. |
| P273 | Avoid release to the environment. |

P280 Wear protective gloves/ eye protection/ face protection.

Response:

P332 + P313 If skin irritation occurs: Get medical advice/attention.

P337 + P313 If eye irritation persists: Get medical advice/attention.

Disposal:

P501 Dispose of contents/container to an approved waste disposal plant.

Product safety labeling primarily follows EU GHS guidance.

Contact phone: all countries: +49-621-7590, USA: 1-800-428-2336

Reagent handling

Ready for use

Storage and stability

ALP2S, ALP2L

Shelf life at 2-8 °C: See expiration date on **cobas c** pack label.

On-board in use and refrigerated on the analyzer: 8 weeks

Diluent NaCl 9 %

Shelf life at 2-8 °C: See expiration date on **cobas c** pack label.

On-board in use and refrigerated on the analyzer: 12 weeks

Specimen collection and preparation

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable.
Serum.

Plasma: Li-heparin plasma.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

Stability:⁷ 7 days at 15-25 °C
7 days at 2-8 °C
2 months at (-15)-(-25) °C

Materials provided

See "Reagents – working solutions" section for reagents.

Materials required (but not provided)

- See "Order information" section
- General laboratory equipment

Assay

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

The performance of applications not validated by Roche is not warranted and must be defined by the user.

Applications for serum and plasma**cobas c 311 test definition**

Assay type Rate A
Reaction time / Assay points 10 / 13-31
Wavelength (sub/main) 480/450 nm

| | | | |
|--------------------|--------------|----------------------------|----------------|
| Reaction direction | Increase | | |
| Units | U/L (µkat/L) | | |
| Reagent pipetting | | Diluent (H ₂ O) | |
| R1 | 75 µL | 25 µL | |
| R2 | 17 µL | 21 µL | |
| Sample volumes | Sample | Sample dilution | |
| | | Sample | Diluent (NaCl) |
| Normal | 2.8 µL | – | – |
| Decreased | 2.8 µL | 20 µL | 80 µL |
| Increased | 2.8 µL | – | – |

cobas c 501 test definition

| | | | |
|------------------------------|--------------|----------------------------|----------------|
| Assay type | Rate A | | |
| Reaction time / Assay points | 10 / 19-48 | | |
| Wavelength (sub/main) | 480/450 nm | | |
| Reaction direction | Increase | | |
| Units | U/L (µkat/L) | | |
| Reagent pipetting | | Diluent (H ₂ O) | |
| R1 | 75 µL | 25 µL | |
| R2 | 17 µL | 21 µL | |
| Sample volumes | Sample | Sample dilution | |
| | | Sample | Diluent (NaCl) |
| Normal | 2.8 µL | – | – |
| Decreased | 2.8 µL | 20 µL | 80 µL |
| Increased | 2.8 µL | – | – |

cobas c 502 test definition

| | | | |
|------------------------------|--------------|----------------------------|----------------|
| Assay type | Rate A | | |
| Reaction time / Assay points | 10 / 19-48 | | |
| Wavelength (sub/main) | 480/450 nm | | |
| Reaction direction | Increase | | |
| Units | U/L (µkat/L) | | |
| Reagent pipetting | | Diluent (H ₂ O) | |
| R1 | 75 µL | 25 µL | |
| R2 | 17 µL | 21 µL | |
| Sample volumes | Sample | Sample dilution | |
| | | Sample | Diluent (NaCl) |
| Normal | 2.8 µL | – | – |
| Decreased | 2.8 µL | 20 µL | 80 µL |
| Increased | 5.6 µL | – | – |

Calibration

| | |
|-----------------------|--|
| Calibrators | S1: H ₂ O S2: C.f.a.s. |
| Calibration mode | Linear |
| Calibration frequency | 2-point calibration <ul style="list-style-type: none"> after reagent lot change as required following quality control procedures |

Traceability: This method has been standardized against the proposed IFCC formulation⁶ using calibrated pipettes together with a manual photometer providing absolute values and the substrate-specific absorptivity, ε.

Quality control

For quality control, use control materials as listed in the "Order information" section.

In addition, other suitable control material can be used.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

Calculation

Roche/Hitachi **cobas c** systems automatically calculate the analyte activity of each sample.

Conversion factor: U/L x 0.0167 = µkat/L

Limitations - interference

Criterion: Recovery within ± 10 % of initial value at an alkaline phosphatase activity of 100 U/L (1.67 µkat/L).

Icterus:⁸ No significant interference up to an I index of 60 for conjugated and unconjugated bilirubin (approximate conjugated and unconjugated bilirubin concentration: 1026 µmol/L or 60 mg/dL).

Hemolysis:⁸ No significant interference up to an H index of 200 (approximate hemoglobin concentration: 124 µmol/L or 200 mg/dL).

Lipemia (Intralipid):⁸ No significant interference up to an L index of 2000. There is poor correlation between the L index (corresponds to turbidity) and triglycerides concentration.

Drugs: No interference was found at therapeutic concentrations using common drug panels.^{9,10}

In very rare cases, gammopathy, in particular type IgM (Waldenström's macroglobulinemia), may cause unreliable results.¹¹

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

ACTION REQUIRED

Special Wash Programming: The use of special wash steps is mandatory when certain test combinations are run together on Roche/Hitachi **cobas c** systems. The latest version of the carry-over evasion list can be found with the NaOHD-SMS-SmpCin1+2-SCCS Method Sheets. For further instructions refer to the operator's manual. **cobas c** 502 analyzer: All special wash programming necessary for avoiding carry-over is available via the **cobas** link, manual input is not required.

Where required, special wash/carry-over evasion programming must be implemented prior to reporting results with this test.

Limits and ranges**Measuring range**

5-1200 U/L (0.084-20.0 µkat/L)

Determine samples having higher activities via the rerun function. Dilution of samples via the rerun function is a 1:5 dilution. Results from samples diluted using the rerun function are automatically multiplied by a factor of 5.

Lower limits of measurement**Lower detection limit of the test**

5 U/L (0.084 µkat/L)

The lower detection limit represents the lowest measurable analyte level that can be distinguished from zero. It is calculated as the value lying 3 standard deviations above that of the lowest standard (standard 1 + 3 SD, repeatability, n = 21).

Expected values

(measured at 37 °C)

Adults¹²

Males (n = 221) 40-129 U/L (0.67-2.15 µkat/L)

Females (n = 229) 35-104 U/L (0.58-1.74 µkat/L)

Consensus values¹³

Males 40-130 U/L (0.67-2.17 µkat/L)

Females 35-105 U/L (0.58-1.75 µkat/L)

Children¹⁴**Males****Age**

0 – 14 days 83-248 U/L (1.39-4.14 µkat/L)

15 days – < 1 year 122-469 U/L (2.04-7.83 µkat/L)

1 – < 10 years 142-335 U/L (2.37-5.59 µkat/L)

10 – < 13 years 129-417 U/L (2.15-6.96 µkat/L)

13 – < 15 years 116-468 U/L (1.94-7.82 µkat/L)

15 – < 17 years 82-331 U/L (1.37-5.53 µkat/L)

17 – < 19 years 55-149 U/L (0.92-2.49 µkat/L)

Females**Age**

0 – 14 days 83-248 U/L (1.39-4.14 µkat/L)

15 days – < 1 year 122-469 U/L (2.04-7.83 µkat/L)

1 – < 10 years 142-335 U/L (2.37-5.59 µkat/L)

10 – < 13 years 129-417 U/L (2.15-6.96 µkat/L)

13 – < 15 years 57-254 U/L (0.95-4.24 µkat/L)

15 – < 17 years 50-117 U/L (0.84-1.95 µkat/L)

17 – < 19 years 45-87 U/L (0.75-1.45 µkat/L)

Roche has not evaluated reference ranges in a pediatric population.

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

Specific performance data

Representative performance data on the analyzers are given below. Results obtained in individual laboratories may differ.

Precision

Precision was determined using human samples and controls in an internal protocol with repeatability (n = 21) and intermediate precision (3 aliquots per run, 1 run per day, 21 days). The following results were obtained:

| Repeatability | Mean | SD | CV |
|------------------------|--------------|--------------|-----|
| | U/L (µkat/L) | U/L (µkat/L) | % |
| Precinorm U | 99.2 (1.65) | 0.7 (0.01) | 0.7 |
| Precipath U | 241 (4.02) | 1 (0.02) | 0.6 |
| Human serum 1 | 54.6 (0.912) | 0.5 (0.008) | 0.9 |
| Human serum 2 | 648 (10.8) | 4 (0.1) | 0.7 |
| Intermediate precision | Mean | SD | CV |
| | U/L (µkat/L) | U/L (µkat/L) | % |
| Precinorm U | 92.8 (1.56) | 2.2 (0.04) | 2.4 |
| Precipath U | 224 (3.74) | 4 (0.06) | 1.7 |
| Human serum 3 | 82.2 (1.37) | 1.8 (0.03) | 2.1 |
| Human serum 4 | 1025 (17.1) | 9 (0.2) | 0.9 |

Method comparison

Alkaline phosphatase values for human serum and plasma samples obtained on a Roche/Hitachi **cobas c** 501 analyzer (y) were compared with those determined using the corresponding reagent on a Roche/Hitachi 917 analyzer (x).

Sample size (n) = 203

| Passing/Bablok ¹⁵ | Linear regression |
|------------------------------|------------------------|
| y = 0.988x + 1.31 U/L | y = 0.991x + 0.799 U/L |
| τ = 0.961 | r = 0.997 |

The sample activities were between 50.0 and 1002 U/L (0.835 and 16.7 µkat/L).

References

- Greiling H, Gressner AM, eds. Lehrbuch der Klinischen Chemie und Pathobiochemie, 3rd ed. Stuttgart/New York: Schattauer Verlag 1995.
- King EJ, Armstrong AR. Can Med Assoc J 1934;31:376
- Ohmori Y. Über die Phosphomomesterase. Enzymologia 1937;4:217-231.
- Bessey OA, Lowry OH, Brock MJ. A method for the rapid determination of alkaline phosphatase with five cubic millimeters of serum. J Biol Chem 1946;164:321-329.
- Hausamen TU, Helger R, Rick W, et al. Optimal conditions for the determination of serum alkaline phosphatase by a new kinetic method. Clin Chim Acta 1967;15:241-245.
- Tietz NW, Rinker AD, Shaw LM. J Clin Chem Clin Biochem 1983;21:731-748.
- Use of Anticoagulants in Diagnostic Laboratory Investigations. WHO Publication WHO/DIL/LAB/99.1 Rev.2.
- Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-475.
- Breuer J. Report on the Symposium "Drug effects in Clinical Chemistry Methods". Eur J Clin Chem Clin Biochem 1996;34:385-386.
- Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385.
- Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med 2007;45(9):1240-1243.
- Abicht K, El-Samalouti V, Junge W, et al. Multicenter evaluation of new GGT and ALP reagents with new reference standardization and determination of 37 °C reference intervals. Clin Chem Lab Med 2001;39:Special Supplement pp S 346.
- Thomas L, Müller M, Schumann G, et al. Consensus of DGKL and VDGH for interim reference intervals on enzymes in serum. J Lab Med 2005;29:301-308.
- Estey MP, Cohen AH, Colantonio DA, et al. CLSI-based transference of the CALIPER database of pediatric reference intervals from Abbott to Beckman, Ortho, Roche and Siemens Clinical Chemistry Assays: Direct validation using reference samples from the CALIPER cohort. Clin Biochem 2013;46:1197-1219.
- Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard:

CONTENT

Contents of kit



Volume after reconstitution or mixing

GTIN

Global Trade Item Number

FOR US CUSTOMERS ONLY: LIMITED WARRANTY

Roche Diagnostics warrants that this product will meet the specifications stated in the labeling when used in accordance with such labeling and will be free from defects in material and workmanship until the expiration date printed on the label. THIS LIMITED WARRANTY IS IN LIEU OF ANY OTHER WARRANTY, EXPRESS OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. IN NO EVENT SHALL ROCHE DIAGNOSTICS BE LIABLE FOR INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES.

COBAS, COBAS C, PRECINORM, PRECIPATH and PRECICONTROL are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2016, Roche Diagnostics



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim
www.roche.com

Distribution in USA by:
Roche Diagnostics, Indianapolis, IN
US Customer Technical Support 1-800-428-2336

