

AMPS2

ONLINE DAT Amphetamines II



Order information

REF	CONTENT	Analyzer(s) on which cobas c pack(s) can be used
04939425 190	ONLINE DAT Amphetamines II (200 tests)	System-ID 07 6980 0 Roche/Hitachi cobas c 501/502
03304671 190	Preciset DAT Plus I calibrator CAL 6	Code 436
07978766 190	Serum DAT Control Low (ACQ Partner Channel*)	
07978740 190	Serum DAT Control High (ACQ Partner Channel*)	

*Roche does not hold the product registration for Partner Channels. The legal manufacturer indicated on the kit is solely responsible for all of the design, legal, and regulatory aspects of the product.

English

System information

For **cobas c** 501 analyzer:

AMQ3S: ACN 746: for qualitative assay, 300 ng/mL

For **cobas c** 502 analyzer:

AMQ3S: ACN 8746: for qualitative assay, 300 ng/mL

Intended use

Amphetamines II (AMPS2) is an in vitro diagnostic test for the qualitative detection of amphetamines and methamphetamines in human serum and plasma on Roche/Hitachi **cobas c** systems. For serum/plasma the cutoff concentration is 300 ng/mL when calibrated with d-methamphetamine.

Amphetamines II provides only a preliminary test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) or Liquid Chromatography coupled with Tandem Mass Spectrometry (LC-MS/MS) is the preferred confirmatory method.¹ Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used.

Summary

The amphetamines are known as the sympathomimetic amines as they mimic the effects of stimulation of the sympathetic nervous system. These small molecules, based on β -phenylethylamine, structurally resemble the body's own catecholamines. A wide variety have been created via substitutions anywhere on the structure. The amphetamines are potent central nervous stimulants. As such they can increase wakefulness, physical activity, and decrease appetite. The amphetamines have some limited indications and approval for use in ADHD, narcolepsy, and obesity. However, because these CNS stimulants convey a sense of self-confidence, well being, and euphoria, they are highly addictive, widely abused, and consequently controlled substances.² Abuse can lead to medical, psychological, and social consequences. Adverse health effects include memory loss, aggression, psychotic behavior, heart damage, malnutrition, and severe dental problems.³ Amphetamine may be self-administered either orally or by intravenous injection in amounts of up to 2000 mg daily by tolerant addicts. It is a metabolite of a number of other drugs including methamphetamine. Normally about 30 % is excreted unchanged in the 24 hour urine, but this may change to as much as 74 % in acid urine and may decrease to 1 % in alkaline urine.⁴

Amphetamines II is calibrated with d-methamphetamine and therefore the sensitivity towards amphetamines is different than d-methamphetamine, as indicated in the "Analytical specificity" section.

Test principle

The assay is based on the kinetic interaction of microparticles in a solution (KIMS)^{5,6} as measured by changes in light transmission. In the absence of sample drug, soluble drug conjugates bind to antibody-bound microparticles, causing the formation of particle aggregates. As the aggregation reaction proceeds in the absence of sample drug, the absorbance increases.

When a serum sample contains the drug in question, this drug competes with the drug derivative conjugate for microparticle-bound antibody. Antibody bound to sample drug is no longer available to promote particle aggregation, and subsequent particle lattice formation is inhibited. The presence of sample drug diminishes the increasing absorbance in proportion to the concentration of drug in the sample. Sample drug content is determined relative to the value obtained for a known cutoff concentration of drug.⁷

Reagents - working solutions

- R1** Conjugated amphetamine and methamphetamine derivatives; buffer; bovine serum albumin; 0.09 % sodium azide
- R2** Microparticles attached to amphetamine and methamphetamine antibodies (mouse monoclonal); buffer; bovine serum albumin; 0.09 % sodium azide

R1 is in position B and R2 is in position C.

Precautions and warnings

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents.

Disposal of all waste material should be in accordance with local guidelines. Safety data sheet available for professional user on request.

Reagent handling

Ready for use

Carefully invert reagent container several times prior to use to ensure that the reagent components are mixed.

Storage and stability

Shelf life at 2-8 °C: See expiration date on **cobas c** pack label

On-board in use and refrigerated on the analyzer: 8 weeks

Do not freeze.

Specimen collection and preparation

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable.

Serum: Serum tubes with and without separating gel.

Plasma: K₂- or K₃-EDTA, lithium heparin.

Stability: 5 days capped at 15-25 °C
14 days capped at 2-8 °C
6 months capped at -20 °C

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

Specimens can be repeatedly frozen and thawed up to 3 times.

Invert thawed specimens several times prior to testing.

CAUTION: Specimen dilutions should only be used to interpret results of Calc.? and Samp.? alarms, or when estimating concentration in preparation for GC/MS or LC-MS/MS. Dilution results are not intended for patient values. Dilution procedures, when used, should be validated.

Materials provided

See "Reagents – working solutions" section for reagents.

Materials required (but not provided)

General laboratory equipment

AMPS2

ONLINE DAT Amphetamines II



See "Order information" section

Assay

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

The performance of applications not validated by Roche is not warranted and must be defined by the user.

Application for serum and plasma

Deselect Automatic Rerun for these applications in the Utility menu, Application screen, Range tab.

cobas c 501/502 test definition

	Qualitative
Assay type	2-Point End
Reaction time / Assay points	10 / 16-46
Wavelength (sub/main)	- /600 nm
Reaction direction	Increase
Unit	mAbs
Reagent pipetting	
R1	90 µL
R2	40 µL
R3	-
Sample volumes	Sample

300 ng/mL cutoff

Normal	6.0 µL
Decreased	6.0 µL
Increased	6.0 µL

Calibration

Calibrators	<i>Qualitative application</i> <i>300 ng/mL cutoff assay</i> S1: Preciset DAT Plus I calibrator - CAL 6, 5000 ng/mL with automatic pre-dilution The drug concentration of the calibrator has been verified by GC/MS.
Calibration K Factor	Enter the K Factor as -1000 into the Calibration menu, Status screen, Calibration Result window.
Calibration mode	<i>Qualitative application</i> Linear
Calibration frequency	Blank calibration - after reagent lot change - as required following quality control procedures

Calibration interval may be extended based on acceptable verification of calibration by the laboratory.

Traceability: This method has been standardized against a primary reference method (GC/MS).

Quality control

For quality control, use control materials as listed in the "Order information" section.

In addition, other suitable control material can be used.

Drug concentrations of the high and low controls have been verified by LC-MS/MS.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

Results

The cutoff calibrator is used as a reference in distinguishing between preliminary positive and negative samples. Samples producing a positive or "0" absorbance value are considered preliminary positive. Preliminary positive samples are flagged with the >Test. Samples producing a negative absorbance value are considered negative. Negative samples are preceded by a minus sign. Results of this assay distinguish preliminary positive (≥ 300 ng/mL) from negative samples only. The amount of drug detected in a preliminary positive sample cannot be estimated.

As with any sensitive test for drugs of abuse on automated clinical chemistry analyzers, the possibility exists for analyte carry-over from a sample with an extremely high concentration to a normal (negative) sample which immediately follows it.

Confirm all preliminary positive results by another method.

Limitations - interference

Criterion: No cross-over at initial values of samples of 150 ng/mL and 450 ng/mL (control levels).

See the "Specific performance data" section of this document for information on substances tested with this assay. There is the possibility that other substances and/or factors may interfere with the test and cause erroneous results (e.g., technical or procedural errors).

A preliminary positive result with this assay indicates the presence of amphetamines/methamphetamines and/or their metabolites in serum. It does not measure the level of intoxication.

Icterus:⁸ No significant interference up to an I index of 60 for conjugated and unconjugated bilirubin (approximate conjugated and unconjugated bilirubin concentration: 1026 µmol/L or 60 mg/dL).

Hemolysis:⁸ No significant interference up to an H index of 1000 (approximate hemoglobin concentration: 622 µmol/L or 1000 mg/dL).

Lipemia (native lipaemic samples):⁸ No significant interference up to an L index of 100. There is poor correlation between the L index (corresponds to turbidity) and triglycerides concentration.

Rheumatoid factors: No significant interference from rheumatoid factors up to 100 IU/mL.

Immunoglobulin: No significant interference from immunoglobulin up to a concentration of 16 g/L (simulated by human immunoglobulin A), up to a concentration of 70 g/L (simulated by human immunoglobulin G) and up to a concentration of 10 g/L (simulated by human immunoglobulin M).

Total protein: No significant interference from total protein up to a concentration of 70 g/L (simulated by human serum albumin).

As with any assay employing mouse antibodies, the possibility exists for interference by human anti-mouse antibodies (HAMA) in the sample, which could cause falsely lowered results.

In very rare cases, gammopathy, in particular type IgM (Waldenström's macroglobulinemia), may cause unreliable results.⁹

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

ACTION REQUIRED

When running Amphetamines II and Tina-quant Hemoglobin A1c II assays, on the same **cobas c** 501 analyzer, avoid processing Amphetamines II as the first test from standby status. If no other testing is pending, a dummy test sample should be processed to prevent the Amphetamines II from being the first test from standby. Order a dummy test for any R1 assay other than HbA1c II.

ACTION REQUIRED

Special Wash Programming: The use of special wash steps is mandatory when certain test combinations are run together on Roche/Hitachi **cobas c** systems. The latest version of the carry-over evasion list can be found with the NaOHD-SMS-SmpCln1+2-SCCS Method Sheets. For further instructions refer to the operator's manual. **cobas c** 502 analyzer: All special wash programming necessary for avoiding carry-over is available via the **cobas** link, manual input is not required.

Where required, special wash/carry-over evasion programming must be implemented prior to reporting results with this test.

AMPS2

ONLINE DAT Amphetamines II



Expected values

Qualitative assay

Results of this assay distinguish preliminary positive (≥ 300 ng/mL) from negative samples only. The amount of drug detected in a preliminary positive sample cannot be estimated.

No drug should be present in individuals that have not ingested amphetamine or methamphetamine.

Specific performance data

Representative performance data on the analyzers are given below. Results obtained in individual laboratories may differ.

Precision

A *d*-methamphetamine (MAMP) solution was added to 9 samples obtained from a human serum sample pool to achieve concentrations at approximately -100 %, -75 %, -50 %, -25 %, ± 0 %, +25 %, +50 %, +75 %, and +100 % of the cutoff value. These samples were tested for precision. Following a CLSI (EP5-A3) precision protocol, samples were tested in 2 replicates per run, 2 runs per day for 21 days, total $n = 84$. The following results were obtained on a Roche/Hitachi **cobas c 501** analyzer.

Drug	Concentration of Sample	Number of Determinations	Results # Neg / # Pos
MAMP	zero drug	84	84 Neg / 0 Pos
MAMP	-75 %	84	84 Neg / 0 Pos
MAMP	-50 %	84	84 Neg / 0 Pos
MAMP	-25 %	84	84 Neg / 0 Pos
MAMP	cutoff	84	42 Neg / 42 Pos
MAMP	+25 %	84	0 Neg / 84 Pos
MAMP	+50 %	84	0 Neg / 84 Pos
MAMP	+75 %	83	0 Neg / 83 Pos
MAMP	+100 %	84	0 Neg / 84 Pos

Accuracy

73 serum samples obtained from a clinical laboratory, where they screened negative in a drug test panel, were evaluated with the Amphetamines II assay. 100 % of these normal serum samples were negative relative to the 300 ng/mL cutoff.

30 samples obtained from a clinical laboratory, where they were screened preliminary positive with a commercially available immunoassay and were subsequently confirmed positive by LC-MS/MS, were evaluated with the Amphetamines II assay. 100 % of these serum samples were positive relative to the 300 ng/mL cutoff.

In addition, 21 samples were found in a concentration of 100-150 % of the cutoff concentration; and 37 samples were found in a concentration of 50-100 % of the cutoff concentration. The following results were obtained with the Amphetamines II assay on the Roche/Hitachi **cobas c 501** analyzer relative to the LC-MS/MS values.

Note: Several LC-MS/MS methods do not distinguish between the *d*- and *l*-enantiomers of amphetamines, whereas the Amphetamines II assay does. AMPS II detects only the *d*-enantiomer, which is the active CNS stimulant (see section "Analytical specificity"). For distribution of *d*- and *l*-enantiomers see reference 10.¹⁰

		n = 161			
		LC-MS/MS			
		neg	neg near cutoff	pos near cutoff	pos
cobas c 501 analyzer	neg	73	26	6	0
	pos	0	11	15	30

Analytical specificity

The specificity of Amphetamines II for various phenethylamines and structurally similar compounds was determined by generating inhibition curves for each of the compounds listed and determining the approximate quantity of each compound that is equivalent in assay reactivity to the

300 ng/mL *d*-methamphetamine assay cutoff. The following results were obtained on a Roche/Hitachi **cobas c 501** analyzer.

Compound	ng/mL Equivalent to 300 ng/mL <i>d</i> -methamphetamine	Approximate % Cross-reactivity
<i>d</i> -Amphetamine	302	99.2
<i>d</i> -Methamphetamine	307	97.9
BDB ^{a)}	915	32.8
MBDB ^{b)}	357	84.1
MDA ^{c)}	248	121
MDEA ^{d)}	554	54.1
MDMA ^{e)}	196	153
<i>l</i> -Methamphetamine	5919	5.07
Phendimetrazine	38281	0.78
Phentermine	88086	0.34
Benfluorex	> 100000	n.d.
Tyramine	93000	0.32
Benzylamine	> 100000	n.d.
Chlorpheniramine	> 100000	n.d.
Chlorpromazine	> 100000	n.d.
Chloramphetamine	135	223
<i>l</i> -Ephedrine	118341	0.25
<i>d</i> -Pseudoephedrine	104227	0.29
<i>d,l</i> -Phenylpropanolamine	299838	0.10
<i>l</i> -Amphetamine	5932	5.06
<i>l</i> -Norpseudoephedrine	99822	0.30
Cathine	30763	0.98
<i>m</i> -CPP	3088	9.72
Labetalol	> 100000	n.d.
PMA ^{f)}	120	250
PMMA ^{g)}	126	238
1-Methyl-3-phenylpropylamine (APB) ^{h)}	1438	20.9

n.d. = not detectable

a) *d,l*-3,4-Methylenedioxyphenyl-2-butamine hydrochloride

b) *d,l*-N-Methyl-1-(3,4-methylenedioxyphenyl)-2-butamine hydrochloride

c) *d,l*-3,4-Methylenedioxyamphetamine

d) *d,l*-3,4-Methylenedioxyethylamphetamine

e) *d,l*-3,4-Methylenedioxymethamphetamine

f) para-Methoxyamphetamine

g) para-Methoxymethamphetamine

h) APB, metabolite of Labetalol

Drug interference

Interfering substances were added to serum containing methamphetamine at -50 % and +50 % of the cutoff level at the concentrations listed below. Samples were tested and the following results were obtained on a Roche/Hitachi **cobas c 501** analyzer.

Compound	Compd. Conc. mg/L	Neg Level	Pos Level
Acetaminophen	200	neg	pos
Acetylcysteine	1660	neg	pos

Acetylsalicylic acid	1000	neg	pos
Amitriptyline	1.00	neg	pos
Ampicillin-Na	1000	neg	pos
Ascorbic acid	300	neg	pos
Caffeine	59.8	neg	pos
Cefoxitin	2500	neg	pos
Cyclosporine	5.00	neg	pos
Dextromethorphan	1.00	neg	pos
Doxycycline	50.0	neg	pos
Erythromycin	59.9	neg	pos
Fenoprofen	6.49	neg	pos
Furosemide	59.9	neg	pos
Gentisic acid	18.0	neg	pos
Heparin	5000 U/L	neg	pos
Hydrochlorothiazide	0.20	neg	pos
Ibuprofen	500	neg	pos
Ketamine	10.0	neg	pos
Levodopa	20.0	neg	pos
Lidocaine	12.0	neg	pos
LSD	2.50	neg	pos
Methyldopa + 1.5 H ₂ O	20.0	neg	pos
Metronidazole	200	neg	pos
Naproxen	499	neg	pos
Phenylbutazone	400	neg	pos
Procaine	2.00	neg	pos
Promethazine	1.20	neg	pos
Quinidine	12.0	neg	pos
Quinine	48.0	neg	pos
Rifampicin	60.0	neg	pos
Tetracycline	15.1	neg	pos
Theophylline	100	neg	pos
Trifluoperazine	1.00	neg	pos

References




- 1 Karch SB, ed. Drug Abuse Handbook. Boca Raton, FL: CRC Press LLC 1998.
- 2 Borenstein M. Central Nervous System Stimulants. In: Troy D. ed. Remington: The Science and Practice of Pharmacy. 21st ed. Baltimore, MD: Lippincott Williams & Wilkins Co 2005;1551.
- 3 NIDA Research Report - Methamphetamine Abuse and Addiction: NIH Publication No. 06-4210. National Institute on Drug Abuse 6001 Executive Blvd. Room 5213, Bethesda, MD 20892-9561.
- 4 Baselt RC. Disposition of Toxic Drugs and Chemicals in Man. 7th ed. Foster City, CA: Biomedical Publications 2004;67.
- 5 Armbruster DA, Schwarzhoff RH, Pierce BL, et al. Method comparison of EMIT II and ONLINE with RIA for drug screening. J Forensic Sci 1993;38:1326-1341.
- 6 Armbruster DA, Schwarzhoff RH, Hubster EC, et al. Enzyme immunoassay, kinetic microparticle immunoassay, radioimmunoassay, and fluorescence polarization immunoassay compared for drugs-of-abuse screening. Clin Chem 1993;39:2137-2146.
- 7 Rouse S, Motter K, McNally A, et al. An Abuscreen OnLine Immunoassay for the Detection of Amphetamine in Urine on the COBAS MIRA Automated Analyzer. Clin Chem 1991;37(6):995. Abstract.

- 8 Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-475.
- 9 Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med 2007;45(9):1240-1243.
- 10 Peters FT, Kraemer T, Maurer HH. Drug testing in blood: validated negative-ion chemical ionization gas chromatographic-mass spectrometric assay for determination of amphetamine and methamphetamine enantiomers and its application to toxicology cases. Clin Chem 2002;48(9):1472-1485.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard (for USA: see <https://usdiagnostics.roche.com> for definition of symbols used):

	Contents of kit
	Volume after reconstitution or mixing
	Global Trade Item Number

COBAS, COBAS C, COBAS MIRA, ONLINE DAT and PRECISET are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2017, Roche Diagnostics



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim
www.roche.com

