

**Order information**

REF	CONTENT	Analyzer(s) on which <b>cobas c</b> pack(s) can be used
04490754 190	ONLINE DAT Barbiturates Plus (200 tests)	System-ID 07 6917 7 Roche/Hitachi <b>cobas c</b> 501/502
03304671 190	Preciset DAT Plus I calibrator CAL 3	Code 433
07978766 190	Serum DAT Control Low (ACQ Partner Channel*)	
07978740 190	Serum DAT Control High (ACQ Partner Channel*)	

\*Roche does not hold the product registration for Partner Channels. The legal manufacturer indicated on the kit is solely responsible for all of the design, legal, and regulatory aspects of the product.

**English****System information**

For **cobas c** 501 analyzer:

**BAQ2S:** ACN 600: for qualitative assay, 200 ng/mL

For **cobas c** 502 analyzer:

**BAQ2S:** ACN 8600: for qualitative assay, 200 ng/mL

**Intended use**

Barbiturates Plus (BARB) is an in vitro diagnostic test for the qualitative detection of barbiturates in human serum and plasma on Roche/Hitachi **cobas c** systems at a cutoff concentration of 200 ng/mL.

The assay provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC-MS) or Liquid Chromatography coupled with Tandem Mass Spectrometry (LC-MS/MS) is the preferred confirmatory method.<sup>1</sup> Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used.

**Summary**

The barbiturates, a class of drugs derived from barbituric acid (malonylurea), are sedative hypnotics with central nervous system (CNS)-depressant activity.<sup>1,2,3,4,5,6</sup> As CNS-depressants, the barbiturates are classified relative to their duration of action (ultra short-, short-, intermediate-, and long-acting). They have been used medically as sedatives to reduce emotional tension and induce sleep, and in certain types of epilepsy to reduce seizure frequency by raising the seizure threshold. Excessive dosages may cause impaired motor coordination (slurred speech, loss of balance), perceptual alterations (faulty judgment, inflated perceptions of performance), and disinhibition euphoria. Overdoses can result in stupor, coma, and death. The combined use of the barbiturates with alcohol, opiates, or other CNS-depressants can result in fatal, additive respiratory depression. Although their utilities as sedative-hypnotic drugs have largely been replaced by the benzodiazepines, the barbiturates still maintain an important role as anesthetic and anticonvulsant drugs.

Oral administration is most common, although the barbiturates may be injected intravenously or intramuscularly. Following ingestion, they are rapidly absorbed from the stomach and enter the circulation. Their resulting distribution and concentration in various tissues is largely dependent on the lipid solubility and protein-binding characteristics of the different barbiturates; fat deposits and protein-rich tissues accumulate the highest concentration. Most of the barbiturates are metabolized by the liver via oxidation and conjugation, nitrogen-dealkylation, nitrogen-hydroxylation, and/or desulfuration of thiobarbiturates. The extent of liver metabolism is drug-dependent; secobarbital, for example, is extensively oxidized to a series of pharmacologically inactive metabolites, while a relatively high percentage of phenobarbital and barbital are excreted unchanged in the urine. As a drug class, the barbiturates are excreted as active drug/metabolite mixes whose ratios and concentrations depend on the specific barbiturate in question.

**Test principle**

The assay is based on the kinetic interaction of microparticles in a solution (KIMS)<sup>7,8</sup> as measured by changes in light transmission. In the absence of sample drug, free antibody binds to drug-microparticle conjugates causing the formation of particle aggregates. As the aggregation reaction proceeds in the absence of sample drug, the absorbance increases.

When a serum sample contains the drug in question, this drug competes with the particle-bound drug derivative for free antibody. Antibody bound to sample drug is no longer available to promote particle aggregation, and subsequent particle lattice formation is inhibited. The presence of sample

drug diminishes the increasing absorbance in proportion to the concentration of drug in the sample. Sample drug content is determined relative to the value obtained for a known cutoff concentration of drug.

**Reagents - working solutions**

**R1** Buffer; 0.09 % sodium azide

**R2** Secobarbital antibody (sheep polyclonal); buffer; bovine serum albumin; 0.09 % sodium azide

**R3** Conjugated secobarbital derivative microparticles; buffer; bovine serum albumin; 0.09 % sodium azide

R1 is in position B, R2 is in position C, and R3 is in position A.

**Precautions and warnings**

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents.

Reagents from different kit lots must not be interchanged. Reagents within kit lots have been matched to ensure optimum test performance. Disposal of all waste material should be in accordance with local guidelines. Safety data sheet available for professional user on request.

**Reagent handling**

Ready for use

Carefully invert reagent container several times prior to use to ensure that the reagent components are mixed.

**Storage and stability**

Shelf life at 2-8 °C: See expiration date on **cobas c** pack label

On-board in use and refrigerated on the analyzer: 8 weeks

**Do not freeze.****Specimen collection and preparation**

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable.

Serum: Serum tubes with and without separating gel.

Plasma: K<sub>2</sub>- or K<sub>3</sub>-EDTA, lithium heparin.

Stability: 5 days capped at 15-25 °C  
14 days capped at 2-8 °C  
6 months capped at -20 °C

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

Specimens can be repeatedly frozen and thawed up to 3 times.

Invert thawed specimens several times prior to testing.

**CAUTION:** Specimen dilutions should only be used to interpret results of Calc.? and Samp.? alarms, or when estimating concentration in preparation

for GC-MS. Dilution results are not intended for patient values. Dilution procedures, when used, should be validated.

#### Materials provided

See "Reagents – working solutions" section for reagents.

#### Materials required (but not provided)

See "Order information" section

General laboratory equipment

#### Assay

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

The performance of applications not validated by Roche is not warranted and must be defined by the user.

#### Application for serum and plasma

Deselect Automatic Rerun for these applications in the Utility menu, Application screen, Range tab.

#### cobas c 501/502 test definition

	Qualitative
Assay type	2-Point End
Reaction time / Assay points	10 / 40-65
Wavelength (sub/main)	– /505 nm
Reaction direction	Increase
Unit	mAbs
Reagent pipetting	
R1	59 µL
R2	59 µL
R3	52 µL
Sample volumes	Sample

#### 200 ng/mL cutoff

Normal	2.3 µL
Decreased	2.3 µL
Increased	2.3 µL

#### Calibration

Calibrators	Qualitative application 200 ng/mL cutoff assay S1: Preciset DAT Plus I calibrator - CAL 3 The drug concentration of the calibrator has been verified by GC-MS.
Calibration K Factor	Enter the K Factor as -1000 into the Calibration menu, Status screen, Calibration Result window.
Calibration mode	Qualitative application Linear
Calibration frequency	Blank calibration - after reagent lot change - as required following quality control procedures

Calibration interval may be extended based on acceptable verification of calibration by the laboratory.

Traceability: This method has been standardized against a primary reference method (GC-MS).

#### Quality control

For quality control, use control materials as listed in the "Order information" section.

In addition, other suitable control material can be used.

Drug concentrations of the high and low controls have been verified by LC-MS/MS.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

#### Results

The cutoff calibrator is used as a reference in distinguishing between preliminary positive and negative samples. Samples producing a positive or "0" absorbance value are considered preliminary positive. Preliminary positive samples are flagged with >Test. Samples producing a negative absorbance value are considered negative. Negative samples are preceded by a minus sign.

As with any sensitive test for drugs of abuse on automated clinical chemistry analyzers, the possibility exists for analyte carry-over from a sample with an extremely high concentration to a normal (negative) sample which immediately follows it.

Confirm all preliminary positive results by another method.

#### Limitations - interference

Criterion: No cross-over at initial values of samples of 100 ng/mL and 300 ng/mL (control levels).

See the "Specific performance data" section of this document for information on substances tested with this assay. There is the possibility that other substances and/or factors may interfere with the test and cause erroneous results (e.g., technical or procedural errors).

A preliminary positive result with this assay indicates the presence of barbiturates and/or their metabolites in serum. It does not measure the level of intoxication.

Icterus:<sup>9</sup> No significant interference up to an I index of 60 for conjugated and unconjugated bilirubin (approximate conjugated and unconjugated bilirubin concentration: 1026 µmol/L or 60 mg/dL).

Hemolysis:<sup>9</sup> No significant interference up to an H index of 1000 (approximate hemoglobin concentration: 622 µmol/L or 1000 mg/dL).

Lipemia (Intralipid):<sup>9</sup> No significant interference up to an L index of 1000. There is poor correlation between the L index (corresponds to turbidity) and triglycerides concentration.

Rheumatoid factors: No significant interference from rheumatoid factors up to 1200 IU/mL.

Immunoglobulin: No significant interference from immunoglobulin up to a concentration of 16 g/L (simulated by human immunoglobulin A), up to a concentration of 70 g/L (simulated by human immunoglobulin G) and up to a concentration of 10 g/L (simulated by human immunoglobulin M).

Total protein: No significant interference from total protein up to a concentration of 70 g/L (simulated by human serum albumin).

In very rare cases, gammopathy, in particular type IgM (Waldenström's macroglobulinemia), may cause unreliable results.<sup>10</sup>

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

#### ACTION REQUIRED

**Special Wash Programming:** The use of special wash steps is mandatory when certain test combinations are run together on Roche/Hitachi cobas c systems. The latest version of the carry-over evasion list can be found with the NaOHD-SMS-SmpCln1+2-SCCS Method Sheets. For further instructions refer to the operator's manual. cobas c 502 analyzer: All special wash programming necessary for avoiding carry-over is available via the cobas link, manual input is not required.

**Where required, special wash/carry-over evasion programming must be implemented prior to reporting results with this test.**

#### Expected values

##### Qualitative assay

Results of this assay distinguish preliminary positive ( $\geq 200$  ng/mL) from negative samples only. The amount of drug detected in a preliminary positive sample cannot be estimated.

**Specific performance data**

Representative performance data on the analyzers are given below. Results obtained in individual laboratories may differ.

**Precision**

A secobarbital solution was added to 9 samples obtained from a human serum sample pool to achieve concentrations at approximately -100 %, -75 %, -50 %, -25 %, ±0 %, +25 %, +50 %, +75 %, and +100 % of the cutoff value. These samples were tested for precision. Following a CLSI (EP5-A3) precision protocol, samples were tested in 2 replicates per run, 2 runs per day for 21 days, total n = 84. The following results were obtained on a Roche/Hitachi **cobas c 501** analyzer:

Drug	Concentration of Sample	Number of Determinations	Results # Neg / # Pos
Secobarbital	zero drug	84	84 Neg / 0 Pos
Secobarbital	-75 %	84	84 Neg / 0 Pos
Secobarbital	-50 %	84	84 Neg / 0 Pos
Secobarbital	-25 %	84	84 Neg / 0 Pos
Secobarbital	Cutoff	84	78 Neg / 6 Pos
Secobarbital	+25 %	84	0 Neg / 84 Pos
Secobarbital	+50 %	84	0 Neg / 84 Pos
Secobarbital	+75 %	83	0 Neg / 83 Pos
Secobarbital	+100 %	84	0 Neg / 84 Pos

**Accuracy**

72 serum samples obtained from a clinical laboratory, where they screened negative in a drug test panel, were evaluated with the Barbiturates Plus assay. 100 % of these normal serum samples were negative relative to the 200 ng/mL cutoff.

63 samples obtained from a clinical laboratory, where they were screened preliminary positive with a commercially available immunoassay and were subsequently confirmed positive by GC-MS, were evaluated with the Barbiturates Plus assay. 100 % of these serum samples were positive relative to the 200 ng/mL cutoff.

In addition, 7 samples were found in a concentration of 100-150 % of the cutoff concentration; and 8 samples were found in a concentration of 50-100 % of the cutoff concentration. The following results were obtained with the Barbiturates Plus assay on the Roche/Hitachi **cobas c 501** analyzer relative to the GC-MS values.

		n = 150			
		GC-MS			
		neg	neg near cutoff	pos near cutoff	pos
<b>cobas c 501</b> analyzer	neg	72	8	0	0
	pos	0	0	7	63

**Analytical specificity**

The specificity of this assay for some common barbiturates and structurally similar compounds was determined by generating inhibition curves for each of the compounds listed and determining the approximate quantity of each compound that is equivalent in assay reactivity to a 200 ng/mL secobarbital assay cutoff. The following results were obtained on a Roche/Hitachi **cobas c 501** analyzer.

Compound	ng/mL Equivalent to 200 ng/mL Secobarbital	Approximate % Cross-reactivity
Allobarbitol	211	94.7
Amobarbital	466	42.9
Aprobarbital	193	104
Barbital	340	58.9
Butabarbital	326	61.4

Butalbital	350	57.2
Cyclopentobarbital	137	146
Pentobarbital	333	60.0
Phenobarbital	312	64.1
p-Hydroxyphenobarbital	590	33.9

**Drug interference**

Interfering substances were added to serum containing secobarbital at -50 % and +50 % of the cutoff level at the concentration listed below. Samples were tested and the following results were obtained on a Roche/Hitachi **cobas c 501** analyzer.

Compound	Comp. Conc. mg/L	Neg Level	Pos Level
Acetaminophen	200	neg	pos
Acetylcysteine	1660	neg	pos
Acetylsalicylic acid	1000	neg	pos
Amitriptyline	1.00	neg	pos
Ampicillin-Na	1000	neg	pos
Ascorbic acid	300	neg	pos
Caffeine	59.8	neg	pos
Cefoxitin	2500	neg	pos
Cyclosporine	5.00	neg	pos
Doxycycline	50.0	neg	pos
d-Pseudoephedrine	9.98	neg	pos
Erythromycin	59.9	neg	pos
Fenopropfen	195	neg	pos
Furosemide	59.9	neg	pos
Gentisic acid	18.0	neg	pos
Heparin	5000 U/L	neg	pos
Hydrochlorothiazide	6.02	neg	pos
Ibuprofen	500	neg	pos
Imipramine	0.70	neg	pos
Ketamine	10.0	neg	pos
Levodopa	20.0	neg	pos
Lidocaine	12.0	neg	pos
Methyldopa + 1.5 H <sub>2</sub> O	20.0	neg	pos
Metronidazole	200	neg	pos
Naproxen	499	neg	pos
Phenylbutazone	400	neg	pos
Procaine	39.9	neg	pos
Promethazine	1.20	neg	pos
Quinidine	12.0	neg	pos
Quinine	48.0	neg	pos
Rifampicin	60.0	neg	pos
Tetracycline	15.1	neg	pos
Theophylline	100	neg	pos
Trifluoperazine	1.00	neg	pos

**References**




- 1 Karch SB, ed. Drug Abuse Handbook. Boca Raton, FL: CRC Press LLC 1998.
- 2 Wesson DR, Smith DE. Barbiturates: Their Use, Misuse, and Abuse. New York, NY: Human Sciences Press 1977.

- 3 Robinson AE, McDowall RD. The distribution of amylobarbitone, butobarbitone, pentobarbitone and quinalbarbitone and the hydroxylated metabolites in man. J Pharm Pharmacol 1979;31:357-365.
- 4 Hardman JG, Limbird LE, Gilman A, eds. Goodman and Gilman's The Pharmacological Basis of Therapeutics. 10th ed. New York, NY: McGraw Hill Pub Co. 2001.
- 5 Baselt RC. Disposition of Toxic Drugs and Chemicals in Man. 7th ed. Foster City, CA: Biomedical Publications 2004.
- 6 Barbiturates - A Medical Dictionary, Bibliography, and Annotated Research Guide to Internet Reference. San Diego, CA: ICON Group International Inc 2004.
- 7 Armbruster DA, Schwarzhoff RH, Pierce BL, et al. Method comparison of EMIT II and ONLINE with RIA for drug screening. J Forensic Sci 1993;38:1326-1341.
- 8 Armbruster DA, Schwarzhoff RH, Hubster EC, et al. Enzyme immunoassay, kinetic microparticle immunoassay, radioimmunoassay, and fluorescence polarization immunoassay compared for drugs-of-abuse screening. Clin Chem 1993;39:2137-2146.
- 9 Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-475.
- 10 Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med 2007;45(9):1240-1243.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

### Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard (for USA: see <https://usdiagnostics.roche.com> for definition of symbols used):

	Contents of kit
	Volume after reconstitution or mixing
	Global Trade Item Number

COBAS, COBAS C, ONLINE DAT and PRECISET are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2017, Roche Diagnostics



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim  
[www.roche.com](http://www.roche.com)

