



07660839001V2.0

# Cell Wash Solution I / NaOH-D

**cobas®**

CellCln1/Detergent 1

## Order information

REF	CONTENT	Analyzer(s) on which kit(s) can be used
07531869 190	Cell Wash Solution I / NaOH-D (0.9 L in a 2 L bottle)	Roche/Hitachi MODULAR Roche/Hitachi <b>cobas c</b>

**NOTE:** Users must order an even number of Cell Wash Solution I / NaOH-D bottles.

## English

### ACTION REQUIRED

Users must combine 2 bottles of Cell Wash Solution I / NaOH-D Cat. no. 07531869 190 (0.9 L each) using the provided funnel to obtain a total volume of 1.8 L before usage on the analyzer.

Bottles containing only 0.9 L of solution must not be used on the analyzer as this could lead to air being drawn into the system.

### Intended use

Cell Wash Solution I / NaOH-D is used as alkaline wash solution for reaction cells on Roche/Hitachi systems.

### Summary

Routine reaction cell (cuvette) washes are required due to potential interference from other reagent/sample mixtures. This procedure maintain the integrity of patient sample results.

### Reagents - working solutions

Sodium hydroxide solution 1 mol/L, 4 %; detergent

### Precautions and warnings

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents.

Disposal of all waste material should be in accordance with local guidelines. Safety data sheet available for professional user on request.

For USA: Caution: Federal law restricts this device to sale by or on the order of a physician.

This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:



Danger

H290 May be corrosive to metals.

H314 Causes severe skin burns and eye damage.

### Prevention:

P280 Wear protective gloves/ protective clothing/ eye protection/ face protection.

### Response:

P301 + P330 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.  
+ P331

P303 + P361 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.  
+ P353

P304 + P340 IF INHALED: Remove person to fresh air and keep comfortable for breathing.  
+ P310 Immediately call a POISON CENTER/ doctor.

P305 + P351 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.  
+ P338 Continue rinsing. Immediately call a POISON CENTER/ doctor.  
+ P310

P390 Absorb spillage to prevent material damage.

Product safety labeling follows EU GHS guidance.

Contact phone: all countries: +49-621-7590, USA: 1-800-428-2336

## Reagent handling

### ACTION REQUIRED

Users must combine 2 bottles of Cell Wash Solution I / NaOH-D Cat. no. 07531869 190 (0.9 L each) using the provided funnel to obtain a total volume of 1.8 L before usage on the analyzer.

Bottles containing only 0.9 L of solution must not be used on the analyzer as this could lead to air being drawn into the system.

### NOTE:

Over time, a slight discoloration or small particles may be observed on the bottom of the bottle. The discoloration or particles do not affect reagent performance.

### Storage and stability

Shelf life at 15-25 °C: See expiration date on labels.

On-board stability: 10 weeks

### Materials provided

See "Reagents – working solutions" section for reagents.

### Funnel

### Materials required (but not provided)

See "Order information" section

Roche/Hitachi analyzers and assay reagents.

See appropriate instructions for use and operator's manuals for additional required materials.

### Assay

For further information, please refer to the appropriate operator's manual for the analyzer concerned, the respective application sheets and method sheets of all necessary components.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

### Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard (for USA: see <https://usdiagnostics.roche.com> for definition of symbols used):

<b>DETERGENT 1</b>	Detergent 1
<b>CONTENT</b>	Contents of kit
<b>REAGENT</b>	Reagent
<b>CALIBRATOR</b>	Calibrator
	Volume after reconstitution or mixing
<b>GTIN</b>	Global Trade Item Number

COBAS and COBAS C are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2017, Roche Diagnostics



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim  
[www.roche.com](http://www.roche.com)

