



12174901001V7.0

ISE Clean. Soln./Elecsys SysClean cobas®

ISE Cleaning Solution / Elecsys SysClean

Order information

REF	CONTENT	Analyzer(s) on which kit(s) can be used
11298500 316	ISE Cleaning Solution / Elecsys SysClean (5 × 100 mL)	Roche/Hitachi Roche/Hitachi cobas c MODULAR ANALYTICS E170 cobas e

English

Intended use

For the cleaning of ISE units on Roche/Hitachi analyzers.
For the cleaning of Elecsys and **cobas e** immunoassay analyzers.

Summary

ISE Cleaning Solution / Elecsys SysClean is an alkaline cleaning solution with antimicrobial properties.

Reagents - working solutions

5 bottles, each containing 100 mL

Components:

Sodium hydroxide 3 mol/L

Sodium hypochlorite solution (< 2 % active chlorine)

Additive

Precautions and warnings

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents.

Disposal of all waste material should be in accordance with local guidelines.

Safety data sheet available for professional user on request.

This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:



Danger

H290 May be corrosive to metals.

H314 Causes severe skin burns and eye damage.

Prevention:

P280 Wear protective gloves/ protective clothing/ eye protection/ face protection.

Response:

P301 + P330 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.
+ P331

P303 + P361 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.
+ P353

P304 + P340 IF INHALED: Remove person to fresh air and keep comfortable for breathing.
+ P310 Immediately call a POISON CENTER/doctor.

P305 + P351 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.
+ P338 Continue rinsing. Immediately call a POISON CENTER/doctor.
+ P310

P390 Absorb spillage to prevent material damage.

Product safety labeling follows EU GHS guidance.

Contact phone: all countries: +49-621-7590

Reagent handling

Ready for use

Storage and stability

Store protected from light.

The solution is stable up to the stated expiration date when stored at 2-8 °C.

Materials provided

See "Reagents – working solutions" section for reagents.

Materials required (but not provided)

- General laboratory equipment
- Roche/Hitachi, Roche/Hitachi **cobas c**, MODULAR ANALYTICS E170 or **cobas e** analyzers

Assay

For further information, please refer to the appropriate operator's manual for the analyzer concerned, the respective application sheets and method sheets of all necessary components.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard (for USA: see <https://usdiagnostics.roche.com> for definition of symbols used):

CONTENT	Contents of kit
SYSTEM	Analyzers/Instruments on which reagents can be used
REAGENT	Reagent
CALIBRATOR	Calibrator
→	Volume after reconstitution or mixing
GTIN	Global Trade Item Number

COBAS, COBAS C, COBAS E and ELECSYS are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2016, Roche Diagnostics



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim
www.roche.com

