

Order information

REF	CONTENT	Analyzer(s) on which cobas c pack(s) can be used
04851013 190	Instrument Check (65 tests)	System-ID 07 6933 9 Roche/Hitachi cobas c 311, cobas c 501/502
20757144 322	Check Solution Sample	

English

System information

For **cobas c** 311/501 analyzers:

CHKS: ACN 473 (check of sample pipetting)

CHKR1: ACN 474 (check of R1 pipetting)

CHKR2: ACN 475 (check of R2 pipetting)

For **cobas c** 502 analyzer:

CHKS: ACN 8473 (check of sample pipetting)

CHKR1: ACN 8474 (check of R1 pipetting)

CHKR2: ACN 8475 (check of R2 pipetting)

Intended use

Test for accuracy and precision of sample and reagent pipetting with dye solution and saline solution on Roche/Hitachi **cobas c** systems.

Summary

Check Test Sample (CHKS):

Sample (diluted Check Solution Sample) and R3 (sodium chloride solution) are pipetted and absorbance is measured by the instrument. Mean and CV of 21 determinations are calculated and compared to the specifications given below (accuracy and precision).

Check Test R1 (CHKR1):

Dye solution (R1) and water are pipetted to sample (water) and absorbance is measured by the instrument. Mean and CV of 21 determinations are calculated and compared to the specifications given below (accuracy and precision).

Check Test R2 (CHKR2):

Dye solution (R2) and water are pipetted to sample (water) and absorbance is measured by the instrument. Mean and CV of 21 determinations are calculated and compared to the specifications given below (accuracy and precision).

Test principle

Sample and reagent are pipetted and mixed, and the absorbance is measured photometrically.

Reagents - working solutions

R1 Orange G: 45 mg/L; TRIS: 3 g/L, pH 7.5; Tween 80: 0.25 g/L; preservatives

R2 Orange G: 45 mg/L; TRIS: 3 g/L, pH 7.5; Tween 80: 0.25 g/L; preservatives

R3 Sodium chloride: 1.54 mol/L

Precautions and warnings

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents.

Disposal of all waste material should be in accordance with local guidelines. Safety data sheet available for professional user on request.

For USA: Caution: Federal law restricts this device to sale by or on the order of a physician.

Reagent handling

Ready for use

Storage and stability

INSTC

Shelf life at 2-8 °C: See expiration date on **cobas c** pack label.

On-board in use and refrigerated on the analyzer: 7 days

Specimen collection and preparation

Sample:

Check Solution Sample

Ingredients: Orange G: 5.2 g/L; TRIS: 3 g/L; Bovine serum albumin: 40 g/L; preservative

Before measurement of the Check Test Sample dilute the Check Solution Sample exactly 1 + 9 with distilled/deionized water.

Stability of dilution: 1 day at 15-25 °C. Store in a closed container.

Materials provided

See "Reagents – working solutions" section for reagents.

Materials required (but not provided)

See "Order information" section

General laboratory equipment

Assay

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

The performance of applications not validated by Roche is not warranted and must be defined by the user.

Applications

cobas c 311 test definition CHKS

Assay type	1 Point		
Reaction time/Assay points	7 / 30		
Wavelength (sub/main)	660/480 nm		
Reaction direction	Increase		
Units	mABS		
Reagent pipetting		Diluent (H ₂ O)	
R3	150 µL	–	
Sample volume	Sample	Sample dilution	
		Sample	Diluent (H ₂ O)
Normal	2.0 µL	–	–

cobas c 501/502 test definition CHKS

Assay type	1 Point		
Reaction time/Assay points	7 / 45		
Wavelength (sub/main)	660/480 nm		
Reaction direction	Increase		
Units	mABS		
Reagent pipetting		Diluent (H ₂ O)	
R3	150 µL	–	
Sample volume	Sample	Sample dilution	
		Sample	Diluent (H ₂ O)
Normal	2.0 µL	–	–

cobas c 311 test definition CHKR1

Assay type	1 Point
Reaction time/Assay points	4 / 14
Wavelength (sub/main)	660/480 nm
Reaction direction	Increase
Units	mABS

Reagent pipetting		Diluent (H ₂ O)
R1	80 µL	100 µL
<i>Sample volume</i>	<i>Sample</i>	<i>Sample dilution</i>
		Sample Diluent (H ₂ O)
Normal	15.0 µL	– –

cobas c 501/502 test definition CHKR1

Assay type	1 Point
Reaction time/Assay points	3 / 20
Wavelength (sub/main)	660/480 nm
Reaction direction	Increase
Units	mABS
Reagent pipetting	Diluent (H ₂ O)
R1	80 µL 100 µL
<i>Sample volume</i>	<i>Sample</i> <i>Sample dilution</i>
	Sample Diluent (H ₂ O)
Normal	15.0 µL – –

cobas c 311 test definition CHKR2

Assay type	1 Point
Reaction time/Assay points	6 / 25
Wavelength (sub/main)	660/480 nm
Reaction direction	Increase
Units	mABS
Reagent pipetting	Diluent (H ₂ O)
R2	40 µL 140 µL
<i>Sample volume</i>	<i>Sample</i> <i>Sample dilution</i>
	Sample Diluent (H ₂ O)
Normal	15.0 µL – –

cobas c 501/502 test definition CHKR2

Assay type	1 Point
Reaction time/Assay points	6 / 36
Wavelength (sub/main)	660/480 nm
Reaction direction	Increase
Units	mABS
Reagent pipetting	Diluent (H ₂ O)
R2	40 µL 140 µL
<i>Sample volume</i>	<i>Sample</i> <i>Sample dilution</i>
	Sample Diluent (H ₂ O)
Normal	15.0 µL – –

Calibration

Calibrator	S1: H ₂ O (concentration 0.0)
Calibration mode	Linear
	Calibration is not performed. Calibration values 0 for S1 ABS and 10000 for K-Factor are entered manually.

Specifications of performance

Refer to target values for accuracy given on the **cobas c** pack (standardized for each new lot).

	Recovery of target value	CV
CHKS	+/- 8 %	≤ 1.5 %

	Recovery of target value	CV
CHKR1	+/- 5 %	≤ 0.5 %
CHKR2	+/- 5 %	≤ 1.0 %

Note: In case the target values are not reached the respective precision run (n = 21) has to be repeated.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard (for USA: see <https://usdiagnostics.roche.com> for definition of symbols used):

	Contents of kit
	Volume after reconstitution or mixing
	Global Trade Item Number

FOR US CUSTOMERS ONLY: LIMITED WARRANTY

Roche Diagnostics warrants that this product will meet the specifications stated in the labeling when used in accordance with such labeling and will be free from defects in material and workmanship until the expiration date printed on the label. THIS LIMITED WARRANTY IS IN LIEU OF ANY OTHER WARRANTY, EXPRESS OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. IN NO EVENT SHALL ROCHE DIAGNOSTICS BE LIABLE FOR INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES.

COBAS and COBAS C are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2017, Roche Diagnostics



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim
www.roche.com

Distribution in USA by:
Roche Diagnostics, Indianapolis, IN
US Customer Technical Support 1-800-428-2336

