

LAMB2

Tina-quant Lambda Gen.2



REF	CONTENT	Analyzer(s) on which kit(s) can be used
06749992 190	Tina-quant Lambda Gen.2 100 tests	System ID 07 6813 8 Roche/Hitachi cobas c 311, cobas c 501/502
11355279 216	Calibrator f.a.s. Proteins (5 x 1 mL)	Code 656
10557897 122	Precinorm Protein (3 x 1 mL)	Code 302
11333127 122	Precipath Protein (3 x 1 mL)	Code 303
04489357 190	Diluent NaCl 9 % (50 mL)	System ID 07 6869 3
05117003 190	PreciControl ClinChem Multi 1 (20 x 5 mL)	Code 391
05947626 190	PreciControl ClinChem Multi 1 (4 x 5 mL)	Code 391
05117216 190	PreciControl ClinChem Multi 2 (20 x 5 mL)	Code 392
05947774 190	PreciControl ClinChem Multi 2 (4 x 5 mL)	Code 392

English

System information

For **cobas c 311/501** analyzers:

LAMB2: ACN 284

For **cobas c 502** analyzer:

LAMB2: ACN 8284

Intended use

Immunoturbidimetric in vitro assay for the quantitative determination of bound and free immunoglobulins of the lambda light chain type in human serum and plasma on Roche/Hitachi **cobas c** systems.

Summary

References^{1,2,3,4,5,6,7}

Measurement of the various amounts of the different types of light chains aids in the diagnosis of multiple myeloma, lymphocytic neoplasms, Waldenström's macroglobulinemia, and connective tissue diseases such as rheumatoid arthritis or systemic lupus erythematosus.

Every plasma cell clone normally produces a uniform immunoglobulin molecule of the kappa or lambda light chain type. The kappa:lambda ratio in serum is normally around 2:1.

Pathological increases of a cell clone lead to elevated formation of monoclonal immunoglobulins or immunoglobulin fragments (free light chains), which bring about a change in the kappa:lambda ratio. A kappa:lambda ratio outside the normal range is indicative of monoclonal gammopathy.

This test encompasses both bound and free immunoglobulins of the light chain type.

It is known that the so-called paraproteins secreted in monoclonal gammopathies (monoclonal immunoglobulinemia) may differ from the respective immunoglobulins of polyclonal origin in amino acid composition and size. This may impair the binding to antibody and consequently cause antigen excess below the limits determined with immunoglobulins of polyclonal origin. Antigen excess may be detected after appropriate dilution of such samples.

Furthermore, the occurrence of two monoclonal gammopathies producing differing light chain types could theoretically lead to kappa:lambda ratios in the normal range.

Accordingly, quantitative determination of the kappa and lambda light chains cannot completely replace high-resolution electrophoresis, immunoelectrophoresis or immunofixation electrophoresis in the diagnosis of monoclonal gammopathy.

Test principle

Immunoturbidimetric assay

Anti-lambda antibodies react with the antigen in the sample to form antigen/antibody complexes that, following agglutination, are measured turbidimetrically.

Reagents - working solutions

R1	TRIS/HCl buffer: 50 mmol/L, pH 8.0; PEG 7 %; stabilizers and preservative
R2	Polyclonal anti-human lambda antibody (goat): dependent on titer; TRIS/HCl buffer: 20 mmol/L, pH 7.5; stabilizers and preservative

R1 is in position B and R2 is in position C.

Precautions and warnings

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents.

Disposal of all waste material should be in accordance with local guidelines. Safety data sheet available for professional user on request.

Storage and stability

LAMB2

Shelf life at 2-8 °C: See expiration date on **cobas c** pack label.

On-board in use and refrigerated on the analyzer: 12 weeks

Diluent NaCl 9 %

Shelf life at 2-8 °C: See expiration date on **cobas c** pack label

On-board in use and refrigerated on the analyzer: 12 weeks

Specimen collection and preparation

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable.

Serum

Plasma: Li-heparin or K₂-EDTA plasma.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

Stability:⁸ 7 days at 15-25 °C
4 weeks at 4-8 °C
2 months at (-15)-(-25) °C

Materials provided

See "Reagents – working solutions" section for reagents.

Materials required (but not provided)

- See "Order information" section
- General laboratory equipment

volume by a factor of 0.25. The results are automatically divided by this factor.

Lower limits of measurement

Lower detection limit of the test

0.2 g/L (20 mg/dL)

The lower detection limit represents the lowest measurable analyte level that can be distinguished from zero. It is calculated as the value lying three standard deviations above that of the lowest standard (standard 1 + 3 SD, repeatability, n = 21).

Expected values

Reference¹

	lambda	kappa/lambda ratio
Serum ^{6,13,14,15}	0.93-2.42 g/L	1.17-2.93

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

Specific performance data

Representative performance data on the analyzers are given below. Results obtained in individual laboratories may differ.

Precision

Precision was determined using human samples and controls in an internal protocol with repeatability (n = 21) and intermediate precision (3 aliquots per run, 1 run per day, 21 days). The following results were obtained:

Repeatability	Mean g/L (mg/dL)	SD g/L (mg/dL)	CV %
Precinorm Protein	1.08 (108)	0.02 (2)	1.7
Precipath Protein	2.26 (226)	0.02 (2)	1.1
Human serum 1	0.88 (88)	0.01 (1)	1.3
Human serum 2	2.39 (239)	0.03 (3)	1.1
Intermediate precision	Mean g/L (mg/dL)	SD g/L (mg/dL)	CV %
Precinorm Protein	1.06 (106)	0.02 (2)	1.9
Precipath Protein	2.30 (230)	0.03 (3)	1.4
Human serum 1	0.88 (88)	0.02 (2)	2.5
Human serum 2	2.41 (241)	0.03 (3)	1.2

Method comparison

Lambda light chain values for human serum samples obtained on a Roche/Hitachi **cobas c 501** analyzer (y) were compared with those determined using the same reagent on a Roche/Hitachi 917 analyzer (x).

Sample size (n) = 268

Passing/Bablok ¹⁶	Linear regression
$y = 0.959x + 0.023 \text{ g/L}$	$y = 0.939x + 0.077 \text{ g/L}$
$r = 0.943$	$r = 0.996$

The sample concentrations were between 0.558 and 6.37 g/L (55.8 and 637 mg/dL).

References

- Skvaril F, Barandum S, Morell A, et al. Imbalances of kappa/lambda immunoglobulin light chain ratios in normal individuals and in immunodeficient patients. In: Proteides of biological fluids, Peeters H, ed 1975;23:415-420.
- Sun T, de Szalay H, Lien YY, et al. Quantitation of kappa and lambda light chains for the detection of monoclonal gammopathy. J Clin Lab Anal 1988;2:84-90.
- Whicher JT, Wallage M, Fifield R. Use of immunoglobulin heavy- and light-chain measurements compared with existing techniques as a means of typing monoclonal immunoglobulins. Clin Chem 1987;33:1771-1773.

- Keren DF, Warren JS, Lowe JB. Strategy to diagnose monoclonal gammopathies in serum: high-resolution electrophoresis, immunofixation and kappa/lambda quantification. Clin Chem 1988;34:2196-2201.
- Duc J, Morel B, Peitrequin R, et al. Identification of monoclonal gammopathies: a comparison of immunofixation, immunoelectrophoresis and measurements of kappa- and lambda-immunoglobulin levels. J Clin Lab Immunol 1988;26:141-146.
- Lievens M. Medical and technical usefulness of measurement of kappa and lambda immunoglobulin light chains in serum with an M-component. J Clin Chem Clin Biochem 1989;27:519-523.
- Whicher JT, Ritchie RF, Johnson AM, et al. New international reference preparation for proteins in human serum (RPPHS). Clin Chem 1994;40:934-938.
- Wu AHB, ed. Tietz Clinical Guide to Laboratory Tests, 4th ed. St. Louis (MO): Saunders Elsevier 2006:674-675.
- Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-475.
- Breuer J. Report on the Symposium "Drug effects in Clinical Chemistry Methods". Eur J Clin Chem Clin Biochem 1996;34:385-386.
- Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385.
- Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med 2007;45(9):1240-1243.
- Hafner G, Endler T, Oppitz M, et al. Effects of standardization with the new international reference preparation for proteins in human serum on method comparability and reference values. Clin Lab 1995;41:743-748.
- Jones RG, Aguzzi F, Bienvenu J, et al. Use of Immunoglobulin Heavy-chain and Light-chain measurement in a multicenter trial to investigate Monoclonal components: I. Detection. Clin Chem 1991;37:1917-1921.
- Jones RG, Aguzzi F, Bienvenu J, et al. Use of Immunoglobulin Heavy-chain and Light-chain measurement in a multicenter trial to investigate Monoclonal components: II. Classification by use of Computer-based algorithms. Clin Chem 1991;37:1922-1926.
- Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard.

CONTENT



Contents of kit

Volume after reconstitution or mixing

COBAS, COBAS C, PRECINORM, PRECIPATH, PRECICONTROL and TINA-QUANT are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Significant additions or changes are indicated by a change bar in the margin.

© 2013, Roche Diagnostics



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim
www.roche.com

