

**Order information**

REF	CONTENT	Analyzer(s) on which <b>cobas c</b> pack(s) can be used
<b>20767123</b> 322	Lactate Dehydrogenase optimized 300 tests	System-ID 07 6712 3 Roche/Hitachi <b>cobas c</b> 311, <b>cobas c</b> 501/502
<b>10759350</b> 190	Calibrator f.a.s. (12 x 3 mL)	Code 401
<b>10759350</b> 360	Calibrator f.a.s. (12 x 3 mL, for USA)	Code 401
<b>12149435</b> 122	Precinorm U plus (10 x 3 mL)	Code 300
<b>12149443</b> 122	Precipath U plus (10 x 3 mL)	Code 301
<b>10171743</b> 122	Precinorm U (20 x 5 mL)	Code 300
<b>10171735</b> 122	Precinorm U (4 x 5 mL)	Code 300
<b>10171778</b> 122	Precipath U (20 x 5 mL)	Code 301
<b>10171760</b> 122	Precipath U (4 x 5 mL)	Code 301
<b>05117003</b> 190	PreciControl ClinChem Multi 1 (20 x 5 mL)	Code 391
<b>05947626</b> 190	PreciControl ClinChem Multi 1 (4 x 5 mL)	Code 391
<b>05947626</b> 160	PreciControl ClinChem Multi 1 (4 x 5 mL, for USA)	Code 391
<b>05117216</b> 190	PreciControl ClinChem Multi 2 (20 x 5 mL)	Code 392
<b>05947774</b> 190	PreciControl ClinChem Multi 2 (4 x 5 mL)	Code 392
<b>05947774</b> 160	PreciControl ClinChem Multi 2 (4 x 5 mL, for USA)	Code 392

**English****System information**For **cobas c** 311/501 analyzers:**LDHL:** ACN 672For **cobas c** 502 analyzer:**LDHL:** ACN 8672**Intended use**In vitro test for the quantitative determination of lactate dehydrogenase in serum and plasma on Roche/Hitachi **cobas c** systems.**Summary**<sup>1,2,3</sup>

The lactate dehydrogenase (LDH) enzyme is widely distributed in tissue, particularly heart, liver, muscle, and kidney. The LDH in serum can be separated into five different isoenzymes based on their electrophoretic mobility. Each isoenzyme is a tetramer composed of two different subunits. These two subunits have been designated heart and muscle, based on their polypeptide chains. There are two homotetramers, LDH-1 (heart) and LDH-5 (muscle), and three hybrid isoenzymes.

Elevated serum levels of LDH have been observed in a variety of disease states. The highest levels are seen in patients with megaloblastic anemia, myocardial infarction, disseminated carcinoma, leukemia, and trauma. Mild increases in LDH activity have been reported in cases of hemolytic anemias, muscular dystrophy, pulmonary infarction, hepatitis, nephrotic syndrome, and cirrhosis.

**Test principle**Optimized standard method according to the Deutsche Gesellschaft für Klinische Chemie (DGKC).<sup>4,5</sup>LDH catalyzes the reaction between pyruvate and NADH to form L-lactate and NAD<sup>+</sup>.

The initial rate of the NADH oxidation is directly proportional to the catalytic LDH activity. It is determined by measuring the decrease in absorbance.

**Reagents - working solutions****R1** Phosphate buffer: 68 mmol/L, pH 7.5; pyruvate:  $\geq 0.73$  mmol/L; preservative**R2** NADH:  $\geq 1.1$  mmol/L; preservative

R1 is in position B and R2 is in position C.

**Precautions and warnings**

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents.

Disposal of all waste material should be in accordance with local guidelines. Safety data sheet available for professional user on request.

**Reagent handling**

Ready for use

**Storage and stability****LDHL**

Shelf life at 2-8 °C:

See expiration date on **cobas c** pack label.

On-board in use and refrigerated on the analyzer: 8 weeks

**Specimen collection and preparation**

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable.

Serum (free from hemolysis)

Plasma: Li-heparin plasma (free from hemolysis)

Plasma may be contaminated with platelets which contain high concentrations of lactate dehydrogenase and should be avoided.<sup>6,7</sup>

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

Stability:<sup>8</sup>

7 days at 15-25 °C

The sample may be stored for 4 days at 2-8 °C or 6 weeks at -20 °C. In connection with certain diseases (e.g. hepatopathy, diseases of skeletal muscle, malignant tumors), the LDH-4 and LDH-5 isoenzyme portions are increased and unstable in cooled and frozen samples; this may lead to an incorrect LDH value in samples collected from patients suffering from such diseases.

**Materials provided**

See "Reagents – working solutions" section for reagents.

**Materials required (but not provided)**

- See "Order information" section
- General laboratory equipment



**Assay**

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

The performance of applications not validated by Roche is not warranted and must be defined by the user.

**Application for serum and plasma****cobas c 311 test definition**

Assay type	Rate A		
Reaction time / Assay points	10 / 18-28		
Wavelength (sub/main)	546/340 nm		
Reaction direction	Decrease		
Units	U/L (μkat/L)		
Reagent pipetting	Diluent (H <sub>2</sub> O)		
R1	100 μL	—	
R2	20 μL	—	
Sample volumes	Sample	Sample dilution	
		Sample	Diluent (H <sub>2</sub> O)
Normal	2.8 μL	—	—
Decreased	1.1 μL	—	—
Increased	2.8 μL	—	—

**cobas c 501 test definition**

Assay type	Rate A		
Reaction time / Assay points	10 / 26-40		
Wavelength (sub/main)	546/340 nm		
Reaction direction	Decrease		
Units	U/L (μkat/L)		
Reagent pipetting	Diluent (H <sub>2</sub> O)		
R1	100 μL	—	
R2	20 μL	—	
Sample volumes	Sample	Sample dilution	
		Sample	Diluent (H <sub>2</sub> O)
Normal	2.8 μL	—	—
Decreased	1.1 μL	—	—
Increased	2.8 μL	—	—

**cobas c 502 test definition**

Assay type	Rate A		
Reaction time / Assay points	10 / 26-40		
Wavelength (sub/main)	546/340 nm		
Reaction direction	Decrease		
Units	U/L (μkat/L)		
Reagent pipetting	Diluent (H <sub>2</sub> O)		
R1	100 μL	—	
R2	20 μL	—	
Sample volumes	Sample	Sample dilution	
		Sample	Diluent (H <sub>2</sub> O)
Normal	2.8 μL	—	—
Decreased	1.1 μL	—	—
Increased	5.6 μL	—	—

**Calibration**

Calibrators	S1: H <sub>2</sub> O S2: C.f.a.s.
Calibration mode	Linear
Calibration frequency	2-point calibration • after reagent lot change • as required following quality control procedures

**Traceability:** This method has been standardized against the Roche system reagent using calibrated pipettes together with a manual photometer providing absolute values and substrate specific absorptivity, ε.

**Quality control**

For quality control, use control materials as listed in the "Order information" section.

In addition, other suitable control material can be used.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

**Calculation**

Roche/Hitachi **cobas c** systems automatically calculate the analyte concentration of each sample.

Conversion factor: U/L x 0.0167 = μkat/L

**Limitations - interference**

**Criterion:** Recovery within ± 10 % of initial value at a lactate dehydrogenase activity of 480 U/L (8.0 μkat/L).

**Icterus:**<sup>9</sup> No significant interference up to an I index of 60 for conjugated and unconjugated bilirubin (approximate conjugated and unconjugated bilirubin concentration: 1026 μmol/L or 60 mg/dL).

**Hemolysis:**<sup>9</sup> No significant interference up to an H index of 15 (approximate hemoglobin concentration: 15 mg/dL or 9.6 μmol/L).

Contamination with erythrocytes will elevate results, because the analyte level in erythrocytes is higher than in normal sera. The level of interference may be variable depending on the content of analyte in the lysed erythrocytes.

**Lipemia (Intralipid):**<sup>9</sup> No significant interference up to an L index of 900. There is poor correlation between the L index (corresponds to turbidity) and triglycerides concentration.

**Drugs:** No interference was found at therapeutic concentrations using common drug panels.<sup>10,11</sup>

In very rare cases, gammopathy, in particular type IgM (Waldenström's macroglobulinemia), may cause unreliable results.<sup>12</sup>

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

**ACTION REQUIRED**

**Special Wash Programming:** The use of special wash steps is mandatory when certain test combinations are run together on Roche/Hitachi **cobas c** systems. The latest version of the carry-over evasion list can be found with the NaOHD/SMS/Multiclean/SCCS or the NaOHD/SMS/SmpCln1+2/SCCS Method Sheets. For further instructions refer to the operator's manual. **cobas c 502** analyzer: All special wash programming necessary for avoiding carry-over is available via the **cobas** link, manual input is not required.

**Where required, special wash/carry-over evasion programming must be implemented prior to reporting results with this test.**

**Limits and ranges****Measuring range**

20-1200 U/L (0.33-20 μkat/L)

Determine samples having higher activities via the rerun function. Dilution of samples via the rerun function is a 1:2.5 dilution. Results from samples diluted using the rerun function are automatically multiplied by a factor of 2.5.



**Lower limits of measurement***Lower detection limit of the test*

20 U/L (0.33 µkat/L)

The lower detection limit represents the lowest measurable analyte level that can be distinguished from zero. It is calculated as the value lying 3 standard deviations above that of the lowest standard (standard 1 + 3 SD, repeatability, n = 21).

**Expected values**Expected values in adults<sup>13</sup> at 37 °C (calculated):

240-480 U/L (4.00-8.00 µkat/L)

A factor of 2.00 was used for converting the reference range from 25 to 37 °C.

Reference ranges for children are given in the brochure "Reference ranges for adults and children; preanalytical considerations" by Heil W, Koberstein R, Zawta B (published by Roche Diagnostics GmbH, 2004).

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

**Specific performance data**

Representative performance data on the analyzers are given below. Results obtained in individual laboratories may differ.

**Precision**

Precision was determined using human samples and controls in an internal protocol with repeatability (n = 21) and intermediate precision (3 aliquots per run, 1 run per day, 21 days). The following results were obtained:

<i>Repeatability</i>	<i>Mean</i>	<i>SD</i>	<i>CV</i>
	<i>U/L (µkat/L)</i>	<i>U/L (µkat/L)</i>	<i>%</i>
Precinorm U	316 (5.28)	2 (0.04)	0.8
Precipath U	513 (8.57)	3 (0.05)	0.6
Human serum 1	281 (4.69)	3 (0.05)	1.1
Human serum 2	731 (12.2)	5 (0.08)	0.7
<i>Intermediate precision</i>	<i>Mean</i>	<i>SD</i>	<i>CV</i>
	<i>U/L (µkat/L)</i>	<i>U/L (µkat/L)</i>	<i>%</i>
Precinorm U	312 (5.22)	3 (0.05)	0.9
Precipath U	502 (8.38)	5 (0.08)	0.9
Human serum 3	239 (3.99)	6 (0.10)	2.6
Human serum 4	632 (10.6)	8 (0.13)	1.2

**Method comparison**

LDH values for human serum and plasma samples obtained on a Roche/Hitachi **cobas c 501** analyzer (y) were compared with those determined using the same reagent on a Roche/Hitachi 917 analyzer (x).

Sample size (n) = 79

Passing/Bablok <sup>14</sup>	Linear regression
$y = 0.995x + 6.03 \text{ U/L}$	$y = 0.991x + 6.36 \text{ U/L}$
$r = 0.971$	$r = 0.999$

The sample activities were between 223 and 1146 U/L (3.72 and 19.1 µkat/L).

**References**

- 1 Dito WR. Lactate dehydrogenase: A brief review. In: Griffiths JC, ed. Clinical Enzymology. New York: Masson Publishing USA 1979:1-8.
- 2 Moss DW, Henderson AR, Kachmar JF. Enzymes. In: Tietz NW, ed. Fundamentals of Clinical Chemistry, 3rd ed. Philadelphia, PA: WB Saunders 1987:346-421.
- 3 Zimmerman HJ, Henry JB In: Henry JB, ed. Clinical Diagnosis and Management by Laboratory Methods. 17th ed. Philadelphia, PA: WB Saunders 1984:251-282.
- 4 Empfehlungen der Deutschen Gesellschaft für Klinische Chemie. Z klin Chem u klin Biochem 1970;8:658-659.

- 5 Empfehlungen der Deutschen Gesellschaft für Klinische Chemie. Z klin Chem u klin Biochem 1972;10:182-190.
- 6 Bais R, Philcox M. Approved recommendations of IFCC methods for the measurement of catalytic concentration of enzymes. Part 8. IFCC method for lactate dehydrogenase. International Federation of Clinical Chemistry (IFCC). Eur J Clin Chem Clin Biochem 1994;32(8):639-655.
- 7 Burtis CA, Ashwood ER. Tietz Textbook of Clinical Chemistry, 3rd ed. Pa: WB Saunders Co 1999:669.
- 8 Use of Anticoagulants in Diagnostic Laboratory Investigations. WHO Publication WHO/DIL/LAB/99.1 Rev. 2. Jan. 2002.
- 9 Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-475.
- 10 Breuer J. Report on the Symposium "Drug effects in Clinical Chemistry Methods". Eur J Clin Chem Clin Biochem 1996;34:385-386.
- 11 Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385.
- 12 Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med 2007;45(9):1240-1243.
- 13 Weißhaar D, Gossrau E, Faderl B. Normalbereiche von α-HBDH, LDH, AP und LAP bei Messung mit substrat-optimierten Testansätzen. Med Welt 1975;26:387-390.
- 14 Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

**Symbols**

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard.

**CONTENT**

Contents of kit



Volume after reconstitution or mixing

**FOR US CUSTOMERS ONLY: LIMITED WARRANTY**

Roche Diagnostics warrants that this product will meet the specifications stated in the labeling when used in accordance with such labeling and will be free from defects in material and workmanship until the expiration date printed on the label. THIS LIMITED WARRANTY IS IN LIEU OF ANY OTHER WARRANTY, EXPRESS OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. IN NO EVENT SHALL ROCHE DIAGNOSTICS BE LIABLE FOR INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES.

COBAS, COBAS C, PRECINORM, PRECIPATH and PRECICONTROL are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Significant additions or changes are indicated by a change bar in the margin.

© 2013, Roche Diagnostics



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim  
www.roche.com

Distribution in USA by:  
Roche Diagnostics, Indianapolis, IN  
US Customer Technical Support 1-800-428-2336

