

Phencyclidine Plus

Order information

REF	CONTENT	Analyzer(s) on which cobas c pack(s) can be used
04490908 190	ONLINE DAT Phencyclidine Plus (200 tests) System-ID 07 6919 3	Roche/Hitachi cobas c 311, cobas c 501/502
03304671 190	Preciset DAT Plus I calibrators CAL 1-6 (6 x 5 mL) Codes 431-436	
03304698 190	C.f.a.s. DAT Qualitative Plus (6 x 5 mL)	
04590856 190	C.f.a.s. DAT Qualitative Plus Clinical (3 x 5 mL) Code 699	
03312950 190	Control Set DAT I PreciPos DAT Set I (2 x 10 mL) PreciNeg DAT Set I (2 x 10 mL)	
04500873 190	Control Set DAT Clinical PreciPos DAT Clinical (2 x 10 mL) PreciNeg DAT Clinical (2 x 10 mL)	

English

System information

For **cobas c** 311/501 analyzers:

PC2QP: ACN 518: for qualitative assay

PC2SP: ACN 519: for semiquantitative assay

PC2QC: ACN 795: for qualitative assay; using C.f.a.s. DAT Qualitative Plus Clinical

For **cobas c** 502 analyzer:

PC2QP: ACN 8518: for qualitative assay

PC2SP: ACN 8519: for semiquantitative assay

PC2QC: ACN 8795: for qualitative assay; using C.f.a.s. DAT Qualitative Plus Clinical

Intended use

Phencyclidine Plus (PCP) is an in vitro diagnostic test for the qualitative and semiquantitative detection of phencyclidine and its metabolites in human urine on Roche/Hitachi **cobas c** systems at a cutoff concentration of 25 ng/mL. Semiquantitative test results may be obtained that permit laboratories to assess assay performance as part of a quality control program. Semiquantitative assays are intended to determine an appropriate dilution of the specimen for confirmation by a confirmatory method such as gas chromatography/mass spectrometry (GC-MS).

Phencyclidine Plus provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. GC-MS is the preferred confirmatory method.¹ Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used.

Summary

Phencyclidine (PCP) is an arylcyclohexylamine with potent analgesic and anesthetic properties.^{1,2,3,4,5,6} Originally developed as an intravenous anesthetic, the occurrence of emergence psychosis side effects negated its potential clinical utility. PCP was never approved for human use because of the post-anesthetic confusion and delirium that arose during clinical studies. Illegally sold on the street, PCP is known by various names such as "angel dust"; whereas, names such as "supergrass" refer to PCP combined with marijuana. PCP possesses hallucinogenic, central nervous system (CNS)-stimulant, and CNS-depressant properties, the expression of which is dose- and species-dependent.⁴ PCP and its structural analog, ketamine, are NMDA (N-methyl-D-aspartate) receptor antagonists.^{2,5} Known as dissociative anesthetics, they produce mind-altering feelings of dissociation from the environment and self. Dextromethorphan, a cough suppressant, can produce similar effects when taken in high doses.

The water-soluble powder of PCP can be ingested, snorted, injected intravenously, or smoked. Typical street doses (1-10 mg) can cause tachycardia, hypertension, hallucinations, stupor, lethargy, sensory isolation, and loss of coordination. Excitation and agitation may also occur, leading to unpredictably violent behavior not usually encountered with other hallucinogens. Repeated use of PCP can result in addiction and higher doses can cause symptoms that mimic schizophrenia and can culminate in convulsions and prolonged or fatal coma.^{2,6}

PCP is metabolized via ring-hydroxylation and oxidation by the cytochrome P450 enzymes.^{3,7} An amino acid metabolite of PCP exists in human urine in

significant quantities.⁸ Significant variations in the PCP elimination half-life have been found in humans; however, phase II metabolism of PCP sulfation and glucuronidation could also contribute to the variation in PCP half-life.⁷

Test principle

The assay is based on the kinetic interaction of microparticles in a solution (KIMS)⁹ as measured by changes in light transmission. In the absence of sample drug, free antibody binds to drug-microparticle conjugates causing the formation of particle aggregates. As the aggregation reaction proceeds in the absence of sample drug, the absorbance increases.

When a urine sample contains the drug in question, this drug competes with the particle-bound drug derivative for free antibody. Antibody bound to sample drug is no longer available to promote particle aggregation, and subsequent particle lattice formation is inhibited. The presence of sample drug diminishes the increasing absorbance in proportion to the concentration of drug in the sample. Sample drug content is determined relative to the value obtained for a known cutoff concentration of drug.

Reagents - working solutions

- R1** Buffer; 0.09 % sodium azide
- R2** PCP antibody (mouse monoclonal); buffer; bovine serum albumin; 0.09 % sodium azide
- R3** Conjugated PCP derivative microparticles; buffer; 0.09 % sodium azide

R1 is in position B, R2 is in position C, and R3 is in position A.

Precautions and warnings

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents.

Disposal of all waste material should be in accordance with local guidelines. Safety data sheet available for professional user on request.

For USA: Caution: Federal law restricts this device to sale by or on the order of a physician.

Reagent handling

Ready for use

Carefully invert reagent container several times prior to use to ensure that the reagent components are mixed.

Storage and stability

Shelf life at 2-8 °C: See expiration date on **cobas c** pack label

On-board in use and refrigerated on the analyzer: 8 weeks

Do not freeze.

Specimen collection and preparation

Only the specimens listed below were tested and found acceptable.

Urine: Collect urine samples in clean glass or plastic containers. Fresh urine specimens do not require any special handling or pretreatment, but an effort should be made to keep pipetted samples free of gross debris.

Phencyclidine Plus

Samples should be within the normal physiological pH range of 5-8. No additives or preservatives are required. It is recommended that urine specimens be stored at 2-8 °C and tested within 5 days of collection.¹⁰

For prolonged storage, freezing of samples is recommended.

Centrifuge highly turbid specimens before testing.

Adulteration or dilution of the sample can cause erroneous results. If adulteration is suspected, another sample should be collected. Specimen validity testing is required for specimens collected under the *Mandatory Guidelines for Federal Workplace Drug Testing Programs*.¹¹

CAUTION: Specimen dilutions should only be used to interpret results of Calc.? and Samp.? alarms, or when estimating concentration in preparation for GC-MS. Dilution results are not intended for patient values. Dilution procedures, when used, should be validated.

Materials provided

See "Reagents – working solutions" section for reagents.

Materials required (but not provided)

See "Order information" section

General laboratory equipment

Assay

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

The performance of applications not validated by Roche is not warranted and must be defined by the user.

Application for urine

Deselect Automatic Rerun for these applications in the Utility menu, Application screen, Range tab.

cobas c 311 test definition

	Semiquantitative	Qualitative	
Assay type	2-Point End	2-Point End	
Reaction time / Assay points	10 / 27-51	10 / 27-51	
Wavelength (sub/main)	– /505 nm	– /505 nm	
Reaction direction	Increase	Increase	
Unit	ng/mL	mAbs	
Reagent pipetting		Diluent (H ₂ O)	
R1	59 µL	–	
R2	59 µL	–	
R3	54 µL	–	
<i>Sample volumes</i>	<i>Sample</i>	<i>Sample dilution</i>	
		Sample	Diluent NaCl
Normal	11.3 µL	–	–
Decreased	11.3 µL	–	–
Increased	11.3 µL	–	–

cobas c 501/502 test definition

	Semiquantitative	Qualitative	
Assay type	2-Point End	2-Point End	
Reaction time / Assay points	10 / 40-58	10 / 40-58	
Wavelength (sub/main)	– /505 nm	– /505 nm	
Reaction direction	Increase	Increase	
Unit	ng/mL	mAbs	
Reagent pipetting		Diluent (H ₂ O)	
R1	59 µL	–	
R2	59 µL	–	
R3	54 µL	–	

<i>Sample volumes</i>	<i>Sample</i>	<i>Sample dilution</i>	
		Sample	Diluent NaCl
Normal	11.3 µL	–	–
Decreased	11.3 µL	–	–
Increased	11.3 µL	–	–

Calibration

Calibrators *Semiquantitative application*

S1-4: Preciset DAT Plus I calibrators, CAL 1-4

0, 12.5, 25, 50 ng/mL

Qualitative application

S1: C.f.a.s. DAT Qualitative Plus, C.f.a.s. DAT Qualitative Plus Clinical, or Preciset DAT Plus I calibrator - CAL 3 25 ng/mL

The drug concentrations of the calibrators have been verified by GC-MS.

Calibration K Factor For the qualitative application, enter the K Factor as -1000 into the Calibration menu, Status screen, Calibration Result window.

Calibration mode *Semiquantitative application*
Result Calculation Mode (RCM)^a

Qualitative application
Linear

Calibration frequency Full (semiquantitative) or blank (qualitative) calibration
- after reagent lot change
- as required following quality control procedures

a) See Results section.

Calibration interval may be extended based on acceptable verification of calibration by the laboratory.

Traceability: This method has been standardized against a primary reference method (GC-MS).

Quality control

For quality control, use control materials as listed in the "Order information" section.

In addition, other suitable control material can be used.

Drug concentrations of the Control Set DAT I and Clinical have been verified by GC-MS.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

Results

For the qualitative assay, the cutoff calibrator is used as a reference in distinguishing between preliminary positive and negative samples. Samples producing a positive or "0" absorbance value are considered preliminary positive. Preliminary positive samples are flagged with >Test. Samples producing a negative absorbance value are considered negative. Negative samples are preceded by a minus sign.

The semiquantitation of preliminary positive results should only be used by laboratories to determine an appropriate dilution of the specimen for confirmation by a confirmatory method such as GC-MS. It also permits the laboratory to establish quality control procedures and assess control performance.

For the semiquantitative assay, the analyzer computer constructs a calibration curve from absorbance measurements of the standards using a 4 parameter logit-log fitting function (RCM). The logit-log function fits a smooth line through the data points. The analyzer computer uses

absorbance measurements of samples to calculate drug or drug metabolite concentration by interpolation of the logit-log fitting function.

NOTE: If a result of Calc.? or Samp.? alarm is obtained, review the Reaction Monitor data for the sample and compare with the Reaction Monitor data for the highest calibrator. The most likely cause is a high concentration of the analyte in the sample, in which case the absorbance value for the sample will be less than that of the highest calibrator. Make an appropriate dilution of the sample using the 0 ng/mL calibrator and rerun the sample. A normal drug-free urine may be substituted for the 0 ng/mL calibrator if the urine and procedure have been validated by the laboratory. To ensure that the sample was not over-diluted, the diluted result, prior to multiplying by the dilution factor, must be at least half the analyte cutoff value. If the diluted result falls below half the analyte cutoff value, repeat the sample with a smaller dilution. A dilution that produces a result closest to the analyte cutoff is the most accurate estimation. To estimate the preliminary positive sample's concentration, multiply the result by the appropriate dilution factor. Dilutions should only be used to interpret results of Calc.? or Samp.? alarms, or when estimating concentration in preparation for GC-MS.

Use caution when reporting results as there are various factors that influence a urine test result, such as fluid intake and other biological factors.

As with any *sensitive* test for drugs of abuse on automated clinical chemistry analyzers, the possibility exists for analyte carry-over from a sample with an extremely high concentration to a normal (negative) sample which immediately follows it.

Confirm all preliminary positive results by another method.

Limitations - interference

See the "Specific performance data" section of this document for information on substances tested with this assay. There is the possibility that other substances and/or factors may interfere with the test and cause erroneous results (e.g., technical or procedural errors).

A preliminary positive result with this assay indicates the presence of PCP and/or its metabolites in urine. It does not measure the level of intoxication.

Interfering substances were added to drug free urine at the concentration listed below. These samples were then spiked to 25 ng/mL using a PCP stock solution. Samples were tested in triplicate (n = 3) on a Roche/Hitachi **cobas c 501** analyzer. The median % recoveries were calculated and are listed below.

Substance	Concentration Tested	% Phencyclidine Recovery
Acetone	1 %	98
Ascorbic acid	1.5 %	105
Bilirubin	0.25 mg/mL	98
Creatinine	5 mg/mL	113
Ethanol	1 %	100
Glucose	2 %	105
Hemoglobin	7.5 g/L	94
Human albumin	0.5 %	102
Oxalic acid	2 mg/mL	98
Sodium chloride	0.5 M	100
Sodium chloride	1 M	102
Urea	6 %	106

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

ACTION REQUIRED

Special Wash Programming: The use of special wash steps is mandatory when certain test combinations are run together on Roche/Hitachi **cobas c** systems. The latest version of the carry-over evasion list can be found with the NaOHD-SMS-SmpCln1+2-SCCS Method Sheets. For further instructions refer to the operator's manual. **cobas c 502** analyzer: All special wash programming necessary for avoiding carry-over is available via the **cobas** link, manual input is not required.

Where required, special wash/carry-over evasion programming must be implemented prior to reporting results with this test.

Expected values

Qualitative assay

Results of this assay distinguish preliminary positive (≥ 25 ng/mL) from negative samples only. The amount of drug detected in a preliminary positive sample cannot be estimated.

Semiquantitative assay

Results of this assay yield only approximate cumulative concentrations of the drug and its metabolites (see "Analytical specificity" section).

Specific performance data

Representative performance data on a Roche/Hitachi analyzer are given below. Results obtained in individual laboratories may differ.

Precision

Precision was determined in an internal protocol by running a series of calibrator and controls (repeatability n = 20, intermediate precision n = 100). The following results were obtained on a Roche/Hitachi **cobas c 501** analyzer.

Semiquantitative precision

Repeatability	Mean ng/mL	SD ng/mL	CV %
Level 1	18.0	0.6	3.6
Level 2	25.1	0.7	2.9
Level 3	30.6	0.6	1.9

Intermediate precision	Mean ng/mL	SD ng/mL	CV %
Level 1	18.1	0.8	4.3
Level 2	24.6	0.8	3.1
Level 3	31.2	0.7	2.2

Qualitative precision

Cutoff (25)	Number tested	Correct results	Confidence level
0.75x	100	100	> 95 % negative reading
1.25x	100	100	> 95 % positive reading

Lower detection limit of the test

1.6 ng/mL

The lower detection limit represents the lowest measurable analyte level that can be distinguished from zero. It is calculated as the value lying 2 standard deviations above that of the lowest standard (standard 1 + 2 SD, repeatability, n = 21).

Accuracy

100 urine samples, obtained from a clinical laboratory where they screened negative in a drug test panel, were evaluated with the Phencyclidine Plus assay. 100 % of these normal urines were negative relative to a 25 ng/mL cutoff. 65 samples obtained from a clinical laboratory, where they screened preliminary positive with a commercially available immunoassay and were subsequently confirmed by GC-MS, were evaluated with the Phencyclidine Plus assay. 99 % of these samples were positive relative to a 25 ng/mL cutoff. In addition, 9 samples with GC-MS values approximately 50-100 % of the cutoff were evaluated with the Phencyclidine Plus assay. Data from the accuracy studies described above were combined with data generated from these samples. The following results were obtained with the Phencyclidine Plus assay on the Roche/Hitachi 917 analyzer relative to the GC-MS values.

Phencyclidine Plus Clinical Correlation (Cutoff = 25 ng/mL)				
		Negative Samples	GC-MS values (ng/mL)	
			Near Cutoff	
			12-23	25-32
				34->1000

Phencyclidine Plus

Phencyclidine Plus Clinical Correlation (Cutoff = 25 ng/mL)					
Roche/Hitachi 917 analyzer	+	0	4	10	54
	-	100	5	1	0

Additional clinical samples were evaluated with this assay on a Roche/Hitachi **cobas c** 501 analyzer and a Roche/Hitachi 917 analyzer. 100 urine samples, obtained from a clinical laboratory where they screened negative in a drug test panel, were evaluated with the Phencyclidine Plus assay. 100 % of these normal urines were negative relative to the Roche/Hitachi 917 analyzer. 54 urine samples, obtained from a clinical laboratory where they screened preliminary positive with a commercially available immunoassay and were subsequently confirmed by GC-MS, were evaluated with the Phencyclidine Plus assay. 100 % of the samples were positive on both the Roche/Hitachi **cobas c** 501 analyzer and the Roche/Hitachi 917 analyzer.

Phencyclidine Plus Correlation (Cutoff = 25 ng/mL)			
		Roche/Hitachi 917 analyzer	
		+	-
cobas c 501 analyzer	+	54	0
	-	0	100

Analytical specificity

The specificity of this assay for structurally similar compounds was determined by generating inhibition curves for each of the compounds listed and determining the approximate quantity of each compound that is equivalent in assay reactivity to a 25 ng/mL phencyclidine assay cutoff. The following results were obtained on a Roche/Hitachi 917 analyzer.

Compound	ng/mL Equivalent to 25 ng/mL Phencyclidine	Approximate % Cross-reactivity
Thienylcyclohexylpiperidine (TCP)	49	51.14
Dextromethorphan	> 100000	0.01
Ketamine	> 100000	0.00

Drug interference

The following compounds were prepared in aliquots of pooled normal human urine to yield a final concentration of 100000 ng/mL. None of these compounds gave values in the assay that were greater than 0.018 % cross-reactivity.

Acetaminophen	Lidocaine
Acetylsalicylic acid	LSD
Aminopyrine	MDA
Amobarbital	MDMA
<i>d</i> -Amphetamine	Melanin
<i>l</i> -Amphetamine	Meperidine
Ampicillin	Methadone
Ascorbic acid	<i>d</i> -Methamphetamine
Aspartame	<i>l</i> -Methamphetamine
Atropine	Methaqualone
Benzocaine	Methylphenidate
Benzoylcegonine (cocaine metabolite)	Methypylon
Benzphetamine	Morphine
Butabarbital	Naloxone
Caffeine	Naltrexone
Calcium hypochlorite	Naproxen
	Niacinamide

Chlordiazepoxide	Norethindrone
Chloroquine	<i>l</i> -Norpseudoephedrine
Chlorpheniramine	Nortriptyline
Chlorpromazine	Oxazepam
Cocaine	Penicillin G
Codeine	Pentobarbital
Dextropropoxyphene	β -Phenethylamine
Diazepam	Phenobarbital
Diphenhydramine	Phenothiazine
Dopamine	Phentermine
Doxepin	Phenylbutazone
Ecgonine	<i>d</i> -Phenylpropanolamine
Ecgonine methyl ester	<i>d,l</i> -Phenylpropanolamine
<i>d</i> -Ephedrine	Procaine
<i>d,l</i> -Ephedrine	Promethazine
<i>l</i> -Ephedrine	<i>d</i> -Pseudoephedrine
Epinephrine	<i>l</i> -Pseudoephedrine
Erythromycin	Quinidine
Estriol	Quinine
Fenoprofen	Secobarbital
Furosemide	Sulindac
Gentisic acid	Tetracycline
Glutethimide	Δ^9 THC-9-carboxylic acid
Guaiacol glycerol ether	Tetrahydrozoline
Hydrochlorothiazide	Trifluoperazine
<i>p</i> -Hydroxyamphetamine	Trimipramine
Ibuprofen	Tyramine
Isoproterenol	Verapamil

The cross-reactivity for Amitriptyline, Desipramine, and Imipramine were tested at a concentration of 100000 ng/mL with the Phencyclidine Plus assay. The results obtained were 0.031 %, 0.022 %, and 0.037 %, respectively.

References




- Karch SB, ed. Drug Abuse Handbook. Boca Raton, FL: CRC Press LLC 1998.
- Hardman JG, Limbird LE, Gilman A, eds. Goodman and Gilman's The Pharmacological Basis of Therapeutics. 10th ed. New York, NY: McGraw Hill Pub Co. 2001.
- Baselt RC. Disposition of Toxic Drugs and Chemicals in Man. 7th ed. Foster City, CA: Biomedical Publications 2004.
- Clouet DH, ed. Phencyclidine: An Update. NIDA Research Monograph 64. National Institute on Drug Abuse 1986.
- De Souza EB, Clouet D, London ED, eds. Sigma, PCP, and NMDA receptors. Research Monograph 133. National Institute on Drug Abuse 1993.
- Leshner AI. Hallucinogens and dissociative, including LSD, PCP, Ketamine, dextromethorphan. NIH Publication Number 01-4209, NIDA Research Report Series 2001.
- Laurenzana EM, Owens SM. Metabolism of phencyclidine by human liver microsomes. Drug Metab Dispos 1997;25:557-663.
- ElSohly MA, Little TL Jr, Mitchell JM, et al. GC/MS analysis of phencyclidine acid metabolite in human urine. J Anal Toxicol 1988;12:180-182.
- Armbruster DA, Schwarzhoff RH, Pierce BL, et al. Method comparison of EMIT II and ONLINE with RIA for drug screening. J Forensic Sci 1993;38:1326-1341.

- 10 Toxicology and Drug Testing in the Clinical Laboratory; Approved Guideline. 2nd ed. (C52-A2). Clinical and Laboratory Standards Institute 2007;27:33.
- 11 Mandatory Guidelines for Federal Workplace Drug Testing Programs. Fed Regist 2008 Nov 25;73:71858-71907.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard (for USA: see <https://usdiagnostics.roche.com> for definition of symbols used):

	Contents of kit
	Volume after reconstitution or mixing
	Global Trade Item Number

FOR US CUSTOMERS ONLY: LIMITED WARRANTY

Roche Diagnostics warrants that this product will meet the specifications stated in the labeling when used in accordance with such labeling and will be free from defects in material and workmanship until the expiration date printed on the label. THIS LIMITED WARRANTY IS IN LIEU OF ANY OTHER WARRANTY, EXPRESS OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. IN NO EVENT SHALL ROCHE DIAGNOSTICS BE LIABLE FOR INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES.

COBAS, COBAS C, ONLINE DAT, and PRECISET are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2017, Roche Diagnostics



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim
www.roche.com



Distribution in USA by:
Roche Diagnostics, Indianapolis, IN
US Customer Technical Support 1-800-428-2336