



05994454001V5.0

SALINE

NaCl Solution

cobas®

Order information

REF	CONTENT	Analyzer(s) on which kit(s) can be used
11557335 216	NaCl solution (12 x 66 mL)	Roche/Hitachi MODULAR analyzers

Some analyzers and kits shown may not be available in all countries. For additional system applications, contact your local Roche Diagnostics representative.

English

System information

Diluent: ACN 951

Serum indices: ACN 850

Intended use

NaCl solution is used as a diluent for automatic and manual dilution of calibrators, patient samples and can also be used for measuring the serum indices on Roche/Hitachi systems.

Summary

NaCl solution is an aqueous preparation used for predilution or dilution of samples with values exceeding the measuring range of assays on Roche/Hitachi systems. When performing serum indices using NaCl solution refer to the appropriate operators manual for detailed instructions. For the dilution of calibrators refer to the instructions for use for the respective Roche/Hitachi assay reagents.

Reagents - working solutions

Sodium chloride 0.9 % (154 mmol/L)

Precautions and warnings

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents.

Disposal of all waste material should be in accordance with local guidelines. Safety data sheet available for professional user on request.

For USA: For prescription use only.

Reagent handling

Ready for use

Storage and stability

Shelf life at 2-8 °C: See expiration date on reagent label.

On-board in use and refrigerated on the analyzer: 12 weeks

Materials provided

See "Reagents – working solutions" section for reagents.

12 additional barcode labels for serum index measurements (ACN 850, bottle sizes 70 mL and 130 mL).

Materials required (but not provided)

See "Order information" section

Roche/Hitachi analyzers and assay reagents.

See appropriate instructions for use and operator's manuals for additional required materials.

Assay

Use NaCl solution as specified in the respective instructions for use of the system reagents.

For further information, please refer to the appropriate operator's manual for the analyzer concerned, the respective application sheets and method sheets of all necessary components.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard:

CONTENT	Contents of kit
REAGENT	Reagent

CALIBRATOR



GTIN

Calibrator

Volume after reconstitution or mixing

Global Trade Item Number

FOR US CUSTOMERS ONLY: LIMITED WARRANTY

Roche Diagnostics warrants that this product will meet the specifications stated in the labeling when used in accordance with such labeling and will be free from defects in material and workmanship until the expiration date printed on the label. THIS LIMITED WARRANTY IS IN LIEU OF ANY OTHER WARRANTY, EXPRESS OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. IN NO EVENT SHALL ROCHE DIAGNOSTICS BE LIABLE FOR INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES.

COBAS is a trademark of Roche.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2016, Roche Diagnostics



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim
www.roche.com

Distribution in USA by:
Roche Diagnostics, Indianapolis, IN
US Customer Technical Support 1-800-428-2336

