

# TOBR2



## Tobramycin

### Order information

REF	CONTENT	Analyzer(s) on which <b>cobas c</b> pack(s) can be used
04491033 190	ONLINE TDM Tobramycin 100 tests	System-ID 07 6985 1 Roche/Hitachi <b>cobas c</b> 311, <b>cobas c</b> 501/502
03375790 190	Preciset TDM I Calibrators CAL A-F (6 x 5 mL) Diluent (1 x 10 mL)	Codes 691-696
04521536 190	TDM Control Set Level I (2 x 5 mL) Level II (2 x 5 mL) Level III (2 x 5 mL)	Code 310 Code 311 Code 312
04708725 190	Sample Cleaner 1 (59 mL)	

### English

#### System information

For **cobas c** 311/501 analyzers:

**TOBR2**: ACN 607

For **cobas c** 502 analyzers:

**TOBR2**: ACN 8607

#### Intended use

In vitro test for the quantitative determination of tobramycin in serum and plasma on Roche/Hitachi **cobas c** systems.

#### Summary

Tobramycin is an aminoglycoside antibiotic used in the treatment of infections caused by *Pseudomonas aeruginosa*, *Proteus species*, *E. coli*, *Klebsiella*, *Serratia*, *Citrobacter*, *Staphylococcus aureus*, *Enterobacter* and other microorganisms. Tobramycin's toxic effect is produced by interfering with ribosomal protein synthesis.<sup>1</sup> Tobramycin undergoes very little, if any, metabolism and is, therefore, eliminated as the parent drug by glomerular filtration. The half-life of tobramycin in serum or plasma correlates closely with renal function and thus is quite variable between individuals and within one individual over time.<sup>2,3</sup> Serum or plasma tobramycin concentration is also impacted by mode of administration, the volume of extracellular fluid, the duration of the treatment and physiological changes during the illness and therapy. The therapeutic range of tobramycin should be measured at peak as well as trough concentrations. In patients with pre-existing renal damage or those to whom tobramycin has been administered for prolonged periods or in doses above the therapeutic range, hearing impairment and/or nephrotoxicity may develop. Therefore, monitoring of peak and trough tobramycin serum or plasma levels is critical in the prevention of these serious complications with the adjustment of dosage administration as indicated.<sup>4</sup>

#### Test principle

The assay is based on a homogeneous enzyme immunoassay technique used for the quantitative analysis of tobramycin in human serum or plasma.<sup>5</sup> The assay is based on competition between drug in the sample and drug labeled with the enzyme glucose-6-phosphate dehydrogenase (G6PDH) for antibody binding sites. Enzyme activity decreases upon binding to the antibody, so the drug concentration in the sample can be measured in terms of enzyme activity. Active enzyme converts oxidized nicotinamide adenine dinucleotide (NAD) to NADH, resulting in an absorbance change that is measured spectrophotometrically. Endogenous serum G6PDH does not interfere because the coenzyme functions only with the bacterial (*Leuconostoc mesenteroides*) enzyme employed in the assay.

#### Reagents - working solutions

- R2** Tobramycin labeled with bacterial G6PDH and bovine serum albumin in buffer
- R3** Anti-tobramycin antibody (sheep polyclonal), G6P, NAD and bovine serum albumin in buffer

R2 is in position A and R3 is in position B.

#### Precautions and warnings

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents.

Disposal of all waste material should be in accordance with local guidelines. Safety data sheet available for professional user on request.

For USA: For prescription use only.

This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:

2-methyl-2H-isothiazol-3-one hydrochloride.

EUH 208 May produce an allergic reaction.

Product safety labeling primarily follows EU GHS guidance.

#### Reagent handling

Ready for use

Mix reagents by gentle inversion numerous times before placing on-board the analyzer.

#### Storage and stability

Shelf life at 2-8 °C : See expiration date on **cobas c** pack label

On-board in use and refrigerated on the analyzer : 12 weeks

#### Do not freeze.

#### Specimen collection and preparation

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable.

Serum: Collect serum using standard sampling tubes.

Plasma: K<sub>2</sub>- or K<sub>3</sub>-EDTA, sodium citrate, fluoride oxalate, or sodium or lithium heparinized plasma.

Stability:<sup>6</sup> 3 days capped at 4-8 °C  
1 month capped at -20 °C

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

Do not induce foaming of specimens. Specimens should not be repeatedly frozen and thawed.

Invert thawed specimens several times prior to testing.

#### Materials provided

See "Reagents – working solutions" section for reagents.

#### Materials required (but not provided)

See "Order information" section

General laboratory equipment

#### Assay

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

The performance of applications not validated by Roche is not warranted and must be defined by the user.

# TOBR2

## Tobramycin



### Application for serum and plasma

Deselect Automatic Rerun for these applications in the Utility menu, Application screen, Range tab.

### cobas c 311 test definition

Assay type	Rate-A assay		
Reaction time /Assay points:	10 / 28-40		
Wavelength (sub/main)	415 /340 nm		
Reaction direction	Increase		
Unit	µg/mL		
Reagent pipetting	Diluent (H <sub>2</sub> O)		
R2	145 µL	–	
R3	73 µL	–	
Sample volumes	Sample	Sample dilution	
		Sample	Diluent (NaCl)
Normal	2.4 µL	–	–
Decreased	2.4 µL	–	–
Increased	2.4 µL	–	–

### cobas c 501/502 test definition

Assay type	Rate-A assay		
Reaction time /Assay points:	10 / 42-54		
Wavelength (sub/main)	415 /340 nm		
Reaction direction	Increase		
Unit	µg/mL		
Reagent pipetting	Diluent (H <sub>2</sub> O)		
R2	145 µL	–	
R3	73 µL	–	
Sample volumes	Sample	Sample dilution	
		Sample	Diluent (NaCl)
Normal	2.4 µL	–	–
Decreased	2.4 µL	–	–
Increased	2.4 µL	–	–

### Calibration

Calibrator	S1-6: Preciset TDM I Calibrators
Calibration mode	RCM
Calibration frequency	6-point calibration <ul style="list-style-type: none"> <li>• after reagent lot change</li> <li>• every 6 weeks</li> <li>• as required following quality control procedures</li> </ul>

Traceability: This method has been standardized against USP reference standards. The calibrators are prepared to contain known quantities of tobramycin in normal human serum.

### ACTION REQUIRED

**NOTE:** Due to potential carryover from the last calibrator (Cal F) into the first quality control sample (Level 1) following calibration, assaying a non-reportable blank quality control sample is required prior to assaying the controls. The blank quality control sample should be programmed in the first position followed by quality control levels 1-3. Use Sample Cleaner 1 (Cat. no. 04708725 190) as the blank quality control sample.

**The blank quality control sample is not required when assaying controls without calibration.**

### Quality control

For quality control, use control materials as listed in the "Order information" section.

In addition, other suitable control material can be used.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

### Calculation

Roche/Hitachi **cobas c** systems automatically calculate the analyte concentration of each sample.

Conversion factor:<sup>7</sup> µg/mL x 2.14 = µmol/L

### Limitations - interference

Criterion: Recovery within ± 10 % of initial value at tobramycin levels of approximately 3.5 and 8 µg/mL (7.5 and 17 µmol/L).

#### Serum/Plasma

Icterus:<sup>8</sup> No significant interference up to an I index of 30 for conjugated bilirubin and unconjugated bilirubin (approximate conjugated and unconjugated bilirubin concentration: 30 mg/dL or 513 µmol/L).

Lipemia (Intralipid):<sup>8</sup> No significant interference up to an L index of 750. There is poor correlation between the L index (corresponds to turbidity) and triglycerides concentration.

Hemolysis:<sup>8</sup> No significant interference up to an H index of 800 (approximate hemoglobin concentration: 800 mg/dL or 497 µmol/L).

Criterion: Recovery within ± 10 % of initial value at a tobramycin level of approximately 3 µg/mL (6.4 µmol/L).

No significant interference from triglycerides up to 750 mg/dL (8.5 mmol/L).

Rheumatoid factors: No interference from rheumatoid factors up to 100 IU/mL.

Total protein: No interference from protein from 2-12 g/dL.

Amikacin cross-reacts with this assay. Kanamycin cross-reacts significantly; however, the assay has not been optimized to quantitate this aminoglycoside. Aminoglycosides are not generally coadministered in clinical practice, although more than one aminoglycoside may be present when switching from treatment with one to another. Samples that contain tobramycin in combination with either amikacin or kanamycin cannot be reliably quantitated by this assay.

In very rare cases, gammopathy, in particular type IgM (Waldenström's macroglobulinemia), may cause unreliable results.<sup>9</sup>

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

### ACTION REQUIRED

**Special Wash Programming:** The use of special wash steps is mandatory when certain test combinations are run together on Roche/Hitachi **cobas c** systems. The latest version of the carry-over evasion list can be found with the NaOHD-SMS-SmpCln1+2-SCCS Method Sheets. For further instructions refer to the operator's manual. **cobas c** 502 analyzer: All special wash programming necessary for avoiding carry-over is available via the **cobas** link, manual input is not required.

**Where required, special wash/carry-over evasion programming must be implemented prior to reporting results with this test.**

### Limits and ranges

#### Measuring range

0.33-10 µg/mL (0.71-21.4 µmol/L)

Manually dilute samples above the measuring range 1 + 1 with the Preciset TDM I Diluent (0 µg/mL) and reassay. Multiply the result by 2 to obtain the specimen value.

#### Lower limits of measurement

##### Lower detection limit of the test

0.33 µg/mL (0.71 µmol/L)

The lower detection limit represents the lowest measurable analyte level that can be distinguished from zero. It is calculated as the value lying 2 standard deviations above that of the lowest standard (standard 1 + 2 SD, repeatability, n = 21).

# TOBR2

## Tobramycin



### Expected values

Investigator	Peak		Trough	
	µg/mL	µmol/L	µg/mL	µmol/L
Baselt and Cravey <sup>10</sup>	6-10	12.8-21.4	0.5-2.0	1.1-4.3
Sande and Mandell <sup>4</sup>	5-8	10.7-17.1	1-2	2.1-4.3
Dipersio <sup>11</sup>	4-8	8.6-17.1	1-2	2.1-4.3
Lew <sup>12</sup>	5-10	10.7-21.4	-	-

Equivalent diagnostic technologies have shown that in most adults, a peak therapeutic response is achieved with tobramycin concentrations in the 6-10 µg/mL (12.8-21.4 µmol/L) range and trough concentrations in the 0.5-2.0 µg/mL (1.1-4.3 µmol/L) range. A peak therapeutic range is suggested for optimal antimicrobial effectiveness. Concentrations above the therapeutic range for a prolonged period of time or in patients with pre-existing renal impairment can cause nephrotoxicity and/or hearing impairment. Elevated or increasing trough levels are an indication of drug accumulation due to renal impairment. Both peak and trough levels should be monitored to ensure prevention of serious complications associated with drug dosage.

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

### Specific performance data

Representative performance data on the analyzers are given below. Results obtained in individual laboratories may differ.

### Precision

Precision was determined using human samples and controls in a modified NCCLS EP5-A protocol (repeatability n = 63, intermediate precision n = 63). The following results were obtained on a Roche/Hitachi **cobas c 501** analyzer.

### Serum/Plasma

Repeatability	Mean		SD		CV
	µg/mL	µmol/L	µg/mL	µmol/L	%
Control 1	1.63	3.49	0.07	0.15	4.2
Control 2	3.64	7.79	0.10	0.21	2.8
Control 3	7.72	16.52	0.22	0.47	2.9
HS 1	4.02	8.60	0.11	0.24	2.7
HS 2	8.72	18.66	0.20	0.43	2.3

Intermediate precision	Mean		SD		CV
	µg/mL	µmol/L	µg/mL	µmol/L	%
Control 1	1.63	3.49	0.07	0.15	4.5
Control 2	3.64	7.79	0.11	0.24	3.1
Control 3	7.72	16.52	0.21	0.45	2.8
HS 1	4.02	8.60	0.12	0.26	2.9
HS 2	8.72	18.66	0.21	0.45	2.4

### Method comparison

#### Serum/plasma

Tobramycin values for human serum and plasma samples obtained on a Roche/Hitachi **cobas c 501** analyzer (y) were compared with those determined using the corresponding reagent on a Roche/Hitachi 917 analyzer (x) and on a COBAS INTEGRA 700 analyzer (x).

Roche/Hitachi 917 analyzer	Sample size (n) = 69
Passing/Bablok <sup>13</sup>	Linear regression
$y = 0.960x + 0.020 \text{ µg/mL}$	$y = 0.973x + 0.007 \text{ µg/mL}$
$r = 0.946$	$r = 0.998$

The sample concentrations were between 0.370 and 9.92 µg/mL (0.792 and 21.2 µmol/L).

COBAS INTEGRA 700 analyzer	Sample size (n) = 68
Passing/Bablok <sup>13</sup>	Linear regression
$y = 0.963x + 0.152 \text{ µg/mL}$	$y = 0.966x + 0.193 \text{ µg/mL}$
$r = 0.916$	$r = 0.995$

The sample concentrations were between 0.3 and 10.0 µg/mL (0.642 and 21.4 µmol/L).

### Analytical specificity

The following compounds were tested for cross-reactivity.

Compound	Concentration Tested (µg/mL)	% Cross-reactivity
Amikacin	100	1.6
Carbenicillin	1000	ND
Cephalothin	1000	ND
Chloramphenicol	1000	ND
Clindamycin	1000	ND
Erythromycin	1000	ND
Gentamicin	100	0.5
Gentamicin	25	ND
Kanamycin	100	5.3
Neomycin	100	ND
Netilmicin	100	ND
Penicillin G	1000	ND
Sisomicin	100	ND
Streptomycin	100	ND
Sulphamethoxazole	600	ND
Tetracycline	1000	ND
Trimethoprim	25	ND
Vancomycin	200	ND
ND = not detectable		

Tests were performed on 16 drugs. No significant interference with the assay was found.

Acetaminophen	Doxycycline (Tetracycline)
Acetyl cysteine	Ibuprofen
Acetylsalicylic acid	Levodopa
Ampicillin-Na	Methyldopa+1.5 H <sub>2</sub> O
Ascorbic acid	Metronidazole
Ca-Dobesilate	Phenylbutazone
Cefoxitin	Rifampicin
Cyclosporine	Theophylline

### References

- Le Goffic F, Capmav ML, Tangy F, et al. Mechanism of action of aminoglycoside antibiotics, binding studies of tobramycin and its 6'-N-acetyl derivatives to the bacterial ribosome and its subunits. Eur J Biochem 1979;102:73-81.
- Cipolle RJ, Seifert RD, Zaske DE, et al. Systematically individualizing tobramycin dosage regimens. J Clin Pharm 1980;20:570-580.
- Naber KG, Westinfelder SR, Madsen PO. Pharmacokinetics of the aminoglycoside antibiotic tobramycin in humans. Antimicrob Agents Chemother 1973;3:469-473.




**Tobramycin**

- 4 Sande MA, Mandell GL. Antimicrobial agents, the aminoglycosides. In: Gilman AG, Goodman LS, Gilman A, eds. The Pharmacological Basis of Therapeutics. New York, NY: MacMillan 1980;1162-1180.
- 5 Hsu P, Ernst R, Levy M. EMIT 2000 tobramycin and vancomycin assays [abstract]. Clin Chem (Abstract 762) 2000;46(suppl 6):A195.
- 6 Guder W, Fonseca-Wollheim W, Heil O, et al. Maximum permissible transport and storage times for analysis of blood (serum, plasma), urine and cerebrospinal fluid. DG Klinische Chemische Mitteilungen 1995;26:205-224.
- 7 Tietz NW, ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia, PA: WB Saunders Company 1995;880.
- 8 Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-475.
- 9 Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med 2007;45(9):1240-1243.
- 10 Baselt RC, Cravey RH. Disposition of Toxic Drugs and Chemicals in Man. 3rd ed 1990:805-807.
- 11 Dipersio JR. Gentamicin and other aminoglycosides. In: Pesce AJ and Kaplan LA eds. Methods in Clinical Chemistry. St.Louis, Mo: CV Mosby Co 1987.
- 12 Lew M. Interpretation of aminoglycoside serum levels. Hosp Pharm 1979;14:465-472.
- 13 Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

**Symbols**

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard.

	Contents of kit
	Volume after reconstitution or mixing
	Global Trade Item Number

**FOR US CUSTOMERS ONLY: LIMITED WARRANTY**

Roche Diagnostics warrants that this product will meet the specifications stated in the labeling when used in accordance with such labeling and will be free from defects in material and workmanship until the expiration date printed on the label. THIS LIMITED WARRANTY IS IN LIEU OF ANY OTHER WARRANTY, EXPRESS OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. IN NO EVENT SHALL ROCHE DIAGNOSTICS BE LIABLE FOR INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES.

COBAS, COBAS C, COBAS INTEGRA, ONLINE TDM and PRECISET are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2015, Roche Diagnostics



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim  
www.roche.com



Distribution in USA by:  
Roche Diagnostics, Indianapolis, IN  
US Customer Technical Support 1-800-428-2336