APOAT

Tina-quant Apolipoproteinas A-1, 2-a versija

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REF	CONTENT	Analizatorius (-iai), su kuriuo (-iais) galima naudoti cobas c pakuotę (-es)
03032566 122	Tina-quant Apolipoprotein A-1 ver.2 100 tyrimų	Sistemos-ID 07 6568 6 Roche/Hitachi cobas c 311, cobas c 501/502
12172623 122	Calibrator f.a.s. Lipids (3 x 1 mL)	Kodas 424
12172623 160	Calibrator f.a.s. Lipids (3 x 1 mL, skirta JAV)	Kodas 424
10781827 122	Precinorm L (4 x 3 mL)	Kodas 304
11285874 122	Precipath L (4 x 3 mL)	Kodas 305
05117003 190	PreciControl ClinChem Multi 1 (20 x 5 mL)	Kodas 391
05947626 190	PreciControl ClinChem Multi 1 (4 x 5 mL)	Kodas 391
05947626 160	PreciControl ClinChem Multi 1 (4 x 5 mL, skirta JAV)	Kodas 391
05117216 190	PreciControl ClinChem Multi 2 (20 x 5 mL)	Kodas 392
05947774 190	PreciControl ClinChem Multi 2 (4 x 5 mL)	Kodas 392
05947774 160	PreciControl ClinChem Multi 2 (4 x 5 mL, skirta JAV)	Kodas 392
04489357 190	Diluent NaCl 9 % (50 mL)	Sistemos-ID 07 6869 3

Lietuvių

Sistemos informacija

Skirta **cobas c** 311/501 analizatoriams:**APOAT**: ACN 168Skirta **cobas c** 502 analizatoriui:**APOAT**: ACN 8168

Paskirtis

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- Becker W, Rapp W, Schwick HG, et al. Methoden zur quantitativen Bestimmung von Plasmaproteinen durch Immunpräzipitation. Z Klin Chem Klin Biochem. 1968;6:113-122.
- Siedel J, Schiefer S, Rosseneu M, et al. Immunturbidimetric Method for Routine Determinations of Apolipoproteins A-I, A-II and B in Normo- and Hyperlipemic Sera compared with Immunonephelometry. Clin Chem. 1988;34(9):1821-1825.
- Rifai N, King ME. Immunturbidimetric Assays of Apolipoproteins A, AI, AII and B in serum. Clin Chem. 1986;32:957-961.
- Naito HK. Reliability of Lipid, Lipoprotein and Apolipoprotein Measurements. Clin Chem. 1988;34:B84-B94.
- Use of Anticoagulants in Diagnostic Laboratory Investigations. WHO Publication WHO/DIL/LAB/99.1 Rev. 2. Jan. 2002.
- Evans K, Mitcheson J, Laker M. Effect of Storage at 4 °C and -20 °C on Lipid, Lipoprotein, and Apolipoprotein Concentrations. Clin Chem. 1995;41:392-396.
- Marcovina SM, Albers JJ, Dati F, et al. International Federation of Clinical Chemistry Standardization Project for Measurements of Apolipoproteins A-I and B. Clin Chem. 1991;37:1676-1682.
- Albers JJ, Marcovina SM, Kennedy H. International Federation of Clinical Chemistry Standardization Project for Measurements of Apolipoproteins A-I and B. II. Evaluation and Selection of Candidate Reference Materials. Clin Chem. 1992;38:658-662.
- Marcovina SM, Albers JJ, Henderson LO, et al. International Federation of Clinical Chemistry Standardization Project for Measurements of Apolipoproteins A-I and B. III. Comparability of Apolipoprotein A-I Values by Use of International Reference Material. Clin Chem. 1993;39:773-781.
- Marcovina SM, Albers JJ, Kennedy H, et al. International Federation of Clinical Chemistry Standardization Project for Measurements of Apolipoproteins A-I and B. IV. Comparability of Apolipoprotein B Values by Use of International Reference Material. Clin Chem. 1994;40:586-592.
- Young DS, et al. Effects of drugs on clinical laboratory tests. American College of Physicians 1998.
- Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem. 1986;32:470-475.
- Breuer J. Report on the Symposium "Drug effects in Clinical Chemistry Methods". Eur J Clin Chem Clin Biochem 1996;34:385-386.
- Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385.
- Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med 2007;45(9):1240-1243.
- Data on file at Roche Diagnostics.
- Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.

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Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim
www.roche.com

JAV platina:

Roche Diagnostics, Indianapolis, IN

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cobas®

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© 2017, Roche Diagnostics