REF	CONTENT	Analizatorius (-iai), su kuriuo (-iais) galima naudoti cobas c pakuotę (-es)
04975723 190	Tina-quant Cystatin C 225 tyrimai	Sistemos-ID 07 6988 6 Roche/Hitachi cobas c 311, cobas c 501/502
04975901 190	C.f.a.s. Cystatin C (4 x 1 mL)	Kodas 404
04975936 190	Cystatin C Control Set (2 x 4 x 1 mL)	Kodas 121 (Control I (low)) Kodas 122 (Control II (high))
04489357 190	Diluent NaCl 9 % (50 mL)	Sistemos-ID 07 6869 3

Lietuvių

Sistemos informacija

Skirta **cobas c** 311/501 analizatoriams:

CYSC: ACN 431

Skirta **cobas c** 502 analizatoriui:

CYSC: ACN 8431

Paskirtis

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Nuorodos

- 1 Levey AS, Coresh J, Balk E, et al. National Kidney Foundation practice guidelines for chronic kidney disease: evaluation, classification, and stratification. *Ann Intern Med* 2003;139:137-147.
- 2 Rule AD, Larson TS, Bergstralh EJ, et al. Using serum creatinine to estimate glomerular filtration rate: accuracy in good health and in chronic kidney disease. *Ann Intern Med* 2004;141:929-937.
- 3 Wasen E, Isoaho R, Mattila K, et al. Estimation of glomerular filtration rate in the elderly: a comparison of creatinine based formulae with serum Cystatin C. *J Intern Med* 2004;256:70-78.
- 4 Levey AS, Coresh JK. DOQI clinical practice guidelines on chronic kidney disease. Guideline 4. Estimation of GFR. *Am J Kidney Dis* 2002;39 (suppl1):76-92.
- 5 Kyhse-Andersen J, Schmidt C, Nordin G, et al. Serum Cystatin C, determined by a rapid, automated particle-enhanced turbidimetric method, is a better marker than serum creatinine for glomerular filtration rate. *Clin Chem* 1994;40(10):1921-1926.
- 6 Grubb A, Simonsen O, Sturfelt G, et al. Serum concentration of Cystatin C, factor D and β_2 -microglobulin as a measure of glomerular filtration rate. *Acta Med Scand* 1985;218(5):499-503.
- 7 Jung K, Jung M. Cystatin C: a promising marker of glomerular filtration rate to replace creatinine. *Nephron* 1995;70(3):370-371.
- 8 Newman DJ, Thakkar H, Edwards RG, et al. Serum Cystatin C measured by automated immunoassay: a more sensitive marker of changes in GFR than serum creatinine. *Kidney Int* 1995;47(1):312-318.
- 9 Mussap M, Dalla Vestra M, Fioretto P, et al. Cystatin C is a more sensitive marker than creatinine for the estimation of GFR in type 2 diabetic patients. *Kidney Int* 2002;61:1453-1461.
- 10 Darnidharka VR, Kwon C, Stevens G. Serum Cystatin C is superior to serum creatinine as a marker of kidney function: a meta analysis. *Am J Kidney Dis* 2002;40:221-226.
- 11 Risch L, Drexel H, Huber AR. Differences in glomerular filtration rates estimated by two cystatin C-based equations. *Clinical Chemistry* 2005;51:2211-2212.
- 12 Grubb A, Nyman U, Bjork J, et al. Simple Cystatin C-Based prediction equations for glomerular filtration rate compared with the modification of diet in renal disease prediction equation for adults and the Schwartz and the Counahan-Barratt prediction equations for children. *Clinical Chemistry* 2005;51:1420-1431.
- 13 Hoek FJ, Kemperman FA, Krediet RT. A comparison between Cystatin C, plasma creatinine and the Cockcroft and Gault formula for the estimation of glomerular filtration rate. *Nephrol Dial Transplant* 2003;18:2034-2031.
- 14 Lassus J, Harjola VP, Sund R, et al. Prognostic value of Cystatin C in acute heart failure in relation to other markers of renal function and NT-proBNP. *European Heart Journal* (February) doi:10.1093/eurheart/ehl507.
- 15 Data on file at Roche Professional Diagnostics.
- 16 Bilrup-Jensen S. Protein standardization II: Dry mass determination procedure for the determination of the dry mass of a pure protein precipitation. *Clin Chem Lab Med* 2001;39(11):1090-1097.

- 17 Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-475.
- 18 Breuer J. Report on the Symposium "Drug effects in Clinical Chemistry Methods". Eur J Clin Chem Clin Biochem 1996;34:385-386.
- 19 Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385.
- 20 Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med 2007;45(9):1240-1243.
- 21 Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.

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Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim
www.roche.com

Platintojas JAV:

Roche Diagnostics, Indianapolis, IN

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© 2013, Roche Diagnostics