REF	CONTENT	Analizatorius (-iai), su kuriuo (-iais) galima naudoti <b>cobas c</b> pakuotę (-es)
03001245 322	Tina-quant D-Dimer (100 tyrimų)	Sistemos-ID 07 6458 2 Roche/Hitachi <b>cobas c</b> 311, <b>cobas c</b> 501/502
11556495 216	D-Dimer Calibrator (1 x 0.5 mL) ir Zero Standard (1 x 2.5 mL)	Kodas 680-685
11556509 216	D-Dimer Control I/II (2 x 0.5 mL)	Kodas 200/Kodas 201

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Skirta **cobas c** 311/501 analizatoriams:

**D-DI:** ACN 123 (citruota plazma)

**D-DIH:** ACN 125 (Li-heparino plazma)

Skirta **cobas c** 502 analizatoriui:

**D-DI:** ACN 8123 (citruota plazma)

**D-DIH:** ACN 8125 (Li-heparino plazma)

**Paskirtis**

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**Nuorodos**

- 1 Gaffney PJ. Fibrinolysis Supplement 2 1993;7:2-8.
- 2 Matsuda M, Yoshida N, Terukina S, et al. Fibrinogen 4. Current basic and clinical aspects. Amsterdam/New York/Oxford: Elsevier Science Publishers 1990;43-48.
- 3 Agency for Healthcare Research and Quality, Evidence Report /Technology Assessment Number 68: Diagnosis and Treatment of Deep Venous Thrombosis and Pulmonary Embolism: Summary. AHRQ Pub No. 03-E012, January, Full report available online at [www.ahrq.com](http://www.ahrq.com) 2003.
- 4 American College of Emergency Physicians Board of Directors. Clinical Policy: Critical Issues in the Evaluation and Management of Adult Patients Presenting with Suspected Lower-Extremity Deep Venous Thrombosis *Ann Em Med* 2003;42(1):124.
- 5 Ramzi DW, Leeper KV. DVT and Pulmonary Embolism: Part 1, Diagnosis. *Am. Fam. Phys* 2004;69(12):2829.
- 6 American College of Emergency Physicians Board of Directors. Clinical Policy: Critical Issues in the Evaluation and Management of Adult Patients Presenting with Suspected Pulmonary Embolism. *Ann Em Med* 2003;41:257.
- 7 Schutgens REG, Ackermack P, Haas FJLM, et al. Combination of a Normal D-Dimer Concentration and a Non-High Pretest Clinical Probability Score is a Safe Strategy to Exclude Deep Venous Thrombosis. *Circulation* 2003;107:593-597.
- 8 Schutgens RE, Haas FJ, Biesma DH. Reduced efficacy of clinical probability score and D-Dimer assay in elderly subjects suspected of having deep vein thrombosis. *Br J Haemat* 2005;129:653-657.
- 9 LeClerq LGL, Lusitan JG, Kooy MvM, et al. Ruling out clinically suspected pulmonary embolism by assessment of clinical probability and D-Dimer levels: a management study. *Thromb Haemost* 2003;89:97-103.
- 10 Van Belle A, Büller HR, Huisman MV, et al. for the Christopher Study Investigators. Effectiveness of managing suspected pulmonary embolism using an algorithm combining clinical probability, D-dimer testing, and computed tomography. *JAMA* 2006. 295(2), 172-179.
- 11 Knecht MF, Heinrich F. Clinical Evaluation of an Immunoturbidimetric D-Dimer Assay in the Diagnostic Procedure of Deep Vein Thrombosis and Pulmonary Embolism. *Thromb Res* 1997;88:413-417.
- 12 Janssen MCH, Heebles AE, deMetz M, et al. Reliability of Five Rapid D-Dimer Assays Compared to ELISA in the Exclusion of Deep Venous Thrombosis. *Thromb Haemost* 1997;77(2):262-266.
- 13 Lindahl TL, Lundahl TH, Fransson SG. Evaluation of an automated micro-latex D-Dimer assay (Tina-quant on Hitachi 911) in symptomatic outpatients. *Thromb Haemost* 1999;82(6):1772-1773.
- 14 Van der Graaf F, van den Borne H, van der Kolk M, et al. Exclusion of deep venous thrombosis with D-dimer testing--comparison of 13 D-dimer methods in 99 outpatients suspected of deep venous thrombosis using venography as reference standard. *Thromb Haemost* 2000;83(2):191-198.
- 15 Fünfsinn N, Caliezi F, Biasiutti FD, et al. Rapid D-Dimer testing and pre-test clinical probability in the exclusion of deep venous thrombosis in symptomatic outpatients. *Blood Coagul Fibrinolysis* 2001;12:165-170.
- 16 Diamond S, Goldweber R, Katz S. Use of D-Dimer to aid in excluding deep venous thrombosis in ambulatory patients. *Am J Surg* 2005;189:23-26.
- 17 Schutgens RE, Haas FJ, Gerritsen WB, et al. The usefulness of five D-Dimer assays in the exclusion of deep venous thrombosis. *J Thromb Haemost* 2003;1:976-981.
- 18 Stolba R, Lenglinger FX, Rezanka E, et al. Diagnostic Value of a new, quantitative D-Dimer assay for the exclusion of pulmonary embolism in symptomatic patients. *J Lab Med* 2000;24(3):153-157.
- 19 De Monyé W, Sanson B-J, Büller HR, et al. ANTELOPE study group. The performance of two rapid quantitative D-Dimer assays in 287 patients with clinically suspected pulmonary embolism. *Thromb Res* 2002;107:283-286.
- 20 Söhne M, Kamphuisen PW, van Mierlo PJWB, et al. Diagnostic strategy using a modified clinical decision rule and D-Dimer test to rule out pulmonary embolism in elderly in- and outpatients. *Thromb Haemost* 2005;94(1):206-210.
- 21 Jennersjö C, Fagerberg I, Karlander S, et al. Normal D-Dimer concentration is a common finding in symptomatic outpatients with distal deep vein thrombosis. *Blood Coagul Fibrinolysis* 2005;16:517-523.
- 22 Angstwurm MW, Reininger AJ, Spannagl M. D-Dimer as marker for microcirculatory failure: correlation with LOD and APACHE II scores. *Thromb Res* 2004;113(6):353-359.
- 23 Wakai A, Gleeson A, Winter D. Role of fibrin D-Dimer testing in emergency medicine. *Emerg Med J* 2003;20:319-325.
- 24 Dempfle CE. Bestimmung des D-Dimer-Antigens in der klinischen Routine. 102, Ausgabe 24 vom 17.06.2005.
- 25 Schutgens REG, Haas FJML, Ruven HJT, et al. No Influence of Heparin Plasma and Other (Pre)analytic variables on D-Dimer Determinations. *Clin Chem* 2002;48(9):1611-1613.
- 26 Guder WG, Narayanan S, Wisser H, et al. List of Analytes; Preanalytical Variables. Brochure in: Samples: From the Patient to the Laboratory. Darmstadt: GIT-Verlag 1996.
- 27 Adema E, Gebert U. Pooled patient samples as reference material for D-Dimer. *Thromb Res* 1995;80(1):85-88.
- 28 Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. *Clin Chem* 1986;32:470-475.
- 29 Breuer J. Report on the Symposium "Drug effects in Clinical Chemistry Methods". *Eur J Clin Chem Clin Biochem* 1996;34:385-386.
- 30 Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. *Ann Clin Biochem* 2001;38:376-385.
- 31 Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. *Clin Chem Lab Med* 2007;45(9):1240-1243.
- 32 Dempfle CE, Hafner G, Lestin HG, et al. Multizentrische Evaluierung von Tina-quant D-Dimer. *J Lab Med* 1996;20:31-37.
- 33 Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. *J Clin Chem Clin Biochem* 1988 Nov;26(11):783-790.
- 34 Wells PS, Ginseberg JF, Anderson DR, et al. Use of a clinical model for safe management of patients with suspected pulmonary embolism. *Ann Intern Med* 1998;129:997-1005.

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JAV platiną:

Roche Diagnostics, Indianapolis, IN

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© 2016, Roche Diagnostics
Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim