REF	CONTENT	Analizatorius (-iai), su kuriuo (-iais) galima naudoti cobas c pakuotę (-es)
04912551 190	Tina-quant D-Dimer Gen.2 100 tyrimų	Sistemos-ID 07 6932 0 Roche/Hitachi cobas c 311, cobas c 501/502
05050901 190	D-Dimer Gen.2 Calibrator Set (6 x 0.5 mL)	Kodai 764-769
05050936 190	D-Dimer Gen.2 Control I/II (2 x 2 x 1 mL)	Kodas 242 Control I Kodas 243 Control II

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Skirta **cobas c** 311/501 analizatoriams:

D-DI2: ACN 102 (citruota plazma)

DDI2H: ACN 403 (Heparino/EDTA plazma)

Skirta **cobas c** 502 analizatoriui:

D-DI2: ACN 8102 (citruota plazma)

DDI2H: ACN 8403 (Heparino/EDTA plazma)

Paskirtis

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Nuorodos

- Gaffney PJ. Fibrinolysis Supplement 2 1993;7:2-8.
- Fibrinogen 4. Current basic and clinical aspects. Matsuda M et al. Amsterdam/New York/Oxford: Elsevier Science Publishers, 1990:43-8.
- Agency for Healthcare Research and Quality, Evidence Report /Technology Assessment Number 68: Diagnosis and Treatment of Deep Venous Thrombosis and Pulmonary Embolism: Summary. AHRQ Pub No. 03-E012, January, Full report available online at www.ahrq.com 2003.
- American College of Emergency Physicians Board of Directors. Clinical Policy: Critical Issues in the Evaluation and Management of Adult Patients Presenting with Suspected Lower-Extremity Deep Venous Thrombosis Ann Em Med 2003;42(1):124.
- Ramzi DW, Leeper KV. DVT and Pulmonary Embolism: Part 1, Diagnosis. Am. Fam. Phys 2004;69(12):2829.
- American College of Emergency Physicians Board of Directors. Clinical Policy: Critical Issues in the Evaluation and Management of Adult Patients Presenting with Suspected Pulmonary Embolism. Ann Em Med 2003;41:257.
- Schutgens REG, Ackermack P, Haas FJLM, et al. Combination of a Normal D-Dimer Concentration and a Non-High Pretest Clinical Probability Score is a Safe Strategy to Exclude Deep Venous Thrombosis. Circulation 2003;107:593-597.
- Schutgens RE, Haas FJ, Biesma DH. Reduced efficacy of clinical probability score and D-Dimer assay in elderly subjects suspected of having deep vein thrombosis. Br J Haemat 2005;129:653-657.
- LeClerq LGL, Lusitan JG, Kooy MvM, et al. Ruling out clinically suspected pulmonary embolism by assessment of clinical probability and D-Dimer levels: a management study. Thromb Haemost 2003;89:97-103.
- Van Belle A, Büller HR, Huisman MV, et al. for the Christopher Study Investigators. Effectiveness of managing suspected pulmonary embolism using an algorithm combining clinical probability, D-dimer testing, and computed tomography. JAMA 2006. 295(2), 172-179.
- Djurabi RK, Klok FA, Nijkeuter M, et al. Comparison of the clinical usefulness of two quantitative D-Dimer tests in patients with a low clinical probability of Pulmonary Embolism. Thromb Res 2009;123:771-774.
- Knecht MF, Heinrich F. Clinical Evaluation of an Immunoturbidimetric D-Dimer Assay in the Diagnostic Procedure of Deep Vein Thrombosis and Pulmonary Embolism. Thromb Res 1997;88:413-417.
- Janssen MCH, Heebles AE, deMetz M, et al. Reliability of Five Rapid D-Dimer Assays Compared to ELISA in the Exclusion of Deep Venous Thrombosis. Thromb Haemost 1997;77(2):262-266.
- Lindahl TL, Lundahl TH, Fransson SG. Evaluation of an automated micro-latex D-Dimer assay (Tina-quant on Hitachi 911) in symptomatic outpatients. Thromb Haemost 1999;82(6):1772-1773.
- Van der Graaf F, van den Borne H, van der Kolk M, et al. Exclusion of deep venous thrombosis with D-Dimer testing: Comparison of 13 D-Dimer methods in 99 outpatients suspected of deep venous thrombosis using venography as reference standard. Thromb Haemost 2000;83:191-8.
- Fünfsinn N, Caliezi F, Biasutti FD, et al. Rapid D-Dimer testing and pre-test clinical probability in the exclusion of deep venous thrombosis in symptomatic outpatients. Blood Coagul Fibrinolysis 2001;12:165-170.
- Diamond S, Goldweber R, Katz S. Use of D-Dimer to aid in excluding deep venous thrombosis in ambulatory patients. Am J Surg 2005;189:23-26.
- Schutgens RE, Haas FJ, Gerritsen WB, et al. The usefulness of five D-Dimer assays in the exclusion of deep venous thrombosis. J Thromb Haemost 2003;1:976-981.
- Stolba R, Lenglinger FX, Rezanka E, et al. Diagnostic Value of a new, quantitative D-Dimer assay for the exclusion of pulmonary embolism in symptomatic patients. J Lab Med 2000;24(3):153-157.
- De Monyé W, Sanson B-J, Büller HR, et al. ANTELOPE study group. The performance of two rapid quantitative D-Dimer assays in 287 patients with clinically suspected pulmonary embolism. Thromb Res 2002;107:283-286.

- 21 Söhne M, Kamphuisen PW, van Mierlo PJWB, et al. Diagnostic strategy using a modified clinical decision rule and D-Dimer test to rule out pulmonary embolism in elderly in- and outpatients. *Thromb Haemost* 2005;94(1):206-210.
- 22 Jennersjö C, Fagerberg I, Karlander S, et al. Normal D-Dimer concentration is a common finding in symptomatic outpatients with distal deep vein thrombosis. *Blood Coagul Fibrinolysis* 2005;16:517-523.
- 23 Angstwurm MW, Reininger AJ, Spannagl M. D-Dimer as marker for microcirculatory failure: correlation with LOD and APACHE II scores. *Thromb Res* 2004;113(6):353-359.
- 24 Wakai A, Gleeson A, Winter D. Role of fibrin D-Dimer testing in emergency medicine. *Emerg Med J* 2003;20:319-325.
- 25 Dempfle CE. Bestimmung des D-Dimer-Antigens in der klinischen Routine. 102, Ausgabe 24 vom 17.06.2005.
- 26 Schutgens REG, Haas FJML, Ruven HJT, et al. No Influence of Heparin Plasma and Other (Pre)analytic variables on D-Dimer Determinations. *Clin Chem* 2002;48(9):1611-1613.
- 27 Guder WG, Narayanan S, Wisser H, et al. List of Analytes; Preanalytical Variables. Brochure in: Samples: From the Patient to the Laboratory. Darmstadt: GIT-Verlag 1996.
- 28 Adema E, Gebert U. Pooled patient samples as reference material for D-Dimer. *Thromb Res* 1995;80(1):85-88.
- 29 Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. *Clin Chem* 1986;32:470-475.
- 30 Breuer J. Report on the Symposium "Drug effects in Clinical Chemistry Methods". *Eur J Clin Chem Clin Biochem* 1996;34:385-386.
- 31 Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. *Ann Clin Biochem* 2001;38:376-385.
- 32 Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. *Clin Chem Lab Med* 2007;45(9):1240-1243.
- 33 Dempfle CE, Hafner G, Lestin HG, et al. Multizentrische Evaluierung von Tina-quant D-Dimer. *J Lab Med* 1996;20:31-37.
- 34 Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. *J Clin Chem Clin Biochem* 1988 Nov;26(11):783-790.
- 35 Wells PS, Ginseberg JF, Anderson DR, et al. Use of a clinical model for safe management of patients with suspected pulmonary embolism. *Ann Intern Med* 1998;129:997-1005.

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© 2015, Roche Diagnostics
Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim