REF	CONTENT		Analizatorius (-iai), su kuriuo (-iais) galima naudoti cobas c pakuotę (-es)
20763454 122	Tina-quant Soluble Transferrin Receptor (80 tyrimų)	Sistemos-ID 07 6345 4	Roche/Hitachi cobas c 311, cobas c 501/502
12148331 122	Preciset sTfR (5 x 1 mL)	Kodai 750-754	
12148340 122	sTfR Control Set		
	Level I (2 x 3 mL)	Level I kodas 211	
	Level II (2 x 3 mL)	Level II kodas 212	
04489357 190	Diluent NaCl 9 % (50 mL)	Sistemos-ID 07 6869 3	

Lietuvių

Sistemos informacija

Skirta **cobas c** 311/501 analizatoriams:**STFR**: ACN 665Skirta **cobas c** 502 analizatoriui:**STFR**: ACN 8665

Paskirtis

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Nuorodos

- Wick M, Pinggera W, Lehmann P, eds. Iron Metabolism, Diagnosis and Therapy of Anemias. 4th ed. Vienna/New York: Springer Verlag 1999.
- Ferguson BJ, Skikne BS, Simpson KM, et al. Serum transferrin receptor distinguishes the anemia of chronic diseases from iron deficiency anemia. J Lab Clin Med 1992;19:385-390.
- Baynes RD, Skikne BS, Cook JD. Circulating transferrin receptors and assessment of iron status. J Nutr Biochem 1994;5:322-330.
- Kuiper-Kramer EPA, Huisman CMS, van Raan J, et al. Analytical and clinical implications of soluble transferrin receptors. Eur J Clin Chem Clin Biochem 1996;34:645-649.
- Baynes RD. Assessment of iron status. Clin Biochem 1996;29:209-215.
- De Paoli Vitali E, Ricci G, Perini L, et al. The determination of plasma transferrin receptor as a good index of erythropoietic activity in renal anemia and after renal transplantation. Nephron 1996;72:552-556.
- Punnonen K, Irjala K, Rajamäki A. Serum transferrin receptor and its ratio to serum ferritin in the diagnosis of iron deficiency. Blood 1997;89:1052-1057.
- Suominen P, Punnonen K, Rajamäki A, et al. Serum transferrin receptor and transferrin receptor-ferritin index identify healthy subjects with subclinical iron deficits. Blood 1998;92:2934-2939.
- Haeckel R, ed. Evaluation methods in laboratory medicine. VCH Verlagsgesellschaft, Weinheim 1993.
- Cook JD, Skikne BS, Baynes RD. Serum transferrin receptor. Annu Rev Med 1993;44:63-74.
- Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-475.
- Breuer J. Report on the Symposium "Drug effects in Clinical Chemistry Methods". Eur J Clin Chem Clin Biochem 1996;34:385-386.
- Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385.
- Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med 2007;45(9):1240-1243.

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Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim
www.roche.com

JAV platina:

Roche Diagnostics, Indianapolis, IN

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© 2016, Roche Diagnostics