

REF



SYSTEM

12017709 122

200

MODULAR ANALYTICS E170

cobas e 411

cobas e 601

cobas e 602

English

System information

For **cobas e 411** analyzer: test number 021

For MODULAR ANALYTICS E170, **cobas e 601** and **cobas e 602** analyzers: Application Code Number 002

Intended use

Immunoassay for the in vitro quantitative determination of thyroxine in human serum and plasma.

The electrochemiluminescence immunoassay "ECLIA" is intended for use on Elecsys and **cobas e** immunoassay analyzers.

Summary

The hormone thyroxine (T4) is the main product secreted by the thyroid gland and is an integral component of the hypothalamus-anterior pituitary-thyroid regulating system. It has the function of anabolically influencing metabolism. Thyroxine is formed in a coupling reaction from two DIT molecules (3,5-diiodotyrosine) in the thyroid gland. It is stored bound to thyroglobulin in the lumina of the thyroid follicles and is secreted as required under the influence of TSH.^{1,2}

The major part (> 99 %) of total thyroxine (T4) in serum is present in protein-bound form. As the concentrations of the transport proteins in serum are subject to exogenous and endogenous effects, the status of the binding proteins must also be taken into account in the assessment of the thyroid hormone concentration in serum. If this is ignored, changes in the binding proteins (e.g. due to estrogen-containing preparations, during pregnancy or in the presence of a nephrotic syndrome etc.) can lead to erroneous assessments of the thyroid metabolic state.^{3,4,5,6,7}

The determination of T4 can be utilized for the following indications: the detection of hyperthyroidism, the detection of primary and secondary hypothyroidism, and the monitoring of TSH-suppression therapy.⁸

The Elecsys T4 assay employs a competitive test principle with an antibody specifically directed against T4. Endogenous T4, released by the action of 8-anilino-1-naphthalene sulfonic acid (ANS), competes with the added biotinylated T4-derivative for the binding sites on the antibodies labeled with the ruthenium complex⁹.

a) Tris(2,2'-bipyridyl)ruthenium(II)-complex (Ru(bpy)₃²⁺)

Test principle

Competition principle. Total duration of assay: 18 minutes.

- 1st incubation: 15 µL of sample and a T4-specific antibody labeled with a ruthenium complex; bound T4 is released from binding proteins in the sample by ANS.
- 2nd incubation: After addition of streptavidin-coated microparticles and biotinylated T4, the still-free binding sites of the labeled antibody become occupied, with formation of an antibody-hapten complex. The entire complex becomes bound to the solid phase via interaction of biotin and streptavidin.
- The reaction mixture is aspirated into the measuring cell where the microparticles are magnetically captured onto the surface of the electrode. Unbound substances are then removed with ProCell/ProCell M. Application of a voltage to the electrode then induces chemiluminescent emission which is measured by a photomultiplier.
- Results are determined via a calibration curve which is instrument-specifically generated by 2-point calibration and a master curve provided via the reagent barcode or e-barcode.

Reagents - working solutions

The reagent rackpack is labeled as T4.

- M Streptavidin-coated microparticles (transparent cap), 1 bottle, 12 mL:
Streptavidin-coated microparticles 0.72 mg/mL; preservative.

- R1 Anti-T4-Ab~Ru(bpy)₃²⁺ (gray cap), 1 bottle, 18 mL:

Polyclonal anti-T4-antibody (sheep) labeled with ruthenium complex 100 ng/mL; ANS 1 mg/mL; phosphate buffer 100 mmol/L, pH 7.4; preservative.

- R2 T4-biotin (black cap), 1 bottle, 18 mL:

Biotinylated T4 20 ng/mL; phosphate buffer 100 mmol/L, pH 7.4; preservative.

Precautions and warnings

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents.

Disposal of all waste material should be in accordance with local guidelines. Safety data sheet available for professional user on request.

Avoid foam formation in all reagents and sample types (specimens, calibrators and controls).

Reagent handling

The reagents in the kit have been assembled into a ready-for-use unit that cannot be separated.

All information required for correct operation is read in from the respective reagent barcodes.

Storage and stability

Store at 2-8 °C.

Do not freeze.

Store the Elecsys reagent kit **upright** in order to ensure complete availability of the microparticles during automatic mixing prior to use.

Stability:	
unopened at 2-8 °C	up to the stated expiration date
after opening at 2-8 °C	12 weeks
on the analyzers	8 weeks

Specimen collection and preparation

Only the specimens listed below were tested and found acceptable.

Serum collected using standard sampling tubes or tubes containing separating gel.

Li-, Na-heparin, K₃-EDTA and sodium citrate plasma.

When sodium citrate is used, the results must be corrected by + 10 %.

Criterion: Recovery within 90-110 % of serum value or slope 0.9-1.1 + intercept within < ± 2x analytical sensitivity (LDL) + coefficient of correlation > 0.95.

When sodium fluoride/potassium oxalate are used, the values found are by approximately 26 % lower than those for serum.

Stable for 7 days at 2-8 °C, 30 days at -20 °C.⁹ Freeze only once.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

Do not use heat-inactivated samples.

Do not use samples and controls stabilized with azide.

Ensure the samples, calibrators and controls are at 20-25 °C prior to measurement.

Due to possible evaporation effects, samples, calibrators and controls on the analyzers should be analyzed/measured within 2 hours.

Materials provided

See "Reagents – working solutions" section for reagents.

Materials required (but not provided)

- [REF] 12017717122, T4 CalSet, 4 x 1.0 mL
- [REF] 11731416190, PreciControl Universal, for 4 x 3.0 mL
- General laboratory equipment
- MODULAR ANALYTICS E170 or **cobas e** analyzer

Accessories for **cobas e** 411 analyzer:

- [REF] 11662988122, ProCell, 6 x 380 mL system buffer
- [REF] 11662970122, CleanCell, 6 x 380 mL measuring cell cleaning solution
- [REF] 11930346122, Elecsys SysWash, 1 x 500 mL washwater additive
- [REF] 11933159001, Adapter for SysClean
- [REF] 11706802001, AssayCup, 60 x 60 reaction cups
- [REF] 11706799001, AssayTip, 30 x 120 pipette tips
- [REF] 11800507001, Clean-Liner

Accessories for MODULAR ANALYTICS E170, **cobas e** 601 and **cobas e** 602 analyzers:

- [REF] 04880340190, ProCell M, 2 x 2 L system buffer
- [REF] 04880293190, CleanCell M, 2 x 2 L measuring cell cleaning solution
- [REF] 03023141001, PC/CC-Cups, 12 cups to prewarm ProCell M and CleanCell M before use
- [REF] 03005712190, ProbeWash M, 12 x 70 mL cleaning solution for run finalization and rinsing during reagent change
- [REF] 12102137001, AssayTip/AssayCup, 48 magazines x 84 reaction cups or pipette tips, waste bags
- [REF] 03023150001, WasteLiner, waste bags
- [REF] 03027651001, SysClean Adapter M

Accessories for all analyzers:

- [REF] 11298500316, ISE Cleaning Solution/Elecsys SysClean, 5 x 100 mL system cleaning solution

Assay

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

Resuspension of the microparticles takes place automatically prior to use. Read in the test-specific parameters via the reagent barcode. If in exceptional cases the barcode cannot be read, enter the 15-digit sequence of numbers (except for the **cobas e** 602 analyzer).

Bring the cooled reagents to approximately 20 °C and place on the reagent disk (20 °C) of the analyzer. Avoid foam formation. The system automatically regulates the temperature of the reagents and the opening/closing of the bottles.

Calibration

Traceability: The Elecsys T4 assay has been checked by ID-GC/MS (isotope dilution gas chromatography mass spectrometry) on various control materials.¹⁰

Every Elecsys reagent set has a barcoded label containing specific information for calibration of the particular reagent lot. The predefined master curve is adapted to the analyzer using the relevant CalSet.

Calibration frequency: Calibration must be performed once per reagent lot using fresh reagent (i.e. not more than 24 hours since the reagent kit was registered on the analyzer).

Calibration interval may be extended based on acceptable verification of calibration by the laboratory.

Renewed calibration is recommended as follows:

- after 8 weeks when using the same reagent lot
- after 7 days when using the same reagent kit on the analyzer

- as required: e.g. quality control findings outside the defined limits

Quality control

For quality control, use PreciControl Universal.

In addition, other suitable control material can be used.

Controls for the various concentration ranges should be run individually at least once every 24 hours when the test is in use, once per reagent kit, and following each calibration.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

If necessary, repeat the measurement of the samples concerned.

Follow the applicable government regulations and local guidelines for quality control.

Calculation

The analyzer automatically calculates the analyte concentration of each sample (either in nmol/L, µg/dL or µg/L).

Conversion factors:

$$\begin{aligned} \text{nmol/L} \times 0.077688 &= \mu\text{g/dL} \\ \mu\text{g/dL} \times 12.872 &= \text{nmol/L} \\ \text{nmol/L} \times 0.77688 &= \mu\text{g/L} \end{aligned}$$

Limitations - interference

The assay is unaffected by icterus (bilirubin < 633 µmol/L or < 37 mg/dL), hemolysis (Hb < 1.4 mmol/L or < 2.3 g/dL), lipemia (triglycerides < 28.5 mmol/L or < 2500 mg/dL) and biotin (< 409 nmol/L or < 100 ng/mL).

Criterion: Recovery within ± 10 % of initial value.

Samples should not be taken from patients receiving therapy with high biotin doses (i.e. > 5 mg/day) until at least 8 hours following the last biotin administration.

No interference was observed from rheumatoid factors up to a concentration of 2400 IU/mL and samples from dialysis patients.

In vitro tests were performed on 15 commonly used pharmaceuticals. No interference with the assay was found.

The test cannot be used in patients receiving treatment with lipid-lowering agents containing D-T4. If the thyroid function is to be checked in such patients, the therapy should first be discontinued for 4-6 weeks to allow the physiological state to become re-established.¹¹

Autoantibodies to thyroid hormones can interfere with the assay.

Binding protein anomalies seen with FDH (familial dysalbuminemic hyperthyroxinemia), for example, may cause values which, while characteristic of the condition, deviate from the expected results.¹²

In rare cases, interference due to extremely high titers of antibodies to analyte-specific antibodies, streptavidin or ruthenium can occur. These effects are minimized by suitable test design.

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

Limits and ranges

Measuring range

5.40-320.0 nmol/L or 0.420-24.86 µg/dL (defined by the lower detection limit and the maximum of the master curve). Values below the lower detection limit are reported as < 5.40 nmol/L or < 0.420 µg/dL. Values above the measuring range are reported as > 320.0 nmol/L or > 24.86 µg/dL.

Lower limits of measurement

Lower detection limit of the test

Lower detection limit: 5.40 nmol/L (0.420 µg/dL)

The lower detection limit represents the lowest measurable analyte level that can be distinguished from zero. It is calculated as the value lying two standard deviations above that of the lowest standard (master calibrator, standard 1 + 2 SD, repeatability study, n = 21).

Dilution

Not necessary due to the broad measuring range.

Expected values

Measurements with the Elecsys T4 assay on 2526 serum samples from euthyroid test subjects in Germany and Japan yielded the following values (2.5th-97.5th percentile):

66-181 nmol/L or 5.1-14.1 µg/dL

FT4 Index (T4/TBI) calculated from 825 serum samples from euthyroid test subjects measured with the Elecsys T4 assay and the Elecsys T-Uptake assay (2.5th-97.5th percentile):

62-164 nmol/L or 4.8-12.7 µg/dL

Following values were determined for the 99 % percentile range from 275 serum and plasma samples from healthy test subjects in USA:

59-154 nmol/L or 4.6-12.0 µg/dL

FT4 Index:

57-147 nmol/L or 4.4-11.4 µg/dL

For detailed information about reference intervals in children, adolescents and pregnant women, refer to the brochure "Reference Intervals for Children and Adults", English: [REF 04640292](#), German: [REF 04625889](#).

This booklet also contains results of a detailed study about influencing factors on thyroid parameters in a well characterized reference group of adults. Different inclusion and exclusion criteria were applied (e.g. sonographic results (thyroid volume and density) as well as criteria according to the guidelines of the National Academy of Clinical Biochemistry - NACB).

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

Specific performance data

Representative performance data on the analyzers are given below. Results obtained in individual laboratories may differ.

Precision

Precision was determined using Elecsys reagents, pooled human sera and controls in a modified protocol (EP5-A) of the CLSI (Clinical and Laboratory Standards Institute): 6 times daily for 10 days (n = 60); repeatability on the MODULAR ANALYTICS E170 analyzer, n = 21. The following results were obtained:

cobas e 411 analyzer								
Sample	Mean		Repeatability			Intermediate precision		
			SD		CV	SD		CV
	nmol/L	µg/dL	nmol/L	µg/dL	%	nmol/L	µg/dL	%
HS ^{b)} 1	33.4	2.59	1.56	0.12	4.7	2.31	0.18	6.9
HS 2	123	9.59	3.38	0.26	2.7	4.56	0.35	3.7
HS 3	237	18.4	5.97	0.46	2.5	6.98	0.54	3.0
PC U ^{c)} 1	113	8.79	2.54	0.20	2.3	3.78	0.29	3.3
PC U2	181	14.0	3.58	0.28	2.0	4.90	0.28	2.7

b) HS = human serum

c) PC U = PreciControl Universal

MODULAR ANALYTICS E170, cobas e 601 and cobas e 602 analyzers					
Sample	Repeatability				
	Mean		SD		CV
	nmol/L	µg/dL	nmol/L	µg/dL	%
HS 1	84.3	6.55	1.13	0.09	1.3
HS 2	63.1	4.90	1.14	0.09	1.8
HS 3	243	18.9	4.07	0.32	1.7
PC U1	90.3	7.01	1.15	0.09	1.3
PC U2	182	14.2	3.07	0.24	1.7

MODULAR ANALYTICS E170, cobas e 601 and cobas e 602 analyzers					
Sample	Intermediate precision				
	Mean		SD		CV
	nmol/L	µg/dL	nmol/L	µg/dL	%
HS 1	65.6	5.09	2.40	0.19	3.7
HS 2	79.1	6.15	2.67	0.21	3.4
HS 3	231	18.0	9.67	0.75	4.2
PC U1	92.9	7.22	3.51	0.27	3.8
PC U2	190	14.7	6.29	0.49	3.3

Method comparison

A comparison of the Elecsys T4 assay (y) with the Enzymun-Test T4 method (x) using clinical samples gave the following correlations (nmol/L):
Number of samples measured: 71

Passing/Bablok¹³ Linear regression

$y = 0.77x + 7.77$ $y = 0.75x + 9.88$

$r = 0.841$ $r = 0.975$

The sample concentrations were between 8 and 250 nmol/L (0.6 and 19 µg/dL).

Analytical specificity

For the antibody derivative used, the following cross-reactivities were found:

L-T4 and D-T4 100 %; L-T3 1.53 %; D-T3 1.38 %; 3-iodo-L-tyrosine 0.002 %; 3,5-diiodo-L-tyrosine 0.01 %; 3,3',5,5'-tetraiodothyroacetic acid 38.5 %.

References


- Wheeler MH, Lazarus JH. Diseases of the Thyroid. London, Glasgow, Weinheim, New York, Tokyo, Melbourne, Madras: Chapman and Hall Medical, 1994:108-115.
- Pfannenstiel P, Saller B. Schilddrüsenkrankheiten Diagnose und Therapie. Berliner Medizinische Verlagsanstalt GmbH, 1995;2:43-62,97-106.
- Wenzel KW. Pharmacological interference with in vitro tests of thyroid function. Metabolism 1981;30(7):717-732.
- Burrow GN. Thyroid status in normal pregnancy. J Clin Endocrinol Metab 1990;71:274-275.
- Lazarus JH, Othman S. Thyroid disease in relation to pregnancy. Clin Endocrinol 1991;34:91-98.
- Fisher DA. Physiological variations in thyroid hormones; physiological and pathophysiological considerations. Clinical Chemistry 1996;42:1,135-139.
- Burrow GN, Fisher DA, Larsen PR. Maternal and fetal thyroid function. N Engl J Med 1994;331:1072-1078.
- Nelson JC, Wilcox RB. Analytical performance of free and total thyroxine assays. Clinical Chemistry 1996;42:1,146-154.
- Tietz NW. Clinical Guide To Laboratory Tests. 3rd ed. Philadelphia, Pa: WB Saunders Co 1995:596.
- Thienpont LM, De Brabandere VI, Stöckl D, et al. Development of a New Method for the Determination of Thyroxine in Serum Based on Isotope Dilution Gas Chromatography Mass Spectrometry. Biological Mass Spectrometry 1994;23:475-482.
- Bantle JP, Hunninghake DB, Frantz ID, et al. Comparison of Effectiveness of Thyrotropin-Suppressive Doses of D- and L-Thyroxine in Treatment of Hypercholesterolemia. Am J Med 1984;3:475-481.
- Wada N, Chiba H, Shimizu C, et al. A Novel Missense Mutation in Codon 218 of the Albumin Gene in a Distinct Phenotype of Familial Dysalbuminemic Hyperthyroxinemia in a Japanese Kindred. Journal of Clinical Endocrinology and Metabolism 1997;82(10):3246-3250.
- Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.

For further information, please refer to the appropriate operator's manual for the analyzer concerned, the respective application sheets, the product information and the Method Sheets of all necessary components (if available in your country).

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard (for USA: see <https://usdiagnostics.roche.com> for definition of symbols used):

CONTENT	Contents of kit
SYSTEM	Analyzers/Instruments on which reagents can be used
REAGENT	Reagent
CALIBRATOR	Calibrator
	Volume after reconstitution or mixing
GTIN	Global Trade Item Number

COBAS, COBAS E, ELECSYS and PRECICONTROL are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2016, Roche Diagnostics



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim
www.roche.com

