# Troponin T hs STAT

cobas®

REF		SYSTEM
05092728 190	100	cobas e 411 cobas e 601 cobas e 602

## Lietuvių

### Sistemos informacija

Skirta **cobas e 411** analizatoriui: tyrimo numeris 090

Skirta **cobas e 601** ir **cobas e 602** analizatoriams: Pritaikymo kodo numeris 107

### Paskirtis

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## Nuorodos

- Katus HA, Remppis A, Looser S, et al. Enzyme linked immunoassay of cardiac troponin T for the detection of acute myocardial infarction in patients. Mol Cell Cardiol 1989;21(12):1349-1353.
- Liebetrau C, Mollmann H, Nef H, et al. Release kinetics of biomarkers in patients undergoing transcatheter ablation of septal hypertrophy. Clin Chem 2012;58(6):1049-1054.
- Katus HA, Scheffold T, Remppis A, et al. Proteins of the troponin complex. Laboratory Medicine 1992;23(5):311-317.
- Reichlin T, Hochholzer W, Bassetti S, et al. Early diagnosis of myocardial infarction with sensitive cardiac troponin assays. NEJM 2009;361(9):858-867.
- Giannitsis E, Becker M, Kurz K, et al. High sensitivity cardiac troponin T for early prediction of evolving non-ST segment elevation myocardial infarction in patients with suspected acute coronary syndrome and negative troponin result on presentation. Clin Chem 2010;56(4):642-650.
- Giannitsis E, Kurz K, Hallermayer K, et al. Analytical validation of a high-sensitivity cardiac troponin T assay. Clin Chem 2010;56(2):254-261.
- Hamm CW, Bassand JP, Agewall S, et al. ESC guidelines for the management of acute coronary syndrome in patients presenting without persistent ST-segment elevation. The Task Force for the management of acute coronary syndrome (ACS) in patient presenting without persistent ST-segment elevation of the European Society of Cardiology (ESC). Eur Heart J 2011;32(23):2999-3054.
- NICE (2014) Myocardial infarction (acute): Early rule out using high-sensitivity troponin tests (Elecsys Troponin T high-sensitive, ARCHITECT STAT High Sensitive Troponin-I and AccuTnI-3 assays). NICE diagnostics guidance DG15. Available at [www.nice.org.uk/dg15](http://www.nice.org.uk/dg15) [NICE guideline]
- Roffi M, Patrono C, Collet JP, et al. 2015 ESC guidelines for the management of acute coronary syndromes in patients presenting without persistent ST-segment elevation. Eur Heart J 2016;37(3):267-315.
- Bandstein N, Ljung R, Johansson M, et al. Undetectable high sensitivity cardiac troponin T level in the emergency department and risk of myocardial infarction. J Am Coll Cardiol 2014;63:2569-2578.
- Body R, Burrows G, Carley S, et al. High-sensitivity cardiac troponin T concentrations below the limit of detection to exclude acute myocardial infarction: A prospective evaluation. Clin Chem 2015;61(7):983-989.
- Rubini Gimenez M, Reichlin T, Zellweger C. Rapid rule out of acute myocardial infarction using undetectable levels of high-sensitivity cardiac troponin. Int J Cardiol 2013;168(4):3896-3901.
- Reichlin T, Schindler C, Drexler B, et al. One-hour rule-out and rule-in of acute myocardial infarction using high-sensitivity cardiac troponin T. Arch Intern Med 2012;172(16):1211-1218.
- Reichlin T, Twerenbold R, Wildi K, et al. Prospective validation of a 1-hour algorithm to rule-out and rule-in acute myocardial infarction using a high-sensitivity cardiac troponin T assay. CMAJ 2015;187(8):E243-252.
- Mueller C, Giannitsis E, Christ M, et al. Multicenter evaluation of a 1h-algorithm in the diagnosis of myocardial infarction using high-sensitivity Troponin T: Manuscript submitted for publication;2015.
- Thygesen K, Alpert JS, Jaffe AS, et al. Third definition of myocardial infarction. Eur Heart J 2012;33(20):2551-2567.
- Aviles RJ, Askari AT, Lindahl B, et al. Troponin T levels in patients with acute coronary syndromes, with or without renal dysfunction. N Engl J Med 2002;346:2047-2052.
- Lindahl B, Venge P, James S. The new high-sensitivity cardiac troponin T assay improves risk assessment in acute coronary syndromes. Am Heart J 2010;160:224-229.
- Haaf P, Reichlin T, Twerenbold R, et al. Risk stratification in patients with acute chest pain using three high-sensitivity cardiac troponin assays. Eur Heart J 2014;35(6):365-375.
- Björman C, Larsson M, Johanson P, et al. Small changes in troponin T levels are common in patients with Non-ST-elevation myocardial infarction and are linked to higher mortality. J Am Coll Cardiol 2014;62(14):1231-1238.
- Lindahl B, Venge P, Wallentin L. Troponin T Identifies Patients With Unstable Coronary Artery Disease Who Benefit From Long-Term Antithrombotic Protection. J Am Coll Cardiol 1997;29(1):43-48.
- Hamm CW, Heeschen C, Goldmann B, et al. Benefit of abciximab in patients with refractory unstable angina in relation to serum troponin T levels. N Engl J Med 1999;340(21):1623-1629.

# Troponin T hs STAT

- 23 Heeschen C, Hamm CW, Goldmann BU, et al. for PRISM study investigators. Troponin concentrations for stratification of patients with acute coronary syndromes in relation to therapeutic efficacy of tirofiban. *Lancet* 1999;354:1757-1762.
- 24 Lindahl B, Diderholm E, Lagerquist B, et al. Effects on mortality of long-term treatment with l.m.w. heparin in relation to troponin T level and ECG findings - a FRISC 2 substudy. *Eur Heart J* 2000;21(Suppl.):521.
- 25 Newby LK, Ohman EM, Christenson RH, et al. Benefit of Glycoprotein IIb/IIIa Inhibition in Patients With Acute Coronary Syndromes and Troponin T-Positive Status: The PARAGON-B Troponin T Substudy. *Circulation* 2001;103:2891-2896.
- 26 Wallentin L, Lindholm D, Siegbahn A, et al. Biomarkers in relation to the effects of Ticagrelor compared with clopidogrel in non-ST-elevation acute coronary syndrome patients managed with or without in-hospital revascularization: a Substudy from the prospective randomized Platelet inhibition and Patient Outcomes (PLATO) Trial. *Circulation* 2014;129(3):293-303.
- 27 Amsterdam EA, Wenger NK, Brindis RG, et al. 2014 AHA/ACC guideline for the management of patients with non-ST elevation acute coronary syndromes: a report of the American College of Cardiology/American Heart Association Task Force on Practice Guidelines. *Circulation* 2014;130(25):e344-426
- 28 Reichlin T, Irfan A, Twerenbold R, et al. Utility of absolute and relative changes in cardiac troponin concentrations in the early diagnosis of acute myocardial infarction. *Circulation* 2011;124(2):136-145.
- 29 Masson S, Anand I, Favero C, et al. Serial measurement of cardiac troponin T using a highly sensitive assay in patients with chronic heart failure: data from two large randomized clinical trials. *Circulation* 2012;125(2):280-288.
- 30 Nambi V, Liu X, Chambless LE, et al. Troponin T and N-terminal pro-B-type natriuretic peptide: a biomarker approach to predict heart failure risk: the Atherosclerosis Risk in Communities Study. *Clin Chem* 2013;59(12):1802-1810.
- 31 Yancy CW, Jessup M, Bozkurt B, et al. 2013 ACCF/AHA guideline for the management of heart failure: a report of the American College of Cardiology Foundation/American Heart Association Task Force on practice guidelines. *Circulation* 2013;128(16):e240-327.
- 32 Cramer G, Bakker J, Gommans F, et al. Relation of highly sensitive cardiac troponin T in hypertrophic cardiomyopathy to left ventricular mass and cardiovascular risk. *Am J Cardiol*. 2014;113(7):1240-1245.
- 33 McGill D, Talaulikar G, Potter JM, et al. Over time, high-sensitivity TnT replaces NT-proBNP as the most powerful predictor of death in patients with dialysis-dependent chronic renal failure. *Clin Chim Acta* 2010;411(13-14):936-939.
- 34 K/DOQI clinical practice guidelines for cardiovascular disease in dialysis patients. *Am J Kidney Dis*. 2005;45 (4 Suppl 3):S1-153.
- 35 Artunc F, Mueller C, Breidhardt T, et al. Sensitive troponins- which suits better for hemodialysis patients? Associated factors and prediction for mortality. *PLoS One*. 2012;7(10):e47610.
- 36 Wolley M, Steward R, Curry E, et al. Variation in and prognostic importance of troponin T measured using a high-sensitivity assay in clinically stable haemodialysis patients. *Clin Kidney J*. 2013;6(4):402-409.
- 37 Honneger Bloch S, Semple D, Sidhu K, et al. Prognostic value and long-term variation of high sensitivity troponin T in clinically stable haemodialysis patients. *N Z Med J*. 2014;127(1402):97-109.
- 38 Twerenbold R, Wildi K, Jeger C, et al. Optimal cutoff levels of more sensitive cardiac troponin assays for the early diagnosis of myocardial infarction in patients with renal dysfunction. *Circulation* 2015;131(23):2041-2050.
- 39 Landeberg G, Jaffe AS, Gilon D, et al. Troponin elevation in severe sepsis and septic shock: the role of left ventricular dysfunction and right ventricular dilatation. *Crit Care Med* 2014;42(4):790-800.
- 40 Hillis GS, Welsh P, Chalmers J, et al. The relative and combined ability of high sensitivity cardiac troponin T and N-terminal pro-BNP to predict cardiovascular events and death in patients with type 2 diabetes mellitus. *Diabetes Care* 2014;37(1):295-303.
- 41 Everett BM, Brooks MM, Vlachos HE, et al. Troponin and cardiac events in stable ischemic heart disease and diabetes. *N Engl J Med*. 2015;373(7):610-620.
- 42 European patent application 890154 by F. Hoffmann-La Roche AG and Roche Diagnostics GmbH. Cardiac Troponin as an indicator of advanced coronary artery disease.
- 43 European patent application 1837659 by F. Hoffmann-La Roche AG and Roche Diagnostics GmbH. Means and methods for the differentiation of acute and chronic myocardial necrosis in symptomatic patients.
- 44 Latini R, Masson S, Anand IS, et al. Prognostic Value of Very Low Plasma Concentrations of Troponin T in Patients with Stable Chronic Heart Failure. *Circulation* 2007;116:1242-1249.
- 45 Omland T, de Lemos JA, Sabatine MS, et al. A sensitive cardiac troponin T assay in stable coronary artery disease. *N Engl J Med* 2009;361:2538-2547.
- 46 Latini R, Masson S, Pirelli S, et al. On the behalf of the GISSI-AF Investigators. Circulating cardiovascular biomarkers in recurrent atrial fibrillation: data from the GISSI-Atrial fibrillation trial. *J Intern Med* 2011;269(2):160-171.
- 47 Lewandowski K. Special topics: cardiac markers in myocarditis. Cardiac transplant rejection and conditions other than acute coronary syndrome. *Clin Lab Med* 2014;34:129-135.
- 48 Swaanenburg JC, Klaase JM, DeJongste M, et al. Troponin I, troponin T, CK-MB-activity and CK-MB mass as markers for the detection of myocardial contusion in patients who experienced blunt trauma. *Clin Chim Acta* 1998;272:171-181.
- 49 Bajaj A, Saleeb M, Rathor P, et al. Prognostic value of troponins in acute noninvasive pulmonary embolism. A meta-analysis. *Heart Lung* 2010;44(4):327-334.
- 50 Newby LK, Rodriguez I, Finkle J, et al. Troponin measurements during drug development-considerations for monitoring and management of potential cardiotoxicity. An educational collaboration among the Cardiac Safety Research Consortium, the Duke Clinical Research Institute, and the US Food and Drug Administration. *Am Heart J* 2011;162(1):64-73.
- 51 European patent 1363128 by F Hoffmann-La Roche AG and Roche Diagnostics GmbH. Diagnosis of myocardial infarction and acute coronary syndrome by combination of markers.
- 52 European patent application 1925943 by F. Hoffmann-La Roche AG and Roche Diagnostics GmbH. Means and methods for optimization of diagnostic and therapeutic approaches in chronic artery disease based on the detection of Troponin T and NT-pro BNP.
- 53 Bosselmann H, Egstrup M, Rossing K, et al. Prognostic significance of cardiovascular biomarkers and renal dysfunction in outpatients with systolic heart failure: a long term follow-up study *Int J Cardiol* 2013;170(2):202-207.
- 54 Dubin RF, Li Y, He J, et al. Predictors of high sensitivity cardiac troponin T in chronic kidney disease patients: a cross-sectional study in the chronic renal insufficiency cohort (CRIC). *BMC Nephrol* 2013;14(1):229.
- 55 Davis GK, Labugger R, Van Eyk JE, et al. Cardiac troponin T is not detected in western blots of diseased renal tissue. *Clin Chem* 2001;47:782-783.
- 56 Ricchiuti V, Voss EM, Ney A, et al. Cardiac Troponin T isoforms expressed in renal diseased skeletal muscle will not cause false positive results by the second generation cTnT assay by Boehringer Mannheim. *Clin Chem* 1998;44(9):1919-1924.
- 57 Hallermayer K, Klenner D, Vogel R. Use of recombinant human cardiac troponin T for standardization of third generation troponin T methods. *Scand J Clin Invest* 1999;59(Suppl 230):128-131.
- 58 Saenger AK, Beyrau R, Braun S, et al. Multicenter analytical evaluation of a high-sensitivity troponin T assay. *Clin Chim Acta* 2011;412(9-10):748-754.
- 59 World Health Organization. Report of the Joint International Society and Federation of Cardiology/World Health Organization Task Force on Standardization of Clinical Nomenclature. Nomenclature and criteria for diagnosis of ischemic heart disease. *Circulation* 1979;59:607-609.

# Troponin T hs STAT

- 60 Müller-Bardorff M, Hallermayer K, Schröder A, et al. Improved troponin T ELISA specific for cardiac troponin T isoform: assay development and analytical validation. Clin Chem 1997;43(3):458-466.
- 61 Klein G, Kampmann M, Baum H, et al. Clinical performance of the new cardiac markers troponin T and CK-MB on the Elecsys 2010. A multicenter evaluation. Wien Klin Wochenschr. 1998;110(Suppl3):40-51.
- 62 Mendis S, Thygesen K, Kuulasmaa K, et al. World Health Organization definition of myocardial infarction: 2008-09 revision. Int J Epidemiol 2011;40(1):139-146.
- 63 Ferjani M, Droc G, Dreux S, et al. Circulating cardiac troponin T in myocardial contusion. Chest 1997;111(2):427-433.
- 64 Erbel C, Taskin R, Doesch A, et al. High-sensitive troponin T measurements early after heart transplantation predict short- and long-term survival. Transpl Int 2013;26(3):267-272.
- 65 Li SF, Zapata J, Tillem E. The prevalence of false-positive cardiac troponin I in ED patients with rhabdomyolysis. Am J Emerg Med 2005;23(7):860-863.
- 66 Zhang L, Wang GC, Ma L, et al. Cardiac involvement in adult polymyositis or dermatomyositis: a systematic review Clin Cardiol 2012;35(11):686-691.
- 67 Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.
- 68 World Health Organization. Report of the Joint International Society and Federation of Cardiology/World Health Organization Task Force on Standardization of Clinical Nomenclature. Nomenclature and criteria for diagnosis of ischemic heart disease. Circulation 1979;59:607-609.
- 69 The Task Force for the diagnosis and Treatment of Non-ST-Segment Elevation Acute Coronary Syndromes of the European Society of Cardiology. Guidelines for the diagnosis and treatment of non-ST elevation acute coronary syndromes. European Heart Journal 2007;28:1598-1660.

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© 2016, Roche Diagnostics
Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim