



AirSpiral™ for less condensate\*



## F&P AIRVO™ 2 | FEATURES

**We've all seen it happen before.** When warm, humid gas touches anything even a fraction of a degree colder than itself, it cools and generates condensate. We see it often in conventional breathing tubes made of thin, dense plastic: their surfaces are easily cooled by the outside air, resulting in unpleasant “rain-out”.

**The new AirSpiral™ tube for AIRVO 2 utilizes an advanced two-spiral design to help protect against condensate:**

**The heating spiral** contains heater wire to warm the breathing gas inside the tube.


**The insulating spiral** uses a bubble of warmed air around the breathing gas to shield it from the outside environment.

Like a double-paned window, the inside and outside surfaces of the tube are physically separated, to reduce cooling and helps maximize the level of humidity remaining in the gas all the way to the patient.


AirSpiral tube produces up to **93% less condensate\*** than a previous AIRVO tube.

Less condensate may reduce interruptions to therapy.

**AirSpiral tube**



**Conventional breathing tubes**



## BENEFITS

### Tubing clip

- The moulded clip attaches to clothes or bedding to allow positioning of the tube during use.



### 2-in-1 tube

- Connects with both Optiflow™ + and Optiflow™ Junior interfaces

### Integrated temperature sensor

- No external probes, cables or adaptors required.

### Ergonomically designed

- Moulded grips at both ends of the tube assist connection and disconnection.



### Product codes

<b>900PT561</b>	AirSpiral tube and chamber kit for AIRVO (10-pack)
<b>900PT560</b>	AirSpiral tube for myAIRVO (10 pack)
<b>900PT560E</b>	AirSpiral tube for myAIRVO (1 pack)

\* When compared with AIRVO tube 900PT501 in internal F&P testing.

F&P, AIRVO, Optiflow and AirSpiral are all trademarks of Fisher & Paykel Healthcare Ltd. For patent information, see [www.fphcare.com/ip](http://www.fphcare.com/ip)

For further information about the AirSpiral tube, including warnings and precautions, please read the user instructions.

[www.fphcare.com/airvo](http://www.fphcare.com/airvo)

**Fisher & Paykel**  
HEALTHCARE