**VIDAS<sup>®</sup> TOXO IgM (TXM)**

VIDAS TOXO IgM is an automated qualitative test for use on the VIDAS family instruments, for the detection of anti-toxoplasma IgM in serum using the ELFA technique (Enzyme Linked Fluorescent Assay).

**SUMMARY AND EXPLANATION**

*Toxoplasma gondii*, an obligate intracellular protozoan parasite, is a significant pathogen in humans. The parasite, whose primary host is the Felidae, is scattered in nature and infects numerous other mammals (1).

Toxoplasmosis is usually benign or asymptomatic in immunocompetent subjects, but can have severe consequences if it occurs in immunodeficient subjects or fetuses. Congenital toxoplasmosis connected with maternal contamination prior to conception is exceptional and most often occurs in immunocompromised women (11). Women who are seronegative are at risk of becoming infected during gestation. The possible ensuing fetal transmission is due to the transplacental migration of toxoplasma during the acute phase of the infection. The frequency and severity of the fetal infection will depend on a number of factors including the virulence of the infecting strain, the size of the original inoculum, the immune response of the patient and when, during gestation, the mother became infected (2-7).

The diagnosis of *T. gondii* infection is most commonly made by biological examination: specific immunoglobulin detection (IgM and IgG) (1, 4, 6).

The diagnosis of an acute acquired infection during pregnancy is established by demonstration of a seroconversion, or by a significant rise in antibody titer in two sequential sera assayed concomitantly (7-10).

VIDAS TOXO IgM is intended for use as an aid in the determination of the immune status of patients.

**PRINCIPLE**

The assay principle combines an enzyme immunoassay method by immunocapture with a final fluorescent detection (ELFA).

The Solid Phase Receptacle (SPR<sup>®</sup>) serves as the solid phase as well as the pipetting device for the assay. Reagents for the assay are ready-to-use and pre-dispensed in the sealed reagent strips.

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After a sample dilution step, the IgM are captured by the polyclonal Ab coating the interior of the SPR. The anti-toxoplasma IgM are specifically detected by inactivated toxoplasma antigen (RH Sabin strain), which is itself revealed by an alkaline phosphatase-labeled murine monoclonal anti-toxoplasma antibody (anti-P30).

During the final detection step, the substrate (4-Methyl-umbelliferyl phosphate) is cycled in and out of the SPR. The conjugate enzyme catalyzes the hydrolysis of this substrate into a fluorescent product (4-Methyl-umbelliferone) the fluorescence of which is measured at 450 nm. The intensity of the fluorescence is proportional to the concentration of antibodies present in the sample. At the end of the assay, an index is automatically calculated by the instrument in relation to the S1 standard stored in memory, and then printed out.

**CONTENT OF THE KIT (60 TESTS):**

60 TXM strips	STR	Ready-to-use.
60 TXM SPRs 2 x 30	SPR	Ready-to-use. Interior of SPRs coated with anti-human $\mu$ chain antibodies (goat).
Positive control TXM 1 x 2 ml (liquid)	C1	Human serum* containing anti-Toxoplasma IgM + protein stabilizer + 1 g/l sodium azide. MLE data indicate the index: confidence interval ("Control C1 (+) Test Value Range").
Negative control TXM 1 x 2 ml (liquid)	C2	Human serum* negative for anti-Toxoplasma IgM + protein stabilizer + 1 g/l sodium azide.
Standard TXM 1 x 1 ml (liquid)	S1	Human serum* containing anti-Toxoplasma IgM + protein stabilizer + 1 g/l sodium azide.
Specifications for the factory master data required to calibrate the test:		
<ul style="list-style-type: none"> <li>MLE data (Master Lot Entry) provided in the kit,</li> </ul> or <ul style="list-style-type: none"> <li>MLE bar codes printed on the box label.</li> </ul>		
1 Package Insert provided in the kit or downloadable from <a href="http://www.biomerieux.com/techlib">www.biomerieux.com/techlib</a> .		

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### The SPR

The interior of the SPR® is coated during production with anti-human  $\mu$  chain antibodies (goat). Each SPR is identified by the TXM code. Only remove the required number of SPRs from the pouch and **carefully reseal the pouch after opening**.

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### Description of the TXM strip

1	Sample well.
2	Serum diluent: TRIS buffer (50 mmol/l) pH 7.4 + protein and chemical stabilizers + 0.9 g/l sodium azide (300 $\mu$ l).
3	Pre-wash buffer: TRIS (50 mmol/l) pH 7.4 + protein and chemical stabilizers + 0.9 g/l sodium azide (600 $\mu$ l).
4 - 5 - 7 - 8	Wash buffer: TRIS (50 mmol/l) pH 7.4 + protein and chemical stabilizers + 0.9 g/l sodium azide (600 $\mu$ l).
6	Conjugate: alkaline phosphatase-labeled immune complex (toxoplasma antigen RH Sabin strain grown in mice (12) - mouse monoclonal anti-P30 antibodies) + 0.9 g/l sodium azide + 0.02% gentamicin (400 $\mu$ l).
9	Empty well
10	Cuvette with substrate: 4-Methyl-umbelliferyl-phosphate (0.6 mmol/l) + diethanolamine (DEA*) (0.62 mol/l or 6.6%, pH 9.2) + 1 g/l sodium azide (300 $\mu$ l).

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### MATERIALS AND DISPOSABLES REQUIRED BUT NOT PROVIDED

- Pipette with disposable tip to dispense 100  $\mu$ l.
- Powderless, disposable gloves.
- For other specific materials and disposables, please refer to the Instrument User's Manual.
- Instrument of the VIDAS family.

### WARNINGS AND PRECAUTIONS

- **For *in vitro* diagnostic use only.**
- **For professional use only.**
- **This kit contains products of human origin. No known analysis method can totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (see Laboratory Biosafety Manual - WHO - Geneva - latest Edition).**
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents.

It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (do not ingest or inhale).

- Do not use the SPRs if the pouch is pierced.
- Do not use visibly deteriorated STRs (damaged foil or plastic).
- Do not use reagents after the expiration date indicated on the label.
- Do not mix reagents (or disposables) from different lots.
- Use **powderless** gloves, as powder has been reported to cause false results for certain enzyme immunoassay tests.
- Kit reagents contain sodium azide which can react with lead or copper plumbing to form explosive metal azides. If any liquid containing sodium azide is disposed of in the plumbing system, drains should be flushed with water to avoid build-up.
- The substrate in well 10 contains an irritant agent (6.6% diethanolamine). Refer to the hazard statements "H" and the precautionary statements "P" above.

- Spills should be wiped up thoroughly after treatment with liquid detergent or a solution of household bleach containing at least 0.5% sodium hypochlorite. See the User's Manual for cleaning spills on or in the instrument. Do not autoclave solutions containing bleach.
- The instrument should be regularly cleaned and decontaminated (see the User's Manual).

### STORAGE CONDITIONS

- Store the VIDAS TXM kit at 2-8°C.
- **Do not freeze reagents.**
- **Store all unused reagents at 2-8°C.**
- After opening the kit, check that the SPR® pouch is correctly sealed and undamaged. If not, do not use the SPRs.
- **Carefully reseal the pouch with the desiccant inside after use to maintain stability of the SPRs and return the complete kit to 2-8°C.**
- If stored according to the recommended conditions, all components are stable until the expiration date indicated on the label.

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The use of cord blood or neonatal sera has not been validated. For these specimens, we recommend the use of the Toxo-ISAGA technique reference 75 361.

**Since the use of hemolyzed, icteric or lipemic samples has not been validated; it is recommended to collect a new sample.**

Heat inactivated sera (56°C for 30 minutes) may be tested with VIDAS TOXO IgM.

#### Specimen stability

Samples can be stored at 2-8°C in stoppered tubes for up to 7 days; if longer storage is required, freeze the sera at -25 ± 6 °C.

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### Procedure

1. **Only remove the required reagents from the refrigerator and allow them to come to room temperature for at least 30 minutes.**
  2. Use one "TXM" strip and one "TXM" SPR for each sample, control or standard to be tested. **Make sure the storage pouch has been carefully resealed after the required SPRs have been removed.**
  3. The test is identified by the "TXM" code on the instrument. The standard must be identified by "S1", and tested in **duplicate**. If the positive control is to be tested, it should be identified by "C1". If the negative control needs to be tested, it should be identified by "C2".
  4. If necessary, clarify samples by centrifugation.
  5. Mix the standard, controls and samples using a vortex-type mixer (for serum separated from the pellet).
6. **For this test, the standard, control, and sample test portion is 100 µl.**
7. Insert the "TXM" SPRs and "TXM" strips into the instrument. Check to make sure the color labels with the assay code on the SPRs and the Reagent Strips match.
  8. Initiate the assay as directed in the User's Manual. All the assay steps are performed automatically by the instrument.
  9. Restopper the vials and return them to 2-8°C after pipetting.
  10. The assay will be completed within approximately **40 minutes**. After the assay is completed, remove the SPRs and strips from the instrument.
  11. Dispose of the used SPRs and strips into an appropriate recipient.

### RESULTS AND INTERPRETATION

Once the assay is completed, results are analyzed automatically by the computer. Fluorescence is measured twice in the Reagent Strip's reading cuvette for each sample tested. The first reading is a background reading of the substrate cuvette before the SPR is introduced into the substrate. The second reading is taken after incubating the substrate with the enzyme remaining on the interior of the SPR. The RFV (Relative Fluorescence Value) is calculated by subtracting the background reading from the final result. This calculation appears on the result sheet.

The instrument calculates a test value for each sample (index), which is the ratio between its RFV and that of the memorized standard.

Interpretation of test results should be made taking into consideration the patient's history, and the results of any other tests performed or other IgM assay methods.

As no international standard is available for the determination of anti-toxoplasma IgM, the VIDAS TOXO IgM reagent is calibrated against collection sera.

**Thresholds and interpretation of results**

Index	Interpretation
$i < 0.55$	Negative
$0.55 \leq i < 0.65$	Equivocal
$i \geq 0.65$	Positive

Equivocal samples (indices between 0.55 and 0.65) should be retested. If the interpretation remains equivocal, a new sample must be collected.

**QUALITY CONTROL**

A positive and a negative control are included in each VIDAS TOXO IgM kit.

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**LIMITATIONS OF THE METHOD**

Interference may be encountered with certain sera containing antibodies directed against reagent components. For this reason, assay results should be interpreted taking into consideration the patient's history, and the results of any other tests performed.

In children under 6 months old, production of IgM is sometimes impossible to detect, in which case a significant increase in specific IgG titer must be noted.

In AIDS infected patients, the presence of IgG associated with results from radiological examinations are useful in diagnosis, since the detection of an increase in IgG or the presence of IgM is unfortunately not very reliable. In nearly 100% of AIDS infected patients who develop toxoplasmosis, IgG are detected in the serum, whereas IgM are seldom detected.

Since the persistence of anti-toxoplasma IgM for several months or several years after seroconversion is a regularly described phenomenon, it is not possible to confirm recent infection, based solely on the presence of anti-toxoplasma IgM, without taking into account the serological history of the patient.

**PERFORMANCE**

Studies performed using VIDAS TOXO IgM gave the following results:

**Precision**Within-run reproducibility:

3 samples were tested 20 times in 3 different runs on the same VIDAS instrument.

Sera	Mean index	Standard deviation	Intra-assay CV %
S1	0.96	0.04	3.8
C1	1.39	0.03	2.0
C2	0.09	0.01	11.5

Between-run reproducibility

3 samples were tested singly in 7 different runs on the same VIDAS instrument.

Sera	Mean index	Standard deviation	Inter-assay CV %
S4	2.09	0.10	4.8
S5	0.76	0.04	5.3
S6	0.11	0.01	8.2

**Relative specificity**

The relative specificity of VIDAS TOXO IgM was calculated in comparison with the Toxo-ISAGA IgM technique Ref. 75361 using unselected fresh sera samples collected from an immunized and non-immunized population, sera collected prior to contamination from pregnant women were in contact with toxoplasmosis, cord blood of neonates not contaminated with toxoplasmosis and sera from children suffering from congenital toxoplasmosis but with no anti-toxoplasma IgM.

Out of the 1,483 sera found to be negative with the Toxo-ISAGA IgM kit, 7 were equivocal and 1,465 were negative with the VIDAS TOXO IgM kit, resulting in a relative specificity of **99.25%\*** [1,465/1,476] (95% confidence interval: 98.66-99.59%).

Out of the 11 samples found to be positive with the VIDAS TOXO IgM kit and negative with the Toxo-ISAGA IgM kit, 4 were positive with another anti-toxoplasma IgM detection technique and 2 corresponded to follow-up sera collected more than 8 weeks after infection (residual IgM).

\*The equivocal results were not used to calculate the performance.

**Relative sensitivity**

The relative sensitivity of VIDAS TOXO IgM was calculated in comparison with the Toxo-ISAGA IgM technique Ref. 75361 using sera from pregnant women who were in contact with toxoplasmosis during their pregnancy and sera from neonates suffering from congenital toxoplasmosis.

A total of 77 biologically and clinically documented cases of seroconversion and 24 cases of congenital infections were studied.

**Acquired toxoplasmosis:**

Out of the 150 seroconversion sera found to be positive with Toxo-ISAGA IgM, 144 were positive with VIDAS TOXO IgM. No equivocal samples were found.

The relative sensitivity of sera from pregnant women who were in contact with toxoplasmosis during their pregnancy is **96.00%** (95% confidence interval: 91.43-98.18%).

Out of the 6 discrepant sera:

- Four corresponded to follow-up sera collected more than 8 weeks after the infection (residual IgM). Two of these sera were also found to be negative with two other anti-toxoplasma IgM detection techniques.

- Two corresponded to specimens collected shortly after contamination: absence or low levels of anti-toxoplasma IgG (negative or weakly positive Lyse Test) and low levels of IgM (weakly positive Toxo-ISAGA IgM result). In both these cases, seroconversion was detected with the specimens collected later.

### Congenital toxoplasmosis:

Out of the 30 congenital toxoplasmosis sera found to be positive with the Toxo-ISAGA IgM technique, 27 were positive with VIDAS TOXO IgM. No equivocal samples were found.

The relative sensitivity in comparison with the Toxo-ISAGA IgM technique using sera from neonates suffering from congenital toxoplasmosis is **90.00%** (95% confidence interval: 73.47-97.89%). The positive predictive value for this population is 100% (95% confidence interval: 87.23-100%).

Out of the 3 discrepant sera, one was collected at birth from cord blood. The other two were collected at one and three months respectively. One serum showed low levels of toxoplasma IgM (weakly positive Toxo-ISAGA IgM result).

**In view of the results obtained for cases of congenital toxoplasmosis, it is recommended that all sera from neonates with congenital toxoplasmosis found to be negative with VIDAS TOXO IgM, be confirmed using the Toxo-ISAGA IgM technique.**

### CROSS REACTIVITY AND RELEVANT INTERFERENTS

During the evaluations, 24 samples containing rheumatoid factor negative for anti-toxoplasma IgM and 10 samples containing anti-EBV IgM negative for anti-toxoplasma IgM were tested. All of these samples were found to be negative with VIDAS TOXO IgM.

### PREVALENCE

*T. gondii* is a strict pathogen whose prevalence differs greatly from one country to another, or even from one region to another.

Contamination by *T. gondii* can vary depending on cultural customs and eating habits, resulting in a prevalence ranging from less than 10% in certain regions of Northern Europe to more than 90% in Africa.

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**Note:***Minor typographical, grammar, and formatting changes are not included in the revision history.*

Release date	Part Number	Change Type	Change Summary
2015/01	05933R	Administrative	INDEX OF SYMBOLS REVISION HISTORY
		Technical	CONTENT OF THE KIT (60 TESTS) WARNINGS AND PRECAUTIONS INSTRUCTIONS FOR USE

BIOMERIEUX, the blue logo, SPR and VIDAS are used, pending and/or registered trademarks belonging to bioMérieux, or one of its subsidiaries, or one of its companies.

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**VIDAS® TOXO IgG II (TXG)**

VIDAS TOXO IgG II is an automated quantitative test for use on the VIDAS family instruments for the quantitative measurement of anti-toxoplasma IgG in human serum or plasma (lithium heparin or EDTA) using the ELFA technique (Enzyme Linked Fluorescent Assay).

**SUMMARY AND EXPLANATION**

*Toxoplasma gondii*, (1) an obligate intracellular protozoan parasite, is a significant pathogen among humans. The parasite, whose primary host is the Felidae, is scattered in nature and invades all orders of mammals.

Toxoplasmosis is usually benign or asymptomatic, but can have severe consequences if it occurs in immunodeficient subjects or fetuses (2). Women who are seropositive before they become pregnant are essentially protected from transmitting the infection to their unborn child (3, 4, 5, 6, 7). Women who are seronegative are at risk of becoming infected during gestation. Transmission to the fetus occurs during the acute stage of the infection in the mother. The frequency with which *Toxoplasma* crosses the placental barrier to cause congenital infection will depend on a number of factors including the virulence of the infecting strain of the parasite, the size of the original inoculum, the immune response of the patient and when, during gestation, the mother became infected (8).

The prevalence of toxoplasmosis varies depending upon geographical location, age and gender of the population studied and other factors. In Europe, the prevalence rate ranges from 20% to 85%. In the United States, the prevalence is lower: 12% to 41%. Prevalence in other countries can vary from 18% to 65%.

The diagnosis of *T.gondii* infection is most commonly made by biological examination: specific immunoglobulin detection (IgM and IgG)(9, 10).

The diagnosis of the acute acquired infection during pregnancy is established by demonstration of a seroconversion, or by a significant rise in antibody titer in sequential sera assayed concomitantly (6, 11).

VIDAS TOXO IgG II is a screening test intended for use as an aid in the determination of the immunity status of patients.

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The Solid Phase Receptacle (SPR®) serves as the solid phase as well as the pipetting device for the assay. Reagents for the assay are ready-to-use and pre-dispensed in the sealed reagent strips.

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After a sample dilution step, the sample is cycled in and out of the SPR. Anti-*T.gondii* IgG antibodies present in the specimen will bind to the *T.gondii* antigen coating the interior of the SPR. Unbound components are eliminated during the washing steps. Mouse monoclonal anti-human IgG conjugated with alkaline phosphatase is cycled through the SPR and will attach to any human IgG bound to the SPR wall. During the final detection step, the substrate (4-Methyl-umbelliferyl phosphate) is cycled in and out of the SPR. The conjugate enzyme catalyzes the hydrolysis of this substrate into a fluorescent product (4-Methyl-umbelliferone) the fluorescence of which is measured at 450 nm. The intensity of the fluorescence is proportional to the concentration of antibodies present in the sample. At the end of the assay, results are automatically calculated by the instrument in relation to the calibration curve stored in memory, and then printed out.

**CONTENT OF THE KIT (60 TESTS):**

60 TXG strips	STR	Ready-to-use.
60 TXG SPRs 2 x 30	SPR®	Ready-to-use. SPRs coated with membrane and cytoplasmic Toxoplasma antigen, RH Sabin strain grown in mice (12).
TXG positive control 1 x 2 mL (liquid)	C1	Ready-to-use. Human serum* containing anti-Toxoplasma IgG + protein stabilizer + 1 g/L sodium azide. MLE data indicate the confidence interval in IU/mL (international unit per milliliter) ("Control C1 (+) Dose Value Range").
TXG negative control 1 x 3 mL (liquid)	C2	Ready-to-use. Human serum* negative for anti-Toxoplasma IgG + protein stabilizer + 1 g/L sodium azide.
TXG calibrator 1 x 1 mL (liquid)	S1	Ready-to-use. Human serum* containing anti-Toxoplasma IgG and calibrated against the 2 <sup>nd</sup> WHO International Standard + protein stabilizer + 1 g/L sodium azide. MLE data indicate the concentration in IU/mL ("Calibrator (S1) Dose Value") and the confidence interval in "Relative Fluorescence Value ("Calibrator (S1) RFV Range").
Specifications for the factory master data required to calibrate the test: • MLE data (Master Lot Entry) provided in the kit, or • MLE bar code printed on the box label.		
1 Package insert provided in the kit or downloadable from <a href="http://www.biomerieux.com/techlib">www.biomerieux.com/techlib</a>		

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**The SPR®**

The interior of the SPR is coated during production with membrane and cytoplasmic Toxoplasma antigen (RH Sabin strain). Each SPR is identified by the TXG code. Only remove the required number of SPRs from the pouch and **carefully reseal the pouch after opening**.

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**Description of the TXG strip**

Wells	Reagents
1	Sample well.
2	Serum diluent: TRIS buffer (50 mmol/l) pH 7.4 + protein and chemical stabilizers + 0.9 g/l sodium azide (600 µl).
3	Pre-wash buffer: TRIS (50 mmol/l) pH 7.4 + protein and chemical stabilizers + 0.9 g/l sodium azide (600 µl).
4 - 5 - 7 - 8	Wash buffer: TRIS (50 mmol/l) pH 7.4 + protein and chemical stabilizers + 0.9 g/l sodium azide (600 µl).
6	Conjugate: Alkaline phosphatase-labeled monoclonal anti-human IgG antibodies (mouse) + 0.9 g/l sodium azide (400 µl).
9	Serum diluent: TRIS buffer (50 mmol/l) pH 7.4 + protein and chemical stabilizers + 0.9 g/l sodium azide (400 µl).
10	Cuvette with substrate: 4-Methyl-umbelliferyl phosphate (0.6 mmol/l) + diethanolamine (DEA*) (0.62 mol/l or 6.6%, pH 9.2) + 1 g/l sodium azide (300 µl).

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- Pipette with disposable tip to dispense 100 µl.
- Powderless, disposable gloves.
- For other specific materials and disposables, please refer to the Instrument User's Manual.
- Instrument of the VIDAS® family.

**WARNINGS AND PRECAUTIONS**

- **For *in vitro* diagnostic use only.**
- **For professional use only.**
- **This kit contains products of human origin. No known analysis method can totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (see Laboratory biosafety manual - WHO - Geneva - latest edition).**
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of

transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (do not ingest or inhale).

- Do not use the SPRs if the pouch is pierced.
- Do not use visibly deteriorated STRs (damaged foil or plastic).
- Do not use reagents after the expiration date indicated on the label.
- Do not mix reagents (or disposables) from different lots.
- Use **powderless** gloves, as powder has been reported to cause false results for certain enzyme immunoassay tests.

- Kit reagents contain sodium azide which can react with lead or copper plumbing to form explosive metal azides. If any liquid containing sodium azide is disposed of in the plumbing system, drains should be flushed with water to avoid build-up.
- The optical cuvette with substrate (well 10) contains an irritant agent (6.6% diethanolamine). Refer to the hazard statements "H" and the precautionary statements "P" above.
- Spills should be wiped up thoroughly after treatment with liquid detergent and a solution of household bleach containing at least 0.5% sodium hypochlorite. See the User's Manual for cleaning spills on or in the instrument. Do not autoclave solutions containing bleach.
- The instrument should be regularly cleaned and decontaminated (see the User's Manual).

### STORAGE CONDITIONS

- Store the VIDAS TOXO IgG II kit at 2-8°C.
- **Do not freeze reagents.**
- **Store all unused reagents at 2-8°C.**
- After opening the kit, check that the SPR pouch is correctly sealed and undamaged. If not, do not use the SPRs.
- **Carefully reseal the pouch with the desiccant inside after use to maintain stability of the SPRs and return the complete kit to 2-8°C.**
- If stored according to the recommended conditions, all components are stable until the expiration date indicated on the label.

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Since the use of hemolyzed, lipemic, or icteric samples has not been validated, it is recommended to collect a new specimen.

Sera inactivated at 56°C for 30 minutes can be tested with VIDAS TOXO IgG II.

#### Specimen stability

Samples can be stored at 2-8°C in stoppered tubes for up to 5 days; if longer storage is required, freeze the sera or plasma at -25 ± 6°C.

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A study performed on frozen samples over a period of one year, showed that the quality of results is not affected.

A study performed on samples frozen for 12 months, showed that the quality of results is not affected.

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### Calibration

Calibration, using the calibrator provided in the kit, must be performed upon receipt of a new lot of reagents after the master lot data have been entered. Calibration should then be performed every 14 days. This operation provides instrument-specific calibration curves and compensates for possible minor variations in assay signal throughout the shelf-life of the kit.

The calibrator, identified by S1, must be tested in **duplicate** (see User's Manual). The calibrator value must be within the set RFV "Relative Fluorescence Value" range. If this is not the case, recalibrate

### Procedure

1. **Only remove the required reagents from the refrigerator and allow them to come to room temperature for at least 30 minutes.**
  2. Use one "TXG" strip and one "TXG" SPR for each sample, control or calibrator to be tested. **Make sure the storage pouch has been carefully resealed after the required SPRs have been removed.**
  3. The test is identified by the "TXG" code on the instrument. The calibrator must be identified by "S1", and tested in duplicate. If the positive control is to be tested, it should be identified by "C1". If the negative control needs to be tested, it should be identified by C2.
  4. Mix the calibrator, controls and samples using a vortex-type mixer (for serum or plasma separated from the pellet).
5. **For this test, the calibrator, control, and sample test portion is 100 µl.**
6. Insert the "TXG" SPRs and "TXG" strips into the instrument. Check to make sure the color labels with the assay code on the SPRs and the Reagent Strips match.
  7. Initiate the assay as directed in the User's Manual. All the assay steps are performed automatically by the instrument.
  8. Restopper the vials and return them to 2-8°C after pipetting.
  9. The assay will be completed within approximately 40 minutes. After the assay is completed, remove the SPRs and strips from the instrument.
  10. Dispose of the used SPRs and strips into an appropriate recipient.

**NOTE: If a comparison needs to be made with a previously collected specimen, test both sera in the same run.**

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The results are automatically calculated by the instrument using calibration curves which are stored by the instrument (4-parameter logistic model) and are expressed in « IU/ml » (2<sup>nd</sup> WHO International Standard). Equivocal samples should be retested. If the interpretation remains equivocal, a new sample must be collected.

Samples with an IgG concentration greater than 300 IU/ml should be reassayed after diluting 1/4 with the negative control.

If the dilution factor was not entered when the Work List was created (see User's Manual), multiply the result by the dilution factor to obtain the sample concentration.

Samples with concentrations greater than 300 IU/ml will be reported as "greater than 300 IU/ml". Samples with high anti-toxoplasma IgM titers may affect dilution test results.

For serological monitoring, it is recommended to test the samples from the different collections, in the same run, using the same lot.

Interpretation of test results should be made taking into consideration the patient's history, and the results of any other tests performed.

## QUALITY CONTROL

A positive and a negative control are included in each VIDAS TOXO IgG II kit.

These controls must be performed immediately after opening a new kit to ensure that reagent performance has not been altered. Each calibration must also be checked using these controls. The instrument will only be able to check the control values if they are identified by C1 and C2. Results cannot be validated if the control values deviate from the expected values.

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## LIMITATIONS OF THE METHOD

1. Sera collected very early in the acute stage of disease may have IgG levels < 4 IU/ml.
2. Positive test results from cord blood should be interpreted with caution. The presence of total or IgG anti-*T. gondii* antibodies in cord blood is usually the result of passive transfer from the mother to the fetus.
3. Positive test results may not be valid in persons who have received blood transfusions or other blood products within the past several months.

4. VIDAS TOXO IgG II results should be used in conjunction with clinical symptoms and laboratory test results, such as anti-*T. gondii* IgM results.

5. Test results from immunosuppressed patients may be difficult to interpret due to diminished immune response.

6. The VIDAS TOXO IgG II assay may only be used with serum or plasma (EDTA, heparin) samples. The use of amniotic fluid or other body fluids has not been established.

7. Interference may be encountered with certain sera containing antibodies directed against reagent components. For this reason, assay results should be interpreted taking into consideration the patient's history, and the results of any other tests performed.

## PERFORMANCE

Studies performed using VIDAS TOXO IgG II gave the following results:

### 1. Sensitivity - Specificity

1940 samples were evaluated on 3 sites in comparison with another immunoenzymatic reference technique.

49 samples were found to be equivocal with one and/or two other techniques, and were not included in the analysis.

Consolidated results for the 3 sites are as follows:

		EIA Technique		
		Positive	Negative	Total
VIDAS	Positive	572	3	575
	Negative	9	1307	1316
Total		581	1310	1891

The 12 discrepant samples were analyzed using a third reference technique: the Sabin-Feldmann test (Dye Test). The definitive interpretation for each of these samples is that obtained with at least two out of the three techniques used (2/3 rule).

Of the 3 VIDAS TOXO IgG II-positive/EIA-negative samples, 2 were positive according to the 2/3 rule.

Of the 9 VIDAS TOXO IgG II-negative/EIA-positive samples, 6 were negative according to the 2/3 rule and 1 could not be tested using the reference technique due to insufficient quantity.

Consolidated results after resolving discrepant samples are as follows:

		2/3 Rule		
		Positive	Negative	Total
VIDAS	Positive	574	1	575
	Negative	2	1313	1315
Total		576	1314	1890

Sensitivity: 99.65 %

(Range at 95 %: 98.75 – 99.96 %)

Specificity: 99.92 %

(Range at 95 %: 99.58 % - 100 %)

## Precision

Precision was tested with a negative sample, a low positive sample and a high positive sample.

Each sample was tested in duplicate two times per day for twenty days.

Within-run reproducibility (intra-assay precision) and between-run reproducibility (total precision) were calculated according to the recommendations of NCCLS Document EP5-T2, volume 12 number 4.

### Within-run reproducibility:

Sample	N	Mean titer	Standard deviation	CV %
Negative	80	1.21	0.25	--*
Low positive	80	17.18	0.88	5.13
High positive	80	168.90	12.17	7.21

### Between-run reproducibility:

Sample	N	Mean titer	Standard deviation	CV %
Negative	80	1.21	0.43	--*
Low positive	80	17.18	1.15	6.70
High positive	80	168.90	18.38	10.88

\* Not significant.

## CROSS REACTIVITY AND RELEVANT INTERFERENTS

No cross-reactivity or interference in the VIDAS TOXO IgG II assay was seen with any of the samples tested positive for rheumatoid factor, antinuclear antibodies, or Epstein Barr virus.

## RANGE OF EXPECTED VALUES

*Toxoplasma gondii* is a strict pathogen whose prevalence differs from one country to another or even one region to another. Contamination by *T. gondii* can vary according to cultural customs and eating habits, resulting in a prevalence ranging from less than 10% in certain regions of Northern Europe to more than 90% in Africa.

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
Release date	Part Number	Change Type	Change Summary
2015/01	09065I	Administrative	INDEX OF SYMBOLS REVISION HISTORY
		Technical	CONTENT OF THE KIT (60 TESTS) WARNINGS AND PRECAUTIONS
2015/06	09065J	Technical	CONTENT OF THE KIT (60 TESTS) INSTRUCTIONS FOR USE

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## VIDAS® HBs Ag Ultra (HBS)

VIDAS HBs Ag Ultra (HBS) is an automated qualitative test for use on the VIDAS family instruments for the detection of hepatitis B surface antigen (HBs Ag) in human serum or plasma, using the ELFA technique (Enzyme Linked Fluorescent Assay).

### SUMMARY AND EXPLANATION

The hepatitis B virus is responsible for acute and chronic hepatitis infections. Acute hepatitis can be asymptomatic or present symptoms of varying severity which may progress to fulminant hepatitis in 0.1 to 0.5% of cases. Chronicity occurs in 5 to 10% of cases in adults, but up to 90% of cases in infants following perinatal transmission. Currently, approximately 350 million people worldwide are chronic carriers of the virus (1). Chronic hepatitis may be asymptomatic and lead to liver lesions of varying severity, possibly evolving to cirrhosis, with an evolution in 5% of cases to hepatocellular carcinoma (2). The hepatitis B virus can be transmitted by parenteral or perinatal pathways or through sexual contact. Persons most at risk are health workers, drug addicts, those with multiple sexual partners, multiple transfusion or hemodialysis patients, as well as close friends and family of an infected subject, and newborns of an infected mother (2).

The discovery of the Australia antigen in 1965 – later known as HBs antigen - combined with viral hepatitis (3, 4) was a major breakthrough in the diagnosis of hepatitis B.

HBs antigen appears several days to several weeks after contact with the virus and can persist for several months. Persistence of HBs antigen for more than 6 months serologically defines chronic HBV infection. Disappearance of the HBs antigen is normally followed by the appearance of anti-HBs antibodies, which is a sign of recovery. The anti-HBs antibody assay is performed to confirm efficacy of the HBV vaccine. Anti-HBc antibodies are normally detected at the onset of the disease (2). Anti-HBc IgM titers are elevated (> 100 UPEI/ml) during acute hepatitis, then their titer decreases or disappears. However, during chronic hepatitis, the appearance of low IgM titers, reflecting hepatic cytolysis, confirms the active phase of the disease (5). Total anti-HBc antibodies, mainly IgG, are detected during acute and chronic hepatitis and persist after recovery.

HBe antigen is a circulating protein with a sequence very similar to that of HBc antigen, but with distinct antigenicity. During acute or chronic hepatitis, the presence of HBe antigen is generally associated with intense viral replication. Its disappearance coincides with the appearance of anti-HBe antibodies. Anti-HBe seroconversion is generally an indicator of a favorable prognosis of recovery (6). Nevertheless, the possible selection of a mutant strain incapable of synthesizing Hbe antigen may also lead to anti-HBe seroconversion, thereby limiting the prognostic value of HBe and anti-HBe markers (7). In this case, only the viral DNA assay is capable of directly detecting viral replication.

### PRINCIPLE

The VIDAS HBs Ag test is an enzyme-linked fluorescent immunoassay (ELFA) that is performed in the automated VIDAS system (see User's Manual). This assay can be performed according to 2 protocols: HBL long protocol (90 minutes), HBS short protocol (60 minutes).

The Solid Phase Receptacle (SPR®) serves as the solid phase as well as the pipetting device for the assay. Reagents for the assay are ready to use and pre-dispensed in the sealed reagent strips.

All of the assay steps are performed automatically by the instrument. The reaction medium is cycled in and out of the SPR several times.

After a preliminary washing step, the antigen present in the sample will bind simultaneously to the monoclonal antibody (8) coating the interior of the SPR and to the antibody conjugated with biotin. Unbound sample components are washed away. The antigen bound to the solid phase and to the biotinylated antibody is in contact with streptavidine conjugated with alkaline phosphatase, which will bind with biotin. Another wash step follows and removes unbound components.

During the final detection step, the substrate (4-Methyl-umbelliferyl phosphate) is cycled in and out of the SPR. The conjugate enzyme catalyzes the hydrolysis of this substrate into a fluorescent product (4-Methyl-umbelliferone), the fluorescence of which is measured at 450 nm. The intensity of the fluorescence is proportional to the concentration of antigen present in the sample.

At the end of the assay, results are analyzed automatically by the instrument and are expressed as an index calculated using a standard.

**CONTENT OF THE KIT (60 TESTS) - RECONSTITUTION OF REAGENTS:**

60 HBS strips	STR	Ready-to-use.
60 HBS SPRs 2 x 30	SPR <sup>®</sup>	Ready-to-use. Interior of SPRs coated with two monoclonal anti-HBs Ag antibodies (mouse) (5).
HBS standard 3 x 1 ml (lyophilized)	S1	Serum base* supplemented with inactivated human plasma HBs antigen + preservatives. The standard must be reconstituted with 1 ml of sterile distilled water (measured exactly). Allow to dissolve for at least 20 mins, then mix using a vortex. After reconstitution, the standard must be stored in aliquots at $-25 \pm 6^{\circ}\text{C}$ for up to 6 months, by freezing within an hour of reconstitution. Avoid successive freezing and thawing.
HBS positive control 1 x 1.5 ml (liquid)	C1	Ready-to-use. Serum base* supplemented with inactivated human plasma HBs antigen + 1 g/l sodium azide. MLE data indicate the index: confidence interval ("Control C1(+) Test Value Range)
Negative control 1 x 1.9 ml (liquid)	C2	Ready-to-use. Phosphate buffer + protein stabilizer of animal origin+ preservatives.
Specifications for the factory master data required to calibrate the test: • MLE data (Master Lot Entry) provided in the kit, or • MLE bar code printed on the box label.		
1 Package insert provided in the kit or downloadable from <a href="http://www.biomerieux.com/techlib">www.biomerieux.com/techlib</a>		

\* This product has been tested and shown to be negative for HBs antigen, antibodies to HIV1, HIV2 and HCV. However, since no existing test method can totally guarantee their absence, this product must be treated as potentially infectious. Therefore, usual safety procedures should be observed when handling.

**The SPR**

The interior of the SPR is coated during production with monoclonal anti-HBs Ag antibody (mouse). Each SPR is identified by the HBS code. Only remove the required number of SPRs from the pouch and **carefully reseal the pouch after opening**.

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Hazard statement

H318 : Causes serious eye damage.

H373 : May cause damage to organs through prolonged or repeated exposure.

H315 : Causes skin irritation.

H302 : Harmful if swallowed.

Precautionary statement

P280 :Wear protective gloves/protective clothing/eye protection/face protection.

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H318 : Causes serious eye damage.

Precautionary statement

P280 :Wear protective gloves/protective clothing/eye protection/face protection.

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- Pipette with disposable tip to dispense 1 ml and 150 µl.
- Powderless, disposable gloves.
- For other specific materials and disposables, please refer to the Instrument User's Manual.
- VIDAS family instrument.

**WARNINGS AND PRECAUTIONS**

- **For *in vitro* diagnostic use only.**
- **For professional use only.**
- **This kit contains products of human origin. No known analysis method can totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (see Laboratory biosafety manual - WHO - Geneva - latest edition).**
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (do not ingest or inhale).
- **Sample to sample contamination through contact with gloves: since high concentrations of HBs Ag may be encountered, it is strongly recommended to keep tips on a stand prior to use. It is also recommended to keep one tube of sample aside for testing HBs antigen.**
- Do not use the SPRs if the pouch is pierced.
- Do not use visibly deteriorated STRs (damaged foil or plastic).
- Do not use reagents after the expiration date indicated on the label.
- Do not mix reagents (or disposables) from different lots.
- Use **powderless** gloves, as powder has been reported to cause false results for certain enzyme immunoassay tests.

- Kit reagents contain sodium azide which can react with lead or copper plumbing to form explosive metal azides. If any liquid containing sodium azide is disposed of in the plumbing system, drains should be flushed with water to avoid build-up.
- The wash buffer (wells 3, 4, 7, 8, 9) contains a harmful agent (diethanolamine). Refer to the hazard statements "H" and the precautionary statements "P" above.
- The optical cuvette with substrate (well 10) contains an irritant agent (diethanolamine). Refer to the hazard statements "H" and the precautionary statements "P" above.
- Spills should be wiped up thoroughly after treatment with liquid detergent or a solution of household bleach containing at least 0.5% sodium hypochlorite. See the User's Manual for cleaning spills on or in the instrument. Do not autoclave solutions containing bleach.
- The instrument should be regularly cleaned and decontaminated (see the User's Manual).

**STORAGE CONDITIONS**

- Store the VIDAS HBs Ag Ultra kit at 2-8°C.
- **Do not freeze reagents, with the exception of S1 standard after reconstitution.**
- Store all unused reagents at 2-8°C, with the exception of S1 standard which must be kept frozen after reconstitution.
- After opening the kit, check that the SPR pouch is correctly sealed and undamaged. If not, do not use the SPRs.
- **Carefully reseal the pouch with the desiccant inside after use to maintain stability of the SPRs and return the complete kit to 2-8°C.**
- If stored according to the recommended conditions, all components are stable until the expiration date indicated on the label. Refer to the kit composition table for special storage conditions.

## SPECIMENS

### Specimen type and collection

Use sera (plain tube, tube with separator gel, plain tube with beads) or plasmas collected in lithium heparin.

Serum and plasma should be stored separated from the pellet. Samples containing impurities must be clarified by centrifugation.

The results obtained were not found to be influenced by: icteric samples (bilirubin concentrations up to 500 µmol/l), hemolyzed samples (hemoglobin concentrations up to 270 µmol/l of monomer) and lipemic samples (up to 30 mg/ml).

Do not inactivate specimens.

### Specimen stability

Samples can be stored for 5 days in stoppered tubes at 2-8°C. If longer storage is required, freeze the sera or plasma at -25 ± 6°C. A study performed on samples frozen for 2 months, showed that the quality of results is not affected.

## INSTRUCTIONS FOR USE

For complete instructions, see the User's Manual.

### Reading Master lot data

Before each new lot of reagents is used, enter the specifications (or factory master data) into the instrument using the master lot entry (MLE) data.

If this operation is not performed **before initiating the tests**, the instrument will not be able to print results.

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It is possible to enter MLE data **manually or automatically** depending on the instrument (refer to the User's Manual).

### Calibration

Calibration, using the standard provided in the kit, must be performed upon receipt of a new lot of reagents after the master lot data have been entered. Calibration should then be performed every 14 days. This operation provides instrument-specific calibration curves and compensates for possible minor variations in assay signal throughout the shelf-life of the kit.

The standard, identified by S1, must be tested **in duplicate** (see User's Manual). The standard value must be within the set RFV "Relative Fluorescence Value" range. If this is not the case, recalibrate.

### Procedure

1. **Only remove the required reagents from the refrigerator and allow them to come to room temperature for at least 30 minutes.**
2. Use one "HBS" strip and one "HBS" SPR for each sample, control or standard to be tested. **Make sure the storage pouch has been carefully resealed after the required SPRs have been removed.**
3. The test is identified by the "HBS" code on the instrument, or by the "HBL" code for the long protocol. The standard must be identified by "S1", and tested **in duplicate**. If the positive control is to be tested, it should be identified by "C1". If the negative control needs to be tested, it should be identified by C2.

4. Mix the standard, controls and samples using a vortex-type mixer (for serum or plasma separated from the pellet).

**5. For this test, the calibrator, control, and sample test portion is 150 µl.**

6. Insert the "HBS" SPRs and "HBS" strips into the instrument. Check to make sure the color labels with the assay code on the SPRs and the Reagent Strips match.
7. Initiate the assay as directed in the User's Manual. All the assay steps are performed automatically by the instrument.
8. Reclose the vials and return them to the required temperature after pipetting.
9. The assay will be completed within approximately 60 to 90 minutes, depending on the protocol selected. After the assay is completed, remove the SPRs and strips from the instrument.
10. Dispose of the used SPRs and strips into an appropriate recipient.

## RESULTS AND INTERPRETATION

Once the assay is completed, results are analyzed automatically by the computer. Fluorescence is measured twice in the Reagent Strip's reading cuvette for each sample tested. The first reading is a background reading of the substrate cuvette before the SPR is introduced into the substrate. The second reading is taken after incubating the substrate with the enzyme remaining on the interior of the SPR. The RFV (Relative Fluorescence Value) is calculated by subtracting the background reading from the final result. This calculation appears on the result sheet. The patient RFV is interpreted by the VIDAS system as follows:

$$i = \text{test value} = \text{patient RFV} / \text{standard RFV}$$

The test value and interpretation are also indicated on the result sheet. Interpretation of results according to test values is as follows:

Test value		
Short protocol	Long protocol	Interpretation
$i < 0.13$	$i < 0.10$	Negative
$i \geq 0.13$	$i \geq 0.10$	Positive

**If a positive result is obtained for a patient with no previous history, the assay must be repeated and confirmed using a neutralization test (VIDAS HBs Ag Ultra Confirmation ref. 30 317) or other tests.**

A positive sample should be confirmed using the same assay protocol (short or long) as was initially used to test the sample.

Before retesting, samples should be centrifuged again so as to eliminate any interference caused by fibrin fragments or cell elements.

Interpretation of test results should be made taking into consideration the patient history, and the results of any other tests performed.

VIDAS HBs Ag Ultra is calibrated against the panel of the Société Française de Transfusion Sanguine (French Blood Transfusion Society) (adw2/ayw3 mixture expressed in ng/ml).

The analytical sensitivity of VIDAS HBs Ag Ultra, which was determined using this panel, is less than 0.20 ng/ml with the short protocol (HBS), and less than 0.15 ng/ml with the long protocol (HBL).

### QUALITY CONTROL

One positive control and one negative control are included in each VIDAS HBs Ag Ultra kit.

These controls must be performed immediately after opening a new kit to ensure that reagent performance has not been altered. Each calibration must also be checked using these controls. The instrument will only be able to check the control values if they are identified by C1 and C2.

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### Note

It is the responsibility of the user to perform Quality Control in accordance with any local applicable regulations.

### LIMITATIONS OF THE METHOD

- As interference (i.e. anti-idiotypic antibody) may be encountered with certain sera, the test should only be declared positive after taking into account the patient's history and the results of the other hepatitis B markers.
- A negative Ag HBs result does not allow infection by the hepatitis B virus to be excluded. The Ag HBs serum concentration may be below the analytical sensitivity of the reagent. The presence of a modified HBs antigen (variant) cannot be excluded; the antigen may, in this case, have been incorrectly recognized or not recognized by the antibodies in the reagent. The results of this test must be interpreted taking into consideration the patient's history and the results of any other tests performed (confirmation by neutralization, HBV DNA etc.).
- In rare cases it is possible to detect both HBs antigen and anti-HBs antibodies.
- This assay has been validated for serum and plasma and should not be used for other biological fluids such as saliva, CSF or urine.
- This assay should not be used with specimens collected post-mortem.
- Do not use serum pools.

### RANGE OF EXPECTED VALUES (9)

The world wide prevalence of hepatitis B varies depending on the countries:

Region of the world	Prevalence of HBs Ag (%)
Europe, North America, Australia	0.2 - 0.5
Eastern Europe, Mediterranean basin, Russia, southwest Asia, South America	2 - 7
Southeast Asia, tropical Africa	8 - 20

### PERFORMANCE

The results of the studies, which demonstrated the conformity of VIDAS HBs Ag Ultra to the Common Technical Specifications of Directive 98/79/EC, are as follows:

### 1. Specificity for blood donor population

5059 blood donor samples from 2 blood transfusion centers were tested using the short protocol and the long protocol. No discrepancies were observed between the two protocols.

VIDAS HBs Ag Ultra (HBS/HBL)	Final interpretation of other EIA techniques	
	Positive	Negative
Positive	0	0
Negative	0	5059

Relative specificity with VIDAS HBs Ag Ultra for this population: 100.00%  
(95% confidence interval: 99.87% - 100.00%).

### 2. Clinical specificity

#### a) for hospitalized patients:

200 samples were tested using the short protocol, the long protocol and another EIA technique. No discrepancies were observed between the two protocols and the technique used as reference.

Relative specificity of VIDAS HBs Ag Ultra, with the two protocols (short and long), for this population: 100.00%.  
(95% confidence interval: 98.04% - 100.00%).

#### b) for outpatients in a screening center:

100 samples were tested using the short protocol, the long protocol and another EIA technique. No discrepancies were observed between the two protocols and the technique used as reference.

Relative specificity of VIDAS HBs Ag Ultra, with the two protocols (short and long), for this population: 100.00%.  
(95% confidence interval: 96.38% - 100.00%).

### 3. Analytical sensitivity

An external study performed using the panel of the French Blood Transfusion Society, showed a sensitivity of 0.12 ng/ml with the short protocol (HBS) and 0.08 ng/ml with the long protocol (HBL).

The sensitivity of VIDAS HBs Ag Ultra, determined using the international standard NIBSC 00/588, was estimated at 0.05 IU/ml with the long protocol (HBL) and 0.075 IU/ml with the short protocol (HBS).

### 4. Diagnostic sensitivity

30 fresh samples with a positive status (collection < 24 hours) were tested and found to be positive using VIDAS HBs Ag Ultra with the short protocol (HBS). 30 fresh samples with a negative status (collection < 24 hours) were tested and found to be negative using VIDAS HBs Ag Ultra with the short protocol (HBS).

31 fresh samples with a positive status (collection < 24 hours) were tested and found to be positive using VIDAS HBs Ag Ultra with the long protocol (HBL). 31 fresh samples with a negative status (collection < 24 hours) were tested and found to be negative using VIDAS HBs Ag Ultra with the long protocol (HBL).

## 5. Detectability

A study was performed using 506 samples characterized as positive, including 28 subtyped or genotyped samples and 46 samples from patients with acute hepatitis. Detectability with the long protocol (threshold of 0.08 ng/ml) and with the short protocol (threshold of 0.12 ng/ml), is 100.00 %.

(95% confidence interval: 99.22% - 100.00%).

## 6. Sensitivity of seroconversion panels

During a study, 32 seroconversion panels were tested using the two VIDAS HBs Ag Ultra protocols.

With the short protocol, HBs Ag was detected earlier by the VIDAS HBs Ag Ultra than with the comparison method in 13 out of 32 seroconversion panels (12 panels one specimen collection earlier, 1 panel two specimen collections earlier). HBs Ag was detected in one panel, one specimen collection later.

With the long protocol, 14 of the 32 panels tested were detected earlier than with the comparison method (9 panels one specimen collection earlier, 5 panels two specimen collections earlier).

## 7. Sensitivity of HBs Ag mutants

A panel of 27 recombinant proteins imitating the most significant mutations in HBs Ag amino acid sequences and 4 native samples harboring HBs Ag mutants were tested successfully (10).

## 8. Precision

Intra-assay and inter-assay reproducibility were determined at two sites and calculated according to the recommendations of the NCCLS EP5-T2 document, volume 12-4.

### Intra-assay reproducibility

The results are expressed as an index:

Sample	N	Short protocol		Long protocol	
		Mean	CV (%)	Mean	CV (%)
Positive	72	3.42	3.41	3.07	2.68
Weak positive	72	0.21	4.98	0.15	3.94
Negative	80	0.01	0.00 *	0.02	0.00 *

\* since the index value is too low to be used for the calculation of the CV, only the standard deviation is indicated.

### Inter-assay reproducibility

The total precision takes into account all the sources of variability.

The results are expressed as an index:

Sample	N	Short protocol		Long protocol	
		Mean	CV (%)	Mean	CV (%)
Positive	72	3.42	11.58	3.07	6.17
Weak positive	72	0.21	12.41	0.15	6.75
Negative	80	0.01	0.01 *	0.02	0.01 *

\* since the index value is too low to be used for the calculation of the CV, only the standard deviation is indicated.

## 9. Cross-reactivity

367 samples from patients whose physiological status is likely to affect the detection of HBs antigen, were tested using both the short and the long protocols. All the samples were found to be negative with another EIA technique. No discrepancies were observed between the two protocols and the technique used as reference.

	VIDAS HBs Ag Ultra positive samples
HCV +	0/20
EBV +	0/10
HIV +	0/10
CMV IgG +	0/10
HAV IgG +	0/10
HSV +	0/10
Syphilis	0/10
Rubella IgG +	0/10
Toxoplasmosis IgG +	0/10
Rheumatoid factor	0/9
Anti-nuclear antibodies	0/9
Non-viral hepatitis	0/5
Dialysis patients	0/10
Children less than 15 years old	0/10
Vaccinated: Anti-HBs +	0/10
Pregnant women*	0/214**

\* including 22 multipara

\*\* including 2 false positive results which were not repeated with the long protocol.

## WASTE DISPOSAL

Dispose of used or unused reagents as well as any other contaminated disposable materials following procedures for infectious or potentially infectious products.

It is the responsibility of each laboratory to handle waste and effluents produced according to their type and degree of hazardousness and to treat and dispose of them (or have them treated and disposed of) in accordance with any applicable regulations.


**LITERATURE REFERENCES**

1. LAVANCHY D., Worldwide epidemiology of HBV infection, disease burden, and vaccine prevention. Journal of Clinical Virology 34 suppl. 1 (2005) S1-S3.
2. HOLLINGER F.B., Hepatitis B virus, in Fields Virology, Third Edition, Lippincott-Raven Publishers, Philadelphia, 1996, 2739-2807.
3. BLUMBERG B.S., ALTER H.J., VISNICH S. JAMA, A "New" Antigen in Leukemia Sera, 1965, 191, 541-546.
4. PRINCE A.M., An antigen detected in blood during the incubation period of serum hepatitis, Proc Natl Acad Sci USA, 1968, 60, 814-821.
5. BRUNETTO, M.R., CERENZIA M.T., OLIVERI F., PIAONTI P., RANDONE A., CALVO P., MANZINI P., ROCCA G., GALLI C. and BONINO F., Monitoring the natural course and response to therapy of chronic hepatitis B with an automated semi-quantitative assay for IgM anti-HBc, Journal of Hepatology, 1993, 19, 431-436.
6. ALDERSHVILE J., FRÖSNER G.G., NEILSEN J.O. et al., Hepatitis B e antigen and antibody measured by radioimmunoassay in acute hepatitis B surface antigen-positive hepatitis, Journal of Infectious Diseases, 1980, 141, 293-298.
7. TONG S.T. and TREPO C., The HBe-minus mutants of hepatitis B virus in The Molecular Medicine of Viral Hepatitis, Harrison T.J. and Zuckerman A.J. ed., 1997, 89-104.
8. JOLIVET-REYNAUD. C., LESENECHAL. M., O' DONNELL. B. et al.- Localization of hepatitis B surface antigen epitopes present on variants and specifically recognised by anti-hepatitis B surface antigen monoclonal antibodies - Journal of Medical Virology - 2001, vol. 65, p.241-249.
9. ZUCKERMAN AJ., Hepatitis Viruses. In: Baron S, eds. Medical Microbiology, 4th ed. The University of Texas M Branch at Galveston, 1996: 849-863.
10. WEBER B., VAN DER TAELEM-BRULE N., BERGER N. et al. Evaluation of a new automated assay for hepatitis B surface antigen (HBs Ag) detection VIDAS HBs Ag Ultra, Journal of Virological Methods 135, 2006, 109-117.

**WARRANTY**

*bioMérieux disclaims all warranties, express or implied, including any implied warranties of MERCHANTABILITY AND FITNESS FOR A PARTICULAR USE. bioMérieux shall not be liable for any incidental or consequential damages. IN NO EVENT SHALL BIOMERIEUX'S LIABILITY TO CUSTOMER UNDER ANY CLAIM EXCEED A REFUND OF THE AMOUNT PAID TO BIOMERIEUX FOR THE PRODUCT OR SERVICE WHICH IS THE SUBJECT OF THE CLAIM.*

**INDEX OF SYMBOLS**

Symbol	Meaning
	Catalog number
	In Vitro Diagnostic Medical Device
	Manufacturer
	Temperature limit
	Use by date
	Batch code
	Consult Instructions for Use
	Contains sufficient for <n> tests
	Date of manufacture

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Release date	Part Number	Change Type	Change Summary
2015/01	11728I	Administrative	INDEX OF SYMBOLS REVISION HISTORY
		Technical	CONTENT OF THE KIT (60 TESTS) - RECONSTITUTION OF REAGENTS WARNINGS AND PRECAUTIONS INSTRUCTIONS FOR USE
2016/05	11728J	Technical	CONTENT OF THE KIT (60 TESTS) – RECONSTITUTION OF REAGENTS

BIOMERIEUX, the BIOMERIEUX logo, SPR and VIDAS are used, pending and/or registered trademarks belonging to bioMérieux, or one of its subsidiaries, or one of its companies.

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
1. LAVANCHY D., Worldwide epidemiology of HBV infection, disease burden, and vaccine prevention. *Journal of Clinical Virology* 34 suppl. 1 (2005) S1-S3.
2. HOLLINGER F.B., Hepatitis B virus, in *Fields Virology*, Third Edition, Lippincott-Raven Publishers, Philadelphia, 1996, 2739-2807.
3. BLUMBERG B.S., ALTER H.J., VISNICH S. JAMA, A "New" Antigen in Leukemia Sera, 1965, 191, 541-546.
4. PRINCE A.M., An antigen detected in blood during the incubation period of serum hepatitis, *Proc Natl Acad Sci USA*, 1968, 60, 814-821.
5. BRUNETTO, M.R., CERENZIA M.T., OLIVERI F., PIANTONI P., RANDONE A., CALVO P., MANZINI P., ROCCA G., GALLI C. and BONINO F., Monitoring the natural course and response to therapy of chronic hepatitis B with an automated semi-quantitative assay for IgM anti-HBc, *Journal of Hepatology*, 1993, 19, 431-436.
6. ALDERSHIVILE J., FRÖSNER G.G., NEILSEN J.O. et al., Hepatitis B e antigen and antibody measured by radioimmunoassay in acute hepatitis B surface antigen-positive hepatitis, *Journal of Infectious Diseases*, 1980, 141, 293-298.
7. TONG S.T. and TREPO C., The HBe-minus mutants of hepatitis B virus in *The Molecular Medicine of Viral Hepatitis*, Harrison T.J. and Zuckerman A.J. ed., 1997, 89-104.
8. JOLIVET-REYNAUD. C., LESENECHAL. M., O' DONNELL. B. et al.- Localization of hepatitis B surface antigen epitopes present on variants and specifically recognised by anti-hepatitis B surface antigen monoclonal antibodies - *Journal of Medical Virology* - 2001, vol. 65, p.241-249.
9. ZUCKERMAN A.J., Hepatitis Viruses. In: Baron S, eds. *Medical Microbiology*, 4th ed. The University of Texas M Branch at Galveston, 1996 : 849-863.
10. WEBER B., VAN DER TAELEM-BRULE N., BERGER N. et al. Evaluation of a new automated assay for hepatitis B surface antigen (HBs Ag) detection VIDAS HBs Ag Ultra, *Journal of Virological Methods* 135, 2006, 109-117.

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## VIDAS<sup>®</sup> HCG (HCG)

VIDAS<sup>®</sup> HCG is an automated quantitative test for use on the VIDAS<sup>®</sup> family instruments, for the quantitative measurement of human Chorionic Gonadotropin in human serum or plasma (lithium heparin or EDTA) using the ELFA technique (Enzyme Linked Fluorescent Assay).

### SUMMARY AND EXPLANATION

Human Chorionic Gonadotropin (hCG) is a glycoprotein with a molecular weight of 40,000 daltons. Two polypeptide subunits, alpha and beta, form hCG. Approximately 30% of the weight of hCG is comprised of glycans (hexoses, sialic acid). Structurally, hCG is similar to the three pituitary hormones LH, FSH and TSH. The alpha subunit (92 amino acids), is practically identical for these four hormones. The beta subunit, composed of 145 amino acids and 5 oligosaccharide groups, contains a unique amino acid sequence which gives biological and immunological specificity to hCG. (1, 2). However, only the alpha-beta dimer is biologically active (3, 4, 5, 6).

hCG is mainly biosynthesized by the placental syncytiotrophoblast. The level of serum hCG will steadily increase and peak during the first three months of pregnancy. The level then decreases until the 16<sup>th</sup> week of pregnancy when it stabilizes until pregnancy term. The physiologic function of hCG in pregnancy is to maintain the gestational yellow body and to foster progesterone and estrogen production during the first three months of pregnancy. hCG also plays a role in the differentiation of the fetal genital tract. Assaying for hCG is an essential element in the diagnosis and follow-up of pregnancy (ectopic pregnancy, hydatidiform mole, etc.) (7, 8).

hCG can also be produced by certain tumor tissues, mainly those of trophoblastic origin (choriocarcinoma, etc.). In this context, hCG represents a biological marker of considerable clinical significance, both for prognosis and therapeutic monitoring (9).

### PRINCIPLE

The assay principle combines an enzyme immunoassay sandwich method with a final fluorescent detection (ELFA). The Solid Phase Receptacle (SPR<sup>®</sup>), serves as the solid phase as well as the pipetting device for the assay. Reagents for the assay are ready-to-use and pre-dispensed in the sealed reagent strips.

All of the assay steps are performed automatically by the instrument.

The reaction medium is cycled in and out of the SPR several times.

Firstly, the sample is taken and transferred into the well containing alkaline-phosphatase labeled anti-hCG antibody. The sample/conjugate mixture is cycled in and out of the SPR several times to increase the reaction speed. The antigen binds to antibodies coated on the SPR and to the conjugate forming a "sandwich".

Then, the remaining free hCG sites are saturated by cycling the conjugate in the fifth well of the strip in and out of the SPR. Unbound components are eliminated during the washing steps.

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At the end of the assay, results are automatically calculated by the instrument in relation to the calibration curve stored in memory, and then printed out.

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60 HCG strips	STR	Ready-to-use.
60 HCG SPRs 2 x 30	SPR	Ready-to-use. Interior of SPRs coated with monoclonal anti-HCG immunoglobulins (mouse).
HCG control 1 x 2 mL (lyophilized)	C1	Reconstitute with 2 mL of distilled water. Wait for 5 to 10 minutes. Mix. Stable after reconstitution for 14 days at 2-8°C or until the expiration date on the kit at -25 ± 6°C. 5 freeze/thaw cycles are possible. Human serum* + human hCG + preservatives. MLE data indicate the confidence interval in "mIU/mL" (milli-international units per milliliter) ("Control C1 Dose Value Range").
HCG calibrator 2 x 2 mL (lyophilized)	S1	Reconstitute with 2 mL of distilled water. Wait for 5 to 10 minutes. Mix. Stable after reconstitution for 14 days at 2-8°C or until the expiration date on the kit at -25 ± 6°C. 5 freeze/thaw cycles are possible. Calf serum + human hCG + preservatives. MLE data indicate the concentration in mIU/mL ("Calibrator (S1) Dose Value") and the confidence interval in "Relative Fluorescence Value" ("Calibrator (S1) RFV Range").
HCG diluent 2 x 25 mL (liquid)	R1	Ready-to-use. Calf serum + 1 g/L sodium azide.
Specifications for the factory master data required to calibrate the test:		
<ul style="list-style-type: none"> <li>• MLE data (Master Lot Entry) provided in the kit,</li> </ul> or <ul style="list-style-type: none"> <li>• MLE bar code printed on the box label.</li> </ul>		
1 Package Insert provided in the kit or downloadable from <a href="http://www.biomerieux.com/techlib">www.biomerieux.com/techlib</a> .		

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**The SPR**

The interior of the SPR<sup>®</sup> is coated during production with monoclonal anti-hCG immunoglobulins (mouse). Each SPR is identified by the HCG code. Only remove the required number of SPRs from the pouch and **carefully reseal the pouch after opening**.

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**Description of the HCG strip**

Wells	Reagents
1	Sample well.
2 - 3 - 4	Empty wells.
5	Conjugate: Alkaline phosphatase-labeled monoclonal anti-hCG immunoglobulins (mouse) + 1 g/L sodium azide (600 µL).
6 - 7	Wash buffer: Sodium phosphate (0.01 mol/L) pH 7.4 + 1 g/L sodium azide (600 µL).
8	Wash buffer: diethanolamine DEA* (1.1 mol/L or 11.5%) pH 9.8 + 1 g/L sodium azide (600 µL).
9	Empty well.
10	Cuvette with substrate: 4 Methyl-umbelliferyl-phosphate (0.6 mmol/L) + DEA** (0.62 mol/L or 6.6%, pH 9.2) + 1 g/L sodium azide (300 µL).

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H318: Causes serious eye damage.

H373: May cause damage to organs through prolonged or repeated exposure.

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P309 + P311: IF exposed or if you feel unwell: Call a POISON CENTER or doctor/physician.

\*\* Signal Word: **DANGER**Hazard statement

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- Pipette with disposable tip to dispense 2 mL and 100 µL.
- Powderless, disposable gloves.
- For other specific materials and disposables, please refer to the Instrument User Manual.
- VIDAS<sup>®</sup> family instrument.

**WARNINGS AND PRECAUTIONS**

- **For *in vitro* diagnostic use only.**
- **For professional use only.**
- **This kit contains products of human origin. No known analysis method can totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (see Laboratory biosafety manual - WHO - Geneva - latest edition).**

- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (do not ingest or inhale).
- Do not use the SPRs if the pouch is pierced.
- Do not use visibly deteriorated STRs (damaged foil or plastic).
- Do not use reagents after the expiration date indicated on the label.
- Do not mix reagents (or disposables) from different lots.
- Use **powderless** gloves, as powder has been reported to cause false results for certain enzyme immunoassay tests.

- Kit reagents contain sodium azide which can react with lead or copper plumbing to form explosive metal azides. If any liquid containing sodium azide is disposed of in the plumbing system, drains should be flushed with water to avoid build-up.
- The wash buffer in well 8 contains a harmful agent (11.5% diethanolamine). Refer to the hazard statements "H" and the precautionary statements "P" above.
- The substrate in well 10 contains an irritant agent (6.6% diethanolamine). Refer to the hazard statements "H" and the precautionary statements "P" above.
- Due to the extremely high concentrations of hCG found in certain samples, it is recommended that the utmost care be taken when pipetting or handling samples, in order to avoid contamination.
- Spills should be wiped up thoroughly after treatment with liquid detergent and a solution of household bleach containing at least 0.5% sodium hypochlorite. See the User Manual for cleaning spills on or in the instrument. Do not autoclave solutions containing bleach.
- The instrument should be regularly cleaned and decontaminated (see the User Manual).

### STORAGE CONDITIONS

- Store the VIDAS<sup>®</sup> HCG kit at 2-8°C.
- **Do not freeze the SPRs, strips and diluent.**
- **Store all unused reagents at 2-8°C.**
- After opening the kit, check that the SPR pouch is correctly sealed and undamaged. If not, do not use the SPRs.
- **Carefully reseal the pouch with the desiccant inside after use to maintain stability of the SPRs and return the complete kit to 2-8°C.**
- If stored according to the recommended conditions, all components are stable until the expiration date indicated on the label. Refer to the kit composition table for special storage conditions.

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- hemolysis (after spiking samples with hemoglobin (monomer): 0 to 300 µmol/L),
- lipemia (after spiking samples with lipids: 0 to 2 g/L equivalent in triglycerides),
- bilirubinemia (after spiking samples with bilirubin: 0 to 513 µmol/L).

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#### Specimen stability:

Samples can be stored at 2-8°C in stoppered tubes for up to 7 days; if longer storage is required, freeze the sera or plasma at -25 ± 6°C. Avoid successive freezing and thawing.

### INSTRUCTIONS FOR USE

**For complete instructions, see the User Manual.**

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#### Calibration

Calibration, using the calibrator provided in the kit, must be performed each time a new lot of reagents is opened, after the master lot data have been entered. Calibration should then be performed every 14 days. This operation provides instrument-specific calibration curves and compensates for possible minor variations in assay signal throughout the shelf-life of the kit.

The calibrator, identified by S1, must be tested in **duplicate** (see User Manual). The calibrator value must be within the set RFV "Relative Fluorescence Value" range. If this is not the case, recalibrate.

#### Procedure

1. **Only remove the required reagents from the refrigerator and allow them to come to room temperature for at least 30 minutes.**
2. Use one "HCG" strip and one "HCG" SPR for each sample, control or calibrator to be tested. **Make sure the storage pouch has been carefully resealed after the required SPRs have been removed.**
3. The test is identified by the "HCG" code on the instrument. The calibrator must be identified by "S1", and tested in **duplicate**. If the control is to be tested, it should be identified by "C1".
4. Mix the controls and samples using a vortex-type mixer (for serum or plasma separated from the pellet).
5. **For this test, the calibrator, control, and sample test portion is 100 µL.**
6. Insert the "HCG" SPRs and "HCG" strips into the instrument. Check to make sure the color labels with the assay code on the SPRs and the Reagent Strips match.
7. Initiate the assay as directed in the User Manual. All the assay steps are performed automatically by the instrument.
8. Reclose the vials and return them to the required temperature after pipetting.
9. The assay will be completed within approximately 30 minutes. After the assay is completed, remove the SPRs and strips from the instrument.
10. Dispose of the used SPRs and strips into an appropriate recipient.

## RESULTS AND INTERPRETATION

Once the assay is completed, results are analyzed automatically by the computer. Fluorescence is measured twice in the Reagent Strip reading cuvette for each sample tested. The first reading is a background reading of the substrate cuvette before the SPR is introduced into the substrate. The second reading is taken after incubating the substrate with the enzyme remaining on the interior of the SPR. The RFV (Relative Fluorescence Value) is calculated by subtracting the background reading from the final result. This calculation appears on the result sheet.

The results are automatically calculated by the instrument using calibration curves which are stored by the instrument (4-parameter logistic model) and the concentrations are expressed in mIU/mL (1<sup>st</sup> IRP 75/537).

Samples with hCG concentrations greater than 1,500 mIU/mL may be reassayed after diluting by 1/20 to 1/200 with the HCG diluent (R1). Any dilution that is not between 1/20 and 1/200 has not been validated and is the customer's responsibility.

If the dilution factor has not been entered when the Work List was created (see User Manual), multiply the result by the dilution factor to obtain the sample concentration.

## QUALITY CONTROL

A control is included in each VIDAS® HCG kit.

This control must be performed immediately after opening a new kit to ensure that reagent performance has not been altered. Each calibration must also be checked using this control. The instrument will only be able to check the control value if it is identified by C1.

Results cannot be validated if the control value deviates from the expected values.

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## LIMITATIONS OF THE METHOD

- This kit has not been validated within the context of Down syndrome screening.
- Low concentrations of hCG may indicate pregnancy, however, the hCG values obtained should be used for diagnosis in association with additional information obtained by the physician (questions, symptoms, clinical observations, other examinations, etc.).
- Since hCG levels in serum increase rapidly during the initial stages of pregnancy, the result can be confirmed by a second assay on a sample collected 48 hours later.
- Certain malignant pathologies, notably disorders of trophoblastic origin and germinal tumors produce similar levels of hCG to those observed in pregnancy. It is therefore recommended to refer to medical information. Likewise, recent injections of hCG used in therapeutic procedures can generate high concentrations.
- Positive results with VIDAS® HCG have been reported in patients treated with retinoic acid derivatives, in particular isotretinoin. In this case, results should be confirmed using another sample, or another technique.
- Interference may be encountered with certain sera containing antibodies directed against reagent components. For this reason, assay results should be interpreted taking into consideration the patient's history, and the results of any other tests performed.

## RANGE OF EXPECTED VALUES

In a clinically healthy population with no tumoral pathologies, the following values were found:

- men: < 3 mIU/mL
- cyclic women: < 4 mIU/mL
- menopausal women: < 13 mIU/mL
- pregnant women:

Week of amenorrhoea	Mean (mIU/mL)	Limits (mIU/mL)
4 - 5	7 400	1 500 - 23 000
5 - 6	32 800	3 400 - 135 300
6 - 7	52 000	10 500 - 161 000
7 - 8	74 000	18 000 - 209 000
8 - 9	100 000	37 500 - 219 000
9 - 10	105 000	42 800 - 218 000
10 - 11	96 000	33 700 - 218 700
11 - 12	75 300	21 800 - 193 200
12 - 13	66 700	20 300 - 166 100
13 - 14	65 900	15 400 - 190 000
2 <sup>nd</sup> trimester (14-26)	26 150	2 800 - 176 100
3 <sup>rd</sup> trimester (26-39)	27 200	2 800 - 144 400

These figures are given as a guide; it is recommended that each laboratory establishes its own reference values from a rigorously selected population.

## PERFORMANCE

Studies performed using VIDAS® HCG gave the following results:

### Measurement range

The measurement range of the VIDAS® HCG kit is: 2-1500 mIU/mL.

### Detection limit

Defined as the smallest concentration of hCG which is significantly different from the zero concentration with a probability of 95%:  $\leq 2$  mIU/mL.

### Hook effect

The hook effect was tested using hCG concentrations up to 1 000 000 mIU/mL. No hook effect was found at the tested concentrations.

### Precision

#### Within-run reproducibility:

Five samples were tested 30 times in the same run.

Sample	1	2	3	4	5
Mean concentration (mIU/mL)	16.6	71.8	257.0	334.0	1211.0
CV %	6.8	4.3	4.6	4.8	5.4

#### Between-run reproducibility

Five samples were tested singly in 24 different runs over a 9-week period.

Sample	1	2	3	4	5
Mean concentration (mIU/mL)	17.4	71.7	253.0	334.0	1152.0
CV %	8.7	5.1	4.3	4.8	5.0

### Specificity

The percentage of cross-reactivity is determined by calculating the ratio of the tested compound concentration over the hCG concentration, for a signal of 1,000 RFV. The results obtained are the following:

Tested compound	Cross-reactivity (%)
hCG	100
hCG $\beta$ - subunit	2.3
hCG $\alpha$ - subunit	0.3
LH	0.1
FSH	0.1
TSH	0.07

Another study performed by addition of 10,000 mIU/mL of LH or 12,500 mIU/mL of FSH to a sample containing 1,356 mIU/mL of hCG did not show any interference.

### Accuracy

#### Dilution test

Three samples were diluted in the HCG diluent (R1) and tested singly in three runs. The ratio of the mean concentration measured over the expected concentration is expressed as a mean recovery percentage.

Sample	Dilution factor	Expected concentration (mIU/mL)	Mean measured concentration (mIU/mL)	Mean recovery percentage
1	1/1	168.19	168.19	100.0
	1/2	84.09	77.44	92.1
	1/4	42.05	36.61	87.1
	1/8	21.02	19.00	90.4
	1/16	10.51	8.30	79.0
	1/32	5.26	4.25	80.8
2	1/1	564.56	564.56	100.0
	1/2	282.28	327.42	116.0
	1/4	141.14	130.12	92.2
	1/8	70.57	64.62	91.6
	1/16	35.29	29.91	84.8
	1/32	17.64	14.33	81.2
3	1/1	1251.04	1251.04	100.0
	1/2	625.52	656.25	104.9
	1/4	312.76	322.58	103.1
	1/8	156.38	165.99	106.1
	1/16	78.19	74.96	95.9
	1/32	39.09	34.51	88.3

### Comparison with other test methods

73 serum samples were assayed in parallel using VIDAS<sup>®</sup> HCG and another commercially available kit. The allometric curve equation obtained is the following:  
 $Y = 0.95 X^{-18}$  with a correlation coefficient of 0.99.

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## INDEX OF SYMBOLS

Symbol	Meaning
	Catalogue number
	In Vitro Diagnostic Medical Device
	Manufacturer
	Temperature limit
	Use by date
	Batch code
	Consult Instructions for Use
	Contains sufficient for <n> tests
	Date of manufacture

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Release date	Part Number	Change Type	Change Summary
2015/01	06038P	Administrative	INDEX OF SYMBOLS REVISION HISTORY
		Technical	CONTENT OF THE KIT (60 TESTS) – RECONSTITUTION OF REAGENTS WARNINGS AND PRECAUTIONS
2015/06	06038Q	Technical	CONTENT OF THE KIT (60 TESTS) – RECONSTITUTION OF REAGENTS INSTRUCTIONS FOR USE
2018-10	06038R	Administrative	LIMITED WARRANTY
		Technical	RESULTS AND INTERPRETATION

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## VIDAS® Ferritin (FER)

**IVD**

VIDAS Ferritin is an automated quantitative test for use on the VIDAS family instruments for the determination of human Ferritin in human serum or plasma (lithium heparin or EDTA) using the ELFA technique (Enzyme Linked Fluorescent Assay).

### SUMMARY AND EXPLANATION

Ferritin is the most common type of iron storage in the human body (1). Its molecules can be found in the cytoplasm of the reticuloendothelial system, specifically in the liver and spleen.

Ferritin has the shape of a hollow spherical protein envelope composed of 24 protein sub-units.

The iron is situated in the center of the molecule in the form of ferric hydroxyphosphate. This molecule may contain as many as 4,500 iron atoms (2).

Iron deficiency anemia is a common ailment which may be caused by insufficient iron intake, pregnancy, hemodialysis or blood donation (3). The drop in serum ferritin levels can indicate an iron deficiency prior to the appearance of anemia. Detection of insufficient ferritin levels therefore enables anticipated treatment (4). Furthermore, iron overload is characteristic of illnesses such as thalassemia and sideroblastic anemia. In such cases, the measurement of serum ferritin levels contributes to diagnosis and patient monitoring (5, 6).

The level of serum ferritin acts as an indicator of the quantities of iron in the human body (7,8). It is also closely correlated to the level of iron in the bone marrow.

Quantitative data are obtained from the determination of serum ferritin concentration, thus avoiding the need for bone marrow biopsy, which is a more invasive technique (9, 10, 11).

### PRINCIPLE

The assay principle combines a one-step enzyme immunoassay sandwich method with a final fluorescent detection (ELFA).

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During the final detection step, the substrate (4-Methyl-umbelliferyl phosphate) is cycled in and out of the SPR. The conjugate enzyme catalyzes the hydrolysis of this substrate into a fluorescent product (4-Methyl-umbelliferone) the fluorescence of which is measured at 450 nm. The intensity of the fluorescence is proportional to the concentration of antigens present in the sample.

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### CONTENT OF THE KIT (60 TESTS):

60 FER strips	STR	Ready-to-use.
60 FER SPRs 2 x 30	SPR	Ready-to-use. Interior of SPRs coated with monoclonal anti-ferritin immunoglobulins (mouse).
FER control 1 x 2 ml (liquid)	C1	Ready-to-use. Tris buffer (0.1 mol/l) pH 7.4 + human ferritin + bovine albumin + protein and chemical stabilizers. MLE data indicate the confidence interval in ng/mL ("Control C1 Dose Value Range").
Calibrator 1 x 2 ml (liquid)	S1	Ready-to-use. Tris buffer (0.1 mol/l) pH 7.4 + human ferritin + bovine albumin + protein and chemical stabilizers. MLE data indicate the concentration in ng/mL ("Calibrator (S1) Dose Value") and the confidence interval in "Relative Fluorescence Value ("Calibrator (S1) RFV Range").
FER dilution buffer 1 x 25 ml (liquid)	R1	Ready-to-use. Tris buffer (0.1 mol/l) pH 7.4 + bovine albumin + protein and chemical stabilizers.
Specifications for the factory master data required to calibrate the test:		
<ul style="list-style-type: none"> <li>MLE data (Master Lot Entry) provided in the kit,</li> </ul> or <ul style="list-style-type: none"> <li>MLE bar code printed on the box label.</li> </ul>		
1 Package insert provided in the kit or downloadable from <a href="http://www.biomerieux.com/techlib">www.biomerieux.com/techlib</a> .		

**The SPR**

The interior of the SPR is coated during production with mouse monoclonal anti-ferritin immunoglobulins. Each SPR is identified by the FER code. Only remove the required number of SPRs from the pouch and **carefully reseal the pouch after opening**.

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**Description of the FER reagent strip**

Wells	Reagents
1	Sample well.
2 - 3 - 4	Empty wells.
5	Conjugate: alkaline phosphatase-labeled monoclonal anti-ferritin immunoglobulins (mouse) + 1 g/l sodium azide (600 µl).
6 - 7	Wash buffer: sodium phosphate (0.01 mol/l) pH 7.4 + 1 g/l sodium azide (600 µl).
8	Wash buffer: diethanolamine* (1.1 mol/l or 11.5 %, pH 9.8) + 1 g/l sodium azide (600 µl).
9	Empty well.
10	Cuvette with substrate: 4-Methyl-umbelliferyl-phosphate (0.6 mmol/l) + diethanolamine (DEA**) (0.62 mol/l or 6.6%, pH 9.2) + 1 g/l sodium azide (300 µl).

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- Pipette with disposable tip to dispense 100 µl.
- Powderless, disposable gloves.
- For other specific materials and disposables, please refer to the Instrument User's Manual
- VIDAS family instrument.

**WARNINGS AND PRECAUTIONS**

- For *in vitro* diagnostic use only.
- For professional use only.
- This kit contains products of human origin. No known analysis method can totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (see Laboratory biosafety manual - WHO - Geneva - latest edition).

- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (do not ingest or inhale).
- Do not use the SPRs if the pouch is pierced.
- Do not use visibly deteriorated STRs (damaged foil or plastic).
- Do not use reagents after the expiration date indicated on the label.
- Do not mix reagents (or disposables) from different lots.
- Use **powderless** gloves, as powder has been reported to cause false results for certain enzyme immunoassay tests.
- Kit reagents contain sodium azide which can react with lead or copper plumbing to form explosive metal azides. If any liquid containing sodium azide is disposed of in the plumbing system, drains should be flushed with water to avoid build-up.
- The wash buffer (well 8) contains a harmful agent (11.5% diethanolamine). Refer to the hazard statements "H" and the precautionary statements "P" above.
- The substrate (well 10) contains an irritant agent (6.6% diethanolamine). Refer to the hazard statements "H" and the precautionary statements "P" above.
- Spills should be wiped up thoroughly after treatment with liquid detergent and a solution of household bleach containing at least 0.5% sodium hypochlorite. See the User's Manual for cleaning spills on or in the instrument. Do not autoclave solutions containing bleach.
- The instrument should be regularly cleaned and decontaminated (see the User's Manual).

### STORAGE CONDITIONS

- Store the FER VIDAS kit at 2-8°C.
- **Do not freeze reagents.**
- **Store all unused reagents at 2-8°C.**
- After opening the kit, check that the SPR pouch is correctly sealed and undamaged. If not, do not use the SPRs.
- **Carefully reseal the pouch with the desiccant inside after use to maintain stability of the SPRs and return the complete kit to 2-8°C.**
- If stored according to the recommended conditions, all components are stable until the expiration date indicated on the label.

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The use of heat inactivated sera has not been validated.

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- hemolysis (after spiking samples with hemoglobin: 0 to 300 µmol/l (monomer)),
- lipemia (after spiking samples with lipids: 0 to 2 g/l equivalent in triglycerides),
- bilirubinemia (after spiking samples with bilirubin: 0 to 513 µmol/l).

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#### Specimen stability

Samples can be stored at 2-8°C in stoppered tubes for up to 7 days; if longer storage is required, freeze the sera or plasma at -25 ± 6° C.

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The calibrator, identified by S1, must be tested **in duplicate** (see User's Manual). The calibrator value must be within the set RFV "Relative Fluorescence Value" range. If this is not the case, recalibrate.

#### Procedure

1. **Only remove the required reagents from the refrigerator and allow them to come to room temperature for at least 30 minutes.**
2. Use one "FER" strip and one "FER" SPR for each sample, control or calibrator to be tested. **Make sure the storage pouch has been carefully resealed after the required SPRs have been removed.**
3. The test is identified by the "FER" code on the instrument. The calibrator must be identified by "S1", and tested **in duplicate**. If the control needs to be tested, it should be identified by "C1".
4. Mix the calibrator, control and samples using a vortex-type mixer (for serum or plasma separated from the pellet).
5. **For this test, the calibrator, control, and sample test portion is 100 µl.**
6. Insert the "FER" SPRs and "FER" strips into the instrument. Check to make sure the color labels with the assay code on the SPRs and the Reagent Strips match.
7. Initiate the assay as directed in the User's Manual. All the assay steps are performed automatically by the instrument.
8. Restopper the vials and return them to 2-8°C after pipetting.
9. The assay will be completed within approximately 30 minutes. After the assay is completed, remove the SPRs and strips from the instrument.
10. Dispose of used SPRs and reagent strips in an appropriate recipient.

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The results are automatically calculated by the instrument using calibration curves which are stored by the instrument (4-parameter logistics model) the concentrations are expressed in ng/ml (preparation NIBSC 80/578).

Samples with ferritin concentrations greater than 1,200 ng/ml must be retested after dilution to 1/10 or 1/100 in FER dilution buffer (R1)

If the dilution factor has not been entered when the Work List was created (see User's Manual), multiply the result by the dilution factor to obtain the sample concentration.

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## RANGE OF EXPECTED VALUES

These figures are given as a guide; it is recommended that each laboratory establish its own reference values from a rigorously selected population.

"Expected values were determined using 206 samples from clinically healthy, hematologically normal subjects with no hepatic disorders".

When the population observed is described using the percentage method, the following results are obtained:

### Men:

Range of values	0 - 68 ng/ml	68 - 208 ng/ml	208 - 434 ng/ml	Mean
Prevalence	5%	45%	45%	236 ng/ml

### Normal menstruating women:

Range of values	0 - 9.3 ng/ml	9.3 - 45 ng/ml	45 - 159 ng/ml	Mean
Prevalence	5%	45%	45%	58 ng/ml

### Menopausal women:

Range of values	0 - 24.4 ng/ml	24.4 - 118 ng/ml	118 - 278 ng/ml	Mean
Prevalence	5%	45%	45%	151 ng/ml

If concentrations are lower than 20 ng/ml in women and 30 ng/ml in men, tests for martial deficiency should be performed. If concentrations are greater than 250 ng/ml in women or 350 ng/ml in men, tests should be performed for inflammatory, infectious, hepatic or tumoral pathologies, or abnormal iron storage (idiopathic or secondary hemochromatosis). "It is recommended that each laboratory establishes its own reference values from a rigorously selected population."

## PERFORMANCE

Studies performed using VIDAS Ferritin gave the following results:

### Measurement range

The measurement range of the VIDAS Ferritin reagent extends from 1.5 to 1200 ng/ml (2<sup>nd</sup>.IS NIBSC 80/578).

### Analytical detection limit

Defined as the smallest concentration of ferritin which is significantly different from the zero concentration with a probability of 95%:  $\leq 1.5$  ng/ml.

## QUALITY CONTROL

A control is included in each VIDAS FER kit.

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### Hook effect

No hook effect was found up to ferritin concentrations of 100,000 ng/ml.

### Precision

#### Within-run reproducibility:

5 samples were tested 30 times in a same run.

Sample	1	2	3	4	5
Dose (ng/ml)	15.3	102	239	466	924
CV%	6.2	4.6	5.0	4.0	4.0

**Between-run reproducibility**

5 samples were tested singly in 24 different runs on the same VIDAS.

Sample	1	2	3	4	5
Dose (ng/ml)	16.5	128	234	537	1121
CV%	4.4	5.9	7.0	4.9	4.6

**Specificity**

Tested compound	Cross-reactivity (%)
Spleen ferritin	106
Liver ferritin	121
Heart ferritin	28
Placental ferritin	137

**Accuracy****Dilution test**

Three samples were diluted in the (R1) diluent and tested singly in 3 series. The ratio of the mean concentration measured over the expected concentration is expressed as a mean recovery percentage.

Sample no.	Dilution factor	Measured concentration (ng/ml)	Expected concentration (ng/ml)	Mean recovery percentage (%)
SC3	1/1	193.7	--	--
	1/2	96.4	96.9	99.5
	1/4	49.8	48.4	102.8
	1/8	25.9	24.2	107.0
	1/16	13.3	12.1	110.0
	1/32	6.5	6.1	108.0
SC4	1/1	422.5	--	--
	1/2	202.6	211.2	96.0
	1/4	106.9	105.6	101.0
	1/8	52.3	52.8	99.0
	1/16	27.2	26.4	102.0
	1/32	13.6	13.2	103.0
SC5	1/1	926.7	--	--
	1/2	473.8	463.4	102.5
	1/4	251.6	231.7	108.0
	1/8	110.2	115.8	95.1
	1/16	55.6	57.9	95.4
	1/32	30.0	29	103.6

**Comparison with other test methods**

Correlation was established between the VIDAS Ferritin kit and another commercially available reagent.

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Release date	Part Number	Change Type	Change Summary
2015/01	06036K	Administrative	INDEX OF SYMBOLS REVISION HISTORY
		Technical	CONTENT OF THE KIT (60 TESTS) WARNINGS AND PRECAUTIONS
2015/06	06036L	Technical	CONTENT OF THE KIT (60 TESTS) INSTRUCTIONS FOR USE

BIOMERIEUX, the BIOMERIEUX logo, SPR and VIDAS are used, pending and/or registered trademarks belonging to bioMérieux, or one of its subsidiaries, or one of its companies.  
Any other name or trademark is the property of its respective owner.

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**VIDAS<sup>®</sup> HIV DUO Quick (HIV6)**

IVD

VIDAS HIV DUO Quick is an automated HIV infection screening test for use on the VIDAS family instruments, for the combined detection of anti-HIV-1 (groups M and O) and anti-HIV-2 total immunoglobulins and HIV-1 p24 antigen in human serum or plasma (lithium heparin or EDTA) using the ELFA technique (Enzyme Linked Fluorescent Assay).

**SUMMARY AND EXPLANATION**

Human immunodeficiency viruses (HIV) are RNA retroviruses transmitted via sexual contact, parenteral and perinatal pathways or the placenta. HIV-1 and HIV-2 were respectively isolated in 1983 and 1985 in patients infected by AIDS (Acquired Immunodeficiency Syndrome). Since then numerous genetic variants have been characterized. These mutations seemed to be without consequence for serological diagnosis until HIV-1 variants of group O (Outlier) were isolated, since they have only 50% homology at the *env* gene level with those of group M (Major). (1, 2, 3, 4, 5).

During 2007, 2.5 million individuals were newly infected by HIV. According to estimates, 33.2 million [30.6 – 36.1 million] people were living with HIV in 2007 (6).

Current diagnosis of HIV infection relies on the detection of anti-HIV serum antibodies using an ELISA method. However, there is a mean period of 3 weeks between contamination and the appearance of the first antibodies (7). During this period p24 antigen is present in most people infected by HIV-1 (8). Thus simultaneous detection of p24 antigenemia, anti-HIV-1 and anti-HIV-2 antibodies enables the time lapse between contamination and diagnosis of the infection to be decreased (9, 10). VIDAS HIV DUO Quick is an automated test based on the combined detection of HIV-1 p24 antigen and anti-HIV-1 and anti-HIV-2 total immunoglobulins, enabling a reduction of the seroconversion window.

Note: The use of VIDAS HIV DUO Quick does not rule out the obligation of using a second technique in countries where legislation requires the use of two different screening tests.

**PRINCIPLE**

The principle of the test combines 2 enzyme immunoassay reactions with a final fluorescent detection (ELFA).

The Solid Phase Receptacle (SPR<sup>®</sup>), serves as the solid phase as well as the pipetting device for the assay. Reagents for the assay are ready-to-use and are pre-dispensed in the sealed reagent strips. All of the assay steps are performed automatically by the instrument. The reaction medium is cycled in and out of the SPR several times.

The upper part of the SPR is coated with monoclonal anti-p24 antibodies for the detection of p24 antigen. The lower part of the SPR enables the detection of anti-HIV-1 and anti-HIV-2 antibodies: it is coated with an HIV-1 gp160 protein, and HIV-1 group O and HIV-2 specific synthetic peptides.

During preliminary incubation, the sample, the biotinylated anti-p24 antibody (rabbit) in the strip and the biotinylated antigen (the same used in the solid phase) are cycled in and out of the SPR. Lysis of the virus occurs during incubation and the released p24 antigen binds to the monoclonal anti-p24 antibody coated on the SPR and is recognized by the biotinylated anti-p24 antibody. At the same time, the anti-HIV-1 and/or anti-HIV-2 antibodies bind to the gp160 and/or the peptides in the lower part of the SPR and are recognized by the biotinylated antigen. Two washing steps remove unbound components.

A second incubation with the biotinylated antigen in the strip (same used in the solid phase) is only performed in the lower part of the SPR. The biotinylated antigen binds to any anti-HIV antibody coated in the lower part of the SPR.

Any excess reagent is eliminated by washing.

A third incubation is performed with alkaline phosphatase-labeled streptavidin. During this step, the streptavidin binds to the biotinylated anti-p24 antibody, if it is present in the upper part of the SPR, and to the biotinylated antigen, if it is present in the lower part of the SPR.

Any excess reagent is eliminated by washing.

During the final detection step, the substrate (4-Methyl-umbelliferyl phosphate) is cycled in and out of the SPR. The conjugate enzyme catalyzes this substrate into a fluorescent product (4-Methyl-umbelliferone), the fluorescence of which is measured at 450 nm. The intensity of the fluorescence is proportional to the presence of anti-HIV antibody and/or p24 antigen in the sample.

At the end of the assay, results are automatically calculated by the instrument in relation to the standard, and then printed out.

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60 HIV6 strips	STR	Ready-to-use.
60 HIV6 SPRs (2 x 30)	SPR	Ready-to-use. Interior of SPRs coated with gp160 (HIV-1) protein, HIV-1 group O and HIV-2 specific synthetic peptides, and monoclonal anti-p24 antibodies.
HIV6 positive antibody control (1 x 2.0 ml) (lyophilized)	C1	Reconstitute with 2 ml of distilled water. Leave for 5 to 10 minutes and then mix. After reconstitution, stable for 2 months at 2-8°C.  Human serum* with anti-HIV-1 antibody + TRIS (0.1 mol/l, pH 7.4) + protein and chemical stabilizers + preservatives.  MLE data indicate the confidence interval for the antibody result (Control C1 (+) Test Value Range).
HIV6 negative control (1 x 1.3 ml) (liquid)	C2	Ready-to-use. Human serum** + TRIS (0.1 mol/l, pH 7.4) + protein and chemical stabilizers + preservatives.
HIV6 positive antigen control (1 x 2.0 ml) (lyophilized)	C3	Reconstitute with 2 ml of distilled water. Leave for 5 to 10 minutes and then mix. After reconstitution, stable for 2 months at 2-8°C.  Human serum** with inactivated HIV-1 viral lysate + TRIS (0.1 mol/l, pH 7.4) + protein and chemical stabilizers + preservatives.  MLE data indicate the confidence interval for the antigen result (Control C3 (+) Test Value Range).-
HIV6 standard (1 x 2.0 ml) (lyophilized)	S1	Reconstitute with 2 ml of distilled water. Leave for 5 to 10 minutes and then mix. After reconstitution, stable for 2 months at 2-8°C.  Human serum* with anti-HIV-1 antibody + TRIS (0.1 mol/l, pH 7.4) + protein and chemical stabilizers + preservatives.  MLE data indicate the confidence interval in "Relative Fluorescence Value (RFV)" ("Standard (S1) RFV Range").
Specifications for the factory master data required to calibrate the test: • MLE data (Master Lot Entry) provided in the kit, or • MLE bar code printed on the box label.		
1 Clip seal		
1 Package insert provided in the kit or downloadable from <a href="http://www.biomerieux.com/techlib">www.biomerieux.com/techlib</a>		

\* This product has been tested and shown to be negative for HBs surface antigen and antibodies to HCV. This product has been heat inactivated (30 minutes at 56°C). However, since no existing test method can totally guarantee their absence, this product must be treated as potentially infectious. Therefore, usual safety procedures should be observed when handling.

\*\* This product has been tested and shown to be negative for HBs surface antigen, and antibodies to HIV1, HIV2 and HCV. However, since no existing test method can totally guarantee their absence, this product must be treated as potentially infectious. Therefore, usual safety procedures should be observed when handling.

**The SPR**

The interior of the SPR is coated during production with a gp160 protein and HIV-1 group O and HIV-2 specific synthetic peptides in the lower part, and monoclonal anti-p24 antibodies in the upper part. The antigens enable the capture of anti-HIV-1 and/or anti-HIV-2 antibodies, and the monoclonal antibodies enable the capture of HIV-1 p24 antigen.

Each SPR is identified by the HIV-6 code. Only remove the required number of SPRs from the pouch and **carefully reseal the pouch after opening using the clip seal provided with the kit.**

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**Description of the VIDAS HIV DUO Quick strip:**

Wells	Reagents
1	Sample well.
2	HEPES buffer + biotin-labeled anti-p24 antibody (rabbit) + triton X100 + goat serum + protein and chemical stabilizers + 0.2 g/l gentamicin sulfate + 0.9 g/l sodium azide (300 µl).
3 - 5 - 6 - 8 - 9	Wash buffer: TRIS + protein and chemical stabilizers + 0.9 g/l sodium azide (600 µl).
4	HEPES buffer + gp160 protein and biotin-labeled synthetic peptides + skimmed milk + protein and chemical stabilizers + 0.2 g/l gentamicin sulfate + 0.9 g/l sodium azide (300 µl).
7	Tracer: alkaline phosphatase-labeled streptavidin + TRIS + skimmed milk + 2.5 g/l bovine albumin + triton X100 + protein and chemical stabilizers + 0.2 g/l gentamicin sulfate + 0.9 g/l sodium azide (400 µl).
10	Reading cuvette with substrate: 4-Methyl-umbelliferyl phosphate (0.6 mmol/l) + diethanolamine (DEA*) (0.62 mol/l or 6.6%) pH 9.2 + 1 g/l sodium azide (300 µl).

\* Signal Word: **DANGER**Hazard statement

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- Pipette with disposable tip to dispense 2 ml and 200 µl.
- Powderless, disposable gloves.
- For other specific materials and disposables, please refer to the Instrument User's Manual.
- VIDAS family instrument.

**WARNINGS AND PRECAUTIONS**

- **For *in vitro* diagnostic use only.**
- **For professional use only.**
- **This kit contains products of human origin. No known analysis method can totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (see Laboratory biosafety manual - WHO - Geneva - Latest edition).**
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (do not ingest or inhale).
- Do not use the SPRs if the pouch is pierced.
- Do not use visibly deteriorated STRs (damaged foil or plastic).
- Do not use reagents after the expiration date indicated on the box label.
- Do not mix reagents (or disposables) from different lots.

- Use **powderless** gloves, as powder has been reported to cause false results for certain enzyme immunoassay tests.
- Kit reagents contain sodium azide which can react with lead or copper plumbing to form explosive metal azides. If any liquid containing sodium azide is disposed of in the plumbing system, drains should be flushed with water to avoid build-up.
- The substrate in well 10 contains an irritant agent (6.6% diethanolamine). Refer to the hazard statements "H" and the precautionary statements "P" above.
- Spills should be wiped up thoroughly after treatment with liquid detergent or a solution of household bleach containing at least 0.5% sodium hypochlorite. See the User's Manual for cleaning spills on or in the instrument. Do not autoclave solutions containing bleach.
- The instrument should be regularly cleaned and decontaminated (see the User's Manual).

**STORAGE CONDITIONS**

- Store the VIDAS HIV DUO Quick kit at 2-8°C.
- **Do not freeze reagents.**
- **Store all unused reagents at 2-8°C.**
- After opening the kit, check that the SPR pouch is correctly sealed and undamaged. If not, do not use the SPRs.
- **To maintain stability of the remaining SPRs, carefully reseal the pouch after use with the dessicant inside using the clip seal provided, and return the complete kit to 2-8°C.**

- If stored according to the recommended conditions, all components are stable until the expiration date indicated on the label. Refer to the kit composition table for special storage conditions.

## SPECIMENS

### Specimen type and collection:

Human serum or plasma (collected in lithium heparin or EDTA).

### **Do not use recalcified plasma.**

It is recommended that each laboratory checks the compatibility of collection tubes used.

Do not inactivate specimens.

None of the following factors have been found to significantly influence this assay.

- hemolysis (after spiking samples with hemoglobin, up to 300 µmol/l (monomer)).
- triglyceride concentration (after spiking samples with up to 30 mg/ml of triglycerides)
- bilirubin concentration (after spiking samples with up to 200 mg/l of bilirubin).

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### Specimen stability

Samples can be stored at 2-8 °C in stoppered tubes for up to 2 days. If longer storage is required, freeze the sera or plasma at -25 ± 6°C.

Avoid successive freezing and thawing.

A study performed on samples frozen for 2 months, showed that the quality of results is not affected.

## INSTRUCTIONS FOR USE

**For complete instructions, see the User's Manual.**

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### Calibration

Calibration, using the standard provided in the kit, must be performed each time a new lot of reagents is opened, after the master lot data have been entered. Calibration should then be performed every 14 days. This operation provides instrument-specific calibration curves and compensates for possible minor variations in assay signal throughout the shelf-life of the kit.

The standard, identified by S1, must be tested in **duplicate** (see User's Manual). The standard value must be within the set RFV "Relative Fluorescence Value" range. If this is not the case, recalibrate.

## Procedure

1. **Only remove the required reagents from the refrigerator and allow them to come to room temperature for at least 30 minutes.**
2. Use one "HIV6" strip and one "HIV6" SPR for each sample, control or standard to be tested. **Make sure the storage pouch has been carefully resealed after the required SPRs have been removed.**
3. The test is identified by the "HIV6" code on the instrument. The standard must be identified by "S1", and tested in **duplicate**. If the positive controls are to be tested, they should be identified by "C1" and "C3". If the negative control needs to be tested, it should be identified by "C2".
4. If necessary, clarify samples by centrifugation.
5. Mix the standard, controls and samples using a vortex-type mixer (for serum or plasma separated from the pellet) in order to improve result reproducibility.
6. **For this test, the calibrator, control, and sample test portion is precisely 200 µl. The test portion is pipetted into the sample well, which is used as the reaction cuvette for this test.**
7. Insert the "HIV6" SPRs and "HIV6" strips into the instrument. Check to make sure the color labels with the assay code on the SPRs and the Reagent Strips match.
8. Initiate the assay as directed in the User's Manual. All the assay steps are performed automatically by the instrument.
9. Restopper the vials and return them to 2–8°C after pipetting.
10. The assay will be completed within approximately 80 minutes. After the assay is completed, remove the SPRs and strips from the instrument.
11. Dispose of the used SPRs and strips into an appropriate recipient.

## RESULTS AND INTERPRETATION

Once the assay is completed, results are analyzed automatically by the computer. Fluorescence is measured twice in each Reagent Strip's reading cuvette for each sample tested. The first reading is a background reading of the substrate cuvette before the SPR is introduced into the substrate. The second reading is taken after the substrate in the SPR has been incubated. The RFV (Relative Fluorescence Value) is calculated by subtracting the background reading from the final result. This calculation appears on the result sheet.

The test value is calculated by the instrument for each sample as follows:

$$\text{Test value} = \text{patient RFV} / \text{standard RFV}$$

This test value and the interpreted result are also included on the result sheet.

The test value is interpreted as follows:

Test value	Interpreted result
< 0.25	Negative
≥ 0.25	Positive

If the message "Invalid" is generated for the sample, the result is invalid. In this case, the sample should be retested. Samples with a negative test value are considered to be negative, within the performance limitations of the reagent. In cases of suspected primary infection, values  $\leq 0.25$  must be interpreted with caution.

Samples with a positive test value should be retested in duplicate:

- Non repeatable positive samples (two negative reactions for three tests) are considered to be negative, within the performance limitations of the reagent.
- Repeatable positive samples (at least two positive reactions for three tests) should be confirmed using additional techniques.

The additional tests can be analyzed by Western Blotting and/or using a second screening test for detection of anti-HIV antibodies, p24 antigen and/or determination of viral load.

It is strongly recommended to test a second sample collected a few days later, particularly in the presence of clinical symptoms and/or risk factors.

**Note:** country-specific requirements for HIV diagnostics must be taken into account if necessary.

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## QUALITY CONTROL

A positive antibody control, a positive antigen control, and a negative control are included in each VIDAS HIV DUO Quick kit. These controls must be performed immediately after opening a new kit to ensure that reagent performance has not been altered. Each calibration must also be checked using these controls. The instrument will only be able to check the control values if they are identified by C1, C2 and C3.

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## LIMITATIONS OF THE METHOD

VIDAS HIV DUO Quick is an HIV infection screening test.

**It should not be used as a specific test for the detection of HIV-1 p24 antigenemia. VIDAS HIV DUO Quick must not be used as a supplementary test for VIDAS HIV DUO Ultra.**

This test has been validated for use with serum or plasma but not with other body fluids such as saliva, CSF or urine.

Since the use of recalcified plasma may generate false positive results, this type of sample must not be used.

Do not use a mixture of samples.

Interference may be encountered with certain sera containing antibodies directed against reagent components. For this reason, assay results should be interpreted taking into consideration the patient history, and the results of any other tests performed.

## PERFORMANCE

Studies performed using VIDAS HIV DUO Quick gave the following results:

### 1. Specificity on a blood donor population:

5027 blood donor samples from 3 blood banks were tested.

VIDAS HIV DUO Quick	Final interpretation	
	Positive	Negative
Positive	0	6
Negative	0	5021

Specificity of the VIDAS HIV DUO Quick reagent on this population: 99.88%

(95% confidence interval: 99.74% - 99.95%).

### 2. Specificity on hospitalized patients:

Two false positive results were found for the 200 samples tested.

Specificity of the VIDAS HIV DUO Quick reagent on this population: 99.00%.

(95% confidence interval: 96.35% - 99.73%)

### 3. Specificity on high-risk behavior patients:

No false positive results were found for the 100 samples tested.

Specificity of the VIDAS HIV DUO Quick reagent on this population: 100.00%.

(95% confidence interval: 96.38% - 100.00%)

### 4. Diagnostic sensitivity:

A study was performed using 440 samples presumed to be HIV-1 positive and 122 HIV-2 positive. All of the samples were found to be positive.

Sensitivity of the VIDAS HIV DUO Quick reagent on this population: 100.00% (95% confidence interval: 99.29% - 100.00%).

In order to check the sensitivity of VIDAS HIV DUO Quick on non-B subtypes, 40 group M samples (4A, 14B, 3C, 3D, 3F, 3G, 3H, 3CRF01-AE and 4 CRF02-AG) and 10 group O samples were tested. All were found to be positive.

31 fresh samples with a positive status (collection < 24 hours) were tested and found to be positive with VIDAS HIV DUO Quick. 31 fresh samples with a negative status (collection < 24 hours) were tested and found to be negative with VIDAS HIV DUO Quick.

In order to check the sensitivity of VIDAS HIV DUO Quick on p24 antigen, at least 50 culture supernatants including different HIV-1 and HIV-2 subtypes and at least 50 p24 antigen-positive samples were tested. All were found to be positive.

## 5. Sensitivity on seroconversion panels

During different studies, 39 seroconversion panels were tested, demonstrating the early detection capabilities of VIDAS HIV DUO Quick.

The dilutions of the 3 HIV-1 samples and the 2 HIV-2 samples of the EFS 96 panel were found to be positive.

The sensitivity of VIDAS HIV DUO Quick, determined on the Ag HIV SFTS panel (subtype B), was estimated at 16.5 pg/ml of HIV-1 antigen. The sensitivity of VIDAS HIV DUO Quick, determined using the international standard NIBSC 90/636 was estimated at 0.5 IU/ml of HIV-1 antigen.

At least 40 early HIV seroconversion samples were tested. The results are conform with the state of the art.

## 6. Precision

Intra-assay and inter-assay reproducibility, determined at two different sites, were calculated according to the recommendations of NCCLS document EP5-T2, volume 12-4.

### Site 1

Sample	n	Mean index	Intra-assay reproducibility %	Inter-assay reproducibility %
Low positive antigen	40	2.48	2.17	3.69
Low positive antibody	40	3.96	3.53	4.14

### Site 2

Sample	n	Mean index	Intra-assay reproducibility %	Inter-assay reproducibility %
Low positive antigen	40	2.23	2.07	9.88
Low positive antibody	40	3.63	3.46	7.11

### Global analysis

Sample	n	Mean index	Inter-assay reproducibility %
Low positive antigen	80	2.36	10.23
Low positive antibody	80	3.79	8.20

In this analysis, the reproducibility takes into account within-run variability, and between-run and between-site dispersion.

## 7. Cross-reactivity

134 HIV negative samples from patients whose physiological status may interfere with the VIDAS HIV DUO Quick assay, were tested. The following results were obtained:

	Positive VIDAS HIV DUO Quick results
anti-HAV + antibody	0/10
anti-HCV + antibody	0/10
anti-HBV + antibody	0/10
anti-Gag + antibody	0/13
anti-EBV + antibody	0/10
anti-HSV + antibody	0/10
anti-CMV + antibody	0/10
anti-nuclear + antibody	0/10
anti- <i>Toxoplasma gondii</i> + antibody	0/6
anti- <i>Treponema pallidum</i> + antibody	0/5
Vaccinated against influenza	0/10
Vaccinated against hepatitis B	0/10
Rheumatoid factor +	0/10
Children under 15 years of age	0/10

Spiking of HIV positive and HIV negative samples with increasing concentrations of biotin (0 to 2 mg/l) did not affect the VIDAS HIV DUO Quick result.

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## LITERATURE REFERENCES

1. F.Barré-Sinoussi, J.C. Chermann, F.Rey et al. Isolation of T-Lymphotropic Retrovirus from a Patient at Risk for acquired Immune Deficiency Syndrome (AIDS). 20 May 1983. Science.Vol.220. page 868-871.
2. F.Clavel, D.Guétard, F.Brun-Vézinet et al. Isolation of a New human Retrovirus from West African Patients with AIDS. 18 July 1986. SCIENCE Vol 233 page 343-346.
3. R.de Leys, B.Vanderborght, M Vanden Haesevelde and al. Isolation and Partial Characterization of an Unusual Human Immunodeficiency Retrovirus from Two Persons of West-Central African Origin. Journal of Virology. Mar.1990.p1207-1216.
4. M.Guyader, M.Emmerman, P.Sonigo et al. Genome organization and transactivation of the human immunodeficiency virus type 2. Nature. 16 april 1987. Vol 326 p.662-669.
5. M.Popovic, M.G.Sarngadharan, E. Read et al. Detection, Isolation, and Continuous Production of Cytopathic Retroviruses (HTLV-III) from Patients with AIDS and Pre-AIDS. Science. 4 May 1984. Vol.224. p 497-500.
6. World Health Organization, AIDS epidemic update. December 2007. p 1-60.

7. S.Kleinman, M.P.Busch, J.J.Korelitz et al. The Incidence/Window Period Model and its Use to Assess the Risk of Transfusion-Transmitted Human Immunodeficiency Virus and Hepatitis C Virus Infection. *Transfusion Medecine Reviews*. July 1997. Vol 11, N°3 p.155-172.
8. M.Von Sydow, H.Gaines, A.Sönnerborg et al. Antigen detection in primary HIV infection. *British Medical Journal* vol.296. 23 January 1988. p.238-240.
9. A.Galetto-Lacour, S.Kinloch-de Loës, b.Hirschel et al. Primo-infection VIH: un diagnostic souvent évoqué mais posé tardivement. *Schwerz.Med Wochenschr* 1995: 125 n°8. p.341-346.
10. T.Meier, E.Knoll, M.Henkes et al. Evidence for a diagnostic window in fourth generation assays for HIV. *Journal of Clinical Virology* vol.23, 2001, p.113-116.

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Not applicable (First publication)

Correction of documentation anomalies

Addition, revision and/or removal of information related to the product

Implementation of non-technical changes noticeable to the user

### Note:

*Minor typographical, grammar, and formatting changes are not included in the revision history.*

Release date	Part Number	Change Type	Change Summary
2015/01	12232E	Administrative	INDEX OF SYMBOLS REVISION HISTORY
		Technical	CONTENT OF THE KIT (60 TESTS) – RECONSTITUTION OF REAGENTS WARNINGS AND PRECAUTIONS
2015/06	12232F	Technical	CONTENT OF THE KIT (60 TESTS) – RECONSTITUTION OF REAGENTS INSTRUCTIONS FOR USE

BIOMERIEUX, the blue logo, SPR and VIDAS are used, pending, and/or registered trademarks belonging to bioMérieux or one of its subsidiaries or one of its companies.

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
**LITERATŪRA**

1. F.Barré-Sinoussi, J.C. Chermann, F.Rey et al. Isolation of T-Lymphotropic Retrovirus from a Patient at Risk for acquired Immune Deficiency Syndrome (AIDS). 20 May 1983. Science.Vol.220. page 868-871.
2. F.Clavel, D.Guébard, F.Brun-Vézinet et al. Isolation of a New human Retrovirus from West African Patients with AIDS. 18 July 1986. SCIENCE Vol 233 page 343-346.
3. R.de Leys, B.Vanderborcht, M Vanden Haesevelde and al. Isolation and Partial Characterization of an Unusual Human Immunodeficiency Retrovirus from Two Persons of West-Central African Origin. Journal of Virology. Mar.1990.p1207-1216.
4. M.Guyader, M.Emerman, P.Sonigo et al. Genome organization and transactivation of the human immunodeficiency virus type 2. Nature. 16 april 1987. Vol 326 p.662-669.
5. M.Popovic, M.G.Sarngadharan, E. Read et al. Detection, Isolation, and Continuous Production of Cytopathic Retroviruses (HTLV-III) from Patients with AIDS and Pre-AIDS. Science. 4 May 1984. Vol.224. p 497-500.
6. World Health Organization, AIDS epidemic update. December 2007. p 1-60.
7. S.Kleinman, M.P.Busch, J.J.Korelitz et al. The Incidence/Window Period Model and its Use to Assess the Risk of Transfusion-Transmitted Human Immunodeficiency Virus and Hepatitis C Virus Infection. Transfusion Medecine Reviews. July 1997. Vol 11, N°3 p.155-172.

8. M.Von Sydow, H.Gaines, A.Sönnernborg et al. Antigen detection in primary HIV infection. British Medical Journal vol.296. 23 January 1988. p.238-240.
9. A.Galetto-Lacour, S.Kinloch-de Loës, b.Hirschel et al. Primo-infection VIH : un diagnostic souvent évoqué mais posé tardivement. Schwerz.Med Wochenschr 1995 : 125 n°8. p.341-346.
10. T.Meier, E.Knoll, M.Henkes et al. Evidence for a diagnostic window in fourth generation assays for HIV. Journal of Clinical Virology vol.23, 2001, p.113-116.

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
## VIDAS® D-Dimer Exclusion II™ (DEX2)

IVD

VIDAS® D-Dimer Exclusion II™ is an automated quantitative test for use on the instruments of the VIDAS® family for the immunoenzymatic determination of fibrin degradation products (FbDP) in human plasma (sodium citrate) using the ELFA technique (Enzyme Linked Fluorescent Assay).

VIDAS® D-Dimer Exclusion II™ is indicated for use in conjunction with a clinical pretest probability assessment model **to exclude deep vein thrombosis (DVT) and pulmonary embolism (PE)** disease in outpatients suspected of DVT or PE.

VIDAS® D-Dimer Exclusion II™ is indicated for use in the HERDOO2 clinical decision rule (CDR) to assess the risk of recurrence of venous thromboembolism (VTE) in women with a first unprovoked VTE. Risk stratification by this CDR is an aid to guide the duration of oral anticoagulant therapy.

### SUMMARY AND EXPLANATION

Fibrin degradation products (FbDP) are a highly heterogeneous group of soluble fragments that appear in the circulation as a result of two simultaneous physiological processes (1):

- Coagulation, resulting in the conversion of soluble fibrinogen into insoluble stabilized fibrin by the enzymes thrombin and factor XIIIa,
- Fibrinolysis, resulting in the dissolution of the fibrin clot by the enzyme plasmin. The D-dimer fragment is the terminal product of this process (2).

Although FbDP's vary in size, they are characterized by the presence of one or more D-dimer epitopes (1, 2). Therefore, FbDP's are collectively referred to as 'D-dimer'. Assays for D-dimer can be performed directly in plasma by virtue of monoclonal antibodies which are able to identify fibrin-specific epitopes without cross-reactivity with fibrinogen or its degradation products (1).

D-dimer reflects the presence of stabilized fibrin and this has made this marker a useful tool in the diagnosis of venous thromboembolism (VTE) (3). Quantitative D-dimer assays based on ELISA techniques have a high sensitivity for the presence of an occluding thrombus and, consequently, are particularly useful in excluding venous thromboembolism (4, 5).

In conjunction with assessment of clinical pretest probability, it is possible to safely rule out the diagnosis of deep vein thrombosis (DVT) and/or pulmonary embolism (PE) in suspected outpatients when the concentration of D-dimer is below a predefined cutoff (determined by rigorous clinical studies) (5-7). The clinical utility of D-dimer ELISA assays in the diagnostic work-up of suspected DVT or PE resides in the significant reduction in the number of imaging tests that are required, with a concomitant reduction in the total cost of diagnosis (8, 9). For patients with suspected pulmonary embolism, the use of an age-adjusted cutoff compared to a single clinical cutoff, enables to reduce the number of imaging tests required (10).

D-dimer is not specific for DVT/PE and elevated levels are also observed in a variety of other conditions where activation of coagulation and fibrinolysis occurs (for example, surgery, trauma, infection, inflammation, pregnancy, cancer) (3). This makes D-dimer less useful for exclusion of DVT or PE in hospitalized patients due to the high proportion of comorbid conditions associated with elevated D-dimer levels (5, 11). Under certain conditions, lower than expected D-dimer results may occur giving rise to false-negatives. Therefore, it is not safe to use D-dimer for exclusion of DVT/PE in patients with high pre-test probability, long duration of DVT/PE symptoms (more than one week) or already under anticoagulant treatment (11, 12).

The VIDAS® D-Dimer Exclusion II™ assay is part of the HERDOO2 clinical decision rule that has been shown to identify women with unprovoked VTE who are at sufficiently low risk of a recurrent VTE event that they could discontinue oral anticoagulants after short-term therapy (13).

The HERDOO2 rule states that this is the case if none or only one of the following four criteria is present: 1) **H**yperpigmentation, **E**dema or **R**edness (HER) on examination of legs; 2) **D**-dimer  $\geq 250$  ng/mL; 3) **O**besity (body mass index  $\geq 30$  kg/m<sup>2</sup>); 4) **O**lder age ( $\geq 65$  years).

The International Society on Thrombosis and Haemostasis (ISTH) suggests that it is safe to discontinue anticoagulants if the risk of recurrent VTE is less than 5% one year after stopping anticoagulant therapy (14).

The REVERSE II study confirmed that the one-year VTE recurrence rate was below the level recommended by the ISTH, since the recurrence rate in women with a first unprovoked VTE and 0 or 1 HERDOO2 criteria was 3.0% in this study (15).

### PRINCIPLE

The assay principle combines a two-step enzyme immunoassay sandwich method with a final fluorescent detection (ELFA).

The Solid Phase Receptacle (SPR®), serves as the solid phase with an anti-FbDP monoclonal antibody adsorbed on its surface, and the pipetting device. The reagents in the single-use reagent strip are ready-to-use and pre-dispensed.

All of the reaction steps are performed automatically by the instrument. The reaction medium is cycled in and out of the SPR® several times.

First the sample is taken by the SPR®, diluted and then cycled in and out of the SPR® several times. The antigen binds to the anti-FbDP immunoglobulins coated on the SPR®. Unbound components are eliminated during a washing step.

During the second step, the conjugate, which contains an alkaline phosphatase-labeled anti-FbDP monoclonal antibody, binds to the antigen coated on the SPR® to form a "sandwich".

Unbound components are eliminated during the washing steps. During the final detection step, the substrate (4-Methyl-umbelliferyl phosphate) is cycled in and out of the SPR®. The conjugate enzyme catalyzes the hydrolysis of this substrate into a fluorescent product (4-Methyl-umbelliferone), the fluorescence of which is measured at 450 nm. The intensity of the fluorescence is proportional to the concentration of antigen present in the sample.

At the end of the assay, results are automatically calculated by the instrument in relation to the calibration curve stored in memory. The results can then be printed.

**CONTENT OF THE KIT (60 TESTS) – RECONSTITUTION OF REAGENTS:**

60 DEX2 Strips	STR	Ready-to-use.
60 DEX2 SPRs 2 x 30	SPR®	Ready-to-use. Interior of SPRs coated with anti-FbDP monoclonal immunoglobulins (mouse).
<u>DEX2 Calibrator:</u> S1 Calibrator 2 x 2 mL (lyophilized)	S1	<b>Reconstitute with 2 mL of distilled water.</b> Wait for 5 minutes and then mix. <b>After reconstitution, the calibrator is stable for 28 days at 2-8°C or until the expiration date of the kit at -25 ± 6°C</b> (freeze immediately after reconstitution). 5 freeze/thaw cycles are possible. FbDP obtained from human plasma* diluted in glycine-albumin bovine buffer + preservatives. MLE data indicate the concentration in ng/mL (FEU) ("Calibrator (S1) Dose Value") and the confidence interval in "Relative Fluorescence Value" ("Calibrator (S1) RFV Range").
<u>DEX2 Controls:</u> C1 Control 2 x 2 mL (lyophilized) C2 Control 2 x 2 mL (lyophilized)	C1  C2	<b>Reconstitute with 2 mL of distilled water.</b> Wait for 5 minutes and then mix. <b>After reconstitution, the controls are stable for 28 days at 2-8°C or until the expiration date of the kit at -25 ± 6°C</b> (freeze immediately after reconstitution). 5 freeze/thaw cycles are possible. FbDP obtained from human plasma* diluted in glycine-albumin bovine buffer + preservatives. MLE data indicate the confidence interval in ng/mL (FEU) ("Control C1 Dose Value Range") or ("Control C2 Dose Value Range").
DEX2 Diluent 1 x 5 mL (liquid)	R1	Ready-to-use. TRIS buffer (0.05 mol/L pH 7.4) + calf serum + preservatives.
Specifications for the factory master data required to calibrate the test: • MLE data (Master Lot Entry) provided in the kit, or • MLE bar code printed on the box label.		
1 Package insert provided in the kit or downloadable from <a href="http://www.biomerieux.com/techlib">www.biomerieux.com/techlib</a> .		

\* This product has been tested and shown to be negative for HBs antigen, and antibodies to HIV1, HIV2 and HCV. However, since no existing test method can totally guarantee their absence, this product must be treated as potentially infectious. Therefore, usual safety procedures should be observed when handling.

**The SPR**

The SPR® is coated during production with anti-FbDP monoclonal immunoglobulins (mouse). Each SPR® is identified by the DEX2 code. Only remove the required number of SPR®s from the pouch and **carefully reseal the pouch after opening**.

**The Reagent Strip**

The strip consists of 10 wells covered with a labeled, foil seal. The label comprises a bar code which mainly indicates the assay code, kit lot number and expiration date. The foil of the first well is perforated to facilitate the introduction of the sample. The last well of each strip is a cuvette in which the fluorometric reading is performed. The wells in the center section of the strip contain the various reagents required for the assay.

**Description of the DEX2 strip**

Wells	Reagents
1	Sample well.
2 - 3 - 4	Empty wells.
5	Conjugate: alkaline phosphatase-labeled anti-FbDP monoclonal immunoglobulins (mouse) in TRIS buffer (0.05 mol/L pH 6.5) + horse serum + preservatives (400 µL).
6 - 7 - 9	Wash buffer: TRIS buffer (0.05 mol/L, pH 7.3) + chemical stabilizers + preservatives (600 µL).
8	Diluent: TRIS buffer (0.05 mol/L, pH 7.4) + calf serum + protein and chemical stabilizers + preservatives (600 µL).
10	Reading cuvette with substrate: 4-Methyl-umbelliferyl phosphate (0.6 mmol/L) + diethanolamine* (DEA) (0.62 mol/L or 6.6%) pH 9.2 + 1 g/L sodium azide (300 µL).

\* Signal Word: **ANGER**

Hazard statement

H318: Causes serious eye damage.

Precautionary statement

P280: Wear protective gloves/protective clothing/eye protection/face protection.

P305 + P351 + P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

For further information, refer to the Material Safety Data Sheet.

**MATERIALS AND DISPOSABLES REQUIRED BUT NOT PROVIDED**

- Single-use pipette and/or micropipettes to dispense the appropriate volumes.
- Powderless, disposable gloves.
- For other specific materials and disposables, please refer to the Instrument User Manual.
- Instrument of the VIDAS® family.

**WARNINGS AND PRECAUTIONS**

- **For *in vitro* diagnostic use only.**
- **For professional use only, by qualified laboratory personnel, in clinical laboratories.**
- This kit contains products of human origin. Source materials from which the controls and calibrator were derived were found negative when tested for HIV1, HIV2, HCV, and HBsAg. However, no known analysis method can totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (see Laboratory Biosafety Manual - WHO - Geneva - Latest Edition).
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (do not ingest; do not inhale).
- Do not use the SPR®s if the pouch is pierced or if the dot sealing a SPR® has come unstuck.
- Do not use visibly deteriorated STRs (damaged foil or plastic).
- Do not use reagents after the expiration date indicated on the box label.
- Do not mix reagents (or disposables) from different lots.
- **Use powderless gloves**, as powder has been reported to cause false results for certain enzyme immunoassay tests.
- Kit reagents contain sodium azide which can react with lead or copper plumbing to form explosive metal azides. If any liquid containing sodium azide is disposed of in the plumbing system, drains should be flushed with water to avoid build-up.
- **The substrate in well 10 contains an irritant agent (6.6% diethanolamine). Refer to the hazard statements "H" and the precautionary statements "P" above.**
- Spills should be wiped up thoroughly after treatment with liquid detergent or a solution of household bleach containing at least 0.5% sodium hypochlorite. See the User Manual for cleaning spills on or in the instrument. Do not autoclave solutions containing bleach.
- The instruments should be regularly maintained (refer to the User Manual for user and preventive maintenance operations).

**STORAGE CONDITIONS**

- Store the VIDAS® D-Dimer Exclusion II™ kit at 2-8°C.
- **Do not freeze reagents, with the exception of the calibrator and controls after reconstitution.**
- **Store all unused reagents at 2-8°C.**
- After opening the kit, check that the SPR® pouch is correctly sealed and undamaged. If not, do not use the SPR®s.
- **Carefully reseal the pouch with the desiccant inside after use to maintain stability of the SPR®s and return the complete kit to 2-8°C.**
- If stored according to the recommended conditions, all components are stable until the expiration date indicated on the label. Refer to the kit composition table for special storage conditions.

**SAMPLES****Sample type and collection**

- Collect blood by clean venipuncture **in trisodium citrate** (0.109 mol/L / 3.2% or 0.129 mol/L / 3.8%) or **CTAD** (sodium citrate, theophylline, adenosine and dipyridamole), observing the correct anticoagulant to blood ratio.
- Refer to the tube manufacturer's recommendations for use.
- Collection of whole blood using a syringe is not recommended in order to avoid formation of microclots in the sample.

**Sample preparation**

The current WHO/DIL/LAB/99.1 document provides recommendations for sample preparation.

For use of sample tubes, refer to the tube manufacturer's recommendations for use.

The pre-analytical step, including the preparation of blood samples, is an essential first step when performing medical analyses.

In accordance with Good Laboratory Practice, this step is performed under the responsibility of the laboratory manager.

**Sample stability**

Plasma samples separated from the pellet can be stored in aliquots at 2-8°C for 3 days. If longer storage is required, plasma can be frozen at  $-25 \pm 6^\circ\text{C}$  for 6 months, with up to 2 freeze/thaw cycles.

Thaw the plasma at 37°C when assaying.

**Sample-related interference**

Interference was studied according to the recommendations of the CLSI® document EP7-A2.

None of the following factors have been found to significantly influence this assay:

- hemolysis (after spiking samples with hemoglobin, up to 300 µmol/L (monomer)),
- lipemia (after spiking samples with lipids, up to 30 g/L equivalent in triglycerides),

- bilirubinemia (after spiking samples with bilirubin, up to 537 µmol/L),
- rheumatoid factor: up to 400 IU/mL (international units per milliliter),
- human albumin: up to 60 g/L.

However, it is recommended not to use samples that are clearly hemolyzed, lipemic or icteric and, if possible, to collect a new sample.

The influence of 51 analytes was also detected *in vitro*: no interference was observed.

Drug	Concentration tested	Drug	Concentration tested	Drug	Concentration tested
Acetaminophen	20 mg/dL	Cyclosporine A	0.4 mg/dL	Lithium chloride	14 mg/dL
Acetylsalicylic acid = aspirin	65 mg/dL	Dabigatran etexilate	18 mg/dL	L-Thyroxine	0.06 mg/dL
Allopurinol	4 mg/dL	Dextran 75	2500 mg/dL	Nicotine	0.1 mg/dL
Amikacin sulfate	10.4 mg/dL	Diazepam	0.5 mg/dL	Nifedipine	0.04 mg/dL
Ampicillin	5 mg/dL	Digoxin	0.0006 mg/dL	Penicillin G sodium salt	2500 U/dL
Apixaban	0.6 mg/dL	D-L methyl dopa hydrochloride	1.8 mg/dL	Pentobarbital	9 mg/dL
Ascorbic Acid (vitamin C)	6 mg/dL	Dopamine hydrochloride	0.1 mg/dL	Phenobarbital	10 mg/dL
Atenolol	1 mg/dL	Edoxaban	3.6 mg/dL	Phenytoine	5 mg/dL
Caffeine	6 mg/dL	Erythromycin	6 mg/dL	Primidone	4 mg/dL
Captopril	0.5 mg/dL	Ethanol	400 mg/dL	Propranolol hydrochloride	0.2 mg/dL
Carbamazepine	3 mg/dL	Ethosuximide	25 mg/dL	Rivaroxaban	1.2 mg/dL
Chloramphenicol	5 mg/dL	Furosemide	6 mg/dL	Theophylline	4 mg/dL
Chlordiazepoxide hydrochloride	1.1 mg/dL	Gentamicin sulfate	1 mg/dL	Urea	500 mg/dL
Chlorpromazine hydrochloride	0.2 mg/dL	Lithium heparin	300 U/dL	Uric acid	24 mg/dL
Cimetidine	2 mg/dL	Sodium heparin	300 U/dL	Valproic acid sodium salt	60 mg/dL
Cinnarizine	3 mg/dL	Ibuprofen	50 mg/dL	Verapamil hydrochloride	0.2 mg/dL
Creatinine	30 mg/dL	Lidocaine	1.2 mg/dL	Warfarin	1 mg/dL

## INSTRUCTIONS FOR USE

For complete instructions, see the User Manual.

### Reading VIDAS® PTC (Protocol Test Change) data and MLE data

#### When using the assay for the first time:

With the external instrument barcode reader, **scan the barcodes (PTC and MLE) in the following order:**

1. According to the instrument used, scan the PTC barcode(s) at the end of the package insert or downloadable from [www.biomerieux.com/techlib](http://www.biomerieux.com/techlib). This reading allows VIDAS® PTC protocol data to be transferred to the instrument software for its update.
2. Scan the MLE data in the kit or on the box label.

#### When opening a new lot of reagents:

With the external instrument barcode reader, scan the MLE data on the box label before performing the test.

**Note: the master lot data need only be entered once for each lot.**

It is possible to enter MLE data **manually or automatically** depending on the instrument (refer to the User Manual).

### Calibration

Calibration, using **the calibrator** provided in the kit, must be performed each time a new lot of reagents is opened, after the MLE data have been entered. Calibration should then be performed **every 28 days**. This operation provides instrument-specific calibration curves and compensates for possible minor variations in assay signal throughout the shelf-life of the kit.

**The calibrator, identified by S1, must be tested in duplicate** (see User Manual).

The calibrator value must be within the set RFV "Relative Fluorescence Value" range. If this is not the case, recalibrate.

### Kit controls

Two controls are included in each VIDAS® D-Dimer Exclusion II™. These controls must be performed immediately after opening a new kit to ensure that reagent performance has not been altered. Each calibration must also be checked using these controls. The instrument will only be able to check the control values if they are identified by C1 and C2.

Results cannot be validated if the control values deviate from the expected values.

### Procedure

1. **Remove the kit from storage at 2-8°C and take out the required reagents. Carefully reseal the SPR® pouch and return the kit to 2-8°C. The reagents can be used immediately.**
2. Use one "DEX2" strip and one "DEX2" SPR® from the kit for each sample, control or calibrator to be tested. Make sure the storage pouch has been carefully resealed after the required SPR®s have been removed.
3. The test is identified by the "DEX2" code on the instrument. The calibrator must be identified by "S1", and tested in duplicate. If the controls are to be tested, they should be identified by "C1" and "C2" and tested singly.
4. If necessary, clarify samples by centrifugation.
5. Mix the calibrator, controls and samples using a vortex-type mixer in order to improve result reproducibility (for plasma separated from the pellet).
6. For optimum results, refer to all the paragraphs in the **SAMPLES** section.
7. Before pipetting ensure that samples, calibrators, controls and diluent are free of bubbles.
8. For this test, the calibrator, control and sample test portion is 200 µL.
9. Insert the "DEX2" SPRs and "DEX2" strips into their appropriate position on the instrument. Check to make sure the color labels with the assay code on the SPR®s and the STRs match.
10. **Initiate the assay immediately** (refer to the User Manual). All the assay steps are performed automatically by the instrument.
11. Reclose the vials and return them to the recommended temperature after pipetting.
12. The assay results are available within 20 minutes. After the assay is completed, remove the SPR®s and STRs from the instrument.
13. Dispose of the used SPR®s and STRs into an appropriate recipient.

### LIMITATIONS OF THE METHOD

- Interference may be encountered with certain plasmas containing antibodies directed against reagent components. For this reason, assay results should be interpreted taking into consideration the patient's history (clinical probability).
- Clinical performance data were determined on an outpatient population. Because D-dimer results are likely to be elevated in an inpatient population due to stasis, chronic illness, post-surgery and other non-specific conditions known to elevate D-dimer levels, the clinical utility of a negative result is not likely to be realized in an inpatient population. Therefore, clinical performance results should not be extrapolated to an inpatient population.
- Results below the detection limit of 45 ng/mL obtained with the VIDAS® D-Dimer Exclusion II™ kit are usually linked to poor pre-analytical conditions or incorrect instrument maintenance, and therefore, the assay must be repeated. Please contact your customer service for further assistance.

### QUALITY CONTROL

Additional quality controls can be performed in accordance with local regulations or requirements related to accreditation, as well as requirements defined in the laboratory's quality control procedure.

### RESULTS AND INTERPRETATION

Once the assay is completed, the results are analyzed automatically by the computer. The results are calculated using the calibration curve that is stored by the instrument (4-parameter logistics model), and then printed. The D-dimer concentrations are expressed in ng/mL of Fibrinogen Equivalent Unit (FEU).

VIDAS® D-Dimer Exclusion II™ has been calibrated against an internal panel of human plasma, the concentrations of which have been determined using the VIDAS D-Dimer Exclusion kit (ref. 30 442). This reference (30 442) is no longer sold.

Samples with D-dimer concentrations greater than 10,000 ng/mL (FEU) can be retested after being diluted by 1/5 in the kit diluent. **If the dilution factor was entered when the worklist was created, the result is calculated automatically. If the dilution factor was not entered, multiply the results by the dilution factor to obtain the sample concentration.**

#### Description of clinical decision cutoffs

- **To exclude DVT and PE in combination with clinical pretest probability assessment.**

For deep vein thrombosis, the **cutoff level is 500 ng/mL.**

For pulmonary embolism, the **cutoff level is either 500 ng/mL, or age-adjusted, as follows:**

- **< 50 years: 500 ng/mL cutoff**
- **≥ 50 years: age x 10 ng/mL** (example: 650 ng/mL cutoff for 65 years)

With a decision cutoff of 500 ng/mL, a D-dimer result ≥ 500 ng/mL (FEU) is considered to be positive and a result < 500 ng/mL (FEU) is considered to be negative.

With an age-adjusted cutoff, a D-dimer result ≥ the age-adjusted cutoff is considered to be positive and a result < the age-adjusted cutoff is considered to be negative.

- **To aid in evaluating the VTE recurrence rate in women with a first unprovoked VTE, in order to guide the duration of oral anticoagulant therapy.**

The cutoff level used in the HERDOO2 clinical decision rule is **250 ng/mL**. Women with a first unprovoked VTE who are on oral anticoagulants are considered positive (1 point in the clinical decision rule) with a VIDAS® D-Dimer Exclusion II™ test result ≥ 250 ng/mL and negative (0 point in the clinical decision rule) with a test result < 250 ng/mL.

- HERDOO2 rule: for postmenopausal women ( $\geq 50$  years) with 0 or 1 HERDOO2 criteria who discontinued oral anticoagulants (n=126), the recurrence rate was 5.7% (95% CI 2.6-10.9%) (15). Bearing in mind that the ISTH suggests that it is safe to discontinue anticoagulants if the recurrence rate is less than 5% one year after stopping the treatment (14), it is all the more important that clinicians take into consideration the clinical history of patients in this group before making any decisions.

## REFERENCE VALUES

In a study carried out using 215 citrated plasma samples from blood donors, 90% of values were below 500 ng/mL (FEU).

It is recommended that each laboratory establish its own reference values from a rigorously selected population.

## PERFORMANCE

Studies performed using VIDAS® D-Dimer Exclusion II™ gave the following results:

### Measurement range

The measurement range of the VIDAS® D-Dimer Exclusion II™ assay is 45 ng/mL (FEU) to 10,000 ng/mL (FEU) (upper limit of quantification).

### Linearity

Linearity was determined on plasmas according to the recommendations of CLSI® document EP6-A. The results obtained have shown that VIDAS® D-Dimer Exclusion II™ is linear from 45 to 10,000 ng/mL (FEU). However, for a few isolated samples, VIDAS® D-Dimer Exclusion II™ may not be linear due to the matrix.

### Detection limit

Defined as the smallest concentration of D-dimer significantly different from the zero concentration with a risk  $\alpha$  of 5%: **< 45 ng/mL (FEU)**. The results of the detection limit were determined according to the standard CLSI® EP17-A.

### Hook effect

No hook effect was found up to D-dimer concentrations of 400,000 ng/mL (FEU).

### Specificity

Interference was studied according to the recommendations of the CLSI® document EP7-A2.

Cross reactant	Concentration (spiked samples)	Cross reactivity
Fibrinogen	$\leq 10$ g/L	No
Fibrinogen degradation products X	$\leq 10$ µg/mL	No
Fibrinogen degradation products Y	$\leq 10$ µg/mL	No
Fibrinogen degradation products D	10 - 100 µg/mL	Yes*

\*: such high levels of fibrinogen degradation products D do not occur in the target population of suspected VTE patients. They may be observed when patients are subjected to therapeutic lysis with thrombolytic agents (16).

**Note:** The specificity of the two antibodies used in this assay has not been tested against fibrinogen degradation products E, therefore cross-reactivity cannot be ruled out.

### Precision

Three samples were tested in duplicate in 40 different runs (2 runs per day) with 2 reagent lots at 3 sites (n=240).

The repeatability (intra-run precision), and between lot reproducibility were calculated using this protocol, based on the recommendations of CLSI® document EP5-A2:

		Repeatability		Between lot reproducibility	
Sample	Mean concentration ng/mL (FEU)	Standard deviation	CV (%)	Standard deviation	CV (%)
Sample 1	277.97	6.88	2.5	18.25	6.6
Sample 2	544.14	11.05	2.0	32.13	5.9
Sample 3	7,788.88	113.83	1.5	468.18	6.0

## Comparison with another test method

### Study 1:

328 fresh citrated plasmas from D-dimer samples coming from the laboratory routine activity of one European site and two North American sites, were tested using two lots of VIDAS® D-Dimer Exclusion II™ and five lots of VIDAS® D-Dimer Exclusion according to the recommendations of CLSI® document EP9-A2. The concentration range of plasma samples investigated is from 58.46 to 9,594.96 ng/mL (FEU). The following results were obtained:

#### • Passing & Bablok regression and correlation coefficient (r):

Sites		n	Slope	Intercept	r
Europe	Site 1	158	1.20	-41.96	0.991
	Site 2	43	1.15	-10.31	0.988
North America		125	1.20	-43.18	0.989
All		326*	1.19	-34.92	0.987

\* 2 statistical outliers were removed from the analysis.

### Study 2:

378 citrated plasmas coming from the laboratory routine activity of one European site, were tested using two lots of VIDAS® D-Dimer Exclusion II™ and two lots of VIDAS® D-Dimer Exclusion according to the recommendations of CLSI® document EP9-A2. The concentration range of plasma samples investigated is from 46.70 to 9,913.53 ng/mL (FEU). The following results were obtained:

#### • Passing & Bablok regression and correlation coefficient (r):

Site 4 Europe	n	Slope	Intercept	r
	374*	1.09	-29.35	0.991

\* 4 statistical outliers were removed from the analysis

In conclusion, these two studies show a **mean slope of 1.14** between VIDAS® D-Dimer Exclusion II™ and VIDAS® D-Dimer Exclusion, with slopes between 1.09 and 1.19 depending on the lots.

## Clinical performance

### Clinical performance using the VIDAS® D-Dimer Exclusion II™ assay in conjunction with a patient's clinical pretest probability to exclude VTE (venous thromboembolism)

A study was performed using 315 frozen **clinically characterized** samples from patients from **previous prospective clinical studies** (17, 18).

The overall prevalence of VTE in the total population studied was 23.5% (74/315). The sensitivity, specificity, negative predictive value (NPV) and positive predictive value (PPV) using a clinical cutoff of 500 ng/mL (FEU) are summarized below with the corresponding 95% confidence intervals (CI).

		Patients with suspected VTE		
		Low and intermediate pretest probability n= 303	High pretest probability n= 12	All probabilities n= 315
% Sensitivity (95% CI)	VIDAS® D-Dimer Exclusion II™	100% (94.2 - 100)	100% (73.5 - 100)	100% (95.1 - 100)
	VIDAS® D-Dimer Exclusion	100% (94.2 - 100)	100% (73.5 - 100)	100% (95.1 - 100)
% Specificity (95% CI)	VIDAS® D-Dimer Exclusion II™	35.7% (29.6 - 42.1)	N/A	35.7% (29.6 - 42.1)
	VIDAS® D-Dimer Exclusion	37.8% (31.6 - 44.2)	N/A	37.8% (31.6 - 44.2)
% NPV (95% CI)	VIDAS® D-Dimer Exclusion II™	100% (95.8 - 100)	N/A	100% (95.8 - 100)
	VIDAS® D-Dimer Exclusion	100% (96.0 - 100)	N/A	100% (96.0 - 100)
% PPV (95% CI)	VIDAS® D-Dimer Exclusion II™	28.6% (22.7 - 35.1)	100% (73.5 - 100)	32.3% (26.3 - 38.8)
	VIDAS® D-Dimer Exclusion	29.2% (23.2 - 35.9)	100% (73.5 - 100)	33.0% (26.9 - 39.6)

\* 95% CI calculated using the exact method

**Agreement study** between the two methods (VIDAS® D-Dimer Exclusion II™ and VIDAS® D-Dimer Exclusion) and the **clinical status** are presented below (n= 315 frozen samples):

VIDAS® D-Dimer Exclusion II™	VIDAS® D-Dimer Exclusion	Total samples	Clinical Status	
			NEGATIVE (NO VTE)	POSITIVE (VTE)
< 500 ng/mL	< 500 ng/mL	85	85	0
	≥ 500 ng/mL	1	1	0
≥ 500 ng/mL	< 500 ng/mL	6	6	0
	≥ 500 ng/mL	223	149	74
<b>Total samples</b>		<b>315</b>	<b>241</b>	<b>74</b>

**There is no significant difference at a risk level of 5% between the sensitivity and the specificity of the two assays.**

A retrospective analysis using the age-adjusted cutoff was performed on patients suspected of pulmonary embolism (PE) with low and intermediate pretest probability. Among these patients, the observed prevalence of PE was 20.5% (62/303). The sensitivity, specificity, negative predictive value (NPV) and positive predictive value (PPV) using the age-adjusted clinical cutoff are summarized below with the corresponding 95% confidence intervals (CI).

	Patients with suspected EP
	Low and intermediate pretest probability n=303
% Sensitivity (95% CI*)	100% (94.2 - 100)
% Specificity (95% CI*)	40.2% (34.3 - 46.5)
% NPV (95% CI*)	100% (96.3 - 100)
% PPV (95% CI*)	30.1% ( 24.2 - 36.7)

\* IC: Confidence Interval

For the VIDAS® D-Dimer Exclusion II™ assay, the specificity using the age adjusted cutoff is significantly higher than the specificity using the clinical cutoff at 500 ng/mL (FEU). There is no significant difference at a risk level of 5%, between the sensitivity using the age adjusted cutoff and the sensitivity using the clinical cutoff at 500 ng/mL (FEU).

#### **Clinical performance using the VIDAS® D-Dimer Exclusion assay in combination with a patient's pretest probability**

##### **❖ To exclude DVT**

Frozen samples collected from patients enrolled in a **multi-center prospective cohort study** (19) were used to validate the diagnostic utility of VIDAS® D-Dimer Exclusion to exclude a diagnosis of deep vein thrombosis (DVT). Consecutive eligible outpatients (n= 556) with a first suspected DVT episode were evaluated at three hospitals during the course of the study. Using the Wells model to estimate the probability of DVT, patients were classified as having a high, intermediate or low pretest probability (PTP) of DVT (20).

The VIDAS® D-Dimer Exclusion assay was performed without knowledge of the PTP assessment results and the clinical outcome of the patients from which the samples were derived. A clinical cutoff level of 500 ng/mL (FEU) was used.

A D-dimer result ≥ 500 ng/mL (FEU) was considered positive and a result < 500 ng/mL (FEU) was considered negative. During the study, patients were grouped according to PTP. Patients with a negative D-dimer test result and a low or intermediate PTP of DVT underwent no further diagnostic testing and were followed up for 3 months for development of DVT. Patients with a positive D-dimer test result or high PTP underwent serial compression ultrasonography.

The overall prevalence of DVT among the total population studied was 10.1% (56/556). One sample from the original study was not tested due to volume limitations. The sensitivity, specificity and negative predictive value (NPV) of the VIDAS® D-Dimer Exclusion assay using a clinical cutoff of 500 ng/mL (FEU) are summarized below with the corresponding 95% confidence intervals (CI).

	Patients with suspected DVT			
	Low pretest probability n= 295	Intermediate pretest probability n= 189	High pretest probability n= 71	All probabilities n= 555
% Sensitivity (95% CI)	100% (81.5% – 100%)	100% (80.5% – 100%)	100% (83.9% – 100%)	100% (93.6% – 100%)
% Specificity (95% CI)	39.7% (33.9% – 45.7%)	26.7% (20.3% – 34.0%)	16.0% (7.2% – 29.1%)	32.9 % (28.8% – 37.2 %)
% NPV (95% CI)	100% (96.7% – 100%)	100% (92.3% – 100.0%)	100% (63.1% – 100%)	100% (97.8% – 100%)

#### ❖ To exclude PE

The diagnostic utility of VIDAS® D-Dimer Exclusion to exclude a diagnosis of PE was validated in a **multi-center prospective cohort study** using samples from 1,290 consecutive patients presenting to the emergency department with suspected PE (21). Of these 1,290 patients enrolled, 325 were eventually excluded for a total of 965 patients included in the final analysis.

All patients enrolled were classified, using the Geneva score, as having a high, intermediate or low pretest probability (PTP) of PE (22).

D-dimer results  $\geq$  the cutoff value of 500 ng/mL (FEU) were considered positive, and results  $<$  500 ng/mL (FEU) were considered negative.

All patients with negative D-dimer test results received no treatment and underwent no further diagnostic testing. Patients with positive D-dimer test results were further evaluated using ultrasound, helical CT scan, and/or angiography. These patients were then treated based upon the results of the additional diagnostic testing. All patients were followed up for a period of three months to evaluate any possible venous thromboembolic events (DVT or PE) and episodes of bleeding.

The overall prevalence of PE among the total population studied was 23.0% (222/965). The sensitivity, specificity, negative predictive value (NPV) and positive predictive value (PPV) using a clinical cutoff value of 500 ng/mL (FEU) are summarized below with the corresponding 95% confidence intervals (CI).

	Patients with suspected PE		
	Low and intermediate pretest probability n= 891	High pretest probability n= 74	All probabilities n= 965
% Sensitivity (95% CI)	100% (97.7% - 100%)	100% (94.3% - 100%)	100% (98.4% - 100%)
% Specificity (95% CI)	37.6% (34.0 % - 41.2%)	45.5% (16.7% - 76.6%)	37.7% (34.2% - 41.3%)
% NPV (95% CI)	100% (98.7% -100%)	100% (47.8% - 100%)	100% (98.7% - 100%)
% PPV (95% CI)	25.8% (22.4% - 29.5%)	91.3% (82.0% - 96.7%)	32.4% (28.9% - 36.1%)

A retrospective analysis using the age-adjusted cutoff was performed on patients suspected of pulmonary embolism (PE) with low and intermediate pretest probability. Among these patients, the observed prevalence of PE was 17.8% (159/891). The sensitivity, specificity, negative predictive value (NPV) and positive predictive value (PPV) using the age-adjusted clinical cutoff, are summarized below with the corresponding 95% confidence intervals (CI)

Patients with suspected EP	
Low and intermediate pretest probability n=891	
% Sensitivity (95% CI*)	98.1% (94.6 - 99.6)
% Specificity (95% CI*)	45.4% (41.8 - 49.0)
% NPV (95% CI*)	99.1% (97.4 - 99.8)
% PPV (95% CI*)	28.1% (24.5 - 31.9)

\* CI: Confidence Interval

For the VIDAS® D-Dimer Exclusion II™ assay, the specificity using the age-adjusted cutoff is significantly higher than the specificity using the clinical cutoff at 500 ng/mL (FEU). There is no significant difference at a risk level of 5%, between the sensitivity using the age-adjusted cutoff and the sensitivity using the clinical cutoff at 500 ng/mL (FEU).

#### **Clinical performance of the VIDAS® D-Dimer Exclusion assay and the VIDAS® D-Dimer Exclusion II™ assay in conjunction with a patient's pretest probability to exclude PE (Pulmonary Embolism)**

A multi-center prospective study (19 hospitals) was conducted, including 3,324 patients with suspected PE (10). The prevalence of PE among the studied population was 19% (631/3324).

D-dimer was measured using the VIDAS® D-Dimer Exclusion and VIDAS® D-Dimer Exclusion II™ assays on 1,345 patients with a non-high or unlikely pretest probability of PE.

Patients with a non-high or unlikely pretest probability of PE and a D-dimer result below their age-adjusted cutoff, received no treatment and underwent no further diagnostic testing. Failure rate in these patients was assessed by a 3-month follow-up period with all suspected venous thromboembolic events and deaths adjudicated by an independent committee.

In this study, for the VIDAS® D-Dimer Exclusion and the VIDAS® D-Dimer Exclusion II™ assays, the PE exclusion rate increased significantly from 31.4% (423/1345) with the cutoff at 500 ng/mL, to 41.1% (553/1345) with the age-adjusted cutoff, which represents a relative increase of 30.7%. For the VIDAS® D-Dimer Exclusion and the VIDAS® D-Dimer Exclusion II™ assays, the 3-month thromboembolic failure rate in patients with D-dimer concentrations > 500 ng/mL but < the age-adjusted cutoff, was 0.0% (95% CI: [ 0.0 – 2.9]).

#### **Clinical performance of the VIDAS® D-Dimer Exclusion assay and the VIDAS® D-Dimer Exclusion II™ assay in the context of the use of the HERDOO2 clinical decision rule**

The HERDOO2 clinical decision rule was validated during an international prospective multi-center clinical study (44 hospitals) (15). 2,785 patients having suffered from unprovoked VTE and having received oral anticoagulants for 5 to 12 months, were recruited (1,572 men and 1,213 women). The HERDOO2 score was calculated for the women in order to differentiate between those with a high or low risk of recurrent VTE. D-dimer measurement was performed using the VIDAS® D-Dimer Exclusion or VIDAS® D-Dimer Exclusion II™ assay while the patients were still taking anticoagulants.

Thanks to the HERDOO2 clinical decision rule, 591 of the 631 women classified as having a low risk of recurrence were able to discontinue anticoagulant therapy. The one-year VTE recurrence rate in women having suffered from unprovoked VTE and having stopped oral anticoagulant therapy after 5 to 12 months was 3.0% (CI 95%: [1.8 ; 4.8]). This recurrence rate is lower than the recommended 5% defined by the ISTH (14).

## WASTE DISPOSAL

Dispose of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.

It is the responsibility of each laboratory to handle waste and effluents produced according to their nature and degree of hazardousness and to treat and dispose of them (or have them treated and disposed of) in accordance with any applicable regulations.

## LITERATURE REFERENCES

1. DEMPFLÉ CE. Validation, calibration, and specificity of quantitative D-dimer assays. *Semin Vasc Med.* 2005;5:315-20.
2. MOSESSON MW. On behalf of the Subcommittee on Fibrinogen of the Scientific and Standardization Committee of the ISTH. Terminology for macromolecular derivatives of crosslinked fibrin. *Thromb Haemost.* 1995;73:725-6.
3. BOCKENSTEDT P. D-dimer in venous thromboembolism. *N Engl J Med.* 2003;349:1203-4.
4. DI NISIO M, SQUIZZATO A, RUTJES AW, BÜLLER HR, ZWINDERMAN AH, BOSSUYT PM. Diagnostic accuracy of D-dimer test for exclusion of venous thromboembolism: a systematic review. *J Thromb Haemost.* 2007;5:296-304.
5. RIGHINI M, PERRIER A, DE MOERLOOSE P, BOUNAMEAUX H. D-Dimer for venous thromboembolism diagnosis: 20 years later. *J Thromb Haemost.* 2008;6:1059-71.
6. TEN CATE-HOEK AJ, PRINS MH. Management studies using a combination of D-dimer test result and clinical probability to rule out venous thromboembolism: a systematic review. *J Thromb Haemost.* 2005;3:2465-70.
7. CARRIER M, RIGHINI M, DJURABI RK, HUISMAN MV, PERRIER A, WELLS PS, RODGER M, WUILLEMIN WA, LE GAL G. VIDAS D-dimer in combination with clinical pre-test probability to rule out pulmonary embolism. A systematic review of management outcome studies. *Thromb Haemost.* 2009;101:886-92.
8. GOODACRE S, SAMPSON F, STEVENSON M, WAILOO A, SUTTON A, THOMAS S, LOCKER T, RYAN A. Measurement of the clinical and cost-effectiveness of non-invasive diagnostic testing strategies for deep vein thrombosis. *Health Technol Assess.* 2006;10:1-168, iii-iv.
9. RIGHINI M, NENDAZ M, LE GAL G, BOUNAMEAUX H, PERRIER A. Influence of age on the cost-effectiveness of diagnostic strategies for suspected pulmonary embolism. *J Thromb Haemost.* 2007;5:1869-77.
10. RIGHINI M, et al. Age-adjusted D-dimer cutoff levels to rule out pulmonary embolism: the ADJUST-PE study. *JAMA.* 2014;311:1117-24.
11. BRUINSTROOP E, VAN DE REE MA, HUISMAN MV. The use of D-dimer in specific clinical conditions: a narrative review. *Eur J Intern Med.* 2009;20:441-6.
12. HOUDIJK W. Proper observation of patient-related factors is an important determinant in the use of the D-dimer test for exclusion of venous thromboembolism in the ED. *Am J Emerg Med.* 2007;25:255-256.
13. RODGER MA, et al. Identifying unprovoked thromboembolism patients at low risk for recurrence who can discontinue anticoagulant therapy. *CMAJ.* 2008;179:417-26.
14. KEARON C, IORIO A, PALARETI G; Subcommittee on Control of Anticoagulation of the SSC of the ISTH. Risk of recurrent venous thromboembolism after stopping treatment in cohort studies: recommendation for acceptable rates and standardized reporting. *J Thromb Haemost.* 2010;8:2313-5.
15. RODGER et al. Validating the HERDOO2 rule to guide treatment duration for women with unprovoked venous thrombosis: multinational prospective cohort management study. *BMJ.* 2017 Mar 17;356:j1065.
16. WALKER JB, NESHEIM ME. The molecular weights, mass distribution, chain composition, and structure of soluble fibrin degradation products released from a fibrin clot perfused with plasmin. *J Biol Chem.* 1999; 274: 5201-12.
17. PERRIER A, et al. Multidetector-row computed tomography in suspected pulmonary embolism. *N Engl J Med.* 2005;352:1760-8.
18. RIGHINI M, et al. Diagnosis of pulmonary embolism by multidetector CT alone or combined with venous ultrasonography of the leg: a randomised non-inferiority trial. *Lancet.* 2008;371:1343-52.
19. BATES S M. et al. A diagnostic strategy involving a quantitative latex D-Dimer assay reliably excludes deep venous thrombosis. *Ann Intern Med.* 2003;138: 787-94.
20. WELLS PS, et al. Value of assessment of pre-test-probability of deep vein thrombosis in clinical management. *Lancet.* 1997; 350: 1795-8.
21. PERRIER A, et al. Diagnosing pulmonary embolism in outpatients with clinical assessment, D-Dimer measurement, venous ultrasound, and helical computed tomography: a multicenter management study. *Am J Med.* 2004; 116:291-9.
22. WICKI J, et al. Assessing clinical probability of pulmonary embolism in the emergency ward: a simple score. *Arch Intern Med.* 2001;161:92-7.

**INDEX OF SYMBOLS**

Symbol	Meaning
	Catalog number
	<i>In Vitro</i> Diagnostic Medical Device
	Manufacturer
	Temperature limit
	Use by date
	Batch code
	Consult Instructions for Use
	Contains sufficient for <n> tests
	Date of manufacture

**LIMITED WARRANTY**

bioMérieux warrants the performance of the product for its stated intended use provided that all procedures for usage, storage and handling, shelf life (when applicable), and precautions are strictly followed as detailed in the instructions for use (IFU).

Except as expressly set forth above, bioMérieux hereby disclaims all warranties, including any implied warranties of merchantability and fitness for a particular purpose or use, and disclaims all liability, whether direct, indirect or consequential, for any use of the reagent, software, instrument and disposables (the "System") other than as set forth in the IFU.

**REVISION HISTORY**Change type categories:

N/A	Not applicable (First publication)
Correction	Correction of documentation anomalies
Technical change	Addition, revision and/or removal of information related to the product
Administrative	Implementation of non-technical changes noticeable to the user

**Note:** *Minor typographical, grammar, and formatting changes are not included in the revision history.*

Release date	Part Number	Change Type	Change Summary
2015/01	14219E	Administrative	INDEX OF SYMBOLS REVISION HISTORY
		Technical	CONTENT OF THE KIT (60 TESTS) – RECONSTITUTION OF REAGENTS WARNINGS AND PRECAUTIONS
2015/06	14219F	Technical	CONTENT OF THE KIT (60 TESTS) – RECONSTITUTION OF REAGENTS INSTRUCTIONS FOR USE
2017/03	21967A	Technical	REFERENCE PRODUCT SUMMARY AND EXPLANATION PRINCIPLE CONTENT OF THE KIT (60 TESTS) - RECONSTITUTION OF REAGENTS MATERIALS AND DISPOSABLES REQUIRED BUT NOT PROVIDED WARNINGS AND PRECAUTIONS SAMPLES INSTRUCTIONS FOR USE QUALITY CONTROL RESULTS AND INTERPRETATION LIMITATIONS OF THE METHOD REFERENCE VALUES PERFORMANCE LITERATURE REFERENCES LIMITED WARRANTY

BIOMERIEUX, the blue logo, VIDAS and SPR are used, pending and/or registered trademarks belonging to bioMérieux, or one of its subsidiaries, or one of its companies.

CLSI is a trademark belonging to Clinical Laboratory and Standards Institute, Inc.

Any other name or trademark is the property of its respective owner.

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1. DEMPFLÉ CE. Validation, calibration, and specificity of quantitative D-dimer assays. *Semin Vasc Med.* 2005;5:315-20.
2. MOSESSON MW. On behalf of the Subcommittee on Fibrinogen of the Scientific and Standardization Committee of the ISTH. Terminology for macromolecular derivatives of crosslinked fibrin. *Thromb Haemost.* 1995;73:725-6.
3. BOCKENSTEDT P. D-dimer in venous thromboembolism. *N Engl J Med.* 2003;349:1203-4.
4. DI NISIO M, SQUIZZATO A, RUTJES AW, BÜLLER HR, ZWINDERMAN AH, BOSSUYT PM. Diagnostic accuracy of D-dimer test for exclusion of venous thromboembolism: a systematic review. *J Thromb Haemost.* 2007;5:296-304.
5. RIGHINI M, PERRIER A, DE MOERLOOSE P, BOUNAMEAUX H. D-Dimer for venous thromboembolism diagnosis: 20 years later. *J Thromb Haemost.* 2008;6:1059-71.
6. TEN CATE-HOEK AJ, PRINS MH. Management studies using a combination of D-dimer test result and clinical probability to rule out venous thromboembolism: a systematic review. *J Thromb Haemost.* 2005;3:2465-70.
7. CARRIER M, RIGHINI M, DJURABI RK, HUISMAN MV, PERRIER A, WELLS PS, RODGER M, WUILLEMIN WA, LE GAL G. VIDAS D-dimer in combination with clinical pre-test probability to rule out pulmonary embolism. A systematic review of management outcome studies. *Thromb Haemost.* 2009;101:886-92.
8. GOODACRE S, SAMPSON F, STEVENSON M, WAILOO A, SUTTON A, THOMAS S, LOCKER T, RYAN A. Measurement of the clinical and cost-effectiveness of non-invasive diagnostic testing strategies for deep vein thrombosis. *Health Technol Assess.* 2006;10:1-168, iii-iv.
9. RIGHINI M, NENDAZ M, LE GAL G, BOUNAMEAUX H, PERRIER A. Influence of age on the cost-effectiveness of diagnostic strategies for suspected pulmonary embolism. *J Thromb Haemost.* 2007;5:1869-77.
10. RIGHINI M, et al. Age-adjusted D-dimer cutoff levels to rule out pulmonary embolism: the ADJUST-PE study. *JAMA.* 2014;311:1117-24.
11. BRUINSTROOP E, VAN DE REE MA, HUISMAN MV. The use of D-dimer in specific clinical conditions: a narrative review. *Eur J Intern Med.* 2009;20:441-6.
12. HOUDIJK W. Proper observation of patient-related factors is an important determinant in the use of the D-dimer test for exclusion of venous thromboembolism in the ED. *Am J Emerg Med.* 2007;25:255-256.
13. RODGER MA, et al. Identifying unprovoked thromboembolism patients at low risk for recurrence who can discontinue anticoagulant therapy. *CMAJ.* 2008;179:417-26.
14. KEARON C, IORIO A, PALARETI G; Subcommittee on Control of Anticoagulation of the SSC of the ISTH. Risk of recurrent venous thromboembolism after stopping treatment in cohort studies: recommendation for acceptable rates and standardized reporting. *J Thromb Haemost.* 2010;8:2313-5.
15. RODGER M, et al. REVERSE II Study: Multi-national validation of the men continue and HERDOO2 rule to identify low risk unprovoked venous thromboembolism patients who can discontinue anticoagulants. *BMJ.* 2017, in press.
16. WALKER JB, NESHEIM ME. The molecular weights, mass distribution, chain composition, and structure of soluble fibrin degradation products released from a fibrin clot perfused with plasmin. *J Biol Chem.* 1999; 274: 5201-12.
17. PERRIER A, et al. Multidetector-row computed tomography in suspected pulmonary embolism. *N Engl J*
18. RIGHINI M, et al. Diagnosis of pulmonary embolism by multidetector CT alone or combined with venous ultrasonography of the leg: a randomised non-inferiority trial. *Lancet.* 2008;371:1343-52.
19. BATES S M, et al. A diagnostic strategy involving a quantitative latex D-Dimer assay reliably excludes deep venous thrombosis. *Ann Intern Med.* 2003;138: 787-94.
20. WELLS PS, et al. Value of assessment of pre-test-probability of deep vein thrombosis in clinical management. *Lancet.* 1997; 350: 1795-8.
21. PERRIER A, et al. Diagnosing pulmonary embolism in outpatients with clinical assessment, D-Dimer measurement, venous ultrasound, and helical computed tomography: a multicenter management study. *Am J Med.* 2004; 116:291-9.
22. WICKI J, et al. Assessing clinical probability of pulmonary embolism in the emergency ward: a simple score. *Arch Intern Med.* 2001;161:92-7.

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# Quality Control VIDAS<sup>®</sup> (QCV)

Automated test for use on the VIDAS<sup>®</sup> system to detect abnormal operation of the VIDAS and mini VIDAS instrument pipette mechanisms and optical systems.

## SUMMARY AND EXPLANATION

Quality Control VIDAS (QCV) is used to detect abnormal operation of pipette mechanisms which may affect the results of biological tests. It is also intended for checking that the optical system is capable of measuring high fluorescence levels.

**The QCV test must be run at least once a month in each position** of the VIDAS and mini VIDAS instruments, or any time that one of the above-mentioned instrument problems is suspected.

## PRINCIPLE

The test corresponds to successive aspirations/dilutions of fluorescent substrate (4-Methyl-umbelliferone) solutions which are standardized and have varying levels of concentration. The aspirations are performed at different speeds to check pump aspiration capability.

The results obtained at the end of the test correspond to:

- a fluorescence ratio "Test Value 1" (TV1 for mini VIDAS) for checking that the pipette mechanism is operating correctly,
- a fluorescence reading (Reading 3 for VIDAS and R3 for mini VIDAS) for checking that the optical system is intact.

The values obtained must be within the acceptable range defined later in this package insert.

## CONTENT OF THE KIT (60 TESTS):

60 QCV Strips	STR	Ready-to-use.
60 QCV SPRs (2 x 30)	SPR	Ready-to-use.
1 Package insert provided in the kit or downloadable from <a href="http://www.biomerieux.com/techlib">www.biomerieux.com/techlib</a>		

## DESCRIPTION

### The SPR<sup>®</sup>:

The SPR is used as a pipetting device and is identified by the code QCV. **Only remove the required number of SPRs from the pouch.**

### The strip:

The strip consists of 10 wells covered with a labeled, foil seal. The label comprises a bar code which mainly indicates the assay code, kit lot number and expiration date. The last well is a cuvette in which the fluorometric reading is performed. The wells in the center of the strip contain the various reagents required for the assay.

### Description of the QCV strip:

Wells	Reagents
1	<b>Empty well: no sample is needed for this test.</b>
2	400 µL of 4-Methyl-umbelliferone solution (100 µmol/L) in CHES buffer (20 mmol/L, pH 9.6) + 0.9 g/L sodium azide.
3	Empty well.
4	400 µL of 4-Methyl-umbelliferone solution (0.6 µmol/L) in CHES buffer (20 mmol/L, pH 9.6) + 0.9 g/L sodium azide.
5 - 6	600 µL of 4-Methyl-umbelliferone solution (0.6 µmol/L) in CHES buffer (20 mmol/L, pH 9.6) + 0.9 g/L sodium azide.
7	Empty well.
8 - 9	400 µL of 4-Methyl-umbelliferone solution (26 µmol/L) in CHES buffer (20 mmol/L, pH 9.6) + 0.9 g/L sodium azide.
10	Optical cuvette with substrate: 300 µL of 4-Methyl-umbelliferone (0.6 µmol/L) in CHES buffer (20 mmol/L, pH 9.6) + 0.9 g/L sodium azide.

## WARNINGS AND PRECAUTIONS

1. **Caution: the strips are photosensitive and must be kept in their box to be protected from light.**
2. Do not use reagents after the expiration date indicated on the label.
3. Kit reagents contain sodium azide which can react with lead or copper plumbing to form explosive metal azides. If any liquid containing sodium azide is disposed of in the plumbing system, drains should be flushed with water to avoid build-up.

4. It is recommended to wear gloves as if handling biological tests.

## STORAGE

- Store the VIDAS QCV kit at 2 - 8°C. Do not freeze reagents. **Store all unused reagents in their box at 2-8°C.**
- If stored according to the recommended conditions, all components are stable until the expiration date indicated on the label. Do not use the reagents after the expiration date.

## INSTRUCTIONS FOR USE

Your instrument should be checked **at least once a month** using the VIDAS QCV test, and each time abnormal operation of the pipette mechanism or optical system is suspected.

**It is essential to perform the VIDAS QCV test on all the positions to ensure that the whole instrument has been checked.**

It is essential to archive patient work lists and their results on a monthly basis in order to track results obtained between two QCV tests.

1. Remove necessary components from the kit and return all unused components to storage at 2-8°C.

2. Create the work list: type QCV and indicate the number of tests to be run.
3. Load the reagent strips and SPRs into the instrument (**Note:** no desiccant is needed in the SPR pouch). Be sure to close all SPR® doors and tray covers.
4. Run the test (see User's Manual). All steps will be performed automatically by the instrument. The test will be completed in approximately 20 minutes.
5. **After the test is completed, print the results and analyze them as instructed later in this package insert. Once interpretation of results has been validated, dispose of the used SPRs and strips into an appropriate recipient.**

## RESULTS AND INTERPRETATION

The first reading is performed at the beginning of the assay to check that the reagent is intact:

reading 1 for VIDAS  
R1 for mini VIDAS

The second reading is taken after successive aspirations which are performed at different speeds; they are used to check the ability of the pump to pipette properly:

reading 2 for VIDAS  
R2 for mini VIDAS

The third reading is taken after a small quantity of highly concentrated substrate is pipetted. This operation is performed to check that the optical system is capable of measuring high fluorescence levels:

reading 3 for VIDAS  
R3 for mini VIDAS


When the test is completed, the computer automatically calculates Test Value 1 and 2 (**TV1** and **TV2** for mini VIDAS) using the 3 fluorescence readings. The calculations appear on a report.

VIDAS:      Test Value (TV1) =  $\frac{\text{reading 2}}{\text{reading 1}}$       Test Value (TV2) =  $\frac{\text{reading 3}}{\text{reading 1}}$

mini VIDAS:      TV1 =  $\frac{R2}{R1}$       TV2 =  $\frac{R3}{R1}$

The values for reading 3 (R3 for mini VIDAS) and Test Value 1 (TV1 for mini VIDAS) must be within the acceptable ranges indicated below (see diagram):


- **Acceptable range for Test Value 1 (TV1 for mini VIDAS):** the Test Value 1 for each position must be **≥ the value indicated on the kit label**. If the result of a particular position is outside the range, two new VIDAS QCV tests must be run successively in all the positions in the section concerned (for VIDAS and mini VIDAS). If at least one other non-compliant result is produced in the same section, independent of the position, take the section offline and contact bioMérieux.
- **Acceptable range for reading 3 (R3 for mini VIDAS):** reading 3 for each position must be **≥ 4100 RFU**. If the result of a particular position is outside the range, take the instrument offline and contact bioMérieux.

**LIMITATIONS OF THE METHOD**

- VIDAS QCV may produce an invalid result when the reading 1 or R1 value is < 200 RFU. In this case, perform the test again using a new QCV strip from a new lot (if possible).
- VIDAS QCV should only be used to diagnose problems with the instrument pipette mechanism or optical system. It is not a calibrator.
- VIDAS QCV does not diagnose temperature problems or leakages.
- VIDAS QCV should only be used to test the instrument.

**INDEX OF SYMBOLS**

Symbol	Meaning
	Catalogue number
	In Vitro Diagnostic Medical Device
	Manufacturer
	Temperature limitation
	Use by
	Batch code
	Consult Instructions for Use
	Contains sufficient for <n> tests
	Protect from light

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