【Name】Urinalysis control

【Volume】8ml/bottle 4ml/bottle

【Intended Use】

Urinalysis control(referred to as control) is used for quality control of urinalysis strips and analyzers. It can conduct quality control for 13 items of strips, such as Glucose, Bilirubin, Ketone, Specific Gravity, Blood, pH, Protein, Urobilinogen, Nitrite, Leucocytes, Microalbumin, Creatinine and Ca.

【Test Principle】

Chemistry reagent contained in the Control reacts with the ingredient within the urinalysis strips, thus, strips color changed.

【Main Composition】

Positive: phosphate buffer 0.2%w/w, glucose 1.0%w/w, sodium chloride 0.5%w/w, hemoglobin 0.1%w/w, albumin 0.7%w/w, ethyl acetoacetate 0.9%w/w, sodium nitrite 0.3%w/w, esterase 0.5%w/w, urea 2%w/w, creatinine 0.2%w/w, Ca 0.1%w/w, bilirubin substitute 0.1%w/w, urobilinogen substitute 0.1%w/w, other non-reactive substances and stabilizers 93.3%w/w.

Negative: urea 2%w/w, sodium chloride 0.5%w/w, phosphate buffer 0.2%w/w, other non-reactive substances and stabilizers 97.3%w/w.

【Storage and Validity】

Its validity period is 12 months when stored at 2°C ~ 8°C, sealed and protected from sunlight. Its validity period is 1 month when stored at 2°C ~ 8°C, sealed and protected from sunlight after the first use.

【Applied Analyzer】Urine Analyzer

【Test Method】

1. Place the Control at room temperature(18°C-25°C) for certain time until it returns to normal temperature. And mix the Control evenly through bottom up and down slightly.
2. Using method on semi-automatic urine analyzer: Adjust the analyzer to QC status;

【Symbols】

drip the Control on the strip, the strip side with control should be placed upward; remove the extra control with absorbent paper, the absorbent paper can not touch the strip pads(avoid cross-contamination); place the strip on the waiting place of analyzer. 3. Using method on automatic urine analyzer: Adjust the analyzer to QC status; pour the Control into a clean tube; place the tube into corresponding position of the analyzer.

【Result Explanation】

All QC results within the target range indicates qualified, otherwise, the result is unqualified.

【Limits】

Bilirubin and Urobilinogen used for the Control are replaced by chemical substances. Therefore, there is a trace difference in color between the result of control and the result of urine.

【Matters Need Attention】

1. The Control follows outside of the control bottle can not be collected back to avoid cross-contamination.
2. Control-skin contact should be avoided. Rinse the contact part with plenty of water if accidental skin contact occur.
3. Tighten the cap immediately after use, and store it at 2°C-8°C.


【Reference】

Refer to the supplied control target list.

The supplied reference is obtained through repeat tests. For the same lab, the results might be different when testing on different dates. The difference might be caused by environment, instrument, reagent or test method difference. Even so, the test results should be within the supplied reference range.

【Documents】ZL01 1 28009.3 Urinalysis Control and Preparation

【Instruction Approved and Modified Date】07/2013

Item	Dirui A10	Bayer 10SG	Dirui E10	Roche Combur10	Dirui M10	Kyoto 10EA	Dirui H series	Dirui H-800
Glucose	1+~3+	1+~3+	3+~4+	3+~4+	±~2+	2+~4+	2+~4+	2+~4+
Bilirubin	2+~3+	2+~3+	1+~3+	1+~3+	2+~4+	2+~4+	1+~3+	1+~3+
Ketone	1+~3+	1+~3+	2+~4+	2+~4+	2+~4+	2+~4+	1+~3+	1+~3+
Blood	1+~3+	1+~3+	3+~5+	3+~5+	1+~3+	1+~3+	1+~3+	1+~3+
Protein	1+~3+	1+~3+	2+~4+	2+~4+	2+~4+	2+~4+	2+~3+	2+~3+
Urobilinogen	1+~3+	1+~3+	3+~4+	3+~4+	2+~4+	1+~3+	1+~3+	1+~3+
Leukocytes	±~2+	±~2+	1+~3+	1+~3+	1+~3+	1+~3+	±~2+	1+~3+
Specific Gravity	1.020~1.030	1.020~1.030	1.010~1.025	1.010~1.025	1.010~1.025	1.005~1.020	1.015~1.030	1.015~1.030
pH	6.0~8.0	6.0~8.0	5.5~7.5	5.5~7.5	5.5~7.5	5.5~7.5	5.5~7.5	6.0~8.0
Microalbumin* <sup>1</sup>	/	/	/	/	/	/	≥0.15	≥0.15
Microalbumin* <sup>2</sup>	/	/	/	/	/	/	80~150	80~150
Creatinine	/	/	/	/	/	/	8.8~26.5mmol/L	8.8~26.5mmol/L
Ca	/	/	/	/	/	/	5.0~10.0mmol/L	5.0~10.0mmol/L
Nitrite	Positive							

Note: Microalbumin\*<sup>1</sup> is the microalbumin item of H11-MA, H11-MA(N), H12-MA, H11-800MA and H12-800MA reagent strips.  
Microalbumin\*<sup>2</sup> is the microalbumin item of H2-Cr, H12-Cr, H13-Cr, H13-Ca, H14-Ca reagent strips.

DIRUI Industrial Co., Ltd.  
95 Yunhe Street, New & High Tech.  
Development Zone  
Changchun, Jilin 130012 P.R.China

Tel :+86(431)85100409  
Fax:+86(431)85172581  
E-mail: dirui@dirui.com.cn  
Http://www.dirui.com.cn

EC REP  
Emergo Europe  
Molenstraat 15  
2513 BH The Hague  
The Netherlands

09/2014

**【Name】** Urinalysis control

**【Volume】** 8ml/bottle- 4ml/bottle

**【Intended Use】**

Urinalysis control(referred to as control) is used for quality control of urinalysis strips and analyzers. It can conduct quality control for 13 items of strips, such as Glucose, Bilirubin, Ketone, Specific Gravity, Blood, pH, Protein, Urobilinogen, Nitrite, Leucocytes, Microalbumin, Creatinine and Ca.

**【Test Principle】**

Chemistry reagent contained in the Control reacts with the ingredient within the urinalysis strips, thus, strips color changed.

**【Main Composition】**

Positive: phosphate buffer 0.2%w/w, glucose 1.0%w/w, sodium chloride 0.5%w/w, hemoglobin 0.1%w/w, albumin 0.7%w/w, ethyl acetoacetate 0.9%w/w, sodium nitrite 0.3%w/w, esterase 0.5%w/w, urea 2%w/w, creatinine 0.2%w/w, Ca 0.1%w/w, bilirubin substitute 0.1%w/w, urobilinogen substitute 0.1%w/w, other non-reactive substances and stabilizers 93.3%w/w.

Negative: urea 2%w/w, sodium chloride 0.5%w/w, phosphate buffer 0.2%w/w, other non-reactive substances and stabilizers 97.3%w/w.

**【Storage and Validity】**

Its validity period is 12 months when stored at 2°C~ 8°C, sealed and protected from sunlight. Its validity period is 1 month when stored at 2°C~ 8°C, sealed and protected from sunlight after the first use.

**【Applied Analyzer】** Urine Analyzer

**【Test Method】**

1. Place the Control at room temperature(18°C-25°C) for certain time until it returns to normal temperature. And mix the Control evenly through bottom up and down slightly.
2. Using method on semi-automatic urine analyzer; Adjust the analyzer to QC status;

**【Symbols】**

drip the Control on the strip, the strip side with control should be placed upward; remove the extra control with absorbent paper, the absorbent paper can not touch the strip pads(avoid cross-contamination); place the strip on the waiting place of analyzer.

**【Result Explanation】**

All QC results within the target range indicates qualified, otherwise, the result is unqualified.

**【Limits】**

Bilirubin and Urobilinogen used for the Control are replaced by chemical substances. Therefore, there is a trace difference in color between the result of control and the result of urine.

**【Matters Need Attention】**

1. The Control follows outside of the control bottle can not be collected back to avoid cross-contamination.
2. Control-skin contact should be avoided. Rinse the contact part with plenty of water if accidental skin contact occur.
3. Tighten the cap immediately after use, and store it at 2°C-8°C.

**【Reference】**

Refer to the supplied control target list.

The supplied reference is obtained through repeat tests. For the same lab, the results might be different when testing on different dates. The difference might be caused by environment, instrument, reagent or test method difference. Even so, the test results should be within the supplied reference range.

**【Documents】** ZL01 1 28009.3 Urinalysis Control and Preparation

**【Instruction Approved and Modified Date】** 11/2014

 Store at 
  Please read package insert 
  Use by 
  LOT Batch code 
  IVD In Vitro Diagnostic Use 
  Manufactured by 
  EC REP Authorised Representative 
  CE European In Vitro Diagnostic Medical Device Directive 98/79/EC(IVDD) 
  REF Catalogue number

Item	Roche	Dirui	Kyoto	Dirui	Dirui	Dirui	Bayer	Dirui
	Combur10	E10	10EA	M10	H series	H-800	10SG	A10
Glucose	Negative	Negative	Negative	Negative	Negative	Negative	Negative	Negative
Bilirubin	Negative	Negative	Negative	Negative	Negative	Negative	Negative	Negative
Ketone	Negative	Negative	Negative	Negative	Negative	Negative	Negative	Negative
Blood	Negative	Negative	Negative	Negative	Negative	Negative	Negative	Negative
Protein	Negative	Negative	Negative	Negative	Negative	Negative	Negative	Negative
Urobilinogen	Normal	Normal	Normal	Normal	3.4~17μmol/L	3.4~17μmol/L	3.2~16μmol/L	3.2~16μmol/L
Leucocytes	Negative	Negative	Negative	Negative	Negative	Negative	Negative	Negative
Specific Gravity	1.000~1.015	1.010~1.020	1.005~1.015	1.010~1.020	1.015~1.025	1.015~1.025	1.015~1.025	1.015~1.025
pH	6.0~7.0	6.0~7.0	6.0~7.0	6.0~7.0	5.5~7.5	5.5~7.5	6.0~7.5	6.0~7.0
Microalbumin*1	/	/	/	/	Negative	Negative	/	/
Microalbumin*2	/	/	/	/	10~30	10~30	/	/
Creatinine	/	/	/	/	0.9~4.4mmol/L	0.9~4.4mmol/L	/	/
Ca	/	/	/	/	1.0~2.5mmol/L	1.0~2.5mmol/L	/	/
Nitrite	Negative							

Note: Microalbumin\* 1 is the microalbumin item of H11-MA, H11-MA(N), H12-MA, H11-800MA and H12-800MA reagent strips.  
Microalbumin\* 2 is the microalbumin item of H2-Cr, H12- Cr, H13- Cr, H13- Ca and H14- Ca reagent strips.

 DIRUI Industrial Co.,Ltd.  
 95 Yunhe Street, New & High Tech.  
 Development Zone  
 Changchun, Jilin 130012 P.R.China

Tel :+86(431)85100409  
 Fax:+86(431)85172581  
 E-mail: dirui@dirui.com.cn  
 Http://www.dirui.com.cn

 EC REP  
 Emergo Europe  
 Molenstraat 15  
 2513 BH The Hague  
 The Netherlands

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Emergo Europe  
Molenstraat 15  
2513 BH The Hague  
The Netherlands

# H-Cr Series

## Urinalysis Reagent Strips User's Guide

Rev:09/2009

**Product Name:** Urinalysis Strips  
**Specification:** 100strips/vial, 50 strips/vial, 25 strips/vial

Products Type	Test Item
H13-Cr	Urobilinogen, Bilirubin, Ketone (acetoacetic acid), Creatinine, Blood, Protein, Microalbumin, Nitrite, Leukocytes, Glucose, Specific Gravity, pH and Ascorbic Acid
H12-Cr	Urobilinogen, Bilirubin, Ketone (acetoacetic acid), Creatinine, Blood, Protein, Microalbumin, Nitrite, Leukocytes, Glucose, Specific Gravity, pH
H2-Cr	Microalbumin and Creatinine

The detection of albumin possesses a significance for a variety of diseases in early detection.

- (1) Value for hypertension. Hypertensive patients' urinary albumin excretion rate is significantly higher than normal urine, elevated urinary albumin is an important indicator of cardiovascular disease progression.
- (2) Micro-albumin can predict the development and early diagnosis of diabetic nephropathy. The earliest clinical signal of diabetic nephropathy is the presence of micro-proteinuria, therefore, measurement of microalbumin will help patients with diabetes to take measures in advance to protect renal function.
- (3) Microalbumin detection is sensitive indicator of diabetic microvascular disease.

Products Type	Applicable Instruments
H13-Cr	DIRUI HC-100/300/500/800, NC -100/300/500/800 Urine Analyzer
H12-Cr	DIRUI HC-100/300/500/800, NC -100/300/500/800 Urine Analyzer
H2-Cr	DIRUI HC-100/300/500/800, NC -100/300/500/800 Urine Analyzer


### Collecting and Preparing Specimen

Collect fresh urine in a clean and dry container. Don't centrifuge the urine. Mix the sample well before taking the test. The urine test must be taken within an hour. All specimens must always be taken and kept under sanitary conditions.

**Note:** Water should not be used as negative control liquid. The preservatives will not prevent the deterioration of ketones, bilirubin or urobilinogen. The growth of bacteria in the long-term storage specimen may affect the test results on glucose, pH, nitrite and blood.

### Visual Reading Technique

1. promptly replace cap after taking out strip
2. Immerse the reagent area of the strip in the urine specimen and remove quickly.
3. Run the edge of the strip against the rim of the container to remove the excess urine.
4. Hold the strip horizontally and compare the result on the strip with the colour chart on the bottle label. Make note of the result. For a semi-quantitative result, read the result according to the time specified on the colour chart. pH and Protein can be read at any time within 60 seconds after dipping. For a qualitative results, the strip should be read between 1-2 minutes after dipping. If a positive result is obtained, repeat the test and compare with the colour chart at the time specified. colour changes beyond 2 minutes are of no diagnostic value.

### Instrumental Reading Technique

Follow the directions given in appropriate instrument-operating manual.

### STORAGE

The strips must be kept in the original bottle. Never use the products after the expiration date. Every strip can be used only once. Do not remove the desiccant(s). Don't take the strips from the bottle if it is for immediate use. Cap the bottle immediately and tightly after taking out the strips. The strips should be stored in dry place at the temperature between 2°C-30°C. Don't store the strips in refrigerator and keep them away from direct sunlight. Do not touch the reagent area of the strip. **PROTECTION AGAINST AMBIENT MOISTURE, LIGHT AND HEAT IS ESSENTIAL TO GUARD AGAINST ALTERED REAGENT REACTIVITY.** Deterioration may result in discoloration or darkening of the reagent area of the strip. If all these happen, or the test results are questionable or inconsistent with the expected results, check and make sure the strips are within the expiration date and also compare with the control urine. Please dispose the used strips as wastes according to Treatment Regulations Of Lab Biohazard Materials.

### Reaction Principles

**Glucose:** The glucose oxidized by glucose oxidase catalyzes the formation of gluconic acid and peroxide hydrogen. Peroxide hydrogen releases neo-ecotyles oxide [O] under the function of peroxidase. [O] oxidizes iodide potassium, which makes the colour changes.

**Bilirubin:** The direct bilirubin and dichlorobenzene diazonium produce azo dyes in a strongly acid medium.

**Ketone:** The acetoacetate and sodium nitroprusside cause reaction in alkaline medium, which produces violet colour.

**Specific Gravity:** Electrolyte (M<sup>+</sup>X<sup>-</sup>) in the form of salt in urine reacts with poly methyl vinyl ether and maleic acid (-COOH) which are weak acid ionic exchanger. The reaction produces hydrogenous ionogen, which reacts with pH indicator that causes the colour change.

**Blood:** Hemoglobin acts as peroxidase. It can cause peroxidase release neo-ecotyles oxide(O<sub>2</sub>). (O<sub>2</sub>) oxidizes the indicator and make the colour change subsequently.

**pH:** The method of pH indicator is applied.

**Protein:** This is based on the protein-error-of-indicator principle. Anion in the specific pH indicator attracted by cation on protein molecule makes the indicator further ionized, which changes its colour.

**Urobilinogen:** Urobilinogen and diazonium produce pink azo dyes under the function of strong acid medium.

**Nitrite:** Nitrite in the urine and aromatic amino sulphanilamide are diazotized to form a diazonium compound. The diazonium compound reacting with tetrahydro benzo(h) quinolin 3-phenol causes the colour change.

**Leukocytes:** Granulocyte leukocytes in urine contains esterases that catalyze the hydrolysis of the pyrrole amino acid ester to liberate 3-hydroxy 5-pheny pyrrole. This pyrrole reacting with diazonium forms a purple colour.

**Microalbumin:** sulfone phthalene dye have high sensitivity to microalbumin by the method of protein error.

**Ascorbic Acid:** Ascorbic acid, with 1,2-dihydroxy alkenes, under the alkaline condition, decolorizes the blue 2,6-dichloroindophenolate into colourless N-(p-pheno)-2,6- dichloro-P-amine phenol.

**Creatinine:** creatinine and 3,5-dinitro-benzoic acid react in strong alkaline condition, and generate colored complex.

### Points for attention

#### Glucose

The test is for specificity of glucose. There is no false positive result occurred in reagent strip, caused by any substance in urine. When the ascorbic acid concentration ≤ 0.28mmol/L, 2.8mmol/L glucose may change the color of reagent strip; When the ascorbic acid concentration ≥ 2.8mmol/L or acetoacetic acid concentration ≥ 1.0mmol/L, the sample of glucose concentration is 3~7mmol/L, may occur false negative result.

#### Bilirubin

Normally, even the most sensitive method can't detect bilirubin in urine. It is abnormal to have little bilirubin in urine, which requires further inspection. Medicines that dyes urine red and anything that shows red itself in an acid medium e.g., phenazopyridine may affect the test result. High concentration of the ascorbic acid may cause false negative result.

#### Ketone

The reagent strip reacts with acetoacetic acid and acetone in urine. It doesn't do with β-hydroxy butyric acid. Normal urine specimens usually conduct negative results in the test. False positive results may occur in highly pigmented urine or those containing a large amount of levodopa metabolites.

#### Specific Gravity

The reagent strip for Specific Gravity allows the urine specimens specific gravity between 1.000 and 1.030. In general, the mean error between the results of the strip test and those from the refractive index method is only 0.005. To make it more accurate, 0.005 may be added to readings from urines with pH equal or greater than 6.5. Urine reading instrument can automatically make these adjustments in strip-readings. The urine nonionic constituents such as glucose or radiopaque dye won't make any changes in the test. Highly buffered alkaline urines may cause the low readings comparing with the other methods. Elevated specific gravity readings may occur in the presence of moderate quantities of protein(1.75g/L).

#### Blood

'Trace' reaction may vary among the patients. Clinical judgments are required for individual cases. The presence of green spots (intact erythrocytes) or green color (haemoglobin/myoglobin) on the reagent area within 60 seconds indicates for further diagnostic check. Blood is often found in the urine of the menstruating females. Haemoglobin 150 μg/L-520 μg/L is approximately equivalent to 5-15 cells/μL intact erythrocytes. The reagent strip is highly sensitive to haemoglobin and thus can be used as a supplementary to the microscopic examination. The sensitivity of the strip might be reduced in urine with a large amount of specific gravity. The strips are equally sensitive to myoglobin as to haemoglobin. Certain oxidizing contaminants, such as hypochlorite, may lead to false positive results. Microbial peroxidase associated with urinary tract infection may also produce a false positive result. Ascorbic acid less than 5.0mmol/L in urine may not influence the result of the test.

#### pH

The strip tests for pH values are generally in the range of 5.0-8.5 visually and 5.0-9.0 instrumentally.

#### Protein and microalbumin

Adult normal urinary albumin concentration is less than 20mg / L. The results within 20 ~ 200mg / L indicate clinical micro-albuminuria; the results greater than 200mg / L indicate clinical proteinuria. When albumin excretion rate within 30 ~ 300mg/24 hours, it is shown as the precursor of micro-albuminuria. Urinary albumin excretion may be increased due to a temporary exercise, urinary tract infection, acute febrile illness.

The Microalbumin test can accurately and specifically determine albumin. A strip result of 15 mg/dL is indicative of clinical albuminuria. The test is not affected by other proteins at concentrations at least nine times greater than the excretion rate considered to be abnormal.

Albumin reagent strip is used for testing urinary volume albumin, and reagent strip has low sensitivity to mucinogen and globulin that it can be tested when more than 60mg/dL. Visible haematuria (≥5mg/dL) may also occur false positive result.

The following possibilities occurs when micro albumin emerges in urine:

(1) Occasional urinary micro-albumin may symbolize physiological proteinuria, such as functional proteinuria, or orthostatic proteinuria caused by diet, exercise, mental stress

(2) The clinical significance is very important if the consecutive urine microalbumin, or micro-albumin and glucose, or micro-albumin and blood emerge in urine at the same time.

The Protein test pad is not specific for a particular protein, and proteins other than albumin can cause a positive response. The test is less sensitive to mucoproteins and globulins, which are generally detected at levels of 60 mg/dL or higher.

A visibly bloody urine (> 5 mg/dL) may cause falsely elevated results.

#### Urobilinogen

The reagent strips can detect urobilinogen in low amount as 3 μmol/L (approximately 0.2 Ehrlich unit/dL) in urine. A result of 33 μmol/L in urine indicate the critical value, representing the transition from normal to abnormal, which requires further check on patients and specimens. The negative results are not final to determine the absence of urobilinogen.

#### Nitrite

Gram-negative bacteria in urine converts nitrate[ derived from foods] into nitrite. The reagent strip is essential to nitrite and won't react with the other substances in urine. Pink spots or edges on the strip should not be interpreted as positive result, but any degrees of uniform pink color development should be taken as positive result. The degrees of colour development the numbers of bacteria are not in direct proportion. The negative result doesn't mean the existence of bacteria in large amount. Negative result may occur (1) when urine doesn't contain organism that caused the conversion from nitrate to nitrite, (2) when urine has not remained in the bladder long enough (four hours up to) let the nitrate covert into nitrite, (3) the nitrate in the foods is absent. Large High volume of specific gravity in urine may reduce the sensitivity of the test. 2.8mmol/L ascorbic acid or less won't interfere the test result.

#### Leukocytes

Test area react with esterase in leucocytes (granulocytic leukocytes). Normal urine specimens generally yield negative result; positive results (+ or greater) are clinically significant. Individually observed 'Trace' results may be of questionable clinical significance; however 'Trace' results observed repeatedly may be clinically significant. 'Positive' results may occasionally be found with random specimens from females due to contamination of the specimen by vaginal discharge. Elevated glucose concentrations (160mmol/L) or high specific gravity may cause decreased test results.

#### Ascorbic Acid

The test area can detect the ascorbic acid in urine. Through the ascorbic acid detection, we will know the level of ascorbic acid in the body and the effect degree that the ascorbic acid bring to the test on glucose, bilirubin, blood and nitrite. It will reduce the sensitivity when the oxidant (such as potassium permanganate, hypochlorite) in the urine.

#### Creatinine

Normal adult urinary creatinine concentration: 0.6 ~ 2.0g/24 hours (strip test result: 4.4 ~ 17.7mmol / L). Random urine creatinine test results vary widely, which may be in the range of from 0.9mmol / L to 26.5mmol / L. High level happens to concentrated urine and morning urine (strip test results may be higher than 17.7mmol / L). Due to excessive drinking and other polyuria to dilute urine, which will lead to a typical low concentration of urine (test results ≤ 50mg / dL) Visible haematuria (≥ 50mg / dL) or cimetidine may cause creatinine false test results elevated.

#### Aulbumin--creatinine rate

Adult normal urine albumin level is less than 30mg albumin / g creatinine (3.4mg albumin / mmol creatinine), sample rate result is in the range of 30 ~ 300mg / g (3.4 ~ 33.9mg/mmol) can be defined as microalbuminuria, and proteinuria when > 300mg / g (> 33.9mg/mmol)

### Result table

Conventional units and international units available

Test item	Abb.	Print or display result			
		Conventional Unit		International Unit	
Microalbumin	MA	10 mg/L	80 mg/L	10 mg/L	80 mg/L
		30 mg/L	150 mg/L	30 mg/L	150 mg/L
Creatinine	Cr	10 mg/dL	200 mg/dL	0.9 mmol/L	17.7 mmol/L
		50 mg/dL	300 mg/dL	4.4 mmol/L	26.5 mmol/L
		100 mg/dL		8.8 mmol/L	
Aulbumin-creatinine rate	A/C	<30 mg/g (Normal)		<3.4 mg/mmol (Normal)	
		30-300 mg/g (Abnormal)		3.4-33.9 mg/mmol (Abnormal)	
		>300 mg/g (Highly abnormal)		>33.9 mg/mmol (Highly abnormal)	

**Remark:** The shaded parts are abnormal results.

### Limitation of Procedures

Like all the other laboratory tests, definitive diagnostic or therapeutic decisions should not be made or based on any single result or method. Specific performance characteristics are based on clinical and analytical studies. In clinical specimens, the sensitivity depends upon several factors; the variability of colour perception, the presence or absence of inhibitory factors typically found in urine, specific gravity, pH, and the lighting conditions when the product is read visually. Each colour block or instrumental display value represents a range of values. Because of specimen and reading variability, specimens with analyte concentrations that fall between nominal levels may give results at either level. Results at levels greater than the second positive level for the Protein, Glucose, Ketone, and Urobilinogen tests will usually be within one level of the true concentration. Exact agreement between visual results and instrumental results might not be found because of the inherent differences between the perception of the human eye and the optical system of the instruments.

### Note

1. In vitro diagnosis for professional only.
2. Reagent strips must be stored in original cylinder; do not take out strips except use them immediately; removing the desiccant is forbidden; tighten the cap immediately after take out strips;
3. The expired should not be used; Please make sure whether the strip is over the validity, and test it by quality control if the deterioration of strips occurs, which leads color shallower or deeper, or any question occurs with the unexpected results.
4. Water can not be used as negative control.
5. Oe time use only for each strip.
6. Please read the instruction carefully before the use.
7. Please refer to the lab biological treatment of harmful materials to dispose the waste strips.
8. Do not place them in refrigerator, however they should be stored in dry circumstance; avoid direct sunlight; do not touch the reactive area of strip; ambient moisture, light and heat must be prevented in order to protect the reactivity of the strips.

### SENSITIVITY AND TEST RANGE OF H-CR SERIES URINALYSIS STRIPS

Item	Sensitivity	Instrumental test range	Visual test range
Glucose(mmol/L)	2.8-5.6		Neg.-56
Protein(g/L)	0.15-0.3	Neg.- 3.0	Neg.-20.0
Microalbumin(mg/L)	30-80	10-150	10-150
Ketone(acetoacetic acid)(mmol/L)	0.5-1.0	Neg.-7.8	Neg.-16
Blood (Ery / μ L)	5-15		Neg.-200
Bilirubin ( μ mol/L)	8.6-17		Neg.-103
Nitrit ( μ mol/L)	13-22		Neg.-Pos.
Leukocytes (leuko / μ L)	5-15		Neg.-500
Urobilinogen ( μ mol/L)	3.4-17		3.4-135
Ascorbic acid(mmol/L)	0.6-1.4		0-5.7
Creatinine(mmol/L)	—	5.0-9.0	0.9-26.5
PH			
Specific Gravity		1.005-1.030	1.000-1.030

### REACTIVE INGREDIENTS (based on dry weight at time of impregnation)

Microalbumin	sulfone phthalein buffer nonreactive ingredients	2.2% w/w 96.0% w/w 1.8% w/w
Protein	tetrabromophenol blue buffer nonreactive ingredients	0.1% w/w 97.4% w/w 2.5% w/w
Blood	diisopropyl/benzene dihydro peroxide tetramethylbenzidine buffer nonreactive ingredients	26.0% w/w 1.5% w/w 35.3% w/w 37.2% w/w
Glucose	glucose oxidase (microbial 123U) peroxidase(horse radish, 203U) potassium iodide buffer nonreactive ingredients	1.7% w/w 0.2% w/w 0.1% w/w 71.8% w/w 26.2% w/w
Ketone	sodium nitroprusside nonreactive ingredients	5.7% w/w 29.9% w/w 64.4% w/w
Leukocytes	pyrrole amino acid ester diazonium salt buffer nonreactive ingredients	4.3% w/w 0.4% w/w 92.6% w/w 2.7% w/w
Nitrite	p-arsamiliacid-N-(1-Naphthol)-ethylendiamine tetrahydroquinoline buffer nonreactive ingredients	1.3% w/w 0.9% w/w 89.1% w/w 8.2% w/w
Specific Gravity	bromthymol blue poly(methyl vinyl ether co maleic anhydride) sodium hydroxide	4.8% w/w 90.2% w/w 5.0% w/w
pH	bromocresol green bromthymol blue nonreactive ingredients	3.3% w/w 55.0% w/w 41.7% w/w
Bilirubin	2,4-dichloroaniline diazonium salt buffer nonreactive ingredients	0.6% w/w 57.3% w/w 42.1% w/w
Urobilinogen	fast B blue buffer nonreactive ingredients	0.2% w/w 98.0% w/w 1.8% w/w
Ascorbic acid	2,6-dichloroindophenolate hydrate buffer nonreactive ingredients	0.8% w/w 40.7% w/w 58.5% w/w
Creatinine	3,5-dinitro-benzoic acid buffer nonreactive ingredients	4.8% w/w 85.2% w/w 10.0% w/w

### Notes on symbols and marks

	Batch code		Use by
	Single use		In Vitro Diagnostic Use
	Manufactured by		Store at
	Please read package insert		These test strips conform to the directive 98/79/EC(IVD-directive)
	Authorized Representative		

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