

OPMI® VARIO 700 on the floor stand

Software Release 1.1



User Manual

G-30-1716-en





Version 4.1

2011-09-14



- About This Manual* The user manual is part of the delivery package.
- Please read the user manual carefully before using the device.
 - Keep it at the place of use.
 - Keep the user manual over the entire service life of the device.
 - Pass it on to any subsequent owners or users of the device.
 - Also observe the user manuals for accessories and other system components.
- Orientation aids*
- The chapter overview at the beginning of the user manual provides a summary of all subjects.
 - The contents of each chapter are specified in detail at the beginning of each chapter.
 - A list of abbreviations, key words and technical terms in the annex facilitates the search for specific terms.
- Scope* The present manual applies to the OPMI® VARIO 700 with the following identification:
- Reference number: 6636
 - Software release 1.1
- Trademark* – OPMI® is a registered trademark of Carl Zeiss.
- Information on the manufacturer:
- | | |
|----------------------------|---|
| Carl Zeiss Meditec AG | Tel: + 49 (0) 7364 - 20-0 |
| Goeschwitzer Strasse 51-52 | Fax: + 49 (0) 7364 - 20 4823 |
| 07745 Jena | Email: surgical@meditec.zeiss.com |
| Germany | Internet: www.meditec.zeiss.com |
- Subject to change in design and scope of delivery and as a result of ongoing technical development. Printed in Germany.
- Copyright* © Carl Zeiss Meditec AG 2011
All rights reserved. Reprints and reproductions, including excerpts, are not permissible. However, this manual may be digitized and archived for internal use. Access to this data by third parties is not permissible and must be prevented.

Chapter Overview

	Chapter: Safety Measures	5
	Chapter: System Overview	43
	Chapter: Preparations for use	75
	Chapter: Operation	115
	Chapter: Procedure in the event of faults	171
	Chapter: Care and Maintenance	187
	Chapter: System Data	203
	Chapter: Annex - (Options)	225
	Chapter: Indexes	267

Safety Measures



Key to symbols	7
Hazard symbols	7
Information symbols.....	7
Target audience.....	8
Field of use	8
Intended use	8
Normal use.....	9
Warranty and liability.....	9
Notes for the operator.....	10
Requirements for operation.....	12
Prior to the very first use.....	12
Risk of burn injuries caused by high illumination intensity	13
Explanation of the visual and acoustic signals	16
Safety devices.....	18
Automatic Iris Control	18
Focus-light control.....	20
180W xenon lighting	22
300-W xenon lighting (optional).....	23
"Light only" mode	24
Symbols and labels on the device	26
Labels on the microscope	26
Labels on the suspension system	30
Labels on the 180W illuminator	38

Labels on the 300-W xenon illumination system (optional)..... 40

Key to symbols

We would like to inform you about safety aspects which must be observed when handling this device. This chapter contains a summary of the most important information concerning matters relevant to instrument safety.

Hazard symbols

The following safety information has been incorporated into the user manual. Please note this information and be particularly careful in these cases.



WARNING

Warning label, which may refer to **fatal injuries** or **severe injuries** if precautions are not followed.



CAUTION

Indicating a hazard, which may lead to **moderately severe injuries** if risks are not avoided.

NOTE

Warning label, which may refer to **minor injuries** and **property damages** if precautions are not followed.

Information symbols

The following information symbols are used in this user manual:

- Listing
- ✓ Prerequisite for an action
- Prompt for action
- Result of an action
-  Additional information and tips

Target audience

**CAUTION****Operation by trained personnel only**

This user manual is intended for physicians, nurses and other medical staff who prepare, operate or maintain the system after appropriate training. It is the duty of the institution operating the system to train and instruct all staff using the system.

Installation and service work that is not described may only be performed by Carl Zeiss Service or specialists authorized by CZS.

Field of use

Intended use

The OPMI VARIO 700 is a surgical microscope intended for the lighting and magnification of the surgical area and for the support of visualization in surgical procedures.

Only use the device in accordance with the intended use.

**CAUTION****Injury to the patient's eye!**

The xenon illumination must not be used for ophthalmic procedures. Make sure that no xenon light enters the patient's eyes or the eyes of those persons in the vicinity of the surgical field. Make sure that no xenon light enters the patient's eyes during paranasal sinus surgery. Cover the eyes of the patient.

**CAUTION****Not for diagnosis purposes**

The live image, video sequences and individual images may not be used for diagnostics purposes. Visualized images may contain variations in shape, contrast and color.

NOTE**Have a suitable replacement system ready at hand!**

Like any technical device used in the OR, this system may fail. We would therefore recommend having a suitable surgical microscope at hand as a replacement during surgery.

Normal use

The OPMI VARIO 700 is a surgical microscope intended for surgical operations and enables the user to visually magnify and illuminate the surgical area.

The system is intended for use in clinics and other human medicine institutions.

Warranty and liability

Warranty and liability depend on the applicable contractual stipulations.

NOTE**Loss of warranty**

No modifications may be made to this system without the manufacturer's approval. If the system is modified, suitable inspection and testing must be completed to ensure that it can still be used safely.

The manufacturer accepts no liability for damage caused to the device as a result of unauthorized modifications. The warranty also becomes null and void.

Notes for the operator

- Observe the legal provisions for accident prevention and occupational health and safety applicable in the relevant country.
- Modifications and repairs of this device or any equipment operated together with this device may only be performed by Carl Zeiss service staff or other persons authorized by Carl Zeiss.
- If required by the regulations and directives applicable in the country of use, connect the system to a special emergency backup power supply.
- Ensure that the installation conditions and the use of the device meet surgical requirements: See ambient conditions page 223.
- Never cover any ventilation openings on the light source and carrier arm as the system switches off in the event of overheating.
- Never attempt to forcefully attach any electrical connections (plugs, bushings). If a plug cannot be connected easily, check again whether the plug and socket are made to fit. If the plug connection is damaged, please call our service department. They will be happy to assist you.
- Do not use multiple sockets!
- Never open the device! The device contains freely accessible live components. If you remove the housing, you run the risk of electric shock.
- Do not touch the system if your body is electrostatically charged and the system is not grounded.
- Please note the information on EMC (electromagnetic compatibility) in the chapter "System data" on page 209.
- To prevent malfunctions, signal lines may not run parallel to power supply lines.
- Switch off the unit at the power switch if you notice any smoke, sparks or unusual noise. Do not use the unit until it has been repaired by our service team.
- Never leave a device unattended with the light source still switched on.
- In order to prevent any impairment of the device's safety due to age, wear, etc., the user must ensure that the device is subjected to the necessary safety checks (see section "Care and maintenance").
- Go through the checklist in the chapter "Operation".
- Never operate the system unattended.

- Always switch off the system
 - before disconnecting it from / connecting it to line power
 - if it will not be used for a prolonged period of time
 - or for cleaning purposes.
- Always use the master switch to turn off the device if it is not in use.
- Insufficient, faulty or wrong methods of cleaning or not disinfecting the device according to this instruction manual may severely increase the risk of infecting the patient and/or medical personnel.

Requirements for operation

Prior to the very first use

Our service representative or an expert authorized by us will install the system. Please make sure that the following requirements continue to be met for further operation:

- ✓ The connecting components have been properly connected. The screw connections have been firmly tightened.
- ✓ All cables and plugs are in perfect condition.
- ✓ The voltage set on the system corresponds to the rated line voltage on the site of installation.
- ✓ The power cord is plugged into a power outlet which has a properly connected protective ground contact.
- ✓ The device is connected with the power cord supplied.

Risk of burn injuries caused by high illumination intensity



CAUTION

If the xenon light is not used as intended, excessive illumination intensity may result in third-degree burns.

General

This system features a powerful xenon light.

- A device with light source switched on may never be left unattended.

The risk of burns is influenced by several different factors.

System-related factors:

- The wavelength range is limited by filters to the visible range between 400 nm and 700 nm. These filters remain stable over a very long period of time and cannot be exchanged by the user.
- With increased aging of the light source, the illumination intensity delivered for a particular setting decreases. When the lamp is replaced, the illumination intensity increases again to the high, original value.
- There are certain accessories (e.g., folding tube $f = 170/260$ mm, or eyepieces with increased magnification) which, when used, may reduce the image brightness at the user's eye. In this case, check whether the user settings are the proper ones for the particular accessories and reduce the light intensity to the minimum required for the application.

Surgery-related factors:

- The selected intensity of the light source is a major factor for the risk of injury. It should always be set to the minimum required for the surgical procedure to be performed.
- The size of the luminous field influences the injury risk in two respects:
 - If the diameter of the luminous field is large, areas of the skin are illuminated which are not monitored as strictly by the surgeon and which are not moistened sufficiently. These areas represent a particular injury risk. Such injuries can be prevented by adjusting the luminous field diameter to the smallest size required for that particular operation.
 - If the luminous field is reduced in size, the intensity increases because the light focus becomes stronger. If the luminous field is reduced in size, the intensity should therefore also be reduced, if possible.

- A long surgical procedure increases the risk of injury, in particular if a standard procedure takes considerably longer than usual.
- Injuries in the peripheral area can be prevented by covering this area with moist, sterile gauze.
The gauze must be moistened at regular intervals to avoid that the area dries out or heats up. The risk is increased if dry drapes are used to cover such areas.
- You should also take into account that some areas of the body may be more sensitive than others.
- Certain preparations of the surgical field, local vasoconstrictive medications and incision drapes may also result in a higher risk of injury (drapes may heat up to varying degrees depending on their color and moisture content).

Patient-related factors:

- The general condition of a patient's health may contribute to the risk of injury.
- The skin type may also play a major role in this respect.
- Certain medications affect the sensitivity to light.
- The interaction of heat and antimicrobial substances in incision foils may increase the patient's reaction to these substances.

Recommendations

Due to the large number of different factors involved and the lack of scientific publications on this topic, Carl Zeiss cannot provide guidance on acceptable intensities and exposure durations. However, this system has several features that can help the user to reduce the risk of burns.

- The initial light intensity can be preset to a low value.
- The automatic Iris Control Function keeps the light intensity in the illuminated area virtually constant, irrespective of the size of the luminous field.
- The Spot function allows reducing the size of the illuminated area to the area of observation during surgery.
- Using the buttons on the hand grip or foot control panel, the surgeon can then set the illumination intensity to the value required for the procedure. Please note that the intensity increases with decreasing luminous field size if the Spot function is used. For this reason, the intensity should only be set after the size of the luminous field has been changed.

- Usually, the magnification factor is increased during surgery. This darkens the image so that the illumination intensity must be increased. This loss in image brightness is automatically compensated if the zoom-dependent brightness control is activated.
- A device with light source switched on may never be left unattended.
- The risk of burns can be reduced by constantly irrigating the illuminated surgical field and by keeping it moist.
- If used, drapes should also be re-moistened at regular intervals in order to prevent heat from accumulating underneath the drape.
- Switch off the light when the microscope is not in use, and ensure it is not directed at unprotected bare skin.

Final remark

Carl Zeiss recommends:

- Reduce the lighting of the surgical field - as far as is possible for the safety of the patient and for a good microscopic image.
The setting of the illumination intensity has been preconfigured at the factory such that an audible signal alerts the user to potential tissue damage caused by excessive light intensity, once the light threshold is exceeded.
- Adjust the light threshold (page 138).

Explanation of the visual and acoustic signals

Device function	Visual signal
Device ready (brakes ready)	Power switch lights permanently green and main menu is shown on the display
"Light only" mode	"Light only" switch lights continually orange
Device off	Power switch does not light up
Indicates that a message is being shown on the display	Warning triangle is shown on the additional display

Situation	Acoustic signal*
For all messages displayed	Sound (ca. 250 ms)
On exceeding the light threshold	Sound (ca. 500 ms)
Speedfocus has focused.	1 sound (ca. 50 ms)
Speedfocus has not focused Focusing not possible.	3 sounds (ca. 3 x 50 ms)

*) the duration of the sound varies but the frequency does not.

Safety devices

Automatic Iris Control

Problem

For a large luminous field diameter, areas of the skin are illuminated which are not monitored as strictly by the surgeon and which are not moistened sufficiently. These areas represent a particular injury risk. These kinds of injuries can be prevented by adjusting the luminous field size of the field of vision (see diagram on the following page).

If the luminous field is reduced in size, the intensity increases because the light is focused more. So, if possible, the intensity should be lowered as soon as the size of the luminous field is reduced.


A long surgical procedure increases the risk of injury, especially if a standard procedure takes considerably longer than usual.


Remedy: Automatic Iris Control

The effect of Automatic Iris Control is that the luminous field size automatically adjusts to the size of the field of vision irrespective of zoom and focus.

Also, the light intensity adjusts to this change, meaning the image of the surgical field the user sees in the eyepiece is more or less uniformly bright (see page 138).

Every user profile has a setting for whether Automatic Iris Control is enabled or disabled.





Focus-light control



The focus-light control is a safety feature intended to prevent possible tissue damage due to excessive illumination intensity.



CAUTION

If xenon illumination is used improperly, excessive illumination intensities may lead to third-degree burns!

The focus-light control regulates the illumination intensity as a function of the working distance. With increasing working distance, the user has accordingly more light available.

The set maximum illumination intensity is shown by a blue bar (1) in the light controller (e.g., max. 24%).

The illumination intensity is reduced to the set limit value (1). It can be set within the range below the blue bar, but not above it.

The focus-light control function can be enabled or disabled by the user as required. It is enabled by default. The factory setting is 25%.

The service technician can set the focus-light control to a value between 5% and 100% in the Service menu.




The service technician can configure the focus-light control in the Service menu such that it cannot be disabled.

For focus-light control, please refer to: Light - Adjusting the light, page 140.


1

Focus-light control enabled

Focus-light control disabled



180W xenon lighting



- The xenon lighting unit has a built-in UV/IR heat protection filter.
 - Display of the light intensity (as a percentage) on the display and additional display.
If the light intensity is greater than or equal to the threshold set, an exclamation mark is displayed after the light intensity display and a warning is sounded on exceeding the threshold.
 - Note on exchanging the lamp
The lamp has a maximum service life of 500 hours.
Once the lamp being used has reached 450 operational hours, the display informs the user on restarting the system that the maximum service life is soon to be reached.
 - Temperature monitoring
The device monitors the temperature. The light source is switched off if too high a temperature is reached.
 - Second lamp
The 180W xenon light source has a lamp module that contains two xenon lamps. The second lamp can be swung into position manually should the first lamp fail.
- 1** "Position second lamp" button
See page 172

Information on the 180W and 300W xenon lighting units (option):



CAUTION

Lamp failure

- The service life of the xenon lamp is limited to 500 hours. If the service life of the lamp is exceeded, the xenon lamp may suddenly fail.
- With increased aging of the light source, the actual illumination intensity at a particular setting decreases.
- Please replace the xenon lamp in good time.
Note the remaining hours display on the screen.


NOTE

Light failure

Never cover any ventilation openings as the light source of the system switches off in the event of overheating.

300-W xenon lighting (optional)

- The 300-W xenon lighting unit has a built-in UV/IR heat protection filter.
- Display of the light intensity (in percent) on the display and additional display.
If the light intensity is greater than or equal to the threshold set, an exclamation mark is displayed after the light intensity display and an audible signal is issued (see page 16).
- Note on lamp change
The lamp has a maximum service life of 500 hours.
Once the lamp in operation has reached 450 operating hours, the display informs the user on restarting the system that the maximum service life will soon be reached.
- Temperature monitoring
The device monitors the temperature. If the temperature is too high, a message is shown on the display and an audible signal is issued. The light source is switched off if the temperature becomes too high.
- Lamp 2 (automatic lamp change, see page 174)
The 300-W xenon lighting unit features two lamps. If the active lamp fails, it is automatically replaced by the other lamp.
However, lamp change can also be triggered via the control panel.
If automatic lamp change fails, the lamp must be replaced manually (see page 176).
A lamp can only be replaced manually (see page 190).



1 Fastener for opening the lamp housing

Opening the lamp housing isolates all accessible parts from the power supply.

"Light only" mode

1 Button for "Light only" mode

If one of the main functions fails (focus, zoom, light control, magnetic brakes) and other functions are affected, press the orange-colored button for "Light only" mode.

The system reacts as follows:

- All electrical components (except the light source) are isolated from the power supply.
- The light source adjusts to a brightness of maximum 75% to permit you to complete the application.
- The magnetic brakes remain engaged and cannot be released.
- Display, additional display and optional TFT screen are switched off.


To be able to continue working:

- manually adjust the magnification on the zoom adjustment knob of the microscope
- manually adjust the working distance on the focus adjustment knob of the microscope
- manually adjust the luminous field diameter
- manually position the stand for focusing and for moving the system in the X-Y directions

Press button (1) once more. The button disengages and the device starts. If you experience problems in rectifying a fault or if faults keep recurring, stop using the system, label it as non-operational and contact Carl Zeiss Service.



To guarantee smooth operation of "Light only" mode, check its function at regular intervals by enabling the mode (without patient).



Symbols and labels on the device



CAUTION

Observe all warnings and information labels!

- If you notice that any label is missing on your device or has become illegible, please contact the Carl Zeiss service department.

Labels on the microscope

F

1 Labeling of the focus button

Z

2 Labeling of the zoom button

Open →



3 Note "Open in direction of the arrow"
on the knurled rings of the left and right image output ports.




4 Symbol for the zoom adjustment knob for manual operation
Identifies the two directions of rotation to increase or reduce the zoom.


**VARIO
700**

5 Labeling of the microscope







6 Autobalance label (option)





7 Manual balancing label: front to back tilt




8 Manual balancing label: lateral tilt



9 Fine adjustment label
Fine adjustment of the microscope's tilt motion



Labels on the suspension system




10 Rating label

The rating label contains the following information:

- Manufacturing symbol 
- Manufacturer (company name) Carl Zeiss Meditec AG
- Manufacturer's address Goeschwitzer Strasse 51-52
07745 Jena, Germany
- Serial number  XXXXXXXXXXXX
- System name OPMI VARIO 700
- Reference number  6636
- Rated voltage (115 V): 100 V-125 V
(230 V): 230 V-240 V
- Current consumption max. 1200 VA
- Line frequency range 50 Hz-60 Hz
- Protection type IP20
- CE mark 



11 NRTL approval



12 SIP label

The SIP label contains the following information:


- Manufacturer (company name)
- Manufacturer's contact data,
- i.e., phone number, fax number and email address of the local contact of the national Carl Zeiss sales organization.
- SIP number
A unique identification number assigned to your system. A product file is maintained by Carl Zeiss for this SIP number.

**13** "Year of manufacture" label

This label indicates the year the system was manufactured.

**14** "Observe disposal regulations" label

Electrical or electronic devices must not be disposed of as normal domestic waste. For more information on the disposal of electrical and electronic devices, please see the chapter "Maintenance and care".







- 15** "Observe instructions for use!"
Observe the user manual or accompanying documents.



- 16** "Risk of crushing" label
Fingers may be crushed. Do not touch this area while moving the surgical microscope or parts of it.



- 17** "Maximum load" label on the microscope body
The maximum load (accessory equipment) on the microscope body may not exceed 7 kg!




- 18** "Moving position" warning label
Indicates the transport position of the system. Before transporting the system, make sure that it is in this position in order to avoid any damage.




- 19** Light Only mode
The OPMI and suspension system functions are deactivated. Illumination continues to operate at a constant intensity.

USB


- 20** USB port (video recording option)
Identifies the connector e.g. for USB sticks or external HDDs.



- 21** "DRAPE" label
The drape is secured here with straps so that it is nearly airtight.



- 22** Label for manual balancing





23 "FCP Gateway WL" label (optional)

This device complies with Part 15 of the FCC Rules and with RSS 210 standard (Canada). Operation is subject to the following conditions:

- The device does not cause harmful interference.
- The device must remain operative even under the influence of unwanted frequencies.




24 Potential equalization

for connection of the system to the potential equalization system.



25 "Connector for foot control panel" label



26 Remote connector

This label identifies a connection on the device with max. 24V/0.5A



max. 500 VA

27 Power outlet socket warning label

Only connect devices with the correct electrical ratings.



28 Defined user group

Identifies the service interface.



max. load
13kg


29 "Maximum load on instrument tray" warning label

Indicates that the maximum load capacity of the instrument tray for accessory equipment is 13 kg.



30 "Videomonitor" label

Identifies the mounted video monitor



NAV

31 Connector for cable for a microscope navigation system (option).



32 "Ethernet connector (LAN)" warning label
Prerequisite: video option

SD connector panel (option)



33 Video signal output/input port (orange)
(e.g. for an SD endoscope camera)



34 Video signal output port: Y/C or S-Video
(e.g. Y/C: video signal output port for an external monitor)



35 Video signal output port: VBS

HD connector panel (option)



36 Video signal output port: YPbPr

37 Video signal output port: DVI-D


HD-SDI

38 Video signal output port: HD-SDI

39 Video signal output port: Y/C or S-Video



40 Video signal output port: VBS



Labels on the 180W illuminator



- 41** Warning label Injury to the patient's eye!
 Xenon illumination must not be used for ophthalmic procedures.
 Make sure that no xenon light enters the patient's eyes.



- 42** Warning label "CAUTION "
 Observe the user manual or accompanying documents.

43 Labels on the inside of the 180W xenon lamp



- CAUTION label, "Hot surfaces" warning label
 Order no. of the lamp module and symbol for the lamp module




- Order no. of the individual lamps
 Lamp symbol: lamp 1 and lamp 2

Xenon-Lamp 180 W





44 Pointing arrow for lamp module change



- 45** Changing the lamps
 The label shows the three steps for the replacement of a lamp.
- A - push button
 - B - pull out lamp module
 - C - replace lamps



- 46** "Observe the user manual" label
 Observe the user manual or accompanying documents.




Labels on the 300-W xenon illumination system (optional)



- 47** Warning label: Injury to the patient's eye!
Xenon illumination may not be used for ophthalmic procedures.
Make sure that no xenon light enters the patient's eyes.




- 48** "CAUTION" label
Observe the user manual or accompanying documents.




- 49** "Lamp change" label
This label marks the cover to be opened for changing the lamp.




- 50** "Hot surfaces" warning label
Caution, risk of burns!
Only change the lamp module after it has cooled down.



- 51** "Lamp change" label
This label marks the cover to be opened for changing the lamp.



- 52** Lamp 1 / Lamp 2





System Overview



Configuration options	44
Overview	46
Options	47
Surgical microscope	48
Handgrips.....	50
Controls for the tubes.....	52
Controls for the widefield eyepieces	53
XY coupling (option)	54
Additional display	55
Components of the floor stand	56
Controls of the xenon lamp housings	62
Controls for manual balancing.....	64
Controls on the control and display panel.....	65
Connector panel on the floor stand	66
Foot control panel (optional)	72
Surgical microscope and laser micromanipulator	74

Configuration options






Overview

Surgical microscope The surgical microscope enables the illumination and visual magnification of the surgical field. For visual observation, the surgical microscope can be equipped with different tubes, eyepieces and coobservation tubes.

XY coupling (option) The optional XY coupling is located on the suspension system of the stand and can be moved horizontally by ± 30 mm in the X and Y directions using the joystick button of the handgrips or the joystick on the foot control panel.

Floor stand The floor stand is used to power and control the microscope functions. It comprises the stand column, stand base, carrier arm, lamp housing, suspension arm and suspension system. All movement axes are provided with electromagnetic spring brakes which can be unlocked and locked by the user via a button on the handgrip. The system can be equipped with various options (see the following page).

Standby position The six axes enable rough positioning of the surgical microscope and moving it into its standby position. In the horizontal position, the suspension arm provides a working range of approx. 1600 mm between the column and the optical axis. The maximum height of the suspension system has been designed in such a way that a user with an eye height of approx. 200 cm can look straight into the eyepiece in an upright, standing position.



The light source in the lamp housing supplies light to the surgical microscope. Depending on the ordered configuration, the system is equipped with a 180 W or 300 W xenon light source.

A handle mounted on the stand column facilitates moving the system. There are two cable brackets on the left and right side of the column that can be used to wind up the cables and to attach the foot control panel. The stand base has four steerable casters, each with a brake tab. The casters must be locked to ensure that the stand is securely positioned. The casters are provided with cable deflectors.

Instrument tray and 17" video monitor The system can optionally be equipped with an instrument tray for Zeiss equipment and a 17" video monitor.

Video The surgical microscope can optionally be equipped with internal or external SD or HD camera systems.

Options

(for ordering data, see page 217)

- AutoBalance
- AutoDrape
- Navigation interface
- XY coupling
- SpeedFocus (requires video option)
- Video focus and Video SpeedFocus
- 300-W Superlux xenon illumination system with automatic lamp change
- Video in (requires 3CCD video standard definition option)
- FCP foot control panel (14 functions), wired
- FCP WL foot control panel (14 functions), wireless via Bluetooth
- 3 CCD HD video (high definition)
- 3 CCD video (standard definition)
- Integrated USB recording
- Ethernet interface
- Instrument support
- Integrated 17" TFT video monitor with carrier arm and tray
- Integrated 17" TFT video monitor with carrier arm

Surgical microscope


- 1 Securing screw
For attaching tubes and accessories (e.g. stereo bridge)
- 2 Handgrip clamp
- 3 Handgrips
- 4 Adjustment knob, zoom
For manual setting of the zoom system
- 5 Adjustment knob, focus
For manual setting of the working distance. The working distance set can be read from the control knob.
- 6 Adjustment knob for luminous field diameter
For manual setting of the luminous field diameter.
- 7 Additional display
- 8 Knob for AutoBalance (option, see page 100 and page 55)
To activate AutoBalance
Manual adjustment not applicable for AutoBalance (9 ,10 and 11).




- 9 Balance setting, microscope tilt axis (manual)
For manually setting the balance of the microscope's swivel motion.
- 10 Balance setting, microscope swivel motion (manual)
For manually setting the balance of the microscope's tilt motion.



- 11 Fine balancing, microscope (manual)
For fine adjustment of the microscope's balance.







Handgrips

- 1 Programmable button A (see page 152)
- 2 Programmable button B (see page 152)
- 3 Joystick
The joysticks on the left and right handgrips are used in the basic settings for motorized fine adjustment in the XY direction (optional).
If the navigation system is connected, the joystick on the right handgrip is used to control a mouse pointer for the functions of the navigation system connected (dependent on system).
- 4 Rocker switch Focus + / Focus - (factory setting)
For continuously variable focusing within the working distance. The current value is shown on the display or additional display.
- 5 Rocker switch Zoom + / Zoom - (factory setting)
To adjust the magnification. The current value is shown on the display or additional display.
- 6 Programmable button C (see page 152)
- 7 Programmable button D (see page 152)
- 8 **Brake release button (AB)**
Brake release button for all stand and microscope axes. Whilst this button remains pressed, all magnetic brakes ("All Brakes") are released and the system can move freely in all directions. On releasing the button, the magnetic brakes lock all axes simultaneously.
- 9 **Brake release button (SB)**
Brake release button for the microscope axes or stand axes (factory setting). Whilst this button remains pressed, only the magnetic brakes of the stand axes or microscope axes are released ("Selected Brakes") and the stand or microscope can move freely. On releasing the button, the magnetic brakes lock the axes simultaneously.




Configuring the handgrip is explained on page 152.



Controls for the tubes



180° tiltable tube

- 1 Eyepiece mount
- 2 Adjustment wheel for pupillary distance
The correct position has been reached when the two eyepiece images merge into one. The interpupillary distance set can be read off the adjustment wheel.



Folding tube f170/f260 mm (optional)

- 3 Eyepiece mount
- 4 Pupillary distance (PD) adjustment knob
- 5 Scale for pupillary distance
- 6 Symbol "normal/larger" for magnification changer
- 7 Dial for integrated magnification changer
- 8 Symbol (released/locked) for manual tube rotation
- 9 Button for manual tube rotation



Controls for the widefield eyepieces



CAUTION

Wide-angled eyepieces with magnetic coupling

Note that the standard rules for handling magnets must be followed for eyepieces removed from the tube:

- Do not place the eyepiece near instruments which may be magnetizable.
- Do not place the eyepiece on sensitive electronic instruments such as infusion pumps, heart pacemakers, measuring instruments or magnetic data carriers such as disks, audio/video tapes and credit cards.
- Always store the eyepiece in its original packaging when not in use.
- Never look at the sun or a light source through the eyepieces or the binocular tube.

1 Eyecup


for adjusting the distance between the eyepiece and the eye.

2 Diopter setting ring

adjustable from -5 D to +5 D

3 Diopter scale

for reading off the prescription set




XY coupling (option)

The X-Y coupling allows motorized positioning of the surgical microscope in a horizontal plane. The range of travel is 30 mm x 30 mm. Use the joystick button (2) on the handgrip or the joystick of the optional foot control panel for positioning. The speed of travel can be set at the display (see page 159). The XY coupling is provided with a reset button. If you press the XY reset button (1), the XY coupling moves to the center position.




1 XY reset button with LED indicator

- The LED is lit as long as you can move the XY coupling in all directions.
- The LED is off when the XY coupling has reached an end stop in one direction.



Additional display


Data in the additional display



You can set the brightness of the additional display to meet your specific needs.

The following current information is displayed to the user on the additional display (see illustration below):


- Light intensity in percent ()
If the light intensity is equal to or higher than the warning threshold, an exclamation mark is displayed after the value.
If the illumination system is switched off, the bulb is shown unlit ()
- Current focus value () (working distance in mm)
- Current magnification factor () (zoom)
- "REC ON" when video recording is in progress, "REC OFF" when no video is being recorded. "REC ON/REC OFF" is not displayed when the recording option is not available.
- Symbol  indicates that a message is shown on the display. If you acknowledge the message on the display, the note on the display and symbol  on the additional display also disappear. The symbol may appear in all menus, with the exception of the standard user menu.
- Autobalance status (Autobalance option, see page 100)
 - 1 Preparation: Do not balance the microscope above the patient
Balancing can be started
 - 2 Autobalance is active
 - 3 Microscope is outside / within the balancing range, see page 102.
 - 4 Balancing procedure cancelled
Repeat the balancing procedure!
 - 5 Autobalance successfully completed



Components of the floor stand

- 1 Extension arm
- 2 Support arm with display
The support arm includes the display for controlling the surgical microscope.
- 3 Stand column
A handle mounted on the stand column facilitates moving the system. There are two cable brackets on the left and right side of the column that can be used to wind up the cables and to attach the foot control panel.
- 4 Stand base
The stand base has four steerable casters (6) enabling positioning of the stand in the operating room.
To prevent the stand from moving or rolling away unintentionally, press down all four locking tabs (7).
Pull up all four locking tabs (7) to release them.
The cable deflectors (5) protect cables lying on the floor from damage after being rolled over.
- 5 Cable deflector
- 6 Steerable casters
- 7 Locking tab

To change the location of the device, see page 95



8 Instrument tray (option)

The instrument tray is designed for holding Zeiss equipment.

**CAUTION****Maximum capacity of instrument tray!**

If the instrument tray is overloaded or accessory tools are insufficiently secured, the risk of injury to the patient is increased.

- Do not exceed the instrument tray's load capacity of 13 kg.
- Ensure all accessory tools are securely placed on the instrument tray. Use tightening belts (incl. in delivery) to securely fix additional accessory tools.

9 17" TFT monitor (option)

The video monitor has a TFT display and is used for dual viewing of the operation for the sterile nurses and other operating room personnel. The flexible support arm can be used to accurately position the video monitor for best viewing.


The video monitor is enabled automatically on switching the system on. During the switch-on process, the video monitor runs through a switch-on sequence in which the signals on the connectors (DVI, VGA, S-Video, Composite and Component) are tested. After detection of the signal present, the correct screen resolution and refresh rate are automatically set.

**CAUTION**

Incorrect diagnosis as a result of video images


The video monitor is not calibrated and not designed for diagnostics purposes. Visualized images may include deviations in scale, shape and color.

- The live image, video sequences and individual images may not be used for diagnostics purposes or treatment.



The xenon lamps generate light whose spectrum resembles that of natural daylight. The color temperature of the light always remains the same regardless of brightness setting. This enables photographic documentation under conditions resembling natural light.


*180W xenon lighting
(standard)*



10 180W xenon lamp housing


The lamp housing with the 180W light source has fiber illumination, a filter wheel (optional) and two xenon lamps.

*300W xenon lighting
(optional)*



11 300W xenon lamp housing

The lamp housing with the 300W light source has fiber illumination, an aperture plate, a filter wheel (optional) and two xenon lamps.




Controls of the xenon lamp housings

180W xenon lamp housing (standard)

- 1 Switch for manual replacement of xenon lamp
- 2 Display: Lamp 2 in use
Lamp 2 is in use if the red segment in the switch (1) is ON.
- 3 Button for opening the xenon lamp module
When this button is pressed, the lamp module is ejected slightly.
For manual lamp replacement, see page 172
For replacing the xenon lamp module, see page 188

Xenon lamp housing
180W



300-W xenon lamp housing (optional)

4 Fastener for opening the lamp housing


5 Housing door

The housing door opens when the fastener (4) is moved in the direction of the arrow.

For automatic lamp change, see page 174

For manual lamp change, see page 176


Xenon lamp housing
300W



Controls for manual balancing

- 1 Balancing of the microscope's lateral tilt movement (manual)
for manual balance setting of the microscope's lateral tilt axis.
- 2 Balance setting knob
for balancing the front to back tilt movement of the surgical microscope
- 3 Fine balancing of the microscope (manual)
for fine balance setting of the microscope
- 4 Crank for weight balancing of the suspension system

For manual balancing, see page 97



Controls on the control and display panel

1 Display with touchscreen function, see page 124

The display is the system's central communication interface. It enables the user to access the settings of the microscope, stand, lighting system, camera (if integrated) and the programmable components of the hand-grips and foot control panel. It is possible to store the selected settings individually for 20 users.

The display is covered by a thin, pressure-sensitive plastic screen so only tip the display with your fingers and not with pointed, hard objects which could damage the screen.

2 Switch for "Light only" mode (see page 24 and page 178)

- All electrical components (except the light source) are isolated from the power supply.
- The light source adjusts to a maximum brightness of 75% to allow you to complete the current operation.

3 USB socket (option)

The USB socket only exists if the "Recording" option has been installed.

NOTE


System malfunctions!

Computer viruses may cause system malfunctions.

- Ensure the USB stick used for the exchange of data is always virus-free.

4 Power switch

for switching the system on and off. The power switch also has an automatic circuit breaker that switches off the device in the event of a short circuit.



Connector panel on the floor stand

- 1 Strain relief device
for securing the power cord and video connection cable
- 2 Potential equalization connector
- 3 Cable connection for the foot control panel
- 4 Remote connector (AUX)
for controlling external devices with a maximum breaking capacity of 24V/0.5A.
- 5 USB connector (for servicing purposes only)
- 6 Ethernet interface (optional)
for networking with OP management systems.
The network connection's cable and connector must comply with Cat-5e EIA/TIA-568A-5 as a minimum, i.e. the Class D values of ISO/IEC 11801:2002 or EN 50173-1:2002.



CAUTION

Risk of patient injury caused by electrical voltage!

- Do not touch the signal interfaces while in contact with the patient.

- 7 Sliding switch for rated voltage
The voltage shown here must correspond to the rated line voltage provided on the site of installation. You can adjust the sliding switch using a suitable tool.





- 8 Power inlet socket
The system must be operated only with 120 V in the USA!
- 9 Power outlet socket
for medical devices of max. 500 VA



CAUTION

Hazard caused by electrical voltage!

Only connect accessories and medical devices intended by Carl Zeiss for use with this system to the power outlet socket (see page 216). When connecting other devices, make sure that safety is guaranteed regarding admissible touch currents and ground leakage currents as per IEC 60601-1:2005.



Video ports with installed 3CCD SD (standard definition) camera (option)



CAUTION

Increased leakage current!

Connecting unapproved or defective accessories may lead to increased leakage current.

- Do not connect any defective or unapproved accessories to the video interfaces.
- Never contact video interfaces while in contact with the patient.



10 Video output "Composite Video" (CVBS)

Analog connector for video components such as recorder, beamer and Zeiss accessories (e.g. MediLive Mindstream).

The video signal is output in accordance with the CVBS standard (shared line for brightness/ color signals). This kind of transmission is suitable when the video signal is to be sent over large distances (e.g. for BNC lines already laid on site).



11 Video output port: Y/C or S-Video (green)

Analog connection for video components, e.g. recorder, projector or Zeiss accessories such as MediLive Mindstream.

Here, the video signal is transmitted according to the S Video or Y/C standard (separate lines for brightness and color signals). This standard provides higher video image quality than VBS.



12 Video input and video output (orange)



Analog connector for video devices such as endoscope camera.

A special cable is required here (10m video in/out connector cable)

CAUTION

Use a galvanic isolation amplifier!

Use video connection cables to other rooms only in combination with a "galvanic isolation amplifier for medical electrical equipment" in accordance with EN 60601-1.



Video ports with installed HD (high definition) camera (option)




CAUTION


Increased leakage current!

Connecting unidentifiable or defective accessories may lead to increased leakage current.

- Do not connect defective or unidentifiable accessories to the video interfaces.
- Never contact video interfaces while touching the patient.




- 13** Video output port: YPbPr
Analog port for high-resolution monitors or TV.



- 14** Video output port: DVI-D
Digital port for LC monitors. DVI (Digital Video Interface) is a technology for the transmission of digital data between the graphics card and the monitor.



When connecting an external HD monitor via the DVI cable 302681-8767-000, make sure to plug the connector with the ferrite core into the stand in order to avoid EMC interference.



- 15** Video output "HD-SDI" (2x)
Digital connector for professional, high-resolution HD monitors. The HD-SDI (High Definition Serial Digital Interface) transfers uncompressed and high-resolution image data in real time. The two connectors are independent of each other and can be configured as required.



Standard BNC cables are not suitable. Only use connector cable "HD-SDI video cable, 75 Ohm, 2 x BNC connector" for these interfaces. Order no.:

305989-8763-000 -	5m long
305989-8764-000 -	10m long
305989-8766-000 -	15m long



- 16** Video output port: Y/C or S-Video (green)
Analog connection for video components, e.g. recorder, projector or Zeiss accessories such as MediLive Mindstream.
Here, the video signal is transmitted according to the S Video or Y/C standard (separate lines for brightness and color signals). This standard provides higher video image quality than VBS.



- 17** Video output "Composite Video" (CVBS)
Analog connector for video components such as recorder, beamer and Zeiss accessories (e.g. MediLive Mindstream).


The video signal is output in accordance with the CVBS standard (shared line for brightness/ color signals). This kind of transmission is suitable when the

video signal is to be sent over large distances (e.g. for BNC lines already laid on site).

CAUTION

Use galvanic isolating amplifier

Only use video connection lines to other rooms in conjunction with a "galvanic isolating amplifier for medical electrical equipment and systems" in accordance with EN 60601-1:2005.



Foot control panel (optional)

The foot control panel enables you to control various functions of the surgical microscope. See page 86 and page 154.




You are only able to activate functions that are configured in the respective system (stand, surgical microscope). Only use the foot control panels mentioned here and in the ordering data:

Follow the correct user manual for the relevant foot control panel.

- Foot control panel with 14 functions (FCP WL), wireless (optional)
- Foot control panel with 14 functions, wired (FCP), 3m (optional)
- Foot control panel with 14 functions, wired (FCP), 6m (optional)

Structure



1 Connector

for connecting the foot control panel with the connector panel of the stand

2 Bracket

for attaching the foot control panel to the stand.

3 Controls


- Buttons A, B, C, D, E, F (freely configurable)
- Joystick (1)
The joystick permits motorized fine movement of the XY coupling.
- Rocker switches G, H, I, J
The rocker switches are used to operate the "Zoom" and "Focus" functions. You can configure the buttons for these functions either horizontally or vertically.

Storage

The bracket on the foot control panel is used to hang it on the column of the floor stand when not in use.

Preconfigured button assignment

Buttons The functions of buttons A, B, C, D, E and F are preconfigured, but they can be re-assigned to meet the user's specific requirements. The factory settings are:




- Button A: DOF aperture -> depth of field low/high
- Button B: light on/off
- Button C: light -
- Button D: light +
- Button E: without function
- Button F: without function

Joystick Joystick (1) is used to control the XY coupling (option).


Rocker switches The focus/zoom rocker switches are preconfigured to the horizontal setting. However, they can be changed to the vertical setting at any time. User-specific settings (see page 154) can be made by yourself at any time. Press button (2) on the display for this purpose.

The horizontal factory settings are:



- Button I: focus rocker switch +
- Button J: focus rocker switch -
- Button G: zoom rocker switch +
- Button H: zoom rocker switch -

The vertical settings are:



- Button I: zoom rocker switch -
- Button J: focus rocker switch -
- Button G: zoom rocker switch +
- Button H: focus rocker switch +

Surgical microscope and laser micromanipulator

A micromanipulator can be attached to the system via the dovetail mount on the bottom of the surgical microscope to permit the use of a laser. Manipulators from other manufacturers can also be attached.

Adjusting the surgical microscope and laser micromanipulator to the same focal plane

The OPMI Vario 700 features a motorized varioscope which is operated via the focus rocker switches of the hand grips, the focus buttons of the foot control panel or an optional autofocus system.

The varioscope is used for motorized adjustment of the working distance and motorized adjustment of the image definition. The focus rocker switches allow you to continuously adjust the working distance between 200 mm and 500 mm.

- Set the working distance to the focus value of the laser micromanipulator. The display and the additional display show the focus value currently set.
- Use the previously described, recommended procedure to check that the focal planes coincide.
- If necessary, correct the focus by appropriate minor adjustment.

The <Fokus Lock> function permits you to deactivate the electrical drive of the focusing system. The focus rocker switches are disabled. This prevents the focal plane setting from being inadvertently changed by motorized movement. If Focus Lock has been activated, no autofocus setting is performed when the brakes are operated, even if the autofocus (optional) has been switched on.

NOTE

If a micromanipulator from another manufacturer is used:

- Check the correct installation and function of the micromanipulator! When using an approved micromanipulator from another manufacturer, activate the <Focus Lock> function (see page 142).
-

Preparations for use



Mounting the microscope components	76
Mounting the tube and the eyepieces.....	76
Mounting documentation / coobservation equipment	78
Connecting navigation systems (option)	80
Connecting the system and accessories	84
Connecting the system to line power.....	84
Connecting wired foot control panel (FCP) to the stand (optional)	86
Pair with wireless foot control panel (FCP WL) - optional	86
Performing foot control pairing (optional).....	88
Fixing cable with strain relief device.....	89
Connecting external video devices to an SD video camera (optional)	90
Connecting external video devices to an HD video camera (optional).....	92
Setting up the system	94
Relocating the system.....	94
Positioning the device in the OR	95
Aligning the video monitor (optional)	96
Manually balancing the surgical microscope	97
Automatically balancing the surgical microscope (optional)	100
Adjusting the surgical microscope	106
Focusing.....	108
Attaching sterile drapes	110

Mounting the microscope components



CAUTION

Risk of injury to the patient caused by lowering of the surgical microscope or by falling parts!

- Never change modules and accessories during a surgical procedure or above the patient!
- After every change in configuration and before every use, make sure that the modules and accessories are securely locked in position. Make sure that all securing screws and locking screws are firmly tightened!

Mounting and removing modules and accessories may lead to unintended tilting of the surgical microscope.

- Always make sure that the maximum permissible load is not exceeded.
- After every change in configuration and before every use, rebalance the completely equipped surgical microscope without the patient and at a safe distance from other persons and equipment.

Mounting the tube and the eyepieces

- Give securing screw(6) a few turns to loosen it.
- Remove cover (5) and keep it in a safe place.

You can attach different tubes (e.g the tiltable tube (2)) directly to the microscope body. First attaching the stereo bridge (1) to the microscope body and then the tubes (2 and 8) is recommended for many applications.

- Place the stereo bridge (1) onto the microscope body and tighten the screw (6) .
- Place the tubes (2 and 8) onto the stereo bridge (1) and tighten the screws (4 and 10) .
- Fully insert the wide-angled eyepieces(3 and 9) into the mounts provided on the binocular tube. The magnetic coupling reliably secures them in position (refer to page 53).



Remove the transportation protective cover (7) from the lens.

If the lighting is switched on, the cover can melt under the heat and cause irreparable damage to optical components.

Mounting documentation / coobservation equipment

The camera adapter and coobservation module shown in the illustration are examples of further accessories that can be mounted on the lateral image exit ports. The method described below can also be used for other accessories. The operating principle of the accessories is described in the relevant user manuals.

- Loosen the respective knurled ring (e.g., 4).



An arrow indicating the "open" direction is located on the shaft next to the knurled ring (e.g., 4).

- Remove dust cover (2) and store it in a safe place.
- Insert the accessory module (1 or 5) into the respective opening of the left or right image exit port.
The mount of the image exit port is fitted with guide projections. The accessory module (1 or 5) has the corresponding grooves. Carefully turn the accessory module until the guide projections fit into the grooves, and slide the accessory module into the opening as far as it will go.
- Screw the corresponding knurled ring (e.g., 4) onto the accessory module.
- Firmly tighten the knurled ring.




Please note that the accessory attached to the documentation port on the left-hand side of the microscope can collide with the microscope attachment for swivel angles greater than -45° .

- If possible, and if an extended swivel range is an absolute requirement, attach the accessory to the right-hand side of the microscope.



If you wish to use documentation equipment, we can supply an eyepiece with a reticle to aid focusing. The retrofitting of a reticle to an eyepiece can only be performed in the factory or by our service staff. Always install the eyepiece with the reticle on that side of the binocular tube where the documentation equipment is located.

Maximum loading
on microscope body



Connecting navigation systems (option)



OPMI VARIO 700 is prepared for the connection of navigation systems.

For the use of a navigation system, the option must be ordered ex factory or installed and activated by Carl Zeiss service staff.

NOTE

The link-up using the navigation interface creates a medical system for which the system supplier (manufacturer of the navigation system) must meet the stipulated requirements (approval, qualifications, etc.). All accompanying papers required will be supplied by the manufacturer of the navigation system.

The navigation interface for OPMI VARIO 700 was developed in cooperation with Carl Zeiss in compliance with the requirements of the communication protocol of the OPMI VARIO 700 navigation interface.

The use of the OPMI VARIO 700 with a connected navigation system requires calibration of the microscope using a technique to be made available by the manufacturer of the connected navigation system. This allows OPMI VARIO 700 in combination with the navigation system to be used like an optical pointer with a variable length (corresponds to the working distance).

- OPMI VARIO 700 can manage various connected navigation systems.
- For safety reasons, the OPMI VARIO 700 navigation interface is not enabled unless a connected system authenticates itself.
- Before every surgery using a connected and authenticated navigation system, the function and accuracy of the navigation system must be verified (e.g., by focusing on a measuring point or comparing the focus with the instrument to be navigated). Please also observe the information contained in the respective user manual for the navigation system used.

**CAUTION****The use of third party drapes may affect navigation systems**

Certain third party drapes (e.g. drapes with curved protective lens covers) may influence the accuracy of navigation systems due to an optical distortion caused by the convex shaped protective lens cover.

Although there may be no perceptible impairment to the visual image quality observed through the microscope eyepieces (oculars), there is an imaging effect which displaces the microscope's focal point displayed by the navigation software. Due to this effect the focal point displayed by the navigation software may not be correctly aligned to the patient's image information from the navigation system (i.e. the patient data set).

Consequently the microscope-related reconstruction of the patient data set by the navigation system could display a wrong focal plane that does not comply with the real focal plane of your surgical microscope.

In systems with such features, depending on the distortion caused by the lens, the object contours of a pre-planned anatomical structure displayed in the microscope video of the navigation software, or in the heads-up-display, may show an incorrect size and/or incorrect position in relation to the actual image in your surgical microscope.

If not noticed by the user, incorrectly displayed information could influence clinical decisions during a navigated surgery. Those clinical decisions could be incorrect and potentially result in ineffective treatment or in serious injury to the patient.

- Always follow the instructions for use of the navigation system manufacturer for calibration, calibration verification and operation.

Connecting and operating navigation systems


The navigation interface of the OPMI VARIO 700 only allows connection and use of systems from authorized manufacturers. Authorized manufacturers are companies or institutions having concluded an appropriate open-interface agreement with Carl Zeiss and having obtained a license for using the integrated navigation interface.

Please see the user's manual provided by the manufacturer of the connected system.

Navigation interface (option)

This navigation interface of the OPMI VARIO 700 comprises the following system components:

- 1 Joystick (right handgrip)
The joystick button on the right handgrip can be used to control software functions of the connected navigation system if this button is supported by the navigation system.
- 2 Connector (example)
for an external microscope navigation system.
- 3 Interface on the surgical microscope
for connecting an antenna module (example)
- 4 Software
for controlling data exchange



Connecting the system and accessories

Connecting the system to line power



In the United States, the device may only be operated at max. 120 V!



CAUTION

- According to EN 60601-1:2005, Section 16.2, the power outlet socket of this system must be a multiple connector that is intended for an ME system. When setting up your ME system (use of the power outlet socket), ensure to comply with the requirements of IEC 60601-1-1:1990 or IEC 60601-1:2005, Section 16. Please also observe the following:
- Never place multiple connectors on the floor.
- Never connect any additional multiple connectors.
- Only connect compatible components to the system.
- Ensure that you do not exceed the max. permissible load capacity of the multiple connectors.
- Only use multiple connectors for components which are part of the system.



CAUTION

Risk of tripping over!

Inappropriately laid cables represent an increased risk of tripping over.

- Always lay cables such that procedures are not hindered.

NOTE

Risk of damage caused by incorrect setting of the line voltage

The line voltage indicated at the window must correspond to the rated voltage available on the site of installation; otherwise, the suspension system will be damaged.


If the line voltage has not been correctly set, adjust the sliding switch in the indicator window using a suitable tool.

- Ensure the power switch (1) is OFF.
- Insert the power cord supplied into power inlet socket (4) on the stand.
- Only connect the power cord to wall outlets having a properly connected protective earth conductor.
- Secure the power cord using the strain relief device (see page 89).

- If required, connect the stand to the potential equalization rail in the operating room. Use the potential equalization bolt (3) for this. This connector can be used for placing other active devices at the same ground potential or for redundant grounding to protective earth.


CAUTION
Risk of patient injury caused by electrical voltage!

Defective or unidentified accessories or Ethernet cables may lead to increased leakage currents and injure the patient. Do not connect defective or unidentified accessories or Ethernet cables and never touch the power outlet socket or other signal interfaces whilst the patient is being treated.



Connecting wired foot control panel (FCP) to the stand (optional)

- Ensure the power switch (1) is OFF.
- Secure the FCP cable with the strain relief device (see page 89).
- Connect the FCP cable to the connector (2) on the stand.

Pair with wireless foot control panel (FCP WL) - optional

Pairing means the fixed relative assignment of the suspension system and foot control panel.

It is required for wireless operation. When pairing the system for the first time, a 20 sec. delay until a signal link between suspension system and foot operator panel has been established may be possible.

Proceed as follows for pairing:





- Turn on the suspension system at power switch (1).
- Put the foot control panel in a vertical position (A) in the immediate vicinity of the suspension system (distance less than 2 m).
- Start the pairing process on the suspension system as described on page 88.

NOTE

Foot control panel non-operational!

If pairing is performed incorrectly or not at all, the foot control panel may be disabled, or activation of a control may trigger functions on a different suspension system not assigned to the foot control panel.

- Perform pairing.
- If the foot control panel continues to be non-operational after pairing, connect it to the system via a 3 m or 6 m cable.



Performing foot control pairing (optional)




Push the <System Config> button in the main menu and press the <Foot control pairing> button to open the "Pairing" submenu. In this submenu, you can assign a wireless foot control panel which is required for wireless operation to the system.

Before pairing, go through the following steps:

- Unless already done, place the foot control panel in its vertical position in the immediate vicinity of the system (max. distance of 1 m) and keep it in this position until pairing is completed.

Performing pairing

- Start the "Pairing" process by pressing the <Start Pairing> button on the touch screen.
 - A message appears: "Activate one of the foot control panel's buttons until the third LED from the top blinks orange. Maintain the vertical position of the foot control panel."



Successful pairing

If pairing was successful, the "Radio link intensity" indicator (1) flashes green for approx. 1 sec and the following message appears:

- "Pairing successfully completed. Put the foot control panel in the horizontal position (working position) and perform a function test. Set the number specified on the suspension system via the wheel on the foot control panel."
- Check whether pairing was successful by simultaneously pressing any two control buttons (2) on the foot control panel. The "Radio link intensity" status indicator will be lit.

Pairing failed

- If pairing failed, the "Radio link intensity" indicator (1) flashes red for approx. 1 s and one of the following message appears:
 - "Pairing failed. No foot control panel was detected."
 - "Pairing failed. Several foot control panels were detected."
- Repeat the pairing process as described above.

Fixing cable with strain relief device




Use the strain relief device to prevent the connected cables from inadvertently being pulled out.

Make sure that the remaining cable length between the strain relief device and the connector (on the connector panel) is approx. 320 mm.

Secure the cable in the strain relief device as follows:

- Form a loop with the cable (2). Do not kink the cable.
- Open the flap (3).
- Feed the cable through the opening (4).
- Close the flap (5).
- Pull the cable until it encloses the flap (6).
Do not kink the cable.

Check the cable length.



Connecting external video devices to an SD video camera (optional)

The OPMI VARIO 700 with integrated SD video camera allows you to connect additional SD video monitors to the system. Parallel SD data recording is possible by connecting the optionally available external record option: MEDIALINK 100 or with the optional record function via the USB port.



CAUTION

Hazard caused by live cables and connectors!

The connection of accessories that have not been approved by Carl Zeiss may lead to injury of the patient.

- When setting up your system, make sure to comply with the requirements of EN 60601-1:2006 Chapter 15 or the requirements of EN 60601-1-1:1990.

Connecting an endoscope camera

- Connect the analog video output port (Video in/out) on the OPMI VARIO 700 and the video input port of your endoscope camera using the following special cable:
 - 10m Video in/out connecting cable


Connecting the SD video monitor (integrated video recording function via USB port)

- Connect the analog video input port (Y/C) on the OPMI VARIO 700 to the video output port (e.g. Y/C) of your SD video monitor using the following video cable:
 - 10m S-Video connecting cable

Connecting the SD video monitor (external MEDIALINK 100)

- Connect the analog video output port (Y/C) on the OPMI VARIO 700 to the video input port (e.g. Y/C) of your MEDIALINK 100 using the following video cable:
 - 10m S-Video connecting cable
- Connect the analog (DVI-D) or digital video output port (Y/C or BNC) of your MEDIALINK 100 and the video input port of your SD video monitor using one of the following video cables (depending on the connector type):
 - 10m S-Video connecting cable
 - 2m BNC cable set
 - 5m DVI-D system cable

*External video devices
on integrated SD camera*



Connecting external video devices to an HD video camera (optional)

The OPMI VARIO 700 system with integrated HD video camera allows you to connect HD video monitors or external recording systems to view high-resolution images.

Parallel SD data recording is possible by connecting the optional MEDIALINK 100 recorder.



CAUTION

Hazard caused by live cables and connectors!

The connection of accessories that have not been approved by Carl Zeiss may lead to injury of the patient.

- When setting up your system, make sure to comply with the requirements of EN 60601-1:2005 Chapter 15 or the requirements of EN 60601-1-1:1990.

Connecting an HD video monitor




When connecting an external HD video monitor with the 5m DVI-D system cable (302681-8767-000), ensure to plug the connector with the ferrite core into the stand to prevent EMC interference.

- Connect the digital video output (DVI-D) on the OPMI VARIO 700 to the video input (e.g. DVI-D or BNC) of your HD video monitor using one of the following video cables (dependent on connector type):
 - HD-SDI connector cable, 75 Ohm, 2x BNC pin
 - 5 m DVI system cable

Connecting an SD video monitor (including MEDIALINK 100)

- Connect the analog video output (Y/C) on the OPMI VARIO 700 to the video input (e.g. Y/C) on your MEDIALINK 100 using the
 - 10 m S-Video connector cable
- Connect the video output (Y/C, BNC or DVI) of your MEDIALINK 100 to the video input on your SD video monitor using one of the following video cables (dependent on connector type):
 - 10 m S-Video connector cable
 - 2 m BNC cable set
 - 5 m DVI system cable

*External video devices
on integrated HD camera*



Setting up the system

Relocating the system



CAUTION

Risk of tipping!

Be careful when moving the device across sills; it may easily tip and/or injure a person next to it.

- When lifting or pulling the system across a doorsill request the help of a second person.



CAUTION




Risk of crushing!

Fingers can become trapped between support arm and extension arm.

- Never touch the area between support arm and extension arm while turning the support arm.

NOTE



Transporting the device

- Turn on the system from the power switch.
- Move the spring arm into the transport position (see figure on left).
- Turn on the system from the power switch.
- Disconnect the power cord from the power socket.
- Wrap the foot control panel and mains cables onto the cable holders.
- Hang the foot control panel on the bracket.
- Release the locking tabs of the casters.
- Hold the device by the bracket and push it carefully. Ensure the device is standing on even ground.
- Be careful of heights when passing through doorways.
- Avoid collisions of any kind.
- Do not go over steps and edges (the stand may topple!)
- Be extremely careful when moving over slopes.
- To lock the device into position, press the locking tabs of all four casters and ensure that the stand cannot be moved unintentionally.


NOTE**Transporting the device over long distances**

- Only transport the device over long distances (e.g. relocation, return for repair) in its original packaging or special return packaging. Please contact your dealer or Carl Zeiss Service.

Positioning the device in the OR

Position the device as follows:

- Release the locking tabs of all casters (2).
- Hold the device by the transport handle (1) and move it carefully into a position convenient for you (such as that shown below).
 - Ensure that the support arm (3) and the extension arm (4) are not fully extended but are at an angle of ca. 90° to each. This provides a greater degree of maneuverability and also a perfect working position.
 - Ensure the device is standing on even ground and is not tilting or wobbling.
- Press the locking tabs of all four casters (2) and ensure that the stand cannot be moved unintentionally.



Aligning the video monitor (optional)



CAUTION



Risk of crushing!

Hands and/or fingers may be crushed between monitor and carrier arm.

- Stay clear of the area between monitor and carrier arm, while moving the monitor.

NOTE


The video monitor can drop downwards

Aging processes can lead to a loss of gas in the gas pressure spring of the spring arm. The video monitor can drop down as a result.

- Compensate for the loss of gas by readjusting the gas pressure spring (see page 192).
- If the video monitor continues to drop downwards, the gas pressure spring is defective. Please inform the Carl Zeiss service department.

How to align the video monitor:

- Pivot the carrier arm (1) and the suspension arm (2) in the desired horizontal position.
- Tilt the suspension arm (2) upward or downward until the proper height is reached.
- Hold on to the upper corner of the video monitor and adjust the proper viewing angle.



Manually balancing the surgical microscope



CAUTION

Risk of injury caused by uncontrolled movement of the system!

It may happen in exceptional cases that the surgical microscope is not correctly balanced although the balancing procedure has been completed. With an incorrectly balanced system, brake release may lead to uncontrolled movements of the suspension system.

- Only use the system after it has been properly balanced!
- To check correct balancing of the system, unlock the brakes. Hold the microscope tightly at both handgrips. If the system has been correctly balanced, the surgical microscope can be moved almost effortlessly. Repeat the balancing procedure, if required.

The balancing procedure and the subsequent test must not be performed above the patient and only at a safe distance from other persons and devices.

Move the stand into position for manual balancing

Turn the surgical microscope in the direction of the stand base and move the stand into a position (see diagram) so that the magnetic brakes and the crank of the balance adjustment can be operated.






Image 1




Image 2



Image 3



Manual balancing through user guidance on the display

Press the <Balance> button (⚖️) on the display in the main menu; images are used to explain the manual balancing procedure.

Press the < > button to access the previous or next image.

Press the <Stop> button to return to the main menu.

i In general, the following applies:

If the microscope moves in the direction of the white or black arrow when the brakes are unlocked, the respective knob must be turned in the direction of the white or black arrow.

1 Adjusting the balance of the suspension arm

- Unlock the stand brakes by pressing the respective handgrip button and turn crank (1) to adjust the balance in such a way that the suspension arm with the microscope remains stationary in any required position.

2 Adjusting the microscope's front-to-back tilt motion

- Unlock the microscope brakes and put the microscope in the position shown in image 2.
- Adjust the balance by turning knob (2) in such a way that the surgical microscope can be moved almost effortlessly through the front-to-back tilt axis.

3 Adjusting the microscope's lateral tilt motion


- Unlock the microscope brakes and put the microscope in the position shown in image 3.
- Adjust the balance by turning knob (3) in such a way that the surgical microscope can be moved almost effortlessly through the lateral tilt axis.

i Procedures 2 and 3 may influence each other.

Therefore, adjust the front-to-back and lateral tilt ranges in such a way that optimum balance is obtained and you can move the surgical microscope almost effortlessly.

i The balancing procedure is completed when the microscope has been balanced to your satisfaction. Perform fine adjustment (see page 104) if the balance is not satisfactory, e.g. when the microscope is used in horizontal viewing direction.



Adjusting the swivel movement of the microscope



Adjusting the tilt movement of the microscope



Adjusting the balance of the extension arm



Automatically balancing the surgical microscope (optional)

The device can be fitted with the AutoBalance option.



CAUTION

Risk of injury caused by uncontrolled movement of the system!

It may happen in exceptional cases that the surgical microscope is not correctly balanced although the balancing procedure has been completed. With an incorrectly balanced system, brake release may lead to uncontrolled movements of the suspension system.

- Only use the system after it has been properly balanced!
- To check correct balancing of the system, unlock the brakes. Hold the microscope tightly at both handgrips. If the system has been correctly balanced, the surgical microscope can be moved almost effortlessly. Repeat the balancing procedure, if required.

The balancing procedure and the subsequent test must not be performed above the patient and only at a safe distance from other persons and devices.




Ensure that any accessories mounted (such as stereo dual viewer) do not collide with the stand during the balancing process.

The balancing process can be started by pressing the <AutoBalance> button (1) on the display or by pressing the balance button (2) on the additional display.

The AutoBalance process can be aborted at any time by pressing the <AutoBalance> button (1) on the display or by pressing the balance button (2) on the additional display. Releasing the magnetic brakes also stops the balancing process.

- 1** Press the <Balance> button on the display → AutoBalance starts
The user is kept up-to-date on the process with messages on the display and symbols on the additional display (see following pages).
- 2** Balancing process by pressing the Balance button → AutoBalance starts
Access to other functions of the device is temporarily disabled during the balancing process.
The user is kept up-to-date on the process with messages on the display and symbols on the additional display (see following pages).




Surgery microscope is within balancing range


The following balancing process applies if the surgical microscope is already in the balancing range (center position):

Display







Additional Display











"Autobalance successfully completed" display automatically disappears from the display/additional display after 3 seconds. If this is too long for you, press the button below the display.

Error message during AutoBalance

Balancing procedure cancelled - Repeat the balancing procedure!

Pay attention to the error messages on the display. The additional display shows a warning triangle in this case, and a beep is heard.






Surgical microscope is outside the balancing range

Surgical microscope is outside the balancing range. Swivel the microscope (with brake released) in the microscope axes until the two marks (1) of the axes are in the balancing range (2). The following balancing process applies:


Display










Additional Display




Bring the microscope into the balancing range (2)




The LED on the additional display blinks.



The LED on the additional display lights up




i "Autobalance successfully completed" display automatically disappears from the display/additional display after 3 seconds. If this is too long for you, press the button below the display.

Surgical microscope is not configured perfectly

The weight distribution of the accessory on the surgical microscope is not perfect. AutoBalance must first perform an approximate balance. This is followed by a fine balance. Swivel the microscope (with brake released) in the microscope axes until the two marks (1) of the axes are in the balancing range (2). The following balancing process applies:

Display

Additional Display

Bring the microscope into the balancing range (2)



"Autobalance successfully completed" display automatically disappears from the display/additional display after 3 seconds. If this is too long for you, press the button below the display.



Perform fine adjustment (see page 104) if the balance is not satisfactory, e.g. when the microscope is used in horizontal viewing direction.


Fine adjustment of microscope

Ideally, adjustment is sensitive and often only half a turn of the knob is sufficient (4).

For large loading or unsymmetrical configuration, greater corrections are necessary with the knob (4).

i Generally, if the microscope body moves in the direction of the black/white arrows after disengaging the brake, the appropriate knob must be turned in the direction of the black/white arrows.

Picture 1




- Release the microscope brakes and position the microscope as shown in Figure 1.
- Adjust the balance by turning the knob(4) so that the surgical microscope can be tilted backwards and forwards with very little force in this position.
 - If the microscope swivels forwards when the brake is released (black arrow), turn the knob (4) in the direction of rotation of the black arrow.
 - If the microscope swivels backwards when the brake is released (white arrow), turn the knob (4) in the direction of rotation of the white arrow.


For the AutoBalance option:

- Re-run an Autobalance in this position.

Picture 2



- Release the microscope brakes and re-position the microscope as shown in Figure 2.
- Correct the balance again by turning the knob (4).




Adjusting the surgical microscope



- Only a correctly adjusted working distance will provide an optimum image at all selectable magnifications.
- For quick setup of the surgical microscope, it is advisable to note the pupillary distance and prescriptions of the different users so that they can be easily preset during system preparation.

Adjusting the interpupillary distance

- Use knob (1) to adjust the distance of eyepieces (2) on the binocular tube in such a way that the two eyepiece images merge into one.
- Set the minimum magnification on the microscope via the handgrips or foot control panel.
- Position the microscope above the surgical field.



- Adjusting the eyepieces*
- Emmetropes:
 - Set diopter setting ring (5) to 0 diopters (D).
 - Eyeglass wearers who perform surgery wearing their eyeglasses:
 - Set diopter setting ring (5) to 0 D.
 - Eyeglass wearers who perform surgery without their eyeglasses (known prescription):
 - Set diopter setting ring (5) to the respective value.
 - Eyeglass wearers who perform surgery without their eyeglasses (unknown prescription):
 - Set diopter setting ring (5) on the eyepiece to + 5 D.

Eyepieces without reticle	Eyepieces with reticle
<ul style="list-style-type: none"> • Remove tube with eyepieces from the microscope body and point it at a distance object • Slowly turn diopter setting ring (5) of one eyepiece clockwise until the object is in focus. • Mount the tube with eyepieces on the microscope body and tighten the screw (3). • Adjust eyecups (4) in such a way that the entire field of view can be seen. 	<ul style="list-style-type: none"> • Slowly turn the diopter setting ring (5) on the eyepieces clockwise until the reticle is in focus. • Focus the main microscope on the object. The reticle and the object must be in focus at the same time.
<ul style="list-style-type: none"> – Viewing with eyeglasses: Screw in eyecups (4) all the way. – Viewing without eyeglasses: Screw out eyecups (4) and adjust them for your field of view. 	

Focusing

Manual focusing

- Move the surgical microscope over the surgical field (1) into an ergonomic position within the working distance.
- Set the minimum magnification (zoom) via the handgrips or foot control panel.
- For coarse focusing, look through the eyepieces and lower the surgical microscope until the surgical field comes into focus.
- Select the highest magnification (zoom).
- Look through the eyepieces and activate the focusing function until the image of the surgical field appears sharply focused.
- Select the magnification required (zoom). Look through the eyepieces of the binocular tube. Adjust the eyepieces in such a way that you can see both the edge of the field of view and the microscope image sharply.

Using the SpeedFocus (autofocus) option for focusing.



The SpeedFocus autofocus function is based on the evaluation of the video image contrast: during focusing of the microscope's Varioskop optics, SpeedFocus continuously determines the video image contrast in the selected ROI. In this way, SpeedFocus controls the microscope's Varioskop optics and adjusts it to the point of maximum contrast.


- Move the surgical microscope over the surgical field into an ergonomic position within the working distance.
- Press the configured <SpeedFocus> button on the handgrip or the foot control panel.
 - SpeedFocus has focused1 beep
 - Focusing could not be performed 3 beeps

Focusing the video image manually on the optional video monitor

- Press the configured <Video focus +/-> buttons on the handgrip or the foot control panel. Or: set the focus using the <+ (4) / - (5)> buttons in the <Focus> menu.

Focusing the video image automatically on the optional video monitor (2)

- Press the configured <Video SpeedFocus> button on the handgrip or the foot control panel.
The video image is automatically in focus, regardless of the microscope image actually seen by the surgeon.



Attaching sterile drapes



CAUTION

The use of third party drapes may affect navigation systems

Certain third party drapes (e.g. drapes with curved protective lens covers) may influence the accuracy of navigation systems due to an optical distortion caused by the convex shaped protective lens cover.

Although there may be no perceptible impairment to the visual image quality observed through the microscope eyepieces (oculars), there is an imaging effect which displaces the microscope's focal point displayed by the navigation software. Due to this effect the focal point displayed by the navigation software may not be correctly aligned to the patient's image information from the navigation system (i.e. the patient data set).

Consequently the microscope-related reconstruction of the patient data set by the navigation system could display a wrong focal plane that does not comply with the real focal plane of your surgical microscope.

In systems with such features, depending on the distortion caused by the lens, the object contours of a pre-planned anatomical structure displayed in the microscope video of the navigation software, or in the heads-up-display, may show an incorrect size and/or incorrect position in relation to the actual image in your surgical microscope.

If not noticed by the user, incorrectly displayed information could influence clinical decisions during a navigated surgery. Those clinical decisions could be incorrect and potentially result in ineffective treatment or in serious injury to the patient.

- Always follow the instructions for use of the navigation system manufacturer for calibration, calibration verification and operation.




Pull the drape until the marked position is reached .
Do not cover the ventilation openings of the light source or carrier arm.

System without Autodrape function


- Unpack the drape and pull it over the surgical microscope and the suspension arm.
Refer to the instructions for use of the sterile drapes!
- Pull the drape over marked position (2) of the suspension arm.
When attaching the drape, please ensure that there is sufficient space for the tilt and rotation movements of the surgical microscope.
- Make the drape airtight at the marked position (2) using one of the straps.



If the AutoDrape option has not been installed, <AutoDrape> button (1) is disabled.



Effective range of drape vacuum system



System with AutoDrape function (optional)

The system features an integrated AutoDrape function which extracts the air from the drape at the suspension arm and surgical microscope.

- Press <AutoDrape> button (1) in the main menu to activate the drape vacuum pump.



The vacuum system operates at maximum power during the first 3 minutes approximately. Then it automatically switches to a defined maintaining power to keep up the vacuum.

- The drape vacuum pump is switched off by pressing <AutoDrape> button (1) in the main menu once again.

We recommend not to switch off the AutoDrape function during surgery.



CAUTION


Contamination with pathogenic germs!

- Observe the instructions for use of the sterile drape and make sure that the drape is attached in such a way that sterility is not impaired (see sterile area in the illustration on the next page and page 200).



1

Effective range of drape vacuum system



Operation



Powering the system up and down	118
Functional test prior to use (checklist)	120
Display	124
Symbols on the display	124
Controls on the display	126
Menu structure	128
Overview	129
User	130
Selecting user and language	130
Adding, deleting, editing a user	132
Light - Adjusting the light	136
How to adjust illumination settings	138
Configuring the video focus function (Options: SpeedFocus, video focus)	142
Configuring the focusing function	144
Setup of record function	146
(Integrated video record function over USB port - option)	146
Perform white balance (option)	148
User settings menu 1	150
Configuring zoom	150
Configuring the handgrips	152
Configuring the foot control panel	154
Configuring the brakes	158


User settings menu 2	159
Configuring the XY coupling (optional)	159
Camera setting 1/2 (option)	160
Camera settings 2/2 - Color setting (option)	162
Configuring video source (camera) - only for SD camera option	163
System settings menu 1	164
Optics	164
Configuring the video settings (optional)	165
Adjusting the brightness of the additional display	168
System settings menu 2	169

Powering the system up and down

- Prerequisites*
- ✓ The correct rated voltage for the stand has been set (see page 66).
 - ✓ The power cable is connected (see page 84).



CAUTION



Risk of injury!

- Do not switch the device on above the patient.
- The device can perform uncontrollable movements or functions (focus, zoom, light control, XY movement(optional)) if the electronics system experiences malfunctions.
- In this case, switch to Light-only mode by pressing the orange-colored Light-only button (1) on the control and display panel.



CAUTION

Malfunction

- After powering on the device, check the system starts correctly. The device is ready for use after the signal tone and display of the start screen:



Proceed as follows to power the system up and down:

- Press the power switch (2) on the stand.
 - The power switch lights up green when the device is switched on.
 - The system performs a self-test after power-up during which a start screen with the Zeiss logo is shown on the display.
 - The power switch does not light up when the device is switched off.

Switching on the surgical microscope for the first time displays the default user profile "Default User". For subsequent log-ins, the system sets the user profile of the user last selected.



The "Default User" profile cannot be edited or changed. It cannot be deleted either.

A functional test (checklist) must be performed before starting to work with the system (see next page). Check the device configuration and the user-specific settings and re-set them if necessary.

Functional test prior to use (checklist)

NOTE

Functional test

- Use the following list to check the function of the device before every operation (with no patient).



CAUTION

Risk of injury!

Incorrect user settings may lead to unexpected behavior of the system and may cause injuries to the patient.

- Always check the settings of the selected user profile before using the system.



CAUTION

The values displayed (SI units) are not calibrated

The values displayed (SI units) are rounded. They are intended to improve the display and do not represent a measurement function.

NOTE

Software faults!

Software faults can impact the behavior of the device.

- In the event of a software fault, switch to Light-only mode and inform our service department.

Microscope

- Check the zoom
 - ✓ The buttons on the handgrips and foot control panel work as per the configuration (see page 152).
 - ✓ The zoom speed has been set on the display and works.
- Check the focus
 - ✓ The buttons on the handgrips and foot control panel work as per the configuration (see page 152).
 - ✓ The focus speed has been set on the display and works.
 - ✓ SpeedFocus (option) has been set on the display and works (see page 108 and page 144).
 - ✓ Video focus (option) has been set on the display and works (see page 108).
- Check the unobstructed movement of the surgical microscope

- ✓ The balance has been set and the microscope is balanced.
- ✓ All accessory parts have been fitted such that the surgical microscope can be positioned safely and without causing damage.
- Check the eyepieces/binocular tube/stereo bridge
 - ✓ The surgical microscope and the tube are in a position convenient for you.
 - ✓ The interpupillary distance has been adjusted (see page 106).
 - ✓ The eyecups have been adjusted such that you can see the full field of view (see page 106).
 - ✓ Check that the correct prescription has been set on the diopter scale (see page 106).
 - ✓ The image quality is correct throughout the entire magnification range.
 - ✓ The lens and eyepieces are clean.
- Check accessories (also third-party accessories)
 - ✓ Other equipment (lighting, video system and third-party accessories such as micromanipulator (see page 74)) has been checked for proper function using the user manuals.

**CAUTION****Loss of light for HDTV components**

Using HDTV components causes a minor loss of light.

- Check and correct the lighting of the surgical field - as far as is possible for the safety of the patient and for a good microscopic image.
- Adjust the luminous field diameter and the lighting intensity to the procedure size required.
- Do not correct the brightness of the camera image by increasing the light intensity of the microscope lighting, instead use the settings on the controller of the HDTV camera.

Suspension system

- Check the balancing
 - ✓ When properly balanced, the surgical microscope must remain stationary in any position of the working range and move neither upward nor downward when the magnetic brakes are unlocked.
- Check the lamp brightness




Please note that, when using the Superlux 180 W light source, the system cannot automatically detect a defective or missing lamp.

- Please ensure that 2 working lamps are always in the system and that defective lamps are replaced immediately.

NOTE

Function test!

- ✓ Perform a function test and, if required, replace the bulbs prior to using the system.
Operate the system only with two functional lamps.
Keep a replacement lamp ready at hand.
- ✓ With increasing age of the light source, the actual illumination intensity delivered at the respective setting decreases. Please replace the xenon lamp in due time. Note the remaining service hours display.
- ✓ The ventilation openings are not covered because the light source of the system switches off in the event of overheating.
 - ✓ The lamp brightness can be varied (see page 138). Brightness variation has an effect on the surgical field illumination.
- Check the XY coupling (optional)
 - ✓ The XY coupling can be positioned over the entire adjustment range. Press the reset button to move the XY coupling to the center position.
- Check the stand for secure positioning
 - ✓ All four brake tabs of the casters have been pressed down and the stand cannot roll away by itself (see page 95).
 - ✓ The stand is in a stable and safe position.
- Check the user settings on the display
 - ✓ The correct user profile has been selected and the settings have been checked.
- Check whether the stand is unobstructed in its movement
 - ✓ The system has been positioned in such a way that it can be moved away from the patient at any time.



Magnetic brakes

- Check proper functioning

On the rear of the hand grips are brake release buttons (1) ("All Brakes") and (2) ("Selected Brakes").

- ✓ The brake release buttons for all stand **and/or** microscope axes work as per the configuration (see).
- ✓ The magnetic brakes are locked when you release the button.



Please note that the magnetic brakes can engage briefly during system start - even when the button is pressed.

- Do not keep the button releasing the magnet brakes pressed for the whole time of the system start.
If the brakes engage nevertheless, please release the button for a few seconds and then press it again.

Foot control panel (optional)

- Check settings

- ✓ The buttons on the foot control panel work as per the configuration (see page 154).

NOTE



Check batteries

- Prior to using the FCP WL wireless foot control panel, ensure its batteries are charged sufficiently. The system may malfunction if the FCP WL is not provided an adequate level of power.

Display


Symbols on the display

Symbol or abbreviation	Definition
	User menu, see page 130
	Light menu, see page 136
	Focus menu, see page 142
	Balance menu, see page 97 (option: AutoBalance, see page 100)
	AutoDrape (option), see page 112
	Integrated Record function via USB port / network, see page 225
	User settings , see page 150 and page 159 Light, zoom, handgrip, focus, brakes, foot control panel (option) camera (option), XY coupling (option, video source (option)
	System settings , see page 164 and page 169 optics, foot control panel pairing (option), camera, additional display service PIN, export of service files, configuration
	Navigation menu (option)
	White balance menu (option), see page 148

Symbol or abbreviation	Definition
 	Shows the light intensity as a % (e.g. 21%) The lamp is shown gray if the lighting is switched off
	Shows the working distance in mm (focus), e.g. 316 mm
	Shows the overall magnification of the system (e.g. 1.6x)
	Displays the charge state of the battery in the wireless foot control panel (option). The number of bars in the symbol gives the charge state.
	Shows that a radio link to the foot control panel (option) is active. The number of bars in the symbol gives the signal strength.
REC	Displays a running video recording if option "Record" is installed.
	The light intensity, focus and zoom values displayed in the header are shown larger if the user touches the header.
	Notifies the user that a message is being shown on the display. This disappears on acknowledging the message.

Controls on the display

OPMI VARIO 700 provides the following controls:

Control element	Definition
	Press the button in the main menu to access a submenu (e.g. Light).
	If the menu comprises two pages, press the tabs on the right-hand side to open the page required.
	An info display appears to inform the user about the current or completed procedure.
	Press this button to close the info display.
	Press the <Previous> button to return to the previous menu, to cancel changes or to retain changes until the user is changed or the system is re-started. Press this button for more than two seconds to access main menu 1.
	On the lower bar in the submenus, you find the <Save> button. Press this button to save the current settings.
	Press the arrow buttons to increase or reduce the numeric values of functions (e.g. light intensity). The current value is displayed in the middle. When the highest or lowest value has been reached, the respective arrow button will be grayed out. The bar on the left graphically indicates the set value.
	Press the - / + buttons to increase or reduce the numeric values of functions. The current value is displayed in the middle. When the highest or lowest value has been reached, the respective arrow button will be grayed out.
	Use the buttons to activate and deactivate functions. <ul style="list-style-type: none"> – Light blue button => function is deactivated – Dark blue button => function is activated – Gray button => function is not available

Control element	Definition
	<p>If text entries need to be made, a virtual keyboard is shown on the display and the following input options are offered:</p> <ol style="list-style-type: none"> 1 <u>Text field</u> The entered text is displayed here (max. 20 characters). 2 <u>ABC keypad</u> Use this keypad to enter letters. 3 <u>Numeric keypad</u> Use this keypad to enter figures and characters. 4 <u><Shift key</u> Use this key to switch between upper case and lower case letters. 5 <u><Space key</u> Use this key to enter spaces between numbers and characters. 6 <u><<Delete> button</u> Use this key to delete characters on the left of the cursor.

Menu structure



The control panel consists of a display with a touchscreen function. All functions can be interactively controlled via menus. The menu comprises a main menu and two configuration menus. The name of the current user (3) is shown in the lower toolbar, the name of the menu (2) above it. If the menu comprises two pages, press the tabs on the right-hand side (1) to open the page required.



The functions available depend on the configuration of your system. Functions which are not available are displayed as disabled.



Start Screen


Main menu

The main menu permits you to set the major functions:

- | | |
|-----------|-------------------------------------|
| 4 - User | 7 - Balance / AutoBalance (option) |
| 5 - Light | 8 - AutoDrape (option) |
| 6 - Focus | 9 - Record - USB / network (option) |



- All settings made by you in the main menu are not saved. The <Settings> menu on page 2 of the main menu enables you to save all settings on a user-specific or system-specific basis only by pressing the  button.
- Changes to the user-specific camera settings are immediately active, but they are reset to the default values unless they are saved ().



User Settings menu


This submenu displays all functions* that can be stored in the system on a user-specific basis.

System Settings menu

This submenu displays all functions* that can be stored in the system on a system-specific basis.

*The functions available depend on the configuration of your system. Functions which are not available are displayed as disabled.

Overview



User



Press button (1) <User > to access the user submenu (2). The following can be set from there:

- User profile selection
- Assign the language of the user interface to the user profile
- Edit user profile
- Delete user profile
- Add user profile

Selecting user and language




Menu "User" enables you to select a default user profile (default user) and 20 freely-configurable user profiles.

- Press the <User> button (1) in the main menu.
 - The "User" menu (2) is displayed.
 - For multi-user configuration, a list (3) is displayed with a scroll bar.
- Click the arrows to scroll up and down.
- Select the required user (4) from the user list.
- Press the <Language> button and select the language (5) you wish to assign to the current user profile.



You must switch to the main menu to enable the user selected.

- Press the  button to return to the main menu.
 - A message is displayed: The current user is being changed.
 - The main menu with the new user profile (6) is displayed.



Adding, deleting, editing a user

Adding a user

Proceed as follows to add a new user profile:

- Press the <User> button in the main menu.
 - The "User" menu (2) is displayed.
- Press the <Add> button in the "User" menu; a virtual keyboard is shown on the display.
- Enter a user name. The maximum text length is 20 characters.
- Press the  button to save the user.
- Answer the prompt with <Yes> if you want to save the new user.
 - This saves the user.
 - Answer the prompt with <Yes> if you want to apply the settings of the existing user.
 - The "User" menu is displayed again.
- If necessary, select a language for the new user in the "Language" selection box.
- Press the  button to return to the main menu.
 - The main menu is displayed again. The new user profile is displayed on the bottom right.






Deleting a user



- The active user cannot be deleted.
- Select a user other than the one you want to delete from the "User" menu.
- The logged-on user can only be activated after the "User" menu has been exited.
- Re-select the "User" menu.
- Select the user to be deleted.
 - The "Default User" cannot be deleted.
- Press the <Delete> button.
- Answer the prompt with <Yes> if you want to delete the user.
 - The "User" menu is displayed again.
- Press the  button to return to the main menu.



Editing the user profile



CAUTION

Risk of injury!

A wrong configuration of the user profile may lead to unexpected behavior of the equipment.

- Always check the user profile before using the equipment.
- Only make changes to your own profile name.

Press the <Edit> button in menu "User" to display a virtual keyboard with the current user name. Proceed as follows to change the user name:

- Delete the existing name with button .
- Enter a new user name (maximum length 20 characters).
- Press the  button to save the changes.
 - The changed user profile is saved. The "User" menu is displayed again.
- Press the  button to return to the main menu.



Light - Adjusting the light



Setting the light intensity

- Switch the light on (1)  →  .
- Adjust the light intensity (2).
The light intensity can be adjusted using the configurable hand grip buttons (see page 152).
→ The light intensity is displayed as a percentage in a window with a scroll bar (2).
- Click on the arrows to increase or reduce the light intensity.
For focus-light control, please refer to page 140.



Remaining time (3)


The maximum service life of the lamp is 500 hours. The remaining operating hours of the two lamps is displayed in the <Light> menu.

- Change the lamp in good time.

Lamp change



- The manual lamp change in the 180W xenon light source is described on page 172.
- Automatic lamp change is only possible in the optional 300W light source and is described on page 174.



How to adjust illumination settings



The software automatically recognizes the light source (Superlux 180 W or Superlux 300 W) installed in your system and displays the appropriate menu.

Press the <User Settings> button in main menu 2 to access the "Light" sub-menu in which you can make the following settings on a user-specific basis:

1 Speed of light intensity control

The different speeds can be set in several steps.

2 Warning threshold for light intensity

The warning threshold for light intensity can be adjusted in a range from 5 to 100 percent.

An acoustic signal is emitted and an exclamation mark shown on the display behind the light intensity value to indicate that the threshold value has been exceeded.

The exclamation mark remains as long as the light intensity value is higher or equal to the threshold value.

Factory setting: 25 percent

3 Light intensity start value

The light intensity start value can be adjusted in a range from 5 to 99 percent.

Factory setting: 24 percent



The light intensity start value can only be adjusted up to a level 1% lower than the setting for the light warning threshold value. If the light warning threshold value is reduced, the intensity start value is automatically set to 1% below the light warning threshold value.

4 Activating and deactivating Automatic Iris Control


If Automatic Iris Control has been activated, the illumination intensity is automatically adapted to the size of the field of view. This means less exposure of the tissue as a result of too much light (see page 18).



A message draws the user's attention to the consequences of working without Automatic Iris Control if the function is deactivated. The message must be acknowledged with YES or NO.

Factory setting: On

- Press the  button to save the changes.
- Press the  button twice to return to the main menu.





Focus-light control

The focus-light control is a safety feature intended to prevent possible tissue damage due to excessive illumination intensity.

The focus-light control regulates the illumination intensity as a function of the working distance. With increasing working distance, the user has accordingly more light available.

The service technician can set the focus-light control to a value between 5% and 100% in the Service menu.

The focus-light control function can be enabled or disabled by the user as required. It is enabled by default. The factory setting is 25%.

Press the <User Settings> button in main menu 2 to access the "Light" submenu (1). In this menu, the focus-light control can be enabled or disabled by the user as required.

Press the Focus-light link button (2) to disable the function.



- Press button  Disable? in the Info window (3) and confirm with <YES> (4).
- Press button  (5) in the "Light" submenu to save your settings.
- Press the  button twice to return to the main menu.

The set maximum illumination intensity is shown by a blue bar (6) in the light controller (e.g., max. 24%).

The illumination intensity is reduced to the set limit value (6). It can be set within the range below the blue bar, but not above it.



The service technician can configure the focus-light control in the Service menu such that it cannot be disabled.



Configuring the video focus function (Options: SpeedFocus, video focus)



- 1 Automatic video focus (optional)**
The "automatic video focus" function ("autofocus") enables you to automatically focus the video image on the monitor.
- 2 Manual video focus (optional)**
This setting permits you to adjust the focus of the image generated by the camera on the monitor (optional).
This focus function on the monitor is independent of the focus of the surgical field actually seen by the surgeon through the microscope.
<Video Focus Adjustment> is only possible if a camera option has been enabled.
- 3 Focus Lock**
The <Focus Lock> function permits you to deactivate the electrical drive of the focusing system. This ensures that the set focal planes are not inadvertently changed when a connected laser micromanipulator is used. The focus rocker switches on the hand grips and the focusing buttons on the foot control panel (optional) are disabled. If Focus Lock has been activated, no autofocus setting is performed when the brakes are operated, even if SpeedFocus (optional) has been switched on.



Configuring the focusing function



Press the <User Settings> button in main menu 2 to access the "Focus" sub-menu. In this submenu, you can make the following settings on a user specific basis:

1 Setting the focus start value

Enter the working distance at which the focusing mechanism starts after the system is powered on.

The working distance can be adjusted and saved in a range from 200 mm to 500 mm.

Factory setting: 300 mm

2 Setting the focusing speed

The focusing speed can be set from minimum to maximum in a range from 5% to 100%.

Factory setting: 50%



Please note that the setting of the focusing speed does not change linearly between 30% and 100%. The consequence is that focusing can still be deemed too quick at 30% of the maximum speed.

In this case, please use speeds slower than 30% or perform fine focusing using the knob on the microscope body.

3 Starting autofocus automatically

(only possible with SpeedFocus option, otherwise the function is disabled)
If autofocus is "On", it automatically focuses on the object (within the working distance) when the brakes are applied.

Factory setting: Off



SpeedFocus (autofocus) option

The SpeedFocus (autofocus) option permits you to select and automatically focus on a feature in the ROI (region of interest, detailed view) in the image of the surgical field.

SpeedFocus (autofocus) is based on the analysis of the video image contrast: Whilst the varioscope lens of the surgery microscope is focusing, the SpeedFocus continually determines the contrast of the video image in the detailed view (ROI) selected. In this way, SpeedFocus (autofocus) controls the varioscope lens of the microscope and adjusts it to the point of maximum contrast which is also the point of maximum image definition.

4 Selecting depth of field

The depth of field type can be set as follows: small / large

- Large depth of field → small aperture → less light
- Small depth of field → large aperture → more light

Factory setting: Small depth of field


5 SpeedFocus ROI (region of interest, detailed view)

(only possible with SpeedFocus (autofocus) option, otherwise the function is disabled)


Set the size of the detailed view to be displayed after focusing.

Settings: small/medium/large/off

Factory setting: large



- Press the  button to save the changes.
→ The setting is saved for the current user.
- Press the  button twice to return to the main menu.



Setup of record function


(Integrated video record function over USB port - option)



The "Record" submenu enables you to record video sequences and single images using the optional "video recording function over USB port".

The recorded video and image data can be stored on an external USB storage medium or in a shared directory within a network. The unambiguous and automatic storage of images and videos ensures reliable assignment to the relevant patient.


For a detailed description of the video record function, please refer to Options annex "Integrated Video Record Function over USB Port" (see page 225).



Perform white balance (option)



(Only if the video camera is installed (optional))





The <White balance> button on page two of the main menu allows you to perform a white balance. The system will adjust the video signal in such a way that white areas in the surgical field are also white on the monitor. In order to proceed with a white balance, follow these steps:

- Press the <White balance> button.
 - The following message will be displayed: Please place a white sheet of paper under the microscope, focus and then press "Continue".
- Align the surgical microscope towards a white object (e.g., a sheet of white paper with a matte surface).
- Focus on the object.
- Press the <Continue> button.
 - The following message is displayed: "Please wait - white balancing in progress!"
 - If the white balancing was successful, the following message will be displayed: "White balancing was successful". If the process was unsuccessful, the following message will be displayed: "White balancing failed".
- If white balancing failed, verify the settings described above and repeat the white balancing procedure.



Please note that a white balance can result in an error when the maximum light intensity is used.

- Please reduce the light intensity or increase the working distance and restart the white balance.
- Press the  button to return to the previous menu.



User settings menu 1

Configuring zoom



Pressing the <User settings> button and then the <Zoom> button accesses submenu "Zoom", from where you are able to configure the following "user-specific" settings:

1 Set initial zoom value

Enter the zoom factor (0.4x... 2.4x) the magnification mechanism starts with when the machine is powered on.

Factory setting: 0.4x

2 Zoom speed

Set the speed for moving the zoom

Setting of the zoom speed from a minimum to a maximum speed (5% to 100%).

Factory setting: 50%


3 Focus zoom link

The higher the magnification (zoom), the lower the focus speed. If the focus zoom link function is enabled, the move speed of the optional XY coupling is controlled dependent on zoom.

This zoom-dependent speed adjustment can be enabled/disabled.

Factory setting: Off

- Press the  button to save changes.
→ The user-specific setting is saved.
- Press the  button twice to return to the main menu.



Configuring the handgrips



Select <User Settings> and press the <Handgrips> button to access the "Handgrips" submenu. This menu permits you to assign different functions to the handgrips on a user-specific basis.

You can configure both handgrips with the same functions by pressing button (2), or you can configure different functions in the right and left handgrip (3, 1).

Differences between the right and left handgrip buttons:

- Light- and video focus- are only available on the left handgrip button
- Light+ and video focus+ are only available on the right handgrip button
- Light+/- and video focus+/- can only be assigned in pairs (A-B or C-D)
- If light+/- or video focus+/- are not assigned in pairs, the second assignment will have no function.

Changing the button assignment

- Press the button whose function you want to change:
 - A list of the selectable functions is displayed. If several selectable functions are available, a scroll bar is displayed. Use the arrow buttons to scroll up and down.
- Tap on the function to be assigned to the button.
- Repeat these two steps until you have assigned the required functions to all configurable handgrip buttons.
- Press button (2) if you want to assign the same functions to both handgrips.



Please note: if you press button (2), the left handgrip automatically adopts the configuration of the right handgrip.

- Successively press buttons (1 and 3) if you want to assign different functions to the right and left handgrips.
- Press the  button to save the changes.
 - The setting is saved for the current user.
- Press the  button twice to return to the main menu.

Joystick The joystick is used to control the XY coupling (option).

Rocker switches The focus/zoom setting of the rocker switches is preconfigured but it can be

changed at any time. User-specific settings can be made by yourself at any time. Press button (4) for this purpose.

Factory settings of buttons A, B, C and D:

A: DOF aperture (depth of field) low/high

B: light on/off


C/D: light +/-

Factory settings of buttons E, F, G and H:	User-specific setting options:
<p>Button E: focus rocker switch +</p> <p>Button G: zoom rocker switch +</p> <p>Button F: focus rocker switch -</p> <p>Button H: zoom rocker switch -</p>	<p>Button I: zoom rocker switch +</p> <p>Button G: zoom rocker switch -</p> <p>Button J: focus rocker switch +</p> <p>Button H: focus rocker switch -</p>

Functions of configurable hand grip buttons

- Left side (buttons A and C):
- D.o.F small/large
 - Light on/off
 - SpeedFocus
 - Video SpeedFocus
 - Video rec. on/off
 - Photo (image acquisition with integrated video camera)
 - Light -
 - VideoFocus -
 - Aux on/off
 - No function

- Right side (buttons B and D):
- D.o.F small/large
 - Light on/off
 - SpeedFocus
 - Video SpeedFocus
 - Video rec. on/off
 - Photo (image acquisition with integrated video camera)
 - Light +
 - VideoFocus +
 - Aux on/off
 - No function



Configuring the foot control panel



Press button <User Settings> and then button <Fußschaltpult> to access the "Foot Control Panel" submenu. This menu allows you to assign different functions to the connected foot control panel on a user-specific basis.

Differences between the right and left buttons:

- Light - and video focus - are only available on the left button
- Light + and video focus + are only available on the right button
- Light +/- and video focus +/- can only be assigned in pairs (A-B, C-D or E-F)
- If light +/- or video focus +/- are not assigned in pairs, the second assignment will have no function.

Changing the button assignment

- Press the button whose function you want to change:
 - A list of the selectable functions is displayed. If several selectable functions are available, a scroll bar is displayed. Use the arrow buttons to scroll up and down.
- Tap on the function to be assigned to the button.
- Repeat the two steps above until you have assigned the desired function to all of the buttons.
- Press the  button to save the changes.
 - The setting is saved for the current user.
- Press the  button twice to return to the main menu.



Please note that the E and F buttons are subject to some limitations if an earlier-generation foot switch is used. This means that functions requiring continuous adjustment can be adjusted in incremental steps only.

- Preferably assign functions which feature only 2 states to these buttons, such as "aperture for depth of field - large or small".




Functions of the configurable foot control panel buttons

Left side (buttons A, C and E):

- D.o.F small/large
- Light on/off
- SpeedFocus
- Video SpeedFocus
- Video rec. on/off
- Photo (image acquisition with integrated video camera)
- Light -
- VideoFocus -
- Aux on/off
- No function

Right side (buttons B, D and F):

- D.o.F small/large
- Light on/off
- SpeedFocus
- Video SpeedFocus
- Video rec. on/off
- Photo (image acquisition with integrated video camera)
- Light +
- VideoFocus +
- Aux on/off
- No function



Preconfigured button assignment




Before using the wireless foot -operated control panel FCP WL ensure its batteries are fully charged.

Buttons

The functions of buttons A, B, C, D, E and F are preconfigured, but they can be re-assigned to meet the user's specific requirements. The factory settings are:

- Button A: DOF aperture -> depth of field low/high
- Button B: light on/off
- Button C: light -
- Button D: light +
- Button E: without function
- Button F: without function




Joystick

Joystick (1) is used to control the XY coupling (option).

Rocker switches


The focus/zoom rocker switches are preconfigured to the horizontal setting. However, they can be changed to the vertical setting at any time. User-specific settings (see page 154) can be made by yourself at any time. Press button (2) on the display for this purpose.

The horizontal factory settings are:



- Button I: focus rocker switch +
- Button J: focus rocker switch -
- Button G: zoom rocker switch +
- Button H: zoom rocker switch -

The vertical settings are:




- Button I: zoom rocker switch -
- Button J: focus rocker switch -
- Button G: zoom rocker switch +
- Button H: focus rocker switch +



Always check the assignment of functions and the functions of the foot control panel before every use.

Configuring the brakes

The brake release buttons are located on the back of the handgrips.



1 AB brake release button (All Brakes)

Brake release button for all stand **and** microscope axes. For as long as you press this button, **all** magnetic brakes are released and the system can be moved in all directions. The magnetic brakes are locked when you release the button.

2 SB brake release button (Selected Brakes)

Brake release button for the microscope axes **or** stand axes. For as long as you press this button, only the magnetic brakes of the stand **or** surgical microscope are released. The magnetic brakes are locked when you release the button.

Select the <User Settings> button and press the <Brakes> button to access the "Brakes" submenu in which you can make the following settings on a user-specific basis:

You can select whether the SB button should release the microscope or the stand axes.

- Select the axes which you want to be released when you press SB button (2) on the handgrip.
 - Stand axes (factory setting)
 - Microscope axes
- Press the  button to save the changes.
- Press the  button twice to return to the main menu.



User settings menu 2

Configuring the XY coupling (optional)




Pressing the <User settings, page 2> button and then <XY Coupling> accesses submenu "XY coupling", from where you are able to set the "user-specific" speed of the XY coupling.

The speed can be set in range 5 to 100%.

Factory setting: 50%

- Press the  button to save changes.
→ The user-specific setting is saved.
- Press the  button twice to return to the main menu.



Camera setting 1/2 (option)




Press button <User Settings, page 2> and then button <Camera> to access the "Camera" submenu. This menu allows you to set the camera parameters on a user-specific basis.

Setting the exposure mode

Select between manual exposure control (button for manual brightness adjustment is activated) and automatic exposure control (button for manual brightness adjustment is deactivated).

Automatic exposure




Automatic exposure control is the preferred method in most applications. In this mode, the camera automatically adjusts the brightness of the video image to the specified value. Exposure control affects the exposure time of the camera. If automatic exposure control is selected, the nominal brightness is set and the camera controls the exposure time according to this nominal value. The nominal brightness can be varied in the range between 0% and 100%.

- Deactivate the <Manual Brightness Setting> button (1).
 - The display of the brightness setting range changes from exposure time (seconds) to percent (%).


Select the metering method to be used in the "Spot Size" field (2).



- Full Field (default)
The exposure is measured and averaged across the full video image. This metering pattern is recommended for surgical fields that are fully and evenly illuminated.




- Large Spot
The exposure is measured in an area in the image center. This metering pattern is ideal for working with an eclipsed surgical field edge (resulting from a reduced illuminated field diameter).



- Small Spot
The exposure is measured in a very small area in the image center. This metering pattern is suitable for working with an extremely small illuminated field diameter.
- However, if the object of interest is not located at the image center, the "Small Spot" setting usually does not provide the desired exposure result. In this case, select either a larger area (Full Field or Large Spot) or use the manual exposure mode.

Manual exposure

Manual exposure control is suitable for difficult light conditions. In this mode, a fixed exposure time is set.



- Activate the <Manual Brightness Setting> button (1).
 - Manual exposure control is active when the "Man. Brightness Setting" button emits dark blue light.
 - The display of the brightness setting range changes from percent (%) to exposure time (seconds).
 - The button under the "Spot Size" field is disabled.


Move the slider to the required exposure value. You can adjust the exposure time between 1/10000 s und 1/8 s.

- Press the  button to save the changes.
 - The setting is saved for the current user.
- Press the  button twice to return to the main menu.

Setting the HD DigiZoom (3)

The HD DigiZoom function allows digital enlargement, but only if an integrated HD camera is installed.

- Press the - / + buttons to zoom out and in, respectively. The current value is displayed in the middle. When the highest or lowest value has been reached, the respective arrow button will be grayed out.



Camera settings 2/2 - Color setting (option)

You are able to configure the following user-specific camera settings on page 2 of submenu "Camera":

1 Chroma (color saturation) of a 3CCD camera or HD camera

Setting range: -10 to +10

Factory setting: 0

2 Hue setting of a 3CCD camera or HD camera

Setting range: -10 to +10

Factory setting: 0



Only the integrated 3 CCD cameras and HD cameras supplied by Carl Zeiss can be adjusted.

- Press the  button to save changes.
→ The user-specific setting is saved.
- Press the  button twice to return to the main menu.



Configuring video source (camera) - only for SD camera option




Select <User Settings, page 2> and press the <Video source> button to access the "Video source" submenu.

This configuration menu permits you to select the video source. It is possible to show the video image of the internal or external video camera (e.g. an endoscope camera) on the (optional) 17" video monitor and record it with the (optional) video recorder. The video signal of the external video camera is sent to the system via the optional ENDO connecting cable (video input and output port, orange).

Settings: integrated / external

Factory setting: integrated

- Press the  button to save the changes.
→ The setting is saved for the current user.
- Press the  button twice to return to the main menu.



System settings menu 1

Optics


Press the <System Settings> button in the main menu to access the "System Settings" submenu. In the "Optics" submenu, the following settings can be assigned to the system:

To ensure correct display of the total magnification, the parameters of the objective lens, tube and eyepieces used must be set here. The settings are applied for all users.

- Select the system components used by pressing the relevant buttons in the "Eyepiece" (2) and "Tube" (1) fields.
 - The parameter is activated when the button is dark blue.
- Press the  button to save the changes.
 - The setting is saved for the current user.
- Press the  button twice to return to the main menu.



When using the folding tube f170/f260 with integrated, manual magnification changer, the overall magnification is not automatically updated for switchover on the folding tube.



Configuring the video settings (optional)

Press the <System Settings> button in the main menu to access the "System Settings" submenu. The "Camera" submenu permits you to make the following camera settings:

- Performing the white balance procedure (see page 148)
- Setting the HD format (only active if a HD video camera is connected)
- Setting the SDTV format (only active if a SD video camera is connected)



Setting the HD format

If a high-resolution display device (e.g. HDTV monitor) is connected, the "HD Format" field permits you to choose between the "720p" and "1080p" resolution.

The maximum resolution of the HD-SDI interface is 1080p. The 1080p signal is only available on the DVI-D port.

- Select the 720p format if your display device is 720p compatible. This mode is similar to the 1080p format but its maximum resolution is only 1280 x 720 pixels.
- Select the 1080p format if your display device is 1080p compatible. The image is set up using real full frames. This prevents any interline flicker, and the image is sharply defined and stable. The maximum resolution is 1920 x 1080 pixels.

Proceed as follows to select a format and to load the respective configuration:

- Press the <1080p> or <720p> button to set the video resolution.
 - The video resolution will be adjusted.
As this process takes some time, the change is not enabled until you press the Save button.
- Press the  button to save the changes.
 - The settings are saved on a system-specific basis.
- Press the  button twice to return to the main menu.
- The main menu is displayed again.

How to set the SDTV format

If a conventional display device (monitor, projector, etc.) with the 4:3 format is connected, the "SDTV Format" field enables you to choose between "Letterbox" and "Sidecut" visualization.

- Use the "Letterbox" format for display devices that are not compatible with the 16:9 widescreen format. A black bar appears at the top and bottom of the screen and the image is displayed in the 16:9 format.
- In the "Sidecut" format, the sides of the image are truncated, leaving the center of the image displayed with a 4:3 aspect ratio.

Proceed as follows to select a format and to load the right configuration:

- Press the currently selected format in the "SDTV Format" field.
 - This shows a list of available formats.
- Select the format to be loaded.
 - Press the  button to save changes.
Because this switchover process takes some time, the setting does not become active until the Save button is pressed.
 - The system-specific settings are saved.
- Press the  button twice to return to the main menu.

The main menu is displayed again.

Adjusting the brightness of the additional display

Select <System Settings> in the main menu and press the <Add. Display> button to access the "Add. Display" submenu in which you can adjust the brightness of the additional display.

The brightness can be adjusted in a range from 5 to 100 percent.

Factory setting: 50%


- Press the  button to save the changes.
→ The setting is saved on a system-specific basis.
- Press the  button twice to return to the main menu.



System settings menu 2

Pressing the <System settings> button in main menu 2 accesses 2-page sub-menu "System settings". Page 2 of this submenu offers the following:

- 1 Service PIN (password protected)
This PIN is only available for persons who have been trained by Carl Zeiss Surgical.
- 2 Export of service files
Send Carl Zeiss Service these files in the event of a fault. The files enable faster troubleshooting.
- 3 Configuration
Menu for displaying the configuration of the device



Procedure in the event of faults



NOTE

Malfunction

If a failure occurs which you cannot rectify using Section "Procedure in the event of faults", label the system as non-functional and inform the Carl Zeiss service department.



CAUTION

Error messages on the display

Watch out for error messages on the display.

If faults occur during operation, an error message is shown on the display and a warning signal is sounded. In this case a warning triangle is also shown on the additional display.

Manual lamp change in the 180W xenon light source..... 172

**Automatic lamp change
300-W xenon light (optional)..... 174**

**Manual lamp change
300-W xenon light (optional)..... 176**

Electrical malfunction (focus, zoom, lighting) 178

Failure of the display 180

Troubleshooting 181

For your safety 181

Malfunctions of the system..... 182

Malfunctions of the wireless foot control panel 184

Malfunctions of the cable-based foot control panel 185

Manual lamp change in the 180W xenon light source



CAUTION

Lamp failure

- The service life of the xenon lamp is limited to 500 hours. If the service life of the lamp is exceeded, the xenon lamp may suddenly fail.
- With increased aging of the light source, the actual illumination intensity at a particular setting decreases.
- Please replace the xenon lamp in good time.
Note the remaining hours display on the screen.


NOTE

Lamp replacement

- A lamp can be replaced within one minute.

Switching to the second lamp


The lamp module contains two xenon lamps. The second lamp can be swung into position manually should the first lamp fail. The lamp module can be opened up if the first xenon lamp fails.



- Before replacing the lamp, switch off the stand from the power switch.
- Press button (2) to open the lamp module.
→ The lamp module is ejected slightly.
- Pull out the lamp module as far as it will go.
- Swing the second xenon lamp by turning knob (1) through 180° until it locks into place.
→ The segment (3) in the knob (1) lights red.
- Push the lamp module all the way back into the lamp housing.
- Switch the stand back on from the power switch.



If the first lamp failed and the second lamp is in use, a warning message is shown on the display on every restart of the device reminding you to keep a new lamp module ready.



Automatic lamp change 300-W xenon light (optional)



CAUTION

Lamp failure

- The service life of the xenon lamp is limited to 500 hours. If the service life of the lamp is exceeded, the xenon lamp may suddenly fail.
- With increased aging of the light source, the actual illumination intensity at a particular setting decreases.
- Please replace the xenon lamp in good time.
Note the remaining hours display on the screen.

NOTE

Lamp replacement

- A lamp can be replaced within one minute.

Switching to the second lamp

The lamp module contains two xenon lamps. In the event of a failure of the active lamp, the system automatically changes to the other lamp.


However, lamp change can also be triggered via the control panel.

- Press button <Light> (1) in the main menu to display submenu "Light".
- Press the <Lamp replacement> (2) button.
- Press the  button twice to return to the main menu.

This should be done if the maximum number of operating hours of the active lamp has been reached and the operator wants to change to the other lamp as a precaution.



- If the first lamp failed and the second lamp is in use, a warning message is displayed reminding the user to keep a lamp module ready as a replacement.
- The operating hours are stored separately in every lamp. The operating hours remaining (3) are displayed when using a used lamp.



Manual lamp change 300-W xenon light (optional)

If automatic lamp change fails, the lamp must be changed manually.

The lamp housing contains two xenon lamps. The slider (2) for manually changing the lamp is always in front of the lamp just having been used. Proceed as follows in the event of a lamp failure and failure of the automatic lamp change function:

- Turn off the suspension system at the power switch.
- Move the fastener in the direction of the arrow to open the housing door (1).
- Push or press the slider (2) all the way in the direction of the second lamp (no longer in use).
- Re-close the housing door (1).
- Switch the suspension system back on again at the power switch.




Be sure to replace the used lamp by a new one after an operation.

- Dispose of the used lamp in line with local guidelines and laws.

NOTE**Lamp replacement**

- A lamp can be replaced within one minute.



Electrical malfunction (focus, zoom, lighting)

NOTE**Failure of major functions!**

A failure of the electronics may affect the focus, zoom and lighting functions. Press the orange-colored button (1) to activate Light-only mode and, e. g., finish an operation that has already been started. Light-only mode ensures continuing illumination at constant and sufficient light intensity.



Working in Light-only mode

In Light-only mode, the lighting is set to an intensity of ca 70% and the power to all motorized functions is switched off (including magnetic brakes).

- Manually move the microscope in the X and Y directions to position it over the surgical field.
- Manually move the microscope up and down to focus the image of the surgical field.
- Use knob (2) to adjust the zoom value (magnification).
- Use knob (3) to adjust the focus (working distance).
- Use knob (4) to adjust the luminous field diameter.

NOTE**Failure of the electronics (magnetic brakes)**

If the magnetic brakes fail (magnetic brakes are closed), you can manually move the axes of the stand by overcoming the locking effect of the brakes.



Failure of the display

**CAUTION****Do not touch the display if it has failed**

- On no account touch the display surface if the display has failed.

Only the lighting of the display may have failed. If so, you might access unwanted menus, trigger functions or unintentionally change values.



If the display no longer responds to input, it is still possible to work with the device with the existing settings.

Proceed as follows to restore the function of the display:

- Switch the stand off and then (after ca. 2 minutes) on again.
After being switched on again, the system automatically enters normal operating mode after the self-test.
- Check all the functions of the surgical microscope and the stand.

You can continue to perform operations in normal mode.

- Please notify your Carl Zeiss Service department should you experience any faults.

Troubleshooting

For your safety

**CAUTION****Error messages on the display!**

Pay attention to the error messages on the display.

In the event of malfunctions occurring during operation, error messages are shown on the display. The additional display shows a warning triangle in this case.

- Check and eliminate the malfunction.
- Acknowledge the message by pressing the relevant button on the display.

Note**Problems in troubleshooting!**

If you have problems eliminating the malfunction or if errors keep recurring, do not continue using the system, attach a sign to the system stating that it is out of order, and contact your Carl Zeiss service partner.

Export of service files

In the event of malfunctions, send these files to your service technician. The files enable faster troubleshooting.

You will find the contact responsible for your country on the following website: www.meditec.zeiss.com

Malfunctions of the system

Problem	Possible cause	Remedy	See
No function at all.	Power plug of suspension system not inserted.	Plug in the power cord.	page 84
	Power switch of suspension system not switched on.	Press the power switch. Green indicator light in the power switch must be lit.	page 118
	Automatic circuit breaker in power switch of suspension system has been activated.	Press power switch again.	page 118
	Power failure.	Contact in-house electrician.	-
Surgical field illumination on microscope not working.			-
	Failure of lamp 1 and lamp 2.	Replace both lamps.	-
	Failure of suspension system electronics.	Illuminate surgical field using an OR illuminator. Contact service dept.	-
	Lamp module in suspension system has no contact.	Insert lamp module as far as it will go.	-
	Light source has not been switched on.	Press light source on/off button on the handgrip / foot control panel.	-
Error message on the display: Lamp overheated - device switches off	Ventilation slots of light source are covered.	Uncover ventilation slots of light source.	-
Insufficient surgical field illumination.	Brightness level set too low.	Adjust brightness on suspension system or handgrip / foot control panel.	-
	Lamp ageing may lead to reduced light intensity.	Replace lamp.	-
	Defective light guide (illumination not uniform).	Replace light guide. Contact service dept.	-

Problem	Possible cause	Remedy	See
Surgical field illumination too bright.	Brightness level set too high.	Adjust brightness on suspension system or handgrip / foot control panel.	-
		Switch off light source on suspension system. Illuminate surgical field using an OR illuminator. Contact service dept.	-
Image on video monitor too dark, or high level of noise.	Defective light control.	Change to the Light Only mode.	-
	Use a camera adapter.	Remove fixed stop from camera adapter. Increase the brightness level of illumination.	User manual for camera -
Motorized focusing and /or zoom function of surgical microscope is inoperative.	Failure of suspension system electronics.	Manually adjust focusing and / or zoom on surgical microscope.	-
No video image.	Connection cable not properly connected.	Check all connections.	-
Stand wobbles.	Floor not level. Stand base not appropriately positioned.	Slightly change the stand base position.	-

Malfunctions of the wireless foot control panel

Problem	Possible cause	Remedy	See
Intermittent malfunction	Batteries are flat	Replace batteries for new ones	User manual for the foot control panel
	Use of rechargeable batteries (accumulators)	Replace batteries for new ones	User manual for the foot control panel
	Failure of individual button functions	Configure button function differently - only possible if functions on the stand are configurable	page 154
	Failure/fault of the radio link for the wireless foot control panel	Establish cable connection (if connector cable present).	User manual for the foot control panel
	Position switch always detects rest position	Establish cable connection (if connector cable present).	User manual for the foot control panel
	Radio link faulty	Establish cable connection (if connector cable present).	User manual for the foot control panel
	Weak radio signal	Establish cable connection (if connector cable present).	User manual for the foot control panel
	No pairing with stand	Perform pairing with the stand	page 88

Problem	Possible cause	Remedy	See
Unintentional triggering of function	Control jams mechanically after being pressed	Move foot switch panel to the rest position	User manual for the foot control panel
		Configure button function differently - only possible if functions on the stand are configurable.	
	Foot switch panel sends erroneous activation signal	Move foot switch panel to the rest position	User manual for the foot control panel
		Configure button function differently - only possible if functions on the stand are configurable	
Swap foot control panels		Check that the label on the stand and on the indicator of the foot control panel match	User manual for the foot control panel
		Perform pairing with the stand	page 88

Malfunctions of the cable-based foot control panel

Malfunction	Possible cause	Remedy	Refer to
Foot control panel without any functions.	Suspension system or foot control panel not plugged in.	Establish connection.	page 86
Sporadic loss of functions.	Loss of individual key functions.	Reconfigure key functions - only possible if functions on suspension system can be configured.	page 154

Care and Maintenance



System maintenance	188
Changing the Superlux 180W xenon lamp module	188
Exchanging xenon lamp for 300W Superlux (option)	190
Counterbalancing the monitor support arm	192
Maintenance intervals.....	194
Technical safety check	198
Care of the device	199
Cleaning.....	199
Sterile use (drapes)	200
Disinfection	201
Environmental protection measures	202
Note on disposal	202

System maintenance

Changing the Superlux 180W xenon lamp module



CAUTION

Risk of injury caused by improper handling of the xenon lamp!

Improper handling of the xenon lamp may lead to damage or injury.

- The lamp module must only be changed by appropriately trained personnel.




CAUTION



Hazard caused by hot lamps!

In the event of malfunction, there is the risk that the lamp may burst due to the high pressure inside the hot lamp. The hot surface of the xenon lamp may also cause burns.

- Only change the lamp module after it has cooled down (this takes approx. 10 min)!




Replacing the lamp module

- Before replacing the lamp, turn off the device at the power switch.
- Press button (1) to open the lamp module.
→ The lamp module is slightly ejected.
- Pull out the lamp module as far as it will go.
- Slide the original transport case (2) over the module, making sure that bolt (3) engages in bore (4). This unlocks the stop.
- Remove the old lamp module and insert the new module by proceeding in the reverse order.
- Turn on the system at its power switch.
- Check the function of the two xenon lamps.
- Pack the old lamp module (5) in the transport package of the new lamp module.
- Fill in the enclosed return card and send the old lamp module to the nearest Zeiss service representative.



Only use the original transport case (2) as it also provides explosion protection should xenon lamps be defective.



Exchanging xenon lamp for 300W Superlux (option)



CAUTION

Risk of injury caused by improper handling of the xenon lamp!

Improper handling of the xenon lamp may lead to damage or injury.

- The lamp module must only be changed by appropriately trained personnel.



CAUTION



Hazard caused by hot lamps!

In the event of malfunction, there is the risk that the lamp may burst due to the high pressure inside the hot lamp. The hot surface of the xenon lamp may also cause burns.

- Only change the lamp module after it has cooled down (this takes approx. 10 min)!






Replacing the 300W lamp

The lamp housing contains two xenon lamps. If a lamp fails, you can open the lamp module to replace the lamp.

- Before replacing the lamp, switch off the suspension system at the power switch.
- Push the lock in the direction of the arrow (1) and open the door (2).
- Pull out the defective lamp (3).
- Push a new lamp all the way into the lamp housing (4).
- Close the door again (5).
- Switch the suspension system back on again at the power switch.



If the first lamp failed and the second lamp is in use, a warning message is shown on the display, reminding you to keep a backup lamp module ready at hand.



Counterbalancing the monitor support arm

If the video monitor does not remain in place in the position required, the following readjustments can be performed.

Increasing the friction of left/right movement of the carrier arm.

- Slightly tighten screw (1) of the carrier arm by turning it clockwise using an M5 hex key.

Increasing the friction of left/right movement of the spring arm.

- Remove the plastic cover (2) on the spring arm joint.
- Loosen securing screw (4) on the support arm by turning it counterclockwise using an M2.5 hex key.
- Tighten the adjustment screw (3) of the spring arm by turning it clockwise until the required friction is attained.
- Firmly retighten the securing screw (4) on the support arm by turning it clockwise using an M2.5 hex key.
- Reattach the plastic cover (2).





Readjusting the gas pressure spring

- For correct adjustment of the gas pressure spring, align the spring and support arms in a horizontal position.
- Tighten setting screw (5) of the gas pressure spring by turning it **counterclockwise** with an M5 hex key until the spring arm with the video monitor no longer moves downwards.
- If the spring arm with video monitor continues to move downward, the gas pressure spring is defective.
- Please inform the Carl Zeiss service department in this case.



Adjusting the movement of the video monitor

- Tighten the securing screw (6) of the monitor bracket by turning it clockwise until the video monitor remains in the required position.



Maintenance intervals

Regular preventive maintenance is required to ensure safe and correct operation and a long service life of the system on an ongoing basis.

The following table lists the relevant maintenance activities and shows at which maximum intervals they must be performed.

For maintenance activities that can be performed by yourself, the relevant procedure is described in the following sections.

All maintenance work not detailed here requires system-specific expert know-how. It is frequently necessary to open the system for this work.

Please contact your local Carl Zeiss service representative in due time for the performance of these maintenance activities.

Component	Inspection features and acceptance criteria	Comments
According to country-specific interval requirements		
Safety check		
Before each use		
Brakes	Electrical and mechanical function test: unlock brakes completely; braking force holds microscope securely in balanced state without any change of position.	
Lamps	Checking: <ul style="list-style-type: none"> – Operating hours within permissible range – Replacement available and ready for use 	
When specified operating hours have been reached		
Lamps	We recommend replacing the light source when the specified operating hours have been reached.	
Every 12 months		
Documentation and identification labels	Visual inspection: Manuals and identification labels available, undamaged, complete and legible	

Component	Inspection features and acceptance criteria	Comments
Controls such as hand-grips and XY coupling (option)	Visual inspection and function test: <ul style="list-style-type: none"> – Optical and acoustic signals functional and react to prompt – Light Only mode is operational – Probes move easily without play 	
Carrier system including microscope mount	Function test: <ul style="list-style-type: none"> – play-free movement of the bearings, stops and suspensions – brakes operational at maximum permissible OPMI configuration <p style="text-align: center;">Maximum admissible load on the microscope body</p> <div style="text-align: right;">  </div>	
Stand column	Function test: Secure without mechanical play in the column connection	
Casters	Visual inspection and function test: <ul style="list-style-type: none"> – play-free attachment, lock-up free turning and rolling, locking tabs and brakes securely positioned – Cable deflectors available 	
Illumination: Fan	Function test: <ul style="list-style-type: none"> – Fan and vacuum capacity noticeable 	
Microscope	Function test: Checking optical imaging and illumination of the field of view	
Zoom, focus and illuminated field diaphragm	Function test: <ul style="list-style-type: none"> – Jolt-free operation of the zoom system – Play-free adjustability of focus, zoom and illuminated field diaphragm when adjusted manually 	
Objective lens, tube and other accessories mounted to the microscope (camera, etc.):	Function test: Attached parts without mechanical play, knurled screws tightened and other removable microscope components securely mounted	
Foot control panel (option)	Visual inspection and function test: <ul style="list-style-type: none"> – Complete, undamaged dust and spray water protection – Replace FCP WL battery 	

Component	Inspection features and acceptance criteria	Comments
Monitor (option)	Visual inspection and function test: <ul style="list-style-type: none"> – No noticeable damage to housing – Checking the image quality – Checking the gas pressure spring of the lift arm 	
Video (option)	Function test: SpeedFocus	

Every 2 years

Illumination: Light guide	Function test: <ul style="list-style-type: none"> – Checking the light guide – Checking the illuminated field 	
Illumination: Optical filters	<ul style="list-style-type: none"> – Heat-absorbing filter and UV filter: – checking the light source filters 	
AutoDrape (option)	Function test: drape pump function	

Every 4 years

Carrier system including microscope mount	Visual inspection: <ul style="list-style-type: none"> – no signs of wear on the supporting structure – No signs of wear or damage on the interfaces and couplings – No visible damage or wear on the connection and mount of the spring packs 	
Cables, plugs and switches	Visual inspection: <ul style="list-style-type: none"> – Undamaged insulation and connectors on system/cable – Checking the cable guide including fasteners 	
Elektronik	We recommend replacing the internal batteries	
Casters	We recommend replacing the casters	
Light guide	We recommend replacing the light guide	

Component	Inspection features and acceptance criteria	Comments
Every 6 years		
Monitor	Checking the gas pressure spring of the lift arm We recommend replacing the monitor carrier arm including the gas pressure spring.	
Navigation	Function test of the interface	By manufacturer of the navigation system
Every 8 years		
AutoBalance (option)	Replacement of wearing parts in the balancing drive	

Technical safety check

**CAUTION****Risk of injury!**

- Make sure that the regular technical safety checks required for this system in accordance with the applicable national regulations are performed on schedule and to the stipulated extent.

To prevent any impairment of the system's safety as a result of ageing, wear, etc., the organization operating the system must ensure, in accordance with the applicable national regulations, that the regular technical safety checks defined for this system are performed on schedule and to the stipulated extent.

The safety checks must only be performed by the manufacturer or qualified personnel.

At a minimum, the scope of the safety checks of the system should comprise the following points:

- Availability of the user manual
- Visual inspection of the system and accessories for damage and legibility of the symbols and labels
- Leakage current test (measuring point: potential equalization bolt)
- Test of the protective ground conductor (measuring point: potential equalization bolt)
- Function and wear test of the steerable casters and locking tabs
- Function test of all switches, buttons, sockets and indicator lamps of the system
- Function test of the Light Only mode without a patient.
- The drive unit of the AutoBalance function (option) is subject to natural wear.
Check the age of the drive unit and replace it if it is older than 8 years.

Care of the device

Cleaning

NOTE**Proper cleaning of the device**

- Follow the cleaning methods below:

Cleaning optical surfaces

The multi-layer T* coating of the optical components (e.g. eyepieces, objective lenses) ensures optimum image quality.

Image quality is impaired by even slight contamination of the optics or by a fingerprint. To protect the internal optics from dust, the system should never be left without the objective lens, binocular tube and eyepieces. After use, cover the system to protect it from dust. Always store objective lenses, eyepieces and accessories in dust-free cases when they are not being used.

Clean the exterior surfaces of the optical components (eyepieces, objective lenses) only when necessary:

- Do not use any chemical cleaning agents.
- Remove dust from the optical surfaces using a squeeze blower or a clean, grease-free brush.

For the regular cleaning of objective lenses and eyepieces of the surgical microscope, we recommend the optics cleaning set available from ZEISS. For the catalog number, please see the section "System data - Ordering data".

Prevention of fogging

To protect the eyepiece optics from fogging, we recommend using an anti-fogging agent. Anti-fogging agents provided by eyecare professionals for use with eyeglass lenses are also suitable for Zeiss eyepieces.



- Please observe the instructions for use supplied with each anti-fogging agent.

Anti-fogging agents do not only ensure fog-free eyepiece optics. They also clean the eyepiece optics and protect them from dirt, grease, dust, fluff and fingerprints.

Cleaning mechanical surfaces

All mechanical surfaces of the system can be cleaned by wiping them with a

damp cloth. Do not use any aggressive or abrasive cleaning agents.

Remove any residue using a mixture of 50% ethyl alcohol and 50% distilled water plus a dash of household dish-washing liquid.



- Apply the cleaning agent on a soft, clean cloth (do not spray or apply it directly on the device) and wipe the device with the damp cloth.
- Do not allow any cleaning agent to seep into the device.

Sterile use (drapes)

Asepsis sets

The asepsis sets supplied by Carl Zeiss contain rubber caps and handgrips that you can sterilize in the autoclave. More detailed information on sterilization is available in the "Preparation of Resterilizable Products" user manual for the particular asepsis sets.

Drapes

Sterile single-use drapes are also available for sterile covering of the system, see page 112.

We recommend the following drape types:

- Model 70, no. **306070**
- Model 25, no. **306025** (recommended with AutoDrape option)
- Model 71, no. **306071** (recommended with AutoDrape option)
- Model 26, no. **306026** (recommended with AutoDrape option)



- When fitting the sterile drape, ensure there is sufficient slack to allow for movement of the microscope carrier and surgical microscope.
- Do not cover the ventilation opens - this ensures the lamps are cooled adequately and no lamp fails.

Disinfection

**CAUTION**

It may be a requirement to disinfect surfaces.

Surface damage on device

- Use a disinfectant based on aldehyde and/or alcohol. The addition of quaternary compounds is acceptable. To prevent damaging surfaces, disinfecting components other than those listed below must not be used.

The maximum concentrations are:

- For alcohol (tested with 2 propanol): 60%
- For aldehyde (tested with glutaraldehyde): 2%
- For quaternary compounds (tested with DDAC): 0,2%

Environmental protection measures

Note on disposal

User information on the disposal of electrical and electronic devices



This symbol means that the product must not be disposed of as normal domestic waste.

The correct disposal of electrical or electronic devices helps to protect the environment and to prevent potential hazards to the environment and/or human health which may occur as a result of improper handling of the devices concerned.

For detailed information on the disposal of the product, please contact your local dealer or the device manufacturer or its legal successor. Please also note the manufacturer's current information on the Internet. In the event of resale of the product or its components, the seller is required to inform the buyer that the product must be disposed of in accordance with the applicable national regulations currently in force.

For customers in the European Union

Please contact your dealer or supplier if you wish to dispose of electrical or electronic devices.


Information on disposal in countries outside the European Union

This symbol is only applicable in the European Union. For the disposal of electrical and electronic devices, please observe the relevant national legislation and other regulations applicable in your country.

Disposal of outer packaging material

- In order to dispose of outer packaging materials please refer to the local guidelines and bylaws.

System Data



Technical data	204
Microscope.....	204
Electrical data	206
Weights.....	208
EMC (electromagnetic compatibility)	209
Electromagnetic interference	210
Electromagnetic immunity for ME equipment and ME systems	211
Electromagnetic immunity for non-life-supporting ME equipment and ME systems.....	213
Recommended safety distances between portable and mobile RF communication equipment and the OPMI® VARIO 700.....	215
System combinations.....	216
Ordering data	217
Further accessories	220
Approval data.....	222
Environmental conditions	223

Technical data

Microscope

Component	Property
Magnification	Magnification factor $\gamma = 0.4x - 2.4x$
Varioscope objective lens	Working distance: $\geq 200 \dots \leq 500$ mm Motorized, adjustable speed. Possibility of manual adjustment (Light Only mode).
Tube / Eyepieces	<ul style="list-style-type: none"> – Tiltable tube, tilt range 180°, focal length $f = 170$ mm – (Optional) tube $f = 170/260$ mm – 10x (12.5x) magnetic widefield eyepieces with integrated eyecups.
XY coupling (optional)	Motorized XY adjustment in the entire microscope tilt and swivel range, adjustment range 30 mm x 30 mm, movement speed can be controlled.
Max. additional load on the microscope	The maximum admissible load (accessory equipment) on the microscope <u>may not exceed 7 kg!</u>

Microscope magnifications

WD in mm	Tube f in mm	10x eyepiece				12.5x eyepiece			
		Zoom	FOV Ø in mm	Magn.	FOV area in mm	Zoom	FOV Ø in mm	Magn.	FOV area in mm
200	170	0.4	97.0	2.2	73.8	0.4	83.1	2.7	54.2
		2.4	16.2	13.0	2.1	2.4	13.9	16.2	1.5
	200	0.4	83.4	2.5	54.6	0.4	71.5	3.1	40.1
		2.4	13.9	15.1	1.5	2.4	11.9	18.8	1.1
350	170	0.4	146.1	1.4	167.5	0.4	125.2	1.8	123.0
		2.4	24.3	8.6	4.7	2.4	20.9	10.8	3.4
	200	0.4	125.6	1.7	123.8	0.4	107.7	2.1	91.0
		2.4	20.9	10.0	3.4	2.4	17.9	12.5	2.5
500	170	0.4	192.8	1.1	291.8	0.4	165.2	1.4	214.3
		2.4	32.1	6.5	8.1	2.4	27.5	8.2	6.0
	200	0.4	165.8	1.3	215.8	0.4	142.1	1.6	158.5
		2.4	27.6	7.6	6.0	2.4	23.7	9.5	4.4


WD = working distance; FOV = field of view; Magn. = magnification; 0.4 = min. / 2.4 = max.

Electrical data

Component	Property
Line connection	<p>Only connect the system to wall outlets which are provided with a properly connected protective ground conductor.</p> <p>To prevent malfunctions, signal lines may not run parallel to power supply lines.</p> <p>In the United States, the device may <u>only</u> be operated at max. 120 V!</p>
Rated voltage	<p>115 VAC (100 - 125 VAC)</p> <p>230 VAC (230 - 240 VAC)</p>
Current consumption	Max. 1200 VA
Rated frequency	50 - 60 Hz
Fuses	Automatic circuit breaker
Electrical outlets	<p>Power socket 100-240 VAC, max. 500 VA</p> <p>Video output YPrPb (only for HD camera)</p> <p>Video output DVI-D (only for HD camera)</p> <p>Video outputs HD-SDI (only for HD camera)</p> <p>Video output Y/C</p> <p>Video output FBAS</p> <p>Remote socket for an external signal with max. 24 V / 0.5 A.</p> <p>Ethernet connector (option)</p>
Xenon light source *	Xenon short-arc reflector lamp, color temperature: approx. 6000 K
Superlux 180	Rated power: approx. 180 W
*) The maximum brightness of the lamp(s) can drop over time (normal device property).	<ul style="list-style-type: none"> - UV/IR heat protection filter - Pinhole stop controlling the light intensity
	Second lamp in lamp housing, swung in manually


Remote socket

View of connector side



Component	Property
Xenon light source * Superlux 300 (option) *) The maximum brightness of the lamp(s) can drop over time (normal device property).	Xenon short-arc reflector lamp, color temperature: approx. 5600 K Rated power: approx. 300 W
	<ul style="list-style-type: none"> - UV/IR heat protection filter - Pinhole stop controlling the light intensity
	Second lamp in lamp housing, swung in manually

Dimension diagram



Weights

Accessories


CAUTION

Do not exceed the maximum weight load.

The maximum load on the microscope body may not exceed 7 kg.

Description	Cat. No.	Weight
Straight tube	303765-0000-000	0.75 kg
180° tiltable tube, f = 170 mm	303791-0000-000	1.11 kg
180° tiltable tube, f = 200 mm	303792-0000-000	1.1 kg
Folding tube f170/f260	303771-9020-000	1.0 kg
10x widefield eyepieces (2 units)	305542-0000-000	0.26 kg
12.5x widefield eyepieces (2 units)	305543-9901-000	0.26 kg
10x widefield eyepieces with reticle (2 units)	000000-1023-184	0.26 kg
12.5x widefield eyepieces with reticle (2 units)	000000-1023-188	0.26 kg
Stereo bridge	1040-085	1.34 kg
Beam splitter 50/50	301513-9902-000	0.4 kg
Beam splitter 20/80	301503-9901-000	0.4 kg
Stereo coobservation tube	1063-869	1.09 kg
Camera adapter f = 340 mm T2	1022-973	0.42 kg
External micromanipulator (for weight data, refer to the user manual of the external micromanipulator)	-	-
Video objective lens f = 85 mm	301677-9085-000	0.4 kg

OPMI VARIO 700

Weight of OPMI VARIO 700 with floor stand (without accessories)	approx. 275 kg
Transport case	approx. 220 kg
Weight of OPMI VARIO 700 with floor stand transportable in the transport case	approx. 495 kg

EMC (electromagnetic compatibility)

The device complies with the EMC requirements of EN 60601-1-2:2007. While operating the device, observe the EMC precautions specified below.

- Only use spare parts approved by Carl Zeiss for this device.
- Do not use any portable or mobile RF communication equipment in the vicinity of the device as this may impair the device's function.
- Do not use a mobile phone in the vicinity of the equipment because the radio interference can cause the equipment to malfunction. The effects of radio interference on medical equipment depend on a number of various factors and are therefore entirely unforeseeable.
- Please note the EMC guidelines on the following pages.

Electromagnetic interference

The OPMI® VARIO 700 is intended for operation in an electromagnetic environment as specified below. The customer or the user of the device OPMI® VARIO 700 is responsible for ensuring that the device is operated in such an environment.

Interference measurements	Compliance	Electromagnetic environment - guidelines
RF emissions as per CISPR11	Group 1	The OPMI® VARIO 700 uses RF energy only for its internal functions. As a result, RF emissions are very low and unlikely to cause any interference in nearby electronic devices.
RF emissions as per CISPR11	Class A	The device OPMI® VARIO 700 is suitable for use in various facilities including locations in residential environments and those directly connected to the public power supply network which also supplies buildings used for residential purposes.
Harmonic emissions as per IEC 61000-3-2	not applicable	
Emission of voltage fluctuations/flicker as per IEC 61000-3-3	not applicable	

Electromagnetic immunity for ME equipment and ME systems

The OPMI® VARIO 700 is intended for operation in an electromagnetic environment as specified below. The customer or the user of the OPMI® VARIO 700 is responsible for ensuring that the device is operated in such an environment.

Immunity tests	IEC 60601 - test level	Compliance level	Electromagnetic environment - guidelines
Electrostatic discharge (ESD) as per IEC 61000-4-2	±6 kV contact discharge ±8 kV air discharge	±6 kV contact discharge ±8 kV air discharge	Floors should be made of wood or concrete or be covered with ceramic tiles. If the flooring contains synthetic materials, the relative humidity must be at least 30%.
Fast transient/burst immunity as per IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	The quality of the supply voltage should be that of a typical business or hospital environment.
Surges as per IEC 61000-4-5	± 1 kV phase-to-neutral voltage ± 2 kV phase/neutral to ground voltage	± 1 kV phase-to-neutral voltage ± 2 kV phase/neutral to ground voltage	The quality of the supply voltage should be that of a typical business or hospital environment.
Voltage dips, short interruptions and voltage variations as per IEC 61000-4-11	< 5 % U_T (> 95 % dip of U_T) for 0.5 cycle 40 % U_T (60 % dip of U_T) for 5 cycles 70 % U_T (30 % dip of U_T) for 25 cycles < 5 % U_T (95 % dip of U_T) for 5s	< 5 % U_T (> 95 % dip of U_T) for 0.5 cycle 40 % U_T (60 % dip of U_T) for 5 cycles 70 % U_T (30 % dip of U_T) for 25 cycles < 5 % U_T (95 % dip of U_T) for 5s	The quality of the supply voltage should be that of a typical business or hospital environment. If the user of OPMI® VARIO 700 requires continued function even in the event of interruptions in the power supply, we recommend to power the device from OPMI® VARIO 700 an uninterruptible power supply or a battery.

Immunity tests	IEC 60601 - test level	Compliance level	Electromagnetic environment - guidelines
Power frequency (50/60Hz) magnetic field as per IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels typical of business and hospital environments.

NOTE: U_T is the AC supply voltage prior to application of the test level.

Electromagnetic immunity for non-life-supporting ME equipment and ME systems

The OPMI® VARIO 700 is intended for operation in an electromagnetic environment as specified below. The customer or the user of the OPMI® VARIO 700 is responsible for ensuring that the device is operated in such an environment.

Immunity tests	IEC 60601 - test level	Compliance level	Electromagnetic environment - guidelines
			<p>Portable and mobile radio communication equipment should not be used closer to the OPMI® VARIO 700, including cables, than the recommended safety distance that is calculated using the equation applicable to the transmission frequency involved.</p> <p>Recommended safety distance:</p>
<p>Conducted RF disturbances as per EN 61000-4-6</p>	<p>3 V_{effective value} 150 kHz to 80 MHz</p>	<p>3V</p>	$d = 1, 17 \sqrt{P}$
<p>Radiated RF disturbances as per EN 61000-4-3</p>	<p>3 V_{effective value} 80 MHz to 2.5 GHz</p>	<p>3 V/m</p>	$d = 1, 17 \sqrt{P} \quad \text{for 80 MHz to 800 MHz}$ $d = 2, 33 \sqrt{P} \quad \text{for 800 MHz to 2.5 GHz}$
			<p>where P is the output power rating of the transmitter in watts (W) according to the transmitter manufacturer's specifications and d is the recommended safety distance in meters (m).</p>
			<p>Field strengths from stationary RF transmitters, as determined by a site survey^a, should be less than the compliance level in all frequency ranges.^b</p>
			<p>Interference may occur in the vicinity of equipment marked with the following symbol:</p>
			

Note 1

At 80 MHz and 800 MHz, the higher frequency range applies.

Immunity tests	IEC 60601 - test level	Compliance level	Electromagnetic environment - guidelines
----------------	------------------------	------------------	--

Note 2 These guidelines may not apply in all situations. Electromagnetic propagation is influenced by absorption and reflection by structures, objects and persons.

^a Field strengths of stationary transmitters such as base stations for mobile telephones and mobile land radio equipment, amateur radio stations, AM and FM radio broadcast and TV broadcast transmitters cannot be theoretically predicted accurately. To assess the electromagnetic environment with respect to stationary RF transmitters, a site study of the electromagnetic phenomena should be considered. If the measured field strength in the location where the device is used exceeds the compliance levels indicated above, the device should be monitored to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the ME equipment or ME system.

^b Field strengths should be less than 3 V/m over the frequency range from 150 kHz to 80 MHz.

Recommended safety distances between portable and mobile RF communication equipment and the OPMI® VARIO 700

The OPMI® VARIO 700 is intended for use in an electromagnetic environment in which RF disturbances are controlled. The customer or the user of OPMI® VARIO 700 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communication equipment (transmitters) and the OPMI® VARIO 700, depending on the output power of the communication equipment as specified below.

Rated output power of transmitter [W]	Separation distance depending on transmission frequency [m]		
	150 KHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d = 1, 17\sqrt{P}$	$d = 1, 17\sqrt{P}$	$d = 2, 33\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.69	3.69	7.38
100	11.67	11.67	23.33

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation indicated for each column, with P being the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer's specifications.

NOTE 1

At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2

These guidelines may not apply in all situations. Electromagnetic propagation is influenced by absorption and reflection by structures, objects and persons.

System combinations

This device can be extended to a tested system with various devices (see accessory equipment in the user manual).

This system combination has been approved by Carl Zeiss.



CAUTION

-
- When setting up an ME system, please make sure to comply with the requirements of EN 60601-1-1: 1990 or EN 60601-1: 2005, Chapter 16, and follow the measures below:
 - Never place multiple connectors on the floor.
 - Never connect any additional multiple connectors.
 - Only connect compatible components to the system.
 - Ensure that you do not exceed the max. permissible load capacity of the multiple connectors.
 - Use multiple connectors only for components which are part of the system.
-
- Any additional equipment connected to medical electrical devices must demonstrably comply with the applicable IEC or ISO standards (e.g. IEC 60950 for data processing equipment). In addition, all configurations must meet the normative requirements for medical systems (see IEC 60601-1-1 or Clause 16 of the 3rd edition of IEC 60601-1). Anyone connecting additional equipment to medical electrical devices is a system configurer and as such responsible for compliance of the system with the normative requirements for systems.

Ordering data

Only operate the system with the accessories included in the delivery package and approved by Carl Zeiss. You will find the contact responsible for orders in your country on this website:

www.meditec.zeiss.com

Eyepieces

Description	Cat. No.
10x widefield eyepieces	305542-0000-000
12.5x widefield eyepieces	305543-9901-000
10x widefield eyepieces with reticle	000000-1023-184
12.5x widefield eyepieces with reticle	000000-1023-188

Tubes

Description	Cat. No.
Straight tube f = 170 mm	303765-0000-000
2 sleeves for tubes with screw thread	305542-0107-000
180° tiltable tube, f = 170 mm	303791-0000-000
180° tiltable tube, f = 200 mm	303792-0000-000
Folding tube f170/f260 (standard) (optional)	303771-9020-000
Rotary dovetail for binocular tubes	301007-0000-000

Components for documentation and coobservation

Description	Cat. No.
Stereo coobservation module with two pivot joints	1063-869
Stereo bridge	1040-085
Camera adapter f=340 mm T2	1022-973
FlexioStill Adapter	301604-9010-000
Still interface	301604-9100-000
FlexioMotion adapter	301604-9050-000

Description	Cat. No.
Motion interface	301604-9500-000
Beam splitter 50/50 (option)	301513-9902-000
Beam splitter 20/80 (option)	301503-9901-000

Video / audio cables

The listed S-VHS cables cannot be connected to the integrated OPMI VARIO 700 cameras because no suitable socket is available.

Video/audio cables (standard definition)

Description	Cat. No.
1 m BNC-BNC connecting cable (option)	301687-9100-000
10m BNC-BNC connecting cable (option)	301687-9101-000
5m S-VHS mini-DIN connecting cable, 2x 4 pins (option)	308203-3040-000
10m RGB cable (option) for YPrPb	302681-8740-000
1.5m audio cable, CINCH-CINCH (option)	301687-9103-000
Cable tie set (20 pcs.) (option)	301687-9106-000
10 m Video Out Lemo, S-Video connecting cable	302681-8756-000
10 m Video In/Out Lemo, connecting cable(option) (only for 3CCD)	302681-8757-000
BNC cable for HD-SDI	
10m Component HD YPbPr connecting cable (option)	308203-3080-000
2m S-VHS mini-DIN connecting cable, 2x 4 pins for external SD monitor (option)	301687-9102-000
10m S-VHS mini-DIN connecting cable, 2x 4 pins for recording (option)	308203-3070-000
2 m BNC cable set (option)	301687-9109-000

Video/audio cables (high definition)

Description	Cat. No.
1 m HD-SDI connecting cable, 75 Ohm, 2x BNC pin (option)	305989-8762-000
5 m HD-SDI connecting cable, 75 Ohm, 2x BNC pin (option)	305989-8763-000
10m HD-SDI connecting cable, 75 Ohm, 2x BNC pin (option)	305989-8764-000
15m HD-SDI connecting cable, 75 Ohm, 2x BNC pin (option)	305989-8766-000
4m TRIO 610 extension cable	308203-2400-000
2m DVI monitor cable	308203-3020-000

Other cables and cable components

Description	Cat. No.
10m Ethernet cable, 2x RJ45	305946-8660-000
Cable clip	302581-8675-000
Cable tie set (20 pcs.) (option)	301687-9106-000


Further accessories

Description	Order number
MEDIALINK 100	308203-2200-000
User manual G-30-1716: "OPMI VARIO 700 on the floor stand"	000000-1586-434


Foot control panels and connecting cables

Description	Order number
Wired 14 function foot control panel (FCP)	304970-9015-000
Foot control panel with 14 functions, wireless (FCP WL)	304970-9020-000
Foot control panel 2, with 14 functions, 6m cable length (optional)	304979-9020-000
Adapter for FSP II (foot control panel 2)	305946-8181-000
3 m cable for 14 function foot control panel (FCP & FCP WL)	304970-8730-000
6 m cable for 14 function foot control panel (FCP & FCP WL)	304970-8760-000

Consumables



1



2

Description	Cat. No.
Replacement lamp module complete with 2 xenon lamps in transport case and with return card (1)	304977-9036-000
300-W xenon lamp (2)	304949-9001-000
Optics cleaning set	000000-1216-071
Sterile drapes type 70	306070-0000-000
Sterile drapes type 71 (only use when AutoDrape option available)	306071-0000-000
Sterile drapes type 25 (only use when AutoDrape option available)	306025-0000-000
Zeiss sterile drapes, type 26, packet with 5 pcs VisionGuard Drape for OPMI ^(r) with 2 observation tubes and camera. Extra-long and -wide version. Dimensions: 132 cm x 391 cm (only use when AutoDrape option available)	306026-0000-000
VisionGuard replacement lenses (sterile, 20/box)	306001-0000-000
Asepsis set - sterilizable caps <ul style="list-style-type: none"> - Internal diameter 12 mm - (pack of 6) 	305810-9002-000
Asepsis set - sterilizable caps for controls on microscope and for PD adjustment knob on 45° inclined tube <ul style="list-style-type: none"> - Internal diameter 22 mm - (pack of 6) 	305810-9001-000
Asepsis set - sterilizable caps <ul style="list-style-type: none"> - Internal diameter 27 mm - (pack of 6) 	305810-9008-000
Asepsis set - sterilizable caps for PD adjustment knob on 180° tiltable tube <ul style="list-style-type: none"> - Internal diameter 51 mm - (pack of 6) 	305810-9003-000

Approval data

Approvals and requirements

Description	Labeling
Electrical standard	EN 60601-1:2005
	CAN/CSA-C22.2 No. 601.1
	Protection class I, degree of protection IP 20
EMC requirements	The system meets the EMC requirements of EN 60601-1-2:2007, Class A
CE mark	The device meets the essential requirements stipulated in Annex I to the Medical Device Directive 93/42/EEC. The system is labeled with 

Approvals and requirements for the optional FCP-WL foot control panel

Description	Labeling
FCC	If a radio module for the optional FCP WL foot control panel is integrated in the device, the device meets the requirements of Part 15 of FCC (Federal Communications Commission, USA)
CE mark	If a radio module for the optional FCP WL foot control panel is integrated in the device, the FCP also meets the requirements of EU Directive 1999/5/EC. The device is labeled with: 

Environmental conditions

**CAUTION**

Observe the storage, transportation and operation conditions.

For operation

Property	Admissible range
Temperature	+ 10°C ... + 35°C
Rel. humidity	30 % ... 75 %


- Do not use power-operated devices included in the delivery package
 - in explosion-risk areas,
 - at a distance of less than 25 cm from inflammable anesthetics or volatile solvents such as alcohol, benzene or similar chemicals.
- Ensure that the installation conditions and the use of the device meet surgical requirements:
 - Low vibration
 - Clean environment
 - Avoidance of extreme mechanical loading
 - The device is not intended for use in environments concentrated with oxygen.
- Do not use or store the system in damp rooms.
Do not expose the system to water splashes, dripping water or sprayed water.

For transportation and storage

Property	Admissible range
Temperature	- 20°C ... + 60°C
Rel. humidity (without condensation)	10 % ... 92 %
Air pressure	500 hPa ... 1060 hPa

The device is not intended for use in environments concentrated with oxygen.

Annex - (Options)



Record (integrated video record function over USB port - optional)

Record (integrated video record function over USB port - optional)	225
Description	228
Normal use.....	228
Storage options.....	229
Video 1/2 button functions.....	230
Video 2/2 button functions.....	232
Menu functions on the video monitor.....	234
Network interface - overview for networking	242
Configuring the device	244
Video option checklist	252
Recording photos or videos	252
Recording photos	252
Recording images.....	253
Recording videos	254
Recording videos	255
Viewing images/ videos	256
Viewing images / videos	257
Importing images/ videos.....	258
Managing patient data.....	260

Importing patient data265

Exporting patient data.....265

Description

Normal use

The recording option is used for easy recording of video sequences and single images in clinical applications.

The recorded video and image data can be stored on an external USB storage medium or in a shared directory within a network. The unambiguous, automatic storage of images and videos ensures reliable assignment to the relevant patient.

The recording option is recommended only in combination with a connected video monitor because all control functions (menus) and image details are shown only on the video monitor.

NOTE

- The user (IT officer) is responsible for data protection and ensuring that no viruses are transferred to the system via the network connection.
- The user is responsible that the USB stick, which is used for the data exchange, is free of any viruses.



The company Carl Zeiss accepts no liability for any loss of patient, image and video data as well as system or user-specific configuration data. If required, arrange for patient, image and video data as well as all system settings to be backed up by your IT administrator on a regular basis.

Storage options

- USB storage medium*
- Network PC**

* Only one USB storage medium at a time is supported.

** Archiving in the network must be configured: Menu / Config / Media / Storage Device : LAN configuration.

Note

USB data security!

USB storage media are not suitable for the permanent storage of patient data. Arrange for data to be backed up on a regular basis by your IT administrator.

- Check the function of USB storage media before use.
 - Remove USB storage media only while the live image is displayed (no menu opened, no action, no recording).
 - Do not remove USB storage media during Record/Play/Photo/Display/Index/Import.
 - Check USB storage media for correct function at regular intervals using a PC.
-

Password protection

The use of password protection is the responsibility of the customer or institution operating the device.

- Configurable password protection on system startup. Password protection can only be deactivated after entry of the password.
- User and administrator passwords comprise 4 figures and can be changed.
The factory setting for the user and administrator passwords is "0000" in both cases.
- User password protection can be activated and deactivated.
Factory setting: user password protection is deactivated.
- If password protection is active, the user or administrator password must be entered after startup to permit any functions to be executed.
Various functions (all network settings and the system settings) are only accessible via admin login.

Video 1/2 button functions



<Record> button

- If you press the button once, video recording is started and the video is saved in accordance with the configuration. REC is displayed at the bottom left in the OSD (on-screen display) and the top right of the display; REC ON or REC OFF is shown on the additional display.
- If you press the button again, the current video recording is stopped. REC disappears from the display (REC OFF is shown on the additional display).



<Live> button

Closes all open dialogs and aborts all actions except <record>.

- Press the button once to return to the live mode.



<Photo> button

Still images can be created both in the live mode and during video recording or playback.

- Pressing the button once triggers the photo function (display at the top left: Photo) and saves the image in accordance with the configuration. Then the live image is displayed again.



<Play> button

- Press the <Play> button or select <Play> in the context menu.
- If you press the button once: "Pause", the current display is shown as a freeze image.
- If you press the button again, the video continues.



<Forward> button (fast forward)

- Press the button once to fast forward the selected video sequence to search it for scenes of interest to you. The sequence is played at a higher speed.
- Press the <Play> button to stop fast forward.



<Reverse> button (fast rewind)


- Press the button once to rewind the selected video sequence to search it for scenes which are of particular importance or interest to you. The sequence is played backward at a higher speed.
- Press the <Play> button to stop fast rewind.



During video recording, the photo function can be used to create additional single images.

These images can also serve as markers. During playback of the video, you can use these marker images to find scenes which are of particular interest to you (see page 254).

You can create further images during playback of the video using the photo function.



Video 2/2 button functions

<Index> button

Press the button once to display the thumbnail index of the current patient.

Use the control buttons to select a video or image. A blue frame marks the selected file.

- Press the <Index> button or <menu> button to end the index display.


<N.Pat.> button

Creates a new patient directory on the external storage medium (patient name: Unknown) and uses it as the **current** patient.

<Menu> button

Opens context menus or closes context menus or dialogs.

- Press the button once to display the menu in the OSD.
- Press the button again or press the <live> button to hide the menu.




<select> control buttons (cursor buttons)

The < > / \blacktriangle control buttons permit you to navigate in the menus displayed in the OSD.

<enter> key (<U346>)

Press this button to confirm your selection or to trigger a function.



Menu functions on the video monitor

You can enter all system-relevant data in a configuration menu where they can be saved.

You have the possibility of writing your saved system settings to a USB medium or to load them from this medium to the system. (This provides an easy method of identically configuring several devices in a hospital).

Config menu



- Use the < > control keys to navigate between the submenus: (Media < > User < > Network < > System).
- Use the ⬆ ⬇ ⬅ control buttons to navigate within the submenus.

<Media>

1 Selection of the supported video and image formats

2 Freeze image function

3 Selecting the storage device

4 <Auto Delete Mode>

This mode is enabled, data cannot be assigned to a patient, but are saved to a neutral folder. Files older than 3 days are automatically deleted. To change to a patient, must be disabled.

5 Audio - Recorder - Player

Setting the microphone sensitivity and playback volume. The audio signal is automatically output when the MPEG video is played back via a player.



Videos recorded without sound can only be played back without sound.



<User>

6 Language selection for OSD text

OSD is always in English, only the language of menus and dialogs is selectable.

7 Sounds

Activating/deactivating the notification sound, confirmation sound

8 Passwords (password protection see page 229)


A password prompt can be activated to restrict system access to authorized persons.

<OSD>


9 OSD position

Selecting the text box position on the optional 17" video monitor:

17" video monitor



– Standard display



– Position selection:

Three OSD elements can be displayed in three lines at six different positions:



Top left - center - right; bottom left - center - right

10 OSD elements

You can select three OSD elements which can be displayed in three lines, one below the other, at six different positions.



The light value is automatically updated on the display when changes are made on the microscope.



<Network>



To configure the network settings, you must log in as "Administrator": (Config/User/Passwords/Login). (Config/User/Passwords/Login).

11 IP settings

MAC address

The MAC address (Media Access Control address) is the hardware address of each individual network adapter that is used for the unambiguous identification of the device in a computer network, if required.

DHCP

DHCP enables the automatic network integration of a new computer into an existing network without the need for its manual configuration. All that normally needs to be done on the client is setting the automatic reference of the IP address. Confirm the IP settings by clicking on "Apply" to make them effective.



When "Apply" is activated under IP Settings, Share is automatically disabled.


- Check the share parameters, as they may need to be adapted. Click on "Apply" again to activate them.

12 Share

The video and image data recorded can be stored in a shared directory on a PC connected to the network.

Requirements for the connection of the archiving PC (client-server)

- You need Admin rights.
- **You need the Ethernet data** of the archiving PC (Run `\cmd\ipconfig`) to configure the record option.
- Create a new user (name/password) in user accounts under which you want to save the recording data.
- Enter the necessary PC attributes in the Share dialog and click "Apply" to enable the connection.



<System>



To perform the system settings, you have to log in as administrator: (Config/User/Passwords/Login).

13 Export/Import

You have the possibility of exporting your saved system settings to a USB medium ("Export Config") and importing these system settings to a different system ("Import Config").

This provides an easy method of identically configuring several devices in a hospital.

All passwords including the administrator password are also overwritten by the new values. Make sure to use current (known) data.

Exporting logfiles (no administration rights required)

Automatically generated logfiles can be saved on a USB stick/HDD as a password-protected ZIP archive and sent to Carl Zeiss Service for evaluation.

14 Special

The installed software version is displayed.


- Software update
A software update can be performed via a USB data medium. Old system configurations are retained in this process.
- Factory settings
The configuration is reset to the factory settings. The passwords are also reset.
- Export license
creates a license file on the USB stick/HDD. A ZIP archive providing information on the licenses is stored on the USB stick/HDD (without password protection).

15 Setting the date and time

- Set the date → yyyy/mm/dd
- Set the time → hh:mm

16 Foot control panel

- Photo - Record
- Record - Capture





Network interface - overview for networking

- Data can be stored in a shared directory on a host computer.
- It is possible to store image and video data in a network and to access this data.
- Patient data can be imported from a shared directory on a host computer.
 - Software programs from other manufacturers can create and read patient-related folders.
- The option can be automatically or manually configured for IP networking.
- The parameters for network access are factory-set.
- The host computer and the shared directory must be manually configured.

Host computer

- Fixed identification is required.
 - The IP address of the host must be known to enable its configuration.
 - The host computer must have a fixed IP address.
- A shared directory is created.
 - The directory must be approved and enabled for shared use.
 - This is done via Microsoft Windows file and printer sharing and the TCP/IP protocol.
- The access parameters used for the shared directory must be the same as those configured in the option.
 - Administrator rights are required for the correct setting of the shared directory.
 - User, password and the full access authorizations for the shared directory must be set.
 - The default parameters of the option can be used.


Network

- Ethernet LAN with 100 Mbit/s must be available.
 - OPMI VARIO 700 is optionally equipped with a RJ45 Ethernet connector.
- A DHCP service can be available.

- The network parameters need not be manually configured.
- Some networks are subject to increased security.
 - The MAC address of a connected device may have to be entered manually.
 - Users may be automatically blocked after incorrect authorization.

USB interface

The read and write speeds of the USB 2.0 storage medium must be higher than 5 MBytes/s.




Configuring the device

For the connection of your recording option to your network, admin rights may be required to be able to create the shared directory and the user account. Make sure that you have admin rights for the computer that you want to configure for networking with the recording option. Otherwise, please contact your system administrator before proceeding.

If your network connection with the host computer is protected by a firewall or another safety feature, please contact your system administrator. If you do not wish to use the factory-configured network settings of the recording option, you have to log in as administrator in the recording configuration mask before proceeding. The preset admin password is "0000" (4x zero, no quotation marks).

Logging in as administrator

- Press the MENU button.
- Select the menu item CONFIG.
- Navigate to the User tab.
- Select Login (1) to log in as administrator and to be able to access the network and system settings.



Configuring your computer for networking

Creating a directory The record function writes video and image data to a directory specified by you on your computer.



- Create a directory in your computer to which the videos and images should be written.

Make sure that the user account has both read and write access to the shared directory on your computer.

- Grant network access rights for the newly created directory so that the record function in the network has access to it.

The record function has been factory-set such that it looks for a shared directory named "medialink" (no quotation marks). If your shared directory has a different name, this name must be changed in the configuration menu for the record function on the display.

Creating a user account To enable the record function to access the shared directory created by you, it must be able to log in to your computer. At the factory, user name "MEDIALINK" and password "MEDIALINK" have been preset for the record function (both in uppercase letters without quotation marks).

- Create a user account on your computer using this name and password, or use a different account. If you use a different user account or password, you must change the default settings in the configuration screen accordingly.
- Always press the "Apply" button to save changes.



Connecting the record function to your network

The recording option has been preset in such a way that DHCP is enabled for networking. As a result, the connection to most networks accepting this configuration type is possible without any problems.

- Connect the device to your network via the LAN connector on the connector panel.
 - If your network has been configured to accept DHCP and if no limitations exist, it should assign a unique IP address to the device.




- If there are any limitations in your network so that only specific MAC addresses are accepted, you have to configure your switch or router in such a way that it accepts the MAC address of the device. The illustration "IP settings" shows where you can find this information. The MAC address of the device (recording option) is indicated at the top left of this mask.
- If your network does not support DHCP or if you want to assign a fixed IP address to the device, remove the checkmark for DHCP and enter the IP address data as shown in the illustration "IP settings". Make sure that the IP address, subnet mask and default gateway are correctly set in accordance with your network requirements.
- Press the "Apply" button in the recording menu on the display for your changes to become effective!

Configuring the connection to your computer

After connecting the OPMI Vario 700 to the network and successfully setting all network parameters, you need to configure the record function so that it knows which computer has been set up for file sharing.

- Enter the IP address of the computer in the Host IP box on the configuration screen for the record function (Fig. "IP setting"). The preset host IP address is 192.168.0.1.
- Test the connection by pressing the Ping button.
 - If the IP address is correct and the computer replies to Ping queries, you will receive a positive confirmation.



Presetting the storage device: USB or LAN

If you want the recording function to store images and videos by default in the network configured by you, you have to set the appropriate configuration in the configuration menu. The USB medium is the default (factory-set) storage device (see illustration "Selecting the storage device").

- Select the USB storage option if the data should be saved to a USB storage medium.



When the device is ready for operation ...

When all changes have been saved and LAN has been preset as your storage device, the recording device tries to log in to the computer which you configured with the user account and password as described above. If the login is successful, a folder for the current patient is generated in the shared directory which you created. All recorded images and videos are saved here at the moment when they are captured.

Saving the configuration

If the transfer of images and videos works perfectly, we recommend saving your configuration to a USB medium. In the unlikely event of a malfunction, you can use this configuration data to restore the system and your settings (including the network configuration) after repair or replacement of the device.

- To save your configuration, simply insert the USB medium and select the "Export Config" option (2).



Video option checklist






- Always check the following points (without the patient) before using the device:
 - The OSD (on-screen display) on the monitor functions correctly.
 - A USB storage medium has been connected and detected, or the network connection has been set up.

Recording photos or videos

Recording photos

- 1 Live image
- 2 Photo function: (<Photo> button)
 - Press the <Photo> button once; the image is displayed for approx. 2 seconds and saved in the graphic index of the current patient (3). Then the live image is displayed again.
 - Active freeze image:
Press the <Photo> button once, the live image is frozen.
Press the <Photo> button once again; the displayed image is saved (3), then the live image appears again.
 - No saving: press the <Live> button, the live image is displayed again.
- 3 <Index>
Graphic index of the current patient.

Recording images



Recording videos

Video function: (<Record> button)


- Press <Record> once to start the video (rec), the video is saved in the graphic index of the current patient (1).
- Press <Record> again to stop the video.

1 <Index>

Graphic index of the current patient.




During video recording, the photo function can be used to create additional single images (e.g. Photo 1 - Photo 4). These images can also serve as markers.






During video playback, the marked positions (e.g. Photo 1 - Photo 4) can successively be accessed using the cursor buttons (\blacktriangledown) on the control panel or the remote control (\blacktriangledown).

Recording videos

Record



Press 1x → recording (>Rec<) is started
Press again → recording is stopped



Viewing images/ videos





<Index> Graphical index of the current patient

- Press once for the graphical index of the current patient.
A video or image is selected with cursor buttons < > . The blue frame denotes the object selected.
- Clicking the object opens a context menu to display, copy or delete the object.

USB or LAN must be configured beforehand as the storage location: Menu / Config / Medium / Storage location : USB or LAN.

Viewing images / videos

Press the Index button



Importing images/ videos

Import (<menu> button)


displays all the files of the non-active storage location (USB or LAN)

- Press the <Menu> button, open <Menu>, open <Import>
- Select the folder with <Enter> to access the graphical index of the current folder.
A video or image is selected with cursor buttons < > . The blue frame denotes the object selected.
- Pressing <Enter> above the object opens a context menu for displaying or importing the object.



Importing images/videos

Press Menu button




Managing patient data



- Use the < > control buttons to navigate between the following actions: (Select, New, Filter, Rename, Delete).
- Use the ⬆ ⬇ control buttons to navigate within the submenus.

Creating a new patient

Using the <New Patient> key on the front panel



- Press <New Patient> key (4).
A new patient directory "Unknown" is created on the external storage medium and used as the **current** patient.
You can change the patient name (see Renaming the patient).

Using the patient menu




- Open <Menu> (1).
- Open <Patient> menu (2) using the ⬅ button or the <Enter> button.
- Select the <New> submenu (5) using the < > control keys and confirm your selection with the ⬅ key or <enter> key.
- Enter the patient's last name (6) and first name (7) via the keyboard displayed on the screen using the cursor buttons (< > or ⬆ ⬇).

Selecting a patient

- Open <Menu> (1).
- Open <Patient> menu (2) using the ⬅ button or the <Enter> button.
- Select <Patient> using the ⬆ ⬇ control buttons (3).
- Confirm your selection by pressing the ⬅ button or <enter> button.

Opening the <Patient> menu via the keyboard (option)

- Open <Menu> by pressing the enter key.
- Open the <Patient> menu using the enter key.
 - Use the < > keys of the keyboard for your selection in the submenu and confirm using the enter key.
 - Or: select <Patient> using the ⬆ ⬇ keys of the keyboard.



Searching for a patient (filter function)

- Press the <Menu> button, open <Menu> (1).
- Open the <Patient> menu using the  button or the <Enter> button.
- Select <Filter> submenu (8) using the < > control keys.
- Confirm your selection by pressing the  button or <Enter> button.
- Select the filter setting using the  control buttons (9).
 - If you select <No Filter>, all patients are displayed.
 - Depending on the selected range, the relevant patients are displayed.

Renaming a patient

- Open <Menu> (1).
- Open the <Patient> menu using the  button or the <Enter> button.
- Select <Rename> submenu (10) using the < > control keys.
- Confirm your selection by pressing the  button or <Enter> button.
- Edit the last name (11) and first name (12) via the respective keyboard in the OSD.

Deleting a patient


- Open <Menu> (1).
- Open the <Patient> menu using the  button or the <Enter> button.
- Select <Delete> (13) using the < > control buttons.
- Confirm your selection by pressing the  button or <Enter> button.
- Select the patient to be deleted using the  control buttons and delete the patient (14).

An inquiry is displayed, asking whether you really want to delete the patient.



The current patient cannot be deleted.


If you try to delete the current patient, the following message appears:
The current patient cannot be deleted!



The image and video files are stored on the USB storage medium according to the following folder structure.

- Under the root directory, a directory named "MEDIALINK" is created.
- Under the "MEDIALINK" directory, the patient directories are created.
- All data collected during the session is filed in a patient directory (structured according to family name/first name) until a new patient is created or the device is switched off and back on again.

The network drive is a shared directory whose name is assigned by the IT administrator of the hospital or office. This directory name is used in the same way as the "MEDIALINK" directory of the USB storage medium.



Importing patient data

The user can import videos and images with the supported formats (MPEG2 Standard, JPG, TIF;no MPEG4, H264, AVI, ...) to an active patient folder. An Explorer is displayed for this purpose, permitting the import of individual files from a USB storage medium.

This data can be copied to a patient folder in the network (LAN).

The contents of the non-active storage medium (if connected) is displayed. Image and video files can be copied from there.

Simultaneous import is possible only from one storage device. The FAT32 and NTFS formats are supported as file systems.

JPG files with up to 8 Mpixels and TIF images with up to 2 Mpixels can be displayed as thumbnail or image. For larger images, a fallback image is displayed.

Exporting patient data




Images and videos can be copied to a connected USB stick using "Copy to USB" (e.g. for patients to take with them).

- Insert a USB stick in one of the USB ports.
- Use the < > control keys to select the required files in the <Index> menu. A blue frame marks the selected file.
- When you click on the file, a context menu is opened.
- Select "Copy to USB", to copy the selected data to the USB stick.



Indexes



List of technical terms

Term	Explanation
Apochromatic optics	Optical system with maximum correction of chromatic aberration
Drapes	Sterile covers for surgical microscopes
Depth of field	Range of depth that appears sharply defined through the microscope
Radiation intensity	Radiation incident or emitted at right angles per unit time and unit area.

List of abbreviations

Term	Explanation
BEV	German acronym for " special emergency backup power supply "
BNC	B ayonet N eill C oncelmann (Coaxial bayonet-locking connector for high frequencies, named after its inventor)
CCD	C harge C oupled D evice Technology (optical sensor)
CE	C ommunauté E uropéenne (European Community) - The manufacturer declares that the device complies with the directives of the European Union.

Term	Explanation
CSA	C anadian S tandards A ssociation - non-governmental organization that tests and certifies product safety
D.o.F.	D epth o f F ield - depth of field management system
DPT	Diopters (refraction power of optical system)
DVI	D igital V ideo I nterface (digital interface for video data)
EMC	E lectromagnetic C ompatib I lity - defines non-interference of electrical and electronic devices with their environment
EN	E uropean N orm (European Standard)
CVBS	C olor V ideo B aseband Signal (television signal for color picture transmission)
FDA	F ood and D rug A dministration - US agency responsible for the safety regulation of foods, drugs and medical products
FCP	Wired 14-function foot control panel
FCP WL	Wired 14-function foot control panel
HD-SDI	H igh D efinition S erial D igital I nterface (digital interface for video and audio data)
HDTV / HD	H igh D efinition T ele V ision
ICC	I nformation and C ontrol C enter
IEC	I nternational E lectrotechnical C ommission
RF	R adio f requency
LED	L ight E mitting D iode
NTSC	N ational T ele V ision S ystems C ommittee (US institution defining the first color television standard)
OP	Operation
OPMI	Surgical microscope
OSD	On Screen Display
PAL	P hase A lternating L ine (color encoding system in analog TV systems)

Term	Explanation
PD	Interpupillary distance
ROI	Region of interest
SCI	S tereo C oaxial I llumination
SDTV / SD	S tandard D efinition T ele V ision (television standard with resolution lower than HDTV, e.g. PAL or NTSC)
UV	U ltra V iolet
Y/C	Analog interface that transmits separate color (luminance Y) and brightness signals (chrominance C).
YPbPr	Color model that transmits brightness information Y and color difference information Pb and Pr separately.

Keyword index

Numerics

180° tiltable tube 52

A

Acoustic signals 16
 Adding a user 132
 Additional display 55
 Additional equipment 216
 Address of manufacturer 2
 Adjusting the OPMI and micromanipulator to the same focal plane 74
 Adjustment wheel for interpupillary distance 52
 Anlog 68, 70
 Approval data 222
 Attaching sterile drapes 111
 Automatic Iris Control 18, 138
 Automatically balancing the surgical microscope (optional) 100

B

Backup lamp 20
 Batteries flat 184
 Brightness of the additional display, adjustment 168
 Burn injuries caused by illumination intensity 13
 Button assignment of foot control panel 73, 156
 Button for "Light only" mode 24

C

Camera setting 1/2 160
 Camera settings 2/2 162
 Changing device location 94
 Checklist 120
 Cleaning agents 200
 Configurations 216
 Configuring the brakes 158
 Configuring the video settings (optional) 165
 Configuring the XY coupling (optional) 159
 Configuring video source (camera) 163
 Configuring zoom 150
 Connecting navigation systems (option) 80
 Connection

Connector panel	66
Power inlet	66
Power outlet	66
Remote	66
Connector	
Ethernet	66
Controls	
Control and display panel	65
Tubes	52
Widefield eyepieces	53
Controls for manual balancing	64
Controls on the display	126
Coobservation equipment, mounting of	78
Creating a new patient	260

D

Delete button	127
Deleting a patient	262
Deleting a user	134
Digital enlargement	161
Diopters	
Scale	53
Setting ring	53
Disinfection	201
Display	124
Display, backup lamp is in use	62
Documentation equipment, mounting of	78
Drapes	200
Drapes, attaching	110

E

Electrical malfunction	178
EMC (electromagnetic compatibility)	209
Emergency backup power supply	10
Environment	
conditions	223
Exchanging xenon lamp for 300W Superlux (option)	190
Exporting patient data	265
Exposure mode	160
Eyecup	53
Eyepiece	
adjusting	107
mount	52
Eyepieces, mounting	76

F

Failure of a function	208
Failure of magnetic brakes	178
Failure of the control panel (touchscreen)	180
Field of use	8
Fine adjustment of microscope	104
First use	12
Fixing cable with strain relief device	89
Focusing	108
Focusing function, configuring	144
Focus-light control	20, 140
Focus-light link	20, 140
Foot control panel	
Button assignment	73, 156
Configuring	154
connecting	86
Foot control panel (optional)	72
Functional test prior to use (checklist)	120

G

Getting started	118
-----------------------	-----

H

Handgrips	
configuring	152
Handgrips on surgical microscope	50
Hazard symbols	7

I

Illumination settings	138
Importing images	258
Importing patient data	265
Importing videos	258
Inclined tube	52
Information symbols	7
Intended use	8
Interface on microscope mount	82
Interpupillary distance, adjustment	106

K

Key to symbols	7
Keypad	127

L

Lamp change of 300-W xenon light (optional)	174
Lamp change, 180W xenon light source	172
Lamp exchange xenon lamp Superlux 300W	190
Light - Adjusting the light	136
Light only	24
Light-only mode	178
Locking the device into position	94

M

Manual balancing through user guidance on the display	98
Manual lamp change in 180W xenon light source	172
Manual lamp change of 300-W xenon light (optional)	174, 176
Manually balancing the surgical microscope	97
Manufacturer's address	2
MediLive	
ImageBox	126
MindStream	68, 70, 126
Micromanipulator	74
Microscope components, mounting of	76
Microscope magnifications	205
Mobile phones	209

N

NAV connector	82
Navigation interfaces	82
Normal use	9

O

Opening the lamp module	62, 63
Optical surfaces, cleaning	199
Optics	164
Ordering data	217

P

Pairing	86
Password protection	229
Performing foot control pairing (optional)	88
Positioning the device in the OR	95
Potential equalization	66, 85
Potential equalization bolt	66
Power switch	65
Powering the system up/down	118

Problem	184
Procedure in the event of faults	171
Protective ground conductor	12

R

Radio link	184
Radio signal	184
Recording photos	252
Recording videos	254
Renaming a patient	262
Reticle	107

S

Safety	
Devices	18
Searching for a patient (filter function)	262
Selecting a patient	260
Selecting depth of field	145
Selecting user and language	130
Service PIN	169
Setting	161
Setting the HD DigiZoom	161
setting the HD format, setting of	166
Setting the SDTV format	167
Shift key	127
Sliding switch for rated voltage	66
Sound	16
Sounds	16
Space key	127
SpeedFocus, focusing	108
Starting for the first time	118
Sterile use	200
Strain relief device	66, 89
Surgical microscope	
Adjustment	106
Suspension system	
Connection	84
Switch	
for manual operation	65
Switching to the backup lamp	188
Switching to the second lamp	172, 174
Symbols and labels on the device	26
Symbols on the display	124
System configuration	216
System overview	43

System settings menu 1	164
System settings menu 2	169

T

Target audience	8
Trademark	2
Transporting the device	94
Tube, mounting	76

U

User	130, 159
User settings menu 1	150
User settings menu 2	159

V

Video 1/2 button functions	230
Video 2/2 button functions	232
Video focus function, configuring	142
Video image, automatic focusing	109
Video image, manual focusing	108
Video option checklist	252
Video output	
Composite - (CVBS)	68, 70
HD-SDI	70
Video input/output	68
Video output port	
DVI-D	70
Y/C (S Video)	68, 70
YPbPr	70
Viewing images	256
Viewing videos	256
Virtual keypad	127
Visual and acoustic signals	16

W

Weight balancing crank	64
Weight of accessories	208
Weights	208
White balance (option)	148
Working in Light-only mode	178, 180



Carl Zeiss Meditec AG

Goeschwitzer Strasse 51-52

07745 Jena

Germany

Fax: + 49 (0)7364 - 20 4823

Email: surgical@meditec.zeiss.com

Internet: www.meditec.zeiss.com