

Technical Data

Vismo

PVM-4761/4751/4731
PVM-4763/4753/4733
Bedside Monitor

Fighting Disease with Electronics



TD. PVM-4000 Op No. 01K Ver.01-04

This technical data may be revised or replaced by Nihon Kohden at any time without notice.
Some products may not be available in your country.

FEATURES

Intended Use

PVM-4000 series bedside monitor is intended for continuous monitoring, recording, and alarming of multiple physiological parameters of adults, pediatrics and neonates in the OR, recovery room, ICU, CCU, HCU, NICU, ER, ward and other areas.

Smart cables and MULTI connectors

Smart cable technology gives you high flexibility in the selection of parameters without traditional modules. When you plug a sensor cable into a MULTI connector, the monitor automatically detects the type of parameter and starts measuring. IBP and CO₂ can be measured with MULTI connector. (PVM-4763/4753/4733 only)

Display

The monitor has a high resolution, wide-viewing angle and 10.4-inch TFT color LCD for easy viewing.

Touch screen operation

The touch screen provides intuitive operation.

Smart guide function

When a technical alarm happens during monitoring, the Guide key appears at the top of the screen. If the Guide key is touched, the guide menu that indicates a possible solution will appear. The guide also shows measurement tips for each parameter.

PI (Pulse-amplitude Index)

The Pulse-amplitude Index indicates the percentage of pulsatile signal in the entire transmitted IR signal. PI supports to assess peripheral perfusion.

(PI is displayed when using Nihon Kohden SpO₂ sensors and Masimo SpO₂ sensors)

SQI bar graph for SpO₂

The SQI (signal quality index) bar graph shows the pulse waveform quality for SpO₂ measurement. (SQI bar graph is available only with Nihon Kohden SpO₂).

iNIBP

iNIBP is Nihon Kohden's unique algorithm to measure NIBP during inflation. It provides fast

and painless measurement of NIBP. Even if patient's blood pressure increase compared to previous measurements, iNIBP provides fast measurement of NIBP.

PWTT-triggered NIBP measurement

PWTT (pulse wave transit time) is continuously and non-invasively measured from the ECG and SpO₂. If a sudden, critical circulation change happens between the periodic NIBP measurements, PWTT may detect it and trigger NIBP measurement to confirm it.

esCCO

esCCO (estimated continuous cardiac output) is new technology to determine the cardiac output using Pulse Wave Transit Time (PWTT) and standard monitoring parameters – ECG, SpO₂ and NIBP. esCCO provides real-time, continuous and non-invasive cardiac output measurement alongside the familiar vitals sign parameters and it is a very effective cost-saving solution because it has no additional running costs or accessories. (QP-470P is required.)

PPV/SPV

PPV (Pulse Pressure Variability) and SPV (Systolic Pressure Variability) is an indicator of intravascular volume. It is useful in guiding fluid therapy for patients receiving mechanical ventilation. (PVM-4763/4753/4733 only)

Mainstream EtCO₂ for both intubated and non-intubated patients

PVM-4763/4753/4733 can use the compact and lightweight cap-ONE mainstream CO₂ sensor kit to obtain smooth and accurate EtCO₂ and respiration rate within 5 or 10 seconds.

Unique nasal adapter effectively catches the nasal and oral expired CO₂ without dilution for non-intubated patient (Available with TG-920P).

cap-ONE mask provides ET/CO₂ measurement and oxygen supply for non-intubated patients (Available with TG-980P).

NICU:

Neonate arrhythmia algorithm

PVM-4000 series monitor has software for monitoring neonates. A dedicated neonate arrhythmia recognition algorithm recognizes narrow QRS.

NIBP inflation pressure for neonate

Connection of a neonate NIBP air hose is detected and the inflation pressure range is automatically changed for patient safety.

OCRG

OCRG (oxycardiorespirogram) combines compressed trends of beat-to-beat heart rate, respiration, and oxygenation levels. OCRG can help doctors detect the cause of apnea attack.

Mainstream CO₂ sensors for neonate

Mainstream CO₂ monitoring for neonate is available with a cap-ONE TG-980P CO₂ sensor kit. The airway adapter YG-213T has a dead space of only 0.5 ml.

Non-adhesive SpO₂ probe

The TL-535U (P206) and TL-260T (P205A) Multi-site Y probe are designed so that the adhesive of the attachment tape does not stick to the delicate skin of neonates or low birth weight infants.

Network:

LS-NET (Life Scope Networking)

You can connect Nihon Kohden central monitors and bedside monitors by LS-NET (Life Scope Networking with Ethernet LAN).

Interbed

When Nihon Kohden bedside monitors are connected in an LS-NET network, you can exchange and view data of other bedside monitors. Both individual display and 8 bed display are available.

Connecting to a network printer without a central monitor

You can connect a network printer to a PVM-4000 series bedside monitor on the LS-NET through the LAN port.

You can print real time and review data on A4 papers.

Other Features:

Sleep mode

Sleep mode prevents the monitor from disturbing the patient during sleep or other times. The screen is darkened and the sync sound and alarm indicator are turned off. Sleep mode is available when the PVM-4000 series bedside monitor is connected to the central monitor.

Graphical and tabular trends

Up to 120 hours of graphic and tabular trends of all parameters can be saved and reviewed.

Full disclosure

Up to 120 hours of full disclosure waveforms for 4 parameters can be reviewed.

Data time is synchronized between review windows

When one review window is changed to another review window, the new review window opens in the same time as the first review window. For example, to check the waveform information for the time of an alarm file in the Alarm History window, select the corresponding file and change the full disclosure window. The full disclosure waveform for the corresponding time is displayed. This makes it easy to compare various vital sign data for the patient at the selected time between the review windows.

Alarm escalation

The alarm escalation function automatically escalates the alarm priority to a selected level if the alarm continues or SpO₂ drops. Alarm escalation applies to vital alarms (SpO₂ limit and apnea) and technical alarms (ECG and SpO₂). Alarm escalation can contribute to decreased medical accidents and improved quality of care.

Specification

◆ Measuring parameters

PVM-4761/4751/4731

ECG(3/6 lead), respiration (impedance method), SpO₂, NIBP, temperature x 2, esCCO*¹

PVM-4763/4753/4733

ECG(3/6 lead), respiration (impedance method), SpO₂, NIBP, temperature x 2, CO₂(mainstream), IBP, esCCO*¹,

*¹: Optional software QP-470P is required.

◆ Display

Display size: 10.4 inch color TFT LCD

Resolution: 800 x 600 dots

Pixel pitch: 0.264 x 0.264

Touch panel operation: available

Waveform display mode: Fixed method

Touch panel operation: available

Sweep speed:

6.25 mm/s, 12.5 mm/s, 25 mm/s, 50 mm/s (when displaying ECG, pulse wave or IBP)

1.56 mm/s, 6.25 mm/s, 12.5 mm/s, 25 mm/s (when displaying respiration waveform or CO₂)

Display colors: 12 colors, selectable

Maximum number of waveform traces:

PVM-4763/4753/4733: Up to 6 traces

PVM-4761/4751/4731: Up to 4 traces

Displayed waveforms on home screen:

ECG (up to 2 waveforms), respiration waveform, IBP (up to 2 waveforms), SpO₂ pulse wave, CO₂ partial pressure curve

Numerical data display:

Heart rate, VPC (per-minute value), ST measurement, SpO₂ value, pulse rate, PI, temperature, NIBP (SYS/DIA/MAP), IBP (SYS/DIA/MAP), PPV, SPV, ETCO₂, FiCO₂, QTc, QRSd, RPP, SI

With QP-470P;

esCCO, esCCI, esSV, esSVI

Synchronization marks: Heart rate sync mark, pulse rate sync mark, respiratory sync mark

◆ Alarm

Alarm items: Upper/lower alarms, arrhythmia alarms, interbed alarms, technical alarms

Alarm levels: Crisis (flashes red), Warning (flashes yellow), Advisory (flashes blue or yellow)

Alarm indication: Message, highlighted numeric value, alarm indicator flashes, alarm sound

Turning alarm off: Silencing alarms, alarm reset, suspend alarms, all alarms off

Alarm escalation: Apnea, SpO₂ lower, technical alarm (check electrodes, cannot analyze, check SpO₂ probe)

Alarm trigger delay:

RR upper/lower: 0 to 30 seconds

HR/PR upper/lower: 0 to 10 seconds

SpO₂: 0 to 30 seconds

Alarm configurable start-up values:

Master Alarm Neonatal, Master Alarm Pediatric, Master Alarm Adult

Number of master setting:

Up to 3 for each of Adult, Pediatric and Neonatal mode.

Alarm limit setting: Displayed for each parameter on the monitoring screen.

Auto setting: Upper/Lower alarm, ST level

Interbed alarm setting: All, Crisis & Warning, Crisis, Off

◆ Power requirement

AC: 100 to 240 V

DC (battery): 6.6 to 8.2 V

Line frequency: 50 or 60 Hz

◆ Battery

Battery pack SB-470P

Type of battery: Lithium ion rechargeable battery

Number of battery: 1

Battery operation time: 6 hours

The above information is a guideline when all the following conditions are fulfilled and ECG is monitored.

- New and fully charged battery pack
- Used at room temperature (around 25 °C)
- Power saving mode on
- No optional components
- No alarm occurs

◆ Interface

Network socket: 1

ZS socket (for ZS-900P): 1

USB socket: 1

Alarm output socket*: 1

Serial communication socket*: 1

RGB output socket*: 1

ECG/BP OUT socket*: 1

*When QI-470P interface is installed

◆ Recording

Recording format: Thermal array recording

Number of traces: Up to 3

Recording width: ≥ 46mm

Recording speed: 12.5, 25, 50 mm/s

◆ Network

Network communication (LS-NET):

Connects to the Nihon Kohden monitor network and communicates with central monitor, other bedside monitors and server.

Network communication (hospital network):

Connects to the hospital network (HL7 output).

◆ Review

Review display time: 120 hours

Trend graph:

Number of graphs: 2

Number of parameters in each graph: Up to 3

Vital sign list:

Number of parameters displayed in each list: Up to 15

List interval: 1, 5, 10, 15, 30 or 60 min

NIBP list:

Number of parameters displayed in the NIBP list:
Up to 15

Full disclosure:

Displayed waveform: Up to 4

Display time for 1 trace: 60, 30, 20, 12 or 5 seconds

Zoom in window (displays a full disclosure a waveform with 1 trace of 5 seconds, and additional information): Available

Alarm history:

Number of files on one screen: 9 items

Arrhythmia recall:

Filters for displayed items: Arrhythmia types

Individual display mode (Arrhythmia event and the waveform for 5 seconds before and after the arrhythmia event): Available

OCRG review: Available

◆ ECG

Leads:

3-electrode cable: I, II, III

6-electrode cable: I, II, III, aVR, aVL, aVF, 2 from V1 to V6

Number of waveforms: Up to 8 (All Leads page)

Frequency response:

DIAG mode: 0.05 to 150 Hz (–3 dB)

MONITOR mode: 0.3 to 40 Hz (–3 dB)

MAXIMUM mode: 1 to 18 Hz (–3 dB)

Heart rate counting range:

0, 15 to 300 beats/min (± 2 beats/min)

Heart rate counting accuracy:

± 2 beats/min

Arrhythmia analysis:

Analysis method: Multi-template matching method

Analysis channels: 2 channels

VPC counting rate: 0 to 99 VPCs/min

Arrhythmia alarm items:

25 items (ASYSTOLE, VF, VT, EXT TACHY, EXT BRADY, VPC RUN, V BRADY, SV TACHY, TACHYCARDIA, BRADYCARDIA, PAUSE, COUPLET, EARLY VPC, MULTIFORM, V RHYTHM, BIGEMINY, TRIGEMINY, FREQ VPC, VPC, A-Fib, End A-Fib, IRREGULAR RR, PROLONGED RR, NO PACER PULSER, PACER NON-CAPTURE)

ST level measurement:

Number of measurement channels:

3-electrode: 1 ch

6-electrode: 2 ch

ST level measuring range: ± 2.5 mV

QTc/QRSD measurement:

QTc/QRSD lead: All, Trace 1, Selected lead

QTc/QRSD value display on the home screen:
Selectable

◆ Respiration (impedance)

Measuring method: Impedance method

Number of channels: Selectable from R-F or R-L

Respiration rate counting range: 0 to 150 counts/min

Respiration rate counting accuracy*: ± 2 counts/min (0 to 150 counts/min)

* Essential performance in EMC standard

Apnea detection time: Off, 5 to 40 s

◆ Non Invasive Blood Pressure, NIBP

Measuring method: Oscillometric

Measuring range: 0 to 300 mmHg

Accuracy: ± 3 mmHg

Measurement accuracy (based on ISO 81060-2)

Standard deviation: $\leq \pm 5$ mmHg

Difference of determination: $\leq \pm 8$ mmHg

Cuff inflation time:

Adult/child: ≤ 11 s (700 cc), 0 to 200 mmHg

Neonate: ≤ 5 s (72 cc), 0 to 200 mmHg

Initial pressurization value (default setting):

Adult: 180 mmHg*

Child: 140 mmHg*

Neonate: 100 mmHg*

*It is settable

Maximum pressurization value:

Adult/Child: 300 mmHg

Neonate: 150 mmHg

Maximum measurement time:

Adult/Child: ≤ 160 s

Neonate: ≤ 80 s

Operation mode: Manual, STAT, periodic, SIM

Measurement mode: Adult/child or neonate is recognized by connected air hose

Air leakage: ≤ 3 mmHg/min

Display items: Systolic (SYS), diastolic (DIA), mean (MAP)

Operation mode: Manual, STAT, periodic, SIM (depends on the site setting)

Other display: Oscillation graph, PR, Cuff pressure (displays cuff pressure during measurement), RPP, SI

Venous puncture mode: Available*

*target pressure is settable

PWTT triggered NIBP measurement: Available

iNIBP (NIBP measurement in inflation method):

Available

RPP: Available

SI:

Old data display: Dim or hide

Time setting until data is recognized as old: 5, 10, 30min, 1, 24 hours

◆ SpO₂

Measuring technology:

Nihon Kohden, Nellcor or Masimo

Alarm item: SpO₂, PR

Sync sound tone setting: 81-100, 40-100%

Nihon Kohden SpO₂, PVM-4761/4763

SpO₂ measurement:Display range: 0 to 100% SpO₂Declared range: 70 to 100 % SpO₂Measuring accuracy *¹ (rms²):±3%SpO₂ (70%SpO₂ ≤ %SpO₂ < 80%SpO₂)±2%SpO₂ (80%SpO₂ ≤ %SpO₂ ≤ 100%SpO₂)**Pulse rate measurement:**

Display range: 30 to 300 beats/min

Declared range: 30 to 300 beats/min

Counting accuracy*¹ (rms²): ±3 beats/min**PI(Pulse-amplitude Index):** available**SQL bar graph:** available**Nellcor SpO₂ , PVM-4751/4753****SpO₂ measurement:**Display range: 1 to 100% SpO₂Declared range: 70 to 100 % SpO₂Measuring accuracy *¹ (rms²):±2%SpO₂ (adult)±3%SpO₂ (neonate)**Pulse rate measurement:**

Display range: 20 to 300 beats/min

Declared range: 20 to 250 beats/min

Counting accuracy*¹ (rms²): ±3 beats/min**Masimo SpO₂ , PVM-4731/4733****SpO₂ measurement:**Display range: 1 to 100% SpO₂Declared range: 70 to 100 % SpO₂Measuring accuracy¹ (rms²):

No motion:

±2%SpO₂ (adult)±3%SpO₂ (neonate)

Motion:

±3%SpO₂ (adult)±3%SpO₂ (neonate)**Pulse rate measurement:**

Display range: 25 to 240 beats/min

Declared range: 25 to 240 beats/min

Counting accuracy*¹ (rms²):

±3 beats/min (No motion)

±5 beats/min (Motion)

PI (Perfusion Index): available**Signal IQ :** available◆ **Temperature**

Complies with ISO 80601-2-56:2009

Thermistor probe: 400 series (YSI)**Number of channels:** Up to 2**Delta TEMP:** Available**Measuring range:** 0 to 45°C, 32 to 113°F**Measuring accuracy*:**

±0.2°C (0°C ≤ TEMP < 25°C)

±0.1°C (25°C ≤ TEMP ≤ 45°C)

* Essential performance in EMC standard

Internal noise: ≤ 0.03°C (at 37°C)**Temperature drift:** within ±0.005°C /1°C**Cable damage detection:** Detects sensor cable damage**Delta TEMP alarm setting:** 0.1 to 45 °C, Off◆ **Invasive blood pressure, IBP (PVM-4763/4753/4733)**

Complies with IEC 60601-2-34: 2011

Calculation: PPV, SPV**Number of channels:** Up to 2**Measuring range:** -50 to 300 mmHg**Complied transducer:**

P23XL-1, P10EZ-1, Merit Medical Systems DX series disposable transducers

Requirements for other equivalent products:

5 μV/V/mmHg, bridge resistor: 200 Ω to 20 kΩ, defibrillation-proof

Volume displacement: 0.04 mm³/100 mmHg (P10EZ-1 or P23XL-1)**Measuring range:** -50 to 300 mmHg**Zero balancing range:** ±200 mmHg**Zero balancing accuracy:** ±1 mmHg**Single measuring accuracy:**

±1 mmHg ±1 digit (-50 mmHg ≤ IBP < 100 mmHg)

±1% ±1 digit (100 mmHg ≤ IBP ≤ 300 mmHg)

Total measuring accuracy*: ±4% or ±4 mmHg,

whichever is greater (When used with ANSI/AAMI BP-22-1994 compliant product)

*Essential performance in EMC standard

Internal noise: within ±1 mmHg**Temperature drift:** ±0.1 mmHg/1°C**Frequency response:** DC to 12 Hz or 20 Hz (selectable)**Cable damage detection:** Available (detects transducer cable damage)**BP sync sound:** Systolic value 20 to 120 mmHg, changes in 20 steps every 5 mmHg**Zero calibration mode:** Zero, Zero All◆ **Carbon Dioxide, CO₂ (Mainstream method, PVM-4763/4753/4733)****Calculation method:**

TG-900P/TG-920P: semi-quantitative

TG-980P quantitative

CO₂ measuring parameter:TG-900P/TG-920P: ETCO₂TG-980P: ETCO₂, CO₂(I)**CO₂ measuring range:**

TG-900P/TG-920P: 0 to 100 mmHg

TG-980P: 0 to 150 mmHg

Warm-up time:

TG-900P/TG-920P: 5 s

TG-980P: 10 s

Total system response time:

TG-900P/GT-920P/TG-980P: ≤ 0.5 seconds

CO₂ measuring accuracy:

TG-900P/TG-920P:

±0.4 kPa (0 ≤ CO₂ ≤ 1.33 kPa)(±3 mmHg (0 ≤ CO₂ ≤ 10 mmHg))±0.53 kPa (1.33 < CO₂ ≤ 5.33 kPa)(±4 mmHg (10 < CO₂ ≤ 40 mmHg))

±10% reading

(5.33 < CO₂ ≤ 13.3 kPa (40 < CO₂ ≤ 100 mmHg))

(At 1 atmospheric pressure, air inspiration, no condensation)

TG-980P:

±0.27 kPa (0 ≤ CO₂ ≤ 5.33 kPa)(±2 mmHg (0 ≤ CO₂ ≤ 40 mmHg))

±5% reading

(5.33 < CO₂ ≤ 9.33 kPa (40 < CO₂ ≤ 70 mmHg))

±7% reading

(9.33 < CO₂ ≤ 13.3 kPa (70 < CO₂ ≤ 100 mmHg))

(When no condensation)

Warm-up time:

TG-900P/TG-920P: 5 s

TG-980P: 10 s

Total system response time:

TG-900P/GT-920P/TG-980P: ≤ 0.5 seconds

CO₂ value display update cycle: Every 3 s or when alarm is generated

◆ **esCCO (Optional QP-470P is required, PVM-4761/4763 only)**

Measuring method: estimated continuous noninvasive cardiac output derived from pulse wave transit time

Measured parameters: esCCO, esCCI, esSV, esSVI

Measuring range:

esCCO: 0.50 to 20.00 L/min

esCCI: 0.50 to 20.00 L/min/m²

esSV: 0 to 300 mL

esSVI: 0 to 200 mL/m²

◆ **ECG/BP Output**

ECG waveform output: Output the ECG on the first trace

BP waveform output: Output the IBP waveform

HT pulse output: Outputs pulse based on the ECG on the first trace

Output impedance:

ECG: $\leq 100 \Omega \pm 20\%$

BP: $\leq 100 \Omega \pm$

Output-waveform:

ECG: ± 5.0 V (at 1 mV/V $\pm 5\%$ sensitivity)

BP: -0.5 to $+3.0$ V (at 100 mmHg/V $\pm 1\%$ sensitivity)

HT: 5.0 to 15.0 V (Open collector output: 0.5 to 50 mA)

Frequency response:

ECG: ≥ 0.5 to 100 Hz (≥ -3 dB) (No reproducibility of pace maker pulse)

BP: DC to 20 Hz ± 3 Hz (-3 dB)

HT pulse width: 15 ms, 100 ms

Offset:

ECG: $\leq \pm 50$ mV

BP: $\leq \pm 10$ mV

Delay:

ECG: Up to 20 ms

BP: Up to 40 ms

HT: Up to 20 ms

◆ **Dimensions and Weight (approximate)**

PVM-4000 series bedside monitor:

276 W x 237 H x 143 D mm

3.3 kg

WS-470P recorder unit:

166.5 W x 67.5 H x 95.2 D mm

0.375 kg

QI-470P Interface:

124.5 W x 123.5 H x 41.5 D mm

0.125 kg

Safety Standard

IEC 60601-1:2005+Amendment 1:2012

IEC 60601-1-2:2014

IEC 60601-1-6:2010+Amendment 1:2013

IEC 60601-1-8:2006+Amendment 1:2012 1, 3

IEC 60601-1-9:2007

IEC 60601-2-27:2011

IEC 80601-2-30:2009+Amendment 1:2013

IEC 60601-2-34:2011 2, 4

IEC 60601-2-49:2011

IEC 62304:2006+Amendment 1:2015

IEC 62366:2007+Amendment 1:2014

EN ISO 14971:2012

ISO 80601-2-55:20112

ISO 80601-2-56:2009

ISO 80601-2-61:2011

1 Except for the interbed alarm sound

2 PVM-4733, PVM-4753, PVM-4763 only

3 To comply with the requirements from international standards for distributed alarm systems, use the monitor with a central monitor that complies with these standards. The system will not comply with some of the requirements from international standards for distributed alarm systems if used with a non-compliant central monitor.

4 Complies with IEC 60601-2-34:2011 only when "Available Alarm Types" in the [Alarm] → [Disply/Sound] is set to "All".

Environment

Transport and Storage Environment

Temperature: -20 to +65°C (-4 to +149°F)

-15 to +55°C (+5 to 131°F) (recording paper)

Humidity: 10 to 95% RH

Atmospheric pressure: 700 to 1060 hPa

Operating Environment and Power

Operating environment

Temperature: 5 to 40°C (41 to 104°F)

SpO₂ accuracy is guaranteed at surrounding temperature of 18 to 40°C (64 to 104°F).

Relative humidity: 15 to 85% RH (noncondensing)

Atmospheric pressure: 700 to 1060 hPa

Power

Line voltage:

AC: 100 to 240 V

DC (battery): 10.8 V

Power input:

AC: 125 VA

DC (battery): 35 W

Line frequency: 50 or 60 Hz

Battery charging current: 2270 mA

Noise

Equipment used in medically used room: ≤ 48 dBA

Cooling system

Natural cooling (no fans)

Electromagnetic Compatibility

IEC 60601-1-2:2014

STANDARD ACCESSORIES

Bedside monitor, PVM-4761/4763/4731/4733

Power cord (220-240V) 1

Bedside monitor, PVM-4751/4753

Power cord (220-240V)
JL-650P, Nellcor SpO2 connection cord 1

Recorder module, WS-470P

Recording paper, FQW50-2-100 1
Thermal head cleaner pen 1

OPTIONS

Recorder module, WS-470P
Battery pack, SB-470P
Interface, QI-470P
Transmitter, ZS-900P
esCCO program, QP-470P
Software upgrade kit, QS-128P
Wireless LAN station, QI-520P
Cart, KC-470P
Adapter for mounting on the KC-470P, DH-470P
Hook for attaching the bedside monitor to the bed board, DZ-470P
Holder for installing the ZS-900P, DI-471P
Holder for installing the QI-520P, DI-470P

OPTIONAL ACCESSORIES

Accessory set

	Y212A	Y212B	Y213A	Y213B
Model	YO-60IY1	YO-60IY2	YO-60AY1	YO-60AY2
Electrode lead, 3 electrode, clip, IEC, K911 (BR-903P), 0.8 m	✓	✓		
Electrode lead, 3 electrode, clip, AHA, K911A (BR-903PA), 0.8 m			✓	✓
ECG connection cord, IEC, K922 (JC-906P), 3 m	✓	✓		
ECG connection cord, AHA, K922A (JC-906PA), 3 m			✓	✓
SpO ₂ connection cord, K931 (JL-900P), 2.5 m	✓	✓	✓	✓
Air hose for NIBP, S902 (YN-901P), 3.5 m	✓	✓	✓	✓
NIBP cuff, 13 cm width, 23-33 cm range, S951D (YP-713T)	✓		✓	
NIBP cuff, 16 cm width, 33-45 cm range, S951E (YP-714T)		✓		✓

When ordering PVM-4751.4753 (Nellcor), PVM-4731/4733 (Masimo), order accessory set below.

Order code	Y212C	Y212D
Model	YO-60IY3	YO-60IY4
Electrode lead, 3 electrode, clip, IEC, K911 (BR-903P), 0.8 m	✓	✓
ECG connection cord, IEC, K922 (JC-906P), 3 m	✓	✓
Air hose for NIBP, S902 (YN-901P), 3.5 m	✓	✓
NIBP cuff, 13 cm width, 23-33 cm range, S951D (YP-713T)	✓	
NIBP cuff, 16 cm width, 33-45 cm range, S951E (YP-714T)		✓

CONSUMABLES

ECG and Respiration (Impedance)

Cable/cord

- K911** Electrode lead, BR-903P, for 3 electrodes, clip, IEC, 0.8 m
K911A Electrode lead, BR-903PA, for 3 electrodes, clip, AHA, 0.8 m
K917 Electrode lead, BR-923P, for 3 electrodes, clip, IEC, 1.5 m
K916 Electrode lead, BR-963P, for 3 electrodes on one lead, clip, IEC, 0.8 m
K905 Electrode lead, BR-973P, for 3 electrodes, clip, IEC, 1.25 m, for low impedance neonate/child
K906 Electrode lead, BR-975P, for 5 electrodes, clip, IEC, 1.25 m, for low impedance neonate/child
K912 Electrode lead, BR-906P, for 6 electrodes, clip, IEC, 0.8 m
K912A Electrode lead, BR-906PA, for 6 electrodes, clip, AHA, 0.8 m
K918 Electrode lead, BR-926P, for 6 electrodes, clip, IEC, 1.5 m
K922 ECG connection cord, JC-906P, 3/6 electrodes, IEC, 3 m
K922A ECG connection cord, JC-906PA, 3/6 electrodes, AHA, 3 m
K925 ECG connection cord, JC-916P, 3/6 electrodes, IEC, 1.5 m
L021 Adapter, JC-990P, for low impedance neonate/child

Disposable electrodes (Vitrode series)

- G207** L-150X, for general use, X-ray (radiolucent), 30 electrodes x 5/set, 35 mm dia, Vitrode L
G210D F-150M, for adult/child, 3 electrodes x 50/set, 25 x 45 mm, Vitrode F
G210C F-150S, for neonate/child, 3 electrodes x 50/set, 18 x 36 mm, Vitrode F
G221 G-600, for ICU/OR, 30 electrodes x 20/set, 47 x 47 mm, Vitrode G
G272A V-090M3, for general use, 3 electrodes x 30/set, DIN, 25 x 45 mm, lead 1.0 m, Vitrode V
G278A V-091O3, for general use, 3 electrodes x 30/set, DIN, 25 x 45 mm, lead 1.5 m, Vitrode V
G272B V-040M4, for general use, 4 electrodes x 10/set, 25 x 45 mm, lead 1.0 m, DIN, Vitrode V
G273A V-041O4, for general use, 4 electrodes x 10/set, DIN, 25 x 45 mm, lead 1.5 m, Vitrode V
G272C V-060M6, for general use, 6 electrodes x 10/set, DIN, 25 x 45 mm, lead 1.0 m, Vitrode V
G274A V-061O6, for general use, 6 electrodes x 10/set, DIN, 25 x 45 mm, lead 1.5 m, Vitrode V
G271A V-120S3, for NICU, 3 electrodes x 40/set, DIN, 20 x 20 mm, Vitrode V
G300A N-03IS3, for NICU, 3 electrodes x 10/set, DIN, 14 x 25 mm, Vitrode N
G300D N-01IS3, for NICU, 3 electrodes x 10/set, DIN, 15 mm dia, Vitrode N

SpO₂

SpO₂ Cord

- K931** SpO₂ connection cord, JL-900P, 2.5 m
K935 Nellcor SpO₂ connection cord, JL-650P, 3 m
K936 Masimo SpO₂ connection cord, for LNOP series, JL-630P, 3.6 m
K937 Masimo SpO₂ connection cord, for LNCS series, JL-631P, 3 m
K943A Masimo SpO₂ connection cord, for RD SET series, JL-632P, 3.6 m

Probes (BluPRO series)

- P225F** Finger probe, TL-201T, finger-clip type
P311C Finger probe, TL-631T3, regular size, attached to finger or toe by tape
P310C Finger probe, TL-630T3, large size, attached to finger or toe by tape
P225G Multi-site probe, TL-220T, attached to adult/child finger/toe or neonate instep and sole by tape
P267 Probe fastener, YS-093P2, 30 pcs/set, for TL-631T3/TL-630T3

Disposable probes (BluPRO series)

- P203A** TL-271T, for adult finger or toe, 24 pcs/set, 0.8 m
P204A TL-271T, for adult finger or toe, 5 pcs/set, 0.8 m
P203E TL-271T3, for adult finger or toe, 24 pcs/set, 1.6 m
P204E TL-271T3, for adult finger or toe, 5 pcs/set, 1.6 m
P203B TL-272T, for child finger or toe, 24 pcs/set, 0.8 m
P204B TL-272T, for child finger or toe, 5 pcs/set, 0.8 m
P203F TL-272T3, for child finger or toe, 24 pcs/set, 1.6 m
P204F TL-272T3, for child finger or toe, 5 pcs/set, 1.6 m
P203C TL-273T, for neonate instep, 24 pcs/set, 0.8 m
P204C TL-273T, for neonate instep, 5 pcs/set, 0.8 m
P203G TL-273T3, for neonate instep, 24 pcs/set, 1.6 m
P204G TL-273T3, for neonate instep, 5 pcs/set, 1.6 m
P203D TL-274T, for infant finger or toe, 24 pcs/set, 0.8 m
P204D TL-274T, for infant finger or toe, 5 pcs/set, 0.8 m
P203H TL-274T3, for infant finger or toe, 24 pcs/set, 1.6 m
P204H TL-274T3, for infant finger or toe, 5 pcs/set, 1.6 m

Multi-site Y probe

- P205A** Multi-site Y probe, TL-260T (for adult, child, neonate), 5 pcs/set
P260E Sponge attachment tape S for TL-260T multi-site Y probe, 24 pcs/set
P260F Sponge attachment tape L for TL-260T multi-site Y probe, 24 pcs/set
P256 Ear clip adapter for TL-260T, reusable

Disposable probes

- P228A** TL-051S, for adult/neonate, 5 pcs/set, 0.8 m cable

P228B TL-052S, for adult/neonate, 5 pcs/set, 1.6 m cable
P229A TL-061S, for child/infant, 5 pcs/set, 0.8 m cable
P229B TL-062S, for child/infant, 5 pcs/set, 1.6 m cable
P260C Foam tape, for TL-051S/052S/061S/062S, 100 pcs
P206 TL-535U, for neonates and preterm infant, 5 pcs/set
P264A Attachment tape, for TL-535U, S: 1kg or less, 24 pcs
P264B Attachment tape, for TL-535U, L: 2.5kg or less, 24 pcs

Nellcor SpO₂ probes are available from MEDTRONIC (www.medtronic.com) or their suppliers

Masimo SpO₂ probes are available from MASIMO (www.masimo.com) or their suppliers

NIBP

Reusable Cuffs (Latex free)

S943A YP-960T: for infant, 5 cm cuff width
S943B YP-961T: for child, 7 cm cuff width
S943C YP-962T: for child, 10 cm cuff width
S944B YP-963T: for adult, 13 cm cuff width
S944C YP-964T: for adult, 15 cm cuff width
S944D YP-965T: for thigh, 19 cm cuff width

Reusable Cuffs (YAWARA-CUFF 2 series, Latex free)

S951A YP-710T: for infant, 5 cm cuff width
S951B YP-711T: for child, 7 cm cuff width
S951C YP-712T: for adult small, 10 cm cuff width
S951D YP-713T: for adult, 13 cm cuff width
S951E YP-714T: for adult large, 16 cm cuff width
S951F YP-715T: for thigh, 19 cm cuff width

Disposable Cuffs

S945C YP-810P: for infant, 6 cm cuff width, 20 pcs
S945D YP-811P: for child, 8 cm cuff width, 20 pcs
S946E YP-812P: for adult, 10 cm cuff width, 20 pcs
S946F YP-813P: for adult, 14 cm cuff width, 20 pcs
S946G YP-814P: for adult, 15 cm cuff width, 20 pcs
S946H YP-815P: for adult, 17 cm cuff width, 20 pcs
S946I YP-816P: for adult, 18 cm cuff width, 20 pcs
S946J YP-817P: for thigh, 20 cm cuff width, 20 pcs

S948A YP-820P: for neonate, 2 cm cuff width, 10 pcs, Luer connector
S948B YP-821P: for neonate, 3 cm cuff width, 10 pcs, Luer connector
S948C YP-822P: for neonate, 4 cm cuff width, 10 pcs, Luer connector
S948D YP-823P: for neonate, 4.5 cm cuff width, 10 pcs, Luer connector
S948E YP-824P: for neonate, 5 cm cuff width, 10 pcs, Luer connector
S948F YP-820S: for neonate, 2 cm cuff width, 10 pcs, SnapQuick connector
S948G YP-821S: for neonate, 3 cm cuff width, 10 pcs, SnapQuick connector
S948H YP-822S: for neonate, 4 cm cuff width, 10 pcs, SnapQuick connector
S948I YP-823S: for neonate, 4.5 cm cuff width, 10 pcs, SnapQuick connector
S948J YP-824S: for neonate, 5 cm cuff width, 10 pcs, SnapQuick connector

Disposable Cuffs (Nihon Kohden original cuff)

S954A YP-840T: for infant, 5cm cuff width, 20 pcs
S954B YP-841T: for child, 7 cm cuff width, 20 pcs
S954C YP-842T: for adult, 10 cm cuff width, 20 pcs
S954D YP-843T: for adult, 13 cm cuff width, 20 pcs
S954E YP-844T: for adult, 16 cm cuff width, 20 pcs
S954F YP-845T: for thigh, 19 cm cuff width, 20 pcs

Air hoses

S901 YN-900P: 1.5 m
S902 YN-901P: 3.5 m
S903 YN-990P: extension hose, 1.5 m
S904 YN-920P: for neonate, 1.5 m, Luer connector
S905 YN-921P: for neonate, 3.5 m, Luer connector for IEC 2nd edition
S904A YN-924P: for neonate, 1.5 m, SnapQuik connector
S905A YN-925P: for neonate, 3.5 m, SnapQuick connector for IEC 3rd edition

IBP

Merit Medical (Argon):

K951 IBP connection cord, JP-900P, Merit Medical blood pressure transducers are available from Merit Medical (www.merit.com) or their supplier.

Edwards Lifesciences:

L901 IBP connection cord, JP-920P, for Edwards Lifesciences transducer, 3.5 m Edwards Lifesciences blood pressure transducers are available from Edwards Lifesciences Corporation (www.edwards.com)

or their suppliers.

Biosensors International:

K957 IBP connection cord, JP-960P, for Biosensor transducer, 3.5 m

Biosensors blood pressure transducers are available from BIOSENSORS

INTERNATIONAL (www.biosensorsintl.com) or their suppliers.

Others:

K952 IBP connection cord, JP-910P

Temperature

K961 Temperature connection cord, JT-900P, 30 cm

Thermistor Probe

P240B For adult rectum/esophagus, 401J

P241B For child rectum/esophagus, 402J

P242D DisG2c type, 409J

P249A Probe cover, for 401J, 10 pcs/set

P252 Temperature insulation pad, for 409J, 60 pcs/set

CO₂ (Mainstream measurement)

Semi-quantitative method

P903 CO₂ sensor kit, TG-900P (TG-101T + JG-900P), 3 m

P922A CO₂ sensor, TG-101T, 1 m

K981 CO₂ adapter, JG-900P, 2 m

R801 Airway adapter, YG-101T, 50 pcs/set

P907 CO₂ sensor kit, TG-920P (TG-121T + JG-920P), 3.5 m

P923 CO₂ sensor, TG-121T, 2 m

K984 CO₂ adapter, JG-920P, 1.5 m

R804 Airway adapter, YG-111T, 30 pcs/set

V921 Nasal adapter, YG-121T, for nasal, 30 pcs/set

V922 Nasal adapter, YG-121T, for nasal/oral, 30 pcs/set

V923 Nasal adapter, YG-122T, for nasal/oral, 30 pcs/set, with O₂ supply tube adapter

V927A Oxygen cannula, 25 pcs/set

Quantitative method

P910A CO₂ sensor kit, TG-980P, 3.5 m

R805 Airway adapter, adult YG-211T, 30 pcs/set

R806 Airway adapter, child/neonate, YG-213T, 30 pcs/set

R807 Airway adapter, pediatric, YG-214T, 30 pcs/set

V939A cap-ONE biteblock for adult, YG-227T, 20 pcs/set

V938A cap-ONE mask, Adult, YG-272T, 30kg or more, 10 pcs/set

V938C cap-ONE mask, Large Adult, YG-282T, 40kg or more, 10 pcs/set

V933 cap-ONE mask, Pediatric, YG-232T, 20-

40kg, 10 pcs/set

V935 cap-ONE mask, Infant, YG-242T, 7-20kg, 10 pcs/set

CO

K962 CO connection cord, JT-950P, 2.0 m

Merit Medical:

Merit Medical injectate supply systems and thermodilution catheters are available from Merit Medical (www.merit.com) or their supplier.

Edwards Lifesciences:

Edwards Lifesciences injectate supply systems and thermodilution catheters are available from Edwards Lifesciences (www.edwards.com).

Recording

A721 Recording paper, FQW50-2-100, 10 pcs/set, 50 mm x 20 mm

Y011 Thermal head cleaning pen, 5 pcs/set

ECG/BP out

K974 ECG/BP output cable, YJ-910P, 5 m

K975 ECG/BP output cable, YJ-920P, 0.3 m

 **NIHON KOHDEN**

NIHON KOHDEN CORPORATION

1-31-4 Nishiiochiai, Shinjuku-ku, Tokyo 161-8560, Japan

Phone +81 (3) 5996-8036 Fax +81 (3) 5996-8100

www.nihonkohden.com